

Red Blood Cell Lysis Buffer

For preferential lysis of red blood cells from human whole blood

Cat. No. 11 814 389 001

100 ml for 50 to 500 reactions, depending on sample size (1 to 500 μ l)

Version 09
Content version: April 2016

Store at +2 to +8°C

1. What this Product Does

Number of Tests

50 to 500 reactions, depending on sample size (1 to 500 μ l)

Contents

100 ml ready-to-use solution.

Storage and Stability

Store the product at +2 to +8°C until the expiration date printed on the label.

Application

Used for both DNA and RNA isolation, the buffer is designed for the preferential lysis of red blood cells from human whole blood, yielding intact white blood cells (free of red blood cells) for further applications. This buffer is not intended for use with whole blood from any other species.

2. How to Use this Product

2.1 Before You Begin

General Handling Recommendations

- Warm the buffer to +15 to +25°C prior to use.
 - ⚠ Do not** use blood that has been frozen and thawed more than three times for any DNA applications.
- Warm blood to +15 to +25°C.
 - ⚠ Do not** use blood that has been stored for longer than 1 month for any DNA applications.
- Perform all centrifugation steps at +15 to +25°C in a variable-speed microfuge.
- Use blood stored in EDTA, citrate, or heparin anticoagulants.
- To obtain RNA, use fresh blood only.
- Use blood that has been stored for ≤ 1 month at +15 to +25°C, +2 to +8°C, or -15 to -25°C. For best results in DNA applications, use fresh blood or blood stored for ≤ 3 days.
- ⚠** In addition, be certain to follow all universal safety precautions governing work with biohazardous materials. For example, wear lab coats, gloves, and safety glasses at all times. Also, properly dispose of all contaminated materials, decontaminate work surfaces, and use a biosafety cabinet whenever aerosols might be generated.

2.2 Procedure

Lysis of Red Blood Cells from 500 μ l Human Whole Blood

The following procedure was optimized for 500 μ l samples of blood. If using < 500 μ l of blood, see "For use with smaller quantities of blood" for necessary modifications to this procedure.

Step Action

- For each blood sample to be processed, add 1 ml Red Blood Cell Lysis Buffer to a sterile 1.5 ml microfuge tube.
 - To each tube, add 500 μ l human whole blood (see "General Handling Recommendations"). Cap the tube, and mix the contents by inversion. Do not vortex.
 - Place the microfuge tube on a rocking platform or gyratory shaker for 10 min at +15 to +25°C. Alternatively, manually invert the sample periodically for 10 min.
 - Centrifuge the tube at 500 $\times g$ for 5 min in a microfuge at +15 to +25°C.
 - With a sterile pipette, **carefully** remove and properly dispose of the clear, red supernatant that is indicative of complete red cell lysis.
 - After removal of the supernatant, a white pellet should be visible at the bottom of the tube. However, if two distinct layers (a cloudy white upper layer containing plasma/leukocytes and a red lower layer containing erythrocytes) are visible in the supernatant, no cell lysis has occurred. If so, repeat steps 1 to 5 with fresh blood and also:
 - verify the Red Blood Cell Lysis Buffer is equilibrated to +15 to +25°C prior to use
 - invert the sample more frequently if mixing by hand
 - use a higher ratio (e.g., 3:1) of Red Blood Cell Lysis Buffer to blood
 - use a 15 min incubation in step 3.
- | If... | Then... |
|---|---|
| the sample will be further purified for downstream applications | disregard steps 7 to 10 |
| no further purification will be performed | continue with steps 7 to 10 to remove residual RBC material |
- Add 1 ml Red Blood Cell Lysis Buffer, cap the tube, and mix by "flicking" the tube until the pellet is resuspended.
 - ⚠ Do not vortex.**
 - Centrifuge the tube at 500 $\times g$ for 3 min in a microfuge at +15 to +25°C.
 - With a sterile pipette, carefully remove and properly dispose of the supernatant, particularly the red ring of blood-cell debris that forms around the outer surface of the white pellet.
 - Resuspend the white pellet in an appropriate buffer.

2.3 For Use with Smaller Quantities of Blood

By slightly adjusting the protocol provided in "For preferential lysis of red blood cells from 500 μ l human whole blood," blood samples of 1 to 500 μ l can be processed. Follow the procedure described above with these modifications:

For a Blood Volume of	Use the following volume of Red Blood Cell Lysis Buffer in steps 1 and 7 of the procedure
400 to 500 μ l	1 ml
300 to < 400 μ l	800 μ l
200 to < 300 μ l	600 μ l
100 to < 200 μ l	400 μ l
25 to < 100 μ l	200 μ l
5 to < 25 μ l	200 μ l
1 to < 5 μ l	200 μ l

3. Additional Information on this Product

How this Product Works

Red Blood Cell Lysis Buffer is designed for the preferential lysis of red blood cells from human whole blood.

Use this buffer on 1 to 500 μ l human whole blood (stored one month or less at +15 to +25°C, +2 to +8°C, or -15 to -25°C) to isolate white blood cells that are free of red blood cells. As most blood cells are red blood cells, which lack nuclei and therefore possess no DNA, the lysis and the centrifugation steps concentrate the nucleated white blood cells. This is important for both DNA and RNA isolation.

Red Blood Cell Lysis Buffer eliminates the need for hazardous organic extractions or chaotropic agents.

The buffer is particularly useful for high-volume research currently requiring Ficoll/Hypaque gradients.

4. Supplementary Information

Changes to Previous Version

- Editorial changes


Text Conventions

To make information consistent and understandable, the following text conventions are used in this document:

Text Convention	Use
Numbered instructions labeled ①, ②, etc.	Steps in a procedure that must be performed in the order listed.

Symbols

The following symbol is used in this document to highlight important information:

Symbol	Description
	Important Note: Information critical to the success of the procedure or use of the product.

Ordering Information

Roche Diagnostics offers a large selection of reagents and systems for life science research. For a complete overview of related products and manuals, please visit and bookmark our home page <http://lifescience.roche.com>.

Trademarks

All brands or product names are trademark of their respective holders.

Regulatory Disclaimer

For life science research only. Not for use in diagnostic procedures.

Disclaimer of License

For patent license limitations for individual products please refer to: [List of biochemical reagent products](#)

Contact and Support

To ask questions, solve problems, suggest enhancements and report new applications, please visit our [Online Technical Support Site](#).

To call, write, fax, or email us, visit sigma-aldrich.com, and select your home country. Country-specific contact information will be displayed.



Roche Diagnostics GmbH
Sandhofer Strasse 116
68305 Mannheim
Germany