

# Annexin-V-FLUOS Staining Kit

Kit for the detection and quantification of apoptosis and differentiation from necrosis at single cell level, based on Annexin-V-labeling

**Cat. No. 1 858 777** 50 tests

**Cat. No. 1 988 549** 250 tests

**Version 3, March 2003**

Store at 2–8°C

## 1. Product overview

### Kit contents

Vial/ Cap	Label	Content/Cat.No.		Use
		1 858 777	1 988 549	
1 green	Annexin-V- Fluorescein	110 µl	500 µl	Ready-to-use
2 red	Propidium iodide	150 µl	500 µl	<ul style="list-style-type: none"> <li>Ready-to-use</li> <li>For the preparation of the Annexin-V-Fluorescein labeling solution</li> </ul>
3 blue	Incubation buffer	50 ml HEPES buffer	4 × 50 ml HEPES buffer	<ul style="list-style-type: none"> <li>Ready-to-use</li> <li>For the dilution of the Annexin-V-Fluorescein solution</li> </ul>


### Introduction

In the early stages of apoptosis, changes occur at the cell surface (1, 2, 3). One of these plasma membrane alterations is the translocation of phosphatidylserine (PS) from the inner part of the plasma membrane to the outer layer, by which PS becomes exposed at the external surface of the cell (4). Fadok et al. showed that macrophages specifically recognize PS exposed on the surface of lymphocytes during the development of apoptosis (2). The recognition and phagocytosis of apoptotic cells and bodies protects organisms from the exposure to cellular compounds leading to inflammation, which mostly accompanies necrosis.

### Assay principle

The analysis of phosphatidylserine on the outer leaflet of apoptotic cell-membranes is performed by using Annexin-V-Fluorescein and Propidium iodide (PI) for the differentiation from necrotic cells or labeling with a cell surface marker for cell characterization. The procedure involves:

Stage	Description
1	Washing the cells in PBS.
2	Incubation of cells with Annexin-V-Fluorescein in a HEPES buffer containing PI or labeling reagent for cell surfaces (e.g., CD-marker).
3	Analysis of the samples under a fluorescence microscope or on a flow cytometer.



Apoptotic cell with phosphatidylserine exposed on the outer leaflet of the membrane

Annexin-V-Fluos

**Fig. 1:** Test principle.

### Application

Annexin V is a  $Ca^{2+}$ -dependent phospholipid-binding protein with high affinity for phosphatidylserine (4). This protein can hence be used as a sensitive probe for PS exposure upon the outer leaflet of the cell membrane and is therefore suited to detect apoptotic cells (4, 5, 6, 7) in cell populations but not on tissue sections. Since necrotic cells also expose PS according to the loss of membrane integrity, apoptotic cells have to be differentiated from these necrotic cells. The simultaneous application of a DNA stain which is used for dye exclusion tests allows the discrimination of necrotic cells from the Annexin V positively stained cell cluster. Any other secondary labeling should be possible, e.g., membrane surface staining with a phycoerythrin or TRITC-labeled monoclonal antibody for further cellular characterization (8).

### Sample material

- Cell lines
- Freshly isolated cells

### Number of tests

For 50 tests (Cat. No. 1 858 777)  
For 250 tests (Cat. No. 1 988 549)

### Preparation

Recombinant Annexin-V is produced in *E. coli* (strain NB42). The GST-tagged protein is purified by standard purification protocols.

### Fluorescence characteristics

Annexin-V-Fluorescein and propidium iodide show the following fluorescence characteristics:

	Fluorescein	Propidium Iodide
<b>Excitation</b>	488 nm	488–540 nm
<b>Emission</b>	518 nm	617 nm

### Specificity

Annexin-V-Fluos binds in a  $Ca^{2+}$ -dependent manner to negatively charged phospholipid surfaces and shows high specificity to phosphatidylserine. Therefore, it stains apoptotic as well as necrotic cells. Propidium iodide stains DNA of leaky necrotic cells only.

### Storage/Stability

Stable at 2–8°C until the expiration date printed on the label.

## 2. Procedures and required material

### Additional solutions required

PBS

### Preparation of Annexin-V-FLUOS labeling solution

Predilute 20 µl Annexin-V-Fluos labeling reagent (vial 1) in 1 ml Incubation buffer (bottle 3) and add 20 µl Propidium iodide solution (vial 2).  
**Note:** 1 ml is enough for 10 samples.

### Staining of cell suspensions

In the following table please find the staining procedure for cell suspensions.

Step	Action
1	Wash $10^6$ cells with PBS and centrifuge cells at $200 \times g$ for 5 min
2	Resuspend the cell pellet in 100 µl of Annexin-V-FLUOS labeling solution. Incubate 10–15 min at 15–25°C.
3	Analyze by fluorescence microscopy or on a flow cytometer (see 3. Analysis).

## Staining of adherent cells

In the following table please find the staining procedure for adherent cells.

Step	Action
1	Remove chambers and silicon borders of cells grown on chamberslides.
2	Remove medium and cover slides with Annexin-V-FLUOS labeling solution (100 µl/chamber).
3	Put coverslips on slides and incubate for 10–15 min at 15–25°C.
4	Analyze by fluorescence microscopy or on a flow cytometer (see 3. Analysis). <b>Note:</b> We do not recommend to analyze adherent cells by flow cytometry, because trypsinization or scraping for monodispersion of the cells results in false positive staining and analysis of non-dispersed cell clusters.

## 3. Analysis

### Fluorescence-microscopy


For evaluation by fluorescence microscopy use an excitation wavelength in the range of 450–500 nm (e.g., 488 nm) and detection in the range of 515–565 nm (green).

### Flowcytometry

Add 0.5 ml Incubation buffer (bottle 3) per 10<sup>6</sup> cells and analyze on a flow cytometer using 488 nm excitation and a 515 nm bandpass filter for fluorescein detection and a filter > 600 nm for PI detection. Electronic compensation of the instrument is required to exclude overlapping of the two emission spectra. Typical histograms of apoptotic versus non-apoptotic and necrotic cells are shown in figure 2.

Figure 2

FACS analysis of apoptotic U937 cells after staining with Annexin-V-FLUOS and propidium iodide. Cultivation for 4 h in the presence (lower row) or absence (upper row) of 4 µg/ml camptothecin: a) single parameter Annexin-V-FLUOS, b) single parameter propidium iodide and c) dual parameter (FL1 = Annexin-V-FLUOS; FL2 = propidium iodide); Cluster R1 = living cells, R2 = apoptotic cells and R3 = necrotic cells.



## 4. References

- Andree, H.A.M. et al. (1990), *J. Biol. Chem.* **265**, 4923.
- Fadok, V. et al. (1992) *J. Immunology* **148**, 2207.
- Creutz, C.E. (1992) *Science* **258**, 924.
- Vermes, I. et al. (1995) *J. Immunol. Methods* **184**, 39.
- Koopman, G. et al (1994) *Blood* **84**, 1415.
- Homburg, C.H.E. et al (1995) *Blood* **85**, 532.
- Verhoven, B. et al (1995) *J. Exp. Med.* **182**, 1597.
- van Engeland, M. et al. (1996) *Cytometrie* **24**, 131.

## 5. Related products

This table only shows a selection of the most important products related to the product described in this pack insert.

For further information please access our web-site address at:

<http://www.roche-applied-science.com>

or the Apoptosis special interest site:

<http://www.roche-applied-science.com/techserv/apoptosis/index.htm>

Apoptosis-specific physiological change	Detection mode/ Product	Pack size	Cat. No.
<b>DNA fragmentation</b>	<b>Gel Electrophoresis</b>		
	Apoptotic DNA-Ladder Kit	20 tests	1 835 246
	<b>In situ assay</b>		
	In Situ Cell Death Detection Kit, TMR red	1 kit (50 tests)	2 156 792
	In Situ Cell Death Detection Kit, Fluorescein	1 kit (50 tests)	1 684 795
	In Situ Cell Death Detection Kit, AP	1 kit (50 tests)	1 684 809
	In Situ Cell Death Detection Kit, POD	1 kit (50 tests)	1 684 817
	<b>Single reagents for TUNEL and supporting reagents</b>		
	TUNEL AP	70 tests (3.5 ml)	1 772 457
	TUNEL POD	70 tests (3.5 ml)	1 772 465
	TUNEL Enzyme	2 × 50 µl (20 tests)	1 767 305
	TUNEL Label	3 × 550 µl (30 tests)	1 767 291
	TUNEL Dilution Buffer	20 ml	1 966 006
	<b>ELISA</b>		
Cell Death Detection ELISA	1 kit	1 544 675	
Cell Death Detection ELISA <sup>PLUS</sup>	1 kit (96 tests)	1 774 425	
Cell Death Detection ELISA <sup>PLUS</sup> , 10×	1 kit	1 920 685	
Cellular DNA Fragmentation ELISA	1 kit (500 tests)	1 585 045	
<b>Cell membrane alterations</b>	<b>Microscopy or FACS</b>		
	Annexin-V-Alexa 568	250 tests	1 985 485
	Annexin-V-Biotin	250 tests	1 828 690
	Annexin-V-FLUOS	250 tests	1 828 681
	Annexin V FLUOS Staining Kit	50 tests 250 tests	1 858 777 1 988 549
<b>Enzymatic activity</b>	<b>Western Blot</b>		
	Anti-Poly (ADP-Ribose) Polymerase	100 µl	1 835 238
	<b>FIENA</b>		
	Caspase 3 Activity Assay	1 kit	2 012 952
	<b>Fluorimetric microplate Assay</b>		
	Homogeneous Caspase Assay, fluorimetric	100 tests 1000 tests	3 005 372 2 236 869
	<b>In situ Assay</b>		
	M30 CytoDEATH (formalin grade)	50 tests 250 tests	2 140 322 2 140 349
	M30 CytoDEATH, Fluorescein	250 tests	2 156 857
	<b>Expression of apoptosis-related proteins</b>	<b>Apoptosis Induction</b>	
Anti-Fas (CD95/APO-1)		1000 tests	1 922 459
<b>In situ Assay/Western Blot</b>			
Anti-bcl-2 Oncoprotein		1 ml	1 624 989
Anti-p53-Protein pan (BMG 1B1)		200 µg	1 810 928
Anti-p53-POD		50 U	1 810 944
<b>ELISA</b>			
p53 pan ELISA	1 kit	1 828 789	

## How to contact Roche Applied Science

[www.roche-applied-science.com](http://www.roche-applied-science.com)

to order, solve technical queries, find product information, or contact your local sales representative.

[www.roche-applied-science.com/pack-insert/1988549a.pdf](http://www.roche-applied-science.com/pack-insert/1988549a.pdf)



Roche Diagnostics GmbH  
Roche Applied Science  
Nonnenwald 2  
82372 Penzberg  
Germany