

Reticulum II Staining Kit

REF

860-024

05279399001


IVD
 75


Figure 1. Reticulum II Staining Kit on liver tissue.

INTENDED USE

The Reticulum II Staining Kit is intended for laboratory use as a qualitative histologic stain to demonstrate reticular fibers by light microscopy in sections of formalin-fixed, paraffin-embedded (FFPE) tissue stained on a BenchMark Special Stains instrument.

This product should be interpreted by a qualified pathologist in conjunction with histological examination, relevant clinical information, and proper controls.

This product is intended for in vitro diagnostic (IVD) use.

SUMMARY AND EXPLANATION

The Reticulum II Staining Kit is a modification of Gordon and Sweets Stain.¹

Reticular fibers are the main component of connective tissue that provide a supporting framework to organs like the liver.² Reticular fibers, also known as reticulin fibers in bone marrow, are loosely packed fibrils made up of type III collagen that provide individual support to cellular organs in a mesh-like network.³ The loss of normal architecture is a key histopathologic feature of chronic liver disease, which can be visualized as replacement of normal tissue with fibrous scar tissue, regeneration and hyperplasia of liver cells, and/or condensation of reticular fibers.⁴⁻⁶

Additionally reticulin fibers are one of the structural fibers that make up bone marrow.^{7,8} Bone marrow fibrosis is characterized by an increase in the deposition of reticulin fibers, and is evaluated based on the density and type of fibrosis.⁸

Reticular fibers are argyrophilic due to their proteoglycan content, which can be visualized by silver stains.^{1,9}

The Reticulum II Staining Kit is used to aid the pathologist in the assessment of liver architecture and diagnosis of liver disease and also assessment of fibrosis in bone marrow biopsies.

PRINCIPLE OF THE PROCEDURE

The staining reaction is based on the affinity of silver for oxidized glycoproteins. The oxidizer, with potassium permanganate, oxidizes the tissue (particularly the carbohydrate component of reticular fibers) to enhance staining of reticular fibers.^{1,2} The Decolorizer, with oxalic acid, removes excess potassium permanganate. The Sensitizer, with ferric ammonium sulfate, is added to form a metal organic compound. The metal organic compound is replaced by the silver in Reticulum II Silver A. The Reticulum II Reducer is applied to develop the deposited silver into visible silver. The Toner contains gold chloride for better contrast and clarity. The Fixer, with thiosulfate, stops the reaction and removes any unreacted silver from the section. The Nuclear Fast Red Counterstain is applied to provide contrasting background.

This kit is optimized for use on BenchMark Special Stains instruments. The reagents are applied to tissue on microscope slides and mixed over the entire specimen.

MATERIAL PROVIDED

The reagent vials are supplied in barcode labeled carriers to insert into the reagent tray of the instrument. Each kit contains sufficient reagent for 75 tests:

One 27 mL vial of Oxidizer contains approximately 1% potassium permanganate.

One 22 mL vial of Decolorizer contains approximately 1% oxalic acid.

One 27 mL vial of Sensitizer contains approximately 2% ferric ammonium sulfate.

One 22 mL vial of Reticulum II Silver A contains approximately 1.5% silver carbonate.

One 22 mL vial of Reticulum II Reducer contains approximately 0.4% formaldehyde.

One 22 mL vial of Toner contains approximately 1% gold chloride

One 22 mL vial of Fixer II contains approximately 2% sodium thiosulfate.

One 22 mL vial of Nuclear Fast Red Counterstain contains approximately 0.2% Nuclear Fast Red and 5.0% aluminum sulfate.

Eight vial inserts with sipping straws.

Reconstitution, Mixing, Dilution, Titration

No reconstitution, mixing, dilution, or titration of kit reagents is required. Further dilution of any of the reagents may result in unsatisfactory staining.

The reagents in this kit have been optimally diluted for use on BenchMark Special Stains instruments.

MATERIALS REQUIRED BUT NOT PROVIDED

Not all products listed in the method sheet may be available in all geographies. Consult your local support representative.

The following reagents and materials may be required for staining but are not provided:

1. Recommended control tissue
2. Microscope slides, positively charged
3. BenchMark Special Stains instrument
4. BenchMark Special Stains Deparaffinization Solution (10X) (Cat. No. 860-036 / 06523102001)
5. BenchMark Special Stains Liquid Coverslip (Cat. No. 860-034 / 06523072001)
6. BenchMark Special Stains Wash II (Cat. No. 860-041 / 08309817001)
7. General purpose laboratory equipment

STORAGE AND STABILITY

The Reticulum II Staining Kit should be stored at 2-8°C. Refrigerated kit components should be brought to room temperature prior to use.

When properly stored, unopened and opened reagents are stable to the expiration date indicated on the vial label. Do not use reagent beyond the expiration date indicated on the kit.

There are no obvious signs to indicate instability of these reagents; therefore, controls should be run simultaneously with unknown specimens. Contact your local support representative if positive control material shows a decrease in staining as it could indicate reagent instability.

SPECIMEN PREPARATION

Routinely processed FFPE tissues are required for use with this product and BenchMark Special Stains instruments. The recommended tissue fixative is 10% neutral buffered formalin.¹

Perform specimen collection and storage according to Histotechnology: A Self Instructional Text.¹ Cut sections to the appropriate thickness, approximately 4 µm, and place the sections on positively charged glass slides.

1. Dry the slides.¹
2. Print appropriate barcode label(s).
3. Apply barcode labels to the frosted end of the slides prior to loading the slides onto the instrument (see the instrument User Guide for correct application of labels).

Refer to the Instructions for Use section for the recommended protocol for the BenchMark Special Stains instrument.







WARNINGS AND PRECAUTIONS

1. For in vitro diagnostic (IVD) use.
2. For professional use only.
3. CAUTION: In the United States, Federal law restricts this device to sale by or on the order of a physician. (Rx Only)
4. Do not use beyond the specified number of tests.
5. Positively charged slides may be susceptible to environmental stresses resulting in inappropriate staining. Ask your Roche representative for more information on how to use these types of slides.
6. Materials of human or animal origin should be handled as biohazardous materials and disposed of with proper precautions. In the event of exposure, the health directives of the responsible authorities should be followed.^{10,11}
7. Avoid contact of reagents with eyes and mucous membranes. If reagents come in contact with sensitive areas, wash with copious amounts of water.

8. Avoid microbial contamination of reagents as it may cause incorrect results.
9. For further information on the use of this device, refer to the BenchMark Special Stains instrument User Guide, and instructions for use of all necessary components located at navifyportal.roche.com.
10. Consult local and/or state authorities with regard to recommended method of disposal.
11. Product safety labeling primarily follows EU GHS guidance. Safety data sheet available for professional user on request.
12. To report suspected serious incidents related to this device, contact the local Roche representative and the competent authority of the Member State or Country in which the user is established.

This product contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

Table 1. Hazard information.

Hazard	Code	Statement
     	H226	Flammable liquid and vapour.
	H290	May be corrosive to metals.
	H314	Causes severe skins burns and eye damage.
	H317	May cause an allergic skin reaction.
	H350	May cause cancer.
	H411	Toxic to aquatic life with long lasting effects.
	P201	Obtain special instructions before use.
	P210	Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
	P273	Avoid release to the environment.
	P280	Wear protective gloves/ protective clothing/ eye protection/ face protection/ hearing protection.
	P303+ P361+ P353	IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.
	P304+ P340+ P310	IF INHALED: Remove person to fresh air and keep comfortable for breathing. Immediately call a POISON CENTER/ doctor.
	P305+ P351+ P338+ P310	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/ doctor.
	P308+ P313	IF exposed or concerned: Get medical advice/ attention.
	P370+ P378	In case of fire: Use dry sand, dry chemical or alcohol-resistant foam to extinguish.
	P391	Collect spillage.

This product contains gold chloride, hydrochloride, trihydrate. May cause an allergic reaction.

INSTRUCTIONS FOR USE

Prepare Reagent Vial

Before first use, a vial insert and sipping straw must be placed in the reagent vial. Remove the shipping cap from the vial and place the insert and straw into the vial. The insert and sipping straw should be left in the vial, once the vial has been opened.

Staining Procedure

1. Load reagents and slides onto the instrument.
2. Place the soft cap into the slot on the reagent holder when the reagent is in use.
3. Perform the staining run according to the recommended protocol in Table 2 and Table 3 and the instructions in the User Guide.
4. When the run is complete, remove the slides from the instrument.
5. Use the soft cap to cover the reagent vial when reagent is not in use.
6. After use, store the reagents according to the recommended storage conditions.

Recommended Protocol

The parameters for the automated procedures can be displayed, printed and edited according to the procedure in the instrument User Guide.

The following procedures allow flexibility to accommodate user preference. This product has been optimized for use with BenchMark Special Stains instruments but the user must validate results obtained with this product.

Table 2. Recommended staining protocol for Reticulum II Staining Kit on a BenchMark Special Stains instrument for other tissues validated (not including bone marrow) listed in Table 4.

Staining Procedure	S Reticulum II
Protocol Step	Method
Deparaffinization	Select to automate paraffin removal.
Baking (optional)	The default is not selected. 75°C for 4 minutes is recommended.
Select Other Tissues or Bone Marrow	Select Other Tissues to run the default protocol.
Optimize Staining	Select Optimize Staining to enable protocol adjustments.*
Incubation Temperature	The default temperature is 40°C. Select temperature from 37-50°C: Adjust temperature in 1°C increments until desired staining intensity of fibers is achieved.* 37°C, lighter staining of fibers 50°C, darker staining of fibers
Oxidizer	The default incubation time is 4 minutes. Select an incubation time from 4-20 minutes.*
Decolorizer	The default incubation time is 8 minutes. Select an incubation time from 4-20 minutes.* The Decolorizer incubation time should remain at least 2 times the incubation of the Oxidizer.
Sensitizer	The default incubation time is 8 minutes. Select an incubation time from 4-20 minutes.*
Optimize Silver Intensity (Reticulum II Silver A)	The default incubation time is 16 minutes. Select an incubation time from 8-20 minutes: <ul style="list-style-type: none"> 8 minutes, lighter staining of fibers 20 minutes, darker staining of fibers
3 Rinse Cycles (optional)	The default is 1 Rinse Cycle.

Staining Procedure	S Reticulum II
Protocol Step	Method
	If necessary, Select 3 Rinse Cycles to reduce non-specific background staining if necessary.
No Toner (optional)	The default is to apply Toner.
Optimize Counterstain Intensity (Nuclear Fast Red Counterstain) (optional)	The default is 4 minutes. Select an incubation time from 4-16 minutes.* <ul style="list-style-type: none"> 4 minutes, lighter counterstain 16 minutes, darker counterstain
No Counterstain (optional)	<ul style="list-style-type: none"> Select No Counterstain if Nuclear Fast Red Counterstain is not desired.

* To adjust staining preferences, increment the stain temperature and incubation time one parameter at a time.

Table 3. Recommended staining protocol for Reticulum II Staining Kit on a BenchMark Special Stains instrument for bone marrow.

Staining Procedure	S Reticulum II
Protocol Step	Method
Deparaffinization	Select to automate paraffin removal.
Baking (optional)	The default is not selected. 75°C for 4 minutes is recommended.
Select Other Tissues or Bone Marrow	Select bone marrow to run the default protocol.
Optimize Staining	Select Optimize Staining to enable protocol adjustments.*
Incubation Temperature	The default temperature is 40°C. Select temperature from 37-50°C: Adjust temperature in 1°C increments until desired staining intensity of fibers is achieved.* 37°C, lighter staining of fibers 50°C, darker staining of fibers
Oxidizer	The default incubation time is 4 minutes. Select an incubation time from 4-20 minutes.*
Decolorizer	The default incubation time is 4 minutes. Select an incubation time from 4-20 minutes.*
Sensitizer	The default incubation time is 20 minutes. Select an incubation time from 4-20 minutes.*
Optimize Silver Intensity (Reticulum II Silver A)	The default incubation time is 20 minutes. Select an incubation time from 8-20 minutes: <ul style="list-style-type: none"> 8 minutes, lighter staining of fibers

Staining Procedure	S Reticulum II
Protocol Step	Method
	<ul style="list-style-type: none"> 20 minutes, darker staining of fibers
3 Rinse Cycles (optional)	The default is 1 Rinse Cycle. Select 3 Rinse Cycles to reduce non-specific background staining if necessary.
No Toner (optional)	The default is to apply Toner.
Optimize Counterstain Intensity (Nuclear Fast Red Counterstain) (optional)	The default is 4 minutes. Select an incubation time from 4-16 minutes.* <ul style="list-style-type: none"> 4 minutes, lighter counterstain 16 minutes, darker counterstain
No Counterstain (optional)	Select No Counterstain, if Nuclear Fast Red Counterstain is not desired.

* To adjust staining preferences, increment the stain temperature and incubation time one parameter at a time.

Recommended Post-Instrument Processing

1. Rinse slides in two changes of 95% ethanol to remove the leftover solution, followed by three changes of 100% ethanol.
2. Dehydrate slides in three changes of 100% xylene.
3. Coverslip with permanent mounting media.

Compatible with the VENTANA HE 600 system coverslipping protocol. For further instructions, refer to the VENTANA HE 600 system User Guide.

QUALITY CONTROL PROCEDURE

An example of a positive control material would be FFPE human tissue with defined reticular fibers (liver, spleen, or lymph nodes).¹ Control tissue should be fresh autopsy, biopsy, or surgical specimen prepared or fixed as soon as possible in a manner identical to test sections. Such tissues should monitor all steps of the analysis, from tissue preparation through staining.

Use of a tissue section fixed or processed differently from the test specimen provides control for all reagents and method steps except fixation and tissue processing. The cellular components of other tissue elements may serve as the negative control.

Optimal laboratory practice is to include a positive control section on the same slide as the test tissue. This helps identify any failures applying reagents to the slide. Control tissue may contain both positive and negative staining elements and serve as both the positive and negative control.

The control tissue must be tested with each run.

Known positive tissue controls should only be utilized for monitoring the correct performance of processed tissues and test reagents, not as an aid in formulating a specific diagnosis of patient samples.

If the positive tissue components fail to demonstrate positive staining, results with the test specimens should be considered invalid. If the negative components demonstrate positive staining, results with patient specimens should also be considered invalid.

Unexplained discrepancies in control results should be referred to the local support representative immediately. If quality control results do not meet specifications, patient results are invalid. The cause must be identified and corrected, and the patient samples repeated.

STAINING INTERPRETATION / EXPECTED RESULTS

The Reticulum II Staining Kit is tested to show reticular fibers.

- Reticular fibers: black
- Background: pink

SPECIFIC LIMITATIONS

Only positively charged microscope slides have been used and validated for this assay. Only hydrochloric acid, EDTA and formic acid decalcification techniques on bone marrow have been used and validated for this assay. Decalcification time was not specified.

PERFORMANCE CHARACTERISTICS

ANALYTICAL PERFORMANCE

Staining tests for sensitivity, specificity, and precision were conducted and the results are listed below.

Sensitivity and Specificity

Analytical sensitivity and specificity for normal and diseased tissue cases was evaluated. All evaluated tissue cases (125/125) passed for acceptable staining as shown in Table 4 and Table 5.

Table 4. Sensitivity/Specificity of Reticulum II Staining Kit was determined by testing the following FFPE normal tissues.

Tissue	# passed / # tested
Liver	10 / 10
Spleen	8 / 8
Lymph node	11 / 11
Bone marrow (unspecified decalcification)	6 / 6
Bone marrow (decalcified with hydrochloric acid)	19 / 19
Bone marrow (decalcified with EDTA)	13 / 13
Bone marrow (decalcified with formic acid)	10 / 10

Table 5. Sensitivity/Specificity of Reticulum II Staining Kit was determined by testing the following FFPE diseased tissues.

Tissue	# passed / # tested
Hepatocellular adenoma (liver)	5 / 5
Hepatocellular carcinoma (liver)	16 / 16
Cirrhosis (liver)	12 / 12
Hyperplasia (pituitary gland)	5 / 5
Adenoma (pituitary gland)	5 / 5
Myelofibrosis (bone marrow, unspecified decalcification)	5 / 5

Precision

Precision of the Reticulum II Staining Kit was determined across multiple runs, days, instruments, and reagent lots using multiple cut slides from 5 normal liver, 5 hepatocellular carcinoma (liver) and 2 normal bone marrow. All acceptance criteria were fully met. Precision studies were performed for the according to Table 6.

Table 6. Precision slide studies for Reticulum II Staining Kit.

Parameters Tested	# of conditions	# passed / # tested
Run to Run	3 runs, same day	54 / 54
Day to Day	5 days	90 / 90
Instrument to Instrument	3 instruments	54 / 54
Intra Run	same day, same instrument	54 / 54
Lot to lot	3 lots	54 / 54

The results demonstrated no significant difference in staining intensity among the slides.

TROUBLESHOOTING

1. Section thickness may affect quality and intensity of staining. If staining is inappropriate, contact your local support representative for assistance.
2. Necrotic or autolyzed tissue may exhibit nonspecific staining.
3. If the positive control is negative, tissue may have been improperly collected, fixed, or deparaffinized. Follow the proper procedure for collection, storage, and fixation.
4. If the positive control is negative, check that the slide has the proper barcode label. If the slide is labeled properly, check the other positive controls from the same run to determine if the controls were properly stained.
5. If excessive background staining occurs: incomplete paraffin removal could cause staining artifacts or no staining. If all paraffin is not removed from the slide, repeat the staining run using the extended deparaffinization option, if available.
6. If excessive background staining occurs with bone marrow section: insufficient rinses after decalcification could cause staining artifacts.
7. If tissue sections wash off the slide, confirm the slides are positively charged.
8. Extended stay of the slides on-board the instrument after run completion may affect quality and intensity of the staining. If the staining is inappropriate, remove slides promptly at the end of the run and proceed to post-instrument processing.
9. For corrective action, refer to the Instructions for Use section, the instrument User Guide or contact your local support representative.

REFERENCES

1. Carson FL, Cappellano C. Histotechnology: A Self-Instructional Text, 5th edition. American Society for Clinical Pathology Press; 2020, 2022.
2. Grizzle WE. Theory and Practice of Silver Staining in Histopathology. Journal of Histotechnology. 1996;19(3):183-195.
3. Bancroft JD, Layton C. Connective and Other Mesenchymal Tissues with Their Stains. In: Bancroft's Theory and Practice of Histological Techniques. 2019:153-175.
4. Ishibashi H, Nakamura M, Komori A, et al. Liver Architecture, Cell Function, and Disease. Semin Immunopathol. 2009;31(3):399-409.
5. Ferrell L. Liver Pathology: Cirrhosis, Hepatitis, and Primary Liver Tumors. Update and Diagnostic Problems. Mod Pathol. 2000;13(6):679-704.
6. Boyd A, Cain O, Chauhan A, et al. Medical Liver Biopsy: Background, Indications, Procedure and Histopathology. Frontline Gastroenterol. 2020;11(1):40-47.
7. Kuter DJ, Bain B, Mufti G, et al. Bone Marrow Fibrosis: Pathophysiology and Clinical Significance of Increased Bone Marrow Stromal Fibres. British Journal of Haematology. 2007;139(3):351-362.
8. Zahr AA, Salama ME, Carreau N, et al. Bone Marrow Fibrosis in Myelofibrosis: Pathogenesis, Prognosis and Targeted Strategies. Haematologica. 2016;101(6):660-671.
9. Ushiki T. Collagen Fibers, Reticular Fibers and Elastic Fibers. A Comprehensive Understanding from a Morphological Viewpoint. Arch Histol Cytol. 2002;65(2):109-126.
10. Occupational Safety and Health Standards: Occupational exposure to hazardous chemicals in laboratories. (29 CFR Part 1910.1450). Fed. Register.
11. Directive 2000/54/EC of the European Parliament and Council of 24 June 2000 on the protection of workers from risks related to exposure to biological agents at work.

NOTE: A point (period/stop) is always used in this document as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Symbols

Ventana uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard (for USA: see elabdoc.roche.com/symbols for definition of symbols used):



Global Trade Item Number



Unique Device Identification



Indicates the entity importing the medical device into the European Union

REVISION HISTORY

Rev	Updates
K	Updates to the Warnings and Precautions, References, and Symbols section.

INTELLECTUAL PROPERTY

VENTANA, BENCHMARK, VENTANA HE and the VENTANA logo are trademarks of Roche.

All other trademarks are the property of their respective owners.


© 2023 Ventana Medical Systems, Inc.

CONTACT INFORMATION



Ventana Medical Systems, Inc.
1910 E. Innovation Park Drive
Tucson, Arizona 85755
USA
+1 520 887 2155
+1 800 227 2155 (USA)

www.roche.com



Roche Diagnostics GmbH
Sandhofer Strasse 116
D-68305 Mannheim
Germany
+800 5505 6606

