# Accutrend Cholesterol cobas®

<span>REF</span>	<span>▽</span>	<span>SYSTEM</span>
11418262012	25	
11418262171	25	Accutrend Plus
11418262370	25	
11418254223	5	

**English**

**Intended use**

The Accutrend Cholesterol test strips are intended for quantitative determination of cholesterol in fresh or lithium-heparinised fresh capillary blood. High cholesterol levels aid to identify clinical conditions including those that are associated with increased risk for cardiovascular disease. Use only with the following meter: Accutrend Plus.

Accutrend Cholesterol is intended for near-patient testing as well as for self-testing.

**Summary**

The determination of cholesterol provides the physician with information on the individual concentration of circulating lipids, where high levels of cholesterol aid to identify conditions that are associated with increased risk for cardiovascular disease. Self-monitoring is not to be seen as a substitute for regular monitoring by a physician. The results of self-monitoring should be noted down and discussed with the physician. The Accutrend Cholesterol test strips can be used in non-critical healthcare environments such as hospital wards and general practitioners' offices.

**How the test works**

Each test strip has a test area containing detection reagents. When blood is applied, a chemical reaction takes place and the test area changes colour. The meter records this change in colour and converts the measurement signal to the displayed result using the data previously entered via the code strip.

**Preparing and performing the measurement**

In order to avoid erroneous measurements, use only the below-specified meters for determinations with Accutrend Cholesterol test strips. Please follow the instructions given in this Method Sheet and in the Operator's Manual of your meter.

**Precautions and warnings**

For *in vitro* diagnostic use.

All components of the pack can be discarded along with household waste or, if used in a laboratory or doctor's office, according to relevant local guidelines. Never allow liquids or disinfectants to get into the vial. The stopper contains a non-toxic silicone-based desiccant. If the desiccant is inadvertently ingested, please drink plenty of water!

**Storage and stability**

Store at 2-30 °C up to the stated expiration date. Do not freeze. Store away from extreme temperatures.

Do not use the test strip after the specified expiration date.

Tightly re-cap the container immediately after removing a test strip.

Perform cholesterol determinations at 18-30 °C.

The test strip must be used within 1 minute after removal from the vial.

Important: If you use more than one type of Accutrend test strip at the same time, make sure that you do not mix up the stoppers. If the container stoppers are mixed up, the specified shelf-life can no longer be guaranteed.

All components of the pack can be discarded along with household waste or, if used in a laboratory or doctor's office, according to relevant local guidelines.

**Materials provided**

- Test strips and 1 code strip

**Materials required (but not provided)**

- [REF] 11418289190, Accutrend Control CH 1 control solution
- [REF] 05050472xxx (mg/dL) Accutrend Plus
- [REF] 05050499xxx (mmol/L) Accutrend Plus
- Lancing device

**Assay**

For optimum performance of the system follow the directions given in this document for the corresponding meter. Refer to the appropriate meter manual for test-specific instructions. The system requires a hanging drop of blood. Use only the above-specified meters for determinations with these test strips to avoid erroneous results.

**Checks required every time you perform a test**

Before taking a measurement, the test strip must be checked to see that it is still fit for use. If the window on the back of the test strip has taken on a bluish colour, there is a risk that the cholesterol values displayed are too high. If this occurs, use a new test strip. Every time you perform a test, make sure the test area is completely covered with blood. Follow also the instructions given in the operating instructions of your meter. Occasionally, the round window on the back of the test strip may be mottled; this uneven coloration is sample-related and does not interfere with the measurement of the Accutrend Cholesterol test strip.

**Coding**

Whenever a new pack of test strips is opened your meter must be recoded with the code strip provided. The procedure is described in the meter manual. If the code number displayed does not match the vial you are using, the meter will not perform a measurement. Always keep the code strip until the last test strip of a package is used.

Store the code strip outside the test strip container. If the code strip is stored in the container, the test strips may be damaged.

**Quality control**

For quality control, use Accutrend Control CH 1.

The acceptable control range is given in the enclosed value sheet.

The target value corresponds to the mean of the lower and upper value of this control range.

***For healthcare professionals:***

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

**Limitations - interference**

The cholesterol determination may be affected by:

- intravenous infusion of ascorbic acid (vitamin C)
- bilirubin conjugated bilirubin > 6.2 mg/dL (105.4 µmol/L) or unconjugated bilirubin > 15.3 mg/dL (260.1 µmol/L, e.g. in hepatitis)
- hemoglobin plasma levels > 800 mg/dL
- haematocrit values below 30 % or above 55 %
- albumin > 40000 mg/L
- Methyl-L-Dopa > 15.0 mg/L

**Limits and ranges**

	mmol/L	mg/dL
Linear measuring range	3.88-7.76	150-300
Detection limit (lowest value displayed)	3.88	150

**Expected values**

Cholesterol values less than 5.2 mmol/L (200 mg/dL) are considered normal.

In some patient groups, the reference ranges may not be within the normal range. You should discuss this with your physician.<sup>1</sup>

**Additional information for healthcare professionals**

**Test principle**

Enzymatic cleavage of the cholesterol ester into fatty acids and cholesterol; oxidation of cholesterol to cholestenone with the simultaneous formation of hydrogen peroxide which oxidises an indicator to its blue radical cation.<sup>2,3</sup>

**Applications**

For the determination of cholesterol in fresh capillary blood in diabetics and non-diabetics; for self-monitoring; for the early detection of a risk of atherosclerosis; for monitoring treatment with lipid-lowering drugs; for screening purposes. If measurements are performed on several persons, as may be the case in a doctor's office or for screening programmes, the fresh capillary blood may also be applied with lithium-heparin-coated capillary pipettes (volume 15 µL to 50 µL).

**User Training**

The operation of the Accutrend Plus meters and the execution of the test is described in the respective meter manuals. No dedicated training is required.

**Precautions and warnings**

For *in vitro* diagnostic use for health care professionals. Exercise the normal precautions required for handling all laboratory reagents.

Infectious or microbial waste:

Warning: handle waste as potentially biohazardous material. Dispose of waste according to accepted laboratory instructions and procedures.

Environmental hazards:

Apply all relevant local disposal regulations to determine the safe disposal.

Safety data sheet available for professional user on request.

Health care professionals performing tests on more than one patient must be aware that there is a potential risk of infection. Any object coming into contact with human blood is a potential source of infection.<sup>4</sup>

**Calibration**

Each test strip lot of the Accutrend Cholesterol test is calibrated against the CHOL2 test (CHOD-PAP) on **cobas c 501** lab analyzer. Calibration has been performed with native human samples. The Accutrend Plus meter automatically reads the lot-specific calibration data from the code strip, eliminating the need for calibration by the user.

**Reagents**

Components per test:

CHE (microorganims) 1.44 U; CHOD (Nocardia erythropolis) 0.12 U; POD (horseradish) 0.20 U; 3,3',5,5'-tetramethylbenzidine 8.5 µg; non-reactive components 1.4 mg.

**Limitations - interference**

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

**Specific performance data**

**Precision**

*Repeatability*

In investigations with EDTA venous blood in a concentration range between 4.11 mmol/L (160 mg/dL) and 7.37 mmol/L (287 mg/dL) typical measurement series showed a CV (coefficient of variation) of 1.35 to 3.75 %.

*Intermediate precision*

The intermediate precision was determined in a test series with a control solution. At concentrations of approximately 4.60 mmol/L (179 mg/dL) and approximately 6.19 mmol/L (241 mg/dL) the CV was typically between 2 % and 3 %.

**Method comparison**

In method comparisons with capillary blood ≥ 95 % of all values measured on 3 lots were found within ±20 % of the reference method (CHOD-PAP, Roche Diagnostics).The mean systematic difference were within 6.3 %.

For further information, please refer to the appropriate Operators Manual for the instrument concerned, and the Method Sheets of all necessary components. A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

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- Expert Panel on Detection, E. and A. Treatment of High Blood Cholesterol in, Executive Summary of The Third Report of The National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, And Treatment of High Blood Cholesterol In Adults (Adult Treatment Panel III). JAMA, 2001, 285(19): p. 2486-97.
- European Atherosclerosis Society. Prevention of coronary heart disease: scientific background and new clinical guideline. Nutr Metab Cardiovasc Disc 1992;2:129.
- Gottschling HD et al. Eur J Clin Chem Clin Biochem 1995;33: 373-383
- Clinical and Laboratory Standards Institute: Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline - Third Edition; CLSI document M29-A3, 2005.

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
	ACCUTREND and COBAS are trademarks of Roche.
---	--

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2024, Roche Diagnostics

	0123
---	------

	Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim www.roche.com
	+800 5505 6606