



cobas[®] eplex

respiratory pathogen panel

IVD

Package Insert

For in vitro Diagnostic Use
For Professional Laboratory Use Only

P/N: 09554998001



GenMark Diagnostics, Inc.
5964 La Place Court
Carlsbad, CA, 92008
USA
+1 760 448 4300

TABLE OF CONTENTS

| | |
|--|----|
| Intended Use..... | 3 |
| Summary and Explanation of Test..... | 3 |
| Summary of Detected Organisms..... | 4 |
| Principles of Technology..... | 6 |
| Materials Provided | 7 |
| Composition of Reagents..... | 7 |
| Reagent Storage, Stability, and Handling..... | 7 |
| Materials Not Provided..... | 7 |
| Equipment | 7 |
| Consumables..... | 8 |
| Warnings and Precautions..... | 8 |
| General..... | 8 |
| Safety | 8 |
| Laboratory | 9 |
| Specimen Collection, Handling, and Storage | 9 |
| Procedure | 9 |
| Procedural Notes..... | 9 |
| Detailed Procedure | 10 |
| Quality Control | 10 |
| Internal Controls | 10 |
| External Controls | 11 |
| Results | 12 |
| Influenza A Results..... | 13 |
| Test Reports | 14 |
| Detection Report..... | 14 |
| External Control Report | 14 |
| Summary Report..... | 14 |
| Limitations of the Procedure | 14 |
| Expected values..... | 16 |
| Performance characteristics | 18 |
| Clinical performance | 18 |
| Comparator Method..... | 18 |
| Prospective Clinical Samples | 19 |
| Prospective Clinical Performance | 19 |
| Retrospective Clinical Samples | 21 |
| Retrospective Clinical Performance..... | 22 |
| Contrived Sample Performance..... | 23 |
| Clinical and Contrived Sample Performance by Target | 23 |
| Co-detections in Prospective Clinical Samples..... | 27 |
| Clinical Study cobas® eplex system Performance | 29 |
| Analytical Performance Characteristics | 29 |
| Limit of Detection..... | 29 |
| Analytical Reactivity (Inclusivity) | 31 |
| Supplemental Analytical Reactivity (Inclusivity) for Influenza A | 34 |
| Analytical Specificity (Cross-Reactivity and Exclusivity) | 38 |
| Reproducibility | 41 |
| Samples with Co-Detected Organisms | 45 |
| Sample Matrix Equivalency | 45 |
| Interfering Substances..... | 46 |
| Carryover and Cross-contamination | 47 |
| Troubleshooting | 48 |
| Technical Support..... | 48 |
| Glossary of Symbols | 49 |
| References..... | 50 |
| Document revision | 52 |
| Trademarks..... | 52 |
| Patent Information..... | 52 |

INTENDED USE

The **cobas® eplex** respiratory pathogen (RP) panel is a multiplexed nucleic acid in vitro diagnostic test intended for use on the **cobas® eplex** system for the simultaneous qualitative detection and identification of multiple respiratory viral and bacterial nucleic acids in nasopharyngeal swabs (NPS) obtained from individuals exhibiting signs and symptoms of respiratory tract infection.

The following virus types, subtypes, and bacteria are identified using the **cobas® eplex** RP panel: adenovirus, coronavirus, human metapneumovirus, human rhinovirus/enterovirus, influenza A, influenza A H1, influenza A H1-2009, influenza A H3, influenza B, parainfluenza virus 1, parainfluenza virus 2, parainfluenza virus 3, parainfluenza virus 4, respiratory syncytial virus (RSV) A, respiratory syncytial virus (RSV) B, *Chlamydia pneumoniae*, and *Mycoplasma pneumoniae*.

The detection and identification of specific viral and bacterial nucleic acids from individuals exhibiting signs and/or symptoms of respiratory tract infection aids in the diagnosis of respiratory infection when used in conjunction with other clinical and epidemiological information. The results of this test should not be used as the sole basis for diagnosis, treatment, or other patient management decisions.

Negative results in the setting of a respiratory illness may be due to infection with pathogens that are not detected by this test, or lower respiratory tract infection that may not be detected by a nasopharyngeal swab specimen. Positive results do not rule out co-infection with other organisms; the organism(s) detected by the **cobas® eplex** RP panel may not be the definite cause of disease. Additional laboratory testing (e.g., bacterial and viral culture, immunofluorescence and radiography) may be necessary when evaluating a patient with possible respiratory tract infection.

SUMMARY AND EXPLANATION OF TEST

The **cobas® eplex** RP panel is an automated qualitative nucleic acid multiplex in vitro diagnostic test for simultaneous detection and identification of multiple respiratory viral and bacterial nucleic acids in nasopharyngeal swabs (NPS) collected in viral transport media (VTM). The test is able to detect 15 respiratory viral targets and 2 bacterial targets as summarized in Table 1. This test is performed on the **cobas® eplex** system.

Respiratory viruses and bacteria are responsible for a wide range of respiratory tract infections including the common cold, influenza, and croup, and represent the most common cause of acute illness. Disease severity can be especially high in the young, the immunocompromised, and elderly patients. Respiratory infections cause more doctor visits and absences from school and work than any other illness.¹ Influenza viruses have a peak season in the winter months in the northern hemisphere and the severity of the flu season varies each year based on the particular strain or strains that are in circulation and how effective the vaccine is for that year.² Globally, seasonal influenza results in about 3-5 million severe cases and 250,000 – 500,000 deaths annually.³

Influenza-like illness is a nonspecific respiratory illness characterized by fever, fatigue, cough, and other symptoms. The majority of influenza-like illnesses are not caused by influenza but by other viruses (e.g., rhinovirus, respiratory syncytial virus, adenovirus, and parainfluenza virus).⁴ Less common causes of influenza-like illness include bacteria such as *Chlamydia pneumoniae* and *Mycoplasma pneumoniae*.⁴

Table 1: Targets Detected by the cobas® eplex RP panel

| Target | Classification (Genome Type) | Seasonal Prevalence* | Most Commonly Infected Demographic |
|--------------------------------------|------------------------------|--|---|
| Adenovirus (A-F) | Adenovirus (DNA) | Late winter to early summer ⁵ | All ages, immunocompromised ⁶ |
| Coronavirus (229E, HKU1, NL63, OC43) | Coronavirus (RNA) | Winter, spring ⁷ | All ages ⁷ |
| Human Metapneumovirus | Paramyxovirus (RNA) | Winter ⁸ | Children, elderly, immunocompromised ⁹ |
| Human Rhinovirus/Enterovirus | Picornavirus (RNA) | Fall, spring ¹⁰ /Summer ¹¹ | All ages, immunocompromised ¹⁰⁻¹² |
| Influenza A | Orthomyxovirus (RNA) | Winter ³ | All ages ³ |
| Influenza A H1 | | | |
| Influenza A H1-2009 | | | |
| Influenza A H3 | | | |
| Influenza B | | | |
| Parainfluenza Virus 1 | Paramyxovirus (RNA) | Fall ¹³ | All ages ¹⁴ |
| Parainfluenza Virus 2 | | Fall, early winter ¹³ | |
| Parainfluenza Virus 3 | | Spring, summer ¹³ | |
| Parainfluenza Virus 4 | | Fall, early winter ¹³ | |
| Respiratory Syncytial Virus A | Paramyxovirus (RNA) | Winter ^{15,16} | Infants, children, older adults ^{15,16} |
| Respiratory Syncytial Virus B | | | |
| <i>Chlamydia pneumoniae</i> | Bacterium (DNA) | No peak season ¹⁷ | All ages, most common in children ¹⁷ |
| <i>Mycoplasma pneumoniae</i> | Bacterium (DNA) | Late summer, fall ¹⁸ | Children, young adults ¹⁸ |

* Based on northern hemisphere seasons

SUMMARY OF DETECTED ORGANISMS

Adenovirus: Adenoviruses are non-enveloped DNA viruses that include seven human species (A - G) and more than 60 serotypes.¹⁹ Adenovirus species B, C, and E are frequently associated with upper respiratory infections; infections are common in children, and outbreaks often occur in crowded environments, such as military barracks.^{6,20} There is no vaccine available to the general public, but the introduction of a live, oral vaccine to the US military in 2011 has reduced the incidence of adenovirus outbreaks in this population.^{6,21} Adenovirus infections generally cause mild illness but can result in severe disease in infants or in immunocompromised patients, particularly in hematopoietic stem cell transplant recipients.^{6,19} In addition to respiratory infections, adenovirus can also cause gastroenteritis, conjunctivitis, and cystitis.^{6,19} Adenovirus species A, D, and F are not typically associated with respiratory infections.

Coronavirus: Human coronaviruses usually cause mild to moderate upper respiratory infections but can cause significant disease in the elderly, young children, and immunocompromised individuals.²²⁻²⁴ Infection with coronaviruses 229E, HKU1, NL63, and OC43 is common worldwide.

Human Metapneumovirus: Human metapneumovirus is a member of the paramyxovirus family and is closely related to RSV.⁹ Metapneumovirus has been identified as an important respiratory pathogen in young children and is the second most common virus identified in pediatric respiratory tract infections.⁹ Illness is more severe in children who are immunocompromised or have underlying conditions, such as

HIV or cardiac disease; it can also cause more severe disease in immunocompromised adults, especially those with COPD (chronic obstructive pulmonary disease), asthma, cancer, or in transplant patients.²⁵

Human Rhinovirus and Enterovirus: Rhinovirus and enterovirus are closely related RNA viruses in the *Picornaviridae* family.¹² There are more than 100 different serotypes that all share high sequence homology.²⁶ Rhinovirus causes up to 80% of all cases of the common cold worldwide and is more common in children than adults. It is the cause of a significant number of mild upper respiratory tract infections throughout the year, especially during the spring and fall seasons.^{10,27} Most infections are mild, but rhinovirus has been associated with severe infections in at-risk populations including young children, the elderly, immunocompromised patients, and those with asthma.^{10,11}

There are 62 non-polio enteroviruses that can cause disease in humans.¹² Enterovirus primarily infects the gastrointestinal tract but can also cause respiratory illness, which is generally mild, like the common cold, but can result in serious complications, especially in infants.¹² A 2014 outbreak of enterovirus D68 (EV-D68) resulted in severe respiratory infections, some of which were fatal.²⁸

Influenza virus: There are three types of influenza viruses: A, B, and C.³ In the northern hemisphere, influenza A and B circulate during the winter months causing seasonal epidemics most years; influenza C infections are less common and not believed to cause epidemics.^{3,29} Both influenza A and B mutate, and the impact of influenza varies from year to year depending on the severity of the changes and effectiveness of influenza vaccines.³⁰ The two most common influenza A subtypes infecting humans are H1N1 (including the 2009 Pandemic H1N1 variant) and H3N2, and prevalence varies annually.²⁹ Other rare influenza A subtypes also known to infect humans, such as H5N1 (avian influenza) and H3N2v, can cause severe illness and, in some cases, death.³¹ Influenza is easily spread from person to person and those most at risk for complications from infection include infants and children, the elderly, and anyone who is immunocompromised or who has co-morbidities such as heart or lung disease.³²

Influenza A 2009 H1N1: During the 2009 - 2010 influenza season, a new strain of influenza A, now known as 2009 H1N1 became the dominant circulating virus, accounting for approximately 95% of reported influenza infections.³³ This strain replaced the H1N1 virus that was previously circulating in humans and is common in both Europe and the U.S.^{3,29}

Parainfluenza Virus: The parainfluenza viruses are members of the paramyxovirus family that commonly cause respiratory infections in children.³⁴ Prevalence of parainfluenza viruses is seasonal and varies by type; most infections are mild and self-limited, but parainfluenza virus can cause life threatening pneumonia in immunocompromised patients, such as those with cystic fibrosis or transplant recipients.³⁵

Respiratory Syncytial Virus: RSV is the most common cause of pediatric viral respiratory infections.⁹ Infection with RSV can occur at any age, and those most at risk for complications and more severe disease are the very young, especially premature infants, the elderly, and anyone with a weakened immune system.³⁶ There are two types of respiratory syncytial virus, RSV A and B. Infections with RSV A are thought to be more severe than infections with RSV B.^{16,37}

***Chlamydia pneumoniae* (formerly known as *Chlamydophila pneumoniae*):** *Chlamydia pneumoniae* is a common cause of upper respiratory infections including atypical pneumonia.³⁸ *C. pneumoniae* is transmitted person-to-person by respiratory secretions and outbreaks are common in close contact settings.¹⁷ Infection severity can be mild or result in more severe disease, particularly in high risk populations such as people with heart or lung disease, diabetes, and the elderly.^{17,39} The true prevalence of *C. pneumoniae* infections is unknown, but the use of molecular diagnostics has improved detection of this organism, as it is difficult to identify using traditional laboratory methods.³⁸

***Mycoplasma pneumoniae*:** *Mycoplasma pneumoniae* is a bacterium lacking a cell wall and is a major cause of respiratory disease.¹⁸ *M. pneumoniae* is transmitted person-to-person by respiratory droplets and is a common cause of atypical, or walking pneumonia.⁴⁰ *M. pneumoniae* is frequently undiagnosed but is estimated to be involved in up to 30% of respiratory infections.¹⁸ Infection often results in mild illness such as tracheobronchitis, or a chest cold, and is most prevalent in young adults and school-aged children.^{18,40} Outbreaks of *M. pneumoniae* occur mostly in crowded environments, like schools, college dormitories, military barracks, and nursing homes.⁴⁰

PRINCIPLES OF TECHNOLOGY


The **cobas® eplex** system automates all aspects of nucleic acid testing including extraction, amplification, and detection, combining electrowetting and eSensor technology in a single-use cartridge. eSensor technology is based on the principles of competitive DNA hybridization and electrochemical detection, which is highly specific and is not based on fluorescent or optical detection.

Electrowetting, or digital microfluidics, uses electrical fields to directly manipulate discrete droplets on the surface of a hydrophobically coated printed circuit board (PCB). Sample and reagents are moved in a programmable fashion in the **cobas® eplex** cartridge to complete all portions of the sample processing from nucleic acid extraction to detection.

A sample is loaded onto the **cobas® eplex** cartridge and nucleic acids are extracted and purified from the specimen via magnetic solid phase extraction. For RNA targets, a reverse transcription step is performed to generate complementary DNA from the RNA, followed by PCR to amplify the targets. Exonuclease digestion creates single-stranded DNA in preparation for eSensor detection.

The target DNA is mixed with ferrocene-labeled signal probes that are complementary to the specific targets on the panel. Target DNA hybridizes to its complementary signal probe and capture probes, which are bound to gold-plated electrodes, as shown below in Figure 1. The presence of each target is determined by voltammetry which generates specific electrical signals from the ferrocene-labeled signal probe.

Figure 1: Hybridization complex. Target-specific capture probes are bound to the gold electrodes in the eSensor microarray on the **cobas® eplex** cartridge. The amplified target DNA hybridizes to the capture probe and to a complementary ferrocene-labeled signal probe. Electrochemical analysis determines the presence or absence of targets using voltammetry.



MATERIALS PROVIDED

Table 2: cobas® eplex respiratory pathogen panel Kit Contents

| Product | Item number | Components (quantity) | Storage |
|---|-------------------|--|---------|
| cobas® eplex respiratory pathogen panel | Roche: 9554998001 | cobas® eplex respiratory pathogen panel cartridge (12) | 2–8 °C |
| | | Sample Delivery Device – RP panel (12) | 2–8 °C |

Safety Data Sheets (SDS) for all reagents provided in this kit may be obtained at: <https://dialog.roche.com>. For paper copies, please reach out to your local affiliate: https://www.roche.com/about/business/roche_worldwide.htm.

COMPOSITION OF REAGENTS

Table 3: Composition of Reagents on the cobas® eplex RP panel Cartridges

| Composition of Reagents on the cobas® eplex RP panel Cartridges | |
|---|--|
| 2-(N-morpholino)ethanesulfonic acid (MES) | NaH ₂ PO ₄ /NaHPO ₄ |
| 6-mercapto-1-hexanol | NaN ₃ |
| Acetonitrile | PEG 8000 |
| Calcium Chloride | Phenol Red |
| Cysteamine HCl | Polydimethylsiloxane |
| Dynol-604 | Ribonuclease inhibitor |
| EDTA | SDS, pH adjusted with HCl |
| EGTA | Sodium perchlorate |
| Ethanol | Sorbitane trioleate |
| Glycerol | Super Q water |
| Guanidinium Hydrochloride | Trehalose |
| Lithium Dodecyl Sulfate | Trimethyl terminated, 5cSt |
| Magnesium Chloride (MgCl ₂) | Tris-HCl |
| MTG, pH adjusted with sodium hydroxide + Tween-20 | Tween-20 |
| NaCl | Urea |

REAGENT STORAGE, STABILITY, AND HANDLING

- Store the cobas® eplex RP panel kit components at 2–8 °C.
- Do not use RP panel kit components beyond the expiration date.
- Do not open a cartridge pouch until you are ready to perform testing.

MATERIALS NOT PROVIDED

Equipment

- cobas® eplex system and Software
- Pipettes calibrated to deliver 200 µL
- Vortex mixer
- Printer (optional) - See cobas® eplex User Assistance Manual for compatibility guidelines

Consumables

- Pipette tips, aerosol resistant, RNase/DNase-free
- Disposable, powder free gloves
- 10% bleach for decontamination of appropriate surfaces
- 70% ethanol or isopropyl alcohol

WARNINGS AND PRECAUTIONS

General

- For in vitro diagnostic use, by laboratory professionals.
- For prescription use only.
- A trained healthcare professional should carefully interpret the results from the **cobas® eplex** RP panel in conjunction with a patient's signs, symptoms, and results from other diagnostic tests.
- Positive results do not rule out co-infection with other viruses or bacteria. The agent detected may not be the definite cause of disease. The use of additional laboratory testing (e.g., bacterial and viral culture, immunofluorescence, and radiography) and clinical presentation must be taken into consideration in the diagnosis of respiratory infection.
- Do not reuse **cobas® eplex** RP panel kit components.
- Do not use reagents beyond the expiration date printed on the labeling.
- Do not use a reagent that is damaged.
- Follow the procedure as described in this package insert. Read all instructions before starting the test.

Safety

- Handle all specimens and waste materials as if they were capable of transmitting infectious agents in accordance with Universal Precautions. Observe safety guidelines such as those outlined in CDC/NIH *Biosafety in Microbiological and Biomedical Laboratories*, CLSI Document M29 *Protection of Laboratory Workers from Occupationally Acquired Infections*, or other appropriate guidelines.
- Follow routine laboratory safety procedures for handling of reagents (e.g., do not pipette by mouth, wear appropriate protective clothing and eye protection).
- Follow your institution's safety procedures for handling biological samples.
- If infection with a novel influenza A virus is suspected based on current clinical and epidemiological screening criteria recommended by public health authorities, specimens should be collected with appropriate infection control precautions for novel virulent Influenza viruses and sent to state or local health department for testing. Viral culture should not be attempted in these cases unless a BSL-3+ facility is available to receive and culture specimens.
- Dispose materials used in this test, including reagents, specimens, and used vials, in accordance with all federal, state, and local regulations.
- Do not stick fingers or other objects inside the **cobas® eplex** system bays.
- Wash hands thoroughly with soap and water after handling reagents. Launder contaminated clothing prior to re-use.
- Do not puncture or pierce reagent blisters on the **cobas® eplex** cartridge. Reagents may cause irritation to skin, eyes, and respiratory tract. Harmful if swallowed or inhaled. Contains oxidizing liquids.
- The **cobas® eplex** RP panel cartridge contains chemicals that are classified as hazardous. Review the Safety Data Sheet (SDS) before use, and in cases of exposure, refer to the SDS for more information. Safety Data Sheets (SDS) are available on request from your local Roche representative or can be accessed via eLabDoc.

Laboratory

- Contamination of the sample may occur if laboratory personnel processing the sample are infected with common respiratory pathogens. To avoid this, specimens should be processed in biosafety cabinets. If a biosafety cabinet is not used, a splash shield or face mask should be used when processing samples.
- A biosafety cabinet that is used for viral or bacterial culture should not be used for sample preparation.
- Samples and cartridges should be handled and/or tested one at a time. To mitigate the risk of sample-to-sample contamination, change gloves after dispensing sample into the cartridge.
- Thoroughly decontaminate the lab and all equipment with 10% bleach followed by 70% ethanol or isopropyl alcohol (or equivalent) prior to processing a specimen.
- Contamination of the sample may occur if the sample is loaded in an area where PCR amplicons for respiratory pathogens are generated. Avoid loading sample in areas that are potentially contaminated with PCR amplicon.

SPECIMEN COLLECTION, HANDLING, AND STORAGE

Nasopharyngeal Swab Collection – Nasopharyngeal swab specimen collection should be performed according to standard technique and placed in viral transport media.

Minimum Sample Volume – 200 µL nasopharyngeal swab specimen in viral transport media is required for testing.

Transport and Storage – Clinical specimens can be stored at room temperature (15–30 °C) for up to 12 hours or refrigerated at 2-8 °C for up to 10 days after collection in transport media. Specimens can also be stored at -20 °C or -80 °C for 12 months with up to 2 freeze/thaw cycles.

PROCEDURE

Procedural Notes

- All frozen samples should be thawed completely before testing.
- Samples should be nasopharyngeal swabs in viral transport media.
- Reagents and cartridge can be used immediately upon removal from 2-8 °C storage. There is no need to equilibrate to room temperature before use.
- Once cartridge is removed from foil pouch, it should be used within 2 hours. Do not open the cartridge pouch until the sample is ready to be tested.
- Once the sample is loaded into the **cobas® eplex** RP panel cartridge, the sample should be tested within 2 hours.
- Do not re-use cartridges or Sample Delivery Devices.
- Do not use a Sample Delivery Device that is empty. Visually verify that the vial contains liquid prior to use by tapping vial on the benchtop. Presence of liquid in the vial indicates that the vial can be used for testing. To prevent damage to the Sample Delivery Device, do not centrifuge the Sample Delivery Device.
- Use a new, sterile pipette tip for loading each sample.
- Do not insert a wet cartridge into the **cobas® eplex** system. If the cartridge or sample has leaked, dispose of cartridge in accordance with all federal, state, and local regulations.
- Samples should be transferred into the **cobas® eplex** RP panel cartridge in an amplicon-free, clean environment.
- Samples, consumables, and lab areas should be protected from aerosol or direct contamination with amplicon. Decontaminate laboratory areas and affected equipment with 10% bleach followed by 70% ethanol or isopropyl alcohol (or equivalent).

- Samples and cartridges should be handled and/or tested one at a time. To mitigate the risk of sample-to-sample contamination, change gloves after dispensing sample into the cartridge.
- Specimens should be processed in biosafety cabinets. If a biosafety cabinet is not used, a splash shield or face mask should be used when processing samples.
- Dispose materials used in this test, including reagents, specimens, and used vials, in accordance with all regulations.

Detailed Procedure

1. Decontaminate the clean area used for setting up the **cobas® eplex** RP panel with 10% bleach followed by 70% ethanol or isopropyl alcohol (or equivalent).
2. Remove one RP panel cartridge pouch and one Sample Delivery Device from kit packaging.
3. Open the RP panel cartridge pouch.
4. Write the accession ID or place a barcode label with accession ID on the RP panel cartridge.
5. Write the accession ID or place a barcode label with accession ID on the Sample Delivery Device.
6. Vortex the sample for 3-5 seconds.
7. Gently tap the Sample Delivery Device on the counter or benchtop surface to collect liquid that may have adhered to the sides of the vial.
NOTE: Contents of vial may adhere to side of vial and inside cap during transit. Visually verify presence of liquid inside vial after tapping vial.
8. Unscrew the purple cap from the Sample Delivery Device.
9. Use a calibrated pipette to aspirate 200 µL of sample and pipette into the Sample Delivery Device.
10. Replace purple cap on Sample Delivery Device. Ensure that cap is securely fastened on the Sample Delivery Device.
11. Vortex the Sample Delivery Device for 10 seconds.
NOTE: This step should be done immediately before loading sample onto cartridge.
12. Remove the white cover from the tip of the Sample Delivery Device cap.
13. Invert the Sample Delivery Device and dispense the entire volume by squeezing the vial and dispensing the drops into the sample loading port of the RP Panel cartridge.
NOTE: Minimize dispensing of bubbles into sample loading port.
14. Close the sample loading port by sliding the cap over the port and firmly pushing down on the cap to securely seal the sample delivery port.
NOTE: Bubbles can be present when closing the cap.
15. Scan the RP panel cartridge using the barcode reader provided with the **cobas® eplex** system.
NOTE: If an accession ID barcode label is not used, manually enter accession ID with the on-screen keyboard and scan the cartridge barcode when prompted by the **cobas® eplex** system.
NOTE: The barcode scanner will read both the accession ID barcode (if placed on the cartridge by the operator) and the 2D barcode printed on the cartridge label; however, the barcode scanner will only beep once to indicate that both barcodes have been read.
16. Insert the RP Panel cartridge into any available bay, indicated by a flashing, white LED light. The test will begin automatically when the cartridge has been inserted into the bay and the pre-run check (cartridge initialization) is completed, indicated by a blue LED light.

QUALITY CONTROL

Internal Controls

Each cartridge includes internal controls that monitor performance of each step of the testing process. A DNA control verifies extraction, amplification and detection of DNA targets, and RNA controls verify amplification and detection of RNA targets.

Each amplification reaction on the cartridge has at least one internal control and in each reaction either the internal control or a target must generate signal above the defined threshold for a valid test result. Internal control results are interpreted by the **cobas® eplex** software and displayed on **cobas® eplex** RP panel Reports as Internal Control with a result of PASS, FAIL, N/A or INVALID. Table 4 includes details on the interpretation of Internal Control results.

Table 4: Internal Control Results

| Internal Control Result | Explanation | Action |
|--------------------------------|---|---|
| PASS | The internal control or a target from each amplification reaction has generated signal above the threshold. The test was completed and internal controls were successful, indicating valid results were generated. | All results are displayed on the RP panel Detection Report. Test is valid, report results. |
| FAIL | Neither the internal control nor any target in at least one amplification reaction generates signal above the threshold. The test was completed but at least one internal control was not detected, indicating that results are not valid. | No results are displayed on the RP panel Detection Report. Test is not valid, repeat the test using a new cartridge. |
| N/A | The internal control in every amplification reaction does not generate signal above the threshold, but a target in every amplification reaction does generate signal above the threshold. The test was completed and internal controls were not successful, however detection of signal above the threshold for a target in every amplification reaction indicates valid results were generated. | All results are displayed on the RP panel Detection Report. Test is valid, report results. |
| INVALID | An error has occurred during processing that prevents analysis of signal data. The test has not successfully completed and results for this test are not valid. This is often due to an instrument or software error. | No results are displayed on the RP panel Detection Report. Test is not valid, repeat the test using a new cartridge. |

External Controls

Positive and negative external controls should be tested as part of good laboratory practice, in accordance with the appropriate accrediting organization as applicable and following the user's laboratory standard quality control procedures. Viral transport medium can be used as the negative control. Previously characterized positive samples or viral transport medium spiked with well characterized organisms can be used as the external positive control. External controls should be run in accordance with laboratory protocols and accrediting organizations, as applicable.

RESULTS

Table 5: Interpretation of Results on the cobas® eplex RP panel Detection Report

| Target Result | Explanation | Action |
|---------------------------|---|--|
| Target Detected | The test was completed successfully, and the target has generated signal above its defined threshold, and the Internal Control was reported as PASS. | All results are displayed on the RP panel Detection Report. Test is valid, report results. |
| Multiple Targets Detected | The test was completed successfully, and multiple targets have generated signal above their defined threshold, and the Internal Control was reported as PASS. | All results are displayed on the RP panel Detection Report. Test is valid, report results. Detection of more than 3 pathogens may indicate contamination. Re-test of the sample is recommended to confirm results. |
| Target Not Detected | The test was completed successfully, and the target did not generate signal above its defined threshold, and the Internal Control was reported as PASS. | All results are displayed on the RP panel Detection Report. Test is valid, report results. |
| Invalid | The test has not successfully completed, and results for this test are not valid. This is often due to an instrument or software error or failure of an internal control. | No results are displayed on the RP panel Detection Report. Test is not valid, repeat test. |

Influenza A Results

The **cobas® eplex** RP panel detects Influenza A and the H1, H1-2009, and H3 subtypes using unique assays for each. Interpretation of results for Influenza A are described in Table 6.

Table 6: Results for Influenza A

| Results for Influenza A and Subtypes | Explanation | Results on Report | Recommended Action |
|---|---|--|---|
| Influenza A detected, at least one subtype (H1, H1-2009, or H3) reported as detected. | This is an expected result. | Result reported as influenza A and influenza A subtype detected. | None |
| Influenza A detected, all subtypes (H1, H1-2009, and H3) reported as not detected | Low virus titers can result in detection of influenza A matrix without a subtype. Detection of influenza A matrix without a subtype can also indicate the presence of a novel strain. | Result reported as influenza A detected. No Influenza A subtype detected. | Re-test to confirm result. If the original result is confirmed, contact the appropriate public health authorities for additional testing. If the re-test provides a different result, test the sample a third time to ensure the accuracy of the result. |
| Influenza A detected and more than one subtype (H1, H1-2009, or H3) reported as detected. | Sample is co-infected with multiple influenza subtypes. Infection with multiple subtypes of influenza are possible but rare. A live intranasal multivalent influenza virus vaccine may cause false positive results for influenza A, A/H1, A/H3, A/H1-2009, and/or influenza B. Contamination has occurred. | Result reported as influenza A and multiple subtypes detected. | Re-test to confirm result. If the re-test result confirms the original result, it is recommended that the sample be further investigated using a different FDA-cleared influenza A subtyping assay. |
| Influenza A not detected, at least one subtype (H1, H1-2009, or H3) reported as detected. | Low virus titers can result in detection of influenza A subtype without the influenza A matrix. Detection of influenza A subtype without the influenza A matrix can also indicate the presence of a novel strain. | Influenza A (subtype) detected. Re-testing of this sample to confirm Influenza A (subtype) is recommended. Refer to package insert for additional information. | Re-test to confirm result. If the re-test result confirms the original result, the influenza A subtype is considered positive. It is recommended that the sample be further investigated using a different FDA-cleared influenza A subtyping assay and/or sending the residual sample to local public health laboratory for further testing. |

TEST REPORTS

There are several different reports that are available on the **cobas® eplex** system. Results are provided in a printable format, may be viewed electronically, or may be exported for additional analysis. Reports can be customized with account specific information such as the address, logo, and institution specific footers on each report. For more information on **cobas® eplex** reports, refer to the **cobas® eplex** User Assistance Manual.

Detection Report

The RP panel Detection Report includes the results for each individual sample run on the **cobas® eplex** system.

The Summary section indicates the overall test result and lists all detected targets in that sample. The Results section includes a list of all targets on the panel with an individual result for each. Results for each target are reported as Detected, Not Detected, or Invalid (displayed as a red **x**); results for the Internal Control are reported as PASS, FAIL, INVALID, or N/A.

External Control Report

The RP panel External Control Report is generated for an external control that has been pre-defined in the **cobas® eplex** RP panel software. For more information on defining external controls on **cobas® eplex** RP panel, refer to the **cobas® eplex** User Assistance Manual.

The Summary section indicates the overall result (Pass or Fail status) and lists all detected targets for that external control. The Results section includes a list of all panel targets with the result, expected result, and Pass/Fail status for each. Results are reported as Detected, Not Detected, or Invalid (displayed as a red **x**). A target is reported as Pass if the actual result matches the expected result (as defined for that control); a target is reported as Fail if the actual result does not match the expected result. If the actual results for each target match the expected result for each target (all targets reported as Pass), the overall result for the external control is reported as Pass in the Summary section. If the actual result for any target does not match the expected result, the overall result for the external control is reported as Fail in the Summary section.

Summary Report

The Summary Report allows the operator to use defined searchable criteria to create customized reports, using specified targets, dates, range of dates, sample, external control, test bay, or operator. For more information on creating Summary Reports, refer to the **cobas® eplex** User Assistance Manual.

LIMITATIONS OF THE PROCEDURE

- This product can be used only with the **cobas® eplex** system.
- Due to the genetic similarity between human rhinovirus/enterovirus and poliovirus, the **cobas® eplex** RP panel cannot reliably differentiate them. If a poliovirus infection is suspected, an **cobas® eplex** RP human rhinovirus/enterovirus result of Detected should be confirmed using an alternate method (e.g. cell culture).
- Due to the genetic similarity between human rhinovirus and enterovirus, this test cannot reliably differentiate them. An **cobas® eplex** RP panel Rhinovirus/Enterovirus positive result should be followed-up using an alternate method (e.g. cell culture or sequence analysis) if differentiation between the viruses is required.
- This test is a qualitative test and does not provide a quantitative value of detected organism present.
- The performance of the test has been evaluated for use with human sample material only.

- This test has not been validated for testing samples other than nasopharyngeal swab samples in viral transport media.
- The performance of this test has not been established for immunocompromised individuals.
- The performance of this test has not been established for patients without signs and symptoms of respiratory infection.
- Results from this test must be correlated with the clinical history, epidemiological data, and other data available to the clinician evaluating the patient.
- The effect of antibiotic treatment on test performance has not been evaluated.
- The performance of this test has not been established for screening of blood or blood products.
- Targets (viral and bacterial nucleic acids) may persist in vivo, independent of viral or bacterial viability. Detection of target(s) does not imply that the corresponding virus(es) or bacteria are infectious, or are the causative agents for clinical symptoms.
- The detection of viral or bacterial nucleic acid is dependent upon proper specimen collection, handling, transportation, storage, and preparation. Failure to observe proper procedures in any one of these steps can lead to incorrect results. There is a risk of false positive or false negative values resulting from improperly collected, transported, or handled samples.
- There is a risk of false negative values due to the presence of sequence variants in the viral or bacterial targets of the test, the presence of inhibitors, technical error, sample mix-up, or an infection caused by an organism not detected by the panel. Test results may be affected by concurrent antibacterial or antiviral therapy or levels of bacteria or virus in the sample that are below the limit of detection for the test. A result of No Targets Detected on the **cobas® eplex** RP panel should not be used as the sole basis for diagnosis, treatment or other patient management decisions.
- A result of No Targets Detected on the **cobas® eplex** RP panel in the setting of a respiratory illness may be due to infection with pathogens that are not detected by this test or lower respiratory tract infection that is not detected by a nasopharyngeal swab sample.
- There is a risk of false positive results due to contamination of the sample with target organisms, their nucleic acids, or amplicons. Particular attention should be given to the Laboratory precautions noted under the *Warnings and Precautions* section.
- There is a risk of false positive results due to non-specific amplification and cross-reactivity with organisms found in the respiratory tract. Erroneous results due to cross-reactivity with organisms that were not specifically evaluated or new variant sequences that emerge are possible.
- If four or more organisms are detected in a sample, retesting is recommended to confirm polymicrobial result.
- The **cobas® eplex** RP panel influenza A subtyping reagents target the influenza A hemagglutinin gene only. The **cobas® eplex** RP panel does not detect or differentiate the influenza A neuraminidase gene.
- The performance of this test has not been established for monitoring treatment of infection with any of the panel organisms.
- Positive and negative predictive values are highly dependent on prevalence. False negative test results are more likely during peak activity when prevalence of disease is high. False positive test results are more likely during periods when prevalence is moderate to low.
- Clinical performance was established when influenza A H3 and influenza A H1-2009 were the predominant influenza A viruses in circulation. When other influenza A viruses emerge, performance may vary.
- Due to the small number of positive samples collected for *Chlamydia pneumoniae* during the prospective and retrospective clinical studies, performance characteristics for *Chlamydia pneumoniae* were established primarily with contrived clinical specimens. Performance characteristics for Influenza A H1 were established using contrived clinical specimens only.
- Clinical evaluation indicates a lower sensitivity for the detection of coronavirus OC43. If infection with coronavirus OC43 is suspected, negative samples should be confirmed using an alternative method.

- The effect of interfering substances has only been evaluated for those listed in this package insert. Interference due to substances other than those described in the “Interfering Substances” section can lead to erroneous results.
- At concentrations greater than 1% weight/volume in the sample, tobramycin was found to inhibit assay performance.
- The performance of this test has not been specifically evaluated for specimens collected from individuals who recently received influenza vaccine. Recent administration of a live intranasal influenza virus vaccine may cause false positive results for influenza A, H1, H3, H1-2009, and/or influenza B.
- The **cobas® eplex** RP panel cannot differentiate variant viruses, such as H3N2v, from seasonal influenza A viruses. If variant virus infection is suspected, clinicians should contact their state or local health department to arrange specimen transport and request a timely diagnosis at a state public health laboratory.

EXPECTED VALUES

A prospective, multicenter clinical study was conducted to evaluate the clinical performance of the **cobas® eplex** RP panel in nasopharyngeal swab samples. 2462 nasopharyngeal swab samples were prospectively-collected at 8 collection sites in 2 phases from patients of all ages and genders presenting with signs and/or symptoms of respiratory infection. In the first phase from March 2013 through August 2014, 1951 samples were prospectively-collected and frozen; from September 2016 through October 2016, 511 samples were prospectively-collected and tested fresh (never frozen). The expected values of individual analytes based on **cobas® eplex** RP panel results in prospective samples for each phase are summarized in Table 7 to Table 10.

Table 7: Expected Value (As Determined by **cobas® eplex** RP panel) Summary By Age Group in the Prospective Clinical Evaluation (Phase 1: March 2013 – August 2014)

| Organism | All Ages (N=1951) n (%) | Age 0-1 (N=315) n (%) | Age >1-5 (N=250) n (%) | Age >5-21 (N=246) n (%) | Age >21-65 (N=745) n (%) | Age >65 (N=395) n (%) |
|------------------------------|-------------------------------|-----------------------------|------------------------------|-------------------------------|--------------------------------|-----------------------------|
| Adenovirus | 72 (3.7) | 31 (9.8) | 24 (9.6) | 7 (2.8) | 7 (0.9) | 3 (0.8) |
| Coronavirus | 102 (5.2) | 19 (6.0) | 18 (7.2) | 16 (6.5) | 32 (4.3) | 17 (4.3) |
| Human Metapneumovirus | 113 (5.8) | 22 (7.0) | 28 (11.2) | 6 (2.4) | 31 (4.2) | 26 (6.6) |
| Human Rhinovirus/Enterovirus | 388 (19.9) | 113 (35.9) | 94 (37.6) | 58 (23.6) | 87 (11.7) | 36 (9.1) |
| Influenza A | 110 (5.6) | 6 (1.9) | 18 (7.2) | 20 (8.1) | 49 (6.6) | 17 (4.3) |
| Influenza A H1 | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) |
| Influenza A H1-2009 | 76 (3.9) | 4 (1.3) | 13 (5.2) | 14 (5.7) | 37 (5.0) | 8 (2.0) |
| Influenza A H3 | 34 (1.7) | 1 (0.3) | 5 (2.0) | 6 (2.4) | 12 (1.6) | 10 (2.5) |
| Influenza B | 62 (3.2) | 4 (1.3) | 9 (3.6) | 10 (4.1) | 24 (3.2) | 15 (3.8) |
| Parainfluenza Virus 1 | 24 (1.2) | 4 (1.3) | 12 (4.8) | 4 (1.6) | 3 (0.4) | 1 (0.3) |
| Parainfluenza Virus 2 | 10 (0.5) | 4 (1.3) | 4 (1.6) | 0 (0.0) | 2 (0.3) | 0 (0.0) |
| Parainfluenza Virus 3 | 99 (5.1) | 31 (9.8) | 20 (8.0) | 3 (1.2) | 27 (3.6) | 18 (4.6) |
| Parainfluenza Virus 4 | 7 (0.4) | 3 (1.0) | 2 (0.8) | 1 (0.4) | 1 (0.1) | 0 (0.0) |
| RSV A | 28 (1.4) | 13 (4.1) | 6 (2.4) | 3 (1.2) | 2 (0.3) | 4 (1.0) |
| RSV B | 83 (4.3) | 33 (10.5) | 19 (7.6) | 6 (2.4) | 15 (2.0) | 10 (2.5) |
| <i>Chlamydia pneumoniae</i> | 3 (0.2) | 0 (0.0) | 0 (0.0) | 1 (0.4) | 1 (0.1) | 1 (0.3) |
| <i>Mycoplasma pneumoniae</i> | 5 (0.3) | 1 (0.3) | 1 (0.4) | 2 (0.8) | 1 (0.1) | 0 (0.0) |

Table 8: Expected Value (As Determined by cobas® eplex RP panel) Summary By Age Group in the Prospective Clinical Evaluation (Phase 2: September 2016 – October 2016)

| Organism | All Ages (N=511) n (%) | Age 0-1 (N=73) n (%) | Age >1-5 (N=75) n (%) | Age >5-21 (N=75) n (%) | Age >21-65 (N=181) n (%) | Age >65 (N=107) n (%) |
|------------------------------|------------------------------|----------------------------|-----------------------------|------------------------------|--------------------------------|-----------------------------|
| Adenovirus | 10 (2.0) | 3 (4.1) | 4 (5.3) | 1 (1.3) | 1 (0.6) | 1 (0.9) |
| Coronavirus | 8 (1.6) | 2 (2.7) | 0 (0.0) | 1 (1.3) | 4 (2.2) | 1 (0.9) |
| Human Metapneumovirus | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) |
| Human Rhinovirus/Enterovirus | 188 (36.8) | 37 (50.7) | 40 (53.3) | 33 (44.0) | 58 (32.0) | 20 (18.7) |
| Influenza A | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) |
| Influenza A H1 | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) |
| Influenza A H1-2009 | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) |
| Influenza A H3 | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) |
| Influenza B | 2 (0.4) | 0 (0.0) | 0 (0.0) | 1 (1.3) | 1 (0.6) | 0 (0.0) |
| Parainfluenza Virus 1 | 1 (0.2) | 0 (0.0) | 1 (1.3) | 0 (0.0) | 0 (0.0) | 0 (0.0) |
| Parainfluenza Virus 2 | 13 (2.5) | 3 (4.1) | 4 (5.3) | 3 (4.0) | 2 (1.1) | 1 (0.9) |
| Parainfluenza Virus 3 | 5 (1.0) | 2 (2.7) | 1 (1.3) | 1 (1.3) | 1 (0.6) | 0 (0.0) |
| Parainfluenza Virus 4 | 8 (1.6) | 1 (1.4) | 4 (5.3) | 2 (2.7) | 1 (0.6) | 0 (0.0) |
| RSV A | 8 (1.6) | 5 (6.8) | 3 (4.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) |
| RSV B | 9 (1.8) | 3 (4.1) | 4 (5.3) | 0 (0.0) | 2 (1.1) | 0 (0.0) |
| <i>Chlamydia pneumoniae</i> | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) |
| <i>Mycoplasma pneumoniae</i> | 4 (0.8) | 0 (0.0) | 1 (1.3) | 2 (2.7) | 1 (0.6) | 0 (0.0) |

Table 9: Expected Value (As Determined by cobas® eplex RP panel) Summary By Sample Collection Site in the Prospective Clinical Evaluation (Phase 1: March 2013 – August 2014)

| Organism | All Sites (N=1951) n (%) | Site 1 (N=165) n (%) | Site 2 (N=248) n (%) | Site 3 (N=350) n (%) | Site 4 (N=892) n (%) | Site 5 (N=296) n (%) |
|------------------------------|--------------------------------|----------------------------|----------------------------|----------------------------|----------------------------|----------------------------|
| Adenovirus | 72 (3.7) | 4 (2.4) | 8 (3.2) | 28 (8.0) | 23 (2.6) | 9 (3.0) |
| Coronavirus | 102 (5.2) | 8 (4.8) | 11 (4.4) | 32 (9.1) | 29 (3.3) | 22 (7.4) |
| Human Metapneumovirus | 113 (5.8) | 10 (6.1) | 23 (9.3) | 27 (7.7) | 30 (3.4) | 23 (7.8) |
| Human Rhinovirus/Enterovirus | 388 (19.9) | 27 (16.4) | 33 (13.3) | 61 (17.4) | 185 (20.7) | 82 (27.7) |
| Influenza A | 110 (5.6) | 5 (3.0) | 21 (8.5) | 48 (13.7) | 19 (2.1) | 17 (5.7) |
| Influenza A H1 | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) |
| Influenza A H1-2009 | 76 (3.9) | 3 (1.8) | 22 (8.9) | 31 (8.9) | 5 (0.6) | 15 (5.1) |
| Influenza A H3 | 34 (1.7) | 2 (1.2) | 0 (0.0) | 18 (5.1) | 12 (1.3) | 2 (0.7) |
| Influenza B | 62 (3.2) | 9 (5.5) | 9 (3.6) | 9 (2.6) | 19 (2.1) | 16 (5.4) |
| Parainfluenza Virus 1 | 24 (1.2) | 0 (0.0) | 0 (0.0) | 5 (1.4) | 2 (0.2) | 17 (5.7) |
| Parainfluenza Virus 2 | 10 (0.5) | 0 (0.0) | 0 (0.0) | 0 (0.0) | 10 (1.1) | 0 (0.0) |
| Parainfluenza Virus 3 | 99 (5.1) | 13 (7.9) | 3 (1.2) | 28 (8.0) | 41 (4.6) | 14 (4.7) |
| Parainfluenza Virus 4 | 7 (0.4) | 0 (0.0) | 0 (0.0) | 1 (0.3) | 4 (0.4) | 2 (0.7) |
| RSV A | 28 (1.4) | 4 (2.4) | 6 (2.4) | 7 (2.0) | 4 (0.4) | 7 (2.4) |
| RSV B | 83 (4.3) | 6 (3.6) | 15 (6.0) | 24 (6.9) | 15 (1.7) | 23 (7.8) |
| <i>Chlamydia pneumoniae</i> | 3 (0.2) | 0 (0.0) | 0 (0.0) | 1 (0.3) | 2 (0.2) | 0 (0.0) |
| <i>Mycoplasma pneumoniae</i> | 5 (0.3) | 1 (0.6) | 0 (0.0) | 3 (0.9) | 0 (0.0) | 1 (0.3) |

Table 10: Expected Value (As Determined by cobas® eplex RP panel) Summary By Sample Collection Site in the Prospective Clinical Evaluation (Phase 2: September 2016 – October 2016)

| Organism | All Sites (N=511) n (%) | Site 5 (N=49) n (%) | Site 6 (N=101) n (%) | Site 7 (N=161) n (%) | Site 8 (N=200) n (%) |
|------------------------------|-------------------------------|---------------------------|----------------------------|----------------------------|----------------------------|
| Adenovirus | 10 (2.0) | 2 (4.1) | 3 (3.0) | 3 (1.9) | 2 (1.0) |
| Coronavirus | 8 (1.6) | 0 (0.0) | 2 (2.0) | 4 (2.5) | 2 (1.0) |
| Human Metapneumovirus | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) |
| Human Rhinovirus/Enterovirus | 188 (36.8) | 24 (49.0) | 49 (48.5) | 62 (38.5) | 53 (26.5) |
| Influenza A | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) |
| Influenza A H1 | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) |
| Influenza A H1-2009 | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) |
| Influenza A H3 | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) |
| Influenza B | 2 (0.4) | 1 (2.0) | 0 (0.0) | 0 (0.0) | 1 (0.5) |
| Parainfluenza Virus 1 | 1 (0.2) | 0 (0.0) | 0 (0.0) | 1 (0.6) | 0 (0.0) |
| Parainfluenza Virus 2 | 13 (2.5) | 2 (4.1) | 4 (4.0) | 3 (1.9) | 4 (2.0) |
| Parainfluenza Virus 3 | 5 (1.0) | 2 (4.1) | 2 (2.0) | 0 (0.0) | 1 (0.5) |
| Parainfluenza Virus 4 | 8 (1.6) | 1 (2.0) | 1 (1.0) | 4 (2.5) | 2 (1.0) |
| RSV A | 8 (1.6) | 0 (0.0) | 8 (7.9) | 0 (0.0) | 0 (0.0) |
| RSV B | 9 (1.8) | 1 (2.0) | 4 (4.0) | 0 (0.0) | 4 (2.0) |
| <i>Chlamydia pneumoniae</i> | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) | 0 (0.0) |
| <i>Mycoplasma pneumoniae</i> | 4 (0.8) | 0 (0.0) | 3 (3.0) | 0 (0.0) | 1 (0.5) |

PERFORMANCE CHARACTERISTICS

CLINICAL PERFORMANCE

Comparator Method

The performance of the cobas® eplex RP panel was compared to an FDA-cleared multiplexed molecular respiratory pathogen panel and analytically validated PCR tests with bi-directional sequencing for confirmation of RSV subtypes. Details of the comparator method are described in Table 11.

Table 11: Comparator Methods Used to Assess cobas® eplex RP panel Clinical Performance

| Target | Comparator Method |
|-------------------------------|---|
| Adenovirus | FDA-cleared multiplexed molecular respiratory pathogen panel |
| Coronavirus | |
| Human Metapneumovirus | |
| Human Rhinovirus/Enterovirus | |
| Influenza A | |
| Influenza A H1 | |
| Influenza A H1-2009 | |
| Influenza A H3 | |
| Influenza B | |
| Parainfluenza Virus 1 | |
| Parainfluenza Virus 2 | |
| Parainfluenza Virus 3 | |
| Parainfluenza Virus 4 | |
| Respiratory Syncytial Virus A | FDA-cleared multiplexed molecular respiratory pathogen panel followed by a PCR test with bi-directional sequencing confirmation |
| Respiratory Syncytial Virus B | |
| <i>Chlamydia pneumoniae</i> | FDA-cleared multiplexed molecular respiratory pathogen panel |
| <i>Mycoplasma pneumoniae</i> | |

Prospective Clinical Samples

Clinical performance was evaluated in clinical nasopharyngeal swab samples in VTM prospectively-collected at 8 clinical sites in 2 phases. From March 2013 through August 2014, 2218 samples were prospectively-collected and frozen; from September 2016 through October 2016, 514 samples were prospectively-collected and tested fresh (never frozen). A total of 2732 samples were collected across the 2 phases. Prior to the start of investigational testing, 263 samples were withdrawn (251 had sample handling deviations, 9 were tested outside of protocol timelines, 2 had insufficient volume, and 1 had incomplete documentation). Of the 2469 prospectively-collected samples eligible for testing, 2462 were evaluable. Samples with final, valid results and a valid comparator result were considered evaluable. Seven prospectively-collected samples were not evaluable because they did not have final, valid cobas® eplex RP panel results and were excluded from performance evaluations. Demographic information for prospectively-collected samples is described in Table 12. Subjects enrolled in this study were from a diverse demographic distribution and represent the intended patient population.

Table 12: Subject Demographic Data for Prospectively-Collected Samples by Collection Site (N=2462)

| | All Sites N=2462 n (%) | Site 1 N=165 n (%) | Site 2 N=248 n (%) | Site 3 N=350 n (%) | Site 4 N=892 n (%) | Site 5 N=345 n (%) | Site 6 N=101 n (%) | Site 7 N=161 n (%) | Site 8 N=200 n (%) |
|--------------------|------------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|
| Sex | | | | | | | | | |
| Male | 1247 (50.6) | 96 (58.2) | 118 (47.6) | 186 (53.1) | 450 (50.4) | 188 (54.5) | 43 (42.6) | 84 (52.2) | 82 (41.0) |
| Female | 1215 (49.4) | 69 (41.8) | 130 (52.4) | 164 (46.9) | 442 (49.6) | 157 (45.5) | 58 (57.4) | 77 (47.8) | 118 (59.0) |
| Age (years) | | | | | | | | | |
| 0–1 | 388 (15.8) | 17 (10.3) | 21 (8.5) | 74 (21.1) | 164 (18.4) | 45 (13.0) | 28 (27.7) | 3 (1.9) | 36 (18.0) |
| > 1–5 | 325 (13.2) | 12 (7.3) | 22 (8.9) | 62 (17.7) | 64 (7.2) | 100 (29.0) | 39 (38.6) | 16 (9.9) | 10 (5.0) |
| > 5–21 | 321 (13.0) | 15 (9.1) | 6 (2.4) | 38 (10.9) | 82 (9.2) | 116 (33.6) | 34 (33.7) | 18 (11.2) | 12 (6.0) |
| > 21–65 | 926 (37.6) | 87 (52.7) | 131 (52.8) | 98 (28.0) | 385 (43.2) | 55 (15.9) | 0 (0.0) | 92 (57.1) | 78 (39.0) |
| > 65 | 502 (20.4) | 34 (20.6) | 68 (27.4) | 78 (22.3) | 197 (22.1) | 29 (8.4) | 0 (0.0) | 32 (19.9) | 64 (32.0) |

Prospective Clinical Performance

Positive percent agreement (PPA) was calculated by dividing the number of true positive (TP) results by the sum of TP and false negative (FN) results, while negative percent agreement (NPA) was calculated by dividing the number of true negative (TN) results by the sum of TN and false positive (FP) results. A TP result was one where the detected cobas® eplex RP panel result matched the detected comparator method result, while a TN result was one where a negative cobas® eplex RP panel result matched a negative comparator method result. The two-sided 95% confidence interval was also calculated.

A total of 2462 prospectively-collected samples (511 tested fresh and 1951 tested after previously frozen) were evaluated for 17 cobas® eplex RP panel organisms. PPA and NPA results are summarized by target in Table 13 and Table 14 below.

Table 13: Positive Percent Agreement (PPA) and Negative Percent Agreement (NPA) in the cobas® eplex RP panel Clinical Study (Fresh)

| Organism | Prevalence | Positive % Agreement | | Negative % Agreement | |
|------------------------------|------------|----------------------|------------------|----------------------|------------------|
| | | TP/TP+FN | PPA (95% CI) | TN/TN+FP | NPA (95% CI) |
| Adenovirus | 1.6% | 6/8 ^a | 75.0 (40.9-92.9) | 499/503 ^a | 99.2 (98.0-99.7) |
| Coronavirus | 1.4% | 7/7 | 100 (64.6-100) | 503/504 | 99.8 (98.9-100) |
| Human Metapneumovirus | 0.0% | 0/0 | --- | 511/511 | 100 (99.3-100) |
| Human Rhinovirus/Enterovirus | 35.8% | 176/183 ^b | 96.2 (92.3-98.1) | 316/328 ^b | 96.3 (93.7-97.9) |
| Influenza A | 0.0% | 0/0 | --- | 511/511 | 100 (99.3-100) |
| Influenza A H1 | 0.0% | 0/0 | --- | 511/511 | 100 (99.3-100) |
| Influenza A H1-2009 | 0.0% | 0/0 | --- | 511/511 | 100 (99.3-100) |
| Influenza A H3 | 0.0% | 0/0 | --- | 511/511 | 100 (99.3-100) |
| Influenza B | 0.2% | 1/1 | 100 (20.7-100) | 509/510 | 99.8 (98.9-100) |
| Parainfluenza Virus 1 | 0.2% | 1/1 | 100 (20.7-100) | 510/510 | 100 (99.3-100) |
| Parainfluenza Virus 2 | 2.5% | 12/13 | 92.3 (66.7-98.6) | 497/498 | 99.8 (98.9-100) |
| Parainfluenza Virus 3 | 1.0% | 5/5 | 100 (56.6-100) | 506/506 | 100 (99.2-100) |
| Parainfluenza Virus 4 | 0.6% | 3/3 | 100 (43.9-100) | 503/508 ^c | 99.0 (97.7-99.6) |
| RSV A | 1.8% | 8/9 | 88.9 (56.5-98.0) | 501/501 | 100 (99.2-100) |
| RSV B | 2.0% | 9/10 | 90.0 (59.6-98.2) | 500/500 | 100 (99.2-100) |
| <i>Chlamydia pneumoniae</i> | 0.0% | 0/0 | --- | 511/511 | 100 (99.3-100) |
| <i>Mycoplasma pneumoniae</i> | 0.6% | 3/3 | 100 (43.9-100) | 507/508 ^d | 99.8 (98.9-100) |

^a Adenovirus was not detected in 2 of 2 FN samples and detected in 4 of 4 FP samples using PCR/sequencing.

^b Human rhinovirus/enterovirus was not detected in 1 of 7 FN samples and detected in 9 of 12 FP samples using PCR/sequencing.

^c Parainfluenza virus 4 was detected in 3 of 5 FP samples using PCR/sequencing.

^d *M. pneumoniae* was detected in the 1 FP sample using PCR/sequencing.

Table 14: Positive Percent Agreement (PPA) and Negative Percent Agreement (NPA) in the cobas® eplex RP panel Clinical Study (After Previously Frozen)

| Organism | Prevalence | Positive % Agreement | | Negative % Agreement | |
|------------------------------|------------|----------------------|------------------|------------------------|------------------|
| | | TP/TP+FN | PPA (95% CI) | TN/TN+FP | NPA (95% CI) |
| Adenovirus | 2.7% | 48/53 ^a | 90.6 (79.7-95.9) | 1874/1898 ^a | 98.7 (98.1-99.1) |
| Coronavirus | 5.6% | 89/110 ^b | 80.9 (72.6-87.2) | 1828/1841 ^b | 99.3 (98.8-99.6) |
| Human Metapneumovirus | 5.8% | 107/113 ^c | 94.7 (88.9-97.5) | 1832/1838 ^c | 99.7 (99.3-99.9) |
| Human Rhinovirus/Enterovirus | 17.2% | 317/336 ^d | 94.3 (91.3-96.4) | 1544/1615 ^d | 95.6 (94.5-96.5) |
| Influenza A ^e | 5.7% | 106/111 ^f | 95.5 (89.9-98.1) | 1836/1840 ^f | 99.8 (99.4-99.9) |
| Influenza A H1 | 0.0% | 0/0 | --- | 1951/1951 | 100 (99.8-100) |
| Influenza A H1-2009 | 3.6% | 70/71 | 98.6 (92.4-99.8) | 1874/1880 ^g | 99.7 (99.3-99.9) |
| Influenza A H3 | 1.9% | 34/37 ^h | 91.9 (78.7-97.2) | 1914/1914 | 100 (99.8-100) |
| Influenza B | 3.3% | 58/65 ⁱ | 89.2 (79.4-94.7) | 1882/1886 ⁱ | 99.8 (99.5-99.9) |
| Parainfluenza Virus 1 | 1.2% | 23/24 | 95.8 (79.8-99.3) | 1926/1927 | 99.9 (99.7-100) |
| Parainfluenza Virus 2 | 0.5% | 9/9 | 100 (70.1-100) | 1941/1942 | 99.9 (99.7-100) |
| Parainfluenza Virus 3 | 5.3% | 94/104 ^j | 90.4 (83.2-94.7) | 1842/1847 ^j | 99.7 (99.4-99.9) |
| Parainfluenza Virus 4 | 0.3% | 5/5 | 100 (56.6-100) | 1944/1946 | 99.9 (99.6-100) |
| RSV A | 1.6% | 27/31 | 87.1 (71.1-94.9) | 1917/1918 | 99.9 (99.7-100) |
| RSV B | 4.4% | 81/86 | 94.2 (87.1-97.5) | 1861/1863 ^k | 99.9 (99.6-100) |
| <i>Chlamydia pneumoniae</i> | 0.3% | 2/5 ^l | 40.0 (11.8-76.9) | 1945/1946 ^l | 99.9 (99.7-100) |
| <i>Mycoplasma pneumoniae</i> | 0.3% | 4/5 ^m | 80.0 (37.6-96.4) | 1945/1946 | 99.9 (99.7-100) |

- ^a Adenovirus was not detected in 1 of 5 FN samples and detected in 9 of 24 FP samples using PCR/sequencing.
- ^b Coronavirus was not detected in 2 of 21 FN samples and detected in 3 of 13 FP samples using PCR/sequencing.
- ^c Human Metapneumovirus was not detected in 1 of 6 FN samples and detected in 4 of 6 FP samples using PCR/sequencing.
- ^d Human rhinovirus/enterovirus was not detected in 6 of 19 FN samples and detected in 33 of 71 FP samples using PCR/sequencing.
- ^e Influenza A comparator results contain 71 samples with A H1-2009, 37 samples with A H3, and 3 samples with no subtype detected.
- ^f Influenza A was not detected in 1 of 3 FN samples (2 samples were not tested by PCR/sequencing) and detected in 1 of 4 FP samples using PCR/sequencing.
- ^g Influenza A H1-2009 was detected in 4 of 6 FP samples using PCR/sequencing.
- ^h Influenza A H3 was not detected in 1 of 3 FN samples using PCR/sequencing.
- ⁱ Influenza B was not detected in 3 of 7 FN samples and detected in 2 of 4 FP samples using PCR/sequencing.
- ^j Parainfluenza virus 3 was not detected in 3 of 10 FN samples and detected in 4 of 5 FP samples using PCR/sequencing.
- ^k RSV B was detected in 1 of 2 FP samples using PCR/sequencing.
- ^l *C. pneumoniae* was not detected in 1 of 3 FN samples and detected in the 1 FP sample using PCR/sequencing.
- ^m *M. pneumoniae* was not detected in the 1 FN sample using PCR/sequencing.

Retrospective Clinical Samples

To supplement the number of positives for targets that were not sufficiently represented in the prospective collection, additional nasopharyngeal swab in VTM samples were retrospectively collected from 6 sites. A total of 535 nasopharyngeal swab samples that had previously tested positive for one or more of the target organisms during standard-of-care (SOC) testing were collected and stored frozen. Prior to the start of investigational testing, 11 samples were withdrawn due to noncompliance with the study protocol, and 52 samples were withdrawn because the organisms present had sufficient representation in other samples. In addition, the composition and integrity of the retrospective samples were confirmed with the same comparator method employed in the prospective clinical study (i.e., an FDA-cleared multiplexed respiratory pathogen panel). As the result of this confirmation testing using the comparator method, 26 additional samples were withdrawn because the original SOC testing positive results for the intended organisms were not confirmed when tested with the comparator method. Of the remaining 446 retrospectively-collected samples eligible for testing, all 446 were evaluable. Demographic information for retrospectively-collected samples is described in Table 15. Subjects enrolled in this study were from a diverse demographic distribution and represent the intended patient population.

Table 15: Subject Demographic Data for Retrospectively-Collected Samples by Collection Site (N=446)

| | All Sites N=446 n (%) | Site 1 N=1 n (%) | Site 2 N=1 n (%) | Site 3 N=129 n (%) | Site 4 N=18 n (%) | Site 5 N=131 n (%) | Site 6 N=166 n (%) |
|--------------------|-----------------------------|------------------------|------------------------|--------------------------|-------------------------|--------------------------|--------------------------|
| Sex | | | | | | | |
| Male | 232 (52.0) | 0 (0.0) | 1 (100) | 76 (58.9) | 11 (61.1) | 68 (51.9) | 76 (45.8) |
| Female | 214 (48.0) | 1 (100) | 0 (0.0) | 53 (41.1) | 7 (38.9) | 63 (48.1) | 90 (54.2) |
| Age (years) | | | | | | | |
| 0 – 1 | 122 (27.4) | 0 (0.0) | 0 (0.0) | 24 (18.6) | 5 (27.8) | 56 (42.7) | 37 (22.3) |
| > 1 – 5 | 107 (24.0) | 0 (0.0) | 1 (100) | 51 (39.5) | 3 (16.7) | 16 (12.2) | 36 (21.7) |
| > 5 – 21 | 59 (13.2) | 0 (0.0) | 0 (0.0) | 9 (7.0) | 2 (11.1) | 19 (14.5) | 29 (17.5) |
| > 21 – 65 | 99 (22.2) | 1 (100) | 0 (0.0) | 11 (8.5) | 8 (44.4) | 31 (23.7) | 48 (28.9) |
| > 65 | 59 (13.2) | 0 (0.0) | 0 (0.0) | 34 (26.4) | 0 (0.0) | 9 (6.9) | 16 (9.6) |

Retrospective Clinical Performance

A total of 446 retrospectively-collected samples were evaluated for 17 cobas® eplex RP panel organisms. The following specimens with the original positive SOC results for the unintended organisms that were not confirmed by the comparator method were excluded from the performance calculation for the respective organism: 1 coronavirus positive specimen, 3 human rhinovirus/enterovirus positive specimens, 1 influenza A positive specimen, 1 influenza A H3 positive specimen, 1 parainfluenza virus positive specimen. In addition, 5 unintended RSV positive specimens by the comparator method were not confirmed by PCR/sequencing with regard to determining RSV subtypes and therefore were excluded from the performance calculations for RSV A and RSV B. PPA and NPA results are summarized by target in Table 16 below.

Table 16: Positive Percent Agreement (PPA) and Negative Percent Agreement (NPA) of the cobas® eplex RP panel With Comparator Methods (Retrospective Collection)

| Organism | Positive % Agreement | | Negative % Agreement | |
|------------------------------|----------------------|------------------|----------------------|------------------|
| | TP/TP+FN | PPA (95% CI) | TN/TN+FP | NPA (95% CI) |
| Adenovirus | 55/56 ^a | 98.2 (90.6-99.7) | 386/390 ^a | 99.0 (97.4-99.6) |
| Coronavirus | 121/138 ^b | 87.7 (81.2-92.2) | 307/307 | 100 (98.8-100) |
| Human Metapneumovirus | 5/7 | 71.4 (35.9-91.8) | 439/439 | 100 (99.1-100) |
| Human Rhinovirus/Enterovirus | 37/41 | 90.2 (77.5-96.1) | 384/402 | 95.5 (93.0-97.1) |
| Influenza A ^c | 75/82 ^d | 91.5 (83.4-95.8) | 363/363 | 100 (99.0-100) |
| Influenza A H1 | 0/0 | --- | 446/446 | 100 (99.1-100) |
| Influenza A H1-2009 | 27/31 ^e | 87.1 (71.1-94.9) | 415/415 | 100 (99.1-100) |
| Influenza A H3 | 45/51 ^f | 88.2 (76.6-94.5) | 394/394 | 100 (99.0-100) |
| Influenza B | 1/1 | 100 (20.7-100) | 445/445 | 100 (99.1-100) |
| Parainfluenza Virus 1 | 43/48 ^g | 89.6 (77.8-95.5) | 396/397 | 99.7 (98.6-100) |
| Parainfluenza Virus 2 | 46/51 | 90.2 (79.0-95.7) | 395/395 | 100 (99.0-100) |
| Parainfluenza Virus 3 | 2/2 | 100 (34.2-100) | 444/444 | 100 (99.1-100) |
| Parainfluenza Virus 4 | 18/20 | 90.0 (69.9-97.2) | 426/426 | 100 (99.1-100) |
| RSV A | 25/27 | 92.6 (76.6-97.9) | 414/414 | 100 (99.1-100) |
| RSV B | 21/22 | 95.5 (78.2-99.2) | 419/419 | 100 (99.1-100) |
| <i>Chlamydia pneumoniae</i> | 1/1 | 100 (20.7-100) | 445/445 | 100 (99.1-100) |
| <i>Mycoplasma pneumoniae</i> | 7/7 | 100 (64.6-100) | 439/439 | 100 (99.1-100) |

^a Adenovirus was not detected in the 1 FN sample and detected in 2 of 4 FP samples using PCR/sequencing.

^b Coronavirus was not detected in 2 of 16 FN samples using PCR/sequencing (1 sample was not tested by PCR/sequencing).

^c Influenza A comparator results contain 31 samples with A H1-2009 and 51 samples with A H3 detected.

^d Influenza A was not detected in 3 of 7 FN samples using PCR/sequencing.

^e Influenza A H1-2009 was not detected in 2 of 4 FN samples using PCR/sequencing.

^f Influenza A H3 was not detected in 1 of 6 FN samples using PCR/sequencing.

^g Parainfluenza virus 1 was not detected in 2 of 5 FN samples using PCR/sequencing.

Contrived Sample Performance

There were 327 contrived samples created and tested to supplement the low prevalence targets on the RP Panel; 104 contained one or more low prevalence organisms and 223 were negative for the contrived organisms. All 327 contrived samples were tested with the **cobas® eplex** RP panel and 326 were evaluable. PPA and NPA results are summarized for these low prevalence organisms in Table 17 below.

Table 17: Positive Percent Agreement (PPA) and Negative Percent Agreement (NPA) of the **cobas® eplex** RP panel With Comparator Method (Contrived Samples)

| Organism | Positive % Agreement | | Negative % Agreement | |
|-----------------------------|----------------------|----------------|----------------------|----------------|
| | TP/TP+FN | PPA (95% CI) | TN/TN+FP | NPA (95% CI) |
| <i>Chlamydia pneumoniae</i> | 52/52 | 100 (93.1-100) | 274/274 | 100 (98.6-100) |
| Influenza A H1 | 51/51 | 100 (93.0-100) | 275/275 | 100 (98.6-100) |

Clinical and Contrived Sample Performance by Target

Table 18 to Table 34 below include the clinical performance by pathogen of prospective samples tested fresh (shown in Table 13), prospective samples tested after previously freezing (shown in Table 14), retrospective samples (shown in Table 16, and contrived samples (shown in Table 17).

Table 18: Positive Percent Agreement (PPA) and Negative Percent Agreement (NPA) for Adenovirus in the **cobas® eplex** RP panel Clinical Study

| Adenovirus | Sample Type | Positive % Agreement | | Negative % Agreement | |
|-----------------------------------|-------------|----------------------|------------------|------------------------|------------------|
| | | TP/TP+FN | PPA (95% CI) | TN/TN+FP | NPA (95% CI) |
| Prospectively-Collected Samples | Fresh | 6/8 ^a | 75.0 (40.9-92.9) | 499/503 ^a | 99.2 (98.0-99.7) |
| | Frozen | 48/53 ^b | 90.6 (79.7-95.9) | 1874/1898 ^b | 98.7 (98.1-99.1) |
| | Total | 54/61 | 88.5 (78.2-94.3) | 2373/2401 | 98.8 (98.3-99.2) |
| Retrospectively-Collected Samples | | 55/56 ^c | 98.2 (90.6-99.7) | 386/390 ^c | 99.0 (97.4-99.6) |

^a Adenovirus was not detected in 2 of 2 FN samples and detected in 4 of 4 FP samples using PCR/sequencing.

^b Adenovirus was not detected in 1 of 5 FN samples and detected in 9 of 24 FP samples using PCR/sequencing.

^c Adenovirus was not detected in the 1 FN sample and detected in 2 of 4 FP samples using PCR/sequencing.

Table 19: Positive Percent Agreement (PPA) and Negative Percent Agreement (NPA) for Coronavirus in the **cobas® eplex** RP panel Clinical Study

| Coronavirus | Sample Type | Positive % Agreement | | Negative % Agreement | |
|--|-------------|----------------------|------------------|------------------------|------------------|
| | | TP/TP+FN | PPA (95% CI) | TN/TN+FP | NPA (95% CI) |
| Prospectively-Collected Samples ^a | Fresh | 7/7 | 100 (64.6-100) | 503/504 | 99.8 (98.9-100) |
| | Frozen | 89/110 ^b | 80.9 (72.6-87.2) | 1828/1841 ^b | 99.3 (98.8-99.6) |
| | Total | 96/117 | 82.1 (74.1-88.0) | 2331/2345 | 99.4 (99.0-99.6) |
| Retrospectively-Collected Samples ^c | | 121/138 ^d | 87.7 (81.2-92.2) | 307/307 | 100 (98.8-100) |

^a 20 FN prospectively-collected frozen samples were repeat tested with the comparator method and 12 had coronavirus detected. Of these 12 samples, 11 were repeat tested with the **cobas® eplex** RP panel and 3 had coronavirus detected.

^b Coronavirus was not detected in 2 of 21 FN samples and detected in 3 of 13 FP samples using PCR/sequencing.

^c 10 FN retrospectively-collected samples were repeat tested with the comparator method and all 10 had coronavirus detected. Of these 10 samples, 9 were repeat tested with the **cobas® eplex** RP panel and 5 had coronavirus detected.

^d Coronavirus was not detected in 2 of 16 FN samples using PCR/sequencing (1 sample was not tested by PCR/sequencing).

Table 20: Positive Percent Agreement (PPA) and Negative Percent Agreement (NPA) for Human Metapneumovirus in the cobas® eplex RP panel Clinical Study

| Human Metapneumovirus | Sample Type | Positive % Agreement | | Negative % Agreement | |
|-----------------------------------|-------------|----------------------|------------------|------------------------|------------------|
| | | TP/TP+FN | PPA (95% CI) | TN/TN+FP | NPA (95% CI) |
| Prospectively-Collected Samples | Fresh | 0/0 | --- | 511/511 | 100 (99.3-100) |
| | Frozen | 107/113 ^a | 94.7 (88.9-97.5) | 1832/1838 ^a | 99.7 (99.3-99.9) |
| | Total | 107/113 | 94.7 (88.9-97.5) | 2343/2349 | 99.7 (99.4-99.9) |
| Retrospectively-Collected Samples | | 5/7 | 71.4 (35.9-91.8) | 439/439 | 100 (99.1-100) |

^a Human Metapneumovirus was not detected in 1 of 6 FN samples and detected in 4 of 6 FP samples using PCR/sequencing.

Table 21: Positive Percent Agreement (PPA) and Negative Percent Agreement (NPA) for Human Rhinovirus/Enterovirus in the cobas® eplex RP panel Clinical Study

| Human Rhinovirus/Enterovirus | Sample Type | Positive % Agreement | | Negative % Agreement | |
|-----------------------------------|-------------|----------------------|------------------|------------------------|------------------|
| | | TP/TP+FN | PPA (95% CI) | TN/TN+FP | NPA (95% CI) |
| Prospectively-Collected Samples | Fresh | 176/183 ^a | 96.2 (92.3-98.1) | 316/328 ^a | 96.3 (93.7-97.9) |
| | Frozen | 317/336 ^b | 94.3 (91.3-96.4) | 1544/1615 ^b | 95.6 (94.5-96.5) |
| | Total | 493/519 | 95.0 (92.8-96.6) | 1860/1943 | 95.7 (94.7-96.5) |
| Retrospectively-Collected Samples | | 37/41 | 90.2 (77.5-96.1) | 384/402 | 95.5 (93.0-97.1) |

^a Human rhinovirus/enterovirus was not detected in 1 of 7 FN samples and detected in 9 of 12 FP samples using PCR/sequencing.

^b Human rhinovirus/enterovirus was not detected in 6 of 19 FN samples and detected in 33 of 71 FP samples using PCR/sequencing.

Table 22: Positive Percent Agreement (PPA) and Negative Percent Agreement (NPA) for Influenza A in the cobas® eplex RP panel Clinical Study

| Influenza A | Sample Type | Positive % Agreement | | Negative % Agreement | |
|--|-------------|----------------------|------------------|------------------------|------------------|
| | | TP/TP+FN | PPA (95% CI) | TN/TN+FP | NPA (95% CI) |
| Prospectively-Collected Samples ^a | Fresh | 0/0 | --- | 511/511 | 100 (99.3-100) |
| | Frozen | 106/111 ^b | 95.5 (89.9-98.1) | 1836/1840 ^b | 99.8 (99.4-99.9) |
| | Total | 106/111 | 95.5 (89.9-98.1) | 2347/2351 | 99.8 (99.6-99.9) |
| Retrospectively-Collected Samples ^c | | 75/82 ^d | 91.5 (83.4-95.8) | 363/363 | 100 (99.0-100) |

^a Influenza A comparator results contain 71 samples with A H1-2009, 37 samples with A H3, and 3 samples with no subtype detected.

^b Influenza A was not detected in 1 of 3 FN samples (2 samples were not tested by PCR/sequencing) and detected in 1 of 4 FP samples using PCR/sequencing.

^c Influenza A comparator results contain 31 samples with A H1-2009 and 51 samples with A H3 detected.

^d Influenza A was not detected in 3 of 7 FN samples using PCR/sequencing.

Table 23: Positive Percent Agreement (PPA) and Negative Percent Agreement (NPA) for Influenza A H1 in the cobas® eplex RP panel Clinical Study

| Influenza A H1 | Sample Type | Positive % Agreement | | Negative % Agreement | |
|-----------------------------------|-------------|----------------------|----------------|----------------------|----------------|
| | | TP/TP+FN | PPA (95% CI) | TN/TN+FP | NPA (95% CI) |
| Prospectively-Collected Samples | Fresh | 0/0 | --- | 511/511 | 100 (99.3-100) |
| | Frozen | 0/0 | --- | 1951/1951 | 100 (99.8-100) |
| | Total | 0/0 | --- | 2462/2462 | 100 (99.8-100) |
| Retrospectively-Collected Samples | | 0/0 | --- | 446/446 | 100 (99.1-100) |
| Contrived Samples | | 51/51 | 100 (93.0-100) | 275/275 | 100 (98.6-100) |

Table 24: Positive Percent Agreement (PPA) and Negative Percent Agreement (NPA) for Influenza A H1-2009 in the cobas® eplex RP panel Clinical Study

| Influenza A H1-2009 | Sample Type | Positive % Agreement | | Negative % Agreement | |
|-----------------------------------|-------------|----------------------|------------------|------------------------|------------------|
| | | TP/TP+FN | PPA (95% CI) | TN/TN+FP | NPA (95% CI) |
| Prospectively-Collected Samples | Fresh | 0/0 | --- | 511/511 | 100 (99.3-100) |
| | Frozen | 70/71 | 98.6 (92.4-99.8) | 1874/1880 ^a | 99.7 (99.3-99.9) |
| | Total | 70/71 | 98.6 (92.4-99.8) | 2385/2391 | 99.7 (99.5-99.9) |
| Retrospectively-Collected Samples | | 27/31 ^b | 87.1 (71.1-94.9) | 415/415 | 100 (99.1-100) |

^a Influenza A H1-2009 was detected in 4 of 6 FP samples using PCR/sequencing.

^b Influenza A H1-2009 was not detected in 2 of 4 FN samples using PCR/sequencing.

Table 25: Positive Percent Agreement (PPA) and Negative Percent Agreement (NPA) for Influenza A H3 in the cobas® eplex RP panel Clinical Study

| Influenza A H3 | Sample Type | Positive % Agreement | | Negative % Agreement | |
|-----------------------------------|-------------|----------------------|------------------|----------------------|----------------|
| | | TP/TP+FN | PPA (95% CI) | TN/TN+FP | NPA (95% CI) |
| Prospectively-Collected Samples | Fresh | 0/0 | --- | 511/511 | 100 (99.3-100) |
| | Frozen | 34/37 ^a | 91.9 (78.7-97.2) | 1914/1914 | 100 (99.8-100) |
| | Total | 34/37 | 91.9 (78.7-97.2) | 2425/2425 | 100 (99.8-100) |
| Retrospectively-Collected Samples | | 45/51 ^b | 88.2 (76.6-94.5) | 394/394 | 100 (99.0-100) |

^a Influenza A H3 was not detected in 1 of 3 FN samples using PCR/sequencing.^b Influenza A H3 was not detected in 1 of 6 FN samples using PCR/sequencing.**Table 26:** Positive Percent Agreement (PPA) and Negative Percent Agreement (NPA) for Influenza B in the cobas® eplex RP panel Clinical Study

| Influenza B | Sample Type | Positive % Agreement | | Negative % Agreement | |
|-----------------------------------|-------------|----------------------|------------------|------------------------|------------------|
| | | TP/TP+FN | PPA (95% CI) | TN/TN+FP | NPA (95% CI) |
| Prospectively-Collected Samples | Fresh | 1/1 | 100 (20.7-100) | 509/510 | 99.8 (98.9-100) |
| | Frozen | 58/65 ^a | 89.2 (79.4-94.7) | 1882/1886 ^a | 99.8 (99.5-99.9) |
| | Total | 59/66 | 89.4 (79.7-94.8) | 2391/2396 | 99.8 (99.5-99.9) |
| Retrospectively-Collected Samples | | 1/1 | 100 (20.7-100) | 445/445 | 100 (99.1-100) |

^a Influenza B was not detected in 3 of 7 FN samples and detected in 2 of 4 FP samples using PCR/sequencing.**Table 27:** Positive Percent Agreement (PPA) and Negative Percent Agreement (NPA) for Parainfluenza Virus 1 in the cobas® eplex RP panel Clinical Study

| Parainfluenza Virus 1 | Sample Type | Positive % Agreement | | Negative % Agreement | |
|-----------------------------------|-------------|----------------------|------------------|----------------------|-----------------|
| | | TP/TP+FN | PPA (95% CI) | TN/TN+FP | NPA (95% CI) |
| Prospectively-Collected Samples | Fresh | 1/1 | 100 (20.7-100) | 510/510 | 100 (99.3-100) |
| | Frozen | 23/24 | 95.8 (79.8-99.3) | 1926/1927 | 99.9 (99.7-100) |
| | Total | 24/25 | 96.0 (80.5-99.3) | 2436/2437 | 100 (99.8-100) |
| Retrospectively-Collected Samples | | 43/48 ^a | 89.6 (77.8-95.5) | 396/397 | 99.7 (98.6-100) |

^a Parainfluenza virus 1 was not detected in 2 of 5 FN samples using PCR/sequencing.**Table 28:** Positive Percent Agreement (PPA) and Negative Percent Agreement (NPA) for Parainfluenza Virus 2 in the cobas® eplex RP panel Clinical Study

| Parainfluenza Virus 2 | Sample Type | Positive % Agreement | | Negative % Agreement | |
|-----------------------------------|-------------|----------------------|------------------|----------------------|-----------------|
| | | TP/TP+FN | PPA (95% CI) | TN/TN+FP | NPA (95% CI) |
| Prospectively-Collected Samples | Fresh | 12/13 | 92.3 (66.7-98.6) | 497/498 | 99.8 (98.9-100) |
| | Frozen | 9/9 | 100 (70.1-100) | 1941/1942 | 99.9 (99.7-100) |
| | Total | 21/22 | 95.5 (78.2-99.2) | 2438/2440 | 99.9 (99.7-100) |
| Retrospectively-Collected Samples | | 46/51 | 90.2 (79.0-95.7) | 395/395 | 100 (99.0-100) |

Table 29: Positive Percent Agreement (PPA) and Negative Percent Agreement (NPA) for Parainfluenza Virus 3 in the cobas® eplex RP panel Clinical Study

| Parainfluenza Virus 3 | Sample Type | Positive % Agreement | | Negative % Agreement | |
|-----------------------------------|-------------|----------------------|------------------|------------------------|------------------|
| | | TP/TP+FN | PPA (95% CI) | TN/TN+FP | NPA (95% CI) |
| Prospectively-Collected Samples | Fresh | 5/5 | 100 (56.6-100) | 506/506 | 100 (99.2-100) |
| | Frozen | 94/104 ^a | 90.4 (83.2-94.7) | 1842/1847 ^a | 99.7 (99.4-99.9) |
| | Total | 99/109 | 90.8 (83.9-94.9) | 2348/2353 | 99.8 (99.5-99.9) |
| Retrospectively-Collected Samples | | 2/2 | 100 (34.2-100) | 444/444 | 100 (99.1-100) |

^a Parainfluenza virus 3 was not detected in 3 of 10 FN samples and detected in 4 of 5 FP samples using PCR/sequencing.

Table 30: Positive Percent Agreement (PPA) and Negative Percent Agreement (NPA) for Parainfluenza Virus 4 in the cobas® eplex RP panel Clinical Study

| Parainfluenza Virus 4 | Sample Type | Positive % Agreement | | Negative % Agreement | |
|-----------------------------------|-------------|----------------------|------------------|----------------------|------------------|
| | | TP/TP+FN | PPA (95% CI) | TN/TN+FP | NPA (95% CI) |
| Prospectively-Collected Samples | Fresh | 3/3 | 100 (43.9-100) | 503/508 ^a | 99.0 (97.7-99.6) |
| | Frozen | 5/5 | 100 (56.6-100) | 1944/1946 | 99.9 (99.6-100) |
| | Total | 8/8 | 100 (67.6-100) | 2447/2454 | 99.7 (99.4-99.9) |
| Retrospectively-Collected Samples | | 18/20 | 90.0 (69.9-97.2) | 426/426 | 100 (99.1-100) |

^a Parainfluenza virus 4 was detected in 3 of 5 FP samples using PCR/sequencing.

Table 31: Positive Percent Agreement (PPA) and Negative Percent Agreement (NPA) for Respiratory Syncytial Virus A in the cobas® eplex RP panel Clinical Study

| Respiratory Syncytial Virus A | Sample Type | Positive % Agreement | | Negative % Agreement | |
|-----------------------------------|-------------|----------------------|------------------|----------------------|-----------------|
| | | TP/TP+FN | PPA (95% CI) | TN/TN+FP | NPA (95% CI) |
| Prospectively-Collected Samples | Fresh | 8/9 | 88.9 (56.5-98.0) | 501/501 | 100 (99.2-100) |
| | Frozen | 27/31 | 87.1 (71.1-94.9) | 1917/1918 | 99.9 (99.7-100) |
| | Total | 35/40 | 87.5 (73.9-94.5) | 2418/2419 | 100 (99.8-100) |
| Retrospectively-Collected Samples | | 25/27 | 92.6 (76.6-97.9) | 414/414 | 100 (99.1-100) |

Table 32: Positive Percent Agreement (PPA) and Negative Percent Agreement (NPA) for Respiratory Syncytial Virus B in the cobas® eplex RP panel Clinical Study

| Respiratory Syncytial Virus B | Sample Type | Positive % Agreement | | Negative % Agreement | |
|-----------------------------------|-------------|----------------------|------------------|------------------------|-----------------|
| | | TP/TP+FN | PPA (95% CI) | TN/TN+FP | NPA (95% CI) |
| Prospectively-Collected Samples | Fresh | 9/10 | 90.0 (59.6-98.2) | 500/500 | 100 (99.2-100) |
| | Frozen | 81/86 | 94.2 (87.1-97.5) | 1861/1863 ^a | 99.9 (99.6-100) |
| | Total | 90/96 | 93.8 (87.0-97.1) | 2361/2363 | 99.9 (99.7-100) |
| Retrospectively-Collected Samples | | 21/22 | 95.5 (78.2-99.2) | 419/419 | 100 (99.1-100) |

^a RSV B was detected in 1 of 2 FP samples using PCR/sequencing.

Table 33: Positive Percent Agreement (PPA) and Negative Percent Agreement (NPA) for Chlamydia pneumoniae in the cobas® eplex RP panel Clinical Study

| <i>Chlamydia pneumoniae</i> | Sample Type | Positive % Agreement | | Negative % Agreement | |
|-----------------------------------|-------------|----------------------|------------------|------------------------|-----------------|
| | | TP/TP+FN | PPA (95% CI) | TN/TN+FP | NPA (95% CI) |
| Prospectively-Collected Samples | Fresh | 0/0 | --- | 511/511 | 100 (99.3-100) |
| | Frozen | 2/5 ^a | 40.0 (11.8-76.9) | 1945/1946 ^a | 99.9 (99.7-100) |
| | Total | 2/5 | 40.0 (11.8-76.9) | 2456/2457 | 100 (99.8-100) |
| Retrospectively-Collected Samples | | 1/1 | 100 (20.7-100) | 445/445 | 100 (99.1-100) |
| Contrived Samples | | 52/52 | 100 (93.1-100) | 274/274 | 100 (98.6-100) |

^a *C. pneumoniae* was not detected in 1 of 3 FN samples and detected in the 1 FP sample using PCR/sequencing.

Table 34: Positive Percent Agreement (PPA) and Negative Percent Agreement (NPA) for Mycoplasma pneumoniae in the cobas® eplex RP panel Clinical Study

| <i>Mycoplasma pneumoniae</i> | Sample Type | Positive % Agreement | | Negative % Agreement | |
|-----------------------------------|-------------|----------------------|------------------|----------------------|-----------------|
| | | TP/TP+FN | PPA (95% CI) | TN/TN+FP | NPA (95% CI) |
| Prospectively-Collected Samples | Fresh | 3/3 | 100 (43.9-100) | 507/508 ^a | 99.8 (98.9-100) |
| | Frozen | 4/5 ^b | 80.0 (37.6-96.4) | 1945/1946 | 99.9 (99.7-100) |
| | Total | 7/8 | 87.5 (52.9-97.8) | 2452/2454 | 99.9 (99.7-100) |
| Retrospectively-Collected Samples | | 7/7 | 100 (64.6-100) | 439/439 | 100 (99.1-100) |

^a *M. pneumoniae* was detected in the 1 FP sample using PCR/sequencing.

^b *M. pneumoniae* was not detected in the 1 FN sample using PCR/sequencing.

Co-detections in Prospective Clinical Samples

The cobas® eplex RP panel identified a total of 135 prospective samples with multiple organisms detected, or 5.5% of all prospectively-collected samples. Of these, 118 (4.8%) had two organisms, 14 (0.6%) had three organisms, and 3 (0.1%) had four organisms detected. Of the 135 co-detected samples, 58 included 1 or more organisms that had not been detected by the comparator method(s). Results are summarized in Table 35 and Table 36.

Table 35: Distinct Co-Detection Combinations Detected by the cobas® eplex RP panel in the Prospective Clinical Samples

| Distinct Co-Detection Combinations Detected by the cobas® eplex RP panel | | | | Total Number Of Co-detections (% of samples) | Number of Discrepant Co-detections | Discrepant Organism(s) ^a |
|--|-------------|------------|------------|--|------------------------------------|---|
| Organism 1 | Organism 2 | Organism 3 | Organism 4 | | | |
| ADV | CoV | | | 2 (0.08%) | 0 | |
| ADV | CoV | HRV/EV | | 2 (0.08%) | 1 | ADV (1) |
| ADV | Flu A (unk) | Flu B | HRV/EV | 1 (0.04%) | 1 | ADV (1), Flu A (unk) (1), Flu B (1), HRV/EV (1) |
| ADV | Flu AH3 | | | 1 (0.04%) | 0 | |
| ADV | Flu B | HRV/EV | RSV B | 1 (0.04%) | 1 | ADV (1), Flu B (1) |
| ADV | FluA09H1 | | | 1 (0.04%) | 1 | ADV (1), FluA09H1 (1) |
| ADV | FluA09H1 | HRV/EV | | 1 (0.04%) | 0 | |
| ADV | FluA09H1 | PIV 3 | | 1 (0.04%) | 1 | PIV 3 (1) |
| ADV | HMPV | | | 3 (0.12%) | 2 | ADV (2) |
| ADV | HMPV | HRV/EV | RSV A | 1 (0.04%) | 1 | RSV A (1) |
| ADV | HRV/EV | | | 18 (0.73%) | 7 | ADV (6), HRV/EV (1) |
| ADV | HRV/EV | Mpneum | | 1 (0.04%) | 0 | |
| ADV | HRV/EV | PIV 1 | | 1 (0.04%) | 1 | PIV 1 (1) |
| ADV | HRV/EV | PIV 4 | | 1 (0.04%) | 1 | ADV (1), PIV 4 (1) |
| ADV | HRV/EV | RSV B | | 1 (0.04%) | 0 | |
| ADV | PIV 2 | | | 2 (0.08%) | 1 | ADV (1) |
| ADV | PIV 3 | | | 2 (0.08%) | 1 | ADV (1) |
| ADV | PIV 4 | | | 1 (0.04%) | 1 | ADV (1) |
| ADV | RSV B | | | 2 (0.08%) | 2 | ADV (2) |
| CPneum | HRV/EV | | | 1 (0.04%) | 0 | |
| CoV | FluA09H1 | | | 1 (0.04%) | 0 | |
| CoV | HMPV | | | 4 (0.16%) | 0 | |
| CoV | HMPV | HRV/EV | | 2 (0.08%) | 0 | |
| CoV | HRV/EV | | | 12 (0.49%) | 4 | CoV (1), HRV/EV (4) |
| CoV | HRV/EV | RSV B | | 1 (0.04%) | 1 | CoV (1) |
| CoV | PIV 1 | | | 1 (0.04%) | 0 | |
| CoV | RSV A | | | 3 (0.12%) | 0 | |
| CoV | RSV B | | | 3 (0.12%) | 2 | CoV (2) |
| Flu A (unk) | HRV/EV | | | 1 (0.04%) | 1 | Flu A (unk) (1) |
| Flu AH3 | HRV/EV | | | 2 (0.08%) | 1 | HRV/EV (1) |
| Flu AH3 | RSV B | | | 1 (0.04%) | 0 | |
| Flu B | HRV/EV | | | 4 (0.16%) | 2 | HRV/EV (2) |
| Flu B | HRV/EV | RSV B | | 1 (0.04%) | 0 | |
| Flu B | PIV 3 | | | 1 (0.04%) | 0 | |

| Distinct Co-Detection Combinations Detected by the cobas® eplex RP panel | | | | Total Number Of Co-detections (% of samples) | Number of Discrepant Co-detections | Discrepant Organism(s) ^a |
|--|------------|------------|------------|--|------------------------------------|-------------------------------------|
| Organism 1 | Organism 2 | Organism 3 | Organism 4 | | | |
| FluA09H1 | HMPV | HRV/EV | | 1 (0.04%) | 1 | HRV/EV (1) |
| FluA09H1 | HRV/EV | | | 2 (0.08%) | 1 | HRV/EV (1) |
| HMPV | HRV/EV | | | 5 (0.20%) | 1 | HRV/EV (1) |
| HMPV | HRV/EV | RSV B | | 1 (0.04%) | 1 | HRV/EV (1) |
| HMPV | PIV 3 | | | 1 (0.04%) | 0 | |
| HRV/EV | PIV 1 | | | 3 (0.12%) | 0 | |
| HRV/EV | PIV 2 | | | 7 (0.28%) | 3 | HRV/EV (1), PIV 2 (2) |
| HRV/EV | PIV 3 | | | 11 (0.45%) | 5 | HRV/EV (5) |
| HRV/EV | PIV 4 | | | 4 (0.16%) | 4 | PIV 4 (4) |
| HRV/EV | RSV A | | | 5 (0.20%) | 0 | |
| HRV/EV | RSV B | | | 11 (0.45%) | 6 | HRV/EV (6) |
| PIV 1 | PIV 4 | | | 1 (0.04%) | 1 | PIV 4 (1) |
| PIV 3 | RSV B | | | 1 (0.04%) | 0 | |
| RSV A | RSV B | | | 1 (0.04%) | 1 | RSV B (1) |
| Total Number of Co-Detections | | | | 135 (5.5%) | 57 | 64/290 ^b |
| Total Number with 2 Organisms Detected | | | | 118 (4.8%) | 47 | 49/236 |
| Total Number with 3 Organisms Detected | | | | 14 (0.6%) | 7 | 8/42 |
| Total Number with 4 Organisms Detected | | | | 3 (0.1%) | 3 | 7/12 |

Note: ADV= adenovirus, CoV= coronavirus, HMPV= human metapneumovirus, HRV/EV= human rhinovirus/enterovirus, Flu= Influenza, (unk)= unknown subtype, PIV= parainfluenza, RSV= respiratory syncytial virus, Cpneum= *C. pneumoniae*, Mpneum= *M. pneumoniae*

^a A discrepant organism is defined as one that was detected by the cobas® eplex RP panel but not by the comparator method(s).

^b 64/64 discrepant organisms were investigated using PCR/sequencing; the discrepant organism was detected in 20/64 cases:

- In 8/18 samples, adenovirus was detected by PCR/sequencing.
- In 1/4 samples, coronavirus was detected by PCR/sequencing.
- In 7/25 samples, human rhinovirus/enterovirus was detected by PCR/sequencing.
- In 1/1 sample, influenza A H1-2009 was detected by PCR/sequencing.
- In 1/1 sample, parainfluenza virus 3 was detected by PCR/sequencing.
- In 2/6 samples, parainfluenza virus 4 was detected by PCR/sequencing.

Table 36: Additional Co-Detection Combinations Detected by the Comparator Method in the Prospective Clinical Samples

| Distinct Co-Detection Combinations Detected by the Comparator Method | | | Total Number Of Co-detections (% of samples) | Number of Discrepant Co-detections | Discrepant Organism(s) ^{a,b} |
|--|------------|------------|--|------------------------------------|---------------------------------------|
| Organism 1 | Organism 2 | Organism 3 | | | |
| ADV | CoV | | 1 (0.04%) | 1 | ADV (1), CoV (1) |
| ADV | HRV/EV | | 4 (0.16%) | 4 | ADV (4) |
| ADV | HRV/EV | PIV 3 | 1 (0.04%) | 1 | HRV/EV (1), PIV 3 (1) |
| ADV | HRV/EV | RSV A | 1 (0.04%) | 1 | ADV (1) |
| CPneum | HRV/EV | | 1 (0.04%) | 1 | CPneum (1) |
| CPneum | PIV 3 | | 1 (0.04%) | 1 | CPneum (1) |
| CoV | FluA09H1 | | 2 (0.08%) | 2 | CoV (2) |
| CoV | HMPV | | 1 (0.04%) | 1 | CoV (1) |
| CoV | HRV/EV | | 6 (0.24%) | 6 | CoV (4), HRV/EV (2) |
| CoV | PIV 3 | | 1 (0.04%) | 1 | CoV (1) |
| CoV | RSV B | | 3 (0.12%) | 3 | CoV (2), RSV B (1) |
| Flu AH3 | HRV/EV | PIV 3 | 1 (0.04%) | 1 | Flu AH3 (1), PIV 3 (1) |
| Flu AH3 | PIV 3 | | 1 (0.04%) | 1 | PIV 3 (1) |
| FluA09H1 | HMPV | HRV/EV | 1 (0.04%) | 1 | HMPV (1), HRV/EV (1) |
| HMPV | HRV/EV | | 1 (0.04%) | 1 | HRV/EV (1) |
| HRV/EV | PIV 1 | | 1 (0.04%) | 1 | HRV/EV (1) |
| HRV/EV | PIV 3 | | 2 (0.08%) | 2 | HRV/EV (2) |
| HRV/EV | PIV 3 | RSV B | 1 (0.04%) | 1 | PIV 3 (1) |
| HRV/EV | RSV A | | 2 (0.08%) | 2 | RSV A (2) |

^a A discrepant organism is defined as one that was detected by the comparator method(s) but not by the cobas® eplex RP panel.

^b 36/36 discrepant organisms were investigated using PCR/sequencing; the discrepant organism was not detected in 10/36 cases:

- In 2/6 samples, adenovirus was not detected by PCR/sequencing.
- In 1/2 samples, *Chlamydia pneumoniae* was not detected by PCR/sequencing.
- In 1/11 samples, coronavirus was not detected by PCR/sequencing.
- In 5/8 samples, human rhinovirus/enterovirus was not detected by PCR/sequencing.
- In 1/1 sample, influenza A H3 was not detected by PCR/sequencing.

Clinical Study cobas® eplex system Performance

A total of 3281 samples (including prospective, retrospective, and contrived samples) were initially tested in the clinical evaluations and 3127/3281 = 95.3% (95% CI: 94.5%-96.0%) generated valid results on the first attempt. After re-test, 8 samples had invalid results; final validity rate was 3273/3281 = 99.8% (95% CI: 99.5%-99.9%).

ANALYTICAL PERFORMANCE CHARACTERISTICS

Limit of Detection

The limit of detection (LoD), or analytical sensitivity was identified and verified for each viral and bacterial target on the cobas® eplex RP panel using quantified reference strains/isolates. Serial dilutions were prepared in a natural clinical matrix (pooled, negative nasopharyngeal swab in VTM samples) with one or more organisms per series, and at least 20 replicates per target were tested in the study. The limit of detection was defined as the lowest concentration at which each target is detected at least 95% of the time. The confirmed LoD for each cobas® eplex RP panel organism is shown in Table 37.

Table 37: LoD Results Summary

| Target | Strain | LoD Concentration |
|--------------------------------|---------------------------------------|---|
| Adenovirus | Type 1 (C) | 1×10^3 TCID ₅₀ /mL |
| | Type 4 (E) | 2×10^0 TCID ₅₀ /mL |
| | Type 7 (B) | 2×10^0 TCID ₅₀ /mL |
| Coronavirus 229E | 229E | 1×10^0 TCID ₅₀ /mL |
| Coronavirus HKU1 | HKU1 ^a | 5×10^4 copies/mL |
| Coronavirus NL63 | NL63 | 7.5×10^0 TCID ₅₀ /mL |
| Coronavirus OC43 | OC43 | 5×10^2 TCID ₅₀ /mL |
| Human Metapneumovirus | A1 IA3-2002 | 2×10^{-1} TCID ₅₀ /mL |
| | A2 IA14-2003 ^b | 2×10^3 TCID ₅₀ /mL |
| | B1 Peru2-2002 | 2×10^2 TCID ₅₀ /mL |
| | B2 Peru1-2002 | 2.25×10^2 TCID ₅₀ /mL |
| Human Rhinovirus/Enterovirus | Enterovirus Type 68 (2007) | 1×10^0 TCID ₅₀ /mL |
| | Rhinovirus 1A | 1.5×10^0 TCID ₅₀ /mL |
| | Rhinovirus B14 | 1×10^0 TCID ₅₀ /mL |
| | Rhinovirus C ^a | 1×10^5 copies/mL |
| Influenza A | H1N1 Brisbane/59/07 | 3×10^{-1} TCID ₅₀ /mL |
| Influenza A H1 | H1N1 Brisbane/59/07 | 3×10^{-1} TCID ₅₀ /mL |
| Influenza A H1-2009 | NY/01/2009 | 1×10^{-1} TCID ₅₀ /mL |
| Influenza A H3 | A/Perth/16/2009 | 1×10^1 TCID ₅₀ /mL |
| | A/Texas/50/2012 | 1×10^0 TCID ₅₀ /mL |
| | A/Victoria/361/2011 | 5×10^{-1} TCID ₅₀ /mL |
| | H3N2 Brisbane/10/07 | 5×10^1 TCID ₅₀ /mL |
| Influenza B (Victoria Lineage) | B/Brisbane/60/2008 | 1×10^0 TCID ₅₀ /mL |
| | B/Montana/5/2012 | 1×10^0 TCID ₅₀ /mL |
| | B/Nevada/03/2011 | 1×10^0 TCID ₅₀ /mL |
| Influenza B (Yamagata Lineage) | Florida/02/06 | 1×10^{-1} TCID ₅₀ /mL |
| | B/Massachusetts/02/2012 | 1×10^2 TCID ₅₀ /mL |
| | B/Texas/06/2011 | 1×10^{-1} TCID ₅₀ /mL |
| | B/Wisconsin/01/2010 | 1×10^0 TCID ₅₀ /mL |
| Parainfluenza Virus 1 | Clinical Isolate | 4×10^{-1} TCID ₅₀ /mL |
| Parainfluenza Virus 2 | Clinical Isolate | 5×10^1 TCID ₅₀ /mL |
| Parainfluenza Virus 3 | Clinical Isolate | 5×10^0 TCID ₅₀ /mL |
| Parainfluenza Virus 4 | 4a | 3×10^1 TCID ₅₀ /mL |
| Respiratory Syncytial Virus A | 2006 Isolate | 1.5×10^0 TCID ₅₀ /mL |
| Respiratory Syncytial Virus B | CH93(18)-18 | 2×10^{-1} TCID ₅₀ /mL |
| <i>Chlamydia pneumoniae</i> | AR-39 | 3×10^2 TCID ₅₀ /mL |
| <i>Mycoplasma pneumoniae</i> | FH strain of Eaton Agent [NCTC 10119] | 3×10^2 CCU/mL |

^a Clinical samples confirmed positive for coronavirus HKU1 and human rhinovirus C by bi-directional sequencing and quantified by real-time RT-PCR were used for determination of LoD.

^b Customer communication from manufacturer dated July 9, 2020 indicated that the human metapneumovirus strain sold as IA14-2003 was actually type B.

Analytical Reactivity (Inclusivity)

A panel of 101 strains/isolates representing the genetic, temporal, and geographic diversity of each target on the cobas® eplex RP panel was evaluated to demonstrate analytical reactivity. Each strain/isolate was tested in triplicate at 3x LoD in natural clinical matrix (pooled, negative nasopharyngeal swab in VTM samples); if the organism was not detected at this concentration, testing of higher concentrations was performed. Additional in silico analysis was also performed on a subset of cobas® eplex RP panel organisms.

All of the 101 strains/isolates tested for inclusivity were detected by the cobas® eplex RP panel. Results of analytical reactivity are shown in Table 38 to Table 48.

Table 38: Analytical Reactivity (Inclusivity) Results for Adenovirus

| Adenovirus Species | Serotype | Concentration | Multiple of LoD Detected |
|--------------------|---------------------|--|--------------------------|
| A | Type 31 | 3×10^3 TCID ₅₀ /mL | 3x |
| B | Type 3 | 6×10^0 TCID ₅₀ /mL | 3x |
| | Type 11 | 6×10^0 TCID ₅₀ /mL | 3x |
| | De Wit Type 14 | 6×10^0 TCID ₅₀ /mL | 3x |
| | Ch.79 Type 16 | 2×10^2 TCID ₅₀ /mL | 100x ^a |
| | Type 21 | 6×10^0 TCID ₅₀ /mL | 3x |
| | Compton Type 34 | 6×10^0 TCID ₅₀ /mL | 3x |
| | Holden Type 35 | 6×10^0 TCID ₅₀ /mL | 3x |
| | Wan Type 50 | 2×10^1 TCID ₅₀ /mL | 10x ^b |
| C | Type 2 | 3×10^3 TCID ₅₀ /mL | 3x |
| | Type 5 | 3×10^3 TCID ₅₀ /mL | 3x |
| | Type 6 | 3×10^3 TCID ₅₀ /mL | 3x |
| D | Type 26 | 3×10^3 TCID ₅₀ /mL | 3x |
| | Type 37 | 3×10^3 TCID ₅₀ /mL | 3x |
| F | Type 40 Dugan | 3×10^3 TCID ₅₀ /mL | 3x |
| | Type 41/ Strain Tak | 3×10^3 TCID ₅₀ /mL | 3x |

^a In silico analysis revealed good homology to primers and probes. Lower sensitivity is likely the result of incorrect estimation of genetic material present in the culture of this or the reference strain (TCID₅₀ value is based only on infectious virus particles).

^b In silico analysis revealed that lower sensitivity may be a result of mismatches in the assay primers and/or probes.

Table 39: Analytical Reactivity (Inclusivity) Results for Human Metapneumovirus

| Metapneumovirus Subtype | Strain | Concentration | Multiple of LoD Detected |
|-------------------------|------------------|---|--------------------------|
| Human Metapneumovirus | Peru6-2003 G, B2 | 6.75×10^2 TCID ₅₀ /mL | 3x |

Table 40: Analytical Reactivity (Inclusivity) Results for Human Rhinovirus/Enterovirus

| Rhinovirus/Enterovirus | Strain | Concentration | Multiple of LoD Detected |
|------------------------|----------|--|--------------------------|
| Human Rhinovirus | Type A2 | 4.5×10^0 TCID ₅₀ /mL | 3x |
| | Type A7 | 1.5×10^1 TCID ₅₀ /mL | 10x ^a |
| | Type A16 | 4.5×10^0 TCID ₅₀ /mL | 3x |
| | Type A18 | 1.5×10^2 TCID ₅₀ /mL | 100x ^a |
| | Type A34 | 4.5×10^0 TCID ₅₀ /mL | 3x |
| | Type A57 | 4.5×10^0 TCID ₅₀ /mL | 3x |
| | Type A77 | 4.5×10^0 TCID ₅₀ /mL | 3x |
| | 277G | 4.5×10^0 TCID ₅₀ /mL | 3x |
| | Type B3 | 1.5×10^1 TCID ₅₀ /mL | 10x ^a |
| | Type B17 | 1.5×10^1 TCID ₅₀ /mL | 10x ^a |
| | Type B42 | 4.5×10^0 TCID ₅₀ /mL | 3x |
| | Type B83 | 4.5×10^0 TCID ₅₀ /mL | 3x |
| | Type B84 | 4.5×10^0 TCID ₅₀ /mL | 3x |
| | FO2-2547 | 4.5×10^0 TCID ₅₀ /mL | 3x |
| Enterovirus | Type 71 | 3×10^0 TCID ₅₀ /mL | 3x |
| Coxsackievirus | A9 | 3×10^0 TCID ₅₀ /mL | 3x |
| | A10 | 3×10^0 TCID ₅₀ /mL | 3x |
| | A21 | 3×10^0 TCID ₅₀ /mL | 3x |
| | A24 | 3×10^0 TCID ₅₀ /mL | 3x |
| | B2 | 1×10^2 TCID ₅₀ /mL | 100x ^a |
| | B3 | 3×10^0 TCID ₅₀ /mL | 3x |
| | B4 | 3×10^0 TCID ₅₀ /mL | 3x |
| | B5 | 1×10^1 TCID ₅₀ /mL | 10x ^a |
| Echovirus | 9 | 3×10^0 TCID ₅₀ /mL | 3x |
| | E6 | 1×10^1 TCID ₅₀ /mL | 10x ^b |
| | 25 | 1×10^1 TCID ₅₀ /mL | 10x ^a |
| | 30 | 3×10^0 TCID ₅₀ /mL | 3x |
| Poliovirus | 1 | 1×10^2 TCID ₅₀ /mL | 100x ^a |

^a In silico analysis revealed that lower sensitivity may be a result of mismatches in the assay primers and/or probes.

^b In silico analysis revealed good homology to primers and probes. Lower sensitivity is likely the result of incorrect estimation of genetic material present in the culture of this or the reference strain (TCID₅₀ value is based only on infectious virus particles).

Table 41: Analytical Reactivity (Inclusivity) Results for Influenza A

Note: Due to different assays for influenza A matrix and influenza A subtypes on the cobas® eplex RP panel, if different LoDs are observed for inclusivity for a Flu A matrix vs. HA subtype, the differences are noted in the Multiple of LoD Detected column.

| Influenza A Subtype | Strain | Concentration | Multiple of LoD Detected |
|-----------------------|--------------------------|--|---|
| Influenza A H1 | A/FM/1/47 | 3 x 10 ⁰ TCID ₅₀ /mL | 10x (Influenza A matrix) ^a 10000x H1 subtype ^b |
| | A/New Caledonia/20/1999 | 9 x 10 ⁻¹ TCID ₅₀ /mL | 3x |
| | A/New Jersey/8/76 | 9 x 10 ⁻¹ TCID ₅₀ /mL | 3x H1 subtype not detected ^c |
| | A/NWS/33 | 3 x 10 ⁰ TCID ₅₀ /mL | 10x (Influenza A matrix) ^a H1 subtype not detected ^d |
| | A/PR/8/34 | 9 x 10 ⁻¹ TCID ₅₀ /mL | 3x (Influenza A matrix) H1 subtype not detected ^e |
| | A/Solomon Islands/3/2006 | 9 x 10 ⁻¹ TCID ₅₀ /mL | 3x |
| | A/Taiwan/42/06 | 9 x 10 ⁰ TCID ₅₀ /mL | 30x ^f |
| Influenza A H3 | A/Hong Kong/8/68 | 1.5 x 10 ² TCID ₅₀ /mL | 3x |
| | A/Port Chalmers/1/73 | | |
| | A/Nanchang/933/95 | | |
| | A/Victoria/3/75 | | |
| | A/Wisconsin/67/05 | | |
| Influenza A 2009 H1N1 | A/California/7/2009 | 1 x 10 ⁰ TCID ₅₀ /mL | 10x ^g |
| | A/Mexico/4108/09 | 3 x 10 ⁻¹ TCID ₅₀ /mL | 3x |
| | A/NY/02/2009 | 1 x 10 ⁰ TCID ₅₀ /mL | 10x ^h |
| | A/Swine NY/03/2009 | 3 x 10 ⁻¹ TCID ₅₀ /mL | 3x |
| | A/Swine/Iowa/15/30 | 3 x 10 ⁻¹ TCID ₅₀ /mL | 3x (Influenza A matrix) 100,000x (H1-2009 subtype) ⁱ |
| | A/Virginia/ATCC1/2009 | 1 x 10 ⁰ TCID ₅₀ /mL | 10x ^j |
| | A/Virginia/ATCC2/2009 | 1 x 10 ¹ TCID ₅₀ /mL | 100x ^j |
| | A/Virginia/ATCC3/2009 | 1 x 10 ² TCID ₅₀ /mL | 1,000x ^j |

^a In silico analysis revealed good homology to primers and probes. Lower sensitivity is likely the result of incorrect estimation of genetic material present in the culture of this or the reference strain (TCID₅₀ value is based only on infectious virus particles).

^b In silico analysis revealed that lower sensitivity may be a result of mismatches in the assay primers and/or probes.

^c H1-2009 subtype was detected in this seasonal influenza A H1 strain at 30x LoD.

^d In silico analysis revealed little homology between this non-contemporary strain sequence and the H1 signal probe/capture probe sequences.

^e In silico analysis revealed little homology between this non-contemporary influenza strain sequence and the H1 primer sequences.

^f For Influenza A matrix, in silico analysis revealed good homology to primers and probes. Lower sensitivity is likely the result of incorrect estimation of genetic material present in the culture of this or the reference strain (TCID₅₀ value is based only on infectious virus particles). For H1 subtype, in silico analysis revealed that lower sensitivity may be a result of mismatches in the assay primers and/or probes.

^g For Influenza A matrix, in silico analysis revealed that lower sensitivity may be a result of mismatches in the assay primers and/or probes. For H1 subtype, in silico analysis revealed good homology to primers and probes. Lower sensitivity is likely the result of incorrect estimation of genetic material present in the culture of this or the reference strain (TCID₅₀ value is based only on infectious virus particles).

^h For Influenza A matrix, in silico analysis revealed good homology to primers and probes. Lower sensitivity is likely the result of incorrect estimation of genetic material present in the culture of this or the reference strain (TCID₅₀ value is based only on infectious virus particles). For H1-2009 subtype, in silico analysis revealed that lower sensitivity may be a result of mismatches in the assay primers and/or probes.

ⁱ In silico analysis revealed little homology between the strain sequence and the H1 or H1-2009 primer, signal probe and capture probe sequences.

^j No sequence data was available to investigate lower sensitivity of the influenza A 2009 H1N1 A/Virginia/ATCC1/2009, A/Virginia/ATCC2/2009 and A/Virginia/ATCC3/2009 strains.

Table 42: Analytical Reactivity (Inclusivity) Results for Influenza A Strains Titered with Methods Different From the Reference Strain

| Influenza A Subtype | Strain | Concentration Detected |
|---------------------|---|---|
| Influenza A H1 | A/Denver/1/57 | 1.6 x 10 ² CEID ₅₀ /mL (Influenza A matrix) 1.6 x 10 ⁸ CEID ₅₀ /mL (H1 subtype) |
| | A/Mal/302/54 | 1.58 x 10 ² CEID ₅₀ /mL (Influenza A matrix) 1.58 x 10 ⁵ CEID ₅₀ /mL (H1 subtype) |
| Influenza A H3 | A/Aichi/2/68 H3N2 | 1.58 x 10 ³ CEID ₅₀ /mL |
| | Alice (vaccine) A/England/42/72 | 5 x 10 ⁰ EID ₅₀ /mL (Influenza A matrix) 5 x 10 ¹ EID ₅₀ /mL (H3 subtype) |
| | MRC-2 Recombinant Strain | 8.89 x 10 ² CEID ₅₀ /mL (Influenza A matrix) 8.89 x 10 ³ CEID ₅₀ /mL (H3 subtype) |
| Influenza A H1N1 | A/Washington/24/2012 (A/H1 pdm09) | 3.16 x 10 ³ EID ₅₀ /mL (Influenza A matrix) 3.16 x 10 ² EID ₅₀ /mL (H1-2009 subtype) |
| Influenza A H1N2 | Kilbourne F63: A/NWS/34 (HA) x A/Rockefeller Institute/5/57 (NA), Reassortant NWS-F- Matrix | 8.89 x 10 ¹ CEID ₅₀ /mL (Influenza A matrix) No subtype detected ^a |
| Influenza A H5N8 | A/Gyrfalcon/Washington/41088- 6/2014 BPL | 1.58 x 10 ³ EID ₅₀ /mL (Influenza A matrix) No subtype detected ^b |
| Influenza A H5N2 | A/Northern Pintail/Washington/40964/2014 BPL | 2.51 x 10 ³ EID ₅₀ /mL (Influenza A matrix) No subtype detected ^b |
| Influenza A H7N9 | A/ANHUI/1/2013 | 7.94 x 10 ³ EID ₅₀ /mL (Influenza A matrix) No subtype detected ^c |
| Influenza A H3N2v | A/Indiana/21/2012 | 2.51 x 10 ⁴ EID ₅₀ /mL (Influenza A matrix and H3 subtype) |

^a In silico analysis revealed little homology between this non-contemporary strain sequence and the H1 Signal Probe/Capture Probe sequences.

^b Detection of the H5 Subtype not expected

^c Detection of the H7 Subtype not expected

NOTE: CEID₅₀/mL= Chick Embryo Infectious Dose; EID₅₀/mL= Egg Infectious Dose

Supplemental Analytical Reactivity (Inclusivity) for Influenza A

For human, avian, and swine influenza strains not available for testing on the **cobas® eplex** RP panel, in silico analysis was performed. Bioinformatics analysis was used to predict a result based on the number and location of mismatches in the primers, capture probes, and signal probes found in the **cobas® eplex** RP panel relative to an alignment of GenBank sequences.

Table 43: Predicted (in silico) Reactivity (Inclusivity) Results for Influenza A

| Influenza A Subtype | Host | Strain | GenBank ID | Predicted cobas® eplex Result |
|---------------------|-------|---|------------|-------------------------------|
| H2N2 | Human | A/Albany/20/1957(H2N2) | CY022014 | Influenza A |
| | | Kilbourne F38: A/Korea/426/68 (HA, NA) x A/Puerto Rico/8/34 | CY037296 | Influenza A |
| | Avian | A/chicken/New York/13828-3/1995(H2N2) | CY014822 | Influenza A |
| | | A/Japan/305/1957(H2N2) | CY014977 | Influenza A |
| | | A/Korea/426/1968(H2N2) | CY031596 | Influenza A |
| H4N6 | Avian | A/Blue-winged teal/Minnesota/Sg-00043/2007(H4N6) | CY063978 | Influenza A |
| H5N1 | | A/Peregrine falcon/Aomori/7/2011 | AB629716 | Influenza A |
| | | A/Chicken/West Bengal/239022/2010 | CY061305 | Influenza A |
| | | A/Chicken/West Bengal/193936/2009 | GU272009 | Influenza A |
| | | A/Chicken/Hunan/1/2009 | HM172150 | Influenza A |

| Influenza A Subtype | Host | Strain | GenBank ID | Predicted cobas® eplex Result |
|---------------------|-------|--|------------|-------------------------------|
| | | A/Chicken/Hunan/8/2008 | GU182162 | Influenza A |
| | | A/Chicken/West Bengal/106181/2008 | GU083632 | Influenza A |
| | | A/Chicken/Primorsky/85/2008 | FJ654298 | Influenza A |
| | | A/Chicken/West Bengal/82613/2008 | GU083648 | Influenza A |
| | | A/Duck/France/080036/2008 | CY046185 | Influenza A |
| | | A/Duck/Vietnam/G12/2008 | AB593450 | Influenza A |
| | | A/Chicken/Thailand/PC-340/2008 | EU620664 | Influenza A |
| | | A/Great egret/Hong Kong/807/2008 | CY036240 | Influenza A |
| | | A/Rook/Rostov-on-Don/26/2007(H5N1) | EU814504 | Influenza A |
| | | A/Turkey/VA/505477-18/2007(H5N1) | GU186510 | Influenza A |
| | | A/Chicken/Bangladesh/1151-10/2010(H5N1) | HQ156766 | Influenza A |
| | Human | A/Bangladesh/3233/2011 | CY088772 | Influenza A |
| | | A/Cambodia/R0405050/2007(H5N1) | HQ200572 | Influenza A |
| | | A/Cambodia/S1211394/2008 | HQ200597 | Influenza A |
| | | A/Hong Kong/486/97(H5N1) | AF255368 | Influenza A |
| | Swine | A/Swine/East Java/UT6010/2007(H5N1) | HM440124 | Influenza A |
| H5N2 | Avian | A/Duck/Pennsylvania/10218/1984(H5N2) | AB286120 | Influenza A |
| | | A/American black duck/Illinois/08OS2688/2008 | CY079453 | Influenza A |
| | | A/American green-winged teal/California/HKWF609/2007 | CY033447 | Influenza A |
| | | A/Canada goose/New York/475813-2/2007 | GQ923358 | Influenza A |
| | | A/Blue-winged teal/Saskatchewan/22542/2007 | CY047705 | Influenza A |
| | | A/Chicken/Taiwan/A703-1/2008 | AB507267 | Influenza A |
| | | A/Duck/France/080032/2008 | CY046177 | Influenza A |
| | | A/Duck/New York/481172/2007 | GQ117202 | Influenza A |
| | | A/Gadwall/Altai/1202/2007 | CY049759 | Influenza A |
| | | A/Mallard/Louisiana/476670-4/2007 | GQ923390 | Influenza A |
| | | A/Waterfowl/Colorado/476466-2/2007 | GQ923374 | Influenza A |
| H5N3 | Avian | A/Duck/Singapore/F119/3/1997(H5N3) | GU052803 | Influenza A |
| H6N1 | | A/Duck/PA/486/1969(H6N1) | EU743287 | Influenza A |
| H6N2 | | A/Mallard/Czech Republic/15902-17K/2009(H6N2) | HQ244433 | Influenza A |
| H7N2 | Avian | A/Chicken/Hebei/1/2002 | AY724263 | Influenza A |
| | | A/Chicken/PA/149092-1/02 | AY241609 | Influenza A |
| | | A/Chicken/NJ/294508-12/2004 | EU743254 | Influenza A |
| | | A/Chicken/New York/23165-6/2005 | CY031077 | Influenza A |
| | | A/Muscovy duck/New York/23165-13/2005 | CY033226 | Influenza A |
| | | A/Muscovy duck/New York/87493-3/2005 | CY034791 | Influenza A |
| | | A/Mallard/Netherlands/29/2006 | CY043833 | Influenza A |
| | | A/Northern shoveler/California/JN1447/2007 | CY076873 | Influenza A |
| | Human | A/New York/107/2003(H7N2) | EU587373 | Influenza A |

| Influenza A Subtype | Host | Strain | GenBank ID | Predicted cobas® eplex Result |
|--------------------------------|----------|---|---------------------|-------------------------------|
| H7N3 | | A/Canada/rv504/2004(H7N3) | CY015007 | Influenza A |
| H7N7 | Avian | A/American green-winged teal/Mississippi/09OS046/2009 | CY079309 | Influenza A |
| | | A/Chicken/Germany/R28/03 | AJ619676 | Influenza A |
| | | A/Chicken/Netherlands/1/03 | AY340091 | Influenza A |
| | | A/Mallard/California/HKWF1971/2007 | CY033383 | Influenza A |
| | | A/Mallard/Korea/GH171/2007 | FJ959087 | Influenza A |
| | | A/Mute swan/Hungary/5973/2007 | GQ240816 | Influenza A |
| | | A/Northern shoveler/Mississippi/ 09OS643/2009 | CY079413 | Influenza A |
| | Human | A/Netherlands/219/03(H7N7) | AY340089 | Influenza A |
| H7N9 | Human | A/Shanghai/1/2013(H7N9) | EPI439493 | Influenza A |
| | Avian | A/Northern shoveler/Mississippi/11OS145/2011(H7N9) | CY133650 | Influenza A |
| | | A/Ruddy turnstone/Delaware Bay/220/1995(H7N9) | CY127254 | Influenza A |
| | | A/Turkey/Minnesota/1/1988(H7N9) | CY014787 | Influenza A |
| | | A/Blue-winged teal/Ohio/566/2006(H7N9) | CY024819 | Influenza A |
| H9N2 | Human | A/Hong Kong/1073/99(H9N2) | AJ278647 | Influenza A |
| | Avian | A/Turkey/Wisconsin/1/1966(H9N2) | CY014664 | Influenza A |
| H10N7 | | A/chicken/Germany/N/1949(H10N7) | GQ176135 | Influenza A |
| H11N9 | | A/Duck/Memphis/546/1974(H11N9) | GQ257441 | Influenza A |
| H1N1 | Swine | A/Swine/Wisconsin/1/1971(H1N1) | CY022414 | Influenza A |
| | Human | A/California/UR06-0393/2007(H1N1) | CY026540 | Influenza A H1 |
| CY026539 | | | | |
| A/New York/297/2003(H1N2) | | CY002664 | Influenza A H1 | |
| | | CY002665 | | |
| A/Aalborg/INS133/2009(H1N1) | | CY063606 | Influenza A H1-2009 | |
| | | CY063607 | | |
| A/South Carolina/02/2010(H1N1) | | KC781370 | Influenza A H1-2009 | |
| | | KC781372 | | |
| H1N2 | Swine | A/Swine/Hong Kong/NS857/2001(H1N2) | GQ229350 | Influenza A |
| | | A/Swine/Sweden/1021/2009(H1N2) | GQ495135 | Influenza A |
| H3N1 | Avian | A/Blue-winged teal/ALB/452/1983(H3N1) | CY004635 | Influenza A |
| H3N2v | Human | A/Iowa/07/2011(H3N2) | JQ070760 | Influenza A H3 |
| | | | JQ290177 | |
| | | A/Iowa/08/2011(H3N2) | JQ070768 | Influenza A H3 |
| | | | JQ290167 | |
| | | A/Iowa/09/2011(H3N2) | JQ070776 | Influenza A H3 |
| | | | JQ290183 | |
| | | A/Indiana/08/2011(H3N2) | JQ070800 | Influenza A H3 |
| | | | JQ070795 | |
| A/Maine/06/2011(H3N2) | JN866181 | Influenza A H3 | | |

| Influenza A Subtype | Host | Strain | GenBank ID | Predicted cobas® eplex Result |
|---------------------|---|---|-------------|-------------------------------|
| | | | JN866186 | |
| | | A/Maine/07/2011(H3N2) | JN992746 | Influenza A |
| | | A/Pennsylvania/09/2011(H3N2) | JN655534 | Influenza A |
| | | A/Pennsylvania/11/2011(H3N2) | JN655540 | Influenza A |
| | | A/Pennsylvania/10/2011(H3N2) | JN655550 | Influenza A |
| | | A/West Virginia/06/2011(H3N2) | JQ290159 | Influenza A H3 |
| | | | JQ290164 | |
| | | A/West Virginia/07/2011(H3N2) | JQ348839 | Influenza A |
| | | A/Indiana/10/2011(H3N2) | KJ942592 | Influenza A H3 |
| | | | JQ070787 | |
| | | A/Boston/38/2008(H3N2) | CY044580 | Influenza A H3 |
| | | | CY044581 | |
| | Swine | A/swine/NY/A01104005/2011(H3N2v) | JN940422 | Influenza A H3 |
| | | A/Maine/06/2011(H3N2) | JN866181 | Influenza A H3 |
| | | | JN866186 | Influenza A H3 |
| | | A/Indiana/08/2011(H3N2) | JN655558 | Influenza A H3 |
| | JN638733 | | | |
| Avian | A/American black duck/North Carolina/675-075/2004(H3N2) | GU051135 | Influenza A | |
| | | GU051136 | Influenza A | |
| | A/Mallard/Netherlands/2/1999(H3N5) | CY060261 | Influenza A | |
| | | CY060264 | Influenza A | |
| | A/American black duck/New Brunswick/25182/2007(H3N6) | CY047696 | Influenza A | |
| | | CY047697 | Influenza A | |
| | A/Northern shoveler/California/HKWF1367/2007(H3N7) | CY033372 | Influenza A | |
| | | CY033375 | Influenza A | |
| | H3N8 | A/American black duck/Washington/699/1978(H3N8) | GU052300 | Influenza A H3 |
| | | | GU052299 | |

Table 44: Analytical Reactivity (Inclusivity) Results for Influenza B

| Influenza B Subtype | Strain | Concentration | Multiple of LoD Detected |
|--------------------------------|--------------------|---|--------------------------|
| Influenza B (Yamagata Lineage) | B/Lee/40 | 3×10^{-1} TCID ₅₀ /mL | 3x |
| | B/Allen/45 | 1×10^0 TCID ₅₀ /mL | 10x ^a |
| | B/Maryland/1/59 | 1×10^1 TCID ₅₀ /mL | 100x ^a |
| | B/Taiwan/2/62 | 1×10^1 TCID ₅₀ /mL | 100x ^a |
| Influenza B (Victoria Lineage) | B/Hong Kong/5/72 | 1×10^1 TCID ₅₀ /mL | 100x ^b |
| | B/Malaysia/2506/04 | 3×10^{-1} TCID ₅₀ /mL | 3x |
| Influenza B (Lineage unknown) | B/GL/1739/54 | 3×10^{-1} TCID ₅₀ /mL | 3x |

^a No sequence data available. Lower sensitivity may be a result of mismatches in the assay primers and/or probes. In addition, the reduced sensitivity may be the result of incorrect estimation of genetic material present in the culture of this or the reference strain (TCID₅₀ value is based only on infectious virus particles).

^b In silico analysis revealed that lower sensitivity may be a result of mismatches in the assay primers and/or probes.

Table 45: Analytical Reactivity (Inclusivity) Results for Parainfluenza Virus

| Parainfluenza Subtype | Strain | Concentration | Multiple of LoD Detected |
|-----------------------|--------|--|--------------------------|
| Parainfluenza Virus 1 | C35 | 1.2×10^0 TCID ₅₀ /mL | 3x |
| Parainfluenza Virus 2 | Greer | 1.5×10^2 TCID ₅₀ /mL | 3x |
| Parainfluenza Virus 3 | C-243 | 5×10^1 TCID ₅₀ /mL | 10x ^a |
| Parainfluenza Virus 4 | 4b | 9×10^1 TCID ₅₀ /mL | 3x |

^a In silico analysis revealed that lower sensitivity may be a result of mismatches in the assay primers and/or probes.

Table 46: Analytical Reactivity (Inclusivity) Results for Respiratory Syncytial Virus

| RSV Subtype | Strain | Concentration | Multiple of LoD Detected |
|-------------------------------|---------------|---|--------------------------|
| Respiratory Syncytial Virus A | A2 | 4.5×10^0 TCID ₅₀ /mL | 3x |
| | Long | 4.5×10^0 TCID ₅₀ /mL | 3x |
| Respiratory Syncytial Virus B | 9320 | 6×10^{-1} TCID ₅₀ /mL | 3x |
| | Wash/18537/62 | 6×10^{-1} TCID ₅₀ /mL | 3x |
| | WV/14617/85 | 6×10^{-1} TCID ₅₀ /mL | 3x |

Table 47: Analytical Reactivity (Inclusivity) Results for Chlamydia pneumoniae

| | Strain | Concentration | Multiple of LoD Detected |
|-----------------------------|------------------|------------------------|--------------------------|
| <i>Chlamydia pneumoniae</i> | CWL-029 | 9×10^2 CFU/mL | 3x |
| | TWAR strain 2043 | 9×10^2 CFU/mL | 3x |

Table 48: Analytical Reactivity (Inclusivity) Results for Mycoplasma pneumoniae

| | Strain | Concentration | Multiple of LoD Detected |
|------------------------------|-----------|------------------------|--------------------------|
| <i>Mycoplasma pneumoniae</i> | [Bru] | 9×10^2 CCU/mL | 3x |
| | M129-B170 | 9×10^2 CCU/mL | 3x |
| | M129-B7 | 9×10^2 CCU/mL | 3x |
| | [M52] | 9×10^2 CCU/mL | 3x |
| | [Mac] | 9×10^2 CCU/mL | 3x |
| | Mutant 22 | 3×10^4 CCU/mL | 100x ^a |
| | PI 1428 | 3×10^4 CCU/mL | 100x ^b |

^a No sequence data available. Lower sensitivity may be a result of mismatches in the assay primers and/or probes. In addition, the reduced sensitivity may be the result of incorrect estimation of genetic material present in the culture of this or the reference strain (CCU/ml value is based only on live bacteria).

^b In silico analysis revealed good homology to primers and probes. The reduced sensitivity is likely the result of incorrect estimation of genetic material present in the culture of this or the reference strain (CCU/ml value is based only on live bacteria).

Analytical Specificity (Cross-Reactivity and Exclusivity)

Cross-reactivity of each viral and bacterial target on the cobas® eplex RP panel was evaluated at high concentrations (1×10^5 TCID₅₀/mL or $>1 \times 10^5$ EID₅₀/mL for viruses, 1×10^6 CFU/mL or CCU/mL for bacterial isolates, or 1×10^6 copies/mL for in vitro transcripts) of quantified strains/isolates diluted in viral transport media. In vitro transcript for coronavirus HKU1 was diluted in PBS. Table 49 summarizes the results of the on-panel viral and bacterial strains/isolates tested. No cross-reactivity was observed between any of the on-panel viruses or bacteria.

Table 49: Cross-reactivity with cobas® eplex RP panel Target Organisms

| Target | Strain | Concentration | Cross-Reactivity Results |
|--------------------------------|--|--|--------------------------|
| Adenovirus A | Type 31 | 1 x 10 ⁵ TCID ₅₀ /mL | Not observed |
| Adenovirus B | Type 7A | 1 x 10 ⁵ TCID ₅₀ /mL | Not observed |
| Adenovirus C | Type 1 | 1 x 10 ⁵ TCID ₅₀ /mL | Not observed |
| Adenovirus D | Type 9 | 1 x 10 ⁵ TCID ₅₀ /mL | Not observed |
| Adenovirus E | Type 4 | 1 x 10 ⁵ TCID ₅₀ /mL | Not observed |
| Adenovirus F | Type 41 | 1 x 10 ⁵ TCID ₅₀ /mL | Not observed |
| Coronavirus | 229E | 1 x 10 ⁵ TCID ₅₀ /mL | Not observed |
| Coronavirus | HKU1 in vitro transcript | 1 x 10 ⁶ copies/mL | Not observed |
| Coronavirus | NL63 | 1 x 10 ⁵ TCID ₅₀ /mL | Not observed |
| Coronavirus | OC43 | 1 x 10 ⁵ TCID ₅₀ /mL | Not observed |
| Enterovirus | Type 68 2007 isolate | 1 x 10 ⁵ TCID ₅₀ /mL | Not observed |
| Human metapneumovirus | B1 | 1 x 10 ⁵ TCID ₅₀ /mL | Not observed |
| Human rhinovirus | 1A | 1 x 10 ⁵ TCID ₅₀ /mL | Not observed |
| Influenza A | A/Brisbane/59/07 | 1 x 10 ⁵ TCID ₅₀ /mL | Not observed |
| Influenza A H1 | A/Brisbane/59/07 | 1 x 10 ⁵ TCID ₅₀ /mL | Not observed |
| Influenza A H1-2009 | A/NY/01/2009 | 1 x 10 ⁵ TCID ₅₀ /mL | Not observed |
| Influenza A H3 | A/Brisbane/10/07 | 1 x 10 ⁵ TCID ₅₀ /mL | Not observed |
| Influenza A H3N2v ^a | A/Indiana/21/2012 | 2.51 x 10 ⁵ EID ₅₀ /mL | Not observed |
| Influenza A H5N2 ^b | A/Northern Pintail Washington/40964/14BPL | 2.51 x 10 ⁵ EID ₅₀ /mL | Not observed |
| Influenza A H5N8 ^c | A/Gyrfalcon/Washington /410886/2014 BPL | 1.58 x 10 ⁵ EID ₅₀ /mL | Not observed |
| Influenza A H7N9 ^d | A/ANHUI/1/2013 | 7.94 x 10 ⁵ EID ₅₀ /mL | Not observed |
| Influenza B | B/Florida/02/06 | 1 x 10 ⁵ TCID ₅₀ /mL | Not observed |
| Parainfluenza Virus 1 | C35 | 1 x 10 ⁵ TCID ₅₀ /mL | Not observed |
| Parainfluenza Virus 2 | Type 2 | 1 x 10 ⁵ TCID ₅₀ /mL | Not observed |
| Parainfluenza Virus 3 | Type 3 | 1 x 10 ⁵ TCID ₅₀ /mL | Not observed |
| Parainfluenza Virus 4 | Type 4a | 1 x 10 ⁵ TCID ₅₀ /mL | Not observed |
| RSV A | 2006 Isolate | 1 x 10 ⁵ TCID ₅₀ /mL | Not observed |
| RSV B | CH93(18)-18 | 1 x 10 ⁵ TCID ₅₀ /mL | Not observed |
| <i>Chlamydia pneumoniae</i> | AR-39 | 1 x 10 ⁶ CFU/mL | Not observed |
| <i>Mycoplasma pneumoniae</i> | FH strain of Eaton Agent [NCTC 10119] | 1 x 10 ⁶ CCU/mL | Not observed |

^a Influenza A H3N2v detected as Influenza A, Influenza A H3^b Influenza A H5N2 detected as Influenza A^c Influenza A H5N8 detected as Influenza A^d Influenza A H7N9 detected as Influenza A

Cross-reactivity of viruses, bacteria, and fungi that are not targets on the cobas® eplex RP panel was evaluated at high concentrations (1 x 10⁵ TCID₅₀/mL or copies/mL for viruses, 1 x 10⁶ CFU/mL for bacterial and yeast isolates, or 1 x 10⁶ copies/mL for plasmid DNA or genomic RNA) by diluting quantified strains in viral transport media. Plasmid for bocavirus and genomic RNA for MERS coronavirus (MERS-CoV) were diluted in PBS. Table 50 summarizes the results of the strains tested. No cross-reactivity was observed between any of the off-panel viruses, bacteria or fungi with the cobas® eplex RP panel targets.

Table 50: Cross-reactivity with Organisms Not Detected by the cobas® eplex RP panel (Exclusivity)

| Target | Strain | Concentration | Cross-Reactivity Results |
|--|--------------------------------|--|--------------------------|
| <i>Acinetobacter baumannii</i> | ATCC® 19606 | 1 x 10 ⁶ CFU/mL | Not observed |
| <i>Bordetella pertussis</i> | 18323 [NCTC 10739] | 1 x 10 ⁶ CFU/mL | Not observed |
| <i>Bordetella parapertussis</i> | ATCC 15311 | 1 x 10 ⁶ CFU/mL | Not observed |
| <i>Burkholderia cepacia</i> | ATCC 25416 | 1 x 10 ⁶ CFU/mL | Not observed |
| <i>Candida albicans</i> | ATCC 10231 | 1 x 10 ⁶ CFU/mL | Not observed |
| <i>Candida glabrata</i> | ATCC 15126 | 1 x 10 ⁶ CFU/mL | Not observed |
| MERS Coronavirus (MERS-CoV) | EMC/2012 ^a | 1 x 10 ⁵ copies/mL | Not observed |
| <i>Corynebacterium diphtheriae</i> | ATCC 13812 | 1 x 10 ⁶ CFU/mL | Not observed |
| Cytomegalovirus | AD 169 | 1 x 10 ⁵ TCID ₅₀ /mL | Not observed |
| Epstein Barr Virus | Strain B95-8 | 1 x 10 ⁵ TCID ₅₀ /mL | Not observed |
| <i>Escherichia coli</i> | ATCC 10279 | 1 x 10 ⁶ CFU/mL | Not observed |
| <i>Haemophilus influenzae</i> | ATCC 43065 | 1 x 10 ⁶ CFU/mL | Not observed |
| Herpes Simplex Virus | Isolate 2 | 1 x 10 ⁵ TCID ₅₀ /mL | Not observed |
| Human bocavirus | Bocavirus plasmid ^b | 1 x 10 ⁶ copies/mL | Not observed |
| <i>Klebsiella pneumoniae</i> | ATCC 51504 | 1 x 10 ⁶ CFU/mL | Not observed |
| <i>Lactobacillus acidophilus</i> | ATCC 314 | 1 x 10 ⁶ CFU/mL | Not observed |
| <i>Lactobacillus plantarum</i> | ATCC 8014 | 1 x 10 ⁶ CFU/mL | Not observed |
| <i>Legionella pneumophila</i> | Philadelphia-1 | 1 x 10 ⁶ CFU/mL | Not observed |
| Measles | N/A | 1 x 10 ⁵ TCID ₅₀ /mL | Not observed |
| <i>Moraxella catarrhalis</i> | ATCC 23246 | 1 x 10 ⁶ CFU/mL | Not observed |
| Mumps | Isolate 2 | 1 x 10 ⁵ TCID ₅₀ /mL | Not observed |
| <i>Mycobacterium tuberculosis</i> | ATCC 25177 | 1 x 10 ⁶ CFU/mL | Not observed |
| <i>Neisseria meningitidis</i> | ATCC 13077 | 1 x 10 ⁶ CFU/mL | Not observed |
| <i>Neisseria sicca</i> | ATCC 29193 | 1 x 10 ⁶ CFU/mL | Not observed |
| <i>Porphyromonas gingivalis</i> | ATCC 33277 | 1 x 10 ⁶ CFU/mL | Not observed |
| <i>Proteus vulgaris</i> | ATCC 33420 | 1 x 10 ⁶ CFU/mL | Not observed |
| <i>Pseudomonas aeruginosa</i> | ATCC 15442 | 1 x 10 ⁶ CFU/mL | Not observed |
| <i>Serratia marcescens</i> | ATCC 13880 | 1 x 10 ⁶ CFU/mL | Not observed |
| <i>Staphylococcus aureus</i> (MRSA) | NRS384 | 1 x 10 ⁶ CFU/mL | Not observed |
| <i>Staphylococcus aureus</i> (MSSA) | ATCC 25923 | 1 x 10 ⁶ CFU/mL | Not observed |
| <i>Staphylococcus epidermidis</i> (MRSE) | ATCC 35983 | 1 x 10 ⁶ CFU/mL | Not observed |
| <i>Staphylococcus epidermidis</i> (MSSE) | ATCC 49134 | 1 x 10 ⁶ CFU/mL | Not observed |
| <i>Staphylococcus haemolyticus</i> | ATCC 29970 | 1 x 10 ⁶ CFU/mL | Not observed |
| <i>Streptococcus agalactiae</i> | ATCC 12401 | 1 x 10 ⁶ CFU/mL | Not observed |
| <i>Streptococcus dysgalactiae</i> | ATCC 35666 | 1 x 10 ⁶ CFU/mL | Not observed |
| <i>Streptococcus mitis</i> | ATCC 15914 | 1 x 10 ⁶ CFU/mL | Not observed |
| <i>Streptococcus pneumoniae</i> | ATCC 49619 | 1 x 10 ⁶ CFU/mL | Not observed |
| <i>Streptococcus pyogenes</i> | ATCC 12384 | 1 x 10 ⁶ CFU/mL | Not observed |
| <i>Streptococcus salivarius</i> | ATCC 13419 | 1 x 10 ⁶ CFU/mL | Not observed |
| Varicella Zoster Virus | 82 | 8.9 x 10 ³ TCID ₅₀ /mL | Not observed |

^a Extracted genomic RNA^b Plasmid does not contain full length viral genome.

Reproducibility

A multisite reproducibility study of the **cobas® eplex** RP panel was performed to evaluate agreement with expected results across major sources of variability, such as site-to-site, lot-to-lot, day-to-day, and operator-to-operator. Testing occurred at 3 sites (2 external, 1 internal) on one **cobas® eplex** system per site with either 3 or 4 towers. Two operators performed testing at each site on 6 days (5 nonconsecutive days) with 3 unique lots of RP Panel cartridges. A reproducibility panel consisting of 3 panel members with 6 organisms (representing 7 RP Panel targets) at 3 concentrations (moderate positive- 3x LoD, low positive- 1x LoD, and negative) was tested in triplicate. The 6 organisms tested included adenovirus, coronavirus, human metapneumovirus, influenza A H3, parainfluenza virus 1, and RSV A; organisms were diluted in natural clinical matrix (pooled, negative nasopharyngeal swab samples). Negative samples consisted of natural clinical matrix only. Each simulated sample was divided into aliquots and stored frozen (-70 °C) prior to testing. Each operator tested 9 samples (3 member reproducibility panel in triplicate) each day; each panel member was tested 108 times (3 replicates x 3 sites x 2 operators x 3 lots x 2 days of testing/operator/lot) for a maximum of 324 tests. After completion of initial and repeat testing for invalid results, 1 low positive sample tested at Site 3 had an invalid result and was excluded from reproducibility performance analyses.

Percent agreement (95% CI) with expected results was 100% for all 7 targets for the moderate positive and negative panel, and 100% for 6 of 7 low positive panel targets (coronavirus, human metapneumovirus, influenza A, influenza A H3, parainfluenza 1, and RSV A); percent agreement was 91.6% for adenovirus. Summary results for the 7 **cobas® eplex** RP panel targets that correspond to the 6 organisms in the reproducibility panel are provided in Table 51 to Table 57. Summary results for the 10 **cobas® eplex** RP panel targets that did not have organisms included in the reproducibility panel are provided in Table 58.

Table 51: Percent Agreement for Adenovirus

| Adenovirus Concentration | Site | Agreement with Expected Results | | |
|---|------|---------------------------------|------|-------------|
| | | Agreed / N | % | 95% CI |
| Moderate Positive 3x LoD 6 x 10 ⁹ TCID ₅₀ /mL | 1 | 36/36 | 100 | (90.4-100) |
| | 2 | 36/36 | 100 | (90.4-100) |
| | 3 | 36/36 | 100 | (90.4-100) |
| | All | 108/108 | 100 | (96.6-100) |
| Low Positive 1x LoD 2 x 10 ⁹ TCID ₅₀ /mL | 1 | 36/36 | 100 | (90.4-100) |
| | 2 | 34/36 | 94.4 | (81.9-98.5) |
| | 3 | 28/35 | 80.0 | (64.1-90.0) |
| | All | 98/107 | 91.6 | (84.8-95.5) |
| Negative | 1 | 36/36 | 100 | (90.4-100) |
| | 2 | 36/36 | 100 | (90.4-100) |
| | 3 | 36/36 | 100 | (90.4-100) |
| | All | 108/108 | 100 | (96.6-100) |

CI=Confidence Interval

Table 52: Percent Agreement for Coronavirus

| Coronavirus Concentration | Site | Agreement with Expected Results | | |
|---|------|---------------------------------|-----|------------|
| | | Agreed / N | % | 95% CI |
| Moderate Positive 3x LoD 1.5 x 10 ³ TCID ₅₀ /mL | 1 | 36/36 | 100 | (90.4-100) |
| | 2 | 36/36 | 100 | (90.4-100) |
| | 3 | 36/36 | 100 | (90.4-100) |
| | All | 108/108 | 100 | (96.6-100) |
| Low Positive 1x LoD 5 x 10 ² TCID ₅₀ /mL | 1 | 36/36 | 100 | (90.4-100) |
| | 2 | 36/36 | 100 | (90.4-100) |
| | 3 | 35/35 | 100 | (90.1-100) |
| | All | 107/107 | 100 | (96.5-100) |
| Negative | 1 | 36/36 | 100 | (90.4-100) |
| | 2 | 36/36 | 100 | (90.4-100) |
| | 3 | 36/36 | 100 | (90.4-100) |
| | All | 108/108 | 100 | (96.6-100) |

Table 53: Percent Agreement for Human Metapneumovirus (hMPV)

| hMPV Concentration | Site | Agreement with Expected Results | | |
|--|------|---------------------------------|-----|------------|
| | | Agreed / N | % | 95% CI |
| Moderate Positive 3x LoD 6.75 x 10 ² TCID ₅₀ /mL | 1 | 36/36 | 100 | (90.4-100) |
| | 2 | 36/36 | 100 | (90.4-100) |
| | 3 | 36/36 | 100 | (90.4-100) |
| | All | 108/108 | 100 | (96.6-100) |
| Low Positive 1x LoD 2.25 x 10 ² TCID ₅₀ /mL | 1 | 36/36 | 100 | (90.4-100) |
| | 2 | 36/36 | 100 | (90.4-100) |
| | 3 | 35/35 | 100 | (90.1-100) |
| | All | 107/107 | 100 | (96.5-100) |
| Negative | 1 | 36/36 | 100 | (90.4-100) |
| | 2 | 36/36 | 100 | (90.4-100) |
| | 3 | 36/36 | 100 | (90.4-100) |
| | All | 108/108 | 100 | (96.6-100) |

Table 54: Percent Agreement for Influenza A

| Influenza A Concentration | Site | Agreement with Expected Results | | |
|---|------|---------------------------------|-----|------------|
| | | Agreed / N | % | 95% CI |
| Moderate Positive 3x LoD 1.5 x 10 ² TCID ₅₀ /mL | 1 | 36/36 | 100 | (90.4-100) |
| | 2 | 36/36 | 100 | (90.4-100) |
| | 3 | 36/36 | 100 | (90.4-100) |
| | All | 108/108 | 100 | (96.6-100) |
| Low Positive 1x LoD 5 x 10 ¹ TCID ₅₀ /mL | 1 | 36/36 | 100 | (90.4-100) |
| | 2 | 36/36 | 100 | (90.4-100) |
| | 3 | 35/35 | 100 | (90.1-100) |
| | All | 107/107 | 100 | (96.5-100) |
| Negative | 1 | 36/36 | 100 | (90.4-100) |
| | 2 | 36/36 | 100 | (90.4-100) |
| | 3 | 36/36 | 100 | (90.4-100) |
| | All | 108/108 | 100 | (96.6-100) |

Table 55: Percent Agreement for Influenza A H3

| Influenza A H3 Concentration | Site | Agreement with Expected Results | | |
|---|------|---------------------------------|-----|------------|
| | | Agreed / N | % | 95% CI |
| Moderate Positive 3x LoD 1.5 x 10 ² TCID ₅₀ /mL | 1 | 36/36 | 100 | (90.4-100) |
| | 2 | 36/36 | 100 | (90.4-100) |
| | 3 | 36/36 | 100 | (90.4-100) |
| | All | 108/108 | 100 | (96.6-100) |
| Low Positive 1x LoD 5 x 10 ¹ TCID ₅₀ /mL | 1 | 36/36 | 100 | (90.4-100) |
| | 2 | 36/36 | 100 | (90.4-100) |
| | 3 | 35/35 | 100 | (90.1-100) |
| | All | 107/107 | 100 | (96.5-100) |
| Negative | 1 | 36/36 | 100 | (90.4-100) |
| | 2 | 36/36 | 100 | (90.4-100) |
| | 3 | 36/36 | 100 | (90.4-100) |
| | All | 108/108 | 100 | (96.6-100) |

Table 56: Percent Agreement for Parainfluenza Virus (PIV) 1

| PIV 1 Concentration | Site | Agreement with Expected Results | | |
|---|------|---------------------------------|-----|------------|
| | | Agreed / N | % | 95% CI |
| Moderate Positive 3x LoD 1.2 x 10 ⁰ TCID ₅₀ /mL | 1 | 36/36 | 100 | (90.4-100) |
| | 2 | 36/36 | 100 | (90.4-100) |
| | 3 | 36/36 | 100 | (90.4-100) |
| | All | 108/108 | 100 | (96.6-100) |
| Low Positive 1x LoD 4 x 10 ⁻¹ TCID ₅₀ /mL | 1 | 36/36 | 100 | (90.4-100) |
| | 2 | 36/36 | 100 | (90.4-100) |
| | 3 | 35/35 | 100 | (90.1-100) |
| | All | 107/107 | 100 | (96.5-100) |
| Negative | 1 | 36/36 | 100 | (90.4-100) |
| | 2 | 36/36 | 100 | (90.4-100) |
| | 3 | 36/36 | 100 | (90.4-100) |
| | All | 108/108 | 100 | (96.6-100) |

Table 57: Percent Agreement for Respiratory Syncytial Virus (RSV) A

| RSV A Concentration | Site | Agreement with Expected Results | | |
|---|------|---------------------------------|-----|------------|
| | | Agreed / N | % | 95% CI |
| Moderate Positive 3x LoD 4.5 x 10 ⁰ TCID ₅₀ /mL | 1 | 36/36 | 100 | (90.4-100) |
| | 2 | 36/36 | 100 | (90.4-100) |
| | 3 | 36/36 | 100 | (90.4-100) |
| | All | 108/108 | 100 | (96.6-100) |
| Low Positive 1x LoD 1.5 x 10 ⁰ TCID ₅₀ /mL | 1 | 36/36 | 100 | (90.4-100) |
| | 2 | 36/36 | 100 | (90.4-100) |
| | 3 | 35/35 | 100 | (90.1-100) |
| | All | 107/107 | 100 | (96.5-100) |
| Negative | 1 | 36/36 | 100 | (90.4-100) |
| | 2 | 36/36 | 100 | (90.4-100) |
| | 3 | 36/36 | 100 | (90.4-100) |
| | All | 108/108 | 100 | (96.6-100) |

Table 58: Negative Percent Agreement with Organisms Not Included in the Reproducibility Panel

| Target | Site | Agreement with Expected Negative Results | | |
|-------------------------------|------|--|------|-------------|
| | | Agreed / N | % | 95% CI |
| Human Rhinovirus/Enterovirus | 1 | 108/108 | 100 | (96.6-100) |
| | 2 | 108/108 | 100 | (96.6-100) |
| | 3 | 104/107 | 97.2 | (92.1-99.0) |
| | All | 320/323 | 99.1 | (97.3-99.7) |
| Influenza A H1 | 1 | 108/108 | 100 | (96.6-100) |
| | 2 | 108/108 | 100 | (96.6-100) |
| | 3 | 107/107 | 100 | (96.5-100) |
| | All | 323/323 | 100 | (98.8-100) |
| Influenza A H1-2009 | 1 | 108/108 | 100 | (96.6-100) |
| | 2 | 108/108 | 100 | (96.6-100) |
| | 3 | 107/107 | 100 | (96.5-100) |
| | All | 323/323 | 100 | (98.8-100) |
| Influenza B | 1 | 108/108 | 100 | (96.6-100) |
| | 2 | 108/108 | 100 | (96.6-100) |
| | 3 | 107/107 | 100 | (96.5-100) |
| | All | 323/323 | 100 | (98.8-100) |
| Parainfluenza Virus 2 | 1 | 108/108 | 100 | (96.6-100) |
| | 2 | 108/108 | 100 | (96.6-100) |
| | 3 | 107/107 | 100 | (96.5-100) |
| | All | 323/323 | 100 | (98.8-100) |
| Parainfluenza Virus 3 | 1 | 108/108 | 100 | (96.6-100) |
| | 2 | 108/108 | 100 | (96.6-100) |
| | 3 | 106/107 | 99.1 | (94.9-99.8) |
| | All | 322/323 | 99.7 | (98.3-99.9) |
| Parainfluenza Virus 4 | 1 | 108/108 | 100 | (96.6-100) |
| | 2 | 108/108 | 100 | (96.6-100) |
| | 3 | 107/107 | 100 | (96.5-100) |
| | All | 323/323 | 100 | (98.8-100) |
| Respiratory Syncytial Virus B | 1 | 108/108 | 100 | (96.6-100) |
| | 2 | 108/108 | 100 | (96.6-100) |
| | 3 | 107/107 | 100 | (96.5-100) |
| | All | 323/323 | 100 | (98.8-100) |
| <i>Chlamydia pneumoniae</i> | 1 | 108/108 | 100 | (96.6-100) |
| | 2 | 108/108 | 100 | (96.6-100) |
| | 3 | 107/107 | 100 | (96.5-100) |
| | All | 323/323 | 100 | (98.8-100) |
| <i>Mycoplasma pneumoniae</i> | 1 | 108/108 | 100 | (96.6-100) |
| | 2 | 107/108 | 99.1 | (94.9-99.8) |
| | 3 | 106/107 | 99.1 | (94.9-99.8) |
| | All | 321/323 | 99.4 | (97.8-99.8) |

Samples with Co-Detected Organisms

Detection of more than one clinically relevant viral organism in a sample was evaluated with the **cobas® eplex** RP panel using a natural clinical matrix (pooled, negative nasopharyngeal swab samples) spiked with two RP Panel organisms: one organism at a low concentration (1-3x LoD) and the second organism at a high concentration (1×10^5 TCID₅₀/mL). Table 59 contains the results of co-detection testing which demonstrated the ability of the **cobas® eplex** RP panel to detect 2 organisms in a sample at both high and low concentrations as indicated in the table.

Table 59: Detection of Co-detections

| Organism 1 | High Titer | Organism 2 | Low Titer | Multiple of LoD |
|-----------------------|--|-----------------------|---|-----------------|
| Influenza A H3 | 1×10^5 TCID ₅₀ /mL | Adenovirus B | 2×10^0 TCID ₅₀ /mL | 1x |
| Adenovirus | 1×10^5 TCID ₅₀ /mL | Influenza A H3 | 5×10^1 TCID ₅₀ /mL | 1x |
| Influenza A H3 | 1×10^5 TCID ₅₀ /mL | RSV A | 1.5×10^0 TCID ₅₀ /mL | 1x |
| RSV A | 1×10^5 TCID ₅₀ /mL | Influenza A H3 | 5×10^1 TCID ₅₀ /mL | 1x |
| Influenza A H1-2009 | 1×10^5 TCID ₅₀ /mL | RSV B | 6×10^{-1} TCID ₅₀ /mL | 3x |
| RSV B | 1×10^5 TCID ₅₀ /mL | Influenza A H1-2009 | 1×10^{-1} TCID ₅₀ /mL | 1x |
| Influenza A H1-2009 | 1×10^5 TCID ₅₀ /mL | Rhinovirus | 1.5×10^0 TCID ₅₀ /mL | 1x |
| Rhinovirus | 1×10^5 TCID ₅₀ /mL | Influenza A H1-2009 | 3×10^{-1} TCID ₅₀ /mL | 3x |
| Influenza A H1-2009 | 1×10^5 TCID ₅₀ /mL | Parainfluenza Virus 3 | 5×10^0 TCID ₅₀ /mL | 1x |
| Parainfluenza Virus 3 | 1×10^5 TCID ₅₀ /mL | Influenza A H1-2009 | 1×10^{-1} TCID ₅₀ /mL | 1x |
| Rhinovirus | 1×10^5 TCID ₅₀ /mL | RSV A | 1.5×10^0 TCID ₅₀ /mL | 1x |
| RSV A | 1×10^5 TCID ₅₀ /mL | Rhinovirus | 1.5×10^0 TCID ₅₀ /mL | 1x |
| Coronavirus | 1×10^5 TCID ₅₀ /mL | RSV A | 1.5×10^0 TCID ₅₀ /mL | 1x |
| RSV A | 1×10^5 TCID ₅₀ /mL | Coronavirus | 7.5×10^0 TCID ₅₀ /mL | 1x |
| Human Metapneumovirus | 1×10^5 TCID ₅₀ /mL | Adenovirus | 2×10^0 TCID ₅₀ /mL | 1x |
| Adenovirus | 1×10^5 TCID ₅₀ /mL | Human Metapneumovirus | 2.25×10^2 TCID ₅₀ /mL | 1x |
| Adenovirus | 1×10^5 TCID ₅₀ /mL | RSV A | 1.5×10^0 TCID ₅₀ /mL | 1x |
| RSV A | 1×10^5 TCID ₅₀ /mL | Adenovirus | 2×10^0 TCID ₅₀ /mL | 1x |

Sample Matrix Equivalency

All analytical studies that utilized viral and bacterial cultures close to LoD were performed by spiking the viral and bacterial cultures into a pool of natural negative NPS in VTM as sample matrix. For analytical studies that used viral and bacterial cultures at a concentration which was at least 10x LoD or higher, the viral and bacterial cultures were spiked into MicroTest™ M5® transport media from Remel instead of negative pooled NPS for ease of use. A sample matrix equivalency study was performed to demonstrate equivalency between natural clinical matrix (pooled, negative nasopharyngeal swab in VTM samples) and viral transport media when spiked with targets at a concentration of approximately 10x LoD. Quantified, representative viral and bacterial strains were diluted in a natural clinical matrix (pooled, negative nasopharyngeal swab in VTM samples) and in viral transport media. All samples were tested in duplicate. There was no difference observed in detection of targets in natural clinical matrix vs. viral transport media.

Interfering Substances

Substances commonly found in respiratory samples, substances that could be introduced during specimen collection, or medications commonly used to treat congestion, allergies, or asthma symptoms that could potentially interfere with the **cobas® eplex** RP panel were individually evaluated. To simulate clinical samples, quantified representative viral and bacterial strains were diluted to 1x LoD in a natural clinical matrix (pooled, negative nasopharyngeal swab specimens) and tested in triplicate for negative and positive interference. Natural clinical matrix (pooled, negative nasopharyngeal swab samples) with no organisms added was used as a control. All substances and organisms tested for interference were shown to be compatible with the **cobas® eplex** RP panel. No potentially interfering substances were found to inhibit the **cobas® eplex** RP panel at the concentrations tested in Table 60.

Table 60: List of Substances for Testing

| Potentially Interfering Substance | Active Ingredient | Testing Concentration |
|--|-------------------------------------|--|
| Control Sample Matrix ^a | Becton Dickinson UVT | N/A |
| Transport Medium ^a | Copan eSwab (Liquid Amies media) | N/A |
| Viral Transport Medium ^a | MicroTest M4 | N/A |
| | MicroTest M4-RT | N/A |
| | MicroTest M5 | N/A |
| | MicroTest M6 | N/A |
| Flocked Swabs | Copan Minitip in UVT | N/A |
| | Copan Regular Tip in UVT | N/A |
| Blood (human) | Blood | 2% v/v |
| | Human gDNA | 50 ng/rxn |
| Throat lozenges, oral anesthetic and analgesic | Benzocaine, menthol | 26% w/v |
| Mucin | Purified mucin protein | 1% w/v |
| Nasal sprays or drops | Phenylephrine HCl (Neo-Synephrine®) | 1.5% v/v |
| | Oxymetazoline HCl (Afrin®) | 1% v/v |
| | Sodium chloride | 0.8% w/v |
| Antibacterial, systemic | Tobramycin ^b | 1% w/v |
| Antibiotic, nasal ointment | Mupirocin | 2% w/v |
| Nasal corticosteroids | Beclomethasone | 1.5% w/v |
| | Dexamethasone | 1.5% w/v |
| | Flunisolide | 1.5% w/v |
| | Budesonide (Rhinocort®) | 0.9% v/v |
| | Triamcinolone (Nasacort®) | 1.5% v/v |
| | Fluticasone (Flonase®) | 1.5% v/v |
| ZICAM® Allergy Relief Nasal Gel | Luffa operculata | 1% v/v |
| | Sulfur | |
| | Galphimia glauca | |
| | Histaminum hydrochloricum | |
| Anti-viral drugs | Zanamivir | 550 ng/mL |
| | Oseltamivir | 142 ng/mL |
| Virus | Cytomegalovirus | 1 x 10 ⁵ TCID ₅₀ /mL |

| Potentially Interfering Substance | Active Ingredient | Testing Concentration |
|-----------------------------------|------------------------------------|----------------------------|
| Bacteria | <i>Streptococcus pneumoniae</i> | 1 x 10 ⁶ CFU/mL |
| | <i>Bordetella parapertussis</i> | |
| | <i>Haemophilus influenza</i> | |
| | <i>Staphylococcus aureus</i> | |
| | <i>Neisseria meningitides</i> | |
| | <i>Corynebacterium diphtheriae</i> | |

^a Testing of media was done by adding a negative NPS collected in the specified media and diluting in the natural clinical matrix.

^b At concentrations greater than 1% weight/volume in the sample, tobramycin was found to inhibit assay performance.

Carryover and Cross-contamination

The carryover/cross-contamination rate of the **cobas® eplex** RP panel and **cobas® eplex** system was tested in a checkerboard approach by running high positive and negative samples interspersed in all bays of a four-tower **cobas® eplex** system (24 bays total) over 5 separate runs on 5 separate days. Quantified parainfluenza virus 3 was prepared in viral transport media at a high concentration (1 x 10⁵ TCID₅₀/mL, 20,000x LoD) to simulate a clinically relevant high positive and was tested as a representative target organism. Transport media was used to represent negative samples. On each round of testing, 24 **cobas® eplex** RP panel cartridges were evaluated. 100% of parainfluenza 3-positive samples generated a result of Detected and 100% of parainfluenza 3-negative samples generated a parainfluenza 3 result of No Target Detected, indicating no carryover or cross-contamination was observed between bays or within bays with the **cobas® eplex** RP panel when testing consecutively or in adjacent bays.

TROUBLESHOOTING

Table 61: Troubleshooting Table



















For a complete list of all **cobas® eplex** error messages and a description of the messages, please refer to the **cobas® eplex** User Assistance Manual.

| Error | Error Messages | Description | Re-test Recommendations |
|---------------------|---|--|--|
| Test did not start | Cartridge failure The cartridge initialization test failed Cartridge not present Bay heater failure Unknown error Bay main / fluid motor failure Bay over pressured Bay temperature out of range The system was unable to read the cartridge Cartridge inserted doesn't match the serial number of the cartridge scanned The system is not ready to accept the cartridge The system failed to prepare the cartridge for processing | <p>An error that occurs during pre-run checks (cartridge initialization) of the cartridge upon insertion into the bay. Cartridge initialization occurs when the cartridge is first inserted into the bay and takes approximately 90 seconds.</p> <p>Upon completion of cartridge initialization, the cartridge cannot be restarted, but prior to this point, the cartridge can be restarted.</p> <p>To verify cartridge initialization has completed, examine the cartridge label upon removal from the bay. If the cartridge label has been pierced, the test has already started and cartridge cannot be reused. If the label has not been pierced, follow the recommendation as stated.</p> | <ol style="list-style-type: none"> 1. Remove cartridge from bay. <ol style="list-style-type: none"> a. Reset bay to clear the error b. Restart cartridge in any available bay 2. If the cartridge is not able to be run on the second try and again generates an error during cartridge initialization, this indicates an issue with the cartridge. This cartridge should be discarded following laboratory procedures and the sample should be repeated using a new cartridge. Bay(s) should be reset to clear the errors. Please contact MAS or GenMark Technical Support to alert them of the issue. <p>If the bay remains in an error state (flashing red) after the cartridge has been removed, then the bay must be reset through the Bay Configuration menu before it can be used to run cartridges.</p> |
| Test did not finish | Bay heater failure Bay main / fluid motor failure Bay voltage failure Bay sub-system communication timeout Cartridge failure Bay over pressured Bay auto-calibration failure Bay temperature out of range The system was unable to eject the cartridge from the bay | <p>This type of error occurs during the run, after pre-run checks (cartridge initialization) have finished, and prevents the cartridge from being processed to completion.</p> | <p>Reagents have been consumed and the cartridge cannot be reused. Contact GenMark Technical Support and proceed with repeat testing of the sample using a new cartridge.</p> <p>If the bay remains in an error state (flashing red) after the cartridge has been removed, then the bay must be reset through the Bay Configuration menu before it can be used to run cartridges.</p> |
| Invalid | - | <p>This is an error that results in no valid results being generated. A test report will be generated, but all targets and the internal control will be invalid.</p> | <p>Reagents have been consumed and the cartridge cannot be reused. Contact GenMark Technical Support and proceed with repeat testing of the sample using a new cartridge.</p> |

Technical Support

For technical support (assistance) please reach out to your local affiliate:
https://www.roche.com/about/business/roche_worldwide.htm.

GLOSSARY OF SYMBOLS

| Symbol | Description | Symbol | Description |
|---|--|---|---|
|  | Batch Code |  | Use by date YYYY-MM-DD |
|  | Caution |  | Serial number |
|  | Contains sufficient for <n> tests |  | Catalog number |
|  | European Union Conformity |  | Biological risks |
|  | In vitro diagnostic medical device |  | Upper limit of temperature |
|  | Consult instructions for use |  | Lower limit of temperature |
|  | Authorized representative in the European Community |  | Temperature range |
|  | Manufacturer |  | Irritant, dermal sensitizer, acute toxicity (harmful), narcotic effects, respiratory tract irritation |
|  | Cartridge Lot |  | Oxidizers |
| Rx Only | For prescription use only | | |

REFERENCES

1. John Hopkins Medicine. Upper respiratory infection (URI or common cold). Accessed: 27 Aug 2024. <https://www.hopkinsmedicine.org/health/conditions-and-diseases/upper-respiratory-infection-uri-or-common-cold>.
2. European Centre for Disease Prevention and Control. Seasonal influenza. Accessed: 27 Aug 2024. <https://www.ecdc.europa.eu/en/seasonal-influenza>.
3. World Health Organization. Influenza (seasonal). Published: 3 Oct 2023; Accessed: 27 Aug 2024. [https://www.who.int/news-room/fact-sheets/detail/influenza-\(seasonal\)](https://www.who.int/news-room/fact-sheets/detail/influenza-(seasonal)).
4. Mossad SB. Upper respiratory tract infections. Published: Aug 2013; Accessed: 27 Aug 2024. <https://teachmeanmedicine.org/cleveland-clinic-upper-respiratory-tract-infections/>.
5. University of Rochester Medical Center. Adenovirus infections in children. Accessed: 27 Aug 2024. <https://www.urmc.rochester.edu/encyclopedia/content.aspx?contenttypeid=90&contentid=p02508>.
6. Centers for Disease Control and Prevention. About adenovirus. Updated: 14 May 2024; Accessed: 27 Aug 2024. <https://www.cdc.gov/adenovirus/about/index.html>.
7. Gaunt ER, Hardie A, Claas EC, Simmonds P, Templeton KE. Epidemiology and clinical presentations of the four human coronaviruses 229E, HKU1, NL63, and OC43 detected over 3 years using a novel multiplex real-time PCR method. *J Clin Microbiol*. 2010;48:2940-7. PMID: 20554810.
8. Centers for Disease Control and Prevention. About human metapneumovirus. Updated: 11 Apr 2024; Accessed: 27 Aug 2024. <https://www.cdc.gov/human-metapneumovirus/about/index.html>.
9. Fabre V, Auwaerter PG. Human metapneumovirus. Last updated: 8 Mar 2020; Accessed: 27 Aug 2024. https://www.hopkinsguides.com/hopkins/view/Johns_Hopkins_ABX_Guide/540614/all/Metapneumovirus.
10. Anzueto A, Niederman MS. Diagnosis and treatment of rhinovirus respiratory infections. *Chest*. 2003;123:1664-72. PMID: 12740288.
11. Auwaerter P. Rhinovirus. Last updated: 18 Mar 2023; Accessed: 27 Aug 2024. https://www.hopkinsguides.com/hopkins/view/Johns_Hopkins_ABX_Guide/540476/all/Rhinovirus?q=rhinovirus&ti=0#0.
12. Fabre V, Auwaerter P. Enterovirus. Last updated: 5 Sept 2020; Accessed: 27 Aug 2024. https://www.hopkinsguides.com/hopkins/view/Johns_Hopkins_ABX_Guide/540204/all/Enterovirus?q=enterovirus&ti=0#0.
13. Henrickson KJ. Parainfluenza viruses. *Clin Microbiol Rev*. 2003;16:242-64. PMID: 12692097.
14. Schomacker H, Schaap-Nutt A, Collins PL, Schmidt AC. Pathogenesis of acute respiratory illness caused by human parainfluenza viruses. *Curr Opin Virol*. 2012;2:294-9. PMID: 22709516.
15. Mahony JB. Detection of respiratory viruses by molecular methods. *Clin Microbiol Rev*. 2008;21:716-47. PMID: 18854489.
16. Resch B, Kurath S, Manzoni P. Epidemiology of respiratory syncytial virus infection in preterm infants. *Open Microbiol J*. 2011;5:135-43. PMID: 22262986.
17. Centers for Disease Control and Prevention. *Chlamydia pneumoniae* infection: Causes and how it spreads. Updated: 30 Jan 2024; Accessed: 27 Aug 2024. <https://www.cdc.gov/cpnemomoniae/causes/index.html>.
18. Auwaerter P. *Mycoplasma pneumoniae*. Last updated: 13 Jun 2024; Accessed: 27 Aug 2024. https://www.hopkinsguides.com/hopkins/view/Johns_Hopkins_ABX_Guide/540373/all/Mycoplasma%20pneumoniae.
19. Spacek LA. Adenovirus. Updated: 5 Feb 2023; Accessed: 28 Aug 2024. https://www.hopkinsguides.com/hopkins/view/Johns_Hopkins_ABX_Guide/540009/all/Adenovirus?q=adenovirus&ti=0#0.
20. Scott MK, Chommanard C, Lu X, et al. Human adenovirus associated with severe respiratory infection, Oregon, USA, 2013-2014. *Emerg Infect Dis*. 2016;22:1044-51. PMID: 27191834.
21. Radin JM, Hawksworth AW, Blair PJ, et al. Dramatic decline of respiratory illness among US military recruits after the renewed use of adenovirus vaccines. *Clin Infect Dis*. 2014;59:962-8. PMID: 24991024.
22. Centers for Disease Control and Prevention. Human coronavirus types. Last updated: 15 Feb 2020; Accessed: 28 Aug 2024. https://archive.cdc.gov/www_cdc_gov/coronavirus/types.html.
23. European Centre for Disease Prevention and Control. Coronaviruses. Accessed: 28 Aug 2024. <https://www.ecdc.europa.eu/en/coronaviruses>.
24. Liu DX, Liang JQ, Fung TS. Human coronavirus-229E, -OC43, -NL63, and -HKU1 (*Coronaviridae*). *Encyclopedia of Virology*. 2021:428-40.
25. Hermos CR, Vargas SO, McAdam AJ. Human metapneumovirus. *Clin Lab Med*. 2010;30:131-48. PMID: 20513544.
26. Tapparel C, Junier T, Gerlach D, et al. New respiratory enterovirus and recombinant rhinoviruses among circulating picornaviruses. *Emerg Infect Dis*. 2009;15:719-26. PMID: 19402957.
27. Jacobs SE, Lamson DM, St George K, Walsh TJ. Human rhinoviruses. *Clin Microbiol Rev*. 2013;26:135-62. PMID: 23297263.

28. European Centre for Disease Prevention and Control. Enterovirus 68 detected in the USA, Canada and Europe – Second update 25 November 2014. Stockholm: ECDC; 2014.
29. Centers for Disease Control and Prevention. Types of influenza viruses. Last reviewed: 30 Mar 2023; Accessed: 28 Aug 2024. <http://www.cdc.gov/flu/about/viruses/types.htm>.
30. European Centre for Disease Prevention and Control. Factsheet about seasonal influenza. Last updated: 12 Apr 2022; Accessed: 28 Aug 2024. <https://www.ecdc.europa.eu/en/seasonal-influenza/facts/factsheet>.
31. Short KR, Richard M, Verhagen JH, et al. One health, multiple challenges: The inter-species transmission of influenza A virus. *One Health*. 2015;1:1-13. PMID: 26309905.
32. Auwaerter PG. Influenza. Last updated: 10 Sept 2023; Accessed: 28 Aug 2024. https://www.hopkinsguides.com/hopkins/view/Johns_Hopkins_ABX_Guide/540285/all/Influenza?q=influenza&ti=0#0.
33. Centers for Disease Control and Prevention. Update: influenza activity - United States, 2009-10 season. *MMWR Morb Mortal Wkly Rep*. 2010;59:901-8. PMID: 20671661.
34. Centers for Disease Control and Prevention. About human parainfluenza viruses (HPIVs). Reviewed: 5 Jun 2024; Accessed: 28 Aug 2024. <https://www.cdc.gov/parainfluenza/about/index.html>.
35. Auwaerter PG. Parainfluenza virus. Last updated: 8 Aug 2022; Accessed: 28 Aug 2024. https://www.hopkinsguides.com/hopkins/view/Johns_Hopkins_ABX_Guide/540415/all/Parainfluenza_virus?q=parainfluenza&ti=0#0.
36. Centers for Disease Control and Prevention. Clinical overview of RSV. Reviewed: 3 Jul 2024; Accessed: 29 Aug 2024. <https://www.cdc.gov/rsv/hcp/clinical-overview/index.html>.
37. Walsh EE, McConnochie KM, Long CE, Hall CB. Severity of respiratory syncytial virus infection is related to virus strain. *J Infect Dis*. 1997;175:814-20. PMID: 9086135.
38. Fabre V. *Chlamydia pneumoniae*. Last updated: 8 Feb 2023; Accessed: 28 Aug 2024. https://www.hopkinsguides.com/hopkins/view/Johns_Hopkins_ABX_Guide/540117/all/Chlamydophila%20pneumoniae.
39. Centers for Disease Control and Prevention. Risk factors for pneumonia. Accessed: 29 Aug 2024. <https://www.cdc.gov/pneumonia/risk-factors/index.html>.
40. Centers for Disease Control and Prevention. *Mycoplasma pneumoniae* infection: Causes and how it spreads. Reviewed: 27 Dec 2023; Accessed: 29 Aug 2024. <https://www.cdc.gov/mycoplasma/causes/index.html>.

DOCUMENT REVISION

| Document Revision Information | |
|-------------------------------|---|
| Doc Rev. 1.0 09/2024 | First publishing for Branchburg based on IFU PI1060-H. Updated SDS website information in Safety section. Moved Rx only text from front page to above legal manufacturer and added “For USA:”. Updated references. Please contact your local Roche Representative if you have any questions. |

TRADEMARKS

COBAS and EPLEX are trademarks of Roche.

Kimwipes™ is a trademark of Kimberly-Clark Corporation.

Flonase® is a registered trademark of GlaxoSmithKline, plc.

Nasacort® is a registered trademark of Sanofi-Aventis Pharmaceuticals.

Rhinocort® is a registered trademark of AstraZeneca AB.

Afrin® is a registered trademark of Bayer.

MicroTest™ M4®, M4RT®, M5®, and M6®, are registered trademarks of Thermo Fisher Scientific.

Neo-Synephrine® is a registered trademark of Foundation Consumer Healthcare, LLC.

Zicam® is a registered trademark of Matrixx Initiatives, Inc.

ATCC® is a registered trademark of the American Type Culture Collection.

PATENT INFORMATION

cobas® eplex respiratory pathogen panel and/or use thereof features technology claimed in one or more of the following United States patents owned or licensed by GenMark Diagnostics, Inc. or its subsidiaries, with multiple additional foreign and domestic patents pending: U.S. Patent Nos. 7,820,391, 8,486,247, 9,222,623, 9,410,663, 9,453,613, 9,498,778, 9,500,663, 9,598,722, 9,957,553, 10,001,476, 10,391,489, 10,495,656, 10,106,847, 10,273,535, 10,352,983, 10,564,211, 10,669,592, 10,670,591, 10,753,986, 10,807,090, 10,864,522, 11,156,605, D881,409, and D900,330, and other international patents and patent applications.

Unless otherwise agreed to in writing, by using a cartridge, Recipient acknowledges that Recipient has read, accepts and agrees to be bound by and comply with the General Terms and Conditions of Sale available on Roche's website which can be amended from time to time by Roche without consent. If Recipient does not accept and agree to be bound by the General Terms and Conditions of Sale, Recipient will immediately cease any further use of the cartridge.

This product is subject to a limited license to use the product in the field of human in vitro diagnostics and research reasonably related thereto. Users are prohibited from using this product for other applications, including in the field of forensics (including human identification testing).

Effective Date: September 2024

©2024 Roche Molecular Systems, Inc.

For USA: Rx only

GenMark Diagnostics, Inc. A Member of the Roche Group

5964 La Place Court, Carlsbad, CA, 92008

+1 800 428 5076

<https://diagnostics.roche.com>