cobas®

REF		Σ	SYSTEM
07007169100	07027168501	100	cobas e 402
07027168190		100	cobas e 801

English

I

	For use	in	the	USA	only	
--	---------	----	-----	-----	------	--

System information

Short name	ACN (application code number)
CPEPTID	10081

Intended use

Immunoassay for the in vitro quantitative determination of C-peptide in human serum, plasma and urine.

The assay is intended for use as an aid in the diagnosis and treatment of patients with abnormal insulin secretion.

The electrochemiluminescence immunoassay "ECLIA" is intended for use on **cobas e** immunoassay analyzers.

Summary

C-peptide is a single chain 31-amino acid (AA 33-63) polypeptide connecting the insulin A chain with the B chain in the proinsulin molecule. It has a molecular weight of approximately 3021 Da.^{1,2}

The proteolytic cleavage of the precursor proinsulin results in the two molecules insulin and C-peptide. Both are secreted in equimolar amounts and released into circulation via the portal vein. As half of the insulin, but almost none of the C-peptide is extracted in the liver, C-peptide has a longer half-life (about 35 minutes) than insulin. 5 to 10 times higher concentration of C-peptide persist in the peripheral circulation, and these levels fluctuate less than insulin. C-peptide is removed from the circulation by the kidneys and degraded, with a fraction excreted unchanged in the urine. The concentration in urine is about 10-20 fold higher than in serum.³

In the past, C-Peptide has been considered biologically inactive. However, recent studies have demonstrated that it is capable of eliciting molecular and physiological effects suggesting that C-peptide is in fact a bioactive peptide.⁴ There is evidence that C-peptide replacement, together with insulin administration, may prevent the development or retard the progression of long-term complications in type 1 diabetes.^{5,6,7,8,9,10}

Measurements of C-peptide, insulin and glucose are used as an aid in the differential diagnosis of hypoglycemia (factitious hypoglycemia and hypoglycemia caused by hyperinsulinism) to ensure an appropriate management and therapy of the patients. To quantify the endogenous insulin secretion, C-peptide is measured basally, after fasting and after stimulation and suppression tests.³ Due to high prevalence of endogenous anti-insulin antibodies C-peptide concentrations reflect the endogenous pancreatic insulin iself.^{2,11} Measurements of C-peptide may therefore be an aid in the assessment of a residual β-cell function in the early stages of type 1 diabetes mellitus and for the differential diagnosis of latent autoimmune diabetes of adults (LADA) and type 2 diabetes.^{2,11,12,13,14}

C-peptide measurements are also used to assess the success of islet transplantation and for monitoring after pancreatectomy. ^2.13,15,16 $\,$

Urinary C-peptide is measured when a continuous assessment of β -cell function is desired, to determine the Urinary C-peptide Creatinine Ratio (UCPCR), in patients with unstable glycemic control, in insulin-dependent diabetes mellitus, or when frequent blood sampling is not practical (e.g. in children).^{2,3,17}

Although testing for C-peptide is not required for the routine monitoring of diabetes, it is a valuable tool for the individual therapeutic decisions which are essential for an optimal long-term metabolic control.^{3,18}

Elevated C-peptide levels may also result from renal insufficiency and obesity. $^{\rm 3,18}$

Test principle

Sandwich principle. Total duration of assay: 18 minutes.

 1st incubation: 12 µL of sample, a biotinylated monoclonal C-peptide-specific antibody, and a monoclonal C-peptide-specific antibody labeled with a ruthenium complex^a) react to form a sandwich complex.

- 2nd incubation: After addition of streptavidin-coated microparticles, the complex becomes bound to the solid phase via interaction of biotin and streptavidin.
- The reaction mixture is aspirated into the measuring cell where the microparticles are magnetically captured onto the surface of the electrode. Unbound substances are then removed with ProCell II M. Application of a voltage to the electrode then induces chemiluminescent emission which is measured by a photomultiplier.
- Results are determined via a calibration curve which is instrumentspecifically generated by 2-point calibration and a master curve provided via the cobas link.

a) Tris(2,2'-bipyridyl)ruthenium(II)-complex (Ru(bpy)₃²⁺)

Reagents - working solutions

The cobas e pack is labeled as CPEPTID.

- M Streptavidin-coated microparticles, 1 bottle, 5.8 mL: Streptavidin-coated microparticles 0.72 mg/mL; preservative.
- R1 Anti-C-peptide-Ab~biotin, 1 bottle, 9.9 mL: Biotinylated monoclonal anti-C-peptide antibody (mouse) 1 mg/L, phosphate buffer 50 mmol/L, pH 6.0; preservative.
- R2 Anti-C-peptide-Ab~Ru(bpy)₃²⁺, 1 bottle, 9.9 mL: Monoclonal anti-C-peptide antibody (mouse) labeled with ruthenium complex 0.4 mg/L; phosphate buffer 50 mmol/L, pH 6.0; preservative.

Precautions and warnings

For in vitro diagnostic use for health care professionals. Exercise the normal precautions required for handling all laboratory reagents. Infectious or microbial waste:

Warning: handle waste as potentially biohazardous material. Dispose of waste according to accepted laboratory instructions and procedures. Environmental hazards:

Apply all relevant local disposal regulations to determine the safe disposal. Safety data sheet available for professional user on request.

For USA: Caution: Federal law restricts this device to sale by or on the order of a physician.

This kit contains components classified as follows in accordance with the Regulation (EC) No. $1272/2008\colon$



Warning

H317 Prevention:	May cause an allergic skin reaction.
P261	Avoid breathing dust/fume/gas/mist/vapours/spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves.
Response:	
P333 + P313	If skin irritation or rash occurs: Get medical advice/attention.
P362 + P364 Disposal:	Take off contaminated clothing and wash it before reuse.

cobas®

I

P501

Dispose of contents/container to an approved waste disposal plant.

Avoid foam formation in all reagents and sample types (specimens, calibrators and controls).

Product safety labeling follows EU GHS guidance.

Contact phone: 1-800-428-2336

Reagent handling

The reagents in the kit have been assembled into a ready-for-use unit that cannot be separated.

All information required for correct operation is available via the cobas link.

Storage and stability

Store at 2-8 °C.

Do not freeze.

Store the **cobas e** pack **upright** in order to ensure complete availability of the microparticles during automatic mixing prior to use.

Stability:

,	
unopened at 2-8 °C	up to the stated expiration date
on the analyzers	16 weeks

Specimen collection and preparation

Only the specimens listed below were tested and found acceptable. Serum collected using standard sampling tubes or tubes containing separating gel.

Li-heparin, K_2 -EDTA and K_3 -EDTA plasma.

Criterion: Slope 0.9-1.1 + coefficient of correlation \ge 0.95.

24-hour urine (must be prediluted 1:10 with Diluent MultiAssay before measurement).

Stability of the 24-hour urine (after collection), serum and plasma samples: 4 hours at 15-25 °C, 24 hours at 2-8 °C, 30 days at -20 °C (\pm 5 °C). Freeze only once.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

Do not use heat-inactivated samples.

Do not use samples and controls stabilized with azide.

Ensure the samples and calibrators are at 20-25 °C prior to measurement.

Due to possible evaporation effects, samples and calibrators on the analyzers should be analyzed/measured within 2 hours.

Sample stability claims were established by experimental data by the manufacturer or based on reference literature and only for the temperatures/time frames as stated in the method sheet. It is the responsibility of the individual laboratory to use all available references and/or its own studies to determine specific stability criteria for its laboratory.

Materials provided

See "Reagents - working solutions" section for reagents.

Materials required (but not provided)

- REF 03184919190, C-Peptide CalSet, for 4 x 1.0 mL
- REF 05341787160, PreciControl Multimarker, for 6 x 2.0 mL
- REF 07299010190, Diluent MultiAssay, 36 mL sample diluent
- General laboratory equipment
- cobas e analyzer

Additional materials for **cobas e** 402 and **cobas e** 801 analyzers:

- REF 06908799190, ProCell II M, 2 x 2 L system solution
- REF 04880293190, CleanCell M, 2 x 2 L measuring cell cleaning solution

- REF 07485409001, Reservoir Cup, 8 cups to supply ProCell II M and CleanCell M
- REF 06908853190, PreClean II M, 2 x 2 L wash solution
- REF 05694302001, Assay Tip/Assay Cup tray, 6 magazines x 6 magazine stacks x 105 assay tips and 105 assay cups, 3 wasteliners
- REF 07485425001, Liquid Flow Cleaning Cup, 2 adaptor cups to supply ISE Cleaning Solution/Elecsys SysClean for Liquid Flow Cleaning Detection Unit
- REF 07485433001, PreWash Liquid Flow Cleaning Cup, 1 adaptor cup to supply ISE Cleaning Solution/Elecsys SysClean for Liquid Flow Cleaning PreWash Unit
- REF 11298500316, ISE Cleaning Solution/Elecsys SysClean, 5 x 100 mL system cleaning solution
- REF 11298500160, ISE Cleaning Solution/Elecsys SysClean, 5 x 100 mL system cleaning solution (for USA)

Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

Resuspension of the microparticles takes place automatically prior to use.

Place the cooled (stored at 2-8 °C) **cobas e** pack on the reagent manager. Avoid foam formation. The system automatically regulates the temperature of the reagents and the opening/closing of the **cobas e** pack.

Calibration

Traceability: This method has been standardized against the WHO International Reference Reagent for C-peptide of human insulin for immunoassay, IRR, code 84/510, established 1986, from the National Institute for Biological Standards and Control (NIBSC).¹⁹

The predefined master curve is adapted to the analyzer using the relevant CalSet.

Calibration frequency: Calibration must be performed once per reagent lot using fresh reagent (i.e. not more than 24 hours since the **cobas e** pack was registered on the analyzer).

Calibration interval may be extended based on acceptable verification of calibration by the laboratory.

Renewed calibration is recommended as follows:

- after 12 weeks when using the same reagent lot
- after 28 days when using the same cobas e pack on the analyzer
- as required: e.g. quality control findings outside the defined limits

Quality control

For quality control, use PreciControl Multimarker.

In addition, other suitable control material can be used.

Controls for the various concentration ranges should be run individually at least once every 24 hours when the test is in use, once per **cobas e** pack, and following each calibration.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

If necessary, repeat the measurement of the samples concerned.

Follow the applicable government regulations and local guidelines for quality control.

Calculation

The analyzer automatically calculates the analyte concentration of each sample in nmol/L, ng/mL or pmol/L (selectable).

Conversion factors:

ng/mL (µg/L) x 0.33333 = nmol/L ng/mL x 333.33 = pmol/L nmol/L x 3.0 = ng/mL pmol/L x 0.003 = ng/mL

cobas®

I

I

I

Limitations - interference

The effect of the following endogenous substances and pharmaceutical compounds on assay performance was tested. Interferences were tested up to the listed concentrations and no impact on results was observed.

Endogenous substances

0	
Compound	Concentration tested
Bilirubin	\leq 855 µmol/L or \leq 50 mg/dL
Hemoglobin	\leq 0.186 mmol/L or \leq 300 mg/dL
Intralipid	≤ 2000 mg/dL
Biotin	\leq 246 nmol/L or \leq 60 ng/mL
Rheumatoid factors	≤ 1200 IU/mL

Criterion: For concentrations \leq 0.5 ng/mL the deviation is \leq 0.2 ng/mL of initial value. For concentrations > 0.5 ng/mL the deviation is \leq 10 % of initial value.

Samples should not be taken from patients receiving therapy with high biotin doses (i.e. > 5 mg/day) until at least 8 hours following the last biotin administration.

There is no high-dose hook effect at C-peptide concentrations up to 60.0 nmol/L (180 ng/mL).

Pharmaceutical substances

In vitro tests were performed on 16 commonly used pharmaceuticals in serum and 12 commonly used pharmaceuticals in urine. No interference with the assay was found.

In rare cases, interference due to extremely high titers of antibodies to analyte-specific antibodies, streptavidin or ruthenium can occur. These effects are minimized by suitable test design.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

Limits and ranges Measuring range

Serum and plasma: 0.007-13.3 nmol/L or 0.02-40 ng/mL (defined by the Limit of Detection and the maximum of the master curve). Values below the Limit of Detection are reported as < 0.007 nmol/L (< 0.02 ng/mL). Values above the measuring range are reported as > 13.3 nmol/L (> 40 ng/mL) (or up to 133 nmol/L or 400 ng/mL for 10-fold diluted samples).

Urine: 0.067-133 nmol/L or 0.2-400 ng/mL (defined by the 10 x Limit of Detection for serum/plasma and the 10 x maximum of the master curve for serum/plasma, thus taking into account the 1:10 predilution of urine samples with Diluent MultiAssay). Values below the Limit of Detection are reported as < 0.067 nmol/L (< 0.2 ng/mL). Values above the measuring range are reported as > 133 nmol/L (> 400 ng/mL) or retested in a higher dilution of the sample.

Lower limits of measurement

Limit of Blank, Limit of Detection and Limit of Quantitation

Serum and plasma:

Limit of Blank = 0.003 nmol/L (0.01 ng/mL)

Limit of Detection = 0.007 nmol/L (0.02 ng/mL)

Limit of Quantitation = 0.050 nmol/L (0.15 ng/mL)

Urine:

Please refer to the values for serum/plasma, taking into account the mandatory 1:10 predilution of urine samples.

The Limit of Blank, Limit of Detection and Limit of Quantitation were determined in accordance with the CLSI (Clinical and Laboratory Standards Institute) EP17-A2 requirements.

The Limit of Blank is the 95th percentile value from n \ge 60 measurements of analyte-free samples over several independent series. The Limit of Blank corresponds to the concentration below which analyte-free samples are found with a probability of 95 %.

The Limit of Detection is determined based on the Limit of Blank and the standard deviation of low concentration samples. The Limit of Detection corresponds to the lowest analyte concentration which can be detected (value above the Limit of Blank with a probability of 95 %).

The Limit of Quantitation is the lowest analyte concentration that can be reproducibly measured with an intermediate precision CV of \leq 20 %.

Dilution

Serum and plasma: Although the necessity for dilutions is unlikely due to the high measuring range, samples with C-peptide concentrations above the measuring range can be diluted with Diluent MultiAssay. The recommended dilution is 1:10 (either automatically by the analyzers or manually). The concentration of the diluted sample must be \geq 1.3 nmol/L (\geq 4 ng/mL).

After manual dilution, multiply the result by the dilution factor.

After dilution by the analyzers, the software automatically takes the dilution into account when calculating the sample concentration.

Urine: All urine samples must be prediluted 1:10 with Diluent MultiAssay before measurement (either automatically by the analyzers or manually).

After dilution by the analyzers, the software automatically takes the dilution into account when calculating the sample concentration.

Urine samples with C-peptide concentrations above the measuring range can be retested using a 1:20 or higher dilution with Diluent MultiAssay (either automatically by the analyzer or manually). The concentration of the diluted sample must be \geq 1.3 nmol/L (\geq 4 ng/mL).

After manual dilution, multiply the result by the dilution factor.

Expected values

Studies with the Elecsys C-Peptide assay were performed using serum samples from apparently healthy fasting males and females, and 24 h urine samples from apparently healthy individuals.

The following results were obtained:

	Ν	Median	5 th -95 th percentile	Unit
C-peptide in	96	1.96	1.1-4.4	ng/mL
serum/plasma	30	0.65	0.37-1.47	nmol/L
C-peptide in	70	54.8	17.2-181	µg/24 h
24-hour urine	15	18.3	5.74-60.3	nmol/24 h

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

Specific performance data

Representative performance data on the analyzers are given below. Results obtained in individual laboratories may differ.

Precision

Serum and plasma:

Precision was determined using Elecsys reagents, pooled human sera and controls in a protocol (EP05-A3) of the CLSI (Clinical and Laboratory Standards Institute): 2 runs per day in duplicate each for 21 days (n = 84). The following results were obtained:

cobas e 402 and cobas e 801 analyzers						
		Repeata	bility	Intermed precisi	diate on	
Sample	Mean nmol/L	SD nmol/L	CV %	SD nmol/L	CV %	
Human serum 1	0.041	0.001	2.9	0.001	3.5	
Human serum 2	0.337	0.003	0.9	0.008	2.3	
Human serum 3	1.34	0.028	2.1	0.044	3.3	
Human serum 4	6.37	0.128	2.0	0.211	3.3	
Human serum 5	12.1	0.343	2.8	0.440	3.6	
PC ^{b)} Multimarker 1	0.670	0.009	1.3	0.017	2.6	
PC Multimarker 2	3.40	0.062	1.8	0.109	3.2	

b) PC = PreciContro

cobas e 402 and cobas e 801 analyzers						
		Repeatability		Intermediate precision		
Sample	Mean ng/mL	SD ng/mL	CV %	SD ng/mL	CV %	
Human serum 1	0.124	0.004	2.9	0.004	3.5	
Human serum 2	1.01	0.009	0.9	0.023	2.3	
Human serum 3	4.01	0.084	2.1	0.133	3.3	
Human serum 4	19.1	0.385	2.0	0.632	3.3	
Human serum 5	36.4	1.03	2.8	1.32	3.6	
PC Multimarker 1	2.01	0.027	1.3	0.052	2.6	
PC Multimarker 2	10.2	0.186	1.8	0.327	3.2	

Urine:

1

Precision was determined using Elecsys reagents and human urine samples in a protocol (EP05-A3) of the CLSI (Clinical and Laboratory Standards Institute): 2 runs per day in duplicate each for 21 days (n = 84). The following results were obtained:

cobas e 402 and cobas e 801 analyzers					
	Repeatability			Intermed precisi	diate on
Sample	Mean	SD	CV	SD	CV
	nmol/L	nmol/L	%	nmol/L	%
Urine 1	0.357	0.020	5.5	0.022	6.2
Urine 2	3.37	0.128	3.8	0.160	4.8
Urine 3	12.8	0.231	1.8	0.277	2.2
Urine 4	63.7	0.917	1.4	2.46	3.9
Urine 5	130	2.55	2.0	3.14	2.4

cobas e 402 and cobas e 801 analyzers						
		Repeatal	bility	Intermed precisi	diate on	
Sample	Mean ng/mL	SD ng/mL	CV %	SD ng/mL	CV %	
Urine 1	1.07	0.059	5.5	0.067	6.2	
Urine 2	10.1	0.384	3.8	0.480	4.8	
Urine 3	38.5	0.692	1.8	0.831	2.2	
Urine 4	191	2.75	1.4	7.38	3.9	
Urine 5	391	7.64	2.0	9.43	2.4	

Method comparison

Serum

A comparison of the Elecsys C-Peptide assay on the **cobas e** 801 analyzer (y) with the **cobas e** 601 analyzer (x) using clinical samples gave the following correlations (ng/mL):

Number of samples measured: 169

Passing/Bablok ²⁰	Linear regression
y = 1.00x - 0.002	y = 1.00x + 0.014
т = 0.994	r = 1.00

The sample concentrations were between 0.091 and 39.0 ng/mL. Urine

A comparison of the Elecsys C-Peptide assay on the **cobas e** 801 analyzer (y) with the **cobas e** 601 analyzer (x) using clinical samples gave the following correlations (ng/mL):

Number of samples measured: 143

Passing/Bablok²⁰

Linear regression

y = 1.02x - 0.091	y = 1.02x
-------------------	-----------

т = 0.992

The sample concentrations were between 0.588 and 381 ng/mL.

Analytical specificity

For the monoclonal antibodies used, the following cross-reactivities were found:

r = 1.00

Substance	Concentration tested µg/mL	Cross-reactivity %
Proinsulin, human ^{c)}	0.10	28.6
Insulin, human ^{d)}	8.66	n. d. ^{e)}
Insulin, porcine ^{f)}	7.50	n. d.
Insulin, bovine ^{g)}	7.69	n. d.
Somatomedin ^{h)} (Insulin-like growth factor 1 - IGF-I)	1.0	n. d.
Human Growth Hormone ⁱ⁾	10.0	n. d.
Glucagon ⁱ⁾	10.0	n. d.

c) WHO preparation 09/296

d) WHO preparation 66/304

e) n. d. = not detectable

f) WHO preparation 83/515

g) WHO preparation 83/511h) NIBSC code 02/254

i) NIBSC code 98/574

i) NIBSC code 69/194

The Elecsys C-Peptide assay uses two monoclonal antibodies specifically directed against human C-peptide. The antibodies show cross-reactivity with the C-chain of human proinsulin and presumably with partially processed proinsulins (split products). The concentrations of proinsulin and split products of fasting healthy subjects are 100 times lower than the C-peptide concentrations and therefore the cross-reactivity is of no clinical significance. In patients with insulinoma, the proinsulin concentrations are reported as up to 60-fold higher than those from fasting healthy subjects.^{21,22}

References

- 1 Clark PM. Assays for insulin, proinsulin(s) and C-peptide. Ann Clin Biochem 1999;36(5):541-564.
- 2 Sacks DB. Chapter 24: Carbohydrates. In: Burtis CA, Ashwood ER (eds). Tietz Textbook of Clinical Chemistry, WB Saunders, Philadelphia, 3rd edition;1999:750-808.
- 3 Jones AG, Hattersley AT. The clinical utility of C-peptide measurement in the care of patients with diabetes. Diabetic Medicine 2013;30:803-817.
- 4 Yosten GLC, Grant KR. The Physiology of Proinsulin C-Peptide: Unanswered Questions and a Proposed Model. Physiology 2015;30(4):327-332.
- 5 Johansson J, Ekberg K, Shafqat J, et al. Molecular effects of proinsulin C-peptide. Biochem Biophys Res Commun 2002;295:1035-1040.
- 6 Kobayashi T, Maruyama T, Shimada A, et al. Insulin Intervention to Preserve β Cells in Slowly Progressive Insulin-Dependent (Type 1) Diabetes Mellitus. Ann N Y Acad Sci 2002;958(4):117-130.
- 7 Forst T, Rave K, Pfuetzner A, et al. Effect of C-Peptide on Glucose Metabolism in Patients With Type 1 Diabetes. Diabetes Care 2002;25(6):1096-1097.
- 8 Shapiro AMJ. Islet Transplants and Impact on Secondary Diabetic Complications: Does C-Peptide Protect the Kidney? J Am Nephrol 2003;14:2214-2216.
- 9 Sima AAF. C-peptide and diabetic neuropathy. Expert Opin Investig Drugs 2003;12(9)1471-1488.
- 10 Wahren J, Jörnvall H. C-peptide makes a comeback. Diabetes Metab Res Rev 2003;19:345-347.



- 0.289

11 Törn C. C-peptide and Autoimmune Markers in Diabetes. Clin Lab 2003;49:1-10.

- 12 Pourmotabbed G, Kitabchi AE. Hypoglycemia. Obst Gynecol Clin North Am 2001;28(2):383-400.
- 13 Batstra MR, Aanstoot H-J, Herbrink P. Prediction and Diagnosis of Type 1 Diabetes Using β-cell Autoantibodies. Clin Lab 2001;47;497-507.
- 14 Meier CH, Ladewig A, Keller U, et al. Clinical Value of the C-Peptide Measurement. Schweiz Rundsch Med Prax 1997;86(34):1289-1295.
- 15 Gottsäter A, Landin-Olsson M, Fernlund P, et al. β-Cell Function in Relation to Islet Cell Antibodies During the First 3 Yr After Clinical Diagnosis of Diabetes in Type II Diabetic Patients. Diabetes Care 1993;16(6):902-910.
- 16 VanBuecken DE, Greenbaum CJ. Residual C-peptide in type 1 diabetes: what do we really know? Pediatric Diabetes 2014;15(2):84-90.
- 17 Cha T, Tahara Y, Ikegami H, et al. Urinary C-peptide as an index of unstable glycemic control in insulin-dependent diabetes mellitus (IDDM). Diabetes Res Clin Pract 1991;13:181-188.
- 18 Sacks DB, Bruns DE, Goldstein DE, et al. Guidelines and Recommendations for Laboratory Analysis in the Diagnosis and Management of Diabetes Mellitus. Clin Chem 2002;48(3):436-472.
- 19 Bristow AF, Gaines-Das RE. WHO international reference reagents for human proinsulin and human insulin C-peptide. J Biol Stand 1988;16:179-186.
- 20 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.
- 21 Houssa P, Dinesen B, Deberg M, et al. First direct assay for intact human proinsulin. Clin Chem 1998;44(7):1514-1519.
- 22 Zilkens TM, Eberle AM, Schmidt-Gayk H. Immunoluminometric assay (ILMA) for intact human proinsulin and its conversion intermediates. Clin Chem Acta 1996;247:23-37.

For further information, please refer to the appropriate operator's manual for the analyzer concerned, the respective application sheets and the Method Sheets of all necessary components (if available in your country).

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard (for USA: see dialog.roche.com for definition of symbols used):

CONTENT	Contents of kit
SYSTEM	Analyzers/Instruments on which reagents can be used
REAGENT	Reagent
CALIBRATOR	Calibrator
\rightarrow	Volume for reconstitution
GTIN	Global Trade Item Number

FOR US CUSTOMERS ONLY: LIMITED WARRANTY

Roche Diagnostics warrants that this product will meet the specifications stated in the labeling when used in accordance with such labeling and will be free from defects in material and workmanship until the expiration date printed on the label. THIS LIMITED WARRANTY IS IN LIEU OF ANY OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. IN NO EVENT SHALL ROCHE DIAGNOSTICS BE LIABLE FOR INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES. COBAS, COBAS E, ELECSYS and PRECICONTROL are trademarks of Roche. INTRALIPID is a trademark of Fresenius Kabi AB.

All other product names and trademarks are the property of their respective owners. Additions, deletions or changes are indicated by a change bar in the margin. © 2022, Roche Diagnostics



Distribution in USA by: Boche Diagnostics, Indianapolis, II

Roche Diagnostics, Indianapolis, IN US Customer Technical Support 1-800-428-2336

