## LightMix<sup>®</sup> Modular RSV-A/B

TIB Molbiol Cat. No. 61-0110-96

610

Roche SAP n° 08 492 018 001

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Total weight (incl. packaging) : 14 g

## 2. Intended Use

**This product is intended to be used to diagnose an infection by detection of the viral genome in a Total-NA nucleic acid extract obtained from respiratory tract specimen from patients with significant respiratory symptoms. This assay has been developed for use in multiplex PCR and can be run combined with a control reaction or combined with other kits as specified in section 5.**

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### 3. Disclaimer

**Please note, a negative test result does not rule out infection, as factors like poor sample quality, low template concentration, or genetic variability may lead to false negatives. Therefore, test results must be interpreted alongside clinical signs, patient history, and other diagnostics. If infection is still suspected or when results are inconsistent with clinical indications, additional testing should be considered.**

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## 4. General Precautions and Warnings

**EN:** This product must be used by qualified personnel. Use the instructions as specified for the respective lot. Do not mix reagents from different lots. Do not use expired reagents. Dispose unused reagents.

The laboratory work-flow must conform to standard practices. Handle reagents at a different place than positive materials. Wear disposable protective gloves, laboratory coats and adequate eye protection. Sample materials and waste must be considered as potentially infectious. Reagent preparation, extraction and PCR should be performed in separated areas. Clean and disinfect work surfaces.

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## 5. Application procedures for this product

FR:	Procédures d'application pour ce produit	PT:	Procedimentos de aplicação para este produto
ES:	Procedimientos de aplicación de este producto	DE:	Anwendungsverfahren für dieses Produkt
IT:	Procedura di applicazione di questo prodotto		

### 5.1. LightCycler® PRO LCAP, LightCycler® 480 Detection Format and Settings:

#### LightCycler® PRO Instrument:

For a matching LightCycler® Analysis Package (LCAP) file for *RSV-A/B* (see Section 6, Table 3) please visit [navifyportal.roche.com](http://navifyportal.roche.com) to download. Please check for the latest version of the LCAP. Save the LCAP file in the **assay** folder of the SFTP or USB device.

**Import and install** the downloaded LCAP onto the LightCycler® PRO Instrument and **Activate** it.

#### LightCycler® 480 Instrument and cobas z 480 Analyzer:

Set Quant Factor 10, Max. Integration time 2 sec:

LightCycler® 480 Instrument:	558-610
LightCycler® 480 II Instrument:	533-610
cobas z 480 Analyzer (open channel):	540-610

Table 1: Cycling condition programming.

Program Step:	RT Step	Denaturation	Cycling			Cooling
<b>Parameter</b>						
Analysis Mode	<b>None</b>	None	Quantification mode			None
Cycles	1	1	45			1
Target °C	55	95	95	60	72	40
Hold hh:mm:ss	00:05:00	00:05:00	00:00:05	00:00:15	00:00:15	00:00:30
Ramp Rate °C/s <b>96</b>	4.4	4.4	4.4	2.2	4.4	1.5
Ramp Rate °C/s <b>384</b>	4.6	4.6	4.6	2.4	4.6	2.0
Acquisition Mode	None	None	None	<b>Single</b>	None	None

### 5.2. Preparation of Parameter-Specific Reagents (PSR, **purple cap**):

Check for the colored pellet, then **add 50 µl** PCR-grade water, mix (vortex or by pipetting up and down 10 times) and spin down. **Use 0.5 µl** PSR per 20 µl PCR reaction.

Store dissolved PSR frozen (-20 °C, until expiry) or cooled (2 to 8 °C, up to 14 days).

### 5.3. Preparation of the Positive Control (**black cap**)

Check for the colored pellet, then **add 160 µl** RNase/DNase-free 10 mM Tris buffer pH 8 - 8.5.

Mix (vortex or by pipetting up and down 10 times) and spin down vial prior to opening. **Use 5 µl** Positive Control per 20 µl PCR reaction.

Store dissolved Positive Control frozen (-20 °C, until expiry).

## 5.4. Preparation of PCR Reaction Mix

In a cooled reaction tube, prepare the reaction mix for one reaction (left) or 100 reactions (right):

Table 2: Preparation of the Reaction Mix.

Reaction Mix to be used with LightMix® Modular RSV-A/B:		
LightCycler® Multiplex RNA Virus Master		Roche Cat. No. 06 754 155 001 Roche Cat. No. 07 083 173 001
One reaction	Component	100 reactions
10.4 µl	Water, PCR-grade	1040 µl
0.5 µl	PSR (purple cap)	50 µl
--	Control Reaction and additional assays (Multiplex PCR) – substitute water	--
4.0 µl	RT-qPCR Reaction Mix, 5x conc. (LightCycler® Multiplex RNA Virus Master)	400 µl
0.1 µl	RT Enzyme Solution, 200x conc. (LightCycler® Multiplex RNA Virus Master)	10 µl
<b>15.0 µl</b>	<b>Volume of Reaction Mix</b>	<b>1500 µl</b>

Mix gently, spin down and transfer 15 µl per well.

Add 5 µl of sample or Positive Control to each well for a final reaction volume of 20 µl.

Seal plate and spin down. **Start run.**

## 6. Multiplex PCR Respiratory Virus

**EN: This assay can be combined with other assays up to 6plex as depicted below. If used in multiplex setup, detection format changes and a color compensation must be applied. Please see IFU for differing detection format and settings setup:**

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LightMix® Universal Color Compensation Hexaplex Plus (TIB Molbiol Cat. No. 40-0320-12; Roche Cat. No. 06 296 971 001) is not necessary for the LightCycler® PRO, but mandatory for all LightCycler® 480 instruments.

Table 3: Options for combining LightMix® Modular RSV-A/B in Multiplex PCR and Instrument Compatibility.

### Multiplex PCR and Instrument Compatibility

500	530	580	610	640	660
			RSV-A/B		
Parainfluenza	Influenza A	Influenza B	RSV-A/B	Picornavirus HRV+EV	EAV or Roche RPC*
Influenza A H1	Influenza A	Influenza B	RSV-A/B	Picornavirus HRV+EV	
Influenza A H1	Influenza A	Influenza B	RSV-A/B	Parechovirus	
Parainfluenza	Influenza A	Influenza B	RSV-A/B	Parechovirus	
Sarbecovirus	Influenza A	Influenza B	RSV-A/B	Picornavirus HRV+EV	EAV

LC PRO	LC 480 II	LC 480 I	z 480	LC96
#	X	X	X	X
1	X	X	§	
2	X	X	§	
3	X	X	§	
4	X	X	§	
5	X	X	§	

# any of the LightCycler® PRO LCAPs mentioned below

1. 0701\_RespiV1-96\_mPCR1 Run profile: ModularDx7plex\_96\_R6G
2. 0702\_RespiV1-96\_mPCR2 Run profile: ModularDx7plex\_96\_R6G
3. 0803\_RespiV2-96\_mPCR3 Run profile: LMM\_7plex\_96
4. 0804\_RespiV2-96\_mPCR4 Run profile: LMM\_7plex\_96
5. 0601\_SARS-CoV2-96\_mPCR1 Run profile: ModularDx7plex\_96\_R6G

\* LightMix® Modular *Influenza B* (TIB Product No. 58-0102-96) cannot be combined with Roche RNA Process Control Kit (RPC, Roche Cat. No. 07 099 592 001).

§ LightMix® Modular *RSV-A/B* may be used with the cobas® z 480 analyzer in a multiplex setup, up to a 5plex, without a product in the channel 500.

## 7. Version History

**EN:** Red notes mark changes in procedures, blue change of components

FR: Notes rouges: changements dans les procédures, bleues: des composants  
 ES: Notas rojas marcan cambios en los procedimientos, en azul de componentes  
 IT: Le note in rosso indicano cambiamenti nelle procedure, in blu di componenti  
 PT: Notas vermelhas marcam mudanças nos procedimentos, a azul nos componentes  
 DE: Geänderte Prozeduren sind rot, bei geänderten Komponenten blau markiert.

V190123	Editorial changes. CoA notes clarified.	2019-03-01
V231220	Addition of the LightCycler® PRO and editorial changes	2023-12-20
V250815	3. Disclaimer added, 5.1 LightCycler® PRO information updated, 6. Table 3 updated and editorial changes	2025-08-15