

ONLINE DAT Benzodiazepines II**Order information**

REF	CONTENT	Analyzer(s) on which cobas c pack(s) can be used
08056927190	ONLINE DAT Benzodiazepines II (850 tests)	System-ID 2028 001 cobas c 303, cobas c 503

Materials required (but not provided):

03304671190	Preciset DAT Plus I CAL 1-6 (1 x 5 mL)	Codes 20431-20436
03304680190	Preciset DAT Plus II CAL 1-6 (1 x 5 mL)	Codes 20437-20442
03304698190	C.f.a.s. DAT Qualitative Plus (6 x 5 mL)	Code 20698
04590856190	C.f.a.s. DAT Qualitative Plus Clinical (3 x 5 mL)	Code 20699
03312950190	Control Set DAT I (for 300 ng/mL assay) PreciPos DAT Set I (2 x 10 mL) PreciNeg DAT Set I (2 x 10 mL)	
03312968190	Control Set DAT II (for 100 ng/mL assay) PreciPos DAT Set II (2 x 10 mL) PreciNeg DAT Set II (2 x 10 mL)	
04500873190	Control Set DAT Clinical (for 100 ng/mL assay) PreciPos DAT Clinical (2 x 10 mL) PreciNeg DAT Clinical (2 x 10 mL)	
03312976190	Control Set DAT III (for 200 ng/mL assay) PreciPos DAT Set III (2 x 10 mL) PreciNeg DAT Set III (2 x 10 mL)	

English**System information****BZ1Q2:** ACN 20280 (Urine): for qualitative assay, 100 ng/mL**BZ2Q2:** ACN 20281 (Urine): for qualitative assay, 200 ng/mL**BZ3Q2:** ACN 20282 (Urine): for qualitative assay, 300 ng/mL**BZ1S2:** ACN 20284 (Urine): for semiquantitative assay, 100 ng/mL**BZ2S2:** ACN 20285 (Urine): for semiquantitative assay, 200 ng/mL**BZ3S2:** ACN 20286 (Urine): for semiquantitative assay, 300 ng/mL**BZQ1C:** ACN 20283 (Urine): for qualitative assay, 100 ng/mL; using C.f.a.s. DAT Qualitative Plus Clinical**BZ3-QP:** ACN 20288 (Urine): for qualitative assay, 300 ng/mL; using C.f.a.s. DAT Qualitative Plus**Intended use**

Benzodiazepines II (BNZ2) is an in vitro diagnostic test for the qualitative and semiquantitative detection of benzodiazepines in human urine on **cobas c** systems at cutoff concentrations of 100 ng/mL, 200 ng/mL, and 300 ng/mL.

Semiquantitative test results may be obtained that permit laboratories to assess assay performance as part of a quality control program.

Benzodiazepines II provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC-MS) or Liquid Chromatography coupled with Tandem Mass Spectrometry (LC-MS/MS) is the preferred confirmatory method.^{1,2} Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

Summary

The benzodiazepines constitute a class of versatile and widely prescribed central nervous system (CNS) depressant drugs with medically useful anxiolytic, sedative, hypnotic, muscle relaxant, and anticonvulsant activities.^{1,2,3,4,5} The absorption rates, distribution, metabolism, and elimination rates differ significantly among the benzodiazepine derivatives. The quantitative differences in their potencies, pharmacodynamic spectra, and pharmacokinetic properties have led to various therapeutic applications. Clinical distinction of short-acting versus long-acting benzodiazepines have been observed in their efficacy, side effect, withdrawal, and dependence potential.^{3,6,7} The extensive and efficacious therapeutic use of the benzodiazepines over the last several decades has inadvertently led to their misuse. Benzodiazepine overdoses are frequently associated with co-administration of drugs of other classes.^{8,9} Acute or chronic alcohol ingestion and benzodiazepines co-administered may lead to various significant toxicological interactions. The net effect may be

influenced by internal, external, and pharmacokinetic factors. Abuse patterns may involve relatively low benzodiazepine doses, as well as high-dose overuse; therefore, urinary drug/metabolite detection requires the proper selection of a cutoff that suits the requirements of the drug testing program.

Following ingestion, the benzodiazepines of the 1,4-substituted class (including the triazolobenzodiazepine derivatives) are absorbed, metabolized, and excreted in the urine at different rates as a variety of structurally related metabolites. Metabolite diversity reflects the different physiochemical properties and metabolic pathways of the individual drugs. Overall metabolic similarities include removal of substituents from the β ring of the 1,4-substituted benzodiazepines, α -hydroxylation of the triazolobenzodiazepines, demethylation, hydroxylation of the three-position carbon of the β ring, and conjugation of hydroxylated metabolites followed by urinary excretion predominantly as glucuronides.^{1,2,3,4,5} The enzymatic hydrolysis of glucuronidated benzodiazepines can increase their cross-reactivities to benzodiazepine immunoassays.^{10,11,12,13,14}

Test principle

The assay is based on the kinetic interaction of microparticles in a solution (KIMS)^{11,15} as measured by changes in light transmission. In the absence of sample drug, free antibody binds to drug-microparticle conjugates causing the formation of particle aggregates that are photometrically detected by turbidity measurements. As the aggregation reaction proceeds in the absence of sample drug, the absorbance increases.

When a urine sample contains the drug in question, this drug competes with the particle-bound drug derivative for free antibody. Antibody bound to sample drug is no longer available to promote particle aggregation, and subsequent particle lattice formation is inhibited. The presence of sample drug diminishes the increasing absorbance in proportion to the concentration of drug in the sample. Sample drug content is determined relative to the value obtained for a known cutoff concentration of drug.

The presence of β -glucuronidase enzyme enhances the Benzodiazepines II assay cross-reactivity to some of the glucuronidated metabolites. Enzymatic cleavage makes the benzodiazepine part of glucuronides more accessible for the antibody.

Reagents - working solutions

R1 Benzodiazepines antibody (sheep polyclonal); buffer; β -glucuronidase enzyme; bovine serum albumin (BSA); 0.09 % sodium azide

R2 Conjugated benzodiazepine derivative microparticles; buffer; 0.09 % sodium azide

R1 is in position B and R2 is in position C.

Precautions and warnings

For in vitro diagnostic use for health care professionals. Exercise the normal precautions required for handling all laboratory reagents.

Infectious or microbial waste:

Warning: handle waste as potentially biohazardous material. Dispose of waste according to accepted laboratory instructions and procedures.

Environmental hazards:

Apply all relevant local disposal regulations to determine the safe disposal.

Safety data sheet available for professional user on request.

Reagent handling

Ready for use

Carefully invert reagent container several times prior to use to ensure that the reagent components are mixed.

Storage and stability

Shelf life at 2-8 °C: See expiration date on **cobas c** pack label

On-board in use and refrigerated on the analyzer: 26 weeks

Do not freeze.

Specimen collection and preparation

Only the specimens listed below were tested and found acceptable.

Urine: Collect urine samples in clean glass or plastic containers. Fresh urine specimens do not require any special handling or pretreatment, but an effort should be made to keep pipetted samples free of gross debris.

Samples should be within the normal physiological pH range of 5-8. No additives or preservatives are required. It is recommended that urine specimens be stored at 2-8 °C and tested within 5 days of collection.¹⁶

For prolonged storage, freezing of samples is recommended.

Centrifuge highly turbid specimens before testing.

See the limitations and interferences section for details about possible sample interferences.

Adulteration or dilution of the sample can cause erroneous results. If adulteration is suspected, another sample should be collected. Specimen validity testing is required for specimens collected under the *Mandatory Guidelines for Federal Workplace Drug Testing Programs*.¹⁷

CAUTION: Specimen dilutions should only be used to interpret results of Calc.? and Samp.? alarms, or when estimating concentration in preparation for GC-MS or LC-MS/MS. Dilution results are not intended for patient values. Dilution procedures, when used, should be validated.

Materials provided

See "Reagents – working solutions" section for reagents.

Materials required (but not provided)

See "Order information" section

General laboratory equipment

Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

The performance of applications not validated by Roche is not warranted and must be defined by the user.

Application for urine

	Semiquantitative	Qualitative
Reporting time	10 min	10 min
Wavelength (sub/main)	– /546 nm	– /546 nm
Reagent pipetting		Diluent (H ₂ O)
R1	63 µL	–
R2	28 µL	–
Sample volumes	Sample	Sample dilution
100 and 200 ng/mL cutoffs	Sample	Diluent (NaCl)

Normal	3.2 µL	–	–
Decreased	3.2 µL	–	–
Increased	3.2 µL	–	–

300 ng/mL cutoff

Normal	1.4 µL	–	–
Decreased	1.4 µL	–	–
Increased	1.4 µL	–	–

For further information about the assay test definitions refer to the application parameters setting screen of the corresponding analyzer and assay.

Calibration

Calibrators	<i>Semiquantitative applications</i> 100 and 200 ng/mL cutoff assays S1-6: Preciset DAT Plus II calibrators, CAL 1-6 0, 50, 100, 200, 400, 1000 ng/mL <i>300 ng/mL cutoff assay</i> S1-6: Preciset DAT Plus I calibrators, CAL 1-6 0, 150, 300, 600, 1000, 3000 ng/mL <i>Qualitative applications</i> <i>100 ng/mL cutoff assay</i> S1: Preciset DAT Plus II calibrator - CAL 3 (<i>Test BZ1Q2</i>), 100 ng/mL, S1: C.f.a.s. DAT Qualitative Plus Clinical (<i>Test BZQ1C</i>), 100 ng/mL <i>200 ng/mL cutoff assay</i> S1: Preciset DAT Plus II calibrator - CAL 4, 200 ng/mL <i>300 ng/mL cutoff assay</i> S1: Preciset DAT Plus I calibrator - CAL 3 (<i>Test BZ3Q2</i>), 300 ng/mL, S1: C.f.a.s. DAT Qualitative Plus (<i>Test BZ3-QP</i>), 300 ng/mL The drug concentrations of the calibrators have been verified by GC-MS.
Calibration K factor	For the qualitative application a K factor of -1000 is predefined in the application settings.
Calibration mode	<i>Semiquantitative applications</i> Non-linear <i>Qualitative applications</i> Linear
Calibration frequency	Full calibration - after reagent lot change - every 13 weeks on-board - as required following quality control procedures

For the cutoff calibrator a value of "0" is encoded in the e-barcode in order to ensure flagging of positive samples with >Test and negative absorbance values for negative samples.

Calibration interval may be extended based on acceptable verification of calibration by the laboratory.

Traceability: This method has been standardized against a primary reference method (GC-MS).

Quality control

For quality control, use control materials as listed in the "Order information" section. In addition, other suitable control material can be used.

Drug concentrations of Control Set DAT I, II, III and Clinical have been verified by GC-MS.

ONLINE DAT Benzodiazepines II

The control intervals and limits should be adapted to each laboratory's individual requirements. It is recommended to perform quality control always after lot calibration and subsequently at least every 26 weeks.

Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

Results

For the qualitative assay, the cutoff calibrator is used as a reference in distinguishing between preliminary positive and negative samples. Samples producing a positive or "0" absorbance value are considered preliminary positive. Preliminary positive samples are flagged with >Test. Samples producing a negative absorbance value are considered negative. Negative samples are preceded by a minus sign.

For the semiquantitative applications **cobas c** systems automatically calculate the drug or metabolite concentration of each sample in the unit ng/mL. Results equal to or greater than the respective cutoff value are considered preliminary positive. Concentration values below the respective cutoff indicate a negative result.

Preliminary positive results should be confirmed by another method.

The semiquantitation of preliminary positive results should only be used by laboratories to determine an appropriate dilution of the specimen for confirmation by a confirmatory method such as GC-MS. It also permits the laboratory to establish quality control procedures and assess control performance.

NOTE: If a result of Calc.? or Samp.? alarm is obtained, review the Reaction Monitor data for the sample and compare with the Reaction Monitor data for the highest calibrator. The most likely cause is a high concentration of the analyte in the sample, in which case the absorbance value for the sample will be less than that of the highest calibrator. Make an appropriate dilution of the sample using the 0 ng/mL calibrator and rerun the sample. A normal drug-free urine may be substituted for the 0 ng/mL calibrator if the urine and procedure have been validated by the laboratory. To ensure that the sample was not over-diluted, the diluted result, prior to multiplying by the dilution factor, must be at least half the analyte cutoff value. If the diluted result falls below half the analyte cutoff value, repeat the sample with a smaller dilution. A dilution that produces a result closest to the analyte cutoff is the most accurate estimation. To estimate the preliminary positive sample's concentration, multiply the result by the appropriate dilution factor. Dilutions should only be used to interpret results of Calc.? or Samp.? alarms, or when estimating concentration in preparation for GC-MS or LC-MS/MS.

As with any sensitive test for drugs of abuse on automated clinical chemistry analyzers, the possibility exists for analyte carry-over from a sample with an extremely high concentration to a normal (negative) sample which immediately follows it.

Use caution when reporting results as there are various factors that influence a urine test result, such as fluid intake and other biological factors.

Limitations - interference

See the "Specific performance data" section of this document for information on substances tested with this assay. There is the possibility that other substances and/or factors may interfere with the test and cause erroneous results (e.g., technical or procedural errors).

A preliminary positive result with this assay indicates the presence of benzodiazepines and/or their metabolites in urine. It does not reflect the degree of intoxication.

Interfering substances were added to urine containing nordiazepam at -25 % and +25 % of the cutoff level at the concentration listed below. Samples were tested and the following results were obtained on a Roche/Hitachi 917 analyzer.

Semiquantitative (ng/mL)		100 ng/mL cutoff		200 ng/mL cutoff		300 ng/mL cutoff	
Compound	Cmpd. conc.	Neg level	Pos level	Neg level	Pos level	Neg level	Pos level
Acetone	1 %	77 (NEG)	134 (POS)	157 (NEG)	260 (POS)	231 (NEG)	402 (POS)

Ascorbic acid	1.5 %	78 (NEG)	132 (POS)	156 (NEG)	262 (POS)	233 (NEG)	399 (POS)
Conjugated bilirubin	0.25 mg/mL	82 (NEG)	129 (POS)	156 (NEG)	247 (POS)	229 (NEG)	392 (POS)
Creatinine	5 mg/mL	81 (NEG)	138 (POS)	158 (NEG)	259 (POS)	230 (NEG)	396 (POS)
Ethanol	1 %	78 (NEG)	136 (POS)	151 (NEG)	261 (POS)	228 (NEG)	395 (POS)
Glucose	20 mg/mL	81 (NEG)	138 (POS)	158 (NEG)	262 (POS)	236 (NEG)	403 (POS)
Hemoglobin	1 mg/mL	76 (NEG)	139 (POS)	159 (NEG)	261 (POS)	228 (NEG)	398 (POS)
Human serum albumin	5 mg/mL	83 (NEG)	140 (POS)	165 (NEG)	273 (POS)	243 (NEG)	422 (POS)
Oxalic acid	2 mg/mL	74 (NEG)	128 (POS)	151 (NEG)	254 (POS)	226 (NEG)	388 (POS)
Sodium chloride	0.5 M	79 (NEG)	139 (POS)	159 (NEG)	262 (POS)	234 (NEG)	389 (POS)
Urea	6 %	80 (NEG)	138 (POS)	157 (NEG)	261 (POS)	233 (NEG)	405 (POS)

The same experiment was performed in the qualitative mode for each cutoff. All negative and positive samples recovered properly in the presence of the interfering substance.

An additional protocol was executed in which samples containing nordiazepam at control levels (± 25 % of cutoff) with specific gravities ranging from 1.006 to 1.034 were tested. As with the other interferences, there were no control cross-overs on any of the 3 assay cutoffs at either extreme specific gravity level.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

ACTION REQUIRED

Special Wash Programming: The use of special wash steps is mandatory when certain test combinations are run together on **cobas c** systems. All special wash programming necessary for avoiding carry-over is available via the **cobas** link. The latest version of the carry-over evasion list can be found with the NaOH/SMS/SCCS Method Sheet. For further instructions, refer to the operator's manual.

Expected values*Qualitative assay*

Results of this assay distinguish preliminary positive (≥ 100 ng/mL, ≥ 200 ng/mL, or ≥ 300 ng/mL depending on the cutoff) from negative samples only. The amount of drug detected in a preliminary positive sample cannot be estimated.

Semiquantitative assay

Results of this assay yield only approximate cumulative concentrations of the drug and its metabolites (see Analytical specificity section).

Specific performance data

Representative performance data on the analyzers are given below. These data represent the performance of the analytical procedure itself.

Results obtained in individual laboratories may differ due to heterogenous sample materials, aging of analyzer components and mixture of reagents running on the analyzer.

Precision

Precision was determined using human samples and controls in accordance with the CLSI (Clinical and Laboratory Standards Institute) EP05-A3 requirements with repeatability (n = 84) and intermediate precision (2 aliquots per run, 2 runs per day, 21 days). Results for repeatability and intermediate precision were obtained on the **cobas c** 503 analyzer.

Semiquantitative precision - 100 ng/mL cutoff

<i>Repeatability</i>	<i>Mean</i>	<i>SD</i>	<i>CV</i>
	<i>ng/mL</i>	<i>ng/mL</i>	<i>%</i>

ONLINE DAT Benzodiazepines II

Urine -50 %	49.5	3.30	6.7
DATCN	79.5	3.10	3.9
DAT2N	80.0	3.51	4.4
Cutoff urine	105	3.68	3.5
DAT2P	129	3.38	2.6
DATCP	128	3.51	2.7
Urine +50 %	153	2.98	1.9

<i>Intermediate precision</i>	<i>Mean ng/mL</i>	<i>SD ng/mL</i>	<i>CV %</i>
Urine -50 %	49.5	4.00	8.1
DATCN	79.5	4.09	5.1
DAT2N	80.0	4.28	5.4
Cutoff urine	105	4.82	4.6
DAT2P	129	4.37	3.4
DATCP	128	4.31	3.4
Urine +50 %	153	4.58	3.0

Qualitative precision - 100 ng/mL cutoff

Cutoff (100)	Number tested	Correct results	Confidence level
Urine -50 %	84	84	> 95 % negative reading
DATCN	84	84	> 95 % negative reading
DAT2N	84	84	> 95 % negative reading
Cutoff urine	84	n.a.*	n.a.*
DAT2P	84	84	> 95 % positive reading
DATCP	84	84	> 95 % positive reading
Urine +50 %	84	84	> 95 % positive reading

*n.a. = not applicable

Semiquantitative precision - 200 ng/mL cutoff

<i>Repeatability</i>	<i>Mean ng/mL</i>	<i>SD ng/mL</i>	<i>CV %</i>
Urine -50 %	104	2.49	2.4
DAT3N	152	3.78	2.5
Cutoff urine	200	2.94	1.5
DAT3P	249	3.61	1.4
Urine +50 %	308	2.52	0.8

<i>Intermediate precision</i>	<i>Mean ng/mL</i>	<i>SD ng/mL</i>	<i>CV %</i>
Urine -50 %	104	4.20	4.0
DAT3N	152	4.31	2.8
Cutoff urine	200	5.85	2.9
DAT3P	249	3.94	1.6
Urine +50 %	308	3.71	1.2

Qualitative precision - 200 ng/mL cutoff

Cutoff (200)	Number tested	Correct results	Confidence level
Urine -50 %	84	84	> 95 % negative reading
DAT3N	84	84	> 95 % negative reading

Cutoff urine	84	n.a.*	n.a.*
DAT3P	84	84	> 95 % positive reading
Urine +50 %	84	84	> 95 % positive reading

*n.a. = not applicable

Semiquantitative precision - 300 ng/mL cutoff

<i>Repeatability</i>	<i>Mean ng/mL</i>	<i>SD ng/mL</i>	<i>CV %</i>
Urine -50 %	150	7.03	4.7
DAT1N	231	8.04	3.5
Cutoff urine	345	7.59	2.2
DAT1P	375	7.66	2.0
Urine +50 %	435	8.89	2.0

<i>Intermediate precision</i>	<i>Mean ng/mL</i>	<i>SD ng/mL</i>	<i>CV %</i>
Urine -50 %	150	10.5	7.0
DAT1N	231	10.0	4.3
Cutoff urine	345	10.6	3.1
DAT1P	375	9.26	2.5
Urine +50 %	435	11.7	2.7

Qualitative precision - 300 ng/mL cutoff

Cutoff (300)	Number tested	Correct results	Confidence level
Urine -50 %	84	84	> 95 % negative reading
DAT1N	84	84	> 95 % negative reading
Cutoff urine	84	n.a.*	n.a.*
DAT1P	84	84	> 95 % positive reading
Urine +50 %	84	84	> 95 % positive reading

*n.a. = not applicable

The data obtained on **cobas c 503** analyzer(s) are representative for **cobas c 303** analyzer(s).

Accuracy

Additional clinical samples were evaluated with this assay on a **cobas c 503** analyzer and on a **cobas c 501** analyzer. 104 urine samples screened negative in a drug test panel, were evaluated with the Benzodiazepines II assay. 100 % of these normal urines were negative for all cutoffs relative to the **cobas c 503** analyzer. 76 urine samples for the 100 ng/mL cutoff, 68 urine samples for the 200 ng/mL cutoff and 51 urine samples for the 300 ng/mL cutoff, screened preliminary positive with a commercially available immunoassay and subsequently confirmed by GC-MS, were evaluated with the Benzodiazepines II assay. At the 100 ng/mL, 200 ng/mL and 300 ng/mL cutoffs, 100 % of the samples were positive on both the **cobas c 501** analyzer and the **cobas c 503** analyzer.

Benzodiazepines II correlation (cutoff = 100 ng/mL)			
		cobas c 501 analyzer	
		+	-
cobas c 503 analyzer	+	76	0
	-	0	104

Benzodiazepines II correlation (cutoff = 200 ng/mL)			
		cobas c 501 analyzer	
		+	-
cobas c 503 analyzer	+	68	0
	-	0	104

Benzodiazepines II correlation (cutoff = 300 ng/mL)			
		cobas c 501 analyzer	
		+	-
cobas c 503 analyzer	+	51	0
	-	0	104

Additional clinical samples were evaluated with this assay on a **cobas c 303** analyzer and on a **cobas c 501** analyzer. 110 urine samples screened negative in a drug test panel, were evaluated with the Benzodiazepines II assay. 100 % of these normal urines were negative for all cutoffs relative to the **cobas c 303** analyzer. 55 urine samples for the 100 ng/mL cutoff, 55 urine samples for the 200 ng/mL cutoff and 55 urine samples for the 300 ng/mL cutoff, screened preliminary positive with a commercially available immunoassay and subsequently confirmed by GC-MS, were evaluated with the Benzodiazepines II assay. At the 100 ng/mL, 200 ng/mL and 300 ng/mL cutoffs, 100 % of the samples were positive on both the **cobas c 501** analyzer and the **cobas c 303** analyzer.

Benzodiazepines II correlation (cutoff = 100 ng/mL)			
		cobas c 501 analyzer	
		+	-
cobas c 303 analyzer	+	55	0
	-	0	110

Benzodiazepines II correlation (cutoff = 200 ng/mL)			
		cobas c 501 analyzer	
		+	-
cobas c 303 analyzer	+	55	0
	-	0	110

Benzodiazepines II correlation (cutoff = 300 ng/mL)			
		cobas c 501 analyzer	
		+	-
cobas c 303 analyzer	+	55	0
	-	0	110

Analytical specificity

The specificity of the Benzodiazepines II assay for various benzodiazepines and benzodiazepine metabolites was determined by generating inhibition curves for each of the compounds listed and determining the approximate quantity of each compound that is equivalent in assay reactivity to a 100 ng/mL, 200 ng/mL, and 300 ng/mL nordiazepam assay cutoff. The following results were obtained on Roche/Hitachi and **cobas c** analyzers.

Compound ^{a)}	ng/mL Equivalent to 100 ng/mL nordiazepam	Approx. % cross- reactivity
Deschloroetizolam	80	125
Flubromazepam	94	107
3-OH-Flubromazepam	126	79
Clonazolam	96	104
Pyrazolam	105	95
Diclazepam	118	85
Flubromazolam	119	84
Etizolam	122	82
Meclonazepam	132	76
Nifoxipam	157	64
Bentazepam	173	58
Estazolam	93	107

Bromazepam	101	99
Nitrazepam	104	96
7-Aminonitrazepam	71	141
7-Acetamidonitrazepam	16909	0.59
Oxazepam	105	95
Oxazepam glucuronide	234	43
Phenazepam	112	89
Alprazolam	113	89
α-Hydroxyalprazolam	115	87
4-Hydroxyalprazolam	117	86
Demoxepam	114	88
Clorazepate	115	87
Clobazam	122	82
Diazepam	128	78
Nordiazepam	101	99
Delorazepam	131	76
Temazepam	133	75
Temazepam glucuronide	302	33
Triazolam	136	74
α-Hydroxytriazolam	145	69
Flunitrazepam	136	73
7-Aminoflunitrazepam	109	92
Desmethylflunitrazepam	114	88
Lormetazepam	138	73
Brotiazolam	144	70
Clonazepam	152	66
7-Aminoclonazepam	107	94
Lorazepam	153	65
Lorazepam glucuronide	275	36
Chlordiazepoxide	156	64
Desmethylchlordiazepoxide	138	73
Norchlordiazepoxide	150	67
Pinazepam	160	63
Flurazepam	164	61
Desalkylflurazepam	106	95
Hydroxyethylflurazepam	127	79
Didesethylflurazepam	144	70
Desmethylmedazepam	168	59
Halazepam	187	53
Midazolam	190	53
α-Hydroxymidazolam	125	80
Prazepam	194	51
Nimetazepam	1045	10
Oxaprozin	2283	4
Zolpidem	106383	0.09

a) Indented compounds are metabolites of the preceding drug.

Compound ^{a)}	ng/mL Equivalent to 200 ng/mL nordiazepam	Approx. % cross- reactivity
------------------------	--	--------------------------------------

ONLINE DAT Benzodiazepines II

Deschloroetizolam	159	126	Chlordiazepoxide	371	54
Flubromazepam	180	111	Desmethylchlordiazepoxide	313	64
3-OH-Flubromazepam	246	81	Norchlordiazepoxide	360	56
Clonazolam	185	108	Prazepam	408	49
Pyrazolam	188	106	Desmethylmedazepam	422	47
Flubromazolam	221	91	Nimetazepam	2191	9
Diclazepam	225	89	Oxaprozin	7500	3
Etizolam	234	86	Zolpidem	206186	0.10
Meclonazepam	329	61	a) Indented compounds are metabolites of the preceding drug.		
Bentazepam	376	53	Compound^{a)}	ng/mL	Approx.
Nifoxipam	391	51		Equivalent to	%
Estazolam	197	101		300 ng/mL	cross-
Bromazepam	208	96		nordiazepam	reactivity
Oxazepam	224	89	Deschloroetizolam	242	124
Oxazepam glucuronide	506	40	Flubromazepam	274	110
Clorazepate	227	88	3-OH-Flubromazepam	358	84
Phenazepam	230	87	Pyrazolam	279	107
Alprazolam	236	85	Clonazolam	290	103
α -Hydroxyalprazolam	241	83	Diclazepam	346	87
4-Hydroxyalprazolam	246	81	Etizolam	343	88
Nitrazepam	243	82	Flubromazolam	351	85
7-Aminonitrazepam	159	126	Meclonazepam	424	71
7-Acetamidonitrazepam	55488	0.36	Bentazepam	504	60
Demoxepam	253	79	Nifoxipam	552	54
Clobazam	256	78	Bromazepam	299	100
Diazepam	258	78	Estazolam	303	99
Nordiazepam	204	98	Oxazepam	325	92
Delorazepam	258	77	Oxazepam glucuronide	684	44
Triazolam	279	72	Phenazepam	346	87
α -Hydroxytriazolam	287	70	Demoxepam	352	85
Temazepam	282	71	Nitrazepam	354	85
Temazepam glucuronide	647	31	7-Aminonitrazepam	218	138
Flunitrazepam	284	70	7-Acetamidonitrazepam	55328	0.54
7-Aminoflunitrazepam	244	82	Alprazolam	372	81
Desmethylflunitrazepam	248	81	4-Hydroxyalprazolam	342	88
Lormetazepam	284	70	α -Hydroxyalprazolam	347	86
Brotiazolam	292	68	Clorazepate	374	80
Clonazepam	318	63	Clobazam	386	78
7-Aminoclonazepam	232	86	Delorazepam	389	77
Flurazepam	333	60	Diazepam	400	75
Desalkylflurazepam	225	89	Nordiazepam	316	95
Hydroxyethylflurazepam	259	77	Lormetazepam	410	73
Didesethylflurazepam	297	67	Temazepam	416	72
Lorazepam	335	60	Temazepam glucuronide	923	33
Lorazepam glucuronide	584	34	Triazolam	425	71
Midazolam	343	58	α -Hydroxytriazolam	440	68
α -Hydroxymidazolam	253	79	Flunitrazepam	439	68
Halazepam	354	56	Desmethylflunitrazepam	338	89
Pinazepam	364	55	7-Aminoflunitrazepam	368	82
			Brotiazolam	464	65

Clonazepam	483	62	Chloroquine	Niacinamide
7-Aminoclonazepam	334	90	Chlorpheniramine	Nicotine
Chlordiazepoxide	499	60	Chlorpromazine	Norethindrone
Desmethylchlordiazepoxide	452	66	Cocaine	<i>l</i> -Norpseudoephedrine
Norchlordiazepoxide	483	62	Codeine	Omeprazole
Lorazepam	506	59	Desipramine HCl	Penicillin G
Lorazepam glucuronide	825	36	Dextromethorphan	Pentazocine
Flurazepam	511	59	Dextropropoxyphene	Pentobarbital
Desalkylflurazepam	336	89	Digoxin	Phencyclidine
Hydroxyethylflurazepam	394	76	Diphenhydramine	Phenobarbital
Didesethylflurazepam	458	65	Diphenylhydantoin	Phenothiazine
Desmethylnedazepam	539	56	Doxepin	Phenylbutazone
Midazolam	564	53	Ecgonine	<i>d,l</i> -Phenylpropanolamine
α -Hydroxymidazolam	428	70	Ecgonine methyl ester	Procaine
Pinazepam	572	52	Enalapril	Promethazine
Halazepam	595	50	<i>d</i> -Ephedrine	<i>d</i> -Pseudoephedrine
Prazepam	637	47	<i>l</i> -Ephedrine	Quinidine
Nimetazepam	3247	9	Epinephrine	Quinine
Oxaprozin	7507	4	Erythromycin	Secobarbital
Zolpidem	200000	0.15	Estriol	Sulindac

a) Indented compounds are metabolites of the preceding drug.

Many benzodiazepines appear in the urine largely as the glucuronidated conjugate. Glucuronidated metabolites may have more or less cross-reactivity than the parent compound. The presence of β -glucuronidase enzyme enhances the Benzodiazepines II assay cross-reactivity to some of the glucuronidated metabolites.

Drug interference

The following compounds were prepared in aliquots of pooled normal human urine to yield a final concentration of 100000 ng/mL. None of these compounds gave values in the assay that were greater than 0.08 % cross-reactivity for the 100 ng/mL and 200 ng/mL cutoffs and 0.13 % cross-reactivity for the 300 ng/mL cutoff.

Acetaminophen	Imipramine
Acetylsalicylic acid	Isoproterenol
Amitriptyline	Ketamine
Amobarbital	Lidocaine
<i>d</i> -Amphetamine	LSD
<i>l</i> -Amphetamine	MDA
Ampicillin	MDMA
Ascorbic acid	Melanin
Aspartame	Meperidine
Atropine	Methadone
Benzocaine	<i>d</i> -Methamphetamine
Benzoylcegonine (cocaine metabolite)	<i>l</i> -Methamphetamine
Benzphetamine	Methaqualone
Buspirone	Methylphenidate
Butabarbital	Methyprylon
Caffeine	Morphine sulfate
Calcium hypochlorite	Naloxone
Cannabidiol	Naltrexone
Captopril	Naproxen

Fenoprofen	Tetracycline
Flumazenil	Δ^9 THC-9-carboxylic acid
Furosemide	Tetrahydrozoline
Gentisic acid	Thioridazine
Glutethimide	Tolmetin
Guaiacol glycerol ether	Trifluoperazine
Hydrochlorothiazide	Trimipramine
Hydroxyindole acetic acid	Tyramine
Hydroxyindole carboxylic acid	Verapamil
Ibuprofen	Zomepirac

In very rare cases, gammopathy, in particular type IgM (Waldenström's macroglobulinemia), may cause unreliable results.¹⁸

References

- 1 Karch SB, ed. Drug Abuse Handbook. Boca Raton, FL: CRC Press LLC 1998.
- 2 Salamone SJ, ed. Benzodiazepines and GHB: Detection and Pharmacology. Totowa, NJ: Humana Press 2001.
- 3 Hardman JG, Limbird LE, Gilman A, eds. Goodman and Gilman's The Pharmacological Basis of Therapeutics. 10th ed. New York, NY: McGraw Hill Pub Co. 2001.
- 4 Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 7th ed. Foster City, CA: Biomedical Publications 2004.
- 5 Laurijsens BE, Greenblatt DJ. Pharmacokinetic-pharmacodynamic relationships for benzodiazepines. Clin Pharmacokinet 1996;30:52-76.
- 6 Hallfors DD, Saxe L. The dependence potential of short half-life benzodiazepines: a meta-analysis. Am J Public Health 1993;83:1300-1304.
- 7 Chouinard G. Issues in the clinical use of benzodiazepines: potency, withdrawal, and rebound. J Clin Psychiatry 2004;65(5):7-12.
- 8 Abernethy DR, Greenblatt DJ, Ochs HR, et al. Benzodiazepine drug-drug interactions commonly occurring in clinical practice. Curr Med Res Opin 1984;8:80-93.
- 9 Tanaka E. Toxicological interactions between alcohol and benzodiazepines. J Toxicol Clin Toxicol 2002;40:69-75.

- 10 Dou C, Bournique JS, Zinda MK, et al. Comparison of the Rates of Hydrolysis of Lorazepam-Glucuronide, Oxazepam-Glucuronide and Temazepam-Glucuronide Catalyzed by E. Coli β -D-Glucuronidase Using the OnLine Benzodiazepine Screening Immunoassay on the Roche/Hitachi 917 Analyzer. J of Forensic Science 2001;46(2):335-340.
- 11 Beck O, Lin Z, Brodin K, et al. The online screening technique for urinary benzodiazepines: comparison with EMIT, FPIA, and GC-MS. J Anal Toxicology 1997;21(7):554-557.
- 12 Salamone SJ, Honasoge S, Brenner C, et al. Flunitrazepam excretion patterns using the Abuscreen OnTrak and OnLine immunoassays: comparison with GC-MS. J Anal Toxicol 1997;21:341-345.
- 13 Klette KL, Wiegand RF, Horn CK, et al. Urine benzodiazepine screening using Roche Online KIMS immunoassay with beta-glucuronidase hydrolysis and confirmation by gas chromatography-mass spectrometry. J Anal Toxicol 2005;29:193-200.
- 14 Valentine JL, Middleton R, Sparks C. Identification of urinary benzodiazepines and their metabolites: comparison of automated HPLC and GC-MS after immunoassay screening of clinical specimens. J Anal Toxicol 1996;20(6):416-424.
- 15 Armbruster DA, Schwarzhoff RH, Hubster EC, et al. Enzyme immunoassay, kinetic microparticle immunoassay, radioimmunoassay, and fluorescence polarization immunoassay compared for drugs-of-abuse screening. Clin Chem 1993;39:2137-2146.
- 16 Toxicology and Drug Testing in the Clinical Laboratory; Approved Guideline. 2nd ed. (C52-A2). Clinical and Laboratory Standards Institute 2007;27:33.
- 17 Mandatory Guidelines for Federal Workplace Drug Testing Programs. Fed Regist 2017 Jan 23;82:7920-7970.
- 18 Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007;45(9):1240-1243.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard (for USA: see dialog. Roche.com for definition of symbols used):

CONTENT	Contents of kit
→	Volume for reconstitution
GTIN	Global Trade Item Number

ABUSCREEN, COBAS, COBAS C, ONLINE DAT and PRECISET are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2022, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.Roche.com

+800 5505 6606

