



LightCycler® Multiplex DNA Master

 **Version 03**

Content version: October 2015

Easy-to-use 5× reaction mix optimized for multiplex qPCR, compatible with the LightCycler® 480, LightCycler® 96, or the LightCycler® 2.0 Real-Time PCR Systems

Cat. No. 07 339 585 001

1 kit
200 reactions of 20 µl final volume each

Cat. No. 07 339 577 001

1 kit
1,000 reactions of 20 µl final volume each

Store the kit at – 15 to – 25°C!

1.	What this Product Does	3
	Number of Tests	3
	Kit Contents	3
	Storage and Stability	3
	Additional Equipment and Reagents Required	3
	Optional	4
	Application	4
	Assay Time	4
2.	How to Use this Product	5
2.1	Before You Begin	5
	Precautions	5
	Sample Material	5
	No Template Controls	5
	Primers	5
	Probes	6
	MgCl ₂	6
2.2	Procedure	6
	LightCycler® 480 and LightCycler® 96 System Protocol	6
	Setup of the qPCR Reaction for the LightCycler® 480 and LightCycler® 96 Instrument	12
	Color Compensation Protocol for the LightCycler® 2.0 Instrument	14
	Protocol for the LightCycler® 2.0 Instrument	15
	Preparation of the Color Compensation Mixes	16
	Setup of the qPCR Reaction for the LightCycler® 2.0 Instrument	17
2.3	Quality Control	18
3.	Results	19
4.	Troubleshooting	21
5.	Supplementary Information	23
5.1	Conventions	23
5.1.1	Text Conventions	23
5.1.2	Symbols	23
5.2	Changes to Previous Version	23
5.3	Ordering Information	24
5.4	Disclaimer of License	25
5.5	Trademarks	26
5.6	Regulatory Disclaimer	26

1. What this Product Does

Number of Tests The kit is designed for 200 or 1,000 reactions with a final reaction volume of 20 μ l each.

Kit Contents

Vial/Cap	Label	Function	Catalog Number
1 red cap	qPCR Reaction Mix, 5 $\times$ conc.	Contains qPCR Reaction Buffer, Aptaq Polymerase, dATP, dCTP, dGTP, and dUTP, MgCl ₂ , and proprietary additives	Cat. No. 07 339 585 001, 1 vial, 880 μ l Cat. No. 07 339 577 001, 5 vials, 880 μ l each
2 colorless cap	Water, PCR grade		Cat. No. 07 339 585 001, 3 vials, 1 ml each Cat. No. 07 339 577 001, 15 vials, 1 ml each

Storage and Stability

This kit is shipped on dry ice.

Store the kit at -15 to -25°C . This kit is stable until the expiration date printed on the label.

The kit is stable at $+2$ to $+8^{\circ}\text{C}$ for 4 weeks.

For Vial 1, avoid repeated freeze/thaw cycles (more than 5 $\times$). Aliquot Vial 1 and freeze or store at $+2$ to $+8^{\circ}\text{C}$ for a maximum of 4 weeks.

Although we recommend working on ice and preparing the reagents right before use, the working solution is stable at $+15$ to $+25^{\circ}\text{C}$ up to 4 hours and is therefore ideal for use in automated workflows.

Additional Equipment and Reagents Required

- Standard laboratory equipment
- Nuclease-free pipette tips
- 1.5 ml RNase-free microcentrifuge tubes to prepare master mixes and dilutions.
- To minimize risk of RNase contamination, autoclave all vessels.

For qPCR:

- Real-Time PCR systems such as LightCycler® 480, LightCycler® 96 or the LightCycler® 2.0 Instrument*
 - LightCycler® 480 Multiwell Plate 96* or 384*
 - LightCycler® Capillaries (20 µl)*
 - Standard swinging-bucket centrifuge with rotor for multiwell plates
-
- For qPCR primer and probe design:
Universal ProbeLibrary Assay Design Center at
www.universalprobelibrary.com
-
- To design and order qPCR assays and panels:
www.realtimeready.roche.com

Optional

For DNA purification

- MagNA Pure 96 Instrument* including disposables
 - MagNA Pure 96 Internal Control Tube*, optional
 - MagNA Pure 96 DNA and Viral NA Kit, Large Volume* or
 - MagNA Pure 96 DNA and Viral NA Kit, Small Volume*
- Alternatively, use other nucleic acid purification instruments and manual sample preparation products such as:
- MagNA Pure LC 2.0 Instrument* with MagNA Pure LC Total Nucleic Acid Isolation Kit*
 - MagNA Pure Compact Instrument* with MagNA Pure Compact Nucleic Acid Isolation Kit I
 - High Pure Viral Nucleic Acid Kit*

Application

The LightCycler® Multiplex DNA Master is designed for fast, highly sensitive and specific real-time PCR analysis of DNA.

The single-vial master mix allows fast and convenient hot start qPCR without the need for upfront polymerase activation incubation. The mix is optimized for use with hydrolysis probes as well as Universal ProbeLibrary (UPL) probes, and does not require optimization of MgCl₂ concentration.

Assay Time

The LightCycler® Multiplex DNA Master can be used for multiplex qPCR protocols. For example, a triplex protocol using 45 cycles requires less than 60 minutes when using the LightCycler® 480 System.

2. How to Use this Product

2.1 Before You Begin

Precautions

Always use nuclease-free techniques. Nuclease contaminated reagents and reaction vessels will degrade template nucleic acids. Please follow these guidelines to minimize risk of contamination:

- Wear disposable gloves and change them frequently.
- Avoid touching surfaces or materials that could cause nuclease carry-over.
- Use only reagents provided in this kit. Substitutions may introduce nucleases.
- Clean and decontaminate work areas and instruments, including pipettes, with commercially available decontamination reagents.
- Use only new nuclease-free aerosol-blocking pipette tips.
- Use a work area specifically designated for nucleic acid work, and if possible use reaction vessels and pipettors dedicated only for work with template nucleic acid.

Sample Material

Use any DNA suitable for qPCR in terms of purity, concentration, and absence of PCR inhibitors. For reproducible isolation of nucleic acids, use either the MagNA Pure LC Instrument*, the MagNA Pure Compact Instrument*, or the MagNA Pure 96 Instrument*, and a dedicated MagNA Pure Nucleic Acid Isolation Kit* (for automated isolation), or a High Pure Nucleic Acid Isolation Kit* (for manual isolation).

For details, see the Roche Catalog or visit the Life Science website: www.lifescience.roche.com

No Template Controls

We highly recommend testing known negative specimens as controls in each run to control for possible contamination. To prepare a no template control, replace the template DNA with PCR grade Water (Vial 2).

Primers

Suitable concentrations of PCR primers range from 0.2 to 0.5 μM final concentration. The recommended starting concentration is 0.5 μM each.

Probes

Suitable concentrations of hydrolysis probes range from 0.05 to 0.5 μM (final concentration in PCR). The recommended starting concentration is 0.25 μM each.

- ③ The optimal probe concentration is the lowest concentration that results in the lowest quantification cycle value (Cq) and adequate fluorescence dynamics for a given target concentration.
- ③ For a hydrolysis probe hybridization complex, the T_m of the hydrolysis probe has to be higher than the T_m of the primers.

MgCl₂

The master mix is optimized with a fixed concentration of MgCl₂, and should not need concentration optimization for most assays.

2.2 Procedure

LightCycler® 480 and LightCycler® 96 System Protocol

The following procedure is optimized for use with the corresponding LightCycler® System you are using.

⚠ Program the LightCycler® Instrument before preparing the reaction mixes.

For details on how to program the experimental protocol, see the current LightCycler® 480 Instrument Operator's Manual, LightCycler® 96 System Operator's Guide, or LightCycler® 2.0 Instrument Operator's Manual B.

The LightCycler® 480 Instrument II protocol contains the following program:

- **Denaturation** of DNA
- **Amplification** of the DNA
- **Cooling** of the thermal block

A) Protocol for use with the LightCycler® 480 Instrument II (Multiwell Plate 96 or 384)

Setup				
Detection Format		Reaction Volume	Block Type	
Hydrolysis Probe or UPL Probe		20 μl (10 μl)	96 (384)	
Programs				
Program Name		Cycles	Analysis Mode	
Initial Denaturation		1	None	
Amplification		45 ¹⁾	Quantification	
Cooling		1	None	
Temperature Targets				
Target [°C]	Acquisition Mode	Hold [hh:mm:ss]	Ramp Rate [°C/s]	Acquisitions [n/°C]
Initial Denaturation				
95	None	00:00:30	4.4 (4.8)	–
Amplification				
95	None	00:00:05	4.4 (4.8)	–
60 ²⁾	Single	00:00:30	2.2 (2.5)	–
Cooling				
40	None	00:00:30	2.2 (2.5)	–

¹⁾ 45 cycles is suitable for most assays. If the assay shows steep amplification curves and early crossing points, 40 cycles should be sufficient. Reducing the number of cycles will improve uniformity of the product and reduce the time required for the assay.

²⁾ The LightCycler® Multiplex DNA Master includes an aptamer for hot start. Anneal / Extend temperatures significantly below 60°C may not be well tolerated because aptamers will inhibit DNA polymerase activity at low temperatures.

Color Compensation Protocol for the LightCycler® 480 Instrument II

For a multicolor, multiplex amplification, the application of a color compensation file is necessary to compensate for optical crosstalk between the different channels.

For the LightCycler® 480 Instrument II, an instrument-specific color compensation file is mandatory for multicolor experiments, and a color compensation object can be generated by performing the following experiment.

The LightCycler® 480 Instrument II protocol contains the following program:

- **Initial Denaturation** of the DNA
- **Amplification** of the DNA
- **Temperature Gradient Step** to create the Color Compensation file
- **Cooling** of the thermal block

For details on how to program the experimental protocol, see the LightCycler® 480 Software Operator's Manual, Version 1.5.

Programming a Customized Detection Format for the LightCycler® 480 System Filter Combination Selection

The Detection Format in the LightCycler® 480 Software, Version 1.5 setup must be customized for the applied multicolor hydrolysis probe format used in the qPCR detection. In the "Tool" module, the "Detection Formats" option allows creating a new detection format specified by the user, including a "Detection Format" list, a "Filter Combination" selection area, and a "Selected Filter Combination List". Different filter settings for the LightCycler® 480 II Instrument are defined. Use the following filter combination:

Example for 3-color Hydrolysis Probes

Detection Formats	Excitation Filter	Emission Filter
FAM	465	510
Red 610	533	610
Cy5	618	660

For the new customized detection format, set the following values for all selected filters in the "Selected Filter Combination List":

Melt Factor	1
Quant Factor	10
Max Integration Time	2

The following table shows the qPCR parameters that must be programmed for a LightCycler® 480 Color Compensation file with a LightCycler® 480 Multiwell Plate 96.

Setup			
Detection Format		Block Type	
Customized (see section above)		96	
Programs			
Program Name	Cycles	Analysis Mode	
Initial Denaturation	1	None	
Amplification	45	Quantification	
Temperature Gradient Step	1	Color Compensation	
Cooling	1	None	
Temperature Targets			
Target [°C]	Acquisition Mode	Hold [hh:mm:ss]	Ramp Rate [°C/s]
Initial Denaturation			
95	None	00:00:30	4.4
Amplification			
95	None	00:00:05	4.4
60	Single	00:00:30	4.4
Temperature Gradient Step			
95	None	00:00:10	4.4
40	None	00:00:10	2.2
95	Continuous		5 Acq. /°C
Cooling			
40	None	00:00:30	2.2

Preparation of the Color Compensation Run

Prepare the calibrator qPCR mix for more than one reaction. Multiply the amount in the volume column by the number of reactions you need (minimum of 3 to 5 replicates) plus additional reactions since there will be a slight loss of liquid during the pipetting steps. In order to ensure accuracy, do not pipette volumes less than 1 µl when adding the individual reagents.

For each dye, set up the following reactions:

Component	1× Buffer	1× for each Dye
qPCR Reaction Mix, 5× conc.(Vial 1)	4.0 µl	4.0 µl
Detection mix for each dye	-	X µl (depending on the assay)
Water, PCR grade (Vial 2)	16.0 µl	Y µl (depending on the assay)
Template, such as DNA or positive samples eluates	-	5.0 µl
Total Volume	20 µl	20 µl

Step	Action
①	Pipette the replicates of each different calibrator mix into a pre-cooled LightCycler® 480 Multiwell Plate 96.
②	Seal the LightCycler® 480 Multiwell Plate using a sealing foil.
③	Place the multiwell plate in a standard swinging-bucket centrifuge containing a rotor for multiwell plates with suitable adaptors. Balance it with a suitable counterweight (e.g., another plate), and centrifuge for 2 minutes at 1,500 × g.
④	Load the multiwell plate into the LightCycler® 480 Instrument II and start the program.

Create Color Compensation Object

When the experiment is finished, open the **Sample Editor**. Select workflow **Color Comp** and designate the used positions of the Multiwell Plate as "Water" for buffer replicates, and the appropriate dyes respectively (e.g., FAM, Red610, Cy5 for the example mentioned above).

Open the **Analysis** module "Color Compensation", click on **Calculate**, and save the color compensation object in the database.

From this point on, you can apply this "CC Object" to each dual-, triple- or quadruple-color experiment performed with FAM, Yellow555, Red610, and Cy5 on the same instrument.

B) Protocol for use with the LightCycler® 96 Instrument

Run Editor			
Detection Format			Reaction Volume
Select dyes used in your assays. If using the DNA Process Control, select Cy5 in Channel 4.			20 µl
Programs			
Temp. [°C]	Ramp Rate [°C/s]	Hold [s]	Acquisition Mode
Preincubation (Initial Denaturation)			
95	4.4	30	None
2-Step Amplification			
No. of Cycles: 45 ¹⁾			
95	4.4	5	None
60 ²⁾	2.2	30	Single

¹⁾ 45 cycles is suitable for most assays. If the assay shows steep amplification curves and early crossing points, 40 cycles should be sufficient. Reducing the number of cycles will improve uniformity of the product and reduce the time required for the assay.

²⁾ The LightCycler® Multiplex DNA Master includes an aptamer for hot start. Anneal / Extend temperatures significantly below 60°C may not be well tolerated because aptamers will inhibit DNA polymerase activity at low temperatures.

Color Compensation Protocol for the LightCycler® 96 Instrument

The LightCycler® 96 Instrument does not require the creation of a color compensation object.

**Setup of the
qPCR Reaction
for the
LightCycler® 480
and
LightCycler® 96
Instrument**

Follow the procedure below to prepare at least ten 20 µl standard reactions. 10 µl amplifications included in brackets to be used for 384-well plate set-ups:

 Do not touch the surface of the LightCycler® 480 Multiwell Plate.

- 1

Thaw the solutions and, to ensure recovery of all the contents, briefly spin vials in a microcentrifuge before opening.
Mix carefully by pipetting up and down or vortexing briefly. Place samples on ice.
- 2

Prepare a 20× conc. solution of your primers and a 20× conc. solution of your probes.
- 3

In a 1.5 ml reaction tube, prepare the qPCR Mix and put on ice. For best results, prepare at least 10 reactions in order to reduce pipetting errors. To prepare more reactions, multiply the amount in the "Volume 1 Reaction" column below by the number of reactions to be run, plus at least one additional reaction.

Component ¹⁾	Volume 1 Reaction	Volume 10 Reactions	Final Conc.
Water, PCR grade (Vial 2)	9 µl (4.5 µl)	90 µl (45 µl)	-
qPCR Reaction Mix, 5× (Vial 1)	4 µl (2 µl)	40 µl (20 µl)	1×
Primer Mix, 20×	1 µl (0.5 µl)	10 µl (5 µl)	1×
Probe Mix, 20×	1 µl (0.5 µl)	10 µl (5 µl)	1×
Total volume	15 µl (7.5 µl)	150 µl (75 µl)	

¹⁾ For eluates derived from stool samples, it is recommended to add 0.2 µg/µl (final) of molecular biology grade Bovine Serum Albumin. Adjust for BSA volume by subtracting from PCR Water volume.

- 4

Mix carefully by pipetting up and down or vortexing briefly. Place on ice.
Although we recommend working on ice and preparing the reagents right before use, the working solution (everything combined except DNA template) is stable at +15 to +25°C up to 4 hours, and is therefore ideal for use in automated workflows.
- 5

Prepare sample concentration of the DNA.
- 6

Pipette 15 µl (7.5 µl) qPCR Mix into a precooled multiwell plate. Add 5 µl (2.5 µl) of the DNA template.
Seal multiwell plate with LightCycler® 480 Sealing Foil.
- 7

Place the Multiwell Plate 96 into a standard swinging-bucket centrifuge with a suitable adapter and balance it with a suitable counterweight (e.g., another multiwell plate).
Centrifuge at 1,500 × g for 0.5 - 2 minutes.
- 8

Load the reaction vessels into the LightCycler® 480 or LightCycler® 96 Instrument.
- 9

Start the PCR program described above.
If you use reaction volumes other than 20 µl, be sure to adapt the right volume in the running protocol. To start, we recommend using the same hold times as for the 20 µl volume.

C) Protocol for Use with the LightCycler® 2.0 Instrument

A LightCycler® 2.0 Instrument protocol that uses LightCycler® Multiplex DNA Master contains the following programs:

- **Denaturation** of the DNA
- **Amplification** of the DNA
- **Cooling** the rotor and thermal chamber

The following table shows the PCR parameters that must be programmed for a LightCycler® 2.0 Instrument PCR run with the LightCycler® Multiplex DNA Master using LightCycler® Capillaries (20 µl).

LightCycler® Software Version 4.1 ¹⁾				
Programs				
Setup ¹⁾	Setting			
Default Channel	Fluorescence Channel			
Seek Temperature	30°C			
Max Seek Pos.	Enter the total number of sample positions for which the instrument should look for.			
Instrument Type	"6 Ch." for LightCycler® 2.0 Instrument			
Capillary Size	Select "20 µl" as the capillary size for the experiment.			
Programs				
Program Name	Cycles	Analysis Mode		
Initial Denaturation	1	None		
Amplification	45 ²⁾	Quantification		
Cooling	1	None		
Temperature Target				
	Target [°C]	Hold [hh:mm:ss]	Ramp Rate [°C/s]	Acquisition Mode [per °C]
Initial Denaturation	95	00:00:30	20	None
Amplification	95	00:00:05	20	None
	60 ³⁾	00:00:30	20	Single
Cooling	40	00:00:30	20	None

¹⁾ The LightCycler® Multiplex DNA Master includes aptamer for hot-start. Anneal / Extend temperatures significantly below 60°C may not be well tolerated because aptamer will inhibit DNA polymerase activity at low temperatures.

²⁾ 45 cycles are suitable for most assays. If the assay shows steep amplification curves and early crossing points, 40 cycles should be sufficient. Reducing the number of cycles will improve uniformity of the product and reduce the time required for the assay.

³⁾ Most available assays are designed for an annealing temperature of +60°C. If the T_m of the primers and probe indicate using different settings, select an annealing temperature +5°C below the calculated primer T_m for initial experiments. For assay optimization, choose the highest annealing temperature which still results in a high yield of amplicon with a low C_q and adequate fluorescence dynamics.

Color Compensation Protocol for the LightCycler® 2.0 Instrument

The following procedure is optimized for use with the LightCycler® 2.0 System. Program the LightCycler® 2.0 Instrument before preparing the reaction mixes.

A LightCycler® 2.0 Instrument color compensation protocol that uses LightCycler® Multiplex DNA Master contains the following programs:

- **Initial Denaturation** of the DNA
- **Amplification** of the DNA
- **Temperature Gradient Step** for Color Compensation
- **Cooling** the rotor and thermal chamber

For details on how to program the experimental protocol, see the LightCycler® 2.0 Instrument Operator's Manual B.

Color Compensation Protocol

The performance of a color compensation is a prerequisite for running a dual color experiment. The generated color compensation file is used to compensate for cross-talk between the individual detection channels when performing multi-color experiments. A color compensation calibration run is performed by running a blank capillary (containing Water, PCR grade), and individual capillaries with one dye each (monocolor PCR reactions), in a RT-PCR program, followed by a color compensation analysis.

**Protocol for the
LightCycler® 2.0
Instrument**

The following tables show the parameters that must be programmed for a LightCycler® Color Compensation calibration run with LightCycler® Multiplex DNA Master.

LightCycler® Software Version 4.1				
Programs				
Setup	Setting			
Default Channel	Fluorescence Channel			
Seek Temperature	30°C			
Max Seek Pos.	"Enter the total number of sample positions for the instrument should look for.			
Instrument Type	"6 Ch." for LightCycler® 2.0 Instrument			
Capillary Size	Select “20 µl” as the capillary size for the experiment.			
Programs				
Program Name	Cycles	Analysis Mode		
Initial Denaturation	1	None		
Amplification	45	Quantification		
Temperature Gradient	1	Color Compensation		
Cooling	1	None		
Temperature Target				
	Target [°C]	Hold [hh:mm:ss]	Ramp Rate [°C/s]	Acquisition Mode [per °C]
Initial Denaturation	95	00:00:05	20	None
Amplification	95	00:00:01	20	None
	60	00:00:15	20	Single
Temperature Gradient	95	00:00:01	20	None
	40	00:01:00	20	None
	95	00:00:00	0.2	Continuous
Cooling	40	00:00:30	20	None

Preparation of the Color Compensation Mixes

 Do not touch the surface of the LightCycler® Capillaries.

- 1

Place three LightCycler® Capillaries (20 µl) into precooled LightCycler® Centrifuge Adapters.
- 2

Prepare the capillaries (20 µl, each), as follows:
- 3

For each dye, set up the following reactions:

Reagent	Volume for each Dye	Capillary with Water
qPCR Reaction Mix, 5× conc (Vial 1)	4 µl	
Detection mix for each dye	X µl (depending on the assay)	
Water, PCR grade (Vial 2)	Y µl (depending on the assay)	20 µl
Bovine Serum Albumin (20 µg/µl) ¹⁾	0.2 µl	
Template, such as DNA or positive sample eluates	5 µl	
Total Volume	20 µl	20 µl

¹⁾ Molecular biology grade Bovine Serum Albumin is recommended.

- 4

Seal each capillary with a stopper using the LightCycler® Capping Tool.
- 5

Place the centrifuge adapters (containing the capillaries) into a standard benchtop microcentrifuge. Place the centrifuge adapters in a balanced arrangement within the centrifuge. Centrifuge at 700 × g for 5 s (3,000 rpm in a standard benchtop microcentrifuge). Alternatively, use the LC Carousel Centrifuge for spinning the capillaries.
- 6

Place the capillaries in the following order in the sample carousel of the LightCycler® 2.0 Instrument:
Carousel rotor position 1: Water
Carousel rotor position 2: monocolour PCR for Dye 1
Carousel rotor position 3: monocolour PCR for Dye 2

Cycle the samples as described above and edit the dominant channel in the "Analysis Type" - "Color Comp" accordingly.

Create Color Compensation Objects

When the experiment is finished, click on the Analysis button and select Color Compensation (Other Methods) from the Analysis Menu. Save the experiment by clicking the Save CC Object button. Place the object in the "Special Data\CCC" folder under your user name.

After doing this, you can apply the specific "CC Object" you created to any dual-color hydrolysis probe experiment that is performed with the same dye combination.

Setup of the qPCR Reaction for the LightCycler® 2.0 Instrument

🕒 This setup can also be used in a qPCR protocol for the LightCycler® 1.x Instrument.

Follow the procedure below to prepare at least ten 20 µl standard reactions:

⚠️ Do not touch the surface of the LightCycler® Capillaries.

- 1 Thaw the solutions and, to ensure recovery of all the contents, briefly spin vials in a microcentrifuge before opening.
Mix carefully by pipetting up and down or vortexing briefly. Place samples on ice.
- 2 Prepare a 20× conc. solution of your primers and a 20× conc. solution of your probes.
- 3 In a 1.5 ml reaction tube, prepare the PCR Mix and put on ice. For best results, prepare at least 10 reactions in order to reduce pipetting errors. To prepare more reactions, multiply the amount in the "Volume 1 Reaction" column below by the number of reactions to be run, plus at least one additional reaction.

Reagent	Volume 1 Reaction	Volume 10 Reactions	Final concentration
Water, PCR grade (Vial 2)	8.7 µl	87 µl	-
qPCR Reaction Mix, 5 x (Vial 1)	4 µl	40 µl	1x
Primer Mix, 20x ¹⁾	1 µl	10 µl	1x
Probe Mix 20x	1 µl	10 µl	1x
Bovine Serum Albumin (20 mg/ml) ²⁾	0.2 µl	2 µl	0.2 µg/µl
Total Volume	15 µl	150 µl	

¹⁾ Due to possible primer/primer interactions that occur during storage, it may be necessary to preheat the PCR primer mix for 1 minute at +95°C before pipetting the complete mixture. This extra step will ensure optimum sensitivity.

²⁾ Molecular biology grade Bovine Serum Albumin is recommended.

- ④ Mix carefully by pipetting up and down or vortexing briefly. Place on ice. Although we recommend working on ice and preparing the reagents right before use, the working solution (everything combined except DNA template) is stable at +15 to +25°C up to 4 hours, and is therefore ideal for use in automated workflows.
- ⑤ Prepare sample concentration of the DNA.
- ⑥ Pipette 15 µl PCR Mix into a LightCycler® Capillary. Add 5 µl of the DNA template.
- ⑦ Seal the LightCycler® Capillaries with a stopper, using the LightCycler® Capping Tool.
- ⑧ If using the LC Carousel Centrifuge 2.0, proceed to step 10.
- ⑨ Alternatively, place the capillaries in cooled adapters in a standard benchtop microcentrifuge in a balanced arrangement, centrifuge at $700 \times g$ (3,000 rpm) for 5 s, and transfer the capillaries to the LightCycler® Sample Carousel.
- ⑩ Place the LightCycler® Sample Carousel in the LightCycler® Carousel-Based Instrument.
- ⑪ Start the PCR program.

2.3 Quality Control

Each lot of LightCycler® Multiplex DNA Master is tested to meet specifications of the qPCR using a duplex qPCR assay on the LightCycler® 480 Instrument II.

3. Results

The following results were obtained using the LightCycler® Multiplex DNA Master on the LightCycler® 480 Instrument. A duplex reaction using β -globin specific assay with BHQ-2 internally quenched probe (FAM) and β 2M specific assay with UPL probe (Yellow555).

Human reference cDNA served as the template.

FAM channel (465-510)




Fig. 1: The FAM channel shows the results for β -globin. Amplification curves shown were generated from a four-fold dilution series [10 ng (far left), 2.5 ng, 625 pg, 156 pg, 40 pg, 10 pg, and 2.4 pg (far right)] of Human Reference cDNA. No template controls are shown in black (flat line). Duplex qPCR was performed in a reaction volume of 20 μ l per well in a LightCycler® 480 Multi-well Plate 96.


Yellow 555 channel (533/580)

Fig. 2: The Yellow555 channel shows the results for $\beta 2M$. Amplification curves shown were generated from a four-fold dilution series [10 ng (far left), 2.5 ng, 625 pg, 156 pg, 39 pg, 10 pg, and 2.4 pg (far right)] of Human Reference cDNA. No template controls are shown in black (flat line). Duplex qPCR was performed in a reaction volume of 20 μ l per well in a LightCycler[®] 480 Multiwell Plate 96.

4. Troubleshooting

	Possible Cause	Recommendation
Fluorescence intensity varies	Some of the reagent is still in the upper part of the microwell or an air bubble is trapped in the microwell.	Repeat centrifugation, but allow sufficient centrifugation time so all reagent is at the bottom of the microwell and air bubbles are expelled.
	Skin oils or dirt on the surface of the microwell plate.	Always wear gloves when handling the multiwell plate.
Fluorescence intensity is very low	Low concentration or deterioration of dyes in the reaction mixtures because dye was not stored properly.	• Keep dye-labeled reagents away from light. • Store the reagents at -15 to -25°C and avoid repeated freezing and thawing.
	Poor PCR efficiency (reaction conditions not optimized)	• Primer concentration should be in the range of 0.2 to 0.5 µM, probe concentration should be in the range of 0.2 to 0.5 µM and half of the primer concentration. • Check annealing temperature of primers and probes. • Check experimental protocol. • Optimize annealing temperature in the reverse transcription step or in the PCR reaction. • Always run a positive control along with your samples.
	Chosen imaging time is too low	• Choose the appropriate detection format in combination with “dynamic” detection mode or • Increase imaging time when using “manual” detection mode. For details, see LightCycler® 480 Software Instrument Operator's Manual.
	qPCR primers and probes are not optimized	• Check sequence and location of the hydrolysis probe on the PCR product. • Check PCR product on an agarose gel.
	PCR has not been optimized	• Check primer design (quality). • Check PCR product on an agarose gel.

	DNA is degraded during isolation or improper storage	• If possible, check DNA quality.
	Pipetting errors and omitted reagents	• Check for missing reagents. • Check the pipetting procedure.
	Impure sample material inhibits reaction	• Dilute sample 1:10 and repeat the analysis. • Repurify the nucleic acids to ensure removal of inhibitory agents.
Negative control sample gives a positive signal	Contamination	Remake all critical reaction mixes. Be sure to use special Pre-PCR setup working areas.

5. Supplementary Information

5.1 Conventions

5.1.1 Text Conventions

To make information consistent and easier to read, the following text conventions are used in this document:

Text Convention	Usage
Numbered instructions labeled 1 , 2 etc.	Steps in a procedure that must be performed in the order listed.
Asterisk *	Denotes a product available from Roche Diagnostics.

5.1.2 Symbols

In this document, the following symbols are used to highlight important information:

Symbol	Description
	Information Note: Additional information about the current topic or procedure.
	Important Note: Information critical to the success of the procedure or use of the product.

5.2 Changes to Previous Version

- Information about the sample material 'template DNA derived from stool samples' has been added.
- LightCycler® 2.0 System Protocol has been added.
- Editorial Changes.

5.3 Ordering Information

Roche offers a large selection of reagents and systems for life science research. For a complete overview of related products and manuals, please visit and bookmark our homepage www.lifescience.roche.com:

- LightCycler® 480 System: www.lightcycler480.com
- LightCycler® 96 System: www.lightcycler96.com
- Automated Sample Preparation (MagNA Pure LC System and MagNA Pure Compact System): www.magnapure.com
- Manual Sample Preparation of Nucleic Acids:
www.lifescience.roche.com/shop/en/us/products/manual-sample-preparation-of-nucleic-acids

Instrument and Accessories

Product	Pack Size	Cat. No.
LightCycler® 480 Instrument II, 96 well	1 instrument (96 well)	05 015 278 001
LightCycler® 480 Instrument II, 384 well	1 instrument (384 well)	05 015 243 001
LightCycler® 480 Block Kit 96 Silver	1 block kit for 96-well PCR Multiwell Plates	05 015 219 001
LightCycler® 480 Block Kit 384 Silver	1 block kit for 384 -well PCR Multiwell Plates	05 015 197 001
LightCycler® 480 Multiwell Plate 96, white	5 × 10 plates with sealing foils	04 729 692 001
LightCycler® 480 Multiwell Plate 384, white	5 × 10 plates with sealing foils	04 729 749 001
LightCycler® 480 Sealing Foil	50 foils	04 729 757 001
LightCycler® 480 Software, Version 1.5	1 software package	04 994 884 001
LightCycler® 96 Instrument	1 instrument	05 815 916 001
LightCycler® 2.0 Instrument	1 instrument plus related products and data station	03 531 414 001
LC Carousel Centrifuge 2.0	1 centrifuge plus rotor and rotor bucket (230 V)	03 709 582 001
	1 centrifuge plus rotor and rotor bucket (115 V)	03 709 507 001

Nucleic Acid Isolation Kits

Product	Pack Size	Cat. No.
LightCycler® Capillaries (20 µl)	1 pack containing 5 boxes, each with 96 capillaries and stoppers	04 929 292 001
LightCycler® Centrifuge Adapters	1 set 32 adapters in a dedicated aluminum cooling block	11 909 312 001
LightCycler® Capping Tool	1 capping tool	03 357 317 001
MagNA Pure 96 Instrument	1 instrument plus accessories	06 541 089 001
MagNA Pure 96 Internal Control Tube	150 tubes (15×10)	05 435 293 001
MagNA Pure LC 2.0 Instrument	1 instrument plus accessories	05 197 686 001
MagNA Pure Compact Instrument	1 instrument plus accessories	03 731 146 001
MagNA Pure 96 DNA and Viral NA Small Volume Kit	3 sets for 192 isolations each	06 543 588 001
MagNA Pure 96 DNA and Viral NA Large Volume Kit	3 sets for 96 isolations each	06 374 891 001
MagNA Pure LC Total Nucleic Acid Isolation Kit	1 kit (192 isolations)	03 038 505 001
MagNA Pure Compact Nucleic Acid Isolation Kit I	1 kit for 32 isolations	03 730 964 001
High Pure Viral Nucleic Acid Kit	1 kit for up to 100 purifications	11 858 874 001
Associated Kits and Reagents	Transcriptor Universal cDNA Master	100 reactions 04 379 012 001
Transcriptor First Strand cDNA Synthesis Kit	1 kit for up to 50 reactions, including 10 control reactions	04 379 012 001
Bovine Serum Albumin	1 ml 20 mg/ml	10 711 454 001

5.4 Disclaimer of License

NOTICE: For patent license limitations for individual products please refer to www.technical-support.roche.com.

5.5 Trademarks

LIGHTCYCLER, MAGNA PURE, and HIGH PURE are trademarks of Roche.
FAM is a trademark of Applied Biosystems Corporation or its subsidiaries in the US and/or certain other countries.
All other product names and trademarks are the property of their respective owners.

5.6 Regulatory Disclaimer

For general laboratory use.

Contact and Support

If you have questions or experience problems with this or any Roche product for Life Science, please contact our Technical Support staff. Our scientists are committed to providing rapid and effective help.

Please also contact us if you have suggestions for enhancing Roche product performance or using our products in new or specialized ways. Such customer information has repeatedly proven invaluable to the research community worldwide.

To ask questions, solve problems, suggest enhancements or report new applications, please visit our [Online Technical Support Site](#).

Visit lifescience.roche.com, to download or request copies of the following Materials:

- Instructions for Use
- Safety Data Sheets
- Certificates of Analysis
- Information Material

To call, write, fax, or email us, visit lifescience.roche.com and select your home country to display country-specific contact information.

