Accutrend Cholesterol cobas® SYSTEM

| | | |
|----------|----|-------------------------------|
| REF | 25 | Accutrend Plus |
| 11418254 | 5 | Accutrend GC Accutrend GCT |

English
Intended use
Test strip for the quantitative determination of cholesterol in fresh or heparinised fresh capillary blood. Use only with the following meters: Accutrend GC, Accutrend GCT or Accutrend Plus. Suitable for self-testing.

Summary
The determination of cholesterol provides the physician with initial information on the state of the patient's lipid metabolism. This information serves as a basis for further diagnostic measures, therapeutic decisions or monitoring the progress of the patient. However, self-monitoring is not to be seen as a substitute for regular monitoring by a physician. The results of self-monitoring should be noted down and discussed with the physician.

How the tests work
Each test strip has a test area containing detection reagents. When blood is applied, a chemical reaction takes place and the test area changes colour. The meter records this change in colour and converts the measurement signal to the displayed result using the data previously entered via the code strip.

Preparing and performing the measurement
In order to avoid erroneous measurements, use only the above specified meters for determinations with Accutrend Cholesterol test strips. Please follow the instructions given in this Method Sheet and the operator's manual of your meter.

Precautions and warnings
For in vitro diagnostic use.
All components of the pack can be discarded along with household waste or, if used in a laboratory or doctor's office, according to relevant local guidelines.
Never allow liquids or disinfectants to get into the vial.
The stopper contains a non-toxic silicate-based desiccant. If the desiccant is inadvertently ingested, please drink plenty of water!

Storage and stability
Store at 2-30 °C up to the stated expiration date. Do not freeze. Store away from extreme temperatures.
Do not use the test strip after the specified expiration date.

Tightly re-cap the container immediately after removing a test strip.
Perform cholesterol determinations at 18-30 °C.

If you perform determinations above or below this temperature range, the result obtained may only be used for orientation as it may deviate considerably from the true value depending on the temperature.

Important: If you use more than one type of Accutrend test strip at the same time, make sure that you do not mix up the stoppers. The stoppers on the Accutrend Cholesterol and Accutrend Triglycerides test strip containers are blue; those on the Accutrend Glucose test strip containers are white. If the container stoppers are mixed up, the specified shelf-life can no longer be guaranteed.

Materials provided
▪ Test strips and 1 code strip

Materials required (but not provided)
▪ [REF 11418289190](#), Accutrend Control CH 1 control solution
▪ Lancing device
▪ One of the above listed meters

Assay
For optimum performance of the system follow the directions given in this document for the corresponding meter. Refer to the appropriate meter manual for test-specific instructions. The system requires a hanging drop of blood. Use only the above specified meters for determinations with these test strips to avoid erroneous results.

Check required every time you perform a test
Before taking a measurement, the test strip must be checked to see that it is still fit for use. If the window on the back of the test strip has taken on a bluish colour, there is a risk that the cholesterol values displayed are too high. If this occurs, use a new test strip. Every time you perform a test, make sure the test area is completely covered with blood. Follow also the instruction given in the operating instructions of your meter. Occasionally, the round window on the back of the test strip may be mottled; this uneven coloration is sample-related and does not interfere with the measurement of the Accutrend Cholesterol test strip.

Coding
Whenever a new pack of test strips is opened your meter must be recoded with the code strip provided. The procedure is described in the meter manual. If the code number displayed does not match the vial you are using, the meter will not perform a measurement. Always keep the code strip until the last test strip of a package is used.

Store the code strip outside the test strip container. If the code strip is stored in the container, the test strips may be damaged.

Quality control
For quality control, use Accutrend Control CH 1.

The acceptable control range is given in the enclosed value sheet.

The target value corresponds to the mean of the lower and upper value of this control range. The control intervals and limits should be adapted to each laboratory's individual requirements.

Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

Limitations - interference
The cholesterol determination may be affected by:

- intravenous infusion of ascorbic acid (vitamin C)
- bilirubin values > 10 mg/dL (171 µmol/L), e.g. in hepatitis
- haematocrit values exceeding 55 %
- methylenoantipyrine
- gentisic acid

Limits and ranges

| | mmol/L | mg/dL |
|--|-----------|---------|
| Linear measuring range | 3.88-7.76 | 150-300 |
| Detection limit (lowest value displayed) | 3.88 | 150 |

Expected values
Cholesterol values less than 5.2 mmol/L (200 mg/dL) are considered normal.

In some patient groups, the reference ranges may be not be within the normal range. You should discuss this with your physician.

Instructions for healthcare professionals

Test principle
Enzymatic cleavage of the cholesterol ester into fatty acids and cholesterol; oxidation of cholesterol to cholestenone with the simultaneous formation of hydrogen peroxide which oxidises an indicator to its blue radical cation.^{1,2}

Applications
For the determination of cholesterol in fresh capillary blood in diabetics and non-diabetics; for self-monitoring; for the early detection of a risk of atherosclerosis; for monitoring treatment with lipid-lowering drugs; for screening purposes. If measurements are performed on several persons, as may be the case in a doctor's office or for screening programmes, the fresh capillary blood may also be applied with heparin-coated capillary pipettes (volume 15 µL to 40 µL).

Precautions and warnings

Exercise the normal precautions required for handling all laboratory reagents. Health care professionals performing tests on more than one patient must be aware that there is a potential risk of infection. Any object coming into contact with human blood is a potential source of infection.³

Safety data sheet available for professional user on request.

Reagents

Components per test:
CHE (microorganisms) 1.44 U; CHOD (Nocardia erythropolis) 0.12 U; POD (horseradish) 0.20 U; 3,3',5,5'-tetramethylbenzidine 8.5 µg; non-reactive components 1.4 mg.

Limitations - interference

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

Specific performance data

Precision

Repeatability (within-run precision)

In investigations with EDTA venous blood in a concentration range between 4.32 mmol/L (167 µg/dL) and 6.75 mmol/L (291 µg/dL) typical measurement series showed a CV (coefficient of variation) of 0.8 to 3.7 %.

Intermediate precision (between-day precision)

The intermediate precision was determined in test series with control solution. At a concentration of 4.96 mmol/L (191 µg/dL) the CV was typically between 1.1 % and 3.8 %, at a concentration of approx. 7.02 mmol/L (270 µg/dL) between 2.4 % and 5.0 %.

Method comparison

In method comparisons with capillary blood in hospitals > 95 % of all values were found within ±15 % of the reference method (CHOD-PAP, Roche Diagnostics). The mean systematic differences were between +2.6 % and -3.6 %.

Point sentence

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Significant additions or changes are indicated by a change bar in the margin.

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