TRAPtest

REF	CONTENT	SYSTEM
08847509190	3 x → 1.0 mL	Multiplate

Česky

Použití

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Odkazy

- Khan ML, Zheng Y-W, Huang W, et al. A dual thrombin receptor system for platelet activation. *Nature*. 1998;394:690–694.
- Kinlough-Rathbone RL, Perry DW, Guccione MA, et al. Degranulation of human platelets by the thrombin receptor peptide SFLLRN: comparison with degranulation by thrombin. *Thromb Haemost*. 1993 Dec 20;70(6):1019-23.
- Desch S, Siegemund A, Scholz U, et al. Platelet inhibition and GP IIb/IIIa receptor occupancy by intracoronary versus intravenous bolus administration of abciximab in patients with ST-elevation myocardial infarction. *Clin Res Cardiol* 2011 Oct 21.
- Halimeh S, Angelis G, Sander A, et al. Multiplate Whole Blood Impedance Point of Care Aggregometry: Preliminary Reference Values in Healthy Infants, Children, and Adolescents. *Klin Padiatr* 2010; 222: 158-163.
- Penz S, Perchuc A, Calatzis A, et al. Effects of GpIIb/IIIa and GpIb blockade on multiple electrode aggregometry (MEA). *GTH 2010; Abstract P10-13*.
- Gryka RJ, Buckley LF, Anderson SM. Vorapaxar: The Current Role and Future Directions of a Novel Protease-Activated Receptor Antagonist for Risk Reduction in Atherosclerotic Disease. *Drugs R D*. 2017 Mar;17(1):65-72
- Jambor C, Weber CF, Gerhardt K, et al. Whole Blood Multiple Electrode Aggregometry Is a Reliable Point-of-Care Test of acetylsalicylic acid-Induced Platelet Dysfunction. *Anesth Analg* 2009;109:25-31.
- von Pape KW, Dzijan-Horn M, Bohner J, et al. Control of acetylsalicylic acid effect in chronic cardiovascular patients using two whole blood platelet function assays. PFA-100 and Multiplate. *Hamostaseologie*. 2007 Aug;27(3):155-60.
- Velik-Salchner C, Maier S, Innerhofer P, et al. Point-of-care whole blood impedance aggregometry versus classical light transmission aggregometry for detecting acetylsalicylic acid and clopidogrel: the results of a pilot study. *Anesth Analg* 2008 Dec;107(6):1798-806.
- Johnson A, Dovlatova N, Heptinstall S. Multiple electrode aggregometry and P2Y(12) antagonists. *Thromb Haemost*. 2008 Jun;99(6):1127-9.
- Sibbing D, Stegheer J, Braun S, et al. A double-blind, randomized study on prevention and occurrence of a rebound phenomenon of platelets after cessation of clopidogrel treatment. *Am Coll Cardiol* 2010 Feb 9;55(6):558-65.
- Djukanovic N, Todorovic Z, Obradovic S, et al. Abrupt Cessation of One-Year Clopidogrel Treatment Is Not Associated With Thrombotic Events. *J Pharmacol Sci* 2011;117(1):12-8.
- Sibbing D, Braun S, Jawansky S, et al. Assessment of ADP-induced platelet aggregation with light transmission aggregometry and multiple electrode platelet aggregometry before and after clopidogrel treatment. *Thromb Haemost* 2008; 99: 121-126.
- CLSI. Platelet Function Testing by Aggregometry; Approved Guideline. CLSI document H58-A. Clinical and Laboratory Standards Institute; 2008.
- Kaiser AF, Neubauer H, Franken CC, et al. Which is the best anticoagulant for whole blood aggregometry platelet function testing? Comparison of six anticoagulants and diverse storage conditions. *Platelets*. 2011 Oct 14.
- Kalb ML, Potura L, Scharbert G et al. The effect of ex vivo anticoagulants on whole blood platelet aggregation. *Platelets*. 2009 Feb;20(1):7-11.
- Tóth O, Calatzis A, Penz S et al. Multiple electrode aggregometry: a new device to measure platelet aggregation in whole blood. *Thromb Haemost*. 2006 Dec;96(6):781-8.
- Jennings LK, White MM. Platelet aggregation. In: Michelson AD, ed. *Platelets*. 2nd ed. Amsterdam: Elsevier; 2007:495-518.
- Velik-Salchner C, Maier S, Innerhofer P et al. Point-of-care whole blood impedance aggregometry versus classical light transmission aggregometry for detecting acetylsalicylic acid and clopidogrel: the results of a pilot study. *Anesth Analg*. 2008 Dec;107(6):1798-806.
- Paniccia R, Antonucci E, Maggini N et al. Assessment of platelet function on whole blood by multiple electrode aggregometry in high-risk patients with coronary artery disease receiving antiplatelet therapy. *Am J Clin Pathol*. 2009 Jun;131(6):834-42.
- CLSI. Platelet Function Testing by Aggregometry; Approved Guideline. CLSI document H58-A. Clinical and Laboratory Standards Institute; 2008, p.31.
- Scharf RE Drugs that Affect Platelet Function *Semin Thromb Hemost*. 2012 Nov 38 (8): 865-83
- Hanke AA, Roberg K, Monaca E et al. Impact of platelet count on results obtained from multiple electrode platelet aggregometry (Multiplate). *Eur J Med Res* 2010 May 18;15(5):214-9.
- Behr T, Behr W & von Scheidt W. Monitoring of clopidogrel treatment by multiplate electrode platelet aggregometry. *J Lab Med* 2010;34(2), 99-107.
- DE MacIntyre, RL Hoover, M Smith, et al. Inhibition of platelet function by cis- unsaturated fatty acids. *Blood* 1984;63:848-57
- F.Driss, E.Vericel, M.Lagarde, et al. Inhibition of platelet aggregation and thromboxane synthesis after intake of small amount of icosapentaenoic acid. *Thromb Res* 1984;36:389-96

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© 2021, Roche Diagnostics
Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim