REF			SYSTEM
09315322190	09315322500	200	<b>cobas e 411</b> <b>cobas e 601</b> <b>cobas e 602</b>

## Svenska

### Systeminformation

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## Referenser

- 1 Katus HA, Remppis A, Looser S, et al. Enzyme linked immunoassay of cardiac troponin T for the detection of acute myocardial infarction in patients. *Mol Cell Cardiol* 1989;21(12):1349-1353.
- 2 Liebetrau C, Mollmann H, Nef H, et al. Release kinetics of biomarkers in patients undergoing transcatheter ablation of septal hypertrophy. *Clin Chem* 2012;58(6):1049-1054.
- 3 Katus HA, Scheffold T, Remppis A, et al. Proteins of the troponin complex. *Laboratory Medicine* 1992;23(5):311-317.
- 4 Reichlin T, Hochholzer W, Bassetti S, et al. Early diagnosis of myocardial infarction with sensitive cardiac troponin assays. *NEJM* 2009;361(9):858-867.
- 5 Giannitsis E, Becker M, Kurz K, et al. High sensitivity cardiac troponin T for early prediction of evolving non-ST segment elevation myocardial infarction in patients with suspected acute coronary syndrome and negative troponin result on presentation. *Clin Chem* 2010;56(4):642-650.
- 6 Giannitsis E, Kurz K, Hallermayer K, et al. Analytical validation of a high-sensitivity cardiac troponin T assay. *Clin Chem* 2010;56(2):254-261.
- 7 Hamm CW, Bassand JP, Agewall S, et al. ESC guidelines for the management of acute coronary syndrome in patients presenting without persistent ST-segment elevation. The Task Force for the management of acute coronary syndrome (ACS) in patient presenting without persistent ST-segment elevation of the European Society of Cardiology (ESC). *Eur Heart J* 2011;32(23):2999-3054.
- 8 NICE (2014) Myocardial infarction (acute): Early rule out using high-sensitivity troponin tests (Elecsys Troponin T high-sensitive, ARCHITECT STAT High Sensitive Troponin-I and AccuTnI+3 assays). NICE diagnostics guidance DG15. Available at [www.nice.org.uk/dg15](http://www.nice.org.uk/dg15) [NICE guideline]
- 9 2020 ESC Guidelines for the management of acute coronary syndromes in patients presenting without persistent ST-segment elevation. *Eur Heart J* 2020. PMID 32860058.
- 10 Bandstein N, Ljung R, Johansson M, et al. Undetectable high sensitivity cardiac troponin T level in the emergency department and risk of myocardial infarction. *J Am Coll Card* 2014;63:2569-2578.
- 11 Body R, Burrows G, Carley S, et al. High-sensitivity cardiac troponin T concentrations below the limit of detection to exclude acute myocardial infarction: A prospective evaluation. *Clin Chem* 2015;61(7):983-989.
- 12 Rubini Gimenez M, Reichlin T, Zellweger C. Rapid rule out of acute myocardial infarction using undetectable levels of high-sensitivity cardiac troponin. *Int J Cardiol* 2013;168(4):3896-3901.
- 13 Reichlin T, Schindler C, Drexler B, et al. One-hour rule-out and rule-in of acute myocardial infarction using high-sensitivity cardiac troponin T. *Arch Intern Med* 2012;172(16):1211-1218.
- 14 Reichlin T, Twerenbold R, Wildi K, et al. Prospective validation of a 1-hour algorithm to rule-out and rule-in acute myocardial infarction using a high-sensitivity cardiac troponin T assay. *CMAJ* 2015;187(8):E243-252.
- 15 Mueller C, Giannitsis E, Christ M, et al. Multicenter evaluation of a 0-hour/1-hour algorithm in the diagnosis of myocardial infarction with high-sensitivity cardiac Troponin T. *Ann Emerg Med* 2016; 68(1):76-87.
- 16 Mokhtari A, Borna C, Gilje P, et al. A 1-h combination algorithm allows fast rule-out and rule-in of major adverse cardiac events *J Am Coll Cardiol* 2016;67(13):1531-1540.
- 17 Shiozaki M, Inoue K, Sua S, et al. Utility of the 0-hour/1-hour high-sensitivity cardiac troponin T algorithm in Asian patients with suspected non-ST elevation myocardial infarction. *Int J Cardiol*. 2017; 249:32-35.
- 18 Chew DP, Lambrakis K, Blyth A, et al. A randomized trial of a 1-hour troponin protocol in suspected acute coronary syndromes: The Rapid Assessment of Possible ACS in the Emergency Department with High Sensitivity Troponin T (RAPID-TnT) study. *Circulation* 2019;140(19):1543-1556.
- 19 Twerenbold R, Costabel JP, Nestelberger T, et al. Outcome of applying the ESC 0/1-hour algorithm in patients with suspected myocardial infarction. *J Am Coll Cardiol* 2019;74(4):483-494.
- 20 Stoyanov KM, Hund H, Biener M et al. RAPID-CPU: a prospective study on implementation of the ESC 0/1-hour algorithm and safety of discharge after rule-out of myocardial infarction. *Eur Heart J Acute Cardiovasc Care*. 2019;9(1):39-51.
- 21 Parsonage W, Greenslade J, Hammett C, et al. Validation of an accelerated high-sensitivity troponin T assay protocol in an Australian cohort with chest pain. *Med J Aust* 2014;200:161-165.
- 22 Reichlin T, Cullen L, Parsonage WA, et al. Two hour algorithm for triage towards rule-out and rule-in of acute myocardial infarction using high-sensitivity cardiac Troponin T. *Am J Med* 2015;128:369-79.
- 23 Meller B, Cullen L, Parsonage WA, et al. Accelerated diagnostic protocol using high-sensitivity cardiac troponin T in acute chest pain patients. *Int J Cardiol*. 2015;184:208-15.
- 24 McRae AD, Innes G, Graham M, et al. Comparative Evaluation of 2-Hour Rapid Diagnostic Algorithms for Acute Myocardial Infarction Using High-Sensitivity Cardiac Troponin T. *Can J Cardiol*. 2017;33(8):1006-1012.
- 25 Wildi K, Cullen L, Twerenbold R, et al. Direct comparison of 2 rule-out strategies for acute myocardial infarction: 2-h accelerated diagnostic protocol vs 2-h algorithm. *Clin Chem*. 2017;63(7):1227-1236.
- 26 Than MP, Pickering JW, Dryden JM, et al. iCare-ACS (Improving Care processes for patients with suspected acute coronary syndrome): a study of cross-system implementation of a national clinical pathway. *Circulation* 2018;137(4):354-363.
- 27 Thygesen K, Alpert JS, Jaffe AS, et al. Fourth definition of myocardial infarction. *J Am Coll Cardiol* 2018;30:72(18):2231-2264.
- 28 Aviles RJ, Askari AT, Lindahl B, et al. Troponin T levels in patients with acute coronary syndromes, with or without renal dysfunction. *N Engl J Med* 2002;346:2047-2052.
- 29 Lindahl B, Venge P, James S. The new high- sensitivity cardiac troponin T assay improves risk assessment in acute coronary syndromes. *Am Heart J* 2010;160:224-229.
- 30 Haaf P, Reichlin T, Twerenbold R, et al. Risk stratification in patients with acute chest pain using three high-sensitivity cardiac troponin assays. *Eur Heart J* 2014;35(6):365-375.
- 31 Bjurman C, Larsson M, Johanson P, et al. Small changes in troponin T levels are common in patients with Non-ST-elevation myocardial infarction and are linked to higher mortality. *J Am Coll Cardiol* 2014;62(14):1231-1238.
- 32 Lindahl B, Venge P, Wallentin L. Troponin T Identifies Patients With Unstable Coronary Artery Disease Who Benefit From Long-Term Antithrombotic Protection. *J Am Coll Cardiol* 1997;29(1):43-48.
- 33 Hamm CW, Heeschen C, Goldmann B, et al. Benefit of abciximab in patients with refractory unstable angina in relation to serum troponin T levels. *N Engl J Med* 1999;340(21):1623-1629.
- 34 Heeschen C, Hamm CW, Goldmann BU, et al. for PRISM study investigators. Troponin concentrations for stratification of patients with acute coronary syndromes in relation to therapeutic efficacy of tirofiban. *Lancet* 1999;354:1757-1762.
- 35 Lindahl B, Diderholm E, Lagerquist B, et al. Effects on mortality of long-term treatment with l.m.w. heparin in relation to troponin T level and ECG findings - a FRISC 2 substudy. *Eur Heart J* 2000;21(Suppl.):521.
- 36 Newby LK, Ohman EM, Christenson RH, et al. Benefit of Glycoprotein IIb/IIIa Inhibition in Patients With Acute Coronary Syndromes and Troponin T-Positive Status: The PARAGON-B Troponin T Substudy. *Circulation* 2001;103:2891-2896.
- 37 Wallentin L, Lindholm D, Siegbahn A, et al. Biomarkers in relation to the effects of Ticagrelor compared with clopidogrel in non-ST-elevation acute coronary syndrome patients managed with or without in-hospital revascularization: a Substudy from the prospective randomized Platelet inhibition and Patient Outcomes (PLATO) Trial. *Circulation* 2014;129(3):293-303.
- 38 Amsterdam EA, Wenger NK, Brindis RG, et al. 2014 AHA/ACC guideline for the management of patients with non-ST elevation acute coronary syndromes: a report of the American College of Cardiology/American Heart Association Task Force on Practice Guidelines. *Circulation* 2014;130(25):e344-426.

- 39 Reichlin T, Irfan A, Twerenbold R, et al. Utility of absolute and relative changes in cardiac troponin concentrations in the early diagnosis of acute myocardial infarction. *Circulation* 2011;124(2):136-145.
- 40 Masson S, Anand I, Favero C, et al. Serial measurement of cardiac troponin T using a highly sensitive assay in patients with chronic heart failure: data from two large randomized clinical trials. *Circulation* 2012;125(2):280-288.
- 41 Nambi V, Liu X, Chambless LE, et al. Troponin T and N-terminal pro-B-type natriuretic peptide: a biomarker approach to predict heart failure risk: the Atherosclerosis Risk in Communities Study. *Clin Chem* 2013;59(12):1802-1810.
- 42 Yancy CW, Jessup M, Bozkurt B, et al. 2013 ACCF/AHA guideline for the management of heart failure: a report of the American College of Cardiology Foundation/American Heart Association Task Force on practice guidelines. *Circulation* 2013;128(16):e240-327.
- 43 Cramer G, Bakker J, Gommans F, et al. Relation of highly sensitive cardiac troponin T in hypertrophic cardiomyopathy to left ventricular mass and cardiovascular risk. *Am J Cardiol*. 2014;113(7):1240-1245.
- 44 McGill D, Talaulikar G, Potter JM, et al. Over time, high-sensitivity TnT replaces NT-proBNP as the most powerful predictor of death in patients with dialysis-dependent chronic renal failure. *Clin Chim Acta* 2010;411(13-14):936-939.
- 45 K/DOQI clinical practice guidelines for cardiovascular disease in dialysis patients. *Am J Kidney Dis*. 2005;45 (4 Suppl 3):S1-153.
- 46 Artunc F, Mueller C, Breidhardt T, et al. Sensitive troponins- which suits better for hemodialysis patients? Associated factors and prediction for mortality. *PLoS One*. 2012;7(10):e47610.
- 47 Wolley M, Steward R, Curry E, et al. Variation in and prognostic importance of troponin T measured using a high-sensitivity assay in clinically stable haemodialysis patients. *Clin Kidney J*. 2013;6(4):402-409.
- 48 Honneger Bloch S, Semple D, Sidhu K, et al. Prognostic value and long-term variation of high sensitivity troponin T in clinically stable haemodialysis patients. *N Z Med J*. 2014;127(1402):97-109.
- 49 Twerenbold R, Wildi K, Jeger C, et al. Optimal cutoff levels of more sensitive cardiac troponin assays for the early diagnosis of myocardial infarction in patients with renal dysfunction. *Circulation* 2015;131(23):2041-2050.
- 50 Landeberg G, Jaffe AS, Gilon D, et al. Troponin elevation in severe sepsis and septic shock: the role of left ventricular dysfunction and right ventricular dilatation. *Crit Care Med* 2014;42(4):790-800.
- 51 Hillis GS, Welsh P, Chalmers J, et al. The relative and combined ability of high sensitivity cardiac troponin T and N-terminal pro-BNP to predict cardiovascular events and death in patients with type 2 diabetes mellitus. *Diabetes Care* 2014;37(1):295-303.
- 52 Everett BM, Brooks MM, Vlachos HE, et al. Troponin and cardiac events in stable ischemic heart disease and diabetes. *N Engl J Med*. 2015;373(7):610-620.
- 53 Latini R, Masson S, Anand IS, et al. Prognostic Value of Very Low Plasma Concentrations of Troponin T in Patients with Stable Chronic Heart Failure. *Circulation* 2007;116:1242-1249.
- 54 Omland T, de Lemos JA, Sabatine MS, et al. A sensitive cardiac troponin T assay in stable coronary artery disease. *N Engl J Med* 2009;361:2538-2547.
- 55 Ponikowski P, Voors AA, Anker SD, et al. 2016 ESC Guidelines for the diagnosis and treatment of acute and chronic heart failure: The Task Force for the diagnosis and treatment of acute and chronic heart failure of the European Society of Cardiology (ESC). Developed with the special contribution of the Heart Failure Association (HFA) of the ESC. *Eur Heart J* 2016;37:2129-2200.
- 56 Latini R, Masson S, Pirelli S, et al. On the behalf of the GISSI-AF Investigators. Circulating cardiovascular biomarkers in recurrent atrial fibrillation: data from the GISSI-Atrial fibrillation trial. *J Intern Med* 2011;269(2):160-171.
- 57 Hijazi Z, Oldgren J, Lindbäck J, et al. The novel biomarker-based ABC (age, biomarkers, clinical history)-bleeding risk score for patients with atrial fibrillation: a derivation and validation study. *Lancet* 2016; 387:2302-2311.
- 58 Berg DD, Ruff C, Jarolim P, et al. Performance of the ABC scores for assessing the risk of stroke or systemic embolism and bleeding in patients with atrial fibrillation in ENGAGE AF-TIMI 48. <https://www.ahajournals.org/doi/10.1161/CIRCULATIONAHA.118.038312>.
- 59 Lewandrowski K. Special topics: cardiac markers in myocarditis. Cardiac transplant rejection and conditions other than acute coronary syndrome. *Clin Lab Med* 2014;34:129-135.
- 60 Swaanenburg JC, Klaase JM, DeJongste M, et al. Troponin I, troponin T, CK-MB-activity and CK-MB mass as markers for the detection of myocardial contusion in patients who experienced blunt trauma. *Clin Chim Acta* 1998;272:171-181.
- 61 Bajaj A, Saleeb M, Rathor P, et al. Prognostic value of troponins in acute noninvasive pulmonary embolism. A meta-analysis. *Heart Lung* 2015;44(4):327-334.
- 62 Dubin RF, Li Y, He J, et al. Predictors of high sensitivity cardiac troponin T in chronic kidney disease patients: a cross-sectional study in the chronic renal insufficiency cohort (CRIC). *BMC Nephrol* 2013;14(1):229.
- 63 Newby LK, Rodriguez I, Finkle J, et al. Troponin measurements during drug development-considerations for monitoring and management of potential cardiotoxicity. An educational collaboration among the Cardiac Safety Research Consortium, the Duke Clinical Research Institute, and the US Food and Drug Administration. *Am Heart J* 2011;162(1):64-73.
- 64 Perrone MA, Spolaore F, Ammirabile M, et al. The assessment of high sensitivity cardiac troponin in patients with COVID-19: A multicenter study. *Int J Cardiol Heart Vasc*. 2021;32:100715.
- 65 Guo T, Fan Y, Chen M, et al. Cardiovascular Implications of Fatal Outcomes of Patients With Coronavirus Disease 2019 (COVID-19). *JAMA Cardiol*. 2020;5(7):811-818.
- 66 Wei JF, Huang FY, Xiong TY, et al. Acute myocardial injury is common in patients with COVID-19 and impairs their prognosis. *Heart*. 2020;106(15):1154-1159.
- 67 De Michieli L, Ola O, Knott JD, et al. High-Sensitivity Cardiac Troponin T for the Detection of Myocardial Injury and Risk Stratification in COVID-19. *Clin Chem*. 2021.
- 68 Poterucha TJ, Elias P, Jain SS, et al. Admission Cardiac Diagnostic Testing with Electrocardiography and Troponin Measurement Prognosticates Increased 30-Day Mortality in COVID-19. *J Am Heart Assoc*. 2021;10(1):e018476.
- 69 Calvo-Fernandez A, Izquierdo A, Subirana I, et al. Markers of myocardial injury in the prediction of short-term COVID-19 prognosis. *Rev Esp Cardiol (Engl Ed)*. 2020.
- 70 Willeit P, Welsh P, Evans JDW, et al. High-Sensitivity Cardiac Troponin Concentration and Risk of First-Ever Cardiovascular Outcomes in 154,052 Participants. *J Am Coll Cardiol* 2017;70(5):558-568.
- 71 deFilippi CR, de Lemos JA, Christenson RH, et al. Association of serial measures of cardiac troponin T using a sensitive assay with incident heart failure and cardiovascular mortality in older adults. *JAMA* 2010;304(22):2494-502.
- 72 Saunders JT, Nambi V, de Lemos JA, et al. Cardiac troponin T measured by a highly sensitive assay predicts coronary heart disease, heart failure, and mortality in the Atherosclerosis Risk in Communities Study. *Circulation* 2011;123 (13):1367-1376.
- 73 Oluleye OW, Folsom AR, Nambi V, Lutsey PL, Ballantyne CM. Troponin T, B-type natriuretic peptide, C-reactive protein, and cause-specific mortality. *Ann Epidemiol* 2013;23 (2):66-73.
- 74 Folsom AR, Nambi V, Bell EJ, et al. Troponin T, N-terminal pro-B-type natriuretic peptide, and incidence of stroke: the Atherosclerosis Risk in Communities Study. *Stroke* 2013;44 (4):961-967.
- 75 Eggers KM, Al-Shakarchi J, Berglund L, et al. High-sensitive cardiac troponin T and its relations to cardiovascular risk factors, morbidity, and mortality in elderly men. *Am Heart J* 2013;166(3):541-548.
- 76 Willeit P, Kaptoge S, Welsh P, et al. Natriuretic peptides and integrated risk assessment for cardiovascular disease: an individual-participant-data meta-analysis. *Lancet Diabetes Endocrinol* 2016;4(10):840-849.

- 77 Bosselmann H, Egstrup M, Rossing K, et al. Prognostic significance of cardiovascular biomarkers and renal dysfunction in outpatients with systolic heart failure: a long term follow-up study *Int J Cardiol* 2013;170(2):202-207.
- 78 Wallentin L, Lindhagen L, Årnström E, et al. Early invasive versus noninvasive treatment in patients with non-ST-elevation acute coronary syndrome (FRISC-II): 15 year follow-up of a prospective, randomised, multicentre study. *Lancet* 2016;388(10054):1903-1911.
- 79 Gillmann HJ, Meinders A, Grohennig A, et al. Perioperative levels and changes of high-sensitivity troponin T are associated with cardiovascular events in vascular surgery patients. *Crit Care Med* 2014; 42(6):1498-1506.
- 80 Botto F, Alonso-Coello P, Chan MT et al. Myocardial injury after noncardiac surgery: a large, international, prospective cohort study establishing diagnostic criteria, characteristics, predictors, and 30-day outcomes. *Anesthesiology* 2014; 120(3):564-78.
- 81 Puelacher C, Lurati Buse G, Seeberger D, et al. Perioperative myocardial injury after noncardiac surgery: incidence, mortality, and characterization. *Circulation*. 2018; 137(12):1221-1232.
- 82 Kristensen SD, Knuuti J, Saraste A, et al. 2014 ESC/ESA Guidelines on non-cardiac surgery: cardiovascular assessment and management: The Joint Task Force on non-cardiac surgery: cardiovascular assessment and management of the European Society of Cardiology (ESC) and the European Society of Anaesthesiology (ESA). *Eur Heart J* 2014; 35(35):2383-431.
- 83 Duceppe E, Parlow J, MacDonald P, et al. Canadian cardiovascular society guidelines on perioperative cardiac risk assessment and management for patients undergoing noncardiac surgery. *Canadian Journal of Cardiology* *Can J Cardiol*. 2017;33(1):17-32.
- 84 Gualandro DM, Yu PC, Caramelli B, et al. 3rd Guideline for perioperative cardiovascular evaluation of the Brazilian society of cardiology. *Arq Bras Cardiol* 2017;109:1-104.
- 85 Davis GK, Labugger R, Van Eyk JE, et al. Cardiac troponin T is not detected in western blots of diseased renal tissue. *Clin Chem* 2001;47(4):782-783.
- 86 Ricchiuti V, Voss EM, Ney A, et al. Cardiac Troponin T isoforms expressed in renal diseased skeletal muscle will not cause false positive results by the second generation cTnT assay by Boehringer Mannheim; *Clin Chem* 1998;44(9):1919-1924.
- 87 Hallermayer K, Klenner D, Vogel R. Use of recombinant human cardiac troponin T for standardization of third generation troponin T methods. *Scand J Clin Invest* 1999;59(Suppl 230):128-131.
- 88 Wu AHB, Christenson RH, Greene DN et al. Clinical laboratory practice recommendations for the use of cardiac troponin in acute coronary syndrome: expert opinion from the Academy of the American Association for Clinical Chemistry and the Task Force on Clinical Applications of Cardiac Bio-Markers of the International Federation of Clinical Chemistry and Laboratory Medicine. *Clin Chem* 2018; 64(4):645-655.
- 89 Giannitsis E, Mueller-Hennessen M, Zeller T et al. Gender-specific reference values for high-sensitivity cardiac troponin T and I in well-phenotyped healthy individuals and validity of high-sensitivity assay designation. *Clin Biochem* 2020; 78:18-24.
- 90 Saenger AK, Beyrau R, Braun S, et al. Multicenter analytical evaluation of a high-sensitivity troponin T assay. *Clin Chim Acta* 2011;412(9-10):748-754.
- 91 Rubini Giménez M, Twerenbold R, Boeddinghaus J, Nestelberger T et al. Clinical Effect of Sex-Specific Cutoff Values of High-Sensitivity Cardiac Troponin T in Suspected Myocardial Infarction. *JAMA Cardiol*. 2016;1(8):912-920.
- 92 Giannitsis E. Sex-specific troponin measures for diagnosis of acute coronary syndrome. *Heart* 2016; 102(2):91-92.
- 93 Mueller-Hennessen M, Lindahl B, Giannitsis E et al. Diagnostic and prognostic implications using age- and gender-specific cut-offs for high-sensitivity cardiac troponin T - Sub-analysis from the TRAPID-AMI study. *Int J Cardiol*. 2016;209:26-33.
- 94 Eggers KM, Jernberg T, Lindahl B. Prognostic importance of sex-specific cardiac troponin T 99(th) percentiles in suspected acute coronary syndrome. *Am J Med* 2016; 129 (8): 880.e1-880.e12.
- 95 Mueller-Hennessen M, Giannitsis E. Do we need to consider age and gender for accurate diagnosis of myocardial infarction?. *Diagnosis (Berl.)* 2016; 3(4):175-181.
- 96 Giannitsis E. Counterpoint: Potential Concerns Regarding the Use of Sex-Specific Cutpoints for High-Sensitivity Troponin Assays. *Clin Chem*. 2017; 63(1):264-266.
- 97 Twerenbold R, Boeddinghaus J, Nestelberger T, et al. Clinical use of high-sensitivity cardiac troponin in patients with suspected myocardial infarction. *J Am Coll Cardiol*. 2017;70(8):996-1012.
- 98 World Health Organization. Report of the Joint International Society and Federation of Cardiology/World Health Organization Task Force on Standardization of Clinical Nomenclature. Nomenclature and criteria for diagnosis of ischemic heart disease. *Circulation* 1979;59:607-609.
- 99 Müller-Bardorff M, Hallermayer K, Schröder A, et al. Improved troponin T ELISA specific for cardiac troponin T isoform: assay development and analytical validation. *Clin Chem* 1997;43(3):458-466.
- 100 Klein G, Kampmann M, Baum H, et al. Clinical performance of the new cardiac markers troponin T and CK-MB on the Elecsys 2010. A multicenter evaluation. *Wien Klin Wochenschr*. 1998;110(Suppl3):40-51.
- 101 Mendis S, Thygesen K, Kuulasmaa K, et al. World Health Organization definition of myocardial infarction: 2008-09 revision. *Int J Epidemiol* 2011;40(1):139-146.
- 102 Ferjani M, Droc G, Dreux S, et al. Circulating cardiac troponin T in myocardial contusion. *Chest* 1997;111(2):427-433.
- 103 Erbel C, Taskin R, Doesch A, et al. High-sensitive troponin T measurements early after heart transplantation predict short- and long-term survival. *Transpl Int* 2013;26(3):267-272.
- 104 Li SF, Zapata J, Tillem E. The prevalence of false-positive cardiac troponin I in ED patients with rhabdomyolysis. *Am J Emerg Med* 2005;23(7):860-863.
- 105 Zhang L, Wang GC, Ma L, et al. Cardiac involvement in adult polymyositis or dermatomyositis: a systematic review *Clin Cardiol* 2012;35(11):686-691.
- 106 Writing committee for the VISION study investigators, Devereaux JP, Biccari BM, Sigamani A, et al. Association of postoperative high-sensitivity troponin levels with myocardial injury and 30-day mortality among patients undergoing noncardiac surgery. *JAMA* 2017; 317(16):1642-1651.
- 107 Devereaux PJ, Szczeklik W. Myocardial injury after non-cardiac surgery: diagnosis and management. *Eur Heart J* 2020;41(32):3083-3091.
- 108 Biccari BM, Rodseth RN. The pathophysiology of peri-operative myocardial infarction. *Anaesthesia* 2010;65(7):733-741.
- 109 DeFilippis AP, Young R, Carrubba CJ et al. An analysis of calibration and discrimination among multiple cardiovascular risk scores in a modern multiethnic cohort *Ann Intern Med* 2015;162:266-275.
- 110 Andersson C, Enserro D, Larson MG et al. Implications of the US cholesterol guidelines on eligibility for statin therapy in the community: comparison of observed and predicted risks in the Framingham Heart Study Offspring Cohort. *J Am Heart Assoc* 2015;4(4):e001888.
- 111 Wolfson J, Vock DM, Bandyopadhyay S et al. Use and customization of risk scores for predicting cardiovascular events using electronic health record data. *J Am Heart Assoc* 2017;6(4):e003670.
- 112 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. *J Clin Chem Clin Biochem* 1988 Nov;26(11):783-790.
- 113 The Task Force for the diagnosis and Treatment of Non-ST-Segment Elevation Acute Coronary Syndromes of the European Society of Cardiology. Guidelines for the diagnosis and treatment of non-ST elevation acute coronary syndromes. *European Heart Journal* 2007;28:1598-1660.

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© 2023, Roche Diagnostics
Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim