

Alkaline Phosphatase acc. to IFCC Gen.2

Order information

| REF | | CONTENT | | Analyzer(s) on which cobas c pack(s) can be used |
|-------------|-------------|---|---------------------|--|
| 03333752190 | 03333752500 | Alkaline Phosphatase acc. to IFCC Gen.2 ALP2S (200 tests) | System-ID 07 6761 1 | cobas c 311 , cobas c 501/502 , COBAS INTEGRA 400 plus |
| 03333701190 | 03333752500 | Alkaline Phosphatase acc. to IFCC Gen.2 ALP2L (400 tests) | System-ID 07 6760 3 | cobas c 311 , cobas c 501/502 , COBAS INTEGRA 400 plus |

Materials required (but not provided):

| | | cobas c 311 , cobas c 501/502 | COBAS INTEGRA 400 plus |
|-------------|---|---|------------------------|
| 10759350190 | Calibrator f.a.s. (12 x 3 mL) | Code 401 | System-ID 07 3718 6 |
| 12149435122 | Precinorm U plus (10 x 3 mL) | Code 300 | System-ID 07 7999 7 |
| 12149443122 | Precipath U plus (10 x 3 mL) | Code 301 | System-ID 07 8000 6 |
| 05117003190 | PreciControl ClinChem Multi 1 (20 x 5 mL) | Code 391 | System-ID 07 7469 3 |
| 05947626190 | PreciControl ClinChem Multi 1 (4 x 5 mL) | Code 391 | System-ID 07 7469 3 |
| 05117216190 | PreciControl ClinChem Multi 2 (20 x 5 mL) | Code 392 | System-ID 07 7470 7 |
| 05947774190 | PreciControl ClinChem Multi 2 (4 x 5 mL) | Code 392 | System-ID 07 7470 7 |
| 04489357190 | Diluent NaCl 9 % (50 mL) | System-ID 07 6869 3 | n.a. |

English

Intended use

In vitro test for the quantitative determination of alkaline phosphatase in human serum and plasma on **cobas c** and COBAS INTEGRA systems.

Summary

Measurement of alkaline phosphatase with this assay in human serum and plasma is used to aid in the diagnosis and monitoring of liver diseases and bone diseases.

Alkaline phosphatases (EC 3.1.3.1) are membrane-bound ectoenzymes that catalyze the hydrolysis of monophosphates from ester linkage under alkaline conditions (pH 8 to 10).¹ Alkaline phosphatase isoforms are encoded by four different genes: the liver-bone-kidney (tissue-nonspecific) variant, the intestinal variant, the placental variant and the variant from the germ cells (placental-like).^{1,2} Alkaline phosphatase activity is present in various tissues, but its concentration varies, and the highest concentrations are typically found in the liver and bone. Although the exact metabolic function of the enzyme is not yet understood, it appears that it is associated with lipid transport in the intestine, with the calcification process in bone, and with host defense through endotoxin dephosphorylation. Minimal amounts of intestinal alkaline phosphatase may also be present and are subjected to increase after a meal.²

Total serum alkaline phosphatase measurement is used extensively as a clinical indicator of liver and bone health.^{1,2,3,4,5,6,7,8,9} Any form of biliary tree obstruction induces the synthesis of alkaline phosphatase by hepatocytes, therefore a rise in the alkaline phosphatase activity in serum occurs with all forms of cholestasis and particularly with obstructive jaundice.^{2,3,4,5} It is also elevated in diseases of the skeletal system associated with increased osteoblastic activity, such as Paget's disease, hyperparathyroidism, rickets and osteomalacia, as well as with fractures and malignant tumors.^{1,6,7,8,9,10} A physiologic rise in the alkaline phosphatase activity is sometimes seen in children and juveniles. It is caused by increased osteoblast activity following accelerated bone growth.^{1,2,10}


Decreased total alkaline phosphatase activity is rarely found in human serum but can occur in hypophosphatasia, in multiple myeloma with osteolytic lesions, secondary to growth hormone deficiency or in hypoparathyroidism.^{1,10}

The assay method was first described by King and Armstrong, modified by Ohmori, Bessey, Lowry and Brock and later improved by Hausamen et al.^{11,12,13,14} In 2011 the International Federation of Clinical Chemistry and Laboratory Medicine (IFCC) Scientific Division, Committee on Reference Systems of Enzymes (C-RSE) recommended a reference procedure for the determination of alkaline phosphatase using an optimized substrate concentration and 2-amino-2-methyl-1-propanol as buffer plus the cations magnesium and zinc at 37 °C.¹⁵ This assay follows the recommendations of the IFCC, but was optimized for performance and stability.

Test principle¹⁵

Colorimetric assay in accordance with a standardized method.

In the presence of magnesium and zinc ions, p-nitrophenyl phosphate is cleaved by phosphatases into phosphate and p-nitrophenol.



The p-nitrophenol released is directly proportional to the catalytic ALP activity. It is determined by measuring the increase in absorbance.

Precautions and warnings

For in vitro diagnostic use for health care professionals. Exercise the normal precautions required for handling all laboratory reagents.

Infectious or microbial waste:

Warning: handle waste as potentially biohazardous material. Dispose of waste according to accepted laboratory instructions and procedures.

Environmental hazards:

Apply all relevant local disposal regulations to determine the safe disposal.

Safety data sheet available for professional user on request.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:



Warning

H315 Causes skin irritation.

H319 Causes serious eye irritation.

Prevention:

P264 Wash skin thoroughly after handling.

P280 Wear protective gloves/ eye protection/ face protection.

Response:

P302 + P352 IF ON SKIN: Wash with plenty of water.

P332 + P313 If skin irritation occurs: Get medical advice/attention.

P337 + P313 If eye irritation persists: Get medical advice/attention.

P362 + P364 Take off contaminated clothing and wash it before reuse.

Product safety labeling follows EU GHS guidance.

Contact phone: all countries: +49-621-7590

Reagent handling

Ready for use

Specimen collection and preparation

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.
Serum.

Plasma: Li-heparin plasma.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

See the limitations and interferences section for details about possible sample interferences.

| | |
|--------------------------|-----------------------------|
| Stability: ¹⁶ | 7 days at 20-25 °C |
| | 7 days at 4-8 °C |
| | 2 months at -20 °C (± 5 °C) |

Freeze only once.

Materials provided

See "Reagents – working solutions" section for reagents.

Materials required (but not provided)

- See "Order information" section
- General laboratory equipment

Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

The performance of applications not validated by Roche is not warranted and must be defined by the user.

Calculation

The systems automatically calculate the analyte concentration of each sample.

Conversion factor: U/L x 0.0167 = µkat/L

Expected values

(measured at 37 °C)

Adults¹⁷

| | | |
|-------------------|------------|--------------------|
| Males (n = 221) | 40-129 U/L | (0.67-2.15 µkat/L) |
| Females (n = 229) | 35-104 U/L | (0.58-1.74 µkat/L) |

Children¹⁸

Males

Age

| | | |
|--------------------|-------------|--------------------|
| 0 – 14 days | 83-248 U/L | (1.39-4.14 µkat/L) |
| 15 days – < 1 year | 122-469 U/L | (2.04-7.83 µkat/L) |
| 1 – < 10 years | 142-335 U/L | (2.37-5.59 µkat/L) |
| 10 – < 13 years | 129-417 U/L | (2.15-6.96 µkat/L) |
| 13 – < 15 years | 116-468 U/L | (1.94-7.82 µkat/L) |
| 15 – < 17 years | 82-331 U/L | (1.37-5.53 µkat/L) |
| 17 – < 19 years | 55-149 U/L | (0.92-2.49 µkat/L) |

Females

Age

| | | |
|--------------------|-------------|--------------------|
| 0 – 14 days | 83-248 U/L | (1.39-4.14 µkat/L) |
| 15 days – < 1 year | 122-469 U/L | (2.04-7.83 µkat/L) |
| 1 – < 10 years | 142-335 U/L | (2.37-5.59 µkat/L) |
| 10 – < 13 years | 129-417 U/L | (2.15-6.96 µkat/L) |
| 13 – < 15 years | 57-254 U/L | (0.95-4.24 µkat/L) |
| 15 – < 17 years | 50-117 U/L | (0.84-1.95 µkat/L) |
| 17 – < 19 years | 45-87 U/L | (0.75-1.45 µkat/L) |

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

cobas c systems

System information

For **cobas c 311/501** analyzers:

ALP2S: ACN 158

ALP2L: ACN 683

For **cobas c 502** analyzer:

ALP2S: ACN 8158

ALP2L: ACN 8683

Reagents - working solutions

R1 2-amino-2-methyl-1-propanol: 1.724 mol/L, pH 10.44 (30 °C);
magnesium acetate: 3.83 mmol/L; zinc sulfate: 0.766 mmol/L;
N-(2-hydroxyethyl)-ethylenediamine triacetic acid: 3.83 mmol/L

R2 p-nitrophenyl phosphate: 132.8 mmol/L, pH 8.50 (25 °C);
preservatives

R1 is in position B and R2 is in position C.

Storage and stability

Shelf life at 2-8 °C: See expiration date on **cobas c** pack label.

On-board in use and refrigerated on the analyzer: 8 weeks

Applications for serum and plasma

cobas c 311 test definition

| | | | |
|------------------------------|----------------------------|-----------------|----------------|
| Assay type | Rate A | | |
| Reaction time / Assay points | 10 / 13-31 | | |
| Wavelength (sub/main) | 480/450 nm | | |
| Reaction direction | Increase | | |
| Units | U/L (µkat/L) | | |
| Reagent pipetting | Diluent (H ₂ O) | | |
| R1 | 75 µL | 25 µL | |
| R2 | 17 µL | 21 µL | |
| Sample volumes | Sample | Sample dilution | |
| | | Sample | Diluent (NaCl) |
| Normal | 2.8 µL | – | – |
| Decreased | 2.8 µL | 20 µL | 80 µL |
| Increased | 2.8 µL | – | – |

cobas c 501 test definition

| | |
|------------------------------|------------|
| Assay type | Rate A |
| Reaction time / Assay points | 10 / 19-48 |
| Wavelength (sub/main) | 480/450 nm |

| | | | |
|--------------------|----------------------------|-----------------|----------------|
| Reaction direction | Increase | | |
| Units | U/L (µkat/L) | | |
| Reagent pipetting | Diluent (H ₂ O) | | |
| R1 | 75 µL | 25 µL | |
| R2 | 17 µL | 21 µL | |
| Sample volumes | Sample | Sample dilution | |
| | | Sample | Diluent (NaCl) |
| Normal | 2.8 µL | – | – |
| Decreased | 2.8 µL | 20 µL | 80 µL |
| Increased | 2.8 µL | – | – |

cobas c 502 test definition

| | | | |
|------------------------------|----------------------------|-----------------|----------------|
| Assay type | Rate A | | |
| Reaction time / Assay points | 10 / 19-48 | | |
| Wavelength (sub/main) | 480/450 nm | | |
| Reaction direction | Increase | | |
| Units | U/L (µkat/L) | | |
| Reagent pipetting | Diluent (H ₂ O) | | |
| R1 | 75 µL | 25 µL | |
| R2 | 17 µL | 21 µL | |
| Sample volumes | Sample | Sample dilution | |
| | | Sample | Diluent (NaCl) |
| Normal | 2.8 µL | – | – |
| Decreased | 2.8 µL | 20 µL | 80 µL |
| Increased | 5.6 µL | – | – |

Calibration

| | |
|-----------------------|--|
| Calibrators | S1: H ₂ O |
| | S2: C.f.a.s. |
| Calibration mode | Linear |
| Calibration frequency | 2-point calibration |
| | - after reagent lot change |
| | - as required following quality control procedures |

Calibration interval may be extended based on acceptable verification of calibration by the laboratory.

Traceability: This method has been standardized against the IFCC procedure (2011).¹⁵

Quality control

For quality control, use control materials as listed in the "Order information" section.

In addition, other suitable control material can be used.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

Limitations - interference

Criterion: Recovery within ± 10 U/L of initial values of samples ≤ 100 U/L and within ± 10 % for samples > 100 U/L.

Icterus:¹⁹ No significant interference up to an I index of 60 for conjugated and unconjugated bilirubin (approximate conjugated and unconjugated bilirubin concentration: 1026 µmol/L or 60 mg/dL).

Hemolysis:¹⁹ No significant interference up to an H index of 200 (approximate hemoglobin concentration: 124 µmol/L or 200 mg/dL).

Lipemia (Intralipid):¹⁹ No significant interference up to an L index of 2000. There is poor correlation between the L index (corresponds to turbidity) and triglycerides concentration.

Drugs: No interference was found at therapeutic concentrations using common drug panels.^{20,21}

In very rare cases, gammopathy, in particular type IgM (Waldenström's macroglobulinemia), may cause unreliable results.²²

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

ACTION REQUIRED

Special Wash Programming: The use of special wash steps is mandatory when certain test combinations are run together on **cobas c** systems. The latest version of the carry-over evasion list can be found with the NaOHD-SMS-SmpCln1+2-SCCS Method Sheets. For further instructions refer to the operator's manual. **cobas c 502** analyzer: All special wash programming necessary for avoiding carry-over is available via the **cobas** link, manual input is required in certain cases.

Where required, special wash/carry-over evasion programming must be implemented prior to reporting results with this test.

Limits and ranges**Measuring range**

5-1200 U/L (0.084-20.0 µkat/L)

Determine samples having higher activities via the rerun function. Dilution of samples via the rerun function is a 1:5 dilution. Results from samples diluted using the rerun function are automatically multiplied by a factor of 5.

Lower limits of measurement

Lower detection limit of the test

5 U/L (0.084 µkat/L)

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying 3 standard deviations above that of the lowest standard (standard 1 + 3 SD, repeatability, n = 21).

Specific performance data

Representative performance data on the analyzers are given below. Results obtained in individual laboratories may differ.

Precision

Precision was determined using human samples and controls in an internal protocol with repeatability (n = 21) and intermediate precision (3 aliquots per run, 1 run per day, 21 days). The following results were obtained on the **cobas c 501** analyzer:

| Repeatability | Mean | SD | CV |
|------------------------|--------------|--------------|-----|
| | U/L (µkat/L) | U/L (µkat/L) | % |
| Precinorm U | 99.2 (1.65) | 0.7 (0.01) | 0.7 |
| Precipath U | 241 (4.02) | 1 (0.02) | 0.6 |
| Human serum 1 | 54.6 (0.912) | 0.5 (0.008) | 0.9 |
| Human serum 2 | 648 (10.8) | 4 (0.1) | 0.7 |
| Intermediate precision | Mean | SD | CV |
| | U/L (µkat/L) | U/L (µkat/L) | % |
| Precinorm U | 92.8 (1.56) | 2.2 (0.04) | 2.4 |
| Precipath U | 224 (3.74) | 4 (0.06) | 1.7 |
| Human serum 3 | 82.2 (1.37) | 1.8 (0.03) | 2.1 |
| Human serum 4 | 1025 (17.1) | 9 (0.2) | 0.9 |

The data obtained on **cobas c 501** analyzer(s) are representative for **cobas c 311** analyzer(s).

Method comparison

Alkaline phosphatase values for human serum and plasma samples obtained on a **cobas c 501** analyzer with the ALP IFCC Gen.2 (ALP2) traceable to IFCC¹⁵ method (y), were compared with those determined on

the same analyzer with the same ALP2 reagent traceable to IFCC²³ method (x).

Sample size (n) = 106

Passing/Bablok²⁴

Linear regression

$y = 1.05x + 0.064 \text{ U/L}$

$y = 1.04x + 0.388 \text{ U/L}$

$\tau = 0.993$

$r = 1.00$

The sample activities were between 16.9 and 1149 U/L (0.282 and 19.2 $\mu\text{kat/L}$).

The data obtained on **cobas c 501** analyzer(s) are representative for **cobas c 311** analyzer(s).

COBAS INTEGRA systems

System information

Serum/plasma

cobas c pack ALP2S, Cat. No. 03333752190:

ALP2S: Test-ID 0-551

cobas c pack ALP2L, Cat. No. 03333701190:

ALP2L: Test-ID 0-550

*Serum/plasma-Primary tube**

cobas c pack ALP2S, Cat. No. 03333752190:

AP2SP: Test-ID 0-554

cobas c pack ALP2L, Cat. No. 03333701190:

AP2LP: Test-ID 0-553

*The applications is intended for customers facing non-valid results due to a contamination of the plasma supernatant in primary tubes with cell aggregates.

Reagents - working solutions

R1 2-amino-2-methyl-1-propanol: 1.724 mol/L, pH 10.44 (30 °C);
magnesium acetate: 3.83 mmol/L; zinc sulfate: 0.766 mmol/L;
N-(2-hydroxyethyl)-ethylenediamine triacetic acid: 3.83 mmol/L

SR p-nitrophenyl phosphate: 132.8 mmol/L, pH 8.50 (25 °C);
preservatives

R1 is in position B and SR is in position C.

Storage and stability

Shelf life at 2-8 °C

See expiration date on
cobas c pack label

On-board in use at 10-15 °C

4 weeks

Application for serum and plasma

Serum/plasma

Test definition

| | |
|-----------------------|------------|
| Measuring mode | Absorbance |
| Abs. calculation mode | Kinetic |
| Reaction mode | R1-S-SR |
| Reaction direction | Increase |
| Wavelength A/B | 409/659 nm |
| Calc. first/last | 41/64 |
| Unit | U/L |

Pipetting parameters

| | | Diluent (H ₂ O) |
|--------------|----------------------|----------------------------|
| R1 | 75 μL | 16 μL |
| Sample | 2.75 μL | 20 μL |
| SR | 17 μL | 10 μL |
| Total volume | 140.75 μL | |

Serum/plasma-Primary tube

Test definition

| | |
|-----------------------|------------|
| Measuring mode | Absorbance |
| Abs. calculation mode | Kinetic |
| Reaction mode | D-R1-S-SR |
| Reaction direction | Increase |
| Wavelength A/B | 409/659 nm |
| Calc. first/last | 41/64 |
| Predilution factor | 10 |
| Unit | U/L |

Pipetting parameters

| | | Diluent (H ₂ O) |
|--------------|---------------------|----------------------------|
| R1 | 75 μL | 11 μL |
| Sample | 27.5 μL | |
| SR | 17 μL | 10 μL |
| Total volume | 140.5 μL | |

Calibration

| | |
|-----------------------|---|
| Calibrator | Calibrator f.a.s. Use deionized water as zero calibrator. |
| Calibration mode | Linear regression |
| Calibration replicate | Duplicate recommended |
| Calibration frequency | 2 point calibration - after reagent lot change - as required following quality control procedures |

Calibration interval may be extended based on acceptable verification of calibration by the laboratory.

Traceability: This method has been standardized against the IFCC procedure (2011).¹⁵

Quality control

| | |
|---------------------------|--|
| Reference range | Precinorm U plus or PreciControl ClinChem Multi 1 |
| Pathological range | Precipath U plus or PreciControl ClinChem Multi 2 |
| Control interval | 24 hours recommended |
| Control sequence | User defined |
| Control after calibration | Recommended |

For quality control, use control materials as listed in the "Order information" section. In addition, other suitable control material can be used.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

Limitations - interference

Criterion: Recovery within $\pm 10 \text{ U/L}$ of initial values of samples $\leq 100 \text{ U/L}$ and within $\pm 10 \%$ for samples $> 100 \text{ U/L}$.

Icterus:¹⁹ No significant interference up to an I index of 42 for conjugated bilirubin and 60 for unconjugated bilirubin (approximate conjugated bilirubin concentration: 718 $\mu\text{mol/L}$ or 42 mg/dL; approximate unconjugated bilirubin concentration: 1026 $\mu\text{mol/L}$ or 60 mg/dL).

Hemolysis:¹⁹ No significant interference up to an H index of 250 (approximate hemoglobin concentration: 155 $\mu\text{mol/L}$ or 250 mg/dL).

Lipemia (Intralipid):¹⁹ No significant interference up to an L index of 2000. There is poor correlation between the L index (corresponds to turbidity) and triglycerides concentration.

Drugs: No interference was found at therapeutic concentrations using common drug panels.^{20,21}

In very rare cases, gammopathy, in particular type IgM (Waldenström's macroglobulinemia), may cause unreliable results.²²

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

ACTION REQUIRED

Special Wash Programming: The use of special wash steps is mandatory when certain test combinations are run together on COBAS INTEGRA analyzers. Refer to the CLEAN Method Sheet for further instructions and for the latest version of the Extra wash cycle list.

Where required, special wash/carry-over evasion programming must be implemented prior to reporting results with this test.

Limits and ranges

Measuring range

3.0-1200 U/L (0.05-20 µkat/L)

Determine samples having higher activities via the rerun function. Dilution of samples via the rerun function is a 1:5 dilution. Results from samples diluted using the rerun function are automatically multiplied by a factor of 5.

Lower limits of measurement

Lower detection limit of the test:

3.0 U/L (0.05 µkat/L)

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying 3 standard deviations above that of a zero sample (zero sample + 3 SD, repeatability, n = 21).

Specific performance data

Representative performance data on the analyzers are given below. Results obtained in individual laboratories may differ.

Precision

Serum/plasma

Precision was determined using human samples and controls in an internal protocol with repeatability (n = 21) and intermediate precision (1 aliquot per run, 1 run per day, 21 days). The following results were obtained on the COBAS INTEGRA 700 analyzer:

| Repeatability | Mean U/L (µkat/L) | SD U/L (µkat/L) | CV % |
|---------------|----------------------|--------------------|---------|
| Precinorm U | 80.1 (1.34) | 1.3 (0.02) | 1.6 |
| Pecipath U | 228 (3.81) | 4 (0.07) | 1.8 |
| Human serum 1 | 72.7 (1.21) | 1.5 (0.03) | 2.0 |
| Human serum 2 | 225 (3.76) | 4 (0.07) | 1.8 |

| Intermediate precision | Mean U/L (µkat/L) | SD U/L (µkat/L) | CV % |
|------------------------|----------------------|--------------------|---------|
| Precinorm U | 81.8 (1.37) | 2.3 (0.04) | 2.8 |
| Pecipath U | 230 (3.84) | 6 (0.10) | 2.8 |
| Human serum 1 | 70.0 (1.17) | 1.9 (0.03) | 2.7 |
| Human serum 2 | 220 (3.67) | 6 (0.10) | 2.7 |

The data obtained on COBAS INTEGRA 700 analyzer(s) are representative for COBAS INTEGRA 400 analyzer(s).

Serum/plasma-Primary tube

Precision was determined using human samples and controls in an internal protocol with repeatability (n = 21) and intermediate precision (1 aliquot per run, 1 run per day, 10 days). The following results were obtained the COBAS INTEGRA 700 analyzer:

| Repeatability | Mean U/L (µkat/L) | SD U/L (µkat/L) | CV % |
|---------------|----------------------|--------------------|---------|
| Precinorm U | 83.7 (1.40) | 0.6 (0.01) | 0.7 |
| Pecipath U | 226 (3.77) | 1 (0.02) | 0.4 |
| Human serum 1 | 84.4 (1.41) | 0.3 (0.01) | 0.4 |
| Human serum 2 | 217 (3.62) | 1 (0.02) | 0.4 |

| Intermediate precision | Mean U/L (µkat/L) | SD U/L (µkat/L) | CV % |
|------------------------|----------------------|--------------------|---------|
| Precinorm U | 84.0 (1.40) | 1.8 (0.03) | 2.1 |
| Pecipath U | 227 (3.80) | 4 (0.07) | 1.7 |
| Human serum 1 | 85.7 (1.43) | 1.7 (0.03) | 2.0 |
| Human serum 2 | 218 (3.64) | 5 (0.08) | 2.2 |

The data obtained on COBAS INTEGRA 700 analyzer(s) are representative for COBAS INTEGRA 400 analyzer(s).

Method comparison

Serum/plasma

ALP values for human serum patient samples obtained on a COBAS INTEGRA 400 plus analyzer with the COBAS INTEGRA ALP IFCC Gen.2 (ALP2) traceable to IFCC¹⁵ method (y), were compared with those determined on the same analyzer with the same ALP2 reagent traceable to IFCC²³ method (x).

COBAS INTEGRA 400 plus analyzer

Sample size (n) = 104

Passing/Bablok²⁴

Linear regression

$y = 1.04x - 0.078$ U/L

$y = 1.04x - 0.039$ U/L

$\tau = 0.997$

$r = 1.00$

The sample activities were between 15.0 and 1036 U/L (0.251 and 17.3 µkat/L).

Serum/plasma -Primary tube

ALP values for human serum and plasma samples obtained on a COBAS INTEGRA 700 analyzer with the COBAS INTEGRA ALP IFCC Gen.2 (ALP2L) reagent and the AP2LP application (y) were compared to those determined using the same reagent but the ALP2L application on a COBAS INTEGRA 700 analyzer (x), and to those determined using the previous reagent using the ALP6P application on a COBAS INTEGRA 700 analyzer (x).

COBAS INTEGRA (ALP2L)

Sample size (n) = 54

Passing/Bablok²⁴

Linear regression

$y = 1.003x - 0.195$ U/L

$y = 0.999x + 2.11$ U/L

$\tau = 0.976$

$r = 1.00$

SD (md 95) = 9.19

Sy.x = 4.06

The sample activities were between 40 and 1010 U/L (0.668 and 16.9 µkat/L).

COBAS INTEGRA (ALP6P)

Sample size (n) = 67

Passing/Bablok²⁴

Linear regression

$y = 0.991x + 0.808$ U/L

$y = 0.997x - 1.35$ U/L

$\tau = 0.976$

$r = 1.00$

SD (md 95) = 7.16

Sy.x = 3.39

The sample activities were between 39 and 862 U/L (0.651 and 15.8 µkat/L).

The data obtained on COBAS INTEGRA 700 analyzer(s) are representative for COBAS INTEGRA 400 analyzer(s).

References

- Fraser WD, Alter DN. Bone and mineral metabolism. In: Rifai N, Chiu RWK, Young I, Burnham CAD, Wittwer CT, editors. Tietz Textbook of Laboratory Medicine, Saunders Elsevier, Philadelphia, 7th edition, 2023, chapter 54, p. 766-766.e85.

- 2 Panteghini M. Serum Enzymes. In: Rifai N, Chiu RWK, Young I, Burnham CAD, Wittwer CT, editors. Tietz Textbook of Laboratory Medicine, Saunders Elsevier, Philadelphia, 7th edition, 2023, chapter 32, p. 350-350.e36.
- 3 Kwo PY, Cohen SM, Lim JK. ACG Clinical Guideline: Evaluation of Abnormal Liver Chemistries. Am J Gastroenterol 2017 Jan;112(1):18-35. doi: 10.1038/ajg.2016.517.
- 4 European Association for the Study of the Liver. EASL Clinical Practice Guidelines: management of cholestatic liver diseases. J Hepatol 2009 Aug;51(2):237-267.
- 5 Newsome PN, Cramb R, Davison SM, et al. Guidelines on the management of abnormal liver blood tests. Gut 2018 Jan;67(1):6-19. doi: 10.1136/gutjnl-2017-314924.
- 6 Singer FR, Bone HG 3rd, Hosking DJ, et al. Endocrine Society. Paget's disease of bone: an endocrine society clinical practice guideline. J Clin Endocrinol Metab 2014 Dec;99(12):4408-4422.
- 7 Vlot MC, den Heijer M, de Jongh RT, et al. Clinical utility of bone markers in various diseases. Bone 2018 Sep;114:215-225.
- 8 Tan A, Goodman K, Walker A, et al. PRISM-EZ Trial Group. Long-Term Randomized Trial of Intensive Versus Symptomatic Management in Paget's Disease of Bone: The PRISM-EZ Study. J Bone Miner Res 2017 Jun;32(6):1165-1173.
- 9 Thacher TD, Smith L, Fischer PR, et al. Optimal Dose of Calcium for Treatment of Nutritional Rickets: A Randomized Controlled Trial. J Bone Miner Res 2016 Nov;31(11):2024-2031.
- 10 Makris K, Mousa C, Cavalier E. Alkaline Phosphatases: Biochemistry, Functions, and Measurement. Calcif Tissue Int 2023 Feb;112(2):233-242.
- 11 King EJ, Armstrong AR. A convenient method for determining serum and bile phosphatase activity. Can Med Assoc J 1934 Oct;31(4):376-381.
- 12 Ohmori Y. Über die phosphomonoesterase. Enzymologia 1937;4:217-231.
- 13 Bessey OA, Lowry OH, Brock MJ. A method for the rapid determination of alkaline phosphatase with five cubic millimeters of serum. J Biol Chem 1946;164:321-329.
- 14 Hausamen TU, Helger R, Rick W, et al. Optimal conditions for the determination of serum alkaline phosphatase by a new kinetic method. Clin Chim Acta 1967;15:241-245.
- 15 Schumann G, Klauke R, Canalias F, et al. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37° C. - Part 9. Reference procedure for the measurement of catalytic concentration of alkaline phosphatase. Clin Chem Lab Med 2011 Sep;49 (9):1439-1446.
- 16 Guder WG, Narayanan S, Wisser H, et al. List of Analytes; Preanalytical Variables. Brochure in: Samples: From the Patient to the Laboratory. Darmstadt: GIT-Verlag 1996.
- 17 Abicht K, El-Samalouti V, Junge W, et al. Multicenter evaluation of new GGT and ALP reagents with new reference standardization and determination of 37 °C reference intervals. Clin Chem Lab Med 2001;39:Special Supplement pp S 346.
- 18 Estey MP, Cohen AH, Colantonio DA, et al. CLSI-based transference of the CALIPER database of pediatric reference intervals from Abbott to Beckman, Ortho, Roche and Siemens Clinical Chemistry Assays: Direct validation using reference samples from the CALIPER cohort. Clin Biochem 2013;46:1197-1219.
- 19 Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-475.
- 20 Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". Eur J Clin Chem Clin Biochem 1996;34:385-386.
- 21 Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.
- 22 Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007;45(9):1240-1243.

- 23 Tietz NW, Rinker AD, Shaw LM. International Federation of Clinical Chemistry. IFCC methods for the measurement of catalytic concentration of enzymes, Part 5. IFCC method for alkaline phosphatase (orthophosphoric-monoester phosphohydrolase, alkaline optimum, EC 3.1.3.1). J Clin Chem Clin Biochem 1983;21:731-748.
- 24 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard:

CONTENT

Contents of kit



Volume for reconstitution

GTIN

Global Trade Item Number

Rx only

For USA: Caution: Federal law restricts this device to sale by or on the order of a physician.

COBAS, NAVIFY, PRECICONTROL, PRECINORM and PRECIPATH are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2024, Roche Diagnostics

CE 0123



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com

+800 5505 6606

