# Combur<sup>5</sup> Test HC

**cobas**<sup>®</sup>

10

REF 11896954176

## English

### Intended use

Five-parameter test strips for graduated determinations of glucose, leukocytes, nitrite, protein and blood in urine.

Suitable for self-testing.

### Summary

Combur<sup>5</sup> Test HC provides a fast and easy means of early detection and clinical monitoring of diabetes mellitus and diseases of the kidneys and urinary tract.

Combur<sup>5</sup> Test HC supports your physician's efforts to detect these diseases in their early stages. It does this by indicating the presence of substances excreted in the urine, thereby providing evidence of disease even when no symptoms are present.

- Test your urine for 2 to 3 days in succession.

- Self-testing cannot, however, be used in place of testing performed by a physician.

- Continue testing every 3 months using two to three test strips, even if your test results are normal or your physician does not detect any abnormalities. This will provide additional security.

- Please contact your physician if you detect + or ++ results on one or more test parameters. Please test your urine again if the test results are unclear. In case all parameters are 0 after the first measuring, please test your urine for 1 or 2 more days and show the results to your physician.

- There is a possibility that test results might be false-positives or false-negatives. In case of doubt please test your urine again.

- Always consult your physician before making any decision that could affect your treatment.

### Precautions and warnings

The stopper of the test strip vial contains a non-toxic silicate-based desiccant which must not be removed. If ingested by accident, drink large quantities of water.

For in vitro diagnostic use.

Disposal of all waste material should be in accordance with local guidelines.

Safety data sheet available for professional user on request.

### Storage and stability

Store the package at 2-30 °C. When stored in the original container the test strips are stable up to the expiration date specified on the box and color label.

Tightly re-cap the container immediately after removing a test strip.

Do not use test strips if they are past the expiration date printed on the container and color label.

Discard expired test strip.

### Specimen collection and preparation

Use only clean, well-rinsed vessels to collect urine.

Do not add preservatives to the urine.

Use fresh urine that has not been centrifuged.<sup>1</sup> The urine specimen should not stand for more than 2 hours before testing. For specimen collection and preparation only use suitable tubes or collection containers, as false-positive readings, particularly for glucose and protein, can result from residues of detergent or strongly oxidizing disinfectants in the specimen collection vessel.<sup>1</sup> Using midstream urine is recommended to avoid contamination by commensal urethral flora in both sexes.<sup>1</sup> Do not expose urine specimens to sunlight as this induces oxidation of bilirubin and urobilinogen and hence leads to artificially low results for these two parameters.<sup>1</sup> Vaginal secretion or menstrual blood may contaminate urine from females.<sup>1</sup>

Diagnosis or therapy should never be based on one test result alone but should be established in the context of all other medical findings. In doubtful cases, it is therefore advisable to repeat the test after discontinuation of the medication.

### Materials provided

10 test strips

### Materials required (but not provided)

- A clean cup for collecting urine

### Directions for use

For optimum performance of the test follow the directions given below:

- When measuring urinary glucose, note that, in some cases, increased glucose values can only be detected after a meal. As such, performing the test 2 hours after a meal, preferably after breakfast, is recommended. If possible, this meal should include 3-4 slices of bread or 1-2 rolls.

- Allowing the urine to remain in the bladder for long periods of time (4-8 hours, ideally overnight) is recommended as a means of ensuring a high level of accuracy.

- In women, this test can yield false results for blood if carried out 3 days before or 3 days after menses; as such, women are advised against performing this test during this period.

### Before you test:

- Always test your first morning urine. Wash your genitals carefully and rinse thoroughly with water (do not dry). To collect urine use only clean, well-rinsed vessels.

- Allow some urine to flow into the toilet. Collect the rest of the stream in the specimen cup. Performing the test:

- Remove test strip. Immediately after you take a test strip out of the vial, seal the vial with the original desiccant stopper. Failure to do this could result in false results, as moisture in the vial could potentially cause color changes in the test parameters.

- Briefly dip the test strip into the urine (approximately 1 sec.), making certain to wet all of the test parameters (fig. 1).

- When removing the test strip, wipe its edge along the rim of the specimen cup in order to remove any excess urine (fig. 2).

- All test parameters should be read at 1 minute (up to 120 seconds for the leukocyte test area for not clearly assignable results). Compare the colors of the test parameters on the strip with the colors on the label (fig. 3).

- Separate sets of color blocks are given for erythrocytes and hemoglobin. Compare the 5th (blood) test parameter area with both rows of colors.


Fig. 1


Fig. 2


Fig. 4

### Result sheet:

- If the color changes in one or more test parameters, place an x in the box next to the line on the results sheet (0, + or ++) that best corresponds to the color of the corresponding test parameter (fig. 4). If a test parameter appears darker than ++, check the box under ++.

- Enter date and your name into results sheet and take it with you to your physician.

- Any color changes appearing only along the edges of the test areas, or developing after more than 2 minutes, do not have any diagnostic significance.

## Additional information for Healthcare Professionals

### Reagents

Each test contains per 1 cm<sup>2</sup> reactive paper area the following:  
**Leukocytes**: Indoxylcarboxylic acid ester 15.5 µg; methoxymorpholinobenzene diazonium salt 5.5 µg  
**Nitrite**: 3-hydroxy-1,2,3,4-tetrahydro-7,8-benzoquinoline 33.5 µg; sulfanilamide 29.1 µg  
**Protein**: 3',3'',5'-trichlorophenol-3,4,5,6-tetrabromosulfophthalimide 13.9 µg  
**Glucose**: 3,3',5,5'-tetramethylbenzidine 103.5 µg; GOD 6 U, POD 35 U  
**Blood**: 3,3',5,5'-tetramethylbenzidine 52.8 µg; 2,5-dimethyl-2,5-dihydroperoxyhexane 297.2 µg

### Test principle

Test parameters on each strip contain indicator reagents. As soon as the test parameters have been wetted with urine, a chemical reaction takes place that causes the color of the test parameters to change in proportion to the concentration of the substance in question. The color change can then be assessed by comparing it with the color scale given on the vial label.

**Glucose (GLU):** Glucose in the urine is almost always an indication of diabetes mellitus. Testing with Combur<sup>5</sup> Test HC is especially important for relatives of diabetics and for those at risk of diabetes (e.g., overweight individuals). Combur<sup>5</sup> Test HC detects substances that are excreted in the urine with a certain degree of variability. Because individual urine samples from diabetics could occasionally be glucose free, for instance, failure of the glucose parameter to change color does not rule out diabetes. Under certain circumstances (such as pregnancy) glucose may be present in the urine, even if blood glucose levels are normal.

**Leukocytes (LEU):** In case of an inflammation or infection of the kidneys or the urinary tract, more white blood cells are excreted than usual. Antibiotics containing imipenem, meropenem and clavulanic acid as active agents can result in false (positive) test results.

**Nitrite (NIT):** The most common cause of urinary tract infections is E. coli, which, along with most other pathogenic organisms of the urinary tract, converts the nitrate consumed in food into nitrite. The presence of nitrite is indicated by a pinkish red color on the test parameter, thereby providing indirect evidence of nitrite-forming bacteria in the urine. Even a faint pink color indicates a significant increase in the bacteria count.

**Protein (PRO):** Excreting protein in the urine can be an indication of kidney and urinary tract disease. Evidence of protein in the urine is not enough by itself to make an unambiguous diagnosis; the presence of protein could also have harmless causes (such as physical strain).

**Blood (ERY/Hb):** Kidney and urinary tract disease are the primary causes of blood in the urine. Hemoglobin, the component of blood that gives red blood cells (erythrocytes) their red color, may be released under certain conditions (e.g., physical strain, burns, ingestion of toxins, infectious processes). The vial label provides separate color scales for erythrocytes and hemoglobin. Individual or clustered green dots on the yellow test parameter are a sign of intact erythrocytes. Hemoglobin and/or damaged erythrocytes or myoglobin (a human protein similar to hemoglobin) are indicated by a uniform green color in the test parameter.

**Specific performance data**  
Representative performance data are given below. Results obtained in individual laboratories may differ. The values for neg. and pos. indicate the proportion of concordant negative or positive results.

The values specified for the limit of detection are defined as the concentration of the analyte which leads to a positive result in ≥ 90 % of the examined urines.

The method comparison data for visual reading are based on the comparison with the instrument cobas i 411 with Combur<sup>5</sup> Test M using at least 232 clinical samples per parameter. All test parameters were covered.

Parameter	Limit of detection	Method comparison <sup>a</sup>
LEU	36 LEU/µL	neg.: 100 % pos.: 62 %
NIT	0.05 mg/dL	neg.: 95 % pos.: 91 %
PRO	14 mg/dL (albumin)	neg.: 85 % pos.: 99 %
GLU	30 mg/dL	neg.: 96 % pos.: 100 %
ERY/Hb	Intact erythrocytes: 6 ERY/µL Hemoglobin: 6 ERY/µL	neg.: 99 % pos.: 94 %

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

**Limitations - interference**  
The following drugs and substances were tested with Combur-Test technology test strips in the latest interference study from November 2013.

Therapeutic drugs	Endogenous substances
Acetaminophen	Hydrochlorothiazide
N-Acetylcysteine	Hydroxylchloroquine
Amoxicillin	Creatinine
Amlodipine besylate	Levodopa
Ascorbic acid	Levothyroxine
Cefoxitin	Hemoglobin
Cetirizine	Methyldopa
Cotrimoxazol	Ofloxacin
Cyclosporine	Phenazopyridine
Furosemide	Urea
Gentamycin sulfate	Salicylic acid
	Tetracycline
	pH 4.5-9

In case of doubt, please check if a repetition is reasonable after discontinuation of the medication.

For more detailed information on interfering substances, please contact our support via the Roche homepage [www.roche.com/contact.htm](http://www.roche.com/contact.htm).

### Common limitations

The test should not be performed within 3 days of treatment with antibiotics and/or chemotherapy. Relatively large quantities of vitamin C (ascorbic acid) can result in false positive results for glucose, particularly after a meal. As such, performing the test 2 hours after breakfast, is recommended. If possible, this meal should include 3-4 slices of bread or 1-2 rolls.

Allowing the urine to remain in the bladder for long periods of time (4-8 hours, ideally overnight) is recommended as a means of ensuring a high level of accuracy.

In women, this test can yield false results for blood if carried out 3 days before or 3 days after menses; as such, women are advised against performing this test during this period.

### Before you test:

- Always test your first morning urine. Wash your genitals carefully and rinse thoroughly with water (do not dry). To collect urine use only clean, well-rinsed vessels.

- Allow some urine to flow into the toilet. Collect the rest of the stream in the specimen cup. Performing the test:

- Remove test strip. Immediately after you take a test strip out of the vial, seal the vial with the original desiccant stopper. Failure to do this could result in false results, as moisture in the vial could potentially cause color changes in the test parameters.

- Briefly dip the test strip into the urine (approximately 1 sec.), making certain to wet all of the test parameters (fig. 1).

- When removing the test strip, wipe its edge along the rim of the specimen cup in order to remove any excess urine (fig. 2).

- All test parameters should be read at 1 minute (up to 120 seconds for the leukocyte test area for not clearly assignable results). Compare the colors of the test parameters on the strip with the colors on the label (fig. 3).

- Separate sets of color blocks are given for erythrocytes and hemoglobin. Compare the 5th (blood) test parameter area with both rows of colors.

The test should not be performed within 3 days of treatment with antibiotics and/or chemotherapy.

**Protein:** False-positive readings may be found after infusion of polyvinylpyrrolidone (blood plasma), or when the urine collection vessel contains chlorhexidine or traces of disinfectants possessing quaternary ammonium groups.<sup>1</sup>

**Glucose:** The effect of ascorbic acid has been largely eliminated so that glucose concentrations of 100 mg/dL and ascorbic acid concentrations up to 400 mg/L are not likely to give false-negative results.<sup>8</sup>

**Blood/ERY:** Ascorbic acid has virtually no effect on the test.<sup>9</sup> In women the test for blood may be falsified from 3 days before to 3 days after a period. It is therefore advisable not to perform the test during that time. After physical activity, e.g. strenuous jogging, raised values for erythrocytes and protein may occur without being signs of disease.<sup>10</sup> For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

### Expected and result values

Based on literature. Current medical guidelines are leading.

Parameter	Expected values	Additional information
LEU	< 10 LEU/µL	

Paramètre	Résultats
PRO	nég., 30, 100 mg/dL nég., 0,3, 1 g/L 0, +, ++
GLU	norm., 50, 100 mg/dL norm., 2,8, 5,6 mmol/L 0, +, ++
ERY	nég., ~ 5-10, ~ 50 ERY/µL 0, +, ++

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