

Order information

REF	CONTENT	Analyzer(s) on which cobas c pack(s) can be used
08057770 190	ONLINE TDM Gentamicin (100 tests)	System-ID 2059 001 Roche/Hitachi cobas c 503
03375790 190	Preciset TDM I calibrators CAL A-F (1 x 5 mL) Preciset TDM I calibrators Diluent (1 x 10 mL)	Codes 20691-20696
04521536 190	TDM Control Set Level I (2 x 5 mL) TDM Control Set Level II (2 x 5 mL) TDM Control Set Level III (2 x 5 mL)	Code 20310 Code 20311 Code 20312

English**System information****GENT2:** ACN 20590**Intended use**

In vitro test for the quantitative determination of gentamicin in serum and plasma on Roche/Hitachi **cobas c** systems.

Summary

Gentamicin is an aminoglycoside antibiotic that displays broad spectrum, high potency, anti-bacterial action for most susceptible organisms.^{1,2,3,4,5,6,7,8,9} At therapeutic serum concentrations ranging from 4 to 10 µg/mL (8.4 to 20.9 µmol/L), gentamicin is capable of inhibiting the growth of many gram positive cocci, especially penicillinase-producing staphylococci. At concentrations of 10 µg/mL (20.9 µmol/L), most strains of *E. coli*, *Proteus spp.*, *Klebsiella*, *Aerobacter*, *Clostridium*, *Brucella spp.*, *Salmonella*, *Serratia*, and *Shigella* are inhibited. At concentrations ranging from 4 to 10 µg/mL (8.4 to 20.9 µmol/L), gentamicin displays activity against most strains of *Pseudomonas aeruginosa*. Because of these characteristics, gentamicin has been most successfully used in the treatment of serious infections, especially those caused by gram-negative bacilli.^{10,11}

The therapeutic range of gentamicin should be measured at peak as well as trough concentrations. In patients with pre-existing renal damage or those to whom gentamicin has been administered for prolonged periods or in doses above the therapeutic range, hearing impairment and/or nephrotoxicity may develop. Therefore, monitoring of peak and trough gentamicin levels is critical in the prevention of these serious complications with the adjustment of dosage administration as indicated.^{12,13}

Test principle

The assay is based on the kinetic interaction of microparticles in a solution (KIMS). Gentamicin antibody is covalently coupled to microparticles and the drug derivative is linked to a macromolecule. The kinetic interaction of microparticles in solutions is induced by binding of drug-conjugate to the antibody on the microparticles and is inhibited by the presence of gentamicin in the sample. A competitive reaction takes place between the drug conjugate and gentamicin in the serum sample for binding to the gentamicin antibody on the microparticles. The resulting kinetic interaction of microparticles is indirectly proportional to the amount of drug present in the sample.

Reagents - working solutions

- R1** Gentamicin conjugate; piperazine-N,N'-bis (ethanesulfonic acid) (PIPES) buffer, pH 7.2; preservative
- R3** Anti-gentamicin antibody (mouse monoclonal); latex microparticle; 3-(N-morpholino) propane sulfonic acid (MOPS) buffer, pH 7.5; stabilizer; preservative

R1 is in position B and R3 is in position C.

Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines. Safety data sheet available for professional user on request.

For USA: Caution: Federal law restricts this device to sale by or on the order of a physician.

Reagent handling

Ready for use

Carefully invert reagent container several times prior to use to ensure that the reagent components are mixed.

Storage and stability

Shelf life at 2-8 °C: See expiration date on **cobas c** pack label

On-board in use and refrigerated on the analyzer: 12 weeks

Do not freeze.**Specimen collection and preparation**

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum: Collect serum using standard sampling tubes.

Plasma: K₂- or K₃-EDTA, sodium citrate, or sodium, lithium, or ammonium heparin plasma.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Stability: 1 week capped at 2-8 °C
4 weeks capped at -20 °C

Centrifuge samples containing precipitates before performing the assay.

See the limitations and interferences section for details about possible sample interferences.

Sample stability claims were established by experimental data by the manufacturer or based on reference literature and only for the temperatures/time frames as stated in the method sheet. It is the responsibility of the individual laboratory to use all available references and/or its own studies to determine specific stability criteria for its laboratory.

Specimens should not be repeatedly frozen and thawed.

Invert thawed specimens several times prior to testing.

Usual sampling time varies dependent upon desired measurement of peak or trough values.¹⁴

Materials provided

See "Reagents – working solutions" section for reagents.

Materials required (but not provided)

See "Order information" section

General laboratory equipment

Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

The performance of applications not validated by Roche is not warranted and must be defined by the user.

Application for serum and plasma**Test definition**

Reporting time 10 min

Wavelength (sub/main)	800/600 nm		
Reagent pipetting		Diluent (H ₂ O)	
R1	65 µL	–	
R3	62 µL	–	
Sample volumes	Sample	Sample dilution	
		Sample	Diluent (H ₂ O)
Normal	1.3 µL	–	–
Decreased	1.3 µL	–	–
Increased	1.3 µL	–	–

For further information about the assay test definitions refer to the application parameters setting screen of the corresponding analyzer and assay.

Calibration

Calibrators	S1-6: Preciset TDM I calibrators
Calibration mode	Non-linear
Calibration frequency	Full calibration - after cobas c pack change - after reagent lot change - as required following quality control procedures

ACTION REQUIRED

Calibration interval may be extended based on acceptable verification of calibration by the laboratory.

Traceability: This method has been standardized against USP reference standards. The calibrators are prepared to contain known quantities of gentamicin in normal human serum.

Quality control

For quality control, use control materials as listed in the "Order information" section. In addition, other suitable control material can be used.

The control intervals and limits should be adapted to each laboratory's individual requirements. It is recommended to perform quality control always after lot calibration and subsequently at least every 12 weeks.

Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

Calculation

Roche/Hitachi **cobas c** systems automatically calculate the analyte concentration of each sample in the unit µg/mL (µmol/L, mg/L).

Conversion factors: ¹⁵	µg/mL x 2.09 = µmol/L
	µg/mL x 1.0 = mg/L

Limitations - interference

Criterion: Recovery within ± 10 % of initial value at gentamicin levels of approximately 2 and 6 µg/mL.

Serum/plasma

Icterus:¹⁶ No significant interference up to an I index of 50 for conjugated and unconjugated bilirubin (approximate conjugated and unconjugated bilirubin concentration: 855 µmol/L or 50 mg/dL).

Hemolysis:¹⁶ No significant interference up to an H index of 1000 (approximate hemoglobin concentration: 621 µmol/L or 1000 mg/dL).

Lipemia (Intralipid):¹⁶ No significant interference up to an L index of 150. There is poor correlation between the L index (corresponds to turbidity) and triglycerides concentration.

Triglycerides: No significant interference from triglycerides up to a concentration of 1000 mg/dL (11.3 mmol/L).

Rheumatoid factors: No significant interference from rheumatoid factors up to a concentration of 100 IU/mL.

Total protein: No significant interference from total protein up to a concentration of 12 g/dL.

Note

A negative bias of up to approximately 20 % has been observed with this assay for some samples artificially spiked with gentamicin sulfate. Patient samples have been verified to recover correctly.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

ACTION REQUIRED

Special Wash Programming: The use of special wash steps is mandatory when certain test combinations are run together on Roche/Hitachi **cobas c** systems. All special wash programming necessary for avoiding carry-over is available via the **cobas** link. The latest version of the carry-over evasion list can be found with the NaOHD/SMS/SCCS Method Sheet for information. For further instructions refer to the operator's manual.

Limits and ranges**Measuring range**

Measuring range: 0.4-10.0 µg/mL (0.84-20.9 µmol/L)

Manually dilute samples having higher concentrations with Preciset TDM I diluent (0 µg/mL) (1 + 1) and reassay. Multiply the result by 2 to obtain the specimen value.

Lower limits of measurement**Limit of Blank, Limit of Detection and Limit of Quantitation**

Limit of Blank = 0.3 µg/mL (0.63 µmol/L)

Limit of Detection = 0.4 µg/mL (0.84 µmol/L)

Limit of Quantitation = 0.8 µg/mL (1.67 µmol/L)

The Limit of Blank, Limit of Detection and Limit of Quantitation were determined in accordance with the CLSI (Clinical and Laboratory Standards Institute) EP17-A2 requirements.

The Limit of Blank is the 95th percentile value from n ≥ 60 measurements of analyte-free samples over several independent series. The Limit of Blank corresponds to the concentration below which analyte-free samples are found with a probability of 95 %.

The Limit of Detection is determined based on the Limit of Blank and the standard deviation of low concentration samples.

The Limit of Detection corresponds to the lowest analyte concentration which can be detected (value above the Limit of Blank with a probability of 95 %).

The Limit of Quantitation is the lowest analyte concentration that can be reproducibly measured with a total error of 20 %. It has been determined using low concentration gentamicin samples.

Expected values

Although optimum values may vary, peak serum values in the range of 6 to 10 µg/mL

(12.5 to 20.9 µmol/L)*
and trough values in the range of 0.5 to 2.0 µg/mL

(1.0 to 4.2 µmol/L)* are generally accepted for therapeutic effectiveness.¹⁷

*calculated by unit conversion factor

The achievement of non-toxic, but therapeutic, serum levels is often difficult, even in patients with normal renal function. Complications encountered with the use of gentamicin are ototoxicity and nephrotoxicity.^{10,18,19,20,21} However, these reactions are predictable, and close patient monitoring is essential for the successful use of this agent. The most serious toxic effect of gentamicin is permanent damage to the vestibular division of the eighth cranial nerve, which has been reported to occur most frequently in patients with renal failure. Since gentamicin is inherently unstable, is not metabolized and is excreted primarily by glomerular filtration, toxic concentrations of the drug may accumulate in the body when the dosage is not adjusted for patients with impaired renal function. While high serum levels can be toxic, indiscriminately low dosages of gentamicin will result in ineffective treatment for many strains of gram-negative bacteria. The indiscriminate use of low dosages of gentamicin may not only engender the emergence of gentamicin-resistant organisms, but also the emergence of aminoglycoside-resistant organisms.^{11,22,23} Current literature reflects increasing interest in once daily dosing versus the conventional administration of drug 2 to 4 times daily.

Adoption of once daily dosing may require a revision of target peak and trough concentrations.^{24,25,26}

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

Specific performance data

Representative performance data on the analyzers are given below. These data represent the performance of the analytical procedure itself.

Results obtained in individual laboratories may differ due to heterogenous sample materials, aging of analyzer components and mixture of reagents running on the analyzer.

Precision

Precision was determined using human samples and controls in accordance with the CLSI (Clinical and Laboratory Standards Institute) EP05-A3 requirements with repeatability (n = 84) and intermediate precision (2 aliquots per run, 2 runs per day, 21 days). The following results were obtained:

Serum/plasma

Repeatability	Mean	SD	CV
	µg/mL	µg/mL	%
TDMC1 ^{a)}	1.81	0.0375	2.1
TDMC2 ^{b)}	4.49	0.0558	1.2
TDMC3 ^{c)}	6.48	0.0786	1.2
Human serum 1	0.663	0.0392	5.9
Human serum 2	1.94	0.0353	1.8
Human serum 3	4.60	0.0518	1.1
Human serum 4	7.76	0.108	1.4
Human serum 5	7.90	0.129	1.6
Intermediate precision	Mean	SD	CV
	µg/mL	µg/mL	%
TDMC1 ^{a)}	1.81	0.0544	3.0
TDMC2 ^{b)}	4.45	0.0720	1.6
TDMC3 ^{c)}	6.48	0.106	1.6
Human serum 1	0.663	0.0525	7.9
Human serum 2	1.89	0.0581	3.1
Human serum 3	4.63	0.0876	1.9
Human serum 4	7.76	0.149	1.9
Human serum 5	7.90	0.167	2.1

a) TDM Control Set Level I

b) TDM Control Set Level II

c) TDM Control Set Level III

Method comparison

Serum/plasma

Gentamicin values for human serum and plasma samples obtained on a Roche/Hitachi **cobas c** 503 analyzer (y) were compared with those determined using the corresponding reagent on a Roche/Hitachi **cobas c** 501 analyzer (x).

Sample size (n) = 82

Passing/Bablok²⁷

$$y = 1.000x - 0.0200 \mu\text{g/mL}$$

$$r = 0.970$$

Linear regression

$$y = 1.002x - 0.0312 \mu\text{g/mL}$$

$$r = 0.999$$

The sample concentrations were between 0.41 and 9.81 µg/mL.

Analytical specificity

The following compounds were tested for cross-reactivity.

Compound	Concentration tested (µg/mL)	% cross-reactivity
Netilmicin	70	9.13
Sisomicin	131	8.16
Methotrexate	23	< 1.0
Tetracycline	40	< 1.0
Amikacin	250	< 0.1
Cephalexin	500	< 0.1
Chloramphenicol	300	< 0.1
Clindamycin	500	< 0.1
Kanamycin	250	< 0.1
Neomycin	100	< 0.1
Spectinomycin	200	< 0.1
Streptomycin	200	< 0.1
Tobramycin	100	< 0.1
Vancomycin	400	< 0.1
Amphotericin B	50	< 0.01
Ampicillin	78	< 0.01
Carbenicillin	500	< 0.01
Cephalosporin C	432	< 0.01
Cephalothin	63	< 0.01
Erythromycin	200	< 0.01
5-Fluorocytosine	700	< 0.01
Furosemide	100	< 0.01
Methylprednisolone	500	< 0.01
Oxytetracycline	37	< 0.01
Prednisolone	500	< 0.01

Tests were performed on 16 drugs. No significant interference with the assay was found.

Acetaminophen	Doxycycline (Tetracycline)
Acetyl cysteine	Ibuprofen
Acetylsalicylic acid	Levodopa
Ampicillin-Na	Methyldopa + 1.5 H ₂ O
Ascorbic acid	Metronidazole
Ca-Dobesilate	Phenylbutazone
Cefoxitin	Rifampicin
Cyclosporine	Theophylline

References

- Black J, Calesnick B, Williams D, et al. Pharmacology of gentamicin, a new broad-spectrum antibiotic. *Antimicrob Agents Chemother* 1963;161:138-147.
- Oden EM, Stander H, Weinstein MJ. Microbiological assay of gentamicin. *Antimicrob Agents Chemother* 1963;161:8-13.
- Waitz JA, Weinstein M. Recent laboratory studies with gentamicin with reference to structure and resistance. *Med J Aust* 1970;1(Suppl):19-25.
- MacMillan BG. Ecology of bacteria colonizing the burned patient given topical and systemic gentamicin therapy: a five-year study. *J Infect Dis* 1971;124(Suppl):278-286.
- Nunnery AW, Riley HD Jr. Gentamicin: clinical and laboratory studies in infants and children. *J Infect Dis* 1969;119(4):460-464.

- 6 Young LS, Louria DB, Armstrong D. Gentamicin in the treatment of severe, hospital-acquired gram-negative infections. *Trans NY Acad Sci* 1967;29(5):579-588.
- 7 Holloway WJ, Taylor WA. Gentamicin and kanamycin in the treatment of gram-negative sepsis: a comparative study. *J Infect Dis* 1971;124(Suppl):180-184.
- 8 Bodey GP, Middleman E, Umsawasdi T, et al. Intravenous gentamicin therapy for infections in patients with cancer. *J Infect Dis* 1971;124(Suppl):174-179.
- 9 Klein JO, Herschel M, Therakan RM, et al. Gentamicin in serious neonatal infections: absorption, excretion, and clinical results in 25 cases. *J Infect Dis* 1971;124(Suppl):224-231.
- 10 Cox CE. Gentamicin. *Med Clin North Am* 1970;54:1305-1315.
- 11 Riff LJ, Jackson GG. Pharmacology of gentamicin in man. *J Infect Dis* 1971;124(Suppl):98-105.
- 12 Kahlmeter G. Gentamicin and tobramycin. Clinical pharmacokinetics and nephrotoxicity. Aspects on assay techniques. *Scand J Infect Dis* 1979;132(Suppl 18):1-40.
- 13 Barza M, Lauermann M. Why monitor serum levels of gentamicin? *Clin Pharm* 1978;3:202-215.
- 14 O'Bey KA, Jim LK, Gee JP, et al. Temperature dependence of the stability of tobramycin mixed with penicillins in human serum. *Am J Hosp Pharm* 1982;39(6):1005-1008.
- 15 Tietz NW. *Fundamentals of Clinical Chemistry*. 3rd ed. Philadelphia, PA: WB Saunders Co 1995:834.
- 16 Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. *Clin Chem* 1986;32:470-475.
- 17 Taylor WJ, Diers-Caviness MH (eds.). *A Textbook for the Clinical Application of Therapeutic Drug Monitoring*. Abbott Laboratories 1986:296.
- 18 Jackson GG, Arcieri GM. Ototoxicity of gentamicin in man: a survey and controlled analysis of clinical experience in the United States. *J Infect Dis* 1971;124(Suppl):130-137.
- 19 Arcieri GM, Falco FG, Smith HM, et al. Clinical research experience with gentamicin. Incidence of adverse reactions. *Med J Aust* 1970;1(Suppl):30-34.
- 20 Falco FG, Smith HM, Arcieri GM. Nephrotoxicity of aminoglycosides and gentamicin. *J Infect Dis* 1969;119(4):406-409.
- 21 Gyselync AM, Forrey A, Cutler R. Pharmacokinetics of gentamicin: distribution and plasma and renal clearance. *J Infect Dis* 1971;124(Suppl):70-76.
- 22 Weinstein MJ, Drube CG, Moss EL Jr, et al. Microbiologic studies related to bacterial resistance to gentamicin. *J Infect Dis* 1971;124(Suppl):11-17.
- 23 Jackson GG, Riff LJ. Pseudomonas bacteremia: pharmacologic and other bases for failure of treatment with gentamicin. *J Infect Dis* 1971;124(Suppl):185-191.
- 24 Barclay ML, Begg EJ, Hickling KG. What is the evidence for once-daily aminoglycoside therapy? *Clin Pharmacokinet* 1994;27(1):32-48.
- 25 Hustinx WN, Hoepelman IM. Aminoglycoside dosage regimens. Is once a day enough? *Clin Pharmacokinet* 1993;25(6):427-432.
- 26 Bertino JS Jr, Rodvold KA, Destache CJ. Cost considerations in therapeutic drug monitoring of aminoglycosides. *Clin Pharmacokinet* 1994;26(1):71-81.
- 27 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. *J Clin Chem Clin Biochem* 1988 Nov;26(11):783-790.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard (for USA: see <https://usdiagnostics.roche.com> for definition of symbols used):

CONTENT



GTIN

Contents of kit

Volume after reconstitution or mixing

Global Trade Item Number

FOR US CUSTOMERS ONLY: LIMITED WARRANTY

Roche Diagnostics warrants that this product will meet the specifications stated in the labeling when used in accordance with such labeling and will be free from defects in material and workmanship until the expiration date printed on the label. THIS LIMITED WARRANTY IS IN LIEU OF ANY OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. IN NO EVENT SHALL ROCHE DIAGNOSTICS BE LIABLE FOR INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES.

COBAS, COBAS C, ONLINE TDM and PRECISET are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2019, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com

Distribution in USA by:
Roche Diagnostics, Indianapolis, IN
US Customer Technical Support 1-800-428-2336

