

anti-Epithelial Related Antigen (MOC-31) Mouse Monoclonal Primary Antibody

REF 790-4561
06374409001

IVD 50




Figure 1. anti-Epithelial Related Antigen (MOC-31) antibody staining of neoplastic cells in lung adenocarcinoma tissue.

INTENDED USE

Anti-Epithelial Related Antigen (MOC-31) Mouse Monoclonal Primary Antibody is intended for laboratory use in the qualitative immunohistochemical detection of epithelial cell adhesion molecule (EpCAM) by light microscopy in sections of formalin-fixed, paraffin-embedded tissue stained on a BenchMark IHC/ISH instrument.

This product should be interpreted by a qualified pathologist in conjunction with histological examination, relevant clinical information, and proper controls. This antibody is intended for in vitro diagnostic (IVD) use.

SUMMARY AND EXPLANATION

Epithelial cell adhesion molecule (EpCAM) is a type 1 transmembrane superficial glycoprotein encoded by the EPCAM gene located on chromosome 2p21.¹ EpCAM mediates calcium-independent cell-cell interaction.² However, it is not structurally related to any of the four major families of adhesion proteins: cadherins, selectins, integrins, and the immunoglobulins.³ The EpCAM protein consists of a large N-terminal extracellular domain, a single transmembrane domain, and a short C-terminal cytoplasmic domain.² The extracellular domain facilitates the adhesions between epithelial cells.⁴ EpCAM also interacts with various adhesion proteins such as CD44, claudins, and E-Cadherin.⁵ In addition to adhesion, EpCAM has numerous roles in various biological processes including signaling, migration, proliferation, and differentiation.^{2,6} EpCAM is an important cell-signaling molecule that controls four independent pathways: nPKC-dependent pathway, Wnt signal pathway, ERas/AKT pathway, and a pathway that regulates cell proliferation.^{6,7} EpCAM also increases cell motility and migration by disrupting the link between α -catenin and actin filaments causing a reduction in cell-cell adhesion.^{2,6,8} Studies show that EpCAM increases the proliferation of tumor cell lines. However, it remains to be determined if EpCAM directly activates components of cell cycle machinery or if EpCAM-activated proliferation is a secondary effect through repression, apoptosis, elevation of cell metabolism, or interruption of anti-proliferative signals.^{9,10,11} Lastly, EpCAM plays an essential role in differentiation during morphogenesis and tissue regeneration.¹²

EpCAM is expressed in epithelial tissues at the lateral and basal membrane of cells as well as neoplasms derived from these tissue types.^{2,3} Expression levels of EpCAM vary between different organs and cell types.^{2,5} Strong immunoreactivity for EpCAM can be found in adult epithelia of the colon, small intestine, pancreas, liver, gall bladder, and endometrium. In most cases, EpCAM expression is positively correlated with proliferating cells and negatively correlated with more differentiated areas.^{2,5}

Detection of EpCAM by immunohistochemistry (IHC) with the Anti-Epithelial Related Antigen (MOC-31) Mouse Monoclonal Primary Antibody (anti-Epithelial Related Antigen (MOC-31) antibody) may be used to aid in the identification of normal and neoplastic epithelial cells and aid in the differentiation of lung adenocarcinoma and mesothelioma. It may be used as part of a panel of IHC studies. The staining pattern is membranous and/or cytoplasmic.

PRINCIPLE OF THE PROCEDURE

The anti-Epithelial Related Antigen (MOC-31) antibody binds to EpCAM (Epithelial Related Antigen) in formalin-fixed, paraffin-embedded (FFPE) tissue sections. This antibody can be visualized using *ultraView* Universal DAB Detection Kit (Cat. No. 760-500 / 05269806001). Refer to the respective method sheet for further information.

MATERIAL PROVIDED

Anti-Epithelial Related Antigen (MOC-31) antibody contains sufficient reagent for 50 tests. One 5 mL dispenser of anti-Epithelial Related Antigen (MOC-31) antibody contains approximately 0.5 μ g of a mouse monoclonal antibody.

The antibody is diluted in Tris-HCl with carrier protein and 0.10% ProClin 300, a preservative.

Specific antibody concentration is approximately 0.1 μ g/mL. There is no known non-specific antibody reactivity observed in this product.

Anti-Epithelial Related Antigen (MOC-31) antibody is a mouse monoclonal antibody produced as cell culture supernatant.

Refer to the appropriate VENTANA detection kit method sheet for detailed descriptions of: Principle of the Procedure, Material and Methods, Specimen Collection and Preparation for Analysis, Quality Control Procedures, Troubleshooting, Interpretation of Results, and Limitations.

MATERIALS REQUIRED BUT NOT PROVIDED

Staining reagents, such as VENTANA detection kits and ancillary components, including negative and positive tissue control slides, are not provided.

Not all products listed in the method sheet may be available in all geographies. Consult your local support representative.

The following reagents and materials may be required for staining but are not provided:

1. Recommended control tissue
2. Microscope slides, positively charged
3. Negative Control (Monoclonal) reagent, (Cat. No. 760-2014 / 05266670001)
4. *ultraView* Universal DAB Detection Kit (Cat. No. 760-500 / 05269806001)
5. EZ Prep Concentrate (10X) (Cat. No. 950-102 / 05279771001)
6. Reaction Buffer Concentrate (10X) (Cat. No. 950-300 / 05353955001)
7. LCS (Predilute) (Cat. No. 650-010 / 05264839001)
8. ULTRA LCS (Predilute) (Cat. No. 650-210 / 05424534001)
9. Cell Conditioning Solution (CC1) (Cat. No. 950-124 / 05279801001)
10. ULTRA Cell Conditioning Solution (ULTRA CC1) (Cat. No. 950-224 / 05424569001)
11. Hematoxylin II (Cat. No. 790-2208 / 05277965001)
12. Bluing Reagent (Cat. No. 760-2037 / 05266769001)
13. Permanent mounting medium
14. Cover glass
15. Automated coverslipper
16. General purpose laboratory equipment
17. BenchMark IHC/ISH instrument

STORAGE AND STABILITY

Upon receipt and when not in use, store at 2-8°C. Do not freeze.

To ensure proper reagent delivery and the stability of the antibody, replace the dispenser cap after every use and immediately place the dispenser in the refrigerator in an upright position.

Every antibody dispenser is expiration dated. When properly stored, the reagent is stable to the date indicated on the label. Do not use reagent beyond the expiration date.

SPECIMEN PREPARATION

Routinely processed FFPE tissues are suitable for use with this primary antibody when used with VENTANA detection kits and BenchMark IHC/ISH instruments. The recommended tissue fixative is 10% neutral buffered formalin.¹³ Sections should be cut at approximately 4 μ m in thickness and mounted on positively charged slides. Slides should be stained immediately, as antigenicity of cut tissue sections may diminish over time. Ask your Roche representative for a copy of "Recommended Slide Storage and Handling" for more information.

It is recommended that positive and negative controls be run simultaneously with unknown specimens.

WARNINGS AND PRECAUTIONS

- For in vitro diagnostic (IVD) use.
- For professional use only.
- Do not use beyond the specified number of tests.
- ProClin 300 solution is used as a preservative in this reagent. It is classified as an irritant and may cause sensitization through skin contact. Take reasonable precautions when handling. Avoid contact of reagents with eyes, skin, and mucous membranes. Use protective clothing and gloves.
- Positively charged slides may be susceptible to environmental stresses resulting in inappropriate staining. Ask your Roche representative for more information on how to use these types of slides.
- Materials of human or animal origin should be handled as biohazardous materials and disposed of with proper precautions. In the event of exposure, the health directives of the responsible authorities should be followed.^{14,15}
- Avoid contact of reagents with eyes and mucous membranes. If reagents come in contact with sensitive areas, wash with copious amounts of water.
- Avoid microbial contamination of reagents as it may cause incorrect results.
- For further information on the use of this device, refer to the BenchMark IHC/ISH instrument User Guide, and instructions for use of all necessary components navifyportal.roche.com.
- Consult local and/or state authorities with regard to recommended method of disposal.
- Product safety labeling primarily follows EU GHS guidance. Safety data sheet available for professional user on request.
- To report suspected serious incidents related to this device, contact the local Roche representative and the competent authority of the Member State or Country in which the user is established.

This product contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

Table 1. Hazard information.

| Hazard | Code | Statement |
|--------|-------------|---|
| | H317 | May cause an allergic skin reaction. |
| | H412 | Harmful to aquatic life with long lasting effects. |
| | P261 | Avoid breathing mist or vapours. |
| | P273 | Avoid release to the environment. |
| | P280 | Wear protective gloves. |
| | P333 + P313 | If skin irritation or rash occurs: Get medical advice/ attention. |
| | P362 + P364 | Take off contaminated clothing and wash it before reuse. |
| | P501 | Dispose of contents/ container to an approved waste disposal plant. |

This product contains CAS # 55965-84-9, reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1).

STAINING PROCEDURE

VENTANA primary antibodies have been developed for use on BenchMark IHC/ISH instruments in combination with VENTANA detection kits and accessories. Refer to Table 2 for recommended staining protocols.

This antibody has been optimized for specific incubation times but the user must validate results obtained with this reagent.

The parameters for the automated procedures can be displayed, printed and edited according to the procedure in the instrument User Guide. Refer to the appropriate VENTANA detection kit method sheet for more details regarding immunohistochemistry staining procedures.

For more details on the proper use of this device, refer to the inline dispenser method sheet associated with P/N 790-4561.

Table 2. Recommended staining protocol for anti-Epithelial Related Antigen (MOC-31) antibody with *ultraView* Universal DAB Detection Kit on BenchMark IHC/ISH instruments.

| Procedure Type | Method | |
|---------------------------------------|---------------------------|----------------------------------|
| | XT | ULTRA or ULTRA PLUS ^a |
| Deparaffinization | Selected | Selected |
| Cell Conditioning (Antigen Unmasking) | CC1, Mild | ULTRA CC1, Mild |
| Enzyme (Protease) | None Required | None Required |
| Antibody (Primary) | 32 minutes, 37°C | 32 minutes, 36°C |
| Counterstain | Hematoxylin II, 4 minutes | |
| Post Counterstain | Bluing, 4 minutes | |

^a Concordance was demonstrated between BenchMark ULTRA and BenchMark ULTRA PLUS instruments using representative assays.

Due to variation in tissue fixation and processing, as well as general lab instrument and environmental conditions, it may be necessary to increase or decrease the primary antibody incubation, cell conditioning or protease pretreatment based on individual specimens, detection used, and reader preference. For further information on fixation variables, refer to "Immunohistochemistry Principles and Advances."¹⁶

NEGATIVE REAGENT CONTROL

In addition to staining with anti-Epithelial Related Antigen (MOC-31) antibody, a second slide should be stained with the appropriate negative control reagent.

POSITIVE TISSUE CONTROL

Optimal laboratory practice is to include a positive control section on the same slide as the test tissue. This helps identify any failures applying reagents to the slide. Tissue with weak positive staining is best suited for quality control. Control tissue may contain both positive and negative staining elements and serve as both the positive and negative control. Control tissue should be fresh autopsy, biopsy, or surgical specimen, prepared or fixed as soon as possible in a manner identical to test sections.

Known positive tissue controls should be utilized only for monitoring performance of reagents and instruments, not as an aid in determining specific diagnosis of test samples. If the positive tissue controls fail to demonstrate positive staining, results of the test specimen should be considered invalid.

Examples of positive control tissues for this antibody are normal kidney, breast, lung and colon. The staining pattern should be a moderate to strong membranous staining of the epithelial cells in these normal tissues.

STAINING INTERPRETATION / EXPECTED RESULTS

The cellular staining pattern for anti-Epithelial Related Antigen (MOC-31) antibody is membranous and/or cytoplasmic.

SPECIFIC LIMITATIONS

All assays might not be registered on every instrument. Please contact your local Roche representative for more information.

PERFORMANCE CHARACTERISTICS

ANALYTICAL PERFORMANCE

Staining tests for specificity, sensitivity, and precision were conducted and the results are listed below.

Sensitivity and Specificity

Table 3. Sensitivity/Specificity of anti-Epithelial Related Antigen (MOC-31) antibody was determined by testing FFPE normal tissues.

| Tissue | # positive / total cases | Tissue | # positive / total cases |
|------------|--------------------------|--------|--------------------------|
| Cerebrum | 0/3 | Lung | 20/33 |
| Cerebellum | 0/3 | Heart | 0/3 |

| Tissue | # positive / total cases | Tissue | # positive / total cases |
|-------------------|--------------------------|--------------------------|--------------------------|
| Adrenal gland | 0/3 | Esophagus | 0/3 |
| Ovary | 0/3 | Stomach | 3/3 |
| Pancreas | 3/3 | Small intestine | 3/3 |
| Parathyroid gland | 3/3 | Colon | 3/3 |
| Pituitary gland | 3/3 | Liver | 2/3 |
| Testis | 0/3 | Salivary gland | 3/3 |
| Thyroid | 3/3 | Kidney | 41/41 |
| Breast | 3/3 | Prostate | 3/3 |
| Spleen | 0/3 | Endometrium | 2/2 |
| Tonsil | 0/3 | Skin | 0/1 |
| Thymus | 2/3 | Mesothelium ^a | 4/10 |
| Bone marrow | 0/3 | | |

^a The normal mesothelium cases listed are from pleura/lung (n=7) and omentum (n=3).

Table 4. Sensitivity/Specificity of anti-Epithelial Related Antigen (MOC-31) antibody was determined by testing a variety of FFPE neoplastic tissues.

| Pathology | # positive / total cases |
|---|--------------------------|
| Glioblastoma (Cerebrum) | 0/1 |
| Meningioma (Cerebrum) | 0/1 |
| Ependymoma (Cerebrum) | 0/1 |
| Oligodendroglioma (Cerebrum) | 0/1 |
| Serous adenocarcinoma (Ovary) | 1/1 |
| Mucinous adenocarcinoma (Ovary) | 1/1 |
| Neuroendocrine neoplasm (Pancreas) | 1/1 |
| Adenocarcinoma (Pancreas) | 1/1 |
| Seminoma (Testis) | 1/1 |
| Embryonal carcinoma (Testis) | 1/1 |
| Medullary carcinoma (Thyroid) | 1/1 |
| Papillary carcinoma (Thyroid) | 1/1 |
| Ductal carcinoma in situ (Breast) | 2/2 |
| Invasive ductal carcinoma (Breast) | 1/1 |
| Small cell carcinoma (Lung) | 7/7 |
| Squamous cell carcinoma (Lung) | 32/36 |
| Adenosquamous carcinoma (Lung) | 2/2 |
| Adenocarcinoma (Lung) | 59/61 |
| Papillary carcinoma (Lung) | 4/4 |
| Mucinous adenocarcinoma (Lung) | 1/2 |
| Adenocarcinoma in situ (Lung) | 4/4 |
| Large cell carcinoma (Lung) | 3/4 |
| Neuroendocrine carcinoma, atypical carcinoid tumor (Lung) | 4/5 |
| Squamous cell carcinoma (Esophagus) | 1/1 |
| Adenocarcinoma (Esophagus) | 1/1 |

| Pathology | # positive / total cases |
|--|--------------------------|
| Mucinous adenocarcinoma (Stomach) | 1/1 |
| Gastrointestinal stromal tumor (GIST) (Small intestine) | 0/1 |
| Adenocarcinoma (Colon) | 2/2 |
| Gastrointestinal stromal tumor (GIST) (Abdominal cavity) | 0/1 |
| Adenocarcinoma (Rectum) | 1/1 |
| Gastrointestinal stromal tumor (GIST) (Rectum) | 0/1 |
| Melanoma (Rectum) | 0/1 |
| Hepatocellular carcinoma (Liver) | 1/1 |
| Hepatoblastoma (Liver) | 0/1 |
| Clear cell carcinoma (Kidney) | 0/1 |
| Adenocarcinoma (Prostate) | 2/2 |
| Leiomyoma (Uterus) | 0/1 |
| Adenocarcinoma (Uterus) | 1/1 |
| Clear cell carcinoma (Uterus) | 1/1 |
| Squamous cell carcinoma (Cervix) | 0/2 |
| Embryonal rhabdomyosarcoma (Striated muscle) | 0/1 |
| Basal cell carcinoma (Skin) | 0/1 |
| Squamous cell carcinoma (Skin) | 1/1 |
| Neurofibroma (Soft tissue) | 0/1 |
| Ganglioneuroblastoma (Retroperitoneum) | 0/1 |
| Mesothelioma | 4/24 |
| Spindle cell rhabdomyosarcoma (Retroperitoneum) | 0/1 |
| Hodgkin lymphoma (Lymph node) | 0/1 |
| Lymphoma, NOS | 0/1 |
| Diffuse large B-cell lymphoma (DLBCL) | 0/2 |
| B-cell lymphoma, NOS | 0/1 |
| Urothelial carcinoma (Bladder) | 1/1 |
| Leiomyosarcoma (Bladder) | 0/1 |
| Osteosarcoma (Bone) | 0/1 |
| Leiomyosarcoma (Smooth muscle) | 0/1 |

Precision

Precision studies for anti-Epithelial Related Antigen (MOC-31) antibody were completed to demonstrate:

- Between lot precision of the antibody.
- Within run and between day precision on a BenchMark XT instrument.
- Between instrument precision on the BenchMark XT and BenchMark ULTRA instruments.
- Between platform precision between the BenchMark XT and BenchMark ULTRA instrument.

All studies met their acceptance criteria.

Precision on the BenchMark ULTRA PLUS instrument was demonstrated using representative assays. Studies included Within-run Repeatability, Between-day and Between-run Intermediate Precision. All studies met their acceptance criteria.

CLINICAL PERFORMANCE

Clinical performance data relevant to the intended purpose of anti-Epithelial Related Antigen (MOC-31) antibody were assessed by systematic review of the literature. The data gathered support the use of the device in accordance with its intended purpose.

REFERENCES

- Herreros-Pomares A, Aguilar-Gallardo C, Calabuig-Farinas S, et al. Epcam Duality Becomes This Molecule in a New Dr. Jekyll and Mr. Hyde Tale. *Crit Rev Oncol Hematol*. 2018;126:52-63.
- Keller L, Werner S, Pantel K. Biology and Clinical Relevance of Epcam. *Cell Stress*. 2019;3(6):165-180.
- Ordóñez NG. Broad-Spectrum Immunohistochemical Epithelial Markers: A Review. *Hum Pathol*. 2013;44(7):1195-1215.
- Appella E, Weber IT, Blasi F. Structure and Function of Epidermal Growth Factor-Like Regions in Proteins. *FEBS Lett*. 1988;231(1):1-4.
- Schnell U, Cirulli V, Giepmans BN. Epcam: Structure and Function in Health and Disease. *Biochim Biophys Acta*. 2013;1828(8):1989-2001.
- Huang L, Yang Y, Yang F, et al. Functions of Epcam in Physiological Processes and Diseases (Review). *Int J Mol Med*. 2018;42(4):1771-1785.
- Maetzel D, Denzel S, Mack B, et al. Nuclear Signalling by Tumour-Associated Antigen Epcam. *Nat Cell Biol*. 2009;11(2):162-171.
- Martowicz A, Seeber A, Untergasser G. The Role of Epcam in Physiology and Pathology of the Epithelium. *Histol Histopathol*. 2016;31(4):349-355.
- Munz M, Kieu C, Mack B, et al. The Carcinoma-Associated Antigen Epcam Upregulates C-Myc and Induces Cell Proliferation. *Oncogene*. 2004;23(34):5748-5758.
- Osta WA, Chen Y, Mikhitarian K, et al. Epcam Is Overexpressed in Breast Cancer and Is a Potential Target for Breast Cancer Gene Therapy. *Cancer Res*. 2004;64(16):5818-5824.
- Chaves-Perez A, Mack B, Maetzel D, et al. Epcam Regulates Cell Cycle Progression Via Control of Cyclin D1 Expression. *Oncogene*. 2013;32(5):641-650.
- Kuechlin S, Schoels M, Slanchev K, et al. Epcam Controls Morphogenetic Programs During Zebrafish Pronephros Development. *Biochem Biophys Res Commun*. 2017;487(2):209-215.
- Carson FL, Cappellano C. *Histotechnology; A Self-Instructional Text*, 5th edition. American Society for Clinical Pathology Press; 2020, 2022.
- Occupational Safety and Health Standards: Occupational exposure to hazardous chemicals in laboratories. (29 CFR Part 1910.1450). *Fed. Register*.
- Directive 2000/54/EC of the European Parliament and Council of 24 June 2020 on the protection of workers from risks related to exposure to biological agents at work.
- Roche PC, Hsi ED. *Immunohistochemistry-Principles and Advances*. Manual of Clinical Laboratory Immunology, 6th edition. In: NR Rose, ed. ASM Press; 2002.

NOTE: A point (period/stop) is always used in this document as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral.

Separators for thousands are not used.

The summary of safety and performance can be found here:

<https://ec.europa.eu/tools/eudamed>

Symbols

Ventana uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard (for USA: see elabdoc.roche.com/symbols for more information).

GTIN Global Trade Item Number

Rx only For USA: Caution: Federal law restricts this device to sale by or on the order of a physician.

REVISION HISTORY

| Rev | Updates |
|-----|---|
| E | Updates to Warnings and Precautions section. Updated to current template. |

INTELLECTUAL PROPERTY

VENTANA, BENCHMARK, and ULTRAVIEW are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

© 2024 Ventana Medical Systems, Inc.

For USA: Rx only

CONTACT INFORMATION



Ventana Medical Systems, Inc.
1910 E. Innovation Park Drive
Tucson, AZ 85755
USA
+1 520 887 2155
+1 800 227 2155 (USA)

www.roche.com



Roche Diagnostics GmbH
Sandhofer Strasse 116
68305 Mannheim
Germany
+800 5505 6606

