


**γ-Glutamyltransferase ver.2 - Standardized against IFCC / Szasz**

### Order information

REF		CONTENT		Analyzer(s) on which <b>cobas c</b> pack(s) can be used
05168775190*	05168775500	γ-Glutamyltransferase ver.2 (1200 tests)	System-ID 03 6598 8	<b>cobas c</b> 701/702
05168775214*	05168775500	γ-Glutamyltransferase ver.2 (1200 tests)	System-ID 03 6598 8	<b>cobas c</b> 701/702

Materials required (but not provided):

10759350190	Calibrator f.a.s. (12 x 3 mL)	Code 401	
12149435122	Precinorm U plus (10 x 3 mL)	Code 300	
12149443122	Precipath U plus (10 x 3 mL)	Code 301	
05117003190	PreciControl ClinChem Multi 1 (20 x 5 mL)	Code 391	
05947626190	PreciControl ClinChem Multi 1 (4 x 5 mL)	Code 391	
05117216190	PreciControl ClinChem Multi 2 (20 x 5 mL)	Code 392	
05947774190	PreciControl ClinChem Multi 2 (4 x 5 mL)	Code 392	
05172152190	Diluent NaCl 9 % (119 mL)	System-ID 08 6869 3	

\* Some kits shown may not be available in all countries.

### English

#### System information

**GGT12:** ACN 8220: assay standardized against IFCC

**GGT2:** ACN 8480: assay standardized against Szasz

#### Intended use

In vitro test for the quantitative determination of γ-glutamyltransferase (GGT) in human serum and plasma on **cobas c** systems.

#### Summary

Measurements of γ-glutamyltransferase (GGT) performed with this assay in human serum and plasma are used in the diagnosis and monitoring of hepatobiliary diseases, as well as a screening test for occult alcoholism.

Mature GGT is a dimeric glycoprotein weighing 68 kDa. It is found in the kidneys, liver, pancreas, and intestine, with the highest abundance in renal tissue. However, the primary source of GGT activity in the serum is the liver.<sup>1</sup>

In clinical practice, GGT serum levels are typically measured alongside a full blood count, bilirubin, albumin, transaminases (ALT and AST), and alkaline phosphatases (ALP) as an initial investigation for potential liver disease.<sup>2</sup> GGT is considered one of the most reliable indicators for the development of liver disease.<sup>3</sup> Multiple guidelines recommend GGT testing as part of the diagnostic workup and monitoring for various liver diseases. Additionally, GGT serves as a well-established marker for alcohol-related liver disease and excessive alcohol consumption.<sup>4,5,6,7,8,9,10</sup> Increased GGT is observed as a result of obesity, excess alcohol consumption or may be induced by drugs, including phenobarbital and phenytoin.<sup>1</sup>

In 1969, Szasz published the first kinetic procedure for GGT in serum using γ-glutamyl-p-nitroanilide as substrate and glycylglycine as acceptor.<sup>11</sup> In order to circumvent the poor solubility of γ-glutamyl-p-nitroanilide, Persijn and van der Slik investigated various derivatives and found the water-soluble substrate L-γ-glutamyl-3-carboxy-4-nitroanilide to be superior in terms of stability and solubility.<sup>12</sup> The results correlate with those derived using the original substrate.

In 2002, the International Federation of Clinical Chemistry (IFCC) recommended the standardized method for determining GGT including optimization of substrate concentrations, employment of NaOH, glycylglycine buffer and sample start.<sup>13,14</sup> The GGT liquid reagent follows the formulation recommendation according to Szasz, but was optimized for performance and stability. The assay is optionally standardized against the original IFCC and Szasz methods. The performance claims and data presented here are independent from the standardization.

#### Test principle<sup>15</sup>

Enzymatic colorimetric assay

γ-glutamyltransferase transfers the γ-glutamyl group of L-γ-glutamyl-3-carboxy-4-nitroanilide to glycylglycine.

L-γ-glutamyl-3-carboxy-4-nitroanilide + glycylglycine

L-γ-glutamyl-glycylglycine + 5-amino-2-nitrobenzoate

The amount of 5-amino-2-nitrobenzoate liberated is proportional to the GGT activity in the sample. It is determined by measuring the increase in absorbance photometrically.

#### Reagents - working solutions

**R1** TRIS: 492 mmol/L, pH 8.25; glycylglycine: 492 mmol/L; preservative; additive

**R3** L-γ-glutamyl-3-carboxy-4-nitroanilide: 22.5 mmol/L; acetate: 10 mmol/L, pH 4.5; stabilizer; preservative

R1 is in position B and R3 is in position C.

#### Precautions and warnings

For in vitro diagnostic use for health care professionals. Exercise the normal precautions required for handling all laboratory reagents.

Infectious or microbial waste:

Warning: handle waste as potentially biohazardous material. Dispose of waste according to accepted laboratory instructions and procedures.

Environmental hazards:

Apply all relevant local disposal regulations to determine the safe disposal.

Safety data sheet available for professional user on request.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:



Warning

H317 May cause an allergic skin reaction.

#### Prevention:

P261 Avoid breathing mist or vapours.

P272 Contaminated work clothing should not be allowed out of the workplace.

P280 Wear protective gloves.

#### Response:

P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.

GGT P362 + P364 Take off contaminated clothing and wash it before reuse.

**Disposal:**

P501 Dispose of contents/container to an approved waste disposal plant.

Product safety labeling follows EU GHS guidance.

Contact phone: all countries: +49-621-7590

**Reagent handling**

Ready for use

**Storage and stability**

Shelf life at 2-8 °C: See expiration date on **cobas c** pack label.

On-board in use and refrigerated on the analyzer: 2 weeks

On-board on the Reagent Manager: 1 hour

**Specimen collection and preparation**

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum: Collect serum using standard sampling tubes.

Plasma: Li-heparin and K<sub>2</sub>-EDTA plasma

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

See the limitations and interferences section for details about possible sample interferences.

Stability:<sup>16,17</sup>  
7 days at 15-25 °C  
7 days at 2-8 °C  
1 year at -20 °C (± 5 °C)

Freeze only once.

**Materials provided**

See "Reagents – working solutions" section for reagents.

**Materials required (but not provided)**

See "Order information" section

General laboratory equipment

**Assay**

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

The performance of applications not validated by Roche is not warranted and must be defined by the user.

**Application for serum and plasma****cobas c 701/702 test definition**

Assay type	Rate A	
Reaction time / Assay points	10/25-38	
Wavelength (sub/main)	700/415 nm	
Reaction direction	Increase	
Units	U/L (μkat/L)	
Reagent pipetting	Diluent (H <sub>2</sub> O)	
R1	25 μL	75 μL
R3	20 μL	–
Sample volumes	Sample	Sample dilution

		Sample	Diluent (NaCl)
Normal	3 μL	–	–
Decreased	3 μL	15 μL	150 μL
Increased	6 μL	–	–

**Calibration**

Calibrators	S1: H <sub>2</sub> O S2: C.f.a.s.
Calibration mode	Linear
Calibration frequency	2-point calibration - after reagent lot change - as required following quality control procedures

Calibration interval may be extended based on acceptable verification of calibration by the laboratory.

Traceability: This method has been standardized against the original IFCC formulation (2002)<sup>13</sup> and against the GGT method published by Persijn and van der Slik (1976)<sup>12</sup>, respectively.

Use the appropriate calibrator value for the corresponding application.

**Quality control**

For quality control, use control materials as listed in the "Order information" section.

In addition, other suitable control material can be used.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

**Calculation**

**cobas c** systems automatically calculate the analyte activity of each sample.

Conversion factor: U/L x 0.0167 = μkat/L

**Limitations - interference**

Criterion: Recovery within ± 4 U/L of initial values of samples ≤ 40 U/L and within ± 10 % for samples > 40 U/L

Icterus:<sup>18</sup> No significant interference up to an I index of 50 for conjugated and 20 for unconjugated bilirubin (approximate conjugated bilirubin concentration: 855 μmol/L or 50 mg/dL) and approximate unconjugated bilirubin concentration: 342 μmol/L or 20 mg/dL).

Hemolysis:<sup>18</sup> No significant interference up to an H index of 200 (approximate hemoglobin concentration: 124 μmol/L or 200 mg/dL).

Lipemia (Intralipid):<sup>18</sup> No significant interference up to an L index of 700. There is poor correlation between the L index (corresponds to turbidity) and triglycerides concentration.

Drugs: No interference was found at therapeutic concentrations using common drug panels.<sup>19,20</sup>

In very rare cases, gammopathy, in particular type IgM (Waldenström's macroglobulinemia), may cause unreliable results.<sup>21</sup>

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

**ACTION REQUIRED**

**Special Wash Programming:** The use of special wash steps is mandatory when certain test combinations are run together on **cobas c** systems. All special wash programming necessary for avoiding carry-over is available via the **cobas** link, manual input is required in certain cases. The latest version of the carry-over evasion list can be found with the NaOH/SMS/SmpCln1+2/SCCS Method Sheet and for further instructions refer to the operator's manual.

**Where required, special wash/carry-over evasion programming must be implemented prior to reporting results with this test.**

**Limits and ranges****Measuring range**3-1200 U/L (0.05-20.0  $\mu$ kat/L)

Determine samples having higher activities via the rerun function. Dilution of samples via the rerun function is a 1:11 dilution. Results from samples diluted using the rerun function are automatically multiplied by 11.

**Lower limits of measurement***Lower detection limit of the test*3 U/L (0.05  $\mu$ kat/L)

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying 3 standard deviations above that of the lowest standard (standard 1 + 3 SD, repeatability, n = 21).

Values below the lower detection limit (< 3 U/L) will not be flagged by the instrument.

**Expected values***Standardized against Szasz (Persijn, van der Slik)<sup>22</sup>*

Men	8-61 U/L	0.13-1.02 $\mu$ kat/L
Women	5-36 U/L	0.08-0.60 $\mu$ kat/L

*Standardized against IFCC*Reference Interval Study at 37 °C (corrected in 2005)<sup>22,23</sup>

Men (n = 216)	10-71 U/L	0.17-1.19 $\mu$ kat/L
Women (n = 228)	6-42 U/L	0.10-0.70 $\mu$ kat/L

Consensus values (IFCC)<sup>24</sup>

Men	< 60 U/L	< 1.00 $\mu$ kat/L
Women	< 40 U/L	< 0.67 $\mu$ kat/L

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

**Specific performance data**

Representative performance data on the analyzers are given below. Results obtained in individual laboratories may differ.

**Precision**

Precision was determined using human samples and controls in an internal protocol with repeatability (n = 21) and intermediate precision (3 aliquots per run, 1 run per day, 21 days). The following results were obtained on the cobas c 701 analyzer:

Repeatability	Mean	SD	CV
	U/L ( $\mu$ kat/L)	U/L ( $\mu$ kat/L)	%
Precinorm U	39.2 (0.655)	0.5 (0.008)	1.4
Precipath U	181 (3.02)	1 (0.02)	0.5
Human serum A	18.0 (0.301)	0.6 (0.010)	3.3
Human serum B	472 (7.88)	2 (0.03)	0.3
Human serum C	867 (14.5)	4 (0.1)	0.5
Intermediate precision	Mean	SD	CV
	U/L ( $\mu$ kat/L)	U/L ( $\mu$ kat/L)	%
Precinorm U	44.1 (0.736)	0.8 (0.013)	1.8
Precipath U	221 (3.69)	4 (0.07)	1.7
Human serum 3	46.8 (0.782)	1.5 (0.025)	3.2
Human serum 4	256 (4.28)	9 (0.15)	3.7

Results for intermediate precision were obtained on the cobas c 501 analyzer.

The data obtained on cobas c 501 analyzer(s) are representative for cobas c 701 analyzer(s).

**Method comparison**

$\gamma$ -glutamyltransferase values for human serum and plasma samples obtained on a cobas c 701 analyzer (y) were compared with those determined using the corresponding reagent on a cobas c 501 analyzer (x).

Sample size (n) = 103

Passing/Bablok <sup>25</sup>	Linear regression
y = 1.014x - 0.39 U/L	y = 1.017x - 0.474 U/L
r = 0.992	r = 1.000

The sample activities were between 14.0 and 1168 U/L (0.234 and 19.5  $\mu$ kat/L).

**References**

- Panteghini M. Serum Enzymes. In: Rifai N, Chiu RWK, Young I, Burnham CAD, Wittwer CT, editors. Tietz Textbook of Laboratory Medicine, Saunders Elsevier, Philadelphia, 7th edition, 2023, chapter 32, p. 350-350.e36.
- Rosenberg WMC, Badrick T, Lo SF, et al. Liver Disease. In: Rifai N, Chiu RWK, Young I, Burnham CAD, Wittwer CT, editors. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. 7th ed. 2023; p. 701-763.e21.
- Dillon JF, Miller MH. Gamma glutamyl transferase 'To be or not to be' a liver function test? Ann Clin Biochem 2016 Nov;53(6):629-631.
- Kwo PY, Cohen SM, Lim JK. ACG Clinical Guideline: Evaluation of Abnormal Liver Chemistries. Am J Gastroenterol 2017 Jan;112(1):18-35. doi: 10.1038/ajg.2016.517.
- Newsome PN, Cramb R, Davison SM, et al. Guidelines on the management of abnormal liver blood tests. Gut 2018 Jan;67(1):6-19. doi: 10.1136/gutjnl-2017-314924.
- O'Shea RS, Dasarathy S, McCullough AJ; Practice Guideline Committee of the American College of Gastroenterology. Alcoholic liver disease. Hepatology 2010 Jan;51(1):307-328.
- European Association for the Study of the Liver. Electronic address: easloffice@easloffice.eu; Clinical practice guidelines panel; Wendon J; Panel members; Cordoba J, Dhawan A, Larsen FS, et al.; EASL Governing Board representative; Bernardi M. EASL Clinical Practical Guidelines on the management of acute (fulminant) liver failure. J Hepatol 2017 May;66(5):1047-1081.
- WHO. Guidelines for the prevention, care and treatment of persons with chronic hepatitis B infection [Internet]. 2015 March. [cited 2024 Apr 15]. Available from: <https://www.who.int/publications/i/item/9789241549059>
- National Institute for Health and Care Excellence (NICE) (2013). Hepatitis B (chronic): diagnosis and management. (Clinical guideline [CG165]) [updated 2017 October; cited 2024 Apr 15]. Available from: <https://www.nice.org.uk/guidance/cg165>
- National Institute for Health and Care Excellence (NICE) (2011). Alcohol-use disorders: diagnosis, assessment and management of harmful drinking (high-risk drinking) and alcohol dependence (Clinical guideline [CG115]) [updated 2022 October; cited 2024 Apr 15]. Available from: <https://www.nice.org.uk/guidance/cg115>
- Szasz G. A kinetic photometric method for serum  $\gamma$ -glutamyltransferase. J Clin Chem 1969;15:124-136.
- Persijn JP, van der Slik W. A new Method for the Determination of  $\gamma$ -Glutamyltransferase. J Clin Chem Clin Biochem 1976;4:421.
- Schumann G, Bonora R, Ceriottti F et al. IFCC Primary Reference Procedures for the Measurement of Catalytic Activity Concentrations of Enzymes at 37 °C – Part 6. Reference Procedure for the Measurement of Catalytic Activity Concentrations of gamma-glutamyltransferase. Clin Chem Lab Med 2002;40(7):734-738.
- Klauke R, Schmidt E, Lorentz K. Recommendations for carrying out standard ECCLS procedures (1988) for the catalytic concentrations of creatine kinase, aspartate aminotransferase, alanine aminotransferase and gamma-glutamyltransferase at 37 degrees. Eur J Clin Chem Clin Biochem 1993;31(12):901-909.

# GGT-2

**γ-Glutamyltransferase ver.2 - Standardized against IFCC / Szasz**




- 15 Szasz G, Weimann G, Stähler F, et al. New Substrates for measuring gamma-glutamyl-transpeptidase activity. Z Klin Chem Klin Biochem 1974;12:228-233.
- 16 Szasz G. Methods of Enzymatic Analysis. 2nd English ed. New York. Academic Press, Inc 1974;717.
- 17 Tietz NW, ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia PA: WB Saunders Company 1995;286.
- 18 Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-475.
- 19 Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". Eur J Clin Chem Clin Biochem 1996;34:385-386.
- 20 Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.
- 21 Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007;45(9):1240-1243.
- 22 Abicht K, El-Samalouti V, Junge W, et al. Multicenter evaluation of new GGT and ALP reagents with new reference standardization and determination of 37 °C reference intervals. Clin Chem Lab Med 2001;39:Special Supplement pp S 346.
- 23 Kytzia H-J. Reference intervals for GGT according to the new IFCC 37°C reference procedure. Clin Chem Lab Med 2005;43:A69 [abstract].
- 24 Thomas L, Müller M, Schumann G, et al. Consensus of DGKL and VDGH for interim reference intervals on enzymes in serum. J Lab Med 2005;29(5):301-308.
- 25 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

**Symbols**

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard:

	Contents of kit
	Volume for reconstitution
	Global Trade Item Number
Rx only	For USA: Caution: Federal law restricts this device to sale by or on the order of a physician.

COBAS, NAVIFY, PRECICONTROL, PRECINORM and PRECIPATH are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2024, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim  
www.roche.com

+800 5505 6606

