


REF		CONTENT		Analyzer(s) on which cobas c pack(s) can be used
05168562190*	05168562500	Creatine Kinase-MB (600 tests)	System-ID 03 5924 4	cobas c 701/702
05168562214*	05168562500	Creatine Kinase-MB (600 tests)	System-ID 03 5924 4	cobas c 701/702

Materials required (but not provided):

11447394216	Calibrator f.a.s. CK-MB (3 x 1 mL)	Code 402	
05117003190	PreciControl ClinChem Multi 1 (20 x 5 mL)	Code 391	
05947626190	PreciControl ClinChem Multi 1 (4 x 5 mL)	Code 391	
05117216190	PreciControl ClinChem Multi 2 (20 x 5 mL)	Code 392	
05947774190	PreciControl ClinChem Multi 2 (4 x 5 mL)	Code 392	
05172152190	Diluent NaCl 9 % (119 mL)	System-ID 08 6869 3	

* Some kits shown may not be available in all countries.

English

System information

CKMB: ACN 8060

CKMB2: ACN 8546

Intended use

In vitro test for the quantitative determination of the catalytic activity of creatine kinase MB subunit (CK-MB) in human serum and plasma on **cobas c** systems.

Summary

Measurements of CK-MB, performed with this assay in human serum and plasma, are used as an aid in diagnosis of myocardial infarction.

Creatine kinase (CK) appears as three isoenzymes which are dimers composed of two types of monomer subunits. The isoenzymes comprise all three combinations of monomers, M (for skeletal muscle derived) and B (for brain derived), as represented by the notations MM, MB, and BB.^{1,2}

Many organs contain CK, but the distribution of isoenzymes is different in each one. Skeletal muscle is very rich in the MM isoenzyme, while brain, stomach, intestine, bladder, and lung contain primarily the BB isoenzyme. The MB isoenzyme has been found in appreciable amounts only in myocardial tissue (15 to 20 percent of the total myocardial CK).³ Therefore, total serum CK activity is elevated in a number of diseases. This lack of specificity limits its diagnostic value. However, the striking difference in the CK isoenzyme patterns from different organs has made CK one of the most useful enzymes for diagnostic purposes in acute myocardial infarction. CK-MB appears in serum reflecting its unique presence in myocardial tissue. It is in supporting the diagnosis of suspected myocardial infarction that serial determinations of CK isoenzymes find their most frequent application in the clinical laboratory.^{1,2}

Because of their higher sensitivity and specificity, cardiac troponins, measured by high-sensitivity assays, are the preferred biomarkers to define myocardial infarction,⁴ and if a troponin assay is not available, the best alternative is CK-MB measured by a mass assay.⁴

After immunoinhibition with antibodies to the CK-M subunit,⁵ the CK-B activity is determined with a standardized method for the determination of CK with activation by NAC as recommended by the German Society for Clinical Chemistry (DGKC)⁶ and the International Federation of Clinical Chemistry (IFCC)^{7,8} in 1977 and 2002 respectively. This assay meets the recommendations of the IFCC and DGKC, but was optimized for performance and stability.

Test principle

Immunological UV assay

- Sample and addition of R1 (buffer/enzymes/coenzyme)
- Addition of R2 (buffer/substrate/antibody) and start of reaction.

Human CK-MB is composed of two subunits, CK-M and CK-B which both have an active site. With the aid of specific antibodies to CK-M, the catalytic activity of CK-M subunits in the sample is inhibited to 99.6 % without affecting the CK-B subunits. The remaining CK-B activity, corresponding to half the CK-MB activity, is determined by the total CK method. As the CK-BB isoenzyme only rarely appears in serum and the catalytic activity of the CK-M and CK-B subunits hardly differ, the catalytic activity of the

CK-MB isoenzyme can be calculated from the measured CK-B activity by multiplying the result by 2.

Reagents - working solutions

- R1** Imidazole buffer: 123 mmol/L, pH 6.5 (37 °C); EDTA: 2.46 mmol/L; Mg²⁺: 12.3 mmol/L; ADP: 2.46 mmol/L; AMP: 6.14 mmol/L; diadenosine pentaphosphate: 19 µmol/L; NADP (yeast): 2.46 mmol/L; N-acetylcysteine: 24.6 mmol/L; HK (yeast): ≥ 36.7 µkat/L; G6P-DH (E. coli): ≥ 23.4 µkat/L; preservative; stabilizers; additives.
- R3** CAPSO* buffer: 20 mmol/L, pH 8.8 (37 °C); glucose: 120 mmol/L; EDTA: 2.46 mmol/L; creatine phosphate: 184 mmol/L; 4 monoclonal anti-CK-M antibodies (mouse), inhibiting capacity: > 99.6 % up to 66.8 µkat/L (4000 U/L) (37 °C) CK-M subunit; preservative; stabilizers; additive.

*CAPSO: 3-(cyclohexylamino)-2-hydroxy-1-propanesulfonic acid

R1 is in position B and R3 is in position C.

Precautions and warnings

For in vitro diagnostic use for health care professionals. Exercise the normal precautions required for handling all laboratory reagents.

Infectious or microbial waste:

Warning: handle waste as potentially biohazardous material. Dispose of waste according to accepted laboratory instructions and procedures.

Environmental hazards:

Apply all relevant local disposal regulations to determine the safe disposal.

Safety data sheet available for professional user on request.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:



Danger

H360D May damage the unborn child.

Prevention:

- P201 Obtain special instructions before use.
- P202 Do not handle until all safety precautions have been read and understood.
- P280 Wear protective gloves/ protective clothing/ eye protection/ face protection/ hearing protection.

Response:

P308 + P313 IF exposed or concerned: Get medical advice/attention.

Storage:

P405 Store locked up.

Disposal:

P501 Dispose of contents/container to an approved waste disposal plant.

Product safety labeling follows EU GHS guidance.

Contact phone: all countries: +49-621-7590

Reagent handling

Ready for use

Storage and stability

Shelf life at 2-8 °C:	See expiration date on cobas c pack label.
On-board in use and refrigerated on the analyzer:	4 weeks
On-board on the Reagent Manager:	24 hours

Specimen collection and preparation

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable. Serum: Nonhemolyzed serum is the specimen of choice and also recommended by IFCC.

Plasma: Li-Heparin, K₂-, K₃-EDTA plasma.Li-heparin in the usual concentration does not interfere with the test, but IFCC warns against its use.⁷

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

See the limitations and interferences section for details about possible sample interferences.

Stability in serum: ⁹	8 hours at 20-24 °C
	8 days at 2-8 °C
	4 weeks at -20 °C (± 5 °C)

Freeze only once.

Stability in heparin plasma: ⁹	8 hours at 20-24 °C
	5 days at 2-8 °C
	8 days at -20 °C (± 5 °C)

Freeze only once.

Stability in EDTA plasma: ¹⁰	2 days at 20-25 °C
	7 days at 4-8 °C
	1 year at -20 °C (± 5 °C)

Freeze only once.

Materials provided

See "Reagents – working solutions" section for reagents.

Materials required (but not provided)

See "Order information" section

General laboratory equipment

Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

The performance of applications not validated by Roche is not warranted and must be defined by the user.

Application for serum and plasma**cobas c 701/702 test definition**

Assay type	Rate A		
Reaction time / Assay points	10 / 30-38		
Wavelength (sub/main)	546/340 nm		
Reaction direction	Increase		
Units	U/L (µkat/L)		
Reagent pipetting		Diluent (H ₂ O)	
R1	100 µL	–	
R3	20 µL	–	
Sample volumes	Sample	Sample dilution	
		Sample	Diluent (NaCl)
Normal	5 µL	–	–
Decreased	15 µL	15 µL	120 µL
Increased	10 µL	–	–

Calibration

Calibrators	S1: H ₂ O		
	S2: C.f.a.s. CK-MB		
Calibration mode	Linear		
Calibration frequency	2-point calibration		
	<ul style="list-style-type: none"> • after reagent lot change • as required following quality control procedures 		

Traceability: This method has been standardized against the IFCC Method for Creatine Kinase⁹ with addition of antibodies.**Quality control**

For quality control, use control materials as listed in the "Order information" section.

In addition, other suitable control material can be used.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

Calculation**cobas c** systems automatically calculate the analyte activity of each sample.

Conversion factor: U/L x 0.0167 = µkat/L

Limitations - interferenceThe total CK activity of the specimen should be determined prior to performing the CK-MB assay. The amount of anti-human CK-M subunit antibody in the CK-MB reagent is sufficient for the complete inhibition of up to 4000 U/L CK-M activity. If the total CK activity exceeds 4000 U/L, the specimen requires dilution because complete inhibition of the CK-M subunit is no longer assured. In patients with a disposition to macro-CK formation, implausibly high CK-MB values may be measured in relation to the total CK, since the macroforms mainly consist of CK-B subunits. As these patients have generally not suffered a myocardial infarction, additional diagnostic measures are necessary.¹¹

Criterion: Recovery within ± 10 % of initial value at a creatine kinase-MB activity of ≥ 25 U/L (≥ 0.42 µkat/L).

Icterus:¹² No significant interference up to an I index of 60 for conjugated and 20 for unconjugated bilirubin (approximate conjugated bilirubin concentration: 1026 µmol/L or 60 mg/dL and approximate unconjugated bilirubin concentration: 342 µmol/L or 20 mg/dL).

Hemolysis:¹² No significant interference up to an H index of 20 (approximate hemoglobin concentration: 12.4 µmol/L or 20 mg/dL).

Lipemia (Intralipid):¹² No significant interference up to an L index of 500. There is poor correlation between the L index (corresponds to turbidity) and triglycerides concentration.

Adenylate kinase: Adenylate kinase (AK) may cause positive interference. Sources of AK in the blood are erythrocytes, muscle, and liver. In order to reduce AK interference to a minimum, AMP and Ap₅A are included in the reagent. The AMP/Ap₅A mixture causes 97 % inhibition of the AK from erythrocytes and muscle, and 95 % inhibition of the AK from liver.⁶ The slight residual AK activity does not influence the assay of total CK, but may affect the low CK-MB activities.

Drugs: No interference was found at therapeutic concentrations using common drug panels.^{13,14}

Cyanokit (hydroxocobalamin) and Cefoxitin at therapeutic concentrations interfere with the test.

In very rare cases, gammopathy, in particular type IgM (Waldenström's macroglobulinemia), may cause unreliable results.¹⁵

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

ACTION REQUIRED

Special Wash Programming: The use of special wash steps is mandatory when certain test combinations are run together on **cobas c** systems. All special wash programming necessary for avoiding carry-over is available via the **cobas** link, manual input is required in certain cases. The latest version of the carry-over evasion list can be found with the NaOHD/SMS/SmpCln1+2/SCCS Method Sheet and for further instructions refer to the operator's manual.

Where required, special wash/carry-over evasion programming must be implemented prior to reporting results with this test.

Limits and ranges

Measuring range

3-2000 U/L (0.050-33.4 µkat/L)

Determine samples having higher activities via the rerun function. Dilution of samples via the rerun function is a 1:3 dilution. Results from samples diluted by the rerun function are automatically multiplied by a factor of 3.

Lower limits of measurement

Lower detection limit of the test

3 U/L (0.05 µkat/L)

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying 3 standard deviations above that of the lowest standard (standard 1 + 3 SD, repeatability, n = 21).

Values below the lower detection limit (< 3 U/L) will not be flagged by the instrument.

Expected values

Reference intervals strongly depend on the patient group regarded and the specific clinical situation.

For healthy people: Reference range (37 °C) according to Klein et al.¹⁶ and consensus values¹⁷:

< 25 U/L ($< 0.418 \mu\text{kat/L}$)

For myocardial infarction diagnosis using the combination CK and CK-MB (activity), and representing a CK consensus value based on long-term experience:^{17,18}

- CK_{men} > 190 U/L (3.17 µkat/L)
CK_{women} > 167 U/L (2.79 µkat/L)
- CK-MB > 24 U/L (0.40 µkat/L)
- The CK-MB activity accounts for 6-25 % of the total CK activity.

When myocardial infarction is suspected the diagnostic strategy proposals in the consensus document of European and American cardiologists should in general be followed.¹⁹

If despite the suspicion of myocardial infarction the values found remain below the stated limits, a fresh infarction may be involved. In such cases the determinations should be repeated after 4 hours.

Maximum diagnostic efficiency of the CK-MB determination will be obtained when a sequential sampling protocol is used and consideration is given to the time pattern of activity over a 6 to 48 hour period. When only CK-MB activity is used, the diagnostic efficiency will be lower and will vary with the sampling time.^{1,11}

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

Specific performance data

Representative performance data on the analyzers are given below. Results obtained in individual laboratories may differ.

Precision

Precision was determined using human samples and controls in an internal protocol with repeatability (n = 21) and intermediate precision (3 aliquots per run, 1 run per day, 21 days). The following results were obtained on the **cobas c 701** analyzer:

Repeatability	Mean	SD	CV
	U/L (µkat/L)	U/L (µkat/L)	%
Precinorm CK-MB	35.0 (0.584)	0.5 (0.008)	1.5
Precipath CK-MB	207 (3.46)	1 (0.01)	0.3
Human serum A	19.5 (0.326)	0.5 (0.008)	2.6
Human serum B	660 (11.0)	1 (0.0)	0.2
Human serum C	1811 (30.2)	6 (0.1)	0.3

Intermediate precision

	Mean	SD	CV
	U/L (µkat/L)	U/L (µkat/L)	%
Precinorm CK-MB	32.1 (0.536)	1.1 (0.018)	3.5
Precipath CK-MB	173 (2.89)	1 (0.02)	0.7
Human serum N	20.4 (0.341)	0.9 (0.015)	4.5
Human serum H	184 (3.07)	2 (0.03)	0.9

Results for intermediate precision were obtained on the Roche/Hitachi 917 analyzer.

The data obtained on Roche/Hitachi 917 analyzer(s) are representative for **cobas c 701** analyzer(s).

Method comparison

Creatine kinase-MB values for human serum and plasma samples obtained on a **cobas c 701** analyzer (y) were compared with those determined using the corresponding reagent on a Roche/Hitachi 917 analyzer (x).

Sample size (n) = 129

Passing/Bablok ²⁰	Linear regression
y = 0.988x + 2.23 U/L	y = 0.993x + 2.56 U/L
τ = 0.924	r = 1.000

The sample activities were between 3.00 and 1861 U/L (0.050 and 31.1 µkat/L).

References

- Moss DW, Henderson AR, Kachmar JF. Enzymes. In: Tietz NW, ed. Fundamentals of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders 1987;346-421.
- Lott JA, Stang JM. Serum enzymes and isoenzymes in the diagnosis and differential diagnosis of myocardial ischemia and necrosis. Clin Chem 1980;26:1241-1250.
- Wu AHB (ed.). General Clinical Tests. In: Tietz clinical guide to laboratory tests, 4th ed. Philadelphia, PA: WB Saunders / Elsevier 2006;313.

- 4 Thygesen K, Alpert JS, Jaffe AS, et al. Executive Group on behalf of the Joint European Society of Cardiology (ESC)/American College of Cardiology (ACC)/American Heart Association (AHA)/World Heart Federation (WHF) Task Force for the Universal Definition of Myocardial Infarction. Fourth Universal Definition of Myocardial Infarction. *Glob Heart* 2018;13(4):305-338.
- 5 Würzburg U, Hennrich N, Lang H, et al. Determination of creatine kinase-MB in serum using inhibiting antibodies. *Klin Wschr* 1976;54(8):357-360.
- 6 Bergmeyer HU, Breuer H, Büttner H, et al. Empfehlungen der Deutschen Gesellschaft für Klinische Chemie. Standard-Methode zur Bestimmung der Aktivität der Creatin-Kinase. *J Clin Chem Clin Biochem* 1977;15:249-254.
- 7 Hørder M, Elser RC, Gerhardt W, et al. IFCC methods for the measurement of catalytic concentration of enzymes. Provisional recommendation IFCC method for creatine kinase Appendix A. *J Int Fed Clin Chem* 1990;2:26-35.
- 8 Schumann G, Bonora R, Ceriotti F, et al. IFCC Primary Reference Procedures for the Measurement of Catalytic Activity Concentrations of Enzymes at 37 °C – Part 2. Reference Procedure for the Measurement of Catalytic Concentration of Creatine Kinase. *Clin Chem Lab Med* 2002;40(6):635-642.
- 9 Braun S, Rösenthaller F, Jarausch J, et al. Analyte Stability of CK-MB Activity and cTnT in ICU Patient Serum and Heparin Plasma. Poster presented at Medica 2004, Düsseldorf. (Roche Diagnostics GmbH No. 04587979990).
- 10 Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO Publication WHO/DIL/LAB/99.1 Rev. 2: Jan 2002.
- 11 Remaley AT, Wilding P. Macroenzymes: Biochemical Characterization, Clinical Significance, and Laboratory Detection. *Clin Chem* 1989;35:2261-2270.
- 12 Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. *Clin Chem* 1986;32:470-475.
- 13 Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". *Eur J Clin Chem Clin Biochem* 1996;34:385-386.
- 14 Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. *Ann Clin Biochem* 2001;38:376-385.
- 15 Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. *Clin Chem Lab Med* 2007;45(9):1240-1243.
- 16 Klein G, Berger A, Bertholf R, et al. Abstract: Multicenter Evaluation of Liquid Reagents for CK, CK-MB and LDH with Determination of Reference Intervals on Hitachi Systems. *Clin Chem* 2001;47:Suppl. A30.
- 17 Thomas L, Müller M, Schumann G, et al. Consensus of DGKL and VDGH for interim reference intervals on enzymes in serum. *J Lab Med* 2005; 29(5):301-308.
- 18 Stein W. Strategie der klinischen-chemischen Diagnostik des frischen Myokardinfarktes. *Med Welt* 1985;36:572-577.
- 19 Myocardial Infarction Redefined - A Consensus Document of the Joint European Society of Cardiology/ American College of Cardiology Committee for the Redefinition of Myocardial Infarction. *Eur Heart J* 2000;21:1502-1513.
- 20 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. *J Clin Chem Clin Biochem* 1988 Nov;26(11):783-790.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

The Summary of Safety & Performance Report can be found here:
<https://ec.europa.eu/tools/eudamed>

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard (for USA: see navifyportal.roche.com for definition of symbols used):

CONTENT

Contents of kit



Volume for reconstitution

GTIN

Global Trade Item Number

Rx only

For USA: Caution: Federal law restricts this device to sale by or on the order of a physician.

COBAS, NAVIFY, PRECICONTROL, PRECINORM and PRECIPATH are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2023, Roche Diagnostics

CE 0123

Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com

+800 5505 6606

