Instructions for Use / Instrucciones de uso / Instructions d'utilisation

CINtec® PLUS Kit

The CINtec® PLUS Kit is an immunocytochemistry assay for the simultaneous qualitative detection of the p16^{INK4a} and Ki-67 proteins in cervical cytology preparations. The kit is intended to be used in medical cytology laboratories. Interpretation of the test results may only be made by a certified professional in conjunction with the patient's clinical history and additional diagnostic tests that have been performed.

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<https://dialog.roche.com>

Distributed by:

Ventana Medical Systems, Inc.
1910 E. Innovation Park Drive
Tucson, Arizona 85755
USA

www.roche.com

	9537
	06595383001
	04015630984374
	
	50
	2 – 8°C

Table of Contents / Índice / Table des matières

ENGLISH

I.	Product Name	4
II.	Intended Use.....	4
III.	Summary and Explanation of the Device	4
	Background	4
	Principle of Procedure	5
IV.	Reagents.....	5
	Materials Provided	5
	Storage	8
	Materials and Reagents Required but not Provided	9
	Equipment Required	10
V.	Warnings and Precautions	10
	Warning	10
	Caution	11
VI.	Procedures	11
	Cytological Specimen Preparations	11
	Heat-induced epitope retrieval treatment of specimens prior to staining	13
	Slide Staining Procedures.....	13
	1. Reagent Preparation.....	13
	1.1 Epitope Retrieval Solution.....	13
	1.2 Wash Buffer	13
	1.3 Substrate-Chromogen Solutions (DAB and Fast Red)	14
	1.4 Counterstain	14
	1.5 Mounting Medium	14
	2. Staining Procedure	15
	2.1 Specimen rehydration	15
	2.2 Epitope Retrieval	16
	2.3 Staining Protocols.....	17
	2.3.1 Staining Protocol for Autostainer Instruments (Dako, LabVision)	17
	2.4 Counterstaining with Hematoxylin	19
	2.5 Mounting.....	19
VII.	Quality Control	20
	Positive Control	20
	Negative Control	20
	Assay Verification	20
VIII.	Interpretation of Results	21
IX.	Limitations	21
X.	Troubleshooting	23
XI.	Symbols	26
XII.	Manufacturer	26
XIII.	Revision Status	26
XIV.	Intellectual Property	26
	Annex 1: References	80

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II. Intended Use

For in vitro diagnostic use.

The CINtec® *PLUS* Kit is an immunocytochemistry assay for the simultaneous qualitative detection of the p16^{INK4a} and Ki-67 proteins in cervical cytology preparations.

The kit is intended to be used in medical cytology laboratories. Interpretation of the test results may only be made by a certified professional in conjunction with the patient's clinical history and additional diagnostic tests that have been performed.

III. Summary and Explanation of the Device

Background

In eukaryotic cells, control of progression of the cell division cycle is effected by a complex pattern of controlled expression and post-translational modifications (e.g., phosphorylation) of cell-cycle regulating proteins. The p16^{INK4a} protein plays a major role in this mechanism of regulation of the eukaryotic cell cycle. It is part of the retinoblastoma protein (pRb)-mediated control of the G₁-S-phase transition, and it triggers cell cycle arrest in the course of cellular differentiation processes. Thus, p16^{INK4a} provides an anti-proliferative effect when expressed during regular cell cycle progression [1]. In terminally differentiated epithelial cells, p16^{INK4a} expression is down-regulated to levels typically not detectable by immunocytochemistry [2].

In cervical dysplasia, overexpression of p16 is regarded as a surrogate biomarker for transforming HPV infections, reflecting the activation of HPV E6/E7 oncoprotein-driven cell proliferation [2-5]. Detection of p16 in cervical cytology preparations has been proposed as a valuable adjunctive marker to triage women with abnormal Pap cytology results as well as with positive HPV test results [3-5]. However, because p16-specific staining may be observed in individual metaplastic or endocervical cells in which p16 may be expressed to exert its physiological normal, growth-suppressive cellular function, interpretation of p16 single-stained cervical cytology preparations requires identification of p16 immunoreactive cells and further classification of these cells regarding signs of morphologic abnormalities [2-4].

The combined simultaneous detection of p16 and the proliferation marker Ki-67 within the same cell by ICC has been shown to be a valuable tool to identify dysplastic cervical cells in cytology preparations without the need for morphologic interpretation [3,6,7]. Ki-67 is a nuclear and nucleolar protein strictly associated with cell proliferation and is undetectable by standard immunostaining methods in resting (G₀) cells [8]. Under normal physiologic conditions, expression of the proliferation-associated protein Ki-67 is mutually exclusive of the anti-proliferative protein p16. In contrast, cells where the retinoblastoma protein (pRB)-mediated pathway controlling the cell-cycle progression is abrogated upstream of the tumor suppressor function of p16 (such as in epithelial cells expressing the high-risk HPV E6/E7 oncoproteins) may proliferate and thus may express Ki-67 in the presence of functional p16 [2,3].

Therefore, the detection of individual cells in cervical cytology preparations that simultaneously co-express p16 and Ki-67 may serve as a morphology-independent indicator of cells with cell cycle dysregulation and thus may be used as an indicator of the presence of transforming HPV infections and underlying cervical intraepithelial neoplasia [2,3]. In the recent past, numerous studies have been performed and published evaluating the potential value and clinical usefulness of p16/Ki-67 dual-stained cytology for the identification of women that may benefit from referral to colposcopy based on various primary cervical cancer screening results, including for the triage of women with Atypical Squamous Cells of Undetermined Significance (ASC-US) or Low grade Squamous Intraepithelial Lesion (LSIL) Pap cytology results, women who are high-risk HPV positive in primary HPV screening, or women who are Negative for Intraepithelial Lesion or Malignancy (NILM)/HPV positive in clinical settings where Pap cytology/HPV co-testing is used for primary screening [6,7;9-25].

Principle of Procedure

The CINtec® *PLUS* Kit contains a set of reagents for the immunocytochemical detection of the p16^{INK4a} and Ki-67 antigens. The kit is designed to perform a two-step immunocytochemical staining procedure for cytological specimens obtained from the cervix uteri. For the detection of the antigens, a primary monoclonal mouse antibody clone E6H4™ directed to human p16^{INK4a} protein and a primary monoclonal rabbit antibody clone 274-11 AC3V1 directed against human Ki-67 protein are used.

Ready-to-use visualization reagents comprising (i) a polymer reagent conjugated to horseradish peroxidase and goat anti-mouse Fab' antibody fragments, and (ii) a polymer reagent conjugated to alkaline phosphatase and goat anti-rabbit Fab' antibody fragments are used. The visualization reagents have been subjected to solid-phase absorption to eliminate cross reactivity with human immunoglobulins.

The chromogen reactions are based on horseradish peroxidase-mediated conversion of a DAB chromogen, and alkaline phosphatase-mediated conversion of Fast Red chromogen to visible reaction products at the respective antigen sites. After counterstaining, a two step mounting protocol must be applied: in a first step, an aqueous mounting of the specimens using an aqueous mounting medium provided with the kit must be used. Subsequently, the slides may be cover-slipped using e.g. a permanent mounting medium. The results can be evaluated by light microscopy inspection.

IV. Reagents

Materials Provided

The materials listed below are included in each kit and are sufficient to perform 50 tests. The number of tests is based on the use of 200 µL of the reagents per slide.

1 Peroxidase Blocking Reagent

Peroxidase Blocking Reagent

11.5 mL, ready-to-use

3% hydrogen peroxide, containing 15 mmol/L sodium azide (NaN₃).

2 Primary Antibodies Solution

Primary Antibodies Solution

11.5 mL, ready-to-use

Monoclonal mouse anti-Human p16^{INK4a} antibody, Clone E6H4™, and monoclonal rabbit anti-Human Ki-67 antibody Clone 274-11 AC3, supplied in 50 mM Tris buffer pH 7.2, containing 15 mmol/L sodium azide (NaN₃) and stabilizing protein.

3 Visualization Reagent HRP

Visualization Reagent HRP

11.5 mL, ready-to-use

Warning

H317 May cause an allergic skin reaction.

P261 Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.

P272 Contaminated work clothing should not be allowed out of the workplace.

P280 Wear protective gloves.

P333 + P313 If skin irritation or rash occurs: Get medical advice/ attention.

P362 + P364 Take off contaminated clothing and wash it before reuse.

P501 Dispose of contents/ container to an approved waste disposal plant.

Contains

26172-54-3 2-methyl-2H-isothiazol-3-one hydrochloride

55965-84-9 reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1)

Polymer reagent conjugated with horseradish peroxidase and affinity purified goat anti-Mouse Fab' antibody fragments, supplied in stabilizing solution comprising preservatives and stabilizing protein.

Contains 5-bromo-5-nitro-1,3-dioxane that may produce an allergic reaction.

4 Visualization Reagent AP

Visualization Reagent AP

11.5 mL, ready-to-use

Warning

H317 May cause an allergic skin reaction.

P261 Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.

P272 Contaminated work clothing should not be allowed out of the work-

place.

P280 Wear protective gloves.

P333 + P313 If skin irritation or rash occurs: Get medical advice/ attention.

P362 + P364 Take off contaminated clothing and wash it before reuse.

P501 Dispose of contents/ container to an approved waste disposal plant.

Contains 26172-54-3 2-methyl-2H-isothiazol-3-one hydrochloride

Polymer reagent conjugated with alkaline phosphatase and affinity purified goat anti-rabbit Fab' antibody fragments, supplied in stabilizing solution comprising preservatives and stabilizing protein.

Contains 5-bromo-5-nitro-1,3-dioxane that may produce an allergic reaction.

5

DAB Buffered Substrate

DAB Buffered Substrate

16.0 mL,

Substrate buffer solution, pH 7.5, containing < 0.1% hydrogen peroxide, stabilizers and enhancers.

6

DAB Chromogen

DAB Chromogen

0.85 mL, 3,3'-diaminobenzidine chromogen solution.

Danger

H314 Causes severe skin burns and eye damage.

H341 Suspected of causing genetic defects.

H350 May cause cancer.

P201 Obtain special instructions before use.

P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.

P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.

P304 + P340 + P310 IF INHALED: Remove person to fresh air and keep comfortable for breathing. Immediately call a POISON CENTER/doctor.

P305 + P351 + P338 + P310 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/doctor.

P308 + P313 IF exposed or concerned: Get medical advice/attention.

Contains 868272-85-9 3,3'-Diaminobenzidine tetrahydrochloride hydrate

NOTE: Consult Federal, State, or local regulations for disposal.

7 **Naphthol Phosphate Substrate**

Naphthol Phosphate Substrate

25.0 mL,

Substrate buffer solution, pH 9.2 containing Naphthol AS-TR Phosphate Substrate, stabilizers and enhancers.

8 **Fast Red Chromogen**

Fast Red Chromogen

1.33 mL, Fast Red chromogen solution.

Danger

H314 Causes severe skin burns and eye damage.

P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.

P301 + P330 + P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.

P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.

P304 + P340 + P310 IF INHALED: Remove person to fresh air and keep comfortable for breathing. Immediately call a POISON CENTER/doctor.

P305 + P351 + P338 + P310 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/doctor.

P501 Dispose of contents/ container to an approved waste disposal plant.

9 **Epitope Retrieval Solution 10 x**

Epitope Retrieval Solution 10 x

500 mL, 100 mM Tris, 10 mM EDTA, pH 9.0, containing 15 mmol/L sodium azide (NaN_3).

10 **CINtec® PLUS Mount**

CINtec® PLUS Mount

18.0 mL, aqueous based permanent mounting medium for the permanent preservation of slide preparations stained with peroxidase and alkaline phosphatase based visualization systems. Contains 7.7 mmol/L sodium azide (NaN_3)

Storage

Store at 2 – 8 °C. Do not use after the expiration date. No data has been generated respective to the storage of the reagents under any conditions other than those stated above.

After opening, kit components are stable for 6 months if stored at 2 – 8 °C. Solutions must be discarded if cloudy in appearance.

The CINtec® *PLUS* Mount vial (Vial 10) should be removed from the kit box after the kit has been opened for the first time and may be further stored at ambient temperature (2 – 30°C).

Diluted Wash Buffer and diluted Epitope Retrieval Solution are stable for up to one month if stored at 2 – 8 °C to reduce viscosity and to facilitate usage. Solutions must not be used if cloudy in appearance.

Materials and Reagents Required but not Provided

CINtec® Wash Buffer 10X to be used with the CINtec® *PLUS* Kit is available under catalog number 06595421001 (8557) from Roche but is not included in the kit. For order details please refer to the website www.roche.com.

500 mmol/L Tris buffer solution with 1.5 mol/L NaCl, pH 7.6, containing detergent and an antimicrobial agent.

Warning

H317 May cause an allergic skin reaction.

P261 Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.

P272 Contaminated work clothing should not be allowed out of the workplace.

P280 Wear protective gloves.

P333 + P313 If skin irritation or rash occurs: Get medical advice/ attention.

P362 + P364 Take off contaminated clothing and wash it before reuse.

P501 Dispose of contents/ container to an approved waste disposal plant.

Contains 55965-84-9 Mixture of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H -isothiazol-3-one (3:1)

Alcohol-free hematoxylin for counterstaining;

Distilled or deionized water;

PreservCyt® Solution; may be required to prepare ThinPrep® slides from residual sample volume remaining in the original ThinPrep® vial. Sample volume must be at least 17 ml to prepare another slide according to the ThinPrep® 2000 System Operator's Manual, Instrument Troubleshooting, pg. 6.27;

SurePath™ Preservative Fluid; is required for preparing SurePath slides intended for CINtec *PLUS* Cytology staining. See section VI. Procedures, Cytological Specimen Preparations for details;

Ethanol, 99%

Xylene-based mounting medium for glass coverslipping;

Slides: SuperFrost® *PLUS*, ThinPrep® Microscopy Slides; BD SurePath™ PreCoat slides

Xylene;

Glass or film coverslips;

Process controls e.g. ThinPrep® slides produced from the residual specimen vial:

- Positive control: ThinPrep® slide from (pooled) specimen(s) with Pap

cytology result confirmed as High-grade squamous intraepithelial lesion (HSIL);

- Negative control: ThinPrep® slide from (pooled) specimen(s) with Pap cytology result confirmed as Negative for intraepithelial lesion or malignancy (NILM).

Equipment Required

Light microscope (10 - 40x objective magnification);

Heat resistant staining jars (plastic);

Measuring cylinders;

Squirt bottle (to be filled with wash buffer);

Timer (capable of 30 seconds - 60 minutes intervals);

Water bath with lid (capable of maintaining Epitope Retrieval Solution at 95 – 99°C);

Thermometer;

Optional: Drying oven capable of maintaining a temperature of 37 or 60°C;

Optional: Drying oven capable of maintaining a temperature of 37 or 60°C;

Optional: Dako or LabVision Autostainer instrument;

V. Warnings and Precautions

Warning

1. Caution! Some of the reagents comprised in this kit contain hazardous chemicals. When handling the components of this kit, adhere to safety precautions for handling hazardous laboratory reagents.
2. Components 1, 2, 9 and 10 of this product contain sodium azide (NaN_3), which is highly toxic in its pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. Components 2, 3 and 4 contain material of animal origin. Adhere to proper handling procedures as is applicable to any product derived from biological sources.
4. Safety Data Sheet for the kit is available upon request.
5. When handling and disposing of cytological specimens, including all specimens before and after fixation, as well as all materials exposed to them, adhere to the safety precautions for handling potentially infectious material as well as applicable waste disposal requirements.
6. Never pipette reagents by mouth. Avoid contacting the skin and mucous membranes with reagents and specimens. In the event that reagents or specimens come in contact with skin or mucous membranes, wash with copious amounts of water.
7. The Visualization reagents, DAB Chromogen and Fast Red Chromogen may be affected adversely if exposed to excessive light levels. Do not store kit components or perform staining in strong light, such as direct sunlight.

8. Wear appropriate personal protective equipment to avoid contact with eyes and skin when handling any of the components included or to be used in conjunction with the CINtec® PLUS Kit. Refer to the Safety Data Sheet for additional information.
9. Product safety labeling primarily follows EU GHS guidance.

Caution

1. For in vitro diagnostic use.
2. For professional use only.
3. Minimize microbial contamination of reagents to avoid non-specific staining.
4. Incubation times, temperatures, or methods other than those specified may give erroneous results.
5. Do not use the kit if the packaging of any of its components is damaged. Should packaging be compromised or components damaged, please notify the manufacturer without delay.
6. Disposal of all waste materials must be in accordance with local guidelines and regulations.
7. All reagents are formulated specifically for use with this test. In order for the test to perform as specified, no substitutions should be made.
8. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the CINtec® PLUS Kit is suspected, immediately refer to the contact information provided in section XII. for further information on technical support.
9. Malfunction of the product due to handling problems or to instability does not result in obvious signs. Therefore, as a quality control measure, positive and negative controls should be run simultaneously with patient specimens.

VI. Procedures

Cytological Specimen Preparations

Cytological specimens must be adequately handled to preserve the specimens for immunocytochemical procedures. All specimens should be subjected to standard methods of cell processing.

ThinPrep® (Hologic™ Inc.) slides prepared on a ThinPrep® 2000 Processor (Hologic™ Inc.) according to the manufacturer's protocol, BD SurePath™ (BD Diagnostics) slides prepared according to the manufacturer's protocol, as well as manually prepared slides (conventional smears) are suitable for use.

Notes:

1. ThinPrep® sample preparation:

Please note that the use of the ThinPrep® 3000 Processor is not recommended, as the spray fixation procedure provided by the instrument may lead to substantial cell loss.

2. BD SurePath™ sample preparation:

As it has been occasionally reported that storing the residual cell material

after the enrichment process in water may have a negative impact on the intensity of the immunocytochemical signal, we strongly recommend to follow the instructions outlined below when preparing slides for CINtec® PLUS testing to avoid any risk of signal loss:

- a. Slide preparation **directly** after processing the Pap slide
 - i. A second slide for each case can be prepared immediately after the Pap staining of the corresponding BD SurePath™ slide is completed.
 - ii. Place a second set of labeled slides into slide racks.
 - iii. On the PrepStain™ using the GYN version 1.1 or 1.2, Slide Preparation, select program „Transfer only“.
- b. Slide preparation **from the enriched cell pellet** after preparation of Pap slide at a later time
 - i. Remove the tube racks from the PrepStain™ System and add approximately 2 mL of BD SurePath™ Preservative Fluid to each tube.
 - ii. Cap each tube and store at room temperature for up to 4 weeks or refrigerated (2 – 8°C) for up to 6 months.
 - iii. If a slide for CINtec® PLUS is to be made on a stored sample, allow the sample to adapt to room temperature for 1 hour. Start with the second centrifugation step of the GYN Enrichment Process and carry through the remainder of the Prep procedure.
- c. Slide preparation **from left-over sample material** remaining in the original sample vial (approximately 2 mL)
 - i. Add 8mL of BD SurePath™ Preservative Fluid to the residual sample in the BD SurePath™ vial (approx. 2 mL).
 - ii. The diluted sample should be processed on the PrepMate™ using the standard technique and on the PrepStain™ using the GYN version 1.1 or 1.2, Slide Preparation, program „Transfer only“.

Immediately after preparation, ThinPrep® or BD SurePath™ slides should be fixed in 99% ethanol for 10 minutes to 1 hour and air-dried for 20 minutes to 16 hours (over night). **ThinPrep® or BD SurePath™ cytology preparation must not be fixed with cytological spray fixation reagent containing polyethylene glycol (e.g. Merckofix®, Merck).**

Conventional smears should be fixed with cytological spray fixation reagent containing polyethylene glycol (e.g. Merckofix®, Merck) immediately after sample collection.

Before beginning the immunostaining procedure, all specimens must be rehydrated following specific protocols described in Section 2.1.

NOTE: Ethanol with minimal impurities is needed for fixation of ThinPrep® liquid based cytology slides to prevent background staining. During the preparation of ThinPrep® slides that are intended to be used for immunostaining renew the ethanol bath after fixation of every 25 ThinPrep® slides.

Heat-induced epitope retrieval treatment of specimens prior to staining

Heat-induced epitope retrieval (HIER) using the Epitope Retrieval Solution included in the kit is necessary for optimal assay performance. Deviation from the described procedure may affect results!

For heat-induced epitope retrieval, the cytological specimen must be heated by immersing into the Epitope Retrieval Solution in a calibrated water bath capable of maintaining the Epitope Retrieval Solution at a temperature of 95 – 99°C. Laboratories located at higher elevations should determine the best method of maintaining the required water bath temperature. Manufacturer does not recommend any deviation from the procedure described herein.

After heat-induced epitope retrieval, the cytological specimen must be kept at ambient temperature for 20 minutes or longer until the temperature has cooled down to 50°C or below before further processing. Thereafter, staining of the cytological specimen must be performed without delay.

Slide Staining Procedures

The kit contains reagent volume which is sufficient to perform 50 tests. The number of tests is based on the use of 200 µL of the reagents per slide.

200 µL per slide is recommended to process ThinPrep® or BD SurePath™ cytological preparations on Dako or LabVision Autostainer instruments (two drop zones per slide with 100 µL each).

To ensure complete coverage of conventional smears on Dako or LabVision Autostainer instruments the use of 300 µL per slide (100 µL for all three available drop zones) is recommended.

1. Reagent Preparation

It is recommended to prepare the following reagents, except the Fast Red working solution, before starting with the staining procedure. All reagents should be equilibrated to room temperature (20 – 25°C) prior to immunostaining.

1.1 Epitope Retrieval Solution

Prepare the amount of Epitope Retrieval Solution sufficient for the staining procedure that is planned by dilution of a quantity of Vial 9 (Epitope Retrieval Solution 10 x) 1:10 using distilled or deionized water.

After dilution, the Epitope Retrieval Solution may be stored at 2 – 8°C for up to one month. The diluted solution must be discarded if it is cloudy in appearance.

NOTE: Use of water with elevated levels of ions for dilution of the epitope retrieval solution may significantly reduce the staining performance of the test. Please make sure that the water used is properly distilled or deionized (i.e. ensure that your ion exchange column for producing deionized water has been checked by routine maintenance). **Do NOT use tap water!**

1.2 Wash Buffer

Use Wash Buffer (10 x), Catalog number 06595413001 (8557), provided by Roche mtm laboratories AG in combination with the CINtec® PLUS Kit. For order details please refer to the website www. Roche.com.

Prepare an amount of Wash Buffer sufficient for the wash steps of the staining procedure that is planned by diluting of a quantity of the Wash Buffer (10 x) 1:10 using distilled or deionized water.

After dilution the Wash Buffer may be stored at 2 – 8°C for up to one month. The diluted solution must be discarded if it is cloudy in appearance.

1.3 Substrate-Chromogen Solutions (DAB and Fast Red)

Ensure that both Chromogenes and Substrates are equilibrated to room temperature (20 – 25°C). Addition of excess Chromogen to the Substrate will result in deterioration of the positive signal.

A) Preparation of the DAB working solution prior to the staining run

The DAB working solution is stable for 8 hours after preparation.

Prepare the DAB working solution as follows:

- i) Transfer 1 mL of the DAB Substrate Solution (Vial 5) into a clean reaction tube;
- ii) Add one drop (25 – 30 µL) of the DAB Chromogen (Vial 6) and mix gently by inverting the tube (do not vortex);
- iii) When using the Autostainer, transfer the DAB working solution into an Autostainer Reagent Vial before starting the staining run and load it on the appropriate Autostainer rack position.

B) Preparation of the Fast Red working solution directly before use

Prepare the Fast Red working solution directly before use, as otherwise decreased staining intensity and thus loss of sensitivity may result. Do not vortex the Fast Red working solution as this may result in formation of precipitates.

Prepare the Fast Red working solution as follows:

- i) Transfer 1 mL of the Naphthol Phosphate Substrate Solution (Vial 7) into a clean reaction tube;
- ii) Add one drop (40 – 45 µL) of the Fast Red Chromogen (Vial 8) and mix gently by inverting the tube (do not vortex);
- iii) When using the Autostainer, prepare the Fast Red working solution when requested by the Autostainer software. Then transfer the working solution into an Autostainer Reagent Vial and load it on the appropriate Autostainer rack position. Avoid delay of the “substrate batch” step during the Autostainer run to minimize the risk of drying artefacts.

1.4 Counterstain

The DAB and Fast Red staining reactions result in water insoluble coloured end-product (DAB: brown; Fast Red: red).

NOTE: Alcohol-free hematoxylin must be used for counterstaining, as the use of alcohol may negatively impact the intensity of the signal generated by the Fast Red, or remove it completely. If used, adhere to the instructions provided by the supplier of the hematoxylin for performing the counterstaining.

1.5 Mounting Medium

For mounting slide specimens after staining a two-step mounting procedure is

needed as described below.

First CINtec® *PLUS* Mount (Vial 10), an aqueous-based permanent mounting medium, is applied manually in a thin coat and allowed to dry (“liquid coverslipping”). In a second step a glass coverslip is applied onto the dried surface of the aqueous mounting media using xylene-based mounting medium.

Equilibrate CINtec® *PLUS* Mount at room temperature before use and store it at ambient temperature for further use.

2. Staining Procedure

The CINtec® *PLUS* Kit has been validated for use on Autostainer Instruments (LabVision Autostainer 480 or Dako Autostainer PLUS) according to the template outlined below (see Section 2.3.1). It may be possible to use other instruments or systems with comparable function after appropriate validation by the user.

Prior to staining, the specimens and reagents should be prepared as stated in Sections 1.1-1.3 and 2.1.

All reagents should be equilibrated to ambient temperature (20 – 25°C) prior to performing the immuno-staining procedure.

Likewise, all steps should be performed at ambient temperature. Do not allow specimens to dry during the staining procedure. Dried specimens may display increased non-specific staining. If prolonged incubations are used, place specimens in a humid environment.

2.1 Specimen rehydration

For all cytological specimens a rehydration step is necessary prior to staining. This step should be performed at ambient temperature (20 - 25°C).

A) ThinPrep® liquid based cytology slides & conventional smears

- Place slides in distilled or deionized water and incubate for 10 (±3) minutes;
- Commence staining procedure as outlined in Section 2.2, Step 1: Epitope retrieval.

B) BD SurePath™ cytology slides

Alcohol fixed BD SurePath™ cytology slides have to undergo a special rehydration procedure to make specially coated glass slides compatible with liquid coverslip medium.

- Make sure that the slides have been completely air-dried;
- Place slides in a fresh bath of 100% xylene (use xylene-resistant jars); dip the slides a couple times and incubate for 2 minutes;
- Transfer the slides into a fresh bath of 99% ethanol and incubate for 2 minutes;
- Transfer the slides into a fresh bath of 70% ethanol and incubate for 2 minutes;
- Transfer the slides in distilled or deionized water and incubate for 2 minutes

Note: Replace the baths used for the dehydration after processing of 50 slides or once a week if a smaller quantity of slides is processed.

2.2 Epitope Retrieval

A) ThinPrep® liquid based cytology slides & conventional smears:

- Fill heat resistant staining jars (plastic) with the diluted Epitope Retrieval Solution (see Procedure, Section 1.1);
- Place staining jars containing Epitope Retrieval Solution in water bath and heat water bath and the Epitope Retrieval Solution to 95 – 99°C. The temperature has to be measured inside the staining jars. It is important to adjust the level of the water in the water bath to make sure that the jars are immersed in the water to a level of 80%. To stabilize the temperature and avoid evaporation, cover jars with lids;
- When the temperature of 95 – 99°C has been reached immerse the re-hydrated cytology slides into the preheated Epitope Retrieval Solution in the staining jars; this step usually will lower the temperature in the jars to less than 90°C. Keep the temperature probe in the jars and close the lid of the staining jars as good as possible.
- Bring the temperature of the water bath and the Epitope Retrieval Solution in the jars back to 95 – 99°C; check the temperature of the Epitope Retrieval Solution in the jars;
- Incubate for 10 (±1) minutes at 95 – 99°C; start counting down only after the temperature of the Epitope Retrieval Solution in the jars has been verified to have reached a temperature of 95 – 99°C;
- Remove the entire jar with slides from the water bath;
- Remove the lid off the staining jars and allow the slides to cool in the Epitope Retrieval Solution for 20 (±1) minutes at ambient temperature until it has reached 50°C or below;
- Transfer the slides into a staining jar filled with Wash Buffer (see Procedure, Section 1.2) and incubate for 5 (±1) minutes prior to loading the slides onto the programmed Autostainer instrument.

B) BD SurePath™ cytology slides

- Fill heat resistant staining jars (plastic) with the diluted Epitope Retrieval Solution (see Procedure, Section 1.1). Please note that when using jars that are not made of heat-resistant plastic, but e.g. metal or heat-resistant glass the time for the epitope retrieval must be modified on an individual basis;
- Immerse slides into the cold Epitope Retrieval Solution; close the lid of the staining jar;
- Place the staining jar in a cold water bath and bring the temperature of the water bath and the Epitope Retrieval Solution to 95 – 99°C and incubate for 15 minutes as soon as the temperature has been reached; please note that the time period for heating up the solution to 95 – 99°C should be no longer than 40 minutes at maximum;
- Remove the entire jar with slides from the water bath;
- Remove the lid off the staining jars and allow the slides to cool in the Epitope Retrieval Solution for 20 (±1) minutes at ambient temperature until it has reached 50°C or below;
- Transfer the slides into a staining jar filled with Wash Buffer (see Procedure,

Section 1.2) and incubate for 5 ($\pm$ 1) minutes prior to loading the slides onto the programmed Autostainer instrument.

NOTE: The Epitope Retrieval Solution is designed for single use application only. Do not re-use.

2.3 Staining Protocols

2.3.1 Staining Protocol for Autostainer Instruments (Dako, LabVision)

Prior to the first application of the CINtec® *PLUS* Kit on an Autostainer Instrument, a new template needs to be set up and the CINtec® *PLUS* Kit reagents need to be added to the Autostainer “Reagent List”. Please refer to the Operator’s Manual for the dedicated Autostainer Instrument.

The use of 200 μ L per slide is recommended to process ThinPrep® or BD SurePath™ cytological preparations. For conventional smears on Dako or LabVision Autostainer instruments the use of 300 μ L of the reagents per slide may be required.

The following is an outline of the program run:

Program Step	ThinPrep® cytological preparations & Conventional Smears	BD Surepath™ cytological preparations
1	rinse*	rinse*
2	Peroxidase-Blocking Reagent 5 minutes	Peroxidase-Blocking Reagent 5 minutes
3	rinse*	rinse*
4	Primary Antibodies Solution 30 minutes	Primary Antibodies Solution 30 minutes
5	rinse*	rinse*
6	Visualization Reagent HRP 15 minutes	Visualization Reagent HRP 15 minutes
7	rinse*	rinse*
8	rinse*	rinse*
9	rinse*	rinse*
10	Visualization Reagent AP 15 minutes	Visualization Reagent AP 15 minutes
11	rinse*	rinse*
12	rinse*	rinse*
13	rinse*	rinse*
14	switch	switch
15	“Substrate” step: DAB 10 minutes	“Substrate” step: DAB 10 minutes
16	rinse with distilled or deionized water	rinse with distilled or deionized water
17	rinse*	rinse*
18	“Substrate-batch” step: Fast Red 15 minutes	“Substrate-batch” step: Fast Red 15 minutes
19	rinse*	rinse*
20	---	For LabVision Autostainer: “Substrate-batch” step: Fast Red 15 minutes For Dako Autostainer: “Substrate” step: Fast Red 15 minutes
21	----	rinse*
22	rinse with distilled or deionized water	rinse with distilled or deionized water
23	switch	switch

*use Wash Buffer for the respective rinsing steps.

NOTE: If the Autostainer Instrument that is used for the staining procedure rinses slides with buffer, the slides must be rinsed with distilled or deionized water after they have been removed from the Autostainer.

Subsequent to programming proceed as follows:

- Transfer the reagents from the kit bottles into Autostainer Reagent Vials. Use the Autostainer-generated map for program times and reagent volumes;
- Place the Autostainer reagent vials into the Autostainer Reagent Rack according to the “Reagent Layout Map”;
- Load the slides onto the Autostainer according to the “Slide Layout Map” and rinse the slides with Wash Buffer to prevent specimen from drying out.

2.4 Counterstaining with Hematoxylin

- Immerse slides in a bath of **alcohol-free hematoxylin**. Incubate for 2 - 10 minutes, depending on the potency of the alcohol-free hematoxylin used;

NOTE: The use of alcohol-free hematoxylin is absolutely mandatory.

- For bluing place slides in a tap water bath or alkaline solution such as a weak ammonia solution (NH₄OH, 0.08% in deionized water) and incubate for 30 seconds - 2 minutes;
- To ensure that all residual hematoxylin has been cleared, subsequently replace the tap water in the jar several times (3 – 5x) with fresh tap water until no stain residue can be observed any more;
- Briefly incubate slides in a bath of distilled or deionized water.

NOTE: Depending on the incubation length and potency of the hematoxylin used, counterstaining will result in a pale to dark blue coloration of the nuclei of the cells. Excessive or incomplete counterstaining may interfere with proper interpretation of results.

2.5 Mounting

To maintain optimal sensitivity and to prevent fading of Fast Red Chromogen a two-step mounting procedure is needed.

After counterstaining & bluing the slides will be mounted following a two-step protocol. Steps A and B have to be performed subsequently:

A) Liquid coverslipping

- Incubate slides in distilled or deionized water for at least 1 min;
- Remove slides from distilled or deionized water without further drying;
- Tap off excess liquid without allowing drying of specimens; carefully wipe the bottom side of the slides with a paper towel to remove water;
- Apply 4 drops of CINtec® PLUS Mount (1 drop corresponds to 35 – 40 µL of this aqueous mounting medium) per LBC slide and 8 drops per conventional smear, respectively. Avoid the generation of air bubbles. To prevent formation of tiny bubbles the first drop can be discarded onto a paper towel before applying CINtec® PLUS Mount on the specimen;

- Gently tilt and rotate the glass slide to generate a thin layer of mounting medium and to fully cover the specimen (Do NOT yet apply a glass or film coverslip!); check the distribution of the mounting medium on the slide by visual inspection;
- For drying place slides in a horizontal position:
 - Incubate ThinPrep® or BD SurePath™ samples at 37-60°C for 1 hour, or alternatively overnight at ambient temperature;
 - Incubate Conventional smears at 37°C for 4 hours, or at 60°C for 1 hour, or alternatively overnight at ambient temperature;

B) Glass or film coverslipping

- After complete drying of the aqueous mount, incubate slides in xylene for a minimum of 1 minute and up to a maximum of 20 minutes. Then, coverslip the slides using a xylene-based mounting medium.

NOTE: Sildes **must not** be dehydrated by ascending series of alcohol before being glass or film coverslipped.

- Let the xylene-based mounting medium dry at room-temperature;

NOTE: To minimize fading, protect slides from light and store at room temperature (20 – 25°C).

VII. Quality Control

Deviations from the recommended procedures for fixation and further processing of the cervical cytological specimens may produce substantial variability in results, necessitating regular performance of in-house controls.

Positive Control

Specimens processed in the same manner as the patient sample(s) should be used as positive controls. Positive controls are indicative of correctly prepared specimens and proper staining techniques. One positive control should be included in each staining run.

Known positive controls should only be utilized for monitoring the correct performance of processed specimens and test reagents, rather than as an aid in formulating a specific diagnosis of patient samples. If the positive controls fail to demonstrate appropriate positive staining, results with the test specimens should be considered invalid.

Negative Control

Use a negative control processed in a manner identical to the patient sample(s) with each staining run to verify the specificity of the staining procedure and to provide an indication of background staining. A variety of different cell types present in representative cervical cytology specimens and that are known to be negative for the expression of the p16^{INK4a} and Ki-67 antigens (such as e.g. superficial cells) may serve as an additional internal negative control to assess any background staining (this should be verified by the user).

Assay Verification

Prior to initial use of an antibody or staining system in a diagnostic procedure, the user should verify the antibody's specificity by testing it on a series of in-house specimens with known immunocytochemical performance characteristics represent-

ing known positive and negative specimens. Refer to the quality control procedures previously outlined in this section of the product insert and to the quality control requirements of the CAP Certification Program for Immunocytochemistry and/or CLSI (NCCLS) Quality Assurance for Immunocytochemistry, Approved Guideline.

VIII. Interpretation of Results

The CINtec® *PLUS* procedure causes two distinct coloured reaction products: a brown precipitate at the p16^{INK4a} antigen sites, and a red precipitate at the Ki-67 antigen sites. Brown staining of cells (cytoplasm and/or nuclei) indicates p16^{INK4a} over-expression. Red staining of cells (nuclei) indicates expression of Ki-67. Cells stained for both antigens exhibit brown cytoplasmic staining with typically pronounced red nuclei. A qualified pathologist/cytologist experienced in immunocytochemical procedures and trained on the interpretation of CINtec® *PLUS* stained slides must evaluate positive and negative controls before interpreting results. Interpretation of results must be made within the context of the patient's history and other diagnostic tests by a certified professional.

For the interpretation of cervical cytology slides stained with the CINtec® *PLUS* Kit the slides should be evaluated with regard to the presence of cervical epithelial cells showing both cytoplasmic brown and nuclear red staining indicative of simultaneous p16 and Ki-67 expression.

The presence of one or more cervical epithelial cells with co-localization of brown cytoplasmic immuno-staining AND red nuclear immuno-staining within the same cell is regarded as a positive CINtec® *PLUS* test result.

In case there is no cervical epithelial cell showing simultaneous brown cytoplasmic immuno-staining AND red nuclear immuno-staining the CINtec® *PLUS* test result is considered negative.

Please note that the presence of cervical epithelial cells that do show a single immuno-reactivity only for one of the two markers (such as e.g. brown staining for p16 only, or red staining for Ki-67 only) is not considered as a positive test result for the CINtec® *PLUS* Kit, even when both potential types of cervical cells showing single immuno-reactivity only are present in the same cytology slide.

In case of the presence of cells showing morphological signs of severe dyskaryosis, but no dual staining for p16 and Ki-67, morphological interpretation criteria should not be disregarded.

IX. Limitations

- For professional use only. Special training is required for the performance of immunocytochemical procedures.
- The clinical interpretation of any positive or negative staining should be evaluated within the context of clinical presentation and cytological criteria. The clinical interpretation of any positive or negative staining should be complemented by morphological studies using proper positive and negative controls as well as other diagnostic tests. It is the responsibility of a qualified pathologist/cytologist who is familiar with the proper use of antibodies, reagents, and methods to interpret all of the steps used to prepare and interpret the final immunocytochemistry preparation.
- The staining results in immunocytochemistry are strongly influenced by the quality of the cells stained. Accordingly, the appropriate handling of the protocol steps for fixation, washing, drying, or heating of the slides, as well as the avoidance of the use of reagents contaminated with bacteria

will significantly contribute to the overall result of the staining procedure. Deviations from the protocol may lead to artefacts, antibody trapping, or false-negative results. Inconsistent results may be due to variations in fixation or non-adequate sampling.

- Excessive or incomplete counterstaining may interfere with proper interpretation of results.
- The manufacturer provides these antibodies/reagents at optimal dilution for use according to the instructions provided herein, for immunocytochemistry testing on prepared liquid-based cytology (LBC) slides or conventional smears. Any deviation from the recommended test procedures may invalidate declared expected results; appropriate controls must be employed and documented. Users who deviate from the recommended test procedures must accept responsibility for interpretation of patient results under these circumstances.
- False-positive results may be seen due to non-immunological binding of proteins or substrate reaction products. They may also be caused by pseudoperoxidase activity (erythrocytes) and endogenous peroxidase activity (e.g. cytochrome C).
- Do not replace kit reagents with reagents carrying different lot numbers or with reagents from other manufacturers.

X. Troubleshooting

Refer to section XII. for contact details in case technical assistance is required.

Problem	Probable Cause	Suggested Action
1. No staining of slides	<u>1a.</u> Deviation from Instructions for Use;	<u>1a.</u> Carefully read Instructions for Use and adhere to the procedures outlined therein;
	<u>1b.</u> Reagent vials were not loaded in the correct locations in the reagent racks;	<u>1b.</u> Check the Reagent Map to verify the proper location of the reagent vials;
	<u>1c.</u> Insufficient reagent in vial;	<u>1c.</u> Ensure that enough reagent is loaded into the vials prior to commencing the run. Refer to Reagent Map for volumes required;
	<u>1d.</u> Missing Wash Buffer for Dako or LabVision Autostainers;	<u>1d.</u> Check that there is sufficient buffer. If not, then fill the container with buffer and prime the pump; Check if the clear plastic tubing around the syringe (placed on the left side of the moving arm) contains air bubbles - if so removal of the bubbles may solve the problem;
2. Weak staining of slides	<u>2a.</u> Inadequate epitope retrieval;	<u>2a.</u> Use freshly prepared Epitope Retrieval Solution and / or make sure that Epitope Retrieval Solution reaches 95 - 99°C for a full 10 minutes and is allowed to cool down for - additional 20 minutes; Ensure that no water slops over into the heat resistant staining jars during epitope retrieval procedure. Review 2.2 for recommendations;
	<u>2b.</u> Inadequate reagent incubation times;	<u>2b.</u> Review 2.3.1 / 2.3.2 Staining protocol recommendations;
	<u>2c.</u> Inappropriate fixation method;	<u>2c.</u> Ensure that cytology preparations are fixed as outlined in Section VII or that no alternative fixative was used;
	<u>2d.</u> Water that has been used to dilute the epitope-retrieval solution contains an ion concentration that is too high;	<u>2d.</u> Ensure that the ion exchange column for producing deionized water has been checked by routine maintenance;
	<u>2e.</u> Missing Wash Buffer for Dako or LabVision Autostainers;	<u>2e.</u> Refer to troubleshooting Section 1d;

3. No or weak Ki-67 staining (red) on positive control slides	3a. Use of alcohol containing Hematoxylin (e.g. Harris');	3a. Fast Red is soluble in EtOH therefore, the use of alcohol-free Hematoxylin is mandatory to prevent fading of Fast Red dye (Ki-67 signal);
	3b. Slides were dehydrated using ascending series of alcohol before being glass or film coverslipped;	3b. Follow the 2-step mounting protocol using CINtec <i>PLUS</i> Mount followed by permanent glass or film coverslipping; Review 2.5 for recommendations;
	3c. Fast Red working solution was not prepared according to the recommendations;	3c. Prepare Fast Red working solution immediately before use; apply the solution to the slides within 15 minutes of preparation, as otherwise a decrease in signal intensity and test sensitivity may occur;
4. Excessive background staining of slides	4a. Incomplete removal of spray fixation reagent (polyethylene glycol) when using conventional smears;	4a. Follow procedure as outlined in Section 2.1.;
	4b. Insufficient rinsing of slides;	4b. Use fresh solution in buffer baths and wash bottles; For Autostainer Instrument, ensure that before running the instrument is primed properly and that sufficient amounts of buffer are available for the whole run; use fresh solutions of buffers and reagents;
	4c. Drying of specimen during staining procedure;	4c. When using the Autostainer Instrument, ensure that the amount of reagent supplied is sufficient and that the hood of the instrument is closed during the run; Do not expose slides to elevated temperatures or to dry environment; Consider whether the programmed reagent volume and dispense location (drop zone) were adequate to cover all areas on a slide which contain cell materials; Check that there is sufficient buffer. If not, then fill the container with buffer and prime the pump; Check if the clear plastic tubing around the syringe (placed on the left side of the moving arm) contains air bubbles - if so removal of the bubbles may solve the problem; Avoid delay of the Fast Red working solution incubation step during autostainer run to minimize the risk of drying artefacts;

	<u>4d.</u> Inappropriate fixation method;	<u>4d.</u> Use only fixatives as recommended herein. Aberrantly fixed cells may exhibit excessive background staining;
	<u>4e.</u> Non-specific binding of reagents to cells;	<u>4e.</u> Check fixation method used for the specimen;
	<u>4f.</u> Drying of specimens during set up of Autostainer Instrument;	<u>4f.</u> Keep specimens covered with reagent (buffer) while initiating the run;
5. Excessively strong specific staining	<u>5a.</u> Inappropriate fixation method;	<u>5a.</u> Ensure proper fixative and fixation method;
	<u>5b.</u> Prolonged reagent incubation times;	<u>5b.</u> Review and adhere to staining protocol given in Sections 2.3.1 / 2.3.2 above;
	<u>5c.</u> Inappropriate wash solution;	<u>5c.</u> Use the Wash Buffer (catalog number 8557);
6. Evaluation of double-immuno-reactive cells not possible	<u>6a.</u> Hematoxylin counterstain too strong;	<u>6a.</u> Reduce incubation time for hematoxylin and/or bluing in tap water or ammonia water;
7. Glass or film coverslipping is not adequate	<u>7a.</u> Glass or film coverslip cannot be applied correctly;	<u>7a.</u> Prevent formation of air-bubbles when applying CIntec® PLUS Mount (aqueous coverslip); Adhere to recommended drying time for liquid coverslip;
8. Crack formation after mounting the slides with xylene-based mounting medium	<u>8a.</u> Insufficient drying of CIntec® PLUS Mount before mounting the slides with xylene-based mounting medium;	<u>8a.</u> Dry slides which have been covered with CIntec® PLUS Mount over night at room temperature or at least for 4h at 37°C. CIntec® PLUS Mount has to be dried completely.

XI. Symbols

Symbol:	Explanation:
	Catalog number
	Batch code
	Global Trade Item Number
	Unique Device Identifier
	In vitro diagnostic medical device
	Manufacturer
	Contains sufficient for <n> tests
	Consult instructions for use
	Use by
	Temperature limitation
	Date of manufacturing
	Do not re-use

XII. Manufacturer

Manufactured by: Roche mtm laboratories AG
Sandhofer Straße 116
68305 Mannheim
Germany
www.roche.com
<https://dialog.roche.com>

Contact for technical support (Telephone): +800 5505 6606

XIII. Revision Status

The current Instructions for use represent Version 1.7 released July 2022.

Changes to previous version (1.6, released May 2021):

- III., Principle of Procedure: antibody clone 274-11 AC3 corrected to 274-11 AC3V1
- IV., Storage, CINtec® PLUS Mount: explanation added regarding storage at ambient temperature
- IV., Equipment Required: Autostainer added (previously already listed in VI., 2. Staining Procedure)
- VI., 1.5: Storage CINtec *PLUS* Mount: harmonised to ambient temperature (previously room temperature)
- VI., 2. Staining Procedure: clarified that only LabVision Autostainer and Dako Autostainer PLUS have been validated by manufacturer, but similar instruments may be validated by user
- Editorial changes

XIV. Intellectual Property

CINtec and E6H4 are trademarks of Roche.

All other trademarks are the property of their respective owners.

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1. Kim WY, Sharpless NE. The regulation of INK4/ARF in cancer and aging. *Cell*. 2006 Oct 20;127(2):265-275.
2. Wentzensen N, von Knebel Doeberitz M. Biomarkers in cervical cancer screening. *Dis Markers*. 2007;23(4):315-330.
3. Cuschieri K, Wentzensen N. Human papillomavirus mRNA and p16 detection as biomarkers for the improved diagnosis of cervical neoplasia. *Cancer Epidemiol Biomarkers Prev*. 2008;17(10):2536-2545.
4. Roelens J, Reuschenbach M, von Knebel Doeberitz M, et al. p16INK4a immunocytochemistry versus human papillomavirus testing for triage of women with minor cytologic abnormalities: a systematic review and meta-analysis. *Cancer Cytopathol*. 2012;120(5):294–307.
5. Carozzi F, Grillo-Tos A, Confortini M, et al. Risk of high-grade cervical intraepithelial neoplasia during follow-up in HPV-positive women according to baseline p16-INK4A results: a prospective analysis of a nested substudy of the NTCC randomised controlled trial. *Lancet Oncol*. 2013;14(2):168–176.
6. Schmidt D, Bergeron C, Denton KJ, et al. p16/ki-67 dual-stain cytology in the triage of ASC-US and LSIL papanicolaou cytology: results from the European equivocal or mildly abnormal Papanicolaou cytology study. *Cancer Cytopathol*. 2011;119(3):158–166.
7. Petry KU, Schmidt D, Scherbring S, et al. Triage Pap cytology negative, HPV positive cervical cancer screening results with p16/Ki-67 Dual-stained cytology. *Gynecol Oncol*. 2011;121(3):505–509.
8. Scholzen T, Gerdes J. The Ki-67 protein: From the known and the unknown. *J. Cell. Physiol*. 2000;182:311–322.
9. Ravarino A, Nemolato S, Macciocu E, et al. CINtec PLUS immunocytochemistry as a tool for the cytologic diagnosis of glandular lesions of the cervix uteri. *Am J Clin Pathol*. 2012;138(5):652–656.
10. Singh M, Mockler D, Akalin A, et al. Immunocytochemical colocalization of p16(INK4a) and Ki-67 predicts CIN2/3 and AIS/adenocarcinoma. *Cancer Cytopathol*. 2012;120(1):26–34.
11. Ikenberg H, Bergeron C, Schmidt D, et al. Screening for cervical cancer precursors with p16/Ki-67 dual-stained cytology: results of the PALMS study. *J Natl Cancer Inst*. 2013;105(20):1550-1557.
12. Wentzensen N, Fetterman B, Castle PE, et al. p16/Ki-67 Dual Stain Cytology for Detection of Cervical Precancer in HPV-Positive Women. *J Natl Cancer Inst*. 2015 Sep 15;107(12):djv257.
13. Uijterwaal MH, Polman NJ, Witte BI, et al. Triage HPV-positive women with normal cytology by p16/Ki-67 dual-stained cytology testing: baseline and longitudinal data. *Int J Cancer*. 2015;136(10):2361-2368.
14. Wentzensen N, Schwartz L, Zuna RE, et al. Performance of p16/Ki-67 immunostaining to detect cervical cancer precursors in a colposcopy referral population. *Clin Cancer Res*. 2012;18(15):4154-4162.
15. Bergeron C, Ikenberg H, Sideri M, et al. Prospective evaluation of p16/Ki-67 dual-stained cytology for managing women with abnormal Papanicolaou cytology: PALMS study results. *Cancer Cytopathol*. 2015;123(6):373-381.
16. Dona MG, Vocaturo A, Giuliani M, et al. p16/Ki-67 dual staining in cervico-vaginal cytology: Correlation with histology, Human Papillomavirus detection and genotyping in women undergoing colposcopy. *Gynecol Oncol*. 2012;126(2):198-202.
17. Allia E, Ronco G, Coccia A, et al. Interpretation of p16(INK4a) /Ki-67 dual immunostaining for the triage of human papillomavirus-positive women by experts and nonexperts in cervical cytology. *Cancer Cytopathol*. 2015;123(4):212-218.

18. Gustinucci D, Giorgi Rossi P, Cesarini E, et al. Use of Cytology, E6/E7 mRNA, and p16INK4a-Ki-67 to Define the Management of Human Papillomavirus (HPV)-Positive Women in Cervical Cancer Screening. *Am J Clin Pathol*. 2016;145(1):35-45.
19. Rossi P, Borghi L, Ferro R, Mencarelli R. A population of 1136 HPV DNA-HR positive women: expression of p16(INK4a)/Ki67 Dual-Stain Cytology and cytological diagnosis. Histological correlations and cytological follow up. *Pathologica*. 2015;107(3-4):185-191.
20. Wright TC, Jr., Behrens CM, Ranger-Moore J, et al. Triaging HPV-Positive Women with p16/Ki-67 Dual-stained Cytology: Results from a Sub-study Nested into the ATHENA Trial. *Gynecol Oncol*. 2017;144:51-56.
21. Clarke MA, Cheung LC, Castle PE, et al. Five-Year Risk of Cervical Precancer Following p16/Ki-67 Dual-Stain Triage of HPV-Positive Women. *JAMA Oncol*. 2019;5(2):181-186.
22. Wentzensen N, Clarke MA, Bremer R, et al. Clinical Evaluation of Human Papillomavirus Screening With p16/Ki-67 Dual Stain Triage in a Large Organized Cervical Cancer Screening Program. *JAMA Intern Med*. 2019;179(7):881-888.
23. Wentzensen N, Fetterman B, Tokugawa D, et al. Interobserver reproducibility and accuracy of p16/Ki-67 dual-stain cytology in cervical cancer screening. *Cancer Cytopathol*. 2014;122(12):914-920.
24. Arean-Cuns C, Mercado-Gutierrez M, Paniello-Alastruey I, et al. Dual staining for p16/Ki67 is a more specific test than cytology for triage of HPV-positive women. *Virchows Arch*. 2018;473(5):599-606.
25. Ebisch RMF, van der Horst J, Hermsen M, et al. Evaluation of p16/Ki-67 dual-stained cytology as triage test for high-risk human papillomavirus-positive women. *Mod Pathol*. 2017;30(7):1021-1031.

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© 2022 Roche