

Elecsys BRAHMS PCT



Materials provided

| REF | | | SYSTEM |
|-------------|-------------|-----|----------------------------|
| 10456299190 | 10456299501 | 100 | cobas e 402 cobas e 801 |

- 2 × 8 bottle labels (calibrators)
- 2 × 14 bottle labels (controls)
- 6 empty labeled snap-cap bottles

For reagents, refer to the "Reagents" section.

Materials required (but not provided)

| REF | Description |
|-------------|-------------------------------------------------|
| 11776576322 | CalSet Vials, 2 × 56 empty snap-cap bottles |
| 03142949122 | ControlSet Vials, 2 × 56 empty snap-cap bottles |
| | General laboratory equipment |
| | Distilled or deionized water |
| | cobas e analyzer |

Additional materials for **cobas e 402** and **cobas e 801** analyzers:

| REF | Description |
|-------------|------------------------------------------------------------------------------------------------------------------------------------------|
| 06908799190 | ProCell II M, 2 × 2 L system solution |
| 04880293190 | CleanCell M, 2 × 2 L measuring cell cleaning solution |
| 07485409001 | Reservoir Cup, 8 cups to supply ProCell II M and CleanCell M |
| 05694302001 | AssayTip/AssayCup tray, 6 magazines × 6 magazine stacks × 105 assay tips and 105 assay cups, 3 wasteliners |
| 07485425001 | Liquid Flow Cleaning Cup, 2 adaptor cups to supply ISE Cleaning Solution / Elecsys SysClean for Liquid Flow Cleaning Detection Unit |
| 07485433001 | PreWash Liquid Flow Cleaning Cup, 1 adaptor cup to supply ISE Cleaning Solution / Elecsys SysClean for Liquid Flow Cleaning PreWash Unit |
| 11298500160 | ISE Cleaning Solution / Elecsys SysClean, 5 × 100 mL system cleaning solution |

For use in the USA only

System information

| | |
|------------|-------------------------------|
| Short name | ACN (application code number) |
| PCT | 10161 |

Indication for use

Immunoassay for the in vitro quantitative determination of PCT (procalcitonin) in human serum and plasma (K2-EDTA, K3-EDTA and Li-Heparin).

The electrochemiluminescence immunoassay "ECLIA" is intended for use on **cobas e** immunoassay analyzers.

Used in conjunction with other laboratory findings and clinical assessments, Elecsys BRAHMS PCT is intended for use as follows:

- to aid in the risk assessment of critically ill patients on their first day of ICU admission for progression to severe sepsis and septic shock,
- to determine the change in PCT level over time as an aid in assessing the cumulative 28-day risk of all-cause mortality for patients diagnosed with severe sepsis or septic shock in the ICU or when obtained in the emergency department or other medical wards prior to ICU admission,
- to aid in decision making on antibiotic therapy, for inpatients or patients in the emergency department with suspected or confirmed lower respiratory tract infections (LRTI) – defined as community-acquired pneumonia (CAP), acute bronchitis, and acute exacerbation of chronic obstructive pulmonary disease (AECOPD),
- to aid in decision making on antibiotic discontinuation for patients with suspected or confirmed sepsis.

Warnings and precautions

- The Elecsys BRAHMS PCT is not indicated to be used as a stand-alone diagnostic assay and should be used in conjunction with clinical signs and symptoms of infection and other diagnostic evidence. In cases where the laboratory results do not agree with the clinical picture or history, additional tests should be performed. Changes in PCT should always be interpreted in the context of the clinical status of the patient and other laboratory results. Decisions regarding antibiotic therapy should NOT be based solely on procalcitonin concentrations.
- There is no uniformly recognized interpretation of the change in PCT concentration levels for the prediction of mortality, and overall mortality is strongly dependent on many factors, including pre-existing patient risk factors and clinical course. The need to continue ICU care at Day 4 and other covariates (e.g., age, SOFA score) are also significant predictors of 28-day cumulative mortality risk. Validation of the Elecsys BRAHMS PCT test as an aid in predicting mortality was performed in a study population with an overall 28-day mortality of 22 %.
- Certain patient characteristics, such as severity of renal failure or insufficiency, may influence procalcitonin values and should be considered as potentially confounding clinical factors when interpreting PCT values. Increased PCT levels may be observed in severe illness such as polytrauma, burns, major surgery, prolonged or cardiogenic shock. PCT levels may not be elevated in patients infected by certain atypical pathogens, such as Chlamydia pneumoniae and Mycoplasma pneumoniae. The safety and performance of PCT-guided therapy for individuals younger than age 17 years, pregnant women, immunocompromised individuals or those on immunomodulatory agents, was not formally analyzed in the supportive clinical trials.

Summary

Sepsis is a daily challenge in the hospital setting. Today, various therapeutic strategies are known to improve survival in patients with sepsis. Early assessment is important for determination of the appropriate treatment.

PCT is the prohormone of the hormone calcitonin, but PCT and calcitonin are distinct proteins. Calcitonin is exclusively produced by C-cells of the thyroid gland in response to hormonal stimuli, whereas PCT can be produced by several cell types and many organs in response to pro-inflammatory stimuli, in particular by bacterial products.¹

In healthy people, plasma PCT concentrations are found to be below 0.1 ng/mL.² Depending on the clinical background, a PCT concentration above 0.1 ng/mL can indicate clinically relevant bacterial infection, requiring antibiotic treatment.³ PCT levels rise rapidly (within 6-12 hours) after a bacterial infectious insult with systemic consequences. The magnitude of the increase in PCT concentration correlates with the severity of the bacterial infection.⁴ At a PCT concentration > 0.5 ng/mL, a patient should be considered at risk of developing severe sepsis or septic shock.^{5,6} On the other hand, the relief of the septic infection is accompanied by a decrease in the PCT concentration which returns to normal with a half-life of 24 hours,^{7,8} i.e., the continuous decline of PCT is indicative of effective source control measures and has been implicated in the safe de-escalation of antibiotic therapy.^{9,10}

By evaluating PCT concentrations, the physician may use the findings to aid in the risk assessment of critically ill patients for progression to severe sepsis and septic shock. In addition, the change of PCT levels over time offers information about the risk of mortality after diagnosis of severe sepsis or septic shock.

Early after multiple traumas, major surgery, severe burns, or in neonates, PCT levels can be elevated independently of an infectious process, but the return to baseline is usually rapid. Viral infections, bacterial colonization, localized infections, allergic disorders, autoimmune diseases, and transplant rejection do not usually induce a significant PCT response (values < 0.5 ng/mL). Therefore, PCT is an important marker enabling specific differentiation between a bacterial infection and other causes of inflammatory reactions.³

The results of the Elecsys BRAHMS PCT assay on the **cobas** e 402 and **cobas** e 801 analyzers should be evaluated in the context of all laboratory findings and the total clinical status of the patient. In cases where laboratory results do not agree with the clinical picture or history, additional tests should be performed.

Test principle

Sandwich principle. Total duration of assay: 18 minutes.

- First incubation: Antigen in the sample (18 µL), a biotinylated monoclonal PCT-specific antibody, and a monoclonal PCT-specific antibody labeled with a ruthenium complex^{a)} react to form a sandwich complex.
- Second incubation: After streptavidin-coated microparticles have been added, the complex becomes bound to the solid phase via interaction of biotin and streptavidin.
- The reaction mixture is aspirated into the measuring cell, where the microparticles are magnetically captured onto the surface of the electrode. Unbound substances are then removed with ProCell II M. Application of a voltage to the electrode then induces chemiluminescent emission, which is measured by a photomultiplier.
- Results are determined via a calibration curve that is instrument-specifically generated by 2-point calibration and a master curve provided via **cobas** link.

a) $\text{Tris}(2,2\text{'-bipyridyl})\text{ruthenium(II)-complex } (\text{Ru}(\text{bpy})_3^{2+})$

Reagents

The **cobas** e pack (M, R1, R2) is labeled as PCT.

Elecsys BRAHMS PCT

| | |
|----------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| M | Streptavidin-coated microparticles, 1 bottle, 6.4 mL: Streptavidin-coated microparticles 0.72 mg/mL; preservative. |
| R1 | Anti-PCT-Ab-biotin, 1 bottle, 8.2 mL: Biotinylated monoclonal anti-PCT antibody (mouse) 2.0 µg/mL; phosphate buffer 95 mmol/L, pH 7.5; preservative. |
| R2 | Anti-PCT-Ab~Ru(bpy) ₃ ²⁺ , 1 bottle, 8.2 mL: Monoclonal anti-PCT antibody (mouse) labeled with ruthenium complex 5.6 µg/mL; Biotin scavenger antibody 1.2 ng/mL; phosphate buffer 95 mmol/L, pH 7.5; preservative. |
| PCT Cal1 | PCT calibrator 1 (lyophilized), 1 bottle for 4 mL: PCT (recombinant) approximately 0.10 ng/mL in a human serum matrix; preservative. |
| PCT Cal2 | PCT calibrator 2 (lyophilized), 1 bottle for 4 mL: PCT (recombinant) approximately 54 ng/mL in a human serum matrix; preservative. |
| PC PCT1 | PreciControl PCT 1 (lyophilized), 2 bottles each for 4 mL: PCT (recombinant) approximately 0.50 ng/mL in a human serum matrix; preservative. |
| PC PCT2 | PreciControl PCT 2 (lyophilized), 2 bottles each for 4 mL: PCT (recombinant) approximately 10 ng/mL in a human serum matrix; preservative. |

Calibrators: The exact lot-specific calibrator values are available via **cobas** link.

Controls: The exact lot-specific target values and ranges are available via **cobas** link.

Warnings and precautions

For in vitro diagnostic use for healthcare professionals. Exercise the normal precautions required for handling all laboratory reagents.

Infectious or microbial waste

Warning: Handle waste as potentially biohazardous material. Dispose of waste according to accepted laboratory instructions and procedures.

Environmental hazards

Apply all relevant local disposal regulations to determine safe disposal.

The Safety Data Sheet is available for professional users on request.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:



Warning

| | |
|------|----------------------------------------------------|
| H317 | May cause an allergic skin reaction. |
| H412 | Harmful to aquatic life with long lasting effects. |

Prevention:

| | |
|------|-----------------------------------|
| P261 | Avoid breathing dust. |
| P273 | Avoid release to the environment. |
| P280 | Wear protective gloves. |

Response:

| | |
|-------------|------------------------------------------------------------------|
| P333 + P313 | If skin irritation or rash occurs: Get medical advice/attention. |
| P362 + P364 | Take off contaminated clothing and wash it before reuse. |

Disposal:

| | |
|------|--------------------------------------------------------------------|
| P501 | Dispose of contents/container to an approved waste disposal plant. |
|------|--------------------------------------------------------------------|

Hazardous components:

- 2-methyl-2H-isothiazol-3-one hydrochloride

Product safety labeling follows EU GHS guidance.

Contact phone: +1-800-428-2336

All human material should be considered potentially infectious. All products derived from human blood are prepared exclusively from the blood of donors tested individually and shown to be free from HBsAg and antibodies to HCV and HIV. The testing methods use assays that have been approved or cleared by the FDA or that are in compliance with the legal rules of the European Union (IVDR 2017/746/EU, IVDD 98/79/EC, Annex II, List A).

However, as no testing method can rule out the potential risk of infection with absolute certainty, the material should be handled with the same level of care as a patient specimen. In the event of exposure, the directives of the responsible health authorities should be followed.^{11,12}

Avoid foam formation in all reagents and sample types (specimens, calibrators, and controls).

Storage and stability

Store at 2-8 °C.

Do not freeze.

Store the **cobas e pack upright** in order to ensure complete availability of the microparticles during automatic mixing prior to use.

| Stability of the cobas e pack | |
|-------------------------------|----------------------------------|
| unopened at 2-8 °C | up to the stated expiration date |
| on the analyzers | 16 weeks |

The lyophilized calibrators and controls are stable up to the stated expiration date.

| Stability of the reconstituted calibrators/controls | |
|-----------------------------------------------------------------|-----------------------------|
| at -20 °C (± 5 °C) | 3 months (freeze only once) |
| reconstituted calibrators/controls on the analyzers at 20-25 °C | 2 hours (use only once) |

Store **upright** in order to prevent the calibrator or control solution from adhering to the snap-cap.

Calibration

Traceability: This method has been standardized against the BRAHMS PCT LIA assay.

The predefined master curve is adapted to the analyzer using the relevant calibrators.

Calibration frequency: Calibration must be performed once per reagent lot using fresh reagent (i.e., not more than 24 hours after the **cobas e pack** was registered on the analyzer).

The calibration interval may be extended based on acceptable calibration verification values determined by the laboratory.

Renewed calibration is recommended as follows:

- every 12 weeks when using the same reagent lot
- every 28 days when using the same **cobas e pack** on the analyzer
- as required, such as when quality control findings are outside the defined limits

Quality control

For routine quality control procedures, use PreciControl PCT1 and PreciControl PCT2 or other suitable controls.

It is recommended to run the controls for the various concentration ranges individually at least once every 24 hours when the test is in use, once per **cobas e pack**, and following each calibration.

Adjust the limits and control intervals based on the laboratory's individual requirements. If values fall outside the limits, each laboratory is advised to establish corrective measures.

If necessary, repeat sample measurement.

Follow the applicable government regulations and local guidelines.

Note: When using two **cobas e packs** with different lots in the same run, the controls will be measured with both reagent lots. Use only control values measured with the corresponding lots.

Specimen collection and preparation

Only the specimens listed below were tested and found acceptable.

Serum collected using standard sampling tubes or tubes containing separator gel.

Li-heparin, K2 EDTA, and K3 EDTA plasma.

Li-heparin plasma tubes containing separator gel can be used.

Stable for 24 hours at 20-25 °C, 48 hours at 2-8 °C, 12 months at -20 °C (± 5 °C). Freeze only once.

The sample types listed were tested with a selection of sample collection tubes or systems that were commercially available at the time of testing. Not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials that could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube / collection system manufacturer.

Centrifuge samples containing precipitates before performing the assay.

Do not use samples and controls stabilized with azide.

Ensure the samples, calibrators, and controls are at 20-25 °C prior to measurement.

Due to possible evaporation effects, analyze and measure samples, calibrators, and controls on the analyzers within 2 hours.

Sample stability claims were established based on experimental data by the manufacturer and only for the temperatures / time frames as stated in the Method Sheet. It is the responsibility of the individual laboratory to use all available references and/or its own studies to determine specific stability criteria for its laboratory.

Test procedure

The reagents (M, R1, R2) in the kit are ready for use and are supplied in **cobas** e packs.

For optimum performance of the assay, follow the instructions given in this document for the corresponding analyzer. For analyzer-specific assay instructions, refer to the corresponding User Guide.

Resuspension of the microparticles takes place automatically prior to use.

Place the cooled (stored at 2-8 °C) **cobas** e pack on the reagent manager.

Avoid foam formation.

The system automatically regulates the temperature of the reagents and the opening/closing of the **cobas** e pack.

Calibrators and controls

Carefully dissolve the contents of 1 bottle by adding exactly 4.0 mL of distilled or deionized water, and allow it to stand closed for 15 minutes to reconstitute.

Mix carefully, avoiding foam formation.

Transfer the reconstituted calibrators / controls into the empty labeled snap-cap bottles supplied.

Unless the entire volume is necessary for calibration and quality control on the analyzers, transfer aliquots of the freshly reconstituted calibrators and controls into empty snap-cap bottles (CalSet Vials / ControlSet Vials). Attach the supplied labels to these additional bottles. Store the aliquots at -20 °C (± 5°C) for later use.

Perform **only 1** calibration or control procedure per aliquot.

Note: Do not combine bottles from different lots. Use only control bottles out of 1 lot with each other. All information required for correct operation is available via **cobas** link.

Calibrators

Place the reconstituted calibrators (in the system-compatible bottles with barcoded labels) in the sample zone.

All the information necessary for calibrating the assay is automatically read into the analyzer.

After calibration has been performed, discard the calibrators.

Controls

Place the reconstituted controls in the sample zone. After the control procedure has been performed, discard the controls.

Calculation

The analyzer automatically calculates the analyte concentration of each sample in ng/mL.

Interpretation of results

This assay is intended for use to determine the change of PCT over time as an aid in assessing the cumulative 28-day risk of all-cause mortality for patients diagnosed with severe sepsis or septic shock in the ICU, or when obtained in the emergency department or other medical wards prior to ICU admission.

SIRS (Systemic Inflammatory Response Syndrome), sepsis, severe sepsis, and septic shock were categorized according to the criteria of the consensus conference of the American College of Chest Physicians/Society of Critical Care Medicine.¹³

PCT should always be interpreted in the clinical context of the patient. Therefore, clinicians should use the PCT results in conjunction with other laboratory findings and clinical signs of the patient.^{5,6}

PCT > 2 ng/mL

A PCT level above 2.0 ng/mL on the first day of ICU admission is associated with a high risk for progression to severe sepsis and/or septic shock.


PCT < 0.5 ng/mL

A PCT level below 0.5 ng/mL on the first day of ICU admission is associated with a low risk for progression to severe sepsis and/or septic shock.

Note: Concentrations < 0.5 ng/mL do not exclude an infection, on account of localized infections (without systemic signs) which can be associated with such low concentrations, or a systemic infection in its initial stages (< 6 hours). Furthermore, increased procalcitonin can occur without infection. PCT concentrations between 0.5 and 2.0 ng/mL should be interpreted taking into account the patient's history. It is recommended to retest PCT within 6-24 hours if any concentrations < 2 ng/mL are obtained.

The change of PCT concentration over time provides prognostic information about the risk of mortality¹⁴ within 28 days for patients diagnosed with severe sepsis or septic shock coming from the emergency department, ICU, other medical wards, or directly from outside the hospital. Data support the use of PCT determinations from the day severe sepsis or septic shock is first diagnosed (Day 0) or the day thereafter (Day 1) and the fourth day after diagnosis (Day 4) for the classification of patients into higher and lower risk for mortality within 28 days according to the workflow below:

$$\Delta PCT = \frac{PCT_{Day0 \text{ (or Day1)}} - PCT_{Day4}}{PCT_{Day0 \text{ (or Day1)}}} \times 100\%$$



ΔPCT ≤ 80 %

A decrease of PCT levels below or equal to 80 % defines a positive ΔPCT test result representing a higher risk for 28-day all-cause mortality of patients diagnosed with severe sepsis or septic shock.

ΔPCT > 80 %

A decrease of PCT levels of more than 80 % defines a negative ΔPCT result representing a lower risk for 28-day all-cause mortality of patients diagnosed with severe sepsis or septic shock.

If Day 0 result is not available, Day 1 result may be used. If more than one PCT value is available on Day 0 (or Day 1), enter the highest value. If more than one PCT value is available on Day 4, enter the most recent value.

For convenience, one can use the Change in Procalcitonin Calculator to determine ΔPCT results from the absolute PCT concentrations of a patient obtained on the day severe sepsis or septic shock was first diagnosed (or 24 hours later) and on Day 4.

Go to www.BRAHMS-PCT-Calculator.com.

Note: By using this website tool, you are visiting a 3rd party website. Roche is not responsible for the content of this website. The upper end of the measuring range for the Elecsys BRAHMS PCT assay is 400 ng/mL when using a 1:4 dilution.

Decision making on antibiotic therapy for patients with suspected or confirmed LRTI

Initiation:

| PCT result | < 0.10 ng/mL | 0.10-0.25 ng/mL | 0.26-0.50 ng/mL | > 0.50 ng/mL |
|-----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------|
| Interpretation | Antibiotic therapy strongly discouraged. Indicates absence of bacterial infection. | Antibiotic therapy discouraged. Bacterial infection unlikely. | Antibiotic therapy encouraged. Bacterial infection possible. | Antibiotic therapy strongly encouraged. Suggestive of presence of bacterial infection. |
| Follow-up | Antibiotic therapy should be considered regardless of PCT result if the patient is clinically unstable, is at high risk for adverse outcome, has strong evidence of bacterial pathogen, or the clinical context indicates antibiotic therapy is warranted. If antibiotics are withheld, reassess if symptoms persist/worsen and/or repeat PCT measurement within 6-24 hours. | | In order to assess treatment success and to support a decision to discontinue antibiotic therapy, follow up samples should be tested once every 1-2 days, ¹⁵ based upon physician discretion taking into account patient's evolution and progress. Antibiotic therapy may be adjusted using the discontinuation table below: | |

Discontinuation:

Antibiotic therapy may be discontinued if the PCT_{Current} is ≤ 0.25 ng/mL or if the ΔPCT > 80 %.

- PCT_{Peak}: Highest observed PCT concentration.
- PCT_{Current}: Most recent PCT concentration.
- ΔPCT: Calculate by using the following equation:

$$\Delta PCT = \frac{PCT_{Peak} \text{ []} - PCT_{Current} \text{ []}}{PCT_{Peak} \text{ []}} \times 100\%$$

Antibiotic therapy may be continued based upon other clinical findings, such as apparent progression on chest x-ray or ongoing/increasing toxicity. If clinical picture has not improved, and PCT remains high, re-evaluate and consider treatment failure or other causes.

Decision making on antibiotic discontinuation for suspected or confirmed septic patients

In order to assess treatment success and to support a decision to discontinue antibiotic therapy, follow up samples should be tested once every 1-2 days,¹⁵ based upon physician discretion taking into account the patients' evolution and progress. Antibiotic therapy may be adjusted using the discontinuation table below:

Antibiotic therapy may be discontinued if the PCT_{Current} is ≤ 0.50 ng/mL or if the ΔPCT > 80 %

- PCT_{Peak}: Highest observed PCT concentration.
- PCT_{Current}: Most recent PCT concentration.
- ΔPCT: Calculate by using the following equation:

$$\Delta PCT = \frac{PCT_{Peak} \text{ []} - PCT_{Current} \text{ []}}{PCT_{Peak} \text{ []}} \times 100\%$$

Antibiotic therapy may be continued based upon other clinical findings, such as failure to control a local infection, or ongoing physiologic instability. **If clinical picture has not improved, and PCT remains high, re-evaluate and consider treatment failure or other causes.**

Recommendations for laboratory reports

It is suggested to report the numerical PCT values (individual or paired). For paired PCT values the report should also indicate if the ΔPCT(%) was ≤ 80 % or > 80 %. The laboratory report should include a reference or a link to the Elecsys BRAHMS PCT assay package insert for a guided interpretation of the test results:

www.BRAHMS-PCT-Calculator.com

Alternatively, the laboratory report may provide such interpretative criteria directly (as found on the website above) together with the absolute PCT concentrations and the "Change in Procalcitonin Result".

High-dose hook effect

There is no high-dose hook effect at PCT concentrations up to 1000 ng/mL.

Limitations and interferences

The effect of the following endogenous substances and pharmaceutical compounds on assay performance was tested. Interferences were tested up to the listed concentrations, and no impact on results was observed.

Endogenous substances

| Compound | Concentration tested |
|--------------------|--------------------------------|
| Bilirubin | ≤ 428 μmol/L or ≤ 25 mg/dL |
| Hemoglobin | ≤ 0.559 mmol/L or ≤ 0.900 g/dL |
| Intralipid | ≤ 1500 mg/dL |
| Rheumatoid factors | ≤ 1200 IU/mL |

Criterion: Recovery within ± 15 % of initial value for PCT concentrations > 0.1 ng/mL and within ± 0.015 ng/mL of initial values for PCT concentrations ≤ 0.1 ng/mL.

Biotin interference

Specimens with biotin concentrations up to 1200 ng/mL did not demonstrate bias in measured PCT values. Specimens with biotin concentrations > 1200 ng/mL and ≤ 2600 ng/mL demonstrated ≤ 10 % negative bias in measured PCT levels. Some studies have shown that serum concentrations of biotin can reach up to 355 ng/mL for subjects consuming supplements of 20 mg biotin per day¹⁶ and up to 1160 ng/mL for subjects after a single dose of 300 mg biotin.¹⁷

Other interferent testing:

- Human Anti-Mouse Antibody (HAMA) interference testing was completed with three PCT analyte concentrations using a high HAMA human serum pool. No interference was detected.
- Samples from patients routinely exposed to animals or animal serum products may contain heterophilic antibodies causing an atypical result. This assay has been formulated to mitigate the risk of this type of interference. However, potential interactions between rare sera and test components can occur.
- In rare cases, interference due to extremely high titers of antibodies to analyte-specific antibodies, streptavidin or ruthenium can occur. These effects are minimized by suitable test design.

In addition, the following 34 common and special drugs were tested. No interference with the assay was found.

Common and special drugs

| Drug | Concentration tested^{A)} |
|----------------------|------------------------------------------|
| Acetylcysteine | 150 mg/L |
| Ampicillin | 75 mg/L |
| Ascorbic acid | 52.5 mg/L |
| Cyclosporine | 1.8 mg/L |
| Cefoxitin | 750 mg/L |
| Heparin | 3300 IU/L |
| Levodopa | 7.5 mg/L |
| Methyl dopa | 22.5 mg/L |
| Metronidazole | 123 mg/L |
| Phenylbutazone | 321 mg/L |
| Doxycycline | 18 mg/L |
| Acetylsalicylic acid | 30 mg/L |
| Rifampicin | 48 mg/L |
| Acetaminophen | 156 mg/L |
| Ibuprofen | 219 mg/L |
| Theophylline | 60 mg/L |
| Imipenem | 1180 mg/L |
| Cefotaxime | 900 mg/L |
| Vancomycin | 3500 mg/L |
| Dopamine | 130 mg/L |
| Noradrenaline | 2 mg/L |
| Dobutamine | 11.2 mg/L |
| Furosemide | 20 mg/L |
| Cromolyn | 24 mg/L |
| Alcohol | 4000 mg/L |
| Azithromycin | 11.5 mg/L |
| Cetirizine HCl | 3.6 mg/L |
| Dextromethorphan | 1.4 mg/L |
| Levofloxacin | 17.5 mg/L |
| Loratadine | 0.3 mg/L |
| Nicotine | 1 mg/L |

| Drug | Concentration tested ^{A)} |
|-------------------|------------------------------------|
| Oxymetazoline HCl | 0.09 mg/L |
| Phenylephrine | 0.18 mg/L |
| Tiotropium | 0.0216 mg/L |

A) Tested on the cobas e 601 analyzer

Recovery was within ± 10 % of the reference value.

Limitations

For diagnostic purposes, always assess the results in conjunction with the patient's medical history, clinical examination, and other findings.

Increased PCT levels may not always be related to systemic infection.^{4, 18, 19, 20}

These include, but are not limited to:

- Patients experiencing major trauma and/or recent surgical procedure including extracorporeal circulation or burns.
- Patients undergoing treatment with OKT3 antibodies, OK-432, interleukins, TNF-alpha and other drugs that stimulate the release of pro-inflammatory cytokines or result in anaphylaxis.
- Patients diagnosed with active medullary C-cell carcinoma, small cell lung carcinoma, or bronchial carcinoid.
- Patients with acute or chronic viral hepatitis and/or decompensated severe liver cirrhosis (Child-Pugh Class C).
- Patients with prolonged or severe cardiogenic shock, prolonged severe organ perfusion anomalies, or after resuscitation from cardiac arrest.
- Patients receiving peritoneal dialysis or hemodialysis treatment.
- Patients with biliary pancreatitis, chemical pneumonitis or heat stroke.
- Patients with invasive fungal infections (e.g., candidiasis, aspergillosis) or acute attacks of plasmodium falciparum malaria.
- Neonates during the first 2 days of life.

The results of the Elecsys BRAHMS PCT assay should be evaluated in the context of all laboratory findings and the total clinical status of the patient. In cases where laboratory results do not agree with the clinical picture or history, additional tests should be performed.

Limits and ranges

Measuring range

0.02-100 ng/mL (defined by the Limit of Detection and the maximum of the master curve). Values below the Limit of Detection are reported as < 0.02 ng/mL. Values above the measuring range are reported as > 100 ng/mL.

Lower limits of measurement

Limit of Blank, Limit of Detection, and Limit of Quantitation

Limit of Blank = 0.015 ng/mL

Limit of Detection = 0.02 ng/mL

Limit of Quantitation = 0.06 ng/mL

The Limit of Blank, the Limit of Detection, and the Limit of Quantitation were determined in accordance with the CLSI (Clinical and Laboratory Standards Institute) EP17-A2 requirements.

The Limit of Blank is the 95th-percentile value from $n \geq 60$ measurements of analyte-free samples over several independent series. The Limit of Blank corresponds to the concentration below which analyte-free samples are found with a probability of 95 %.

The Limit of Detection is determined based on the Limit of Blank and the standard deviation of low-concentration samples. The Limit of Detection corresponds to the lowest analyte concentration that can be detected (value above the Limit of Blank with a probability of 95 %).

The Limit of Quantitation is defined as the lowest amount of analyte in a sample that can be accurately quantitated with a total allowable relative error of 20 %.

| Expected value ng/mL | Elecsys BRAHMS PCT (data from cobas e 601 analyzer) | | |
|-------------------------|-----------------------------------------------------|--------|-------|
| | % CV | % BIAS | % TE |
| 2.00 | 0.19 | 0.12 | 0.43 |
| 0.50 | 0.77 | 0.54 | 1.81 |
| 0.30 | 1.30 | 0.94 | 3.08 |
| 0.25 | 1.56 | 1.15 | 3.72 |
| 0.15 | 2.64 | 2.00 | 6.34 |
| 0.10 | 3.99 | 3.10 | 9.69 |
| 0.05 | 8.12 | 6.56 | 19.97 |

Dilution

Samples with PCT concentrations above the measuring range can be diluted manually with PCT-negative human serum or plasma. The recommended dilution is 1:4. The concentration of the diluted sample must be ≥ 20 ng/mL.

After manual dilution, multiply the result by the dilution factor.

Specific performance data

Representative performance data is given below. The precision data was generated on the **cobas e 402** and **cobas e 801** analyzers. However, since the **cobas e 402** and **cobas e 801** analyzers are members of the Elecsys instrument family of analyzers, some of the data below may have been generated on other members of the Elecsys instrument family. Results obtained in individual laboratories may differ.

Clinical performance

The Elecsys BRAHMS PCT assay was evaluated for the prediction of cumulative 28-day all-cause mortality in a prospective clinical trial²¹ study of 858 adult patients diagnosed with severe sepsis or septic shock admitted to ICU care in which PCT levels were measured on Days 0, 1, and 4 across 13 investigational sites in the US. The per protocol analysis population (598 subjects) was comprised of 44 % female and 56 % male patients with a mean age of 64 years. About half of the patients had severe sepsis (51 %) versus septic shock (49 %). Infections were mainly community acquired (91 %).


The binary test result (Δ PCT > 80 % or ≤ 80 %) was significantly associated with 28-day cumulative mortality (vital status on day 28) (Two-sided Fisher's Exact Test p-value = 0.002). Adjusted for ICU vs. non-ICU patient subgroups (based on hospital location at Day 4 after initial diagnosis), the association remained significant (Cochran-Mantel-Haenszel Test p-value = 0.020). In each binary Δ PCT subgroup, the 28-day cumulative mortality rate was stratified by need to continue ICU care on Day 4 and/or the selection of Day 0 vs. Day 1 as the baseline measurement day for the Δ PCT calculation. The data are as follows:

Prediction performances of binary Δ PCT stratified by ICU care on Day 4

| 28-Day mortality risk stratified by patient location on Day 4: Δ PCT > 80 % = Test Negative; Δ PCT ≤ 80 % = Test Positive per-protocol population | | | | | |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------|---------------------------|--------------------------|-----------------------------------|------------------|
| Day 4 patient location | Measurement interval | 28-Day mortality risk (%) | | Prognostic accuracy ^{A)} | |
| | | Δ PCT > 80 % | Δ PCT ≤ 80 % | Sensitivity | Specificity |
| ICU | Δ PCT Day 0-4 | 22.1 (13.3-31.0) | 29.6 (22.9-36.4) | 73.4 (62.9-83.8) | 35.0 (28.2-41.8) |
| | Δ PCT Day 1-4 | 21.5 (13.0-29.9) | 30.4 (23.4-37.3) | 71.6 (60.8-82.3) | 38.7 (31.7-45.7) |
| Non-ICU | Δ PCT Day 0-4 | 5.6 (1.8-9.4) | 11.0 (6.6-15.5) | 72.3 (55.9-88.6) | 44.4 (38.4-50.3) |
| | Δ PCT Day 1-4 | 7.1 (2.8-11.3) | 9.9 (5.7-14.2) | 65.4 (48.0-82.7) | 43.3 (37.3-49.2) |

A) Prognostic accuracy refers to how accurate the Δ PCT (> 80 % vs. ≤ 80 %) can predict mortality risk using 28-day mortality as the clinical reference.

Kaplan-Meier survival curves show that patients with a positive Δ PCT result (≤ 80 %) had a clearly lower survival probability from study Day 4 till the end of follow-up time compared to Δ PCT-negative (> 80 %) patients.



Additional stratification of patients based on absolute initial PCT levels ($>$ or ≤ 2.0 ng/mL) at Day 0 (or Day 1) revealed subgroups with particularly reduced or elevated mortality risk considering their hospital disposition on Day 4. Mortality rates and prognostic performance are given for the following subgroups in the table below:

1. Patients with PCT ≤ 2.0 ng/mL at Day 0 (or Day 1) receiving ICU care on Day 4.
2. Patients with PCT > 2.0 ng/mL at Day 0 (or Day 1) receiving ICU care on Day 4.
3. Patients with PCT ≤ 2.0 ng/mL at Day 0 (or Day 1) without ICU care on Day 4.
4. Patients with PCT > 2.0 ng/mL at Day 0 (or Day 1) without ICU care on Day 4.

| 28-Day mortality risk stratified by patient location on Day 4, absolute PCT value on Day 0: ΔPCT > 80 % = Test Negative; ΔPCT ≤ 80 % = Test Positive | | | | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------|--------------|---------------------------|---------------------|-----------------------------------|------------------|
| Per-Protocol Population: ΔPCT 0-4 stratified PCT at Day 0 | | | | | |
| Day 4 patient location | PCT at Day 0 | 28-Day mortality risk (%) | | Prognostic accuracy ^{A)} | |
| | | ΔPCT Day 0-4 > 80 % | ΔPCT Day 0-4 ≤ 80 % | Sensitivity (%) | Specificity (%) |
| ICU | ≤ 2.0 ng/mL | 10.4 (0.0-29.7) | 24.9 (15.2-34.6) | 94.9 (85.2-100) | 13.3 (4.9-21.8) |
| | > 2.0 ng/mL | 23.6 (14.0-33.2) | 33.2 (24.0-42.4) | 65.1 (51.8-78.3) | 46.3 (37.5-55.0) |
| Non-ICU | ≤ 2.0 ng/mL | 5.6 (0.0-16.3) | 8.3 (3.6-12.9) | 91.7 (76.0-100) | 12.3 (6.2-18.4) |
| | > 2.0 ng/mL | 5.6 (1.6-9.7) | 17.5 (7.5-27.5) | 58.6 (35.1-82.1) | 71.4 (64.1-78.8) |


A) Prognostic accuracy refers to how accurate the ΔPCT (> 80 % vs. ≤ 80 %) can predict mortality risk.


| 28-Day mortality risk stratified by patient location on Day 4, absolute PCT value on Day 1: ΔPCT > 80 % = Test Negative; ΔPCT ≤ 80 % = Test Positive | | | | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------|--------------|------------------------------|---------------------|-----------------------------------|------------------|
| Per-Protocol Population: ΔPCT 1-4 stratified PCT at Day 1 | | | | | |
| Day 4 patient location | PCT at Day 1 | 28-Day mortality risk (%) | | Prognostic accuracy ^{A)} | |
| | | ΔPCT Day 1-4 > 80 % | ΔPCT Day 1-4 ≤ 80 % | Sensitivity (%) | Specificity (%) |
| ICU | ≤ 2.0 ng/mL | 11.8 (0.0-33.6) | 25.5 (15.4-35.6) | 94.7 (84.6-100.0) | 12.7 (3.8-21.7) |
| | > 2.0 ng/mL | 22.5 (13.5-31.5) | 34.0 (24.4-43.5) | 62.9 (49.4-76.5) | 51.3 (42.4-59.7) |
| Non-ICU | ≤ 2.0 ng/mL | 0.0 (0.0-17.6) ^{B)} | 7.5 (3.0-12.0) | 100.0 (69.2-100.0) ^{B)} | 13.4 (7.3-19.5) |
| | > 2.0 ng/mL | 8.2 (3.3-13.1) | 15.3 (6.0-24.7) | 47.2 (24.6-69.7) | 69.4 (61.6-77.2) |

A) Prognostic accuracy refers to how accurate the ΔPCT (> 80 % vs. ≤ 80 %) can predict mortality risk.

B) Normality approximation of within-imputation variance not valid, therefore the estimate corresponds to within-imputation variation based on exact confidence intervals [Clopper & Pearson, 1934]

Kaplan-Meier Plots are depicted below to illustrate the time-to-event structure in the extended patient subgroups according to hospital location on Day 4 and initial PCT level.





Time-to-event analysis reveals that patients in the ICU or with an initial PCT value > 2.0 ng/mL had a lower survival probability from study Day 4 until the end of follow-up time (28 days) when the ΔPCT test result was positive compared to ΔPCT-negative patients as illustrated by the Kaplan-Meier curves above (patient subgroups according to hospital location on Day 4 and initial PCT level).

A generally lower mortality rate was observed in the non-ICU subgroup. The mortality rates for ΔPCT > 80 % vs. ΔPCT ≤ 80 % patient subgroups were:

| | |
|--------------------------------------------------------------------|-------------------|
| Patients with PCT ≤ 2.0 ng/mL at Day 0 receiving ICU care on Day 4 | 10.4 % vs. 24.9 % |
| Patients with PCT > 2.0 ng/mL at Day 0 receiving ICU care on Day 4 | 23.6 % vs. 33.2 % |
| Patients with PCT ≤ 2.0 ng/mL at Day 0 without ICU care on Day 4 | 5.6 % vs. 8.3 % |
| Patients with PCT > 2.0 ng/mL at Day 0 without ICU care on Day 4 | 5.6 % vs. 17.5 % |

Based on relative mortality ratios a decrease by more than 80 % from Day 0 (or Day 1) to Day 4 constitutes a lower risk for mortality within 28 days compared to smaller declines in each subgroup. For the prediction of absolute mortality risks ICU disposition at Day 4 and initial PCT concentrations should be considered:

1. An initial PCT level ≤ 2.0 ng/mL on Day 0 followed by a PCT decline of more than 80 % until Day 4 indicates an almost 2-fold lower mortality (10.4 %) for patients with severe sepsis or septic shock who are still in the ICU by Day 4 compared to those patients with an initial PCT level > 2.0 ng/mL (23.6 %). Regardless of the initial PCT level, patients in the ICU on Day 4 that do not decline by more than 80 % in PCT plasma concentration from Day 0 to Day 4 have an even higher mortality risk of 24.9-33.2 %.
2. An initial PCT level > 2.0 ng/mL that does not decline by more than 80 % until Day 4 signals that such patients remain at high mortality risk (17.5 %) even when they are no longer receiving ICU care on Day 4. Mortality was otherwise observed between 5.6 to 8.3 % for patients discharged from the ICU by Day 4.

In conclusion, hazard ratios for binary ΔPCT for both subgroups enrolled show an increase in mortality with ΔPCT. The mortality risk stratification using ΔPCT is valid when calculated from the day severe sepsis or septic shock is first diagnosed (Day 0) or the day thereafter (Day 1) and compared to the fourth day (Day 4) after diagnosis. The ΔPCT at Day 4 combined with the patient’s clinical course provides important information for the 28-day all-cause mortality risk prediction after a diagnosis of severe sepsis or septic shock.

The prognostic value of ΔPCT was quantified by pooled p-values of Wald statistics (Rubin’s rule). The hazard ratio of binary ΔPCT at Day 4 in the univariate model is 1.80 (p = 0.011) for patients of the Per-Protocol population. That means **the risk of death is increased 1.8-fold if an individual has a positive test result for ΔPCT.**

As a comparison, the table below lists the univariate hazard ratios for other clinical factors evaluated as separate predictors of mortality in the per protocol study population.

| Univariate hazard ratios for 28-day all-cause mortality of ΔPCT and clinical covariates | | | |
|-----------------------------------------------------------------------------------------|--------------------------------|------------------|---------|
| Factor | Comparison | Hazard ratio | p-Value |
| ΔPCT Day 0-4 | ≤ 80 % vs. > 80 % | 1.80 (1.15-2.82) | 0.011 |
| ΔPCT Day 1-4 | ≤ 80 % vs. > 80 % | 1.61 (1.04-2.49) | 0.034 |
| APACHE | Difference of 5 units | 1.36 (1.22-1.53) | < 0.001 |
| Max SOFA | Difference of 3 units | 1.73 (1.50-2.00) | < 0.001 |
| Antibiotic adequacy | No vs. yes | 1.59 (1.00-2.53) | 0.051 |
| Sepsis severity | Septic shock vs. severe sepsis | 1.19 (0.80-1.76) | 0.386 |

| Univariate hazard ratios for 28-day all-cause mortality of Δ PCT and clinical covariates | | | |
|-------------------------------------------------------------------------------------------------|------------------------------|------------------|---------|
| Factor | Comparison | Hazard ratio | p-Value |
| Biologic infection type | Gram pos vs. gram neg | 0.83 (0.48-1.45) | 0.522 |
| Biologic infection type | Other vs. gram neg | 0.99 (0.63-1.54) | 0.960 |
| Biologic infection type | Fungal vs. gram neg | 2.44 (0.87-6.84) | 0.090 |
| Clinical infection type | Nosocomial vs. community | 0.76 (0.35-1.64) | 0.481 |
| Positive blood culture | Yes vs. no | 1.05 (0.69-1.58) | 0.834 |
| Baseline PCT | > 2 ng/mL vs. \leq 2 ng/mL | 1.43 (0.94-2.17) | 0.095 |
| Age | Difference of 5 years | 1.16 (1.08-1.24) | < 0.001 |
| Gender | Male vs. female | 0.95 (0.64-1.40) | 0.782 |
| ICU care on Day 4 | Yes vs. no | 3.45 (2.24-5.31) | < 0.001 |

Δ PCT from Day 0 (or Day 1) to Day 4 remains a prognostic parameter for the risk of cumulative 28-day mortality in patients diagnosed with severe sepsis or septic shock even when hazard ratios are adjusted for other mortality predictors in multivariate models. The relative mortality risk estimates for Δ PCT and selected predictors are given below with 95 % confidence intervals. For continuous predictors, the hazard ratio (HR) was calculated for one standard deviation (SD) change in the predictor. For binary predictors, the risk estimate compares the hazards for the two binary results.

Hazard ratios for Δ PCT were calculated for each group as univariate and when applying certain covariates and are tabulated below.

| Hazard ratios for Δ PCT per group | |
|---------------------------------------------------------------------|-----------------------------------|
| Hazard ratio for Δ PCT \leq 80 % | Per-protocol |
| Δ PCT Day 0-4 (univariate) | 1.80 [1.15-2.82]; p=0.0106 |
| Δ PCT Day 0-4 with binary APACHE + covariates ^{A)} | 1.94 [1.14-3.31]; p=0.0140 |
| Δ PCT Day 0-4 with numeric APACHE + covariates ^{A)} | 1.72 [1.00-2.95]; p=0.0487 |
| Δ PCT Day 0-4 with binary SOFA + covariates ^{A)} | 1.76 [1.05-2.96]; p=0.0320 |
| Δ PCT Day 0-4 with numeric SOFA + covariates ^{A)} | 1.46 [0.86-2.48]; p=0.1595 |
| Δ PCT Day 1-4 (univariate) | 1.61 [1.04-2.49]; p=0.0345 |
| Δ PCT Day 1-4 with binary APACHE + covariates ^{A)} | 1.68 [1.03-2.76]; p=0.0392 |
| Δ PCT Day 1-4 with numeric APACHE + covariates ^{A)} | 1.61 [0.98-2.65]; p=0.0625 |
| Δ PCT Day 1-4 with binary SOFA + covariates ^{A)} | 1.64 [1.00-2.69]; p=0.0483 |
| Δ PCT Day 1-4 with numeric SOFA + covariates ^{A)} | 1.47 [0.89-2.42]; p=0.1300 |

A) Antibiotic adequacy, Sepsis severity, ICU Care on Day 4, Biological infection type, Clinical infection type, Positive blood culture, PCT on Day 0, Age, Gender

| Hazard ratios for Δ PCT and selected predictors from multivariate Cox Regression Models | | | |
|------------------------------------------------------------------------------------------------|----------------------------------|-----------------------------------------|-----------------------------------------|
| Model | | Hazard ratio (95 % confidence interval) | |
| | | Binary predictors ^{A)} | |
| Δ PCT interval | Score + Covariates ^{B)} | Δ PCT (\leq 80 % vs > 80 %) | Day 4 patient location (ICU vs. no ICU) |
| Day 0 until Day 4 | APACHE | 1.72 (1.00-2.95) | 2.61 (1.63-4.19) |
| | max SOFA | 1.46 (0.86-2.48) | 1.71 (1.04-2.81) |
| Day 1 until Day 4 | APACHE | 1.61 (0.98-2.65) | 2.63 (1.64-4.21) |
| | max SOFA | 1.47 (0.89-2.42) | 1.73 (1.06-2.84) |

A) In the analysis, missing values for predictors were multiple imputed assuming they were Missing at Random (MAR), with the multiple imputations combined according to Rubin's rules (Rubin D.B., Wiley New York 1987; Multiple Imputation for Nonresponse in Surveys).

B) The models also included the following predictors (HR results not shown): Antibiotic adequacy, Sepsis severity, Biological infection type, Clinical infection type, Positive blood culture, PCT on Day 0, Gender.

| Hazard ratios for ΔPCT and selected predictors from multivariate Cox Regression Models | | | | |
|----------------------------------------------------------------------------------------|----------------------------------|-------------------------------------------------------------|------------------------|--------------------|
| Model | | Hazard ratio (95 % confidence interval) | | |
| | | Continuous predictors ^{A)} (Hazard ratio per 1 SD) | | |
| ΔPCT interval | Score + Covariates ^{B)} | APACHE (1 SD = 8.13) | max SOFA (1 SD = 3.98) | Age (1 SD = 16.18) |
| Day 0 until Day 4 | APACHE | 1.25 (0.99-1.57) | ---- | 1.59 (1.27-1.99) |
| | max SOFA | ---- | 1.97 (1.53-2.53) | 1.69 (1.35-2.10) |
| Day 1 until Day 4 | APACHE | 1.29 (1.04-1.62) | ---- | 1.57 (1.25-1.96) |
| | max SOFA | ---- | 2.00 (1.56-2.56) | 1.67 (1.33-2.08) |

A) In the analysis, missing values for predictors were multiple imputed assuming they were Missing at Random (MAR), with the multiple imputations combined according to Rubin's rules (Rubin D.B., Wiley New York 1987; Multiple Imputation for Nonresponse in Surveys).

B) The models also included the following predictors (HR results not shown): Antibiotic adequacy, Sepsis severity, Biological infection type, Clinical infection type, Positive blood culture, PCT on Day 0, Gender.

The change of PCT over time can also be described by the ratio of PCT values from Day 4 and Day 0 (or Day 1):

$$PCT_{ratio} = \frac{PCT_{Day\ 4}}{PCT_{Day\ 0\ (or\ Day\ 1)}}$$

A decline of ΔPCT = 80 % translates into a PCT ratio of 0.2. The PCT ratio has values larger than 0.2 when the ΔPCT decline is below 80 % which is associated with a higher risk for cumulative 28-day all-cause mortality in patients diagnosed with severe sepsis or septic shock. Likewise, a PCT ratio below 0.2 indicates a lower risk for mortality within 28 days. On a continuous scale, the relative mortality risk for patients diagnosed with severe sepsis or septic shock is higher the larger the PCT ratio. The following tables list the hazard ratios for an increase by the factor 2 in PCT ratio, i.e. the relative increase in mortality risk for a patient with any given PCT ratio compared to a patient with a 2-fold lower PCT ratio. For comparison, selected predictors are indicated with corresponding equivalents in standard deviation. For the patient location at Day 4, the risk estimate compares the hazards for patients with vs. without ICU care on Day 4.

| Hazard ratios for ΔPCT and selected predictors from multivariate Cox Regression Models | | |
|----------------------------------------------------------------------------------------|----------------------------------|-----------------------------------------|
| Model | | Hazard ratio (95 % confidence interval) |
| | | Binary predictors ^{A)} |
| ΔPCT interval | Score + Covariates ^{B)} | Day 4 patient location (ICU vs. no ICU) |
| Day 0 until Day 4 | APACHE | 2.57 (1.59-4.13) |
| | max SOFA | 1.70 (1.03-2.80) |
| Day 1 until Day 4 | APACHE | 2.57 (1.60-4.11) |
| | max SOFA | 1.74 (1.06-2.86) |

A) In the analysis, missing values for predictors were multiple imputed assuming they were Missing at Random (MAR), with the multiple imputations combined according to Rubin's rules (Rubin D.B., Wiley New York 1987; Multiple Imputation for Nonresponse in Surveys).

B) The models also included the following predictors (HR results not shown): Antibiotic adequacy, Sepsis severity, Biological infection type, Clinical infection type, Positive blood culture, PCT on Day 0, Gender.

| Hazard ratios for ΔPCT and selected predictors from multivariate Cox Regression Models | | | | | |
|----------------------------------------------------------------------------------------|----------------------------------|-------------------------------------------------------------------------------------------------------------|--------------------------------------|--------------------------|---------------------|
| Model | | Hazard ratio (95 % confidence interval) | | | |
| | | Continuous predictors ^{A)} (Hazard ratio per 2-fold increase in PCT ratio or per equivalent in SD) | | | |
| ΔPCT interval | Score + Covariates ^{B)} | PCT ratio (2-fold increase) | APACHE (SD equivalent) ^{C)} | Max SOFA (SD equivalent) | Age (SD equivalent) |
| Day 0 until Day 4 | APACHE | 1.26 (1.12-1.42) | 1.08 (0.96-1.22) | ---- | 1.29 (1.15-1.45) |
| | max SOFA | 1.20 (1.07-1.35) | ---- | 1.37 (1.20-1.57) | 1.32 (1.18-1.49) |
| Day 1 until Day 4 | APACHE | 1.29 (1.11-1.49) | 1.19 (1.02-1.39) | ---- | 1.37 (1.18-1.60) |
| | max SOFA | 1.23 (1.06-1.44) | ---- | 1.58 (1.33-1.87) | 1.43 (1.23-1.67) |

A) In the analysis, missing values for predictors were multiple imputed assuming they were Missing at Random (MAR), with the multiple imputations combined according to Rubin's rules (Rubin D.B., Wiley New York 1987; Multiple Imputation for Nonresponse in Surveys).

B) The models also included the following predictors (HR results not shown): Antibiotic adequacy, Sepsis severity, Biological infection type, Clinical infection type, Positive blood culture, PCT on Day 0, Gender.

C) A unit change of ΔPCT on log-2-scale corresponded to 0.52 SD of ΔPCT from Day 0 until Day 4 (0.69 SD for ΔPCT from Day 1 until Day 4). Accordingly, the reported ΔPCT hazard ratios refer to an increase of ΔPCT by a factor of 2. For comparability, hazard ratios of the other continuous predictors were estimated for the same fractional SDs, i.e. 0.52 or 0.69, respectively.

Cumulative 28-day all-cause mortality did not differ significantly for male vs. female patients (χ^2 p-value = 0.84).

Demographics with outcome information are shown below:

| Variable | Class | Per-protocol population (N = 598) | | | |
|--------------------------|------------------|-----------------------------------|----------|-----------|---------------|
| | | All (N) | Dead (N) | Alive (N) | Mortality (%) |
| Gender | Female | 264 | 46 | 218 | 17.4 |
| | Male | 334 | 55 | 279 | 16.5 |
| Age, years (categorized) | ≤ 30 | 39 | 1 | 38 | 2.6 |
| | > 30, ≤ 45 | 45 | 4 | 41 | 8.9 |
| | > 45, ≤ 55 | 74 | 8 | 66 | 10.8 |
| | > 55, ≤ 65 | 149 | 26 | 123 | 17.4 |
| | > 65, ≤ 75 | 125 | 21 | 104 | 16.8 |
| | > 75 | 166 | 41 | 125 | 24.7 |
| Ethnicity | African-American | 202 | 32 | 170 | 15.8 |
| | Asian | 7 | 0 | 7 | 0.0 |
| | Caucasian | 362 | 64 | 298 | 17.7 |
| | Hispanic | 23 | 5 | 18 | 21.7 |
| | Other | 4 | 0 | 4 | 0.0 |
| Baseline PCT, ng/mL | < 0.5 | 125 | 19 | 106 | 15.2 |
| | ≥ 0.5, ≤ 2.0 | 104 | 13 | 91 | 12.5 |
| | > 2.0 | 353 | 69 | 284 | 19.5 |
| | Missing | 16 | 0 | 16 | 0.0 |

Reference range

In a population of 282 self-reported healthy individuals, the 95th percentile, upper reference range limit was calculated at 0.08 ng/mL.

Analytical performance data

Representative performance data on the analyzers are given below. Results obtained in individual laboratories may differ.

Precision

Precision was determined using Elecsys reagents, samples, and controls based on a protocol (EP05-A3) of the CLSI (Clinical and Laboratory Standards Institute): 2 runs per day, in duplicate each, for 21 days (n = 84). The following results were obtained:

| cobas e 801 analyzer ^{A)} | | | | |
|------------------------------------|--------------|---------------|------------------------|---------------|
| Sample | Mean (ng/mL) | Repeatability | Intermediate precision | % Total error |
| | | CV % | CV % | |
| Human serum 1 | 0.036 | 6.9 | 10.1 | 24.6 |
| Human serum 2 | 0.043 | 7.4 | 8.1 | 19.7 |
| Human serum 3 | 1.26 | 1.5 | 1.9 | 4.78 |
| Human serum 4 | 22.1 | 2.1 | 2.7 | 6.37 |
| Human serum 5 | 39.3 | 1.5 | 1.8 | 4.43 |
| Human serum 6 | 75.0 | 1.6 | 2.3 | 5.45 |
| Human serum 7 | 0.082 | 3.9 | 4.2 | 10.3 |
| Human serum 8 | 99.7 | 1.4 | 1.9 | 4.51 |
| PreciControl PCT1 | 0.488 | 1.5 | 1.8 | 4.24 |
| PreciControl PCT2 | 9.25 | 1.3 | 1.7 | 4.28 |

A) The precision data generated on the cobas e 402 analyzer was equivalent to that of the cobas e 801 analyzer.

Method comparison

A comparison of the Elecsys BRAHMS PCT assay on the cobas e 801 analyzer (y) with the Elecsys BRAHMS PCT assay on the cobas e 411 analyzer (x) using clinical samples gave the following correlation (ng/mL):

Elecsys BRAHMS PCT



Number of samples measured: 146

Passing/Bablok²²

$y = 1.037x - 0.00133$

$r = 0.985$

The sample concentrations were between 0.025 and 92.7 ng/mL.

Analytical specificity

The Elecsys BRAHMS PCT assay on the **cobas** e 402 and **cobas** e 801 analyzers does not show any significant cross-reactions with the following substances, tested with PCT concentrations of approximately 0.05 ng/mL and 2.0 ng/mL (maximum tested concentration):

| Substances | Non-interfering concentrations |
|--------------------------------|--------------------------------|
| Human katacalcin | 30 ng/mL |
| Human calcitonin | 10 ng/mL |
| Human alpha-CGRP ^{A)} | 10000 ng/mL |
| Human beta-CGRP | 10000 ng/mL |
| Calcitonin eel | 30 mg/L |
| Calcitonin salmon | 30 mg/L |

A) Calcitonin Gene-Related Peptide

This product may not be used by purchaser to conduct any Point-of-Care testing including but not limited to any near patient testing utilizing analyzers with a throughput equal to or less than 50 tests per hour. No general patent or other license of any kind other than this specific right of use from purchase is granted hereby.

Reagent developed in collaboration with B·R·A·H·M·S.

B·R·A·H·M·S and B·R·A·H·M·S PCT are registered trademarks of BRAHMS GmbH.

B · R · A · H · M · S



Additional information

Additions, deletions, or changes are indicated by a change bar in the margin.

For further information, refer to the User Guide for the corresponding analyzer, to the corresponding application sheets, and to the Method Sheets of all necessary components.

Report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the member state in which the user and/or patient is established.

FOR US CUSTOMERS ONLY: LIMITED WARRANTY

Roche Diagnostics warrants that this product meets the specifications stated in the labeling when used in accordance with the labeling and is free from defects in material and workmanship until the expiration date printed on the label. THIS LIMITED WARRANTY IS IN LIEU OF ANY OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. IN NO EVENT SHALL ROCHE DIAGNOSTICS BE LIABLE FOR INCIDENTAL, INDIRECT, SPECIAL, OR CONSEQUENTIAL DAMAGES.

Symbols

For definition of symbols used, refer to navifyportal.roche.com.

In addition to the ISO 15223-1 standard, Roche Diagnostics uses the following symbols and signs:

| | |
|-------------------|------------------------------------------------------------------------------------------------|
| CONTENT | Contents of kit |
| SYSTEM | Analyzers/Instruments on which reagents can be used |
| REAGENT | Reagent |
| CALIBRATOR | Calibrator |
| | Volume for reconstitution |
| GTIN | Global Trade Item Number |
| Rx only | For USA: Caution: Federal law restricts this device to sale by or on the order of a physician. |

References

- 1 Christ-Crain M, Müller B. Procalcitonin in bacterial infections-hype, hope, more or less? *Swiss Med Wkly* 2005 Aug 6;135(31-32):451-460.
- 2 Morgenthaler NG, Struck J, Fischer-Schulz C, et al. Detection of procalcitonin (PCT) in healthy controls and patients with local infection by a sensitive ILMA. *Clin Lab* 2002;48(5-6):263-270.
- 3 Christ-Crain M, Jaccard-Stolz D, Bingisser R, et al. Effect of procalcitonin-guided treatment on antibiotic use and outcome in lower respiratory tract infections: cluster-randomised, single-blinded intervention trial. *Lancet* 2004 Feb 21;363(9409):600-607.
- 4 Meisner M. ISBN:978-3-8374-1241-3 UNI-MED Science, 2010. Procalcitonin-Biochemistry and Clinical Diagnosis.
- 5 Müller B, Becker KL, Schächinger H, et al. Calcitonin precursors are reliable markers of sepsis in medical intensive care unit. *Crit Care Med* 2000 Apr;28(4):977-983.
- 6 Harbarth S, Holeckova K, Froidevaux C, et al. Diagnostic value of procalcitonin, interleukin-6, and interleukin-8 in critically ill patients admitted with suspected sepsis. *Am J Resp Crit Care Med* 2001 Aug 1;164(3):396-402.
- 7 Luyt CE, Guérin V, Combes A, et al. Procalcitonin kinetics as a prognostic marker of ventilator-associated pneumonia. *Am J Respir Crit Care Med*. 2005 Jan1;171(1):48-53.
- 8 Brunkhorst FM, Heinz U, Forycki ZF. Kinetics of procalcitonin in iatrogenic sepsis. *Intensive Care Med*. 1998 Aug;24(8):888-889.
- 9 Schuetz P, Christ-Crain M, Thomann R, et al. Effect of procalcitonin-based guidelines vs standard guidelines on antibiotic use in lower respiratory tract infections: the ProHOSP randomized controlled trial. *JAMA* 2009 Sep 9;302(10):1059-1066.
- 10 Bouadma L, Luyt CE, Tubach F, et al. Use of procalcitonin to reduce patients' exposure to antibiotics in intensive care units (PRORATA trial): a multicentre randomised controlled trial. *Lancet* 2010 Feb 6;375(9713):463-474.
- 11 Occupational Safety and Health Standards: Bloodborne pathogens. (29 CFR Part 1910.1030). Fed. Register.
- 12 Directive 2000/54/EC of the European Parliament and Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work.
- 13 American College of Chest Physicians/Society of Critical Care Medicine Definitions for sepsis and organ failure and guidelines for the use of innovative therapies in sepsis. *Crit Care Med* 1992;20:864-874.
- 14 Schuetz P, Maurer P, Punjabi V, et al. Procalcitonin decrease over 72 hours in US critical care units predicts fatal outcome in sepsis patients. *Crit Care* 2013 Jun 20;17(3):R115.
- 15 Schuetz P, Raad I, Devendra NA. Using procalcitonin-guided algorithms to improve antimicrobial therapy in ICU patients with respiratory infections and sepsis. *Curr Opin Crit Care* 2013 Oct;19(5):453-60.
- 16 Grimsey P, Frey N, Bendig G, et al. Population pharmacokinetics of exogenous biotin and the relationship between biotin serum levels and in vitro immunoassay interference. *Int J Pharmacokinet* 2017;2(4):247-256.
- 17 Piketty ML, Prie D, Sedel F, et al. High-dose biotin therapy leading to false biochemical endocrine profiles: validation of a simple method to overcome biotin interference. *Clin Chem Lab Med* 2017;55(6):817-825.
- 18 Meisner M, Tschaikowsky K, Hutzler A, et al. Postoperative plasma concentrations of procalcitonin after different types of surgery. *Int Care med* 1998 Jul;24(7):680-684.
- 19 Chiesa C, Panero A, Rossi N, et al. Reliability of Procalcitonin Concentrations for the Diagnosis of Sepsis in Critically ill neonates. *Clin Infect Dis* 1998;26:664-672.
- 20 Reith HB, Mittelkötter U, Debus ES, et al. Procalcitonin in early detection of postoperative complications. *Dig Surg* 1998;15(3):260-265.
- 21 Schuetz P, Birkhahn R, Sherwin R, et al. Serial procalcitonin predicts mortality in severe sepsis patients: Results from the multicenter procalcitonin MONitoring SEpsis (MOSES) Study, *Critical Care Medicine* 2017 May;45(5):781-789.
- 22 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. *J Clin Chem Clin Biochem*. 1988 Nov;26(11):783-790.

COBAS, ELECSYS, NAVIFY, and PRECICONTROL are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

© 2026 Roche

Rx only



Roche Diagnostics GmbH
 Sandhofer Strasse 116
 68305 Mannheim, Germany
www.roche.com
 ☎ +800 5505 6606



Elecsys BRAHMS PCT



Roche Diagnostics
9115 Hague Road
Indianapolis, IN 46256, USA
☎ +1 800 428 2336