



CoaguChek[®] XS Plus System

Supplement to the Operator's Manual v9.0

Document version	Date	Content
1.0	2019-11	New document, created as an amendment to Operator's Manual Version 9.0

Purpose of document

This document provides a description of changes to the Operator's Manual version 9.0 for SW 03.05.xx. This addendum is valid for and intended to be used in conjunction with the CoaguChek XS Plus Operator's Manual version 9.0. It is not a replacement for the complete CoaguChek XS Plus Operator's Manual version 9.0. **Changes or additions to the information** in the CoaguChek XS Plus Operator's Manual version 9.0 are shown in **color**.



Be sure to read the complete CoaguChek XS Plus Operator's Manual version 9.0.

The contents of this document, including all graphics, are the property of Roche Diagnostics. No part of this document may be reproduced or transmitted in any form or by any means, electronic or mechanical, for any purpose, without the express written permission of Roche Diagnostics.

Roche Diagnostics has made every reasonable effort to ensure that all the information contained in this manual is correct at the time of printing. However, Roche Diagnostics reserves the right to make any changes necessary without notice as part of ongoing product development.

Caution: Federal law restricts this device to sale by or on the order of a physician.

This manual was created by the Roche Diagnostics Engineering Operations department. Direct questions or concerns regarding the contents of this document to:

Roche Diagnostics
Engineering Operations Department
9115 Hague Road
PO. Box 50457
Indianapolis, IN 46250-0457

ACCU-CHEK, COAGUCHEK and SAFE-T-PRO are trademarks of Roche.

All other trademarks are the property of their respective owners.

© 2019, Roche Diagnostics.

US Order Number 09099093001 (01) 2019-11 USA

This document is available electronically at www.coaguchek-usa.com

1	What is new in SW 03.05?	5
	Options setup	5
	Database	5
	Result Retention	6
	Delete Database	7
	Diagnostics	8
	Error History	8
2	Amendments to Operator's Manual version 9.0	9
	Revised sections	9
	Important safety instructions and additional information	9
	Local Area Network: protection from unauthorized access	9
	Wired network connection	10
	Checking the software version	10
	Settings summary	11


This page intentionally left blank.

1 What is new in SW 03.05?

Options setup

Database

In the *Database* menu¹ you can set up your CoaguChek XS Plus meter to automatically delete test results after a certain time period, and/or delete stored information (test results and code chip data) manually.




- 1 From the *Main Menu*, touch **Setup** to open the meter settings.
- 2 From the *Setup Menu*, touch **Options**.
- 3 From the *Options* menu, touch **Database**.

1. If an *Admin.* password has been set on the meter, the *Connection* and *Database* setup areas are only available to the system administrator.

Result Retention

In the *Result Retention* menu you can enable or disable automatic deletion of test results after a selectable time period.




- 1 From the *Database* menu, touch **Result Retention**.

You may select from the following deletion options:

- **Off** (test results will never be deleted automatically)
 - **On** (test results will be deleted automatically after the time period set below)
- 2 Touch or to select the time period in days, after which data will be deleted automatically.
 - 3 Touch to save this setting, or touch to exit this menu without saving any changes. The display automatically returns to the previous screen.

Delete Database

In the Delete Database menu you can manually delete all test results, code chip data, or the entire database.



1 From the *Database* menu, select the type of data you want to delete:


- **Results** (all stored test results)
- **Code Chip Data** (all stored code chip data)
- **Entire Database** (all test results, comments, code chip data, operator and patient IDs)

2 Touch ✓ to confirm the deletion, or touch ✕ to exit this function without deleting data.

Diagnostics

Error History

Use the *Diagnostics* menu¹ to access the *Error History* of the meter. All errors occurring during internal self tests or regular operation are logged in this list.



- 1 From the *Main Menu*, touch **Setup** to open the meter settings.
- 2 From the *Setup Menu*, touch **Diagnostics**.
- 3 From the *Diagnostics* menu, touch **Error History**.
- 4 Touch and to scroll to the entry of choice on the screen.

1. If an *Admin.* password has been set on the meter, the *Diagnostics* setup area is only available to the system administrator.

2 Amendments to Operator's Manual version 9.0

Revised sections

Important safety instructions and additional information

The following information has been **added** on page 16:



Avoidance of electrical shock, fire, and explosions

- Only use Roche original accessories (cables, power supply units, battery packs, and spare parts). Third-party cables, power supply units, and battery packs can cause the battery pack to explode or the meter to become damaged.
 - Do not use loose power sockets or damaged power supply units, cables, plugs, or battery packs **and batteries**.
 - Do not short circuit the power supply unit, the handheld base unit contacts, or the battery pack.
 - Do not drop the CoaguChek XS Plus meter, the power supply unit, or the battery pack and protect these against shaking and vibrations.
-

The following information has been **added** on page 20:

Local Area Network: protection from unauthorized access

- If this meter is connected to a local area network, this network must be protected against unauthorized access. In particular, it must not be linked directly to any other network or the Internet. Customers are responsible for the security of their local area network, especially in protecting it against malicious software and attacks. This protection might include measures, such as a firewall, to separate the device from uncontrolled networks as well as measures that ensure that the connected network is free of malicious code. If user data and/or passwords are altered by means of a compromised DMS the meter might become unusable.
- If you use a customized data management system solution, ensure that sensitive data transmitted via the POCT1-A interface is protected by appropriate security measures.
- Ensure that the instrument is protected against unauthorized physical access and theft.
- Do not use shared user or operator accounts on meter, DMS, and network.
- Use a strong password for user or operator accounts on the meter, DMS, and network. Observe your own facility guidelines on password management where available, or apply the rules for strong passwords, see "Characteristics of strong passwords" below.

Wired network connection

If the Handheld Base Unit from Roche is used to connect this meter to a local area network, the Handheld Base Unit must be protected against unauthorized access by means of a **strong password management**. Observe your own facility guidelines on password management where available, or apply the following rules:

Characteristics of strong passwords


- Passwords must not contain the user's account name or parts of the user's full name that exceed two consecutive characters.
- Passwords must be at least eight characters in length.
- Passwords must contain characters from at least three of the following four categories:
 - English **uppercase alphabetic** characters (A through Z)
 - English **lowercase alphabetic** characters (a through z)
 - **Numeric** characters (0 through 9)
 - **Non-alphabetic** characters (for example, !, \$, #, %)

Examples of weak passwords

- **uhxwze11** contains no upper case letter.
- **UHXW13SF** contains no lower case letter.
- **uxxxxx7F** contains the same character more than four times.
- **x12useridF** contains a substring of the user ID longer than four characters.

Checking the software version


The following information has been **revised** on page 30:



After displaying the Roche logo, the meter briefly displays the *Init* (for "initialization") screen. Here you can check which software version is currently running on your meter. (The Init screen shown here is for illustration purposes only. Version numbers on your meter may differ.)

Settings summary

The following illustration has been revised on page 32:



The following information has been added on page 34:

Group	Subgroup	Setting	Values *
Options	Sort		Date/Time *
			Patient ID
			Patient Name
	Beeper	Beeper	Off
			Low
			Medium *
			High
		Key Click	Off *
			On
	Auto Off	[minutes]	Off
			1 ... (5*) ... 10
			15
			20
			25
			30
			40
			50
			60
	Connection		Off *
			Computer
	Database	Result Retention	Automatic result deletion On/Off*
			Retention time in days: 1 ... 30* ... 1000
		Delete Database	Results
			Code Chip Data
			Entire Database

* Default settings are labelled with an asterisk (*).

The following information has been added on page 35:

Group	Subgroup	Setting	Values *
ID Setup	Admin. (Administrator)		Blank (Off) *
	Operator (Operator List is optional)		Active
			Inactive *
	Patient		No *
			Optional
			Required
QC Settings	QC Range	Default Range*	
		Custom Range	Display Target Value On/Off *
			Deviation from Target Value Percentage (%)
	Operator Lockout (only if the <i>Operator</i> option is set to <i>Active</i> and an operator list is available)		Off *
			Weekly
			Monthly
			Every 3 months
			Every 6 months
			Yearly
	QC Lockout	New code	Yes
			No *
		General	No *
			Daily
			Weekly
			Monthly
Diagnostics	Error History		

* Default settings are labelled with an asterisk (*).

This page intentionally left blank.

ACCU-CHEK, COAGUCHEK and
SAFE-T-PRO are trademarks of Roche.



Roche Diagnostics GmbH
Sandhofer Strasse 116
68305 Mannheim, Germany

Made in Germany

Manufactured for and distributed in the U.S.A. by:
Roche Diagnostics
9115 Hague Road
Indianapolis, IN 46250-0457

Rx only

www.roche.com
www.coaguchek-usa.com
diagnostics.roche.com