

## Instructions for Use

# CINtec<sup>®</sup> Histology Kit

The CINtec<sup>®</sup> Histology Kit is an immunohistochemistry assay for the qualitative detection of the p16<sup>INK4a</sup> antigen on tissue sections prepared from formalin-fixed, paraffin-embedded cervical biopsies. It is indicated to be used in conjunction with H&E-stained slides prepared from the same cervical tissue specimen as an aid to increase diagnostic accuracy and inter-observer agreement in the diagnosis of high-grade cervical intraepithelial neoplasia.



Roche Diagnostics GmbH  
Sandhofer Strasse 116  
68305 Mannheim  
Germany

<https://navifyportal.roche.com>

REF

10213370001

GTIN

07613336230350



50



2 – 8 °C



0123

IVD



## Table of Contents / Inhalt / Table des matières / Indice / Índice

|   |    |
|---|----|
| ENGLISH .....   | 2  |
| I. Product Name.....                                    | 2  |
| II. Intended Use .....                                  | 2  |
| III. Summary and Explanation of the Device .....        | 2  |
| Summary and Explanation .....                           | 2  |
| Clinical Significance .....                             | 3  |
| Principle of Procedure.....                             | 4  |
| IV. Reagents.....                                       | 4  |
| Materials Provided .....                                | 4  |
| Storage .....   | 6  |
| Materials and Reagents Required but not Provided .....  | 6  |
| Equipment Required .....                                | 7  |
| V. Warnings and Precautions.....                        | 7  |
| Warning .....   | 7  |
| Caution .....   | 8  |
| VI. Procedure.....                                      | 8  |
| Specimen Preparation.....                               | 8  |
| Paraffin-embedded Tissue Specimens.....                 | 9  |
| Heat-Induced Epitope Retrieval .....                    | 9  |
| Staining Procedure .....                                | 9  |
| 1. Reagent Preparation.....                             | 9  |
| 1.1 Epitope Retrieval Solution.....                     | 9  |
| 1.2 Wash Buffer.....                                    | 10 |
| 1.3 Substrate-Chromogen Solution (DAB) .....            | 10 |
| 1.4 Counterstain .....                                  | 10 |
| 1.5 Mounting Medium .....                               | 10 |
| 2. Staining Procedure for Autostainer Instruments ..... | 10 |
| 2.1 Deparaffinization and Rehydration.....              | 10 |
| 2.2 Staining Protocol for Autostainer Instruments ..... | 11 |
| 3. Staining Procedure for Manual Use .....              | 13 |
| 3.1 Deparaffinization and Rehydration.....              | 13 |
| 3.2 Staining Protocol for Manual Use.....               | 14 |
| VII. Quality Control.....                               | 16 |
| VIII. Interpretation of Results.....                    | 16 |
| IX. Limitations .....                                   | 17 |
| X. Performance characteristics.....                     | 18 |
| XI. Troubleshooting.....                                | 23 |
| XII. Symbols .....                                      | 24 |
| XIII. Manufacturer .....                                | 25 |
| XIV. Revision Status.....                               | 25 |
| XV. Intellectual Property.....                          | 25 |
| Annex 1 References .....                                | 26 |
| Annex 2 .....   | 29 |

## ENGLISH

### I. Product Name

CINtec® Histology Kit

### II. Intended Use

For in-vitro diagnostic use.

The CINtec® Histology Kit is an immunohistochemistry assay for the qualitative detection of the p16<sup>INK4a</sup> antigen on tissue sections prepared from formalin-fixed, paraffin-embedded cervical biopsies.

It is indicated to be used in conjunction with H&E-stained slides prepared from the same cervical tissue specimen as an aid to increase diagnostic accuracy and inter-observer agreement in the diagnosis of high-grade cervical intraepithelial neoplasia.

The test is intended for manual use, or for use on Autostainer instruments.

This assay should be interpreted by a qualified pathologist in conjunction with histological examination, relevant clinical information, and proper controls.

### III. Summary and Explanation of the Device

#### Summary and Explanation

The CINtec® Histology Kit is based on a mouse monoclonal antibody (clone E6H4®) directed against human p16<sup>INK4a</sup> protein.

The p16<sup>INK4a</sup> protein is a cyclin-dependent kinase inhibitor which plays a major role in the regulation of the eukaryotic cell cycle. It is part of the retinoblastoma protein (pRB)-mediated control of the G1-S-phase transition, and it triggers cell cycle arrest in the course of cellular differentiation processes. In terminally differentiated epithelial cells, p16<sup>INK4a</sup> is expressed at levels typically not detectable by immunohistochemistry.

In various tumor entities, the p16<sup>INK4a</sup> gene was found to be functionally inactivated by gene mutation or promoter hypermethylation. This inactivation of the p16<sup>INK4a</sup> tumor suppressor gene has been shown to contribute to cell-cycle dysregulation and to loss of control of cell proliferation.

In replication competent cervical epithelial cells however, where high-risk type Human Papillomavirus (HR-HPV) oncoproteins have initiated the cellular transformation process, p16<sup>INK4a</sup> expression has been shown to be highly up-regulated [1; 2]. This strong over-expression of p16<sup>INK4a</sup> has been closely linked at the molecular level to the activity of the E7 oncoproteins from HR-HPV. The p16<sup>INK4a</sup> over-expression was shown to reflect the E7 oncoprotein-mediated inactivation of the functional complex between pRB and transcription factor E2F, which is one of the key events during the HR-HPV induced cellular transformation [3].

There are numerous studies in the published literature reporting that over-expression of the p16<sup>INK4a</sup> protein was immunohistochemically observed in a very high proportion of cases of high-grade pre-cancerous cervical dysplasia (i.e. 80 – 100% of CIN2 lesions, and virtually all CIN3 lesions) and invasive cancers. Low-grade cervical intra-epithelial lesions (CIN1) have been demonstrated to show p16<sup>INK4a</sup> over-expression at variable rates, typically within a range of 30 – 60% [1; 2; 4-16].

Interpretation of the results must take into consideration the fact that p16<sup>INK4a</sup> is a

cellular protein that may be expressed at detectable levels in cervical high-grade dysplastic lesions and cervical cancers as well as in some conditions not associated with cervical dysplasia, albeit at differing levels and with different patterns of expression. Histological tissue preparations have intact tissue morphology to aid in the interpretation of p16<sup>INK4a</sup> positivity of the cervical lesion. A diffuse staining pattern, i.e. a continuous staining of cells of the basal and parabasal cell layers, with or without staining of cells of superficial cell layers, has been suggested to be rated as a positive test result for p16<sup>INK4a</sup> over-expression. This staining pattern was shown to provide the highest level of both sensitivity and specificity for high-grade CIN [1; 4; 5]. In contrast, a focal staining pattern (staining of isolated cells or small cell clusters, i.e. non-continuous staining, particularly not of the basal and parabasal cells) as well as lack of any immunoreactivity is regarded as a negative test result for p16<sup>INK4a</sup> over-expression [1; 4; 5].

The interpretation of slides stained for p16<sup>INK4a</sup> using the CINtec® Histology Kit should be performed in conjunction with H&E-stained slides prepared from the same cervical tissue specimen. The additional information provided by the CINtec®-stained slides should be combined with the preliminary morphology-based diagnosis established on the H&E-stained slides in order to create a final diagnosis.

## **Clinical Significance**

The conjunctive reading of H&E-stained slides comprising cervical biopsy sections together with consecutive slides from the same tissue specimen and immuno-stained for p16<sup>INK4a</sup> has been shown to improve the diagnostic accuracy and the inter-observer agreement in diagnosing high-grade cervical intraepithelial neoplasia (CIN2+).

The pathologic diagnosis made on the H&E-stained cervical tissue sections establishes the basis for the decision to pursue further treatment. Consequently, the impact of an inaccurate diagnosis is significant. Inaccurate diagnosis may lead to inappropriate management of the patient, i.e. over-treatment of essentially healthy women, or under-treatment of women with established high-grade dysplastic lesions.

The diagnostic interpretation of histological H&E-stained tissue sections is subject to high rates of discordance among pathologists. Low rates of inter- and intra-observer agreement in histology of the cervix have been reported in various publications [17-20].

A large multicenter trial in the United States assessed the interpretation of histological cervical specimen surveys (2237 colposcopic biopsies and 535 LEEP conization biopsy specimens) by multiple well-trained pathologists; the reproducibility of the histopathological interpretations was only moderate (kappa=0.46 for punch biopsies and kappa=0.49 for LEEP biopsy specimens) [17]. Using the WHO and a modified Bethesda grading system, the interobserver agreement of six histopathologists assessing 125 colposcopic biopsy specimens has been shown to be poor for both grading systems [18]. Likewise, a study conducted in the UK showed poor interobserver agreement among eight expert histopathologists examining 100 colposcopic biopsy specimens (unweighted kappa value of 0.358) [19].

By adding slides stained with the CINtec® Histology Kit to the conventional H&E-stained slides used for establishing the diagnoses, the overall accuracy of this histomorphological diagnostic procedure is improved [21-39].

## **Principle of Procedure**

The CINtec® Histology Kit contains a set of reagents for the immunohistochemical detection of the p16<sup>INK4a</sup> antigen. The Kit is designed to perform a two-step immunohistochemical staining procedure for formalin-fixed, paraffin-embedded tissue specimens obtained from cervical biopsies. For the detection of the antigen, a primary monoclonal mouse antibody clone E6H4® directed to human p16<sup>INK4a</sup> protein is used.

A ready-to-use visualization reagent comprising a polymer reagent conjugated with horseradish peroxidase and goat anti-mouse Fab' antibody fragments is used. The Visualization Reagent has been subjected to solid-phase absorption to eliminate cross reactivity with human immunoglobulins. The chromogen reaction is based on horseradish peroxidase-mediated conversion of a DAB chromogen to a visible reaction product at the antigen site. After counterstaining, the specimen may be cover-slipped and the results may be evaluated by light microscopy inspection.

## **IV. Reagents**

### **Materials Provided**

The materials listed below are included in each kit and are sufficient to perform 50 tests and 50 negative control reactions. The number of tests is based on the use of 200 µL of the reagents per slide.

#### **1 Peroxidase Blocking Reagent**

##### **Peroxidase Blocking Reagent**

2 x 11.5 mL, ready-to-use

3% hydrogen peroxide, containing 15 mmol/L sodium azide (NaN<sub>3</sub>).

EUH210: Safety data sheet available on request.

#### **2 Mouse anti-Human p16<sup>INK4a</sup> Antibody**

##### **Mouse anti-Human p16<sup>INK4a</sup> Antibody**

11.5 mL, ready-to-use

Monoclonal mouse anti-Human p16<sup>INK4a</sup> antibody (~1 µg/mL), Clone E6H4™, supplied in 50 mmol/L Tris buffer pH 7.2, containing 15 mmol/L sodium azide (NaN<sub>3</sub>) and stabilizing protein.

#### **3 Visualization Reagent**

##### **Visualization Reagent**

2 x 11.5 mL, ready-to-use



Warning

H317 May cause an allergic skin reaction.

P261 Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.

P272 Contaminated work clothing should not be allowed out of the workplace.

P280 Wear protective gloves.

P333 + P313 If skin irritation or rash occurs: Get medical advice/ attention.

P362 + P364 Take off contaminated clothing and wash it before reuse.

P501 Dispose of contents/ container to an approved waste disposal plant.

Contains:

26172-54-3 2-methyl-2H-isothiazol-3-one hydrochloride

55965-84-9 reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H -isothiazol-3-one [EC no. 220-239-6] (3:1)

Polymer reagent conjugated with horseradish peroxidase and affinity purified goat anti-Mouse Fab' antibody fragments, supplied in stabilizing solution comprising preservatives and stabilizing protein.

#### 4 Negative Reagent Control

##### **Negative Reagent Control**

11.5 mL, ready-to-use

Monoclonal mouse anti-Rat oxytocin-related neurophysin antibody (~1 µg/mL), supplied in 50 mmol/L Tris buffer pH 7.2, containing 15 mmol/L sodium azide (NaN<sub>3</sub>) and stabilizing protein. For verification of the specificity of the staining. Rat oxytocin-related neurophysin is not present in human tissues.

#### 5 DAB Buffered Substrate

##### **DAB Buffered Substrate**

31 mL

Substrate buffer solution, pH 7.5, containing < 0.1 % hydrogen peroxide, stabilizers and enhancers.

#### 6 DAB Chromogen

##### **DAB Chromogen**

0.85 mL, 3,3'-diaminobenzidine chromogen solution.



Danger

H314 Causes severe skin burns and eye damage.

H341 Suspected of causing genetic defects.

H350 May cause cancer.

P201 Obtain special instructions before use.

P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.

P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.

P304 + P340 + P310 IF INHALED: Remove person to fresh air and keep comfortable for breathing. Immediately call a POISON CENTER/doctor.

P305 + P351 + P338 + P310 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/doctor.

P308 + P313 IF exposed or concerned: Get medical advice/ attention.

Contains 868272-85-9 3,3'-Diaminobenzidine tetrahydrochloride hydrate

**NOTE:** Consult Federal, State or local regulations for disposal.

## 7 Epitope Retrieval Solution 10X

### **Epitope Retrieval Solution 10X**

500 mL, 100 mmol/L Tris buffer pH 9 containing 10 mmol/L EDTA and 15 mmol/L sodium azide (NaN<sub>3</sub>).

### **Storage**

Store at 2 – 8 °C. Do not use after the expiration date. No data has been generated respective to the storage of the reagents under any conditions other than those stated above.

After opening, kit components are stable for 6 months if stored at 2 – 8 °C. Solutions must be discarded if cloudy in appearance.

Diluted Wash Buffer and diluted Epitope Retrieval Solution are stable for up to one month if stored at 2 – 8 °C. Solutions must not be used if cloudy in appearance.

### **Materials and Reagents Required but not Provided**

**CINtec® Wash Buffer 10X** to be used with the CINtec® Histology Kit is available under catalog number 10215364001 from Roche but is not included in the kit. For order details please refer to the website [www.roche.com](http://www.roche.com).

500 mmol/L Tris buffer solution with 1.5 mol/L NaCl, pH 7.6, containing detergent and an antimicrobial agent.



Warning

H317 May cause an allergic skin reaction.

P261 Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.

P272 Contaminated work clothing should not be allowed out of the workplace.

P280 Wear protective gloves.

P333 + P313 If skin irritation or rash occurs: Get medical advice/ attention.

P362 + P364 Take off contaminated clothing and wash it before reuse.

P501 Dispose of contents/ container to an approved waste disposal plant.

Contains 55965-84-9 reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H -isothiazol-3-one [EC no. 220-239-6] (3:1)

Absorbent wipes;  
Hematoxylin counterstain;  
Distilled or deionized water (Washing Water);  
Ethanol, 95% and 70%;  
Mounting medium;  
Positive and Negative Tissues to use as process controls;  
Slides (SuperFrost® Plus or equivalent);  
Xylene;  
Coverslips.

### **Equipment Required**

Optional: Drying oven, capable of maintaining 60 °C or less;  
Optional: Dako or LabVision Autostainer instrument;  
Humid chamber (optional);  
Light microscope (4 – 40x objective magnification);  
Staining jars or baths;  
Wash bottles;  
Timer (capable of 2 – 60 minute intervals);  
Water bath with lid (capable of maintaining Epitope Retrieval Solution at 95 – 99 °C).

## **V. Warnings and Precautions**

### **Warning**

1. Caution! Some of the reagents comprised in this kit contain hazardous chemicals. When handling the components of this kit, adhere to safety precautions for handling hazardous laboratory reagents.
2. Components 1, 2, 4, and 7 of this product contain sodium azide ( $\text{NaN}_3$ ), which is highly toxic in its pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. Components 2, 3 and 4 contain material of animal origin. Adhere to proper handling procedures as is applicable to any product derived from biological sources.
4. Safety Data Sheet for the kit is available upon request.
5. When handling and disposing of histology specimens, including all specimens before and after fixation, as well as all materials exposed to them, adhere to the safety precautions for handling potentially infectious material as well as applicable waste disposal requirements.
6. Never pipette reagents by mouth. Avoid contacting the skin and mucous



membranes with reagents and specimens. In the event that reagents or specimens come in contact with skin or mucous membranes, wash with copious amounts of water.

7. The Visualization Reagent and DAB Chromogen may be affected adversely if exposed to excessive light levels. Do not store kit components or perform staining in strong light, such as direct sunlight.
8. Wear appropriate personal protective equipment to avoid contact with eyes and skin when handling any of the components included or to be used in conjunction with the CINtec® Histology Kit. Refer to the Safety Data Sheet (SDS) for additional information.
9. Product safety labeling primarily follows EU GHS guidance.

## **Caution**

1. For in-vitro diagnostic use.
2. For professional use only.
3. Minimize microbial contamination of reagents to avoid non-specific staining.
4. Incubation times, temperatures, or methods other than those specified may give erroneous results.
5. Do not use the kit if the packaging of any of its components is damaged. Should packaging be compromised or components damaged, please notify the manufacturer without delay.
6. Disposal of all waste materials must be in accordance with local guidelines and regulations.
7. All reagents are formulated specifically for use with this test. In order for the test to perform as specified, no substitutions should be made.
8. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the CINtec® Histology Kit is suspected, immediately refer to the contact information provided in section XIII. for further information on technical support.
9. Malfunction of the product due to handling problems or to instability does not result in obvious signs. Therefore, as a quality control measure, positive and negative controls should be run simultaneously with patient specimens.
10. To report suspected serious incidents related to this device, contact the local Roche representative and the competent authority of the member State or Country in which the user is established.

## **VI. Procedure**

### **Specimen Preparation**

The CINtec® Histology Kit is designed to be used with tissue specimens preserved for immunohistochemistry procedures. Specimens should be prepared in accordance to standard methods of tissue processing.

Positively charged slides such as SuperFrost® Plus slides are recommended for optimal performance.

## **Paraffin-embedded Tissue Specimens**

Neutral buffered formalin-fixed, paraffin-embedded tissue specimens processed using routine standard methods are appropriate for use with this kit. If specimens are prepared using a different preservation method, the user must verify the appropriateness of the method.

Specimens from the biopsy should be fixed for 18 – 24 hours in neutral-buffered formalin (10% recommended) and blocked into a thickness of 3 or 4 mm. The tissue blocks are then dehydrated in a series of alcohols and xylene, followed by infiltration by melted paraffin held at no more than 60 °C. Each tissue block will be sectioned at 4 – 5 µm and mounted on SuperFrost® Plus microscopy slides by a histopathology laboratory. Slides should be stained immediately, as antigenicity of cut tissue sections may diminish over time.

## **Heat-Induced Epitope Retrieval**

For heat-induced epitope retrieval, the tissue sections mounted on slides must be heated by immersing into the Epitope Retrieval Solution in a calibrated water bath capable of maintaining the Epitope Retrieval Solution at a temperature of 95 – 99 °C. Laboratories located at higher elevations should determine the best method of maintaining the required water bath temperature. Manufacturer does not recommend any deviation from the procedure described herein.

After heat-induced epitope retrieval, tissue sections must be cooled at room temperature for 20 minutes before further processing. Thereafter staining of the tissue sections must be performed without delay.

## **Staining Procedure**

### **1. Reagent Preparation**

All reagents should be brought to ambient temperature (20 – 25 °C) before use in immunostaining. Accordingly all subsequent steps are to be performed at ambient temperature.

Care should be taken to avoid drying of the specimens during the immunostaining procedure as drying may lead to staining artefacts.

The following reagents should be prepared before starting with the staining procedure:

#### **1.1 Epitope Retrieval Solution**

Prepare the amount of Epitope Retrieval Solution sufficient for the staining procedure that is planned by dilution of a quantity of Vial 7 (Epitope Retrieval Solution 10X) 1:10 using distilled or deionized water.

After dilution, the Epitope Retrieval Solution may be stored at 2 – 8 °C for up to one month. The diluted solution must be discarded if it is cloudy in appearance.

**NOTE:** Use of water with elevated levels of ions for dilution of the epitope retrieval solution may significantly reduce the staining performance of the test. Please make sure that the water used is properly deionized (i.e. ensure that your ion exchange column for producing deionized water has been checked by routine maintenance). Do not use Tap water!

## 1.2 Wash Buffer

Use CINtec® Wash Buffer 10X, catalog number 10215364001, provided by Roche in combination with the CINtec® Histology Kit. For order details please refer to the website [www.roche.com](http://www.roche.com).

Prepare an amount of Wash Buffer sufficient for the wash steps of the staining procedure that is planned by diluting of a quantity of the Wash Buffer 10X, 1:10 using distilled or deionized water.

After dilution the Wash Buffer may be stored at 2 – 8 °C for up to one month. The diluted solution must be discarded if it is cloudy in appearance.

## 1.3 Substrate-Chromogen Solution (DAB)

For preparation of Substrate-Chromogen Solution, one drop of DAB Chromogen must be added to 2 mL of DAB Buffered Substrate. Proceed as follows:

- i) transfer 2 mL of DAB Buffered Substrate from Vial 5 to a test tube;
- ii) add one drop (25 – 30 µL) of DAB Chromogen from Vial 6. Mix and apply to tissue sections with a pipette.

2 mL of the Substrate-Chromogen Solution (DAB) prepared according to the instruction above is usually sufficient for staining five tissue sections including the corresponding five control specimens.

**NOTE:** Use the prepared Substrate-Chromogen Solution (DAB) within the same day.

**NOTE:** Addition of excess DAB Chromogen to the DAB Buffered Substrate will result in deterioration of the positive signal.

## 1.4 Counterstain

The DAB staining reaction results in a water insoluble coloured end product. Alcohol or water-based hematoxylin may be used for counterstaining. If used, adhere to the instructions provided by the supplier of the hematoxylin for performing the counterstaining.

## 1.5 Mounting Medium

For mounting slide specimens after staining, a non-aqueous, permanent mounting medium is recommended. However, aqueous mounting is also acceptable.

Eukitt Mounting Medium is recommended for non-aqueous mounting. Aquatex Merck is recommended for aqueous mounting.

## 2. Staining Procedure for Autostainer Instruments

The CINtec® Histology Kit has been adapted for use on the Autostainer Instruments (the Lab Vision Autostainer 480 or Dako Autostainer Plus) according to the template outlined below. It may be possible to use other instruments or systems with comparable function after appropriate validation by the user. Prior to staining on the Autostainer Instrument, the specimens and reagents should be prepared as stated in sections 1.1 – 1.5 and 2.1

### 2.1 Deparaffinization and Rehydration

Prior to deparaffinization, place slides in a drying oven at a temperature of no more

than 60 °C for at least 20 minutes but not more than one hour to quantitatively remove water thereby improving adherence of the tissue to the glass slide (“baking”) and to melt the paraffin. Tissue slides must be deparaffinized to remove embedding medium and must then be rehydrated before the staining procedure can be performed. It is crucial to avoid incomplete removal of paraffin as residual embedding medium will result in increased non-specific staining. Incubate the slides at ambient temperature (20 – 25 °C) according to the following steps.

- 5 (±1) minutes in a xylene bath;
- repeat this step once with a fresh bath;
- remove excess liquid;
- 3 (±1) minutes in 95% ethanol;
- repeat this step once with a fresh bath;
- remove excess liquid;
- 3 (±1) minutes in 70% ethanol;
- repeat this step once with a fresh bath;
- remove excess liquid;
- minimum of 30 seconds in distilled or deionized water.

Start staining procedure as described in Section 2.2, Step 1: Epitope retrieval.

Xylene and alcohol solutions should not be used for more than 40 slides.

**NOTE:** Users should note that variations in equipment temperature or exposure times during preanalytical sample preparation may lead to incomplete removal of paraffin from tissue slides. Residual paraffin can lead to incomplete staining with any histology stain, including the CINtec® Histology stain. Histopathology laboratories should include regular monitoring of the equipment to reduce variation in sample preparation before staining. Observation of sharply defined borders within immuno-reactive tissue areas or other staining inconsistencies within a slide may, with any immuno-histochemical stain, be an indicator of non-optimal or incomplete pre-analytical processing of the specimen. Users should consider checking the equipment and the preanalytical sample preparation methods if inconsistent staining is observed.

## 2.2 Staining Protocol for Autostainer Instruments

### Step 1: Epitope Retrieval

- fill staining jars, *e.g.*, plastic Coplin jars, with the diluted Epitope Retrieval Solution (see Procedure, Section 1.1);
- place staining jars containing Epitope Retrieval Solution in water bath and heat water bath and the Epitope Retrieval Solution to 95 – 99 °C. It is important to adjust the level of the water in the water bath to make sure that the jars are immersed in the water to a level of 80%. To stabilize the temperature and avoid evaporation, cover jars with lids;
- immerse deparaffinized sections into the preheated Epitope Retrieval Solution in the staining jars; this step usually will lower the temperature in the jars to less than 90 °C;
- bring the temperature of the water bath **and** the Epitope Retrieval Solution

in the jars back to 95 – 99 °C; check the temperature of the Epitope Retrieval Solution in the jars;

- incubate for 10 ( $\pm$ 1) minutes at 95 – 99 °C; start count down only after the temperature of the Epitope Retrieval Solution in the jars has been verified to have reached a temperature of 95 – 99 °C;
- remove the entire jar with slides from the water bath;
- allow the slides to cool in the Epitope Retrieval Solution for 20 ( $\pm$ 1) minutes at room temperature;
- decant the Epitope Retrieval Solution and rinse sections in the Wash Buffer (see Procedure, Section 1.2);
- for optimal performance, soak sections in Wash Buffer for 5 minutes after epitope retrieval and prior to staining.

**NOTE:** The Epitope Retrieval Solution is designed for single use application only. Do not re-use.

## Step 2: Programming of the Instrument

Prior to the first application of the CIntec® Histology Kit on an Autostainer Instrument, a new template needs to be set up. Please refer to the Operator's Manual for the dedicated Autostainer Instrument.

## Step 3: Autostainer Procedure

- transfer the reagents from the kit bottles into graduated Autostainer Reagent Vials. Use the Autostainer-generated map for program times and reagent volumes (see point 4 for specific times and volumes);
- place the Autostainer reagent vials in the Autostainer Reagent Rack according to the computer-generated Reagent Layout Map;
- load the slides onto the Autostainer according to the computer-generated Slide Layout Map;
- to prevent drying the specimens should be sprinkled with wash buffer after loading onto the Autostainer;
- the following is an outline of the program run:
  - rinse\*;
  - 200  $\mu$ L Peroxidase-Blocking Reagent - 5 minutes;
  - rinse\*;
  - 200  $\mu$ L Primary p16<sup>INK4a</sup> antibody or Negative Reagent Control - 30 minutes;
  - rinse\*;
  - 200  $\mu$ L Visualization Reagent - 30 minutes;
  - rinse\*;
  - rinse\*;
  - rinse\*;
  - switch;
  - 200  $\mu$ L Substrate-Chromogen Solution (DAB) - 10 minutes;
  - rinse\*;

- rinse slides in deionized water after the substrate-chromogen step.
- \*use Wash Buffer for the respective rinsing steps.

**NOTE:** If the Autostainer Instrument used rinses slides in buffer, the slides must be rinsed with deionized water after they have been removed from the Autostainer.

#### **Step 4: Counterstain (Instructions are for Hematoxylin)**

- immerse slides in a bath of hematoxylin. Incubate for 2 – 5 minutes, depending on the strength of hematoxylin used;
- place slides in a tap water bath and rinse gently with running tap water. Ensure all residual hematoxylin has been cleared;
- briefly rinse slides gently in a bath of distilled or deionized water;
- Counterstain may be performed directly on the Autostainer Instrument.

**NOTE:** Depending on the incubation length and potency of the hematoxylin used, counterstaining will result in a pale to dark blue colouration of the cell nuclei. Excessive or incomplete counterstaining may interfere with proper interpretation of results.

#### **Step 5: Mounting**

Non-aqueous, permanent mounting medium is recommended. Otherwise, aqueous mounting medium is also acceptable. Adhere to the instructions for use of the supplier for mounting medium.

**NOTE:** To minimize fading, protect slides from light and store at ambient temperature (20 – 25 °C).

### **3. Staining Procedure for Manual Use**

**NOTE:** Prevent drying of tissue sections during the staining procedure. Dried tissue sections may result in increased non-specific staining. For prolonged incubations keep tissues in a humid environment.

Adhere to standard procedures used in manual immunohistochemistry staining when using the CINtec® Histology Kit for manual staining.

#### **3.1 Deparaffinization and Rehydration**

Prior to deparaffinization, place slides in a drying oven at a temperature of no more than 60 °C for at least 20 minutes but not more than one hour to quantitatively remove water thereby improving adherence of the tissue to the glass slide (“baking”) and to melt the paraffin. Tissue slides must be deparaffinized to remove embedding medium and must then be rehydrated before the staining procedure can be performed. It is crucial to avoid incomplete removal of paraffin as residual embedding medium will result in increased non-specific staining. Incubate the slides at ambient temperature (20 – 25 °C) according to the following steps.

- 5 (±1) minutes in a xylene bath;
- Repeat this step once with a fresh bath;
- Remove excess liquid;
- 3 (±1) minutes in 95% ethanol;
- Repeat this step once with a fresh bath;

- Remove excess liquid;
- 3 ( $\pm$ 1) minutes in 70% ethanol;
- Repeat this step once with a fresh bath;
- Remove excess liquid;
- Minimum of 30 seconds in distilled or deionized water.

Start staining procedure as described in Section 3.2, Step 1: Epitope retrieval. Xylene and alcohol solutions should not be used for more than 40 slides.

## 3.2 Staining Protocol for Manual Use

### Step 1: Epitope Retrieval

- fill staining jars, e.g., plastic Coplin jars, with the diluted Epitope Retrieval Solution (see Procedure, Section 1.1);
- place staining jars containing Epitope Retrieval Solution in water bath and heat water bath and the Epitope Retrieval Solution to 95 – 99 °C. At that step it is important to adjust the level of the water in the water bath to make sure that the jars are immersed in the water to a level of 80%. To stabilize the temperature and avoid evaporation, cover jars with lids;
- immerse deparaffinized sections into the preheated Epitope Retrieval Solution in the staining jars; this step usually will lower the temperature in the jars to less than 90°C;
- bring the temperature of the water bath **and** the Epitope Retrieval Solution in the jars back to 95 – 99 °C; check the temperature of the Epitope Retrieval Solution in the jars;
- incubate for 10 ( $\pm$ 1) minutes at 95 – 99 °C; start count down only after the temperature of the Epitope Retrieval Solution in the jars has been verified to have reached a temperature of 95 – 99 °C;
- remove the entire jar with slides from the water bath;
- allow the slides to cool in the Epitope Retrieval Solution for 20 ( $\pm$ 1) minutes at room temperature;
- decant the Epitope Retrieval Solution and rinse sections in the diluted Wash Buffer (see Procedure, Section 1.2);
- for optimal performance, soak sections in Wash Buffer for 5 ( $\pm$ 1) minutes after epitope retrieval and prior to staining.

**NOTE:** The Epitope Retrieval Solution is designed for single use application only. Do not re-use.

### Step 2: Peroxidase-Blocking Reagent

- apply 200  $\mu$ L Peroxidase-Blocking Reagent to cover specimen;
- incubate for 5 ( $\pm$ 1) minutes;
- tap off excess liquid and place slides in a fresh Wash Buffer bath for 5 ( $\pm$ 1) minutes.

### Step 3: Primary Antibody or Negative Reagent Control

- remove excess buffer;
- cover specimen with 200 µL of primary antibody (Mouse Anti-Human p16<sup>INK4a</sup> or Negative Reagent Control);
- incubate for 30 (±1) minutes;
- tap off excess liquid and place slides in a fresh Wash Buffer bath for 5 (±1) minutes.

### Step 4: Visualization Reagent

- remove excess buffer;
- cover specimen with 200 µL of Visualization Reagent;
- incubate for 30 (±1) minutes;
- tap off excess liquid and place slides in a fresh buffer bath for 5 (±1) minutes;
- repeat this step twice with a fresh Wash Buffer bath.

### Step 5: Substrate-Chromogen Solution (DAB)

- cover specimen with 200 µL of Substrate-Chromogen Solution (DAB) that has been prepared according to the procedure described in 1.3 above;
- incubate for 10 (±1) minutes;
- tap off excess liquid and rinse gently with distilled or deionized water.

Collect Substrate-Chromogen Solution (DAB) waste in a hazardous materials container for proper disposal.

### Step 6: Counterstain (Instructions are for Hematoxylin)

- immerse slides in a bath of hematoxylin. Incubate for 2 – 5 minutes, depending on the strength of hematoxylin used;
- place slides in a tap water bath and rinse gently with running tap water. Ensure all residual hematoxylin has been cleared;
- briefly rinse slides gently in a bath of distilled or deionized water.

**NOTE:** Depending on the incubation length and potency of the hematoxylin used, counterstaining will result in a pale to dark blue colouration of the cell nuclei. Excessive or incomplete counterstaining may interfere with proper interpretation of results.

### Step 7: Mounting

Non-aqueous, permanent mounting medium is recommended. For xylene-based permanent mounting media a dehydration procedure is necessary, e.g.

- distilled or deionized water
- 3 min 70% Ethanol



- 3 min 70% Ethanol
- 3 min 96% Ethanol
- 3 min 99% Ethanol
- 5 min Xylene
- 5 min Xylene

Otherwise, aqueous mounting medium is also acceptable. Adhere to instructions of use of the supplier for mounting medium.

**NOTE:** To minimize fading, protect slides from light and store at ambient temperature (20 – 25 °C).

## VII. Quality Control

Deviations from the recommended procedures for specimen fixation and processing in the user's laboratory may produce significant variability in results, necessitating regular performance of in-house controls.

### Positive Tissue Control

External Positive control materials should be fresh autopsy/biopsy/surgical specimens fixed, processed and embedded as soon as possible in the same manner as the patient sample(s). Positive tissue controls are indicative of correctly prepared tissues and proper staining techniques. One positive external tissue control for each set of test conditions should be included in each staining run. The tissues used for the external positive control materials should be selected from patient specimens with known positive staining for p16<sup>INK4a</sup>. If the positive controls fail to demonstrate appropriate positive staining, results with the test specimens should be considered invalid.

### Negative Tissue Control

Cell types known to be negative for p16<sup>INK4a</sup> present in most tissue sections can be used by the laboratorian as internal negative control sites to verify the IHC's performance specifications. If specific staining (false positive staining) occurs in the negative tissue control, results with the patient specimens should be considered invalid.

### Non-specific Negative Reagent Control

Use the nonspecific Negative Reagent Control in place of the primary antibody with a section of each patient specimen to evaluate non-specific staining and allow better interpretation of specific staining at the antigen site.

If specific staining (false positive staining) occurs with the non-specific Negative Reagent Control, results with the patient specimens should be considered invalid.

## VIII. Interpretation of Results

Control Specimens stained with the Negative Reagent Control as a primary reagent must not show specific staining.

Positive staining using the monoclonal mouse anti-Human p16<sup>INK4a</sup> antibody, Clone E6H4™, should be assessed within the context of any non-specific background staining of the Negative Reagent Control. As with any immunohistochemical test, a negative result means that the antigen was not detected, not that the antigen was

absent in the cells/tissue assayed.

Interpretation of the results must take into consideration the fact that p16<sup>INK4a</sup> is a cellular protein that may be expressed at detectable levels in cervical high-grade dysplastic lesions and cervical cancers as well as in some conditions not associated with cervical dysplasia, albeit at differing levels and with different patterns of expression.

The stained slide specimens are evaluated according to a binary rating system composed of the ratings “positive” and “negative”.

The rating “positive” is assigned if the p16<sup>INK4a</sup>-stained slide specimen shows a continuous staining of cells of the basal and parabasal cell layers of the squamous cervical epithelium, with or without staining of cells of superficial cell layers (“diffuse staining pattern”). An example for a slide rated as “positive” (“diffuse staining pattern”) is shown in Annex 2, Figure 1.

The rating “negative” is assigned if the p16<sup>INK4a</sup>-stained slide specimen shows either a negative staining reaction in the squamous epithelium (“negative staining pattern”) or a staining of isolated cells or small cell clusters; *i.e.*, a non-continuous staining, particularly not of the basal and parabasal cells (“focal staining pattern”). An example for a slide rated as “negative” (“focal staining pattern”) is shown in Annex 2, Figure 2.

The interpretation of slides stained for p16<sup>INK4a</sup> using the CINtec® Histology Kit should be performed in conjunction with H&E-stained slides prepared from the same cervical tissue specimen. The additional information provided by the CINtec®-stained slides should be combined with the preliminary morphology-based diagnosis established on the H&E-stained slides in order to create a final diagnosis.

## IX. Limitations

- For professional use only. Special training is required for the performance of immunohistochemistry procedures.
- The clinical interpretation of any positive or negative staining should be evaluated within the context of clinical presentation, morphology, and other histopathological criteria. The clinical interpretation of any positive or negative staining should be complemented by morphological studies using proper positive and negative internal and external controls as well as other diagnostic tests. It is the responsibility of a qualified pathologist who is familiar with the proper use of IHC antibodies, reagents, and methods to interpret all of the steps used to prepare and interpret the final IHC preparation.
- The staining results in immunohistochemistry are strongly influenced by the quality of the tissue stained. Accordingly, the steps of fixation, washing, drying, heating, sectioning, or contamination with other tissues significantly contribute to the overall result of the staining and may lead to artefacts, antibody trapping, or false-negative results. Inconsistent results may be due to variations in fixation and embedding methods or to inherent irregularities within the tissue.
- Excessive or incomplete counterstaining may interfere with proper interpretation of results.

- The manufacturer provides these antibodies/reagents at optimal dilution for use according to the instructions provided herein, for IHC testing on prepared tissue sections. Any deviation from the recommended test procedures may invalidate declared expected results; appropriate controls must be employed and documented. Users who deviate from the recommended test procedures must accept responsibility for interpretation of patient results under these circumstances.
- False-positive results may be seen due to non-immunological binding of proteins or substrate reaction products. They may also be caused by pseudoperoxidase activity (erythrocytes) and endogenous peroxidase activity (cytochrome C).
- Reagents may demonstrate unexpected reactions in previously untested tissues. The possibility of unexpected reactions even in tested tissue groups cannot be completely eliminated due to biological variability of antigen expression in neoplasms, or other pathological tissues. Contact Roche mtm laboratories AG with documented unexpected reaction(s). For information on technical support, please refer to the contact information provided in section XIII.
- Do not replace kit reagents with reagents carrying other lot numbers or with reagents from other manufacturers.

## **X. Performance characteristics**

### Clinical performance

The clinical performance of the CINtec® Histology Kit has been evaluated in a controlled clinical study using formalin-fixed, paraffin-embedded cervical tissue specimens [21]. The study was designed to demonstrate the appropriateness of the CINtec® Histology Kit as an aid to increase the diagnostic accuracy and inter-observer agreement for the detection of high-grade cervical intra-epithelial neoplasia (CIN2+).

The clinical study was performed on retrospectively collected cervical punch and conization biopsies. A total of 500 cervical specimens retrieved from two different European pathology laboratories and enriched for high-grade dysplasia based on their original diagnoses were used for preparation of H&E-stained slides and consecutive slides stained with the CINtec® Histology Kit according to the manufacturer's instructions.

Three European expert gynecopathologists established their independent diagnosis on each case based on the H&E-stained slide. Cases with discrepant results were subjected to a joint second review process during an adjudication meeting, and majority diagnoses (consensus diagnoses, two out of three) served as the reference diagnoses for the study.

Twelve investigators (certified pathologists reading cervical pathology on a regular basis) from 4 European countries (France, Italy, Spain, and Germany) participated as panel pathologists in the study. In a first round of slide interpretations, all panel pathologists established their individual diagnoses for each case based on H&E-stained slides only. The pathologists were blinded to the original diagnoses and the reference diagnoses at all times. After a wash-out period of more than 4 weeks, the

same set of (re-labeled) H&E-stained slides were reviewed again by all twelve panel pathologists, conjunctively with the matched slides of each case stained with the CINtec® Histology Kit.

For the assessment of improvements in diagnostic accuracy for the CINtec® Histology-stained cervical slides conjunctively read with the respective H&E-stained slides in comparison to the H&E-stained cervical slides alone, the results from each panel pathologist’s readings for both methods were compared to the consensus diagnoses established as reference diagnoses by the three expert gynecopathologists.

Results:

A total of 482 cases with complete diagnoses by all study pathologists were included in the data analysis. The frequencies of the various diagnostic categories based on the consensus diagnoses by the three expert gynecopathologists was Negative for Dysplasia (n=194), CIN1 (n=96), CIN2 (69), and CIN3 (n=123).

Diagnostic accuracy for the detection of CIN2+

The overall sensitivity for the identification of CIN2+ was increased from 1787 (H&E) to 2018 (H&E plus CINtec® Histology) true positive CIN2+ results, with only a slight reduction in overall specificity from 3088 (H&E) to 3051 (H&E plus CINtec® Histology) true negative ≤CIN1 results.

Table 1

Improvement of diagnostic accuracy for high-grade CIN (CIN2+) by the conjunctive reading of H&E plus CINtec® Histology-stained slides versus H&E-stained slides alone; numbers of true-positive, false-negative, false- - positive and true-negative reading results as compared to the consensus diagnoses by the expert pathologists are given (Note: total agreement with the consensus diagnoses of the expert pathologists would have resulted in 2304 (192 CIN2+ cases, x 12 panel pathologists) true-positive results)

|                          | True<br>Positives | False<br>Negatives | False<br>Positives | True<br>Negatives |
|--------------------------|-------------------|--------------------|--------------------|-------------------|
| H&E                      | 1787              | 517                | 392                | 3088              |
| H&E+CINtec®<br>Histology | 2018              | 286                | 429                | 3051              |

Mixed-Effects ANOVA Model for Jackknife Pseudovalues, commonly known as Dorfman-Berbaum-Metz (DBM) method, was used to test for the effects of the CINtec® Histology test on accuracy.

The null hypothesis that the diagnostic accuracy based on H&E versus H&E plus CINtec® Histology-stained slides for CIN2+ is equal was rejected with p=0.0004 (area under the curve (AUC) for H&E: 0.877; AUC for H&E plus CINtec® Histology: 0.925).

Table 2

Performance characteristics for the conjunctive reading of H&E plus CINtec® Histology-stained slides versus H&E-stained slides alone for the identification of high-grade CIN (CIN2+); consensus diagnoses established

by three expert gynecopathologists on H&E-stained slides used as reference diagnoses

|                    | <b>H&amp;E</b><br>% (95% CI) | <b>H&amp;E plus CINtec® Histology</b><br>% (95% CI) |
|--------------------|------------------------------|---|
| <b>Sensitivity</b> | 77.6% (75.8, 79.3)           | 87.6% (86.2, 88.9)                                  |
| <b>Specificity</b> | 88.7% (87.6, 89.8)           | 87.7% (86.5, 88.8)                                  |
| <b>PPV</b>         | 82.0% (80.3, 83.6)           | 82.5% (80.9, 84.0)                                  |
| <b>NPV</b>         | 85.7% (84.5, 86.8)           | 91.4% (90.4, 92.4)                                  |
| <b>DLR+</b>        | 6.885 (6.256, 7.578)         | 7.105 (6.494, 7.773)                                |
| <b>DLR–</b>        | 0.253 (0.234, 0.273)         | 0.142 (0.127, 0.158)                                |

95% CI, 95% confidence intervals; PPV, positive predictive value; NPV, negative predictive value; DLR+, positive diagnostic likelihood ratio; DLR–, negative diagnostic likelihood ratio

Using the consensus diagnoses of the three expert gynecopathologists as the reference, the sensitivity for the detection of CIN2+ was improved by relative 13% (sensitivity for H&E diagnoses of 77.6% increased to 87.6% for diagnoses established on H&E plus CINtec® Histology-stained slides).

The increase in diagnostic accuracy for the detection of CIN2+ was independently proven with statistical significance for the subgroups of cervical punch biopsies (n=249; AUC for H&E: 0.895; AUC for H&E plus CINtec® Histology: 0.929; p=0.0053) and cervical conization specimens (n=233; AUC for H&E: 0.887; AUC for H&E plus CINtec® Histology: 0.948; p=0.009).

#### Inter-observer agreement for the detection of CIN2+

For the assessment of the improvement of the inter-observer agreement for the detection of CIN2+ between panel pathologists, a multiple-rater form of the Kappa statistics between diagnoses performed on H&E-stained cervical slides and diagnoses performed on H&E-stained slides conjunctively read with CINtec® Histology-stained cervical slides was used.

*Table 3*

Improvement of inter-observer agreement for the detection of CIN2+ by the conjunctive reading of H&E plus CINtec® Histology-stained slides versus H&E-stained slides alone. Kappa values as a chance-corrected measure of agreement are given

| <b>Diagnostic Category</b>      | <b>Kappa H&amp;E</b> | <b>Kappa H&amp;E plus CINtec®</b> | <b>Statistical significance</b> |
|---------------------------------|----------------------|-----------------------------------|---------------------------------|
| CIN2+, all cases                | 0.580                | 0.756                             | p<0.0001                        |
| CIN2+, punch biopsies only      | 0.598                | 0.748                             | p<0.0001                        |
| CIN2+, conization biopsies only | 0.548                | 0.765                             | p<0.0001                        |

The kappa statistic as a chance-corrected measure of agreement between the

panel pathologists for the detection of CIN2+ was improved significantly for the H&E plus CINtec® Histology-stained slide readings, versus H&E slide reading alone (Increase in kappa values for all cases from 0.580 to 0.756;  $p < 0.0001$ ).

### Reproducibility in rating the p16<sup>INK4a</sup> staining pattern

The reliability among the panel pathologists to assess the p16<sup>INK4a</sup> staining pattern in cervical tissue specimens as diffusely p16<sup>INK4a</sup>-positive, focally p16<sup>INK4a</sup>-positive, or negative for p16<sup>INK4a</sup> was evaluated.

There was a high level of reproducibility between pathologists to rate the p16<sup>INK4a</sup> -staining pattern as either positive (diffuse staining pattern), or negative (focal staining pattern, or no immuno-reactivity). The Kappa value (chance corrected measure of agreement) for reproducibility between the twelve panel pathologists for rating the p16<sup>INK4a</sup> staining pattern as either positive (diffuse staining pattern) or negative (focal staining pattern, and no immunoreactivity) was found to be excellent (mean kappa = 0.899; median kappa = 0.903).

### Analytical performance

#### Analytical Sensitivity

The sensitivity study was conducted with CINtec® Histology Kits from 3 lots. Thirty-five cervical biopsies with CIN3+ were analyzed in the Sensitivity Study, including 6 biopsies classified as CIN3 and 29 biopsies classified as squamous cell carcinoma.

Five out of 6 (83.3%) of cases classified as CIN3 showed strong diffuse p16 staining. One CIN3 biopsy case showed weak p16 expression but missed a clear p16-positive signal.

Twenty-two out of 29 (75.9%) of cases classified as squamous cell carcinoma showed strong diffuse p16 staining. Five cases of squamous cell carcinoma showed weak p16 expression but missed a clear p16-positive signal. Only 2 out of 29 (6.9%) of the cases were completely negative for p16. This can be explained as it was previously reported that 3-8% of squamous cell carcinoma cases were negative for p16 staining [1; 6]. This may reflect a small percentage of cases where dedifferentiation and chromosomal rearrangements led to inactivation or deletion of the p16 gene locus.

For all 3 batches staining on a well-characterized panel of rat tissue sections using the Negative Reagent Control resulted in strong and specific staining of single neurons in rat brain. In addition, some canalicular cells of the kidney, single macrophages in the spleen, and single macrophages and plasma cells in the small intestine showed weak staining. Other tissues of the rat tissue array were completely negative.

The anti-p16 antibody E6H4™ is able to render continuous staining of cells of the basal and parabasal cell layers, with or without staining of cells of intermediate or intermediate and superficial layers of the squamous epithelium in cervical biopsies with high-grade CIN (CIN2, CIN3). The Negative Reagent Control (NRC) is able to detect specifically neurons of the rat.

#### Analytical Specificity

The specificity of the mouse anti-human p16<sup>INK4a</sup> antibody clone E6H4™ has been verified by Western blot analysis (positive for HeLa cell line lysate; also refer to [1]).

The specificity study was conducted with CINtec® Histology Kits from 3 lots on a well-characterized panel of 90 NBF (neutral buffered formalin) fixed normal tissue samples (30 different types of tissue) and 54 tumor tissues other than cervix uteri (arranged in Multi Tissue Arrays = MTAs).

For staining using the p16-specific antibody (clone E6H4™) on normal tissues, 16 different p16-negative tissues [cerebrum, cerebellum, adrenal gland, thyroid, bone marrow, heart, esophagus, stomach, intestine, colon, liver, kidney, striated muscle, skin, mesothelium, cervix] and 14 different p16-positive tissues [Weak staining: hypophysis, lung, thymus gland, prostate; positive: nerve, intestine, tonsil, pancreas, spleen; strongly positive: uterus, ovary, breast, testis, parathyroid] were observed. For staining using the p16-specific antibody (clone E6H4™) on tumor tissues other than cervix uteri, 22 different p16-negative tumor cases and 32 different p16-positive tumor cases were observed.

For all 3 lots of CINtec® Histology the results were negative for all tissues (normal tissues and tumor tissues) tested when using the Negative Reagent Control. The mouse monoclonal anti-rat oxytocin-related neurophysin antibody does not react significantly with human normal tissue samples and human tumor tissues.

### Reproducibility

#### Inter-Run Reproducibility

Inter-Run reproducibility was determined for the CINtec® Histology Kit using the manual protocol by staining 36 slides prepared from tissue blocks diagnosed as CIN2+. The dysplastic areas on all sections were stained with comparable intensity ( $\pm 0.5$  on a scale of 0 to 3) across all runs. Normal areas on all slides did not show specific staining.

#### Intra-Run Reproducibility

Intra-Run reproducibility was determined with the CINtec® Histology Kit with the manual protocol and on the autostainer instrument at three different days. In total three tissue blocks diagnosed as CIN2+ were included. From each block one serial section was included at each day. The dysplastic area on the sections from one block was stained with comparable intensity ( $\pm 0.5$  on a scale of 0 to 3) on all days and also in manual and automated staining. Normal areas on all slides did not show specific staining.

#### Inter-Lot Reproducibility

For determining inter-lot reproducibility CINtec® Histology Kits from 3 different lots were used for staining of sections from CIN2+ cases. Manual and automated staining was performed in accordance to the protocol given in the instructions for use. The dysplastic area on the sections from one block was stained with comparable intensity ( $\pm 0.5$  on a scale of 0 to 3) with the reagents from all three Lots in manual and automated staining. Normal areas on all slides did not show specific staining.

Note that the staining intensity scoring method with a scale of 0 to 3 was solely used for analytical performance evaluation purposes and shall not be used for the interpretation of the staining of tissue sections in clinical practice. Instead, the qualitative interpretation of stained slides as described in Section VIII shall be used for the routine interpretation.

## XI. Troubleshooting

Refer to section XIII. for contact details in case technical assistance is required.

| Problem                           | Probable Cause  | Suggested Action  |
|-----------------------------------|---|---|
| <b>1. No staining of slides</b>   | <b><u>1a.</u></b> Deviation from Instructions for Use;  | <b><u>1a.</u></b> Carefully read Instructions for Use and adhere to the procedures outlined therein;  |
| <b>2. Weak staining of slides</b> | <b><u>2a.</u></b> Inadequate epitope retrieval;   | <b><u>2a.</u></b> Use freshly prepared Epitope Retrieval Solution and / or make sure that Epitope Retrieval Solution reaches 95 – 99 °C for a full 10 minutes and is allowed to cool for an additional 20 minutes;  |
|                                   | <b><u>2b.</u></b> Inadequate reagent incubation times;  | <b><u>2b.</u></b> Review 2.2. / 3.2. Staining protocol recommendations;   |
|                                   | <b><u>2c.</u></b> Inappropriate fixation method;  | <b><u>2c.</u></b> Ensure that patient tissue is not over-fixed or that no alternative fixative was used;  |
|                                   | <b><u>2d.</u></b> Water that has been used to dilute the epitope-retrieval solution contains an ion concentration that is too high; | <b><u>2d.</u></b> Ensure that your ion exchange column for producing deionized water has been checked by routine maintenance;   |
|                                   | <b><u>2e.</u></b> Inappropriate deparaffinization;  | <b><u>2e.</u></b> Users should note that variations in equipment temperature or exposure times during preanalytical sample preparation may lead to incomplete removal of paraffin from tissue slides. Residual paraffin can lead to incomplete staining with any histology stain, including the CINtec® Histology stain.<br><br>Histopathology laboratories should include regular monitoring of the equipment to reduce variation in sample preparation before staining. Observation of sharply defined borders within immuno-reactive tissue areas or other staining inconsistencies within a slide may, with any immuno-histochemical stain, be an indicator of non-optimal or incomplete pre-analytical processing of the specimen. Users should consider checking the equipment and the preanalytical sample preparation methods if inconsistent staining is observed; |



|   |   |  |
|---|---|--|
| <b>3. Excessive background staining of slides</b> | <b>3a.</b> Incomplete removal of paraffin;                                  | <b>3a.</b> Use fresh xylene baths and follow procedure as outlined in Section 2.1. / 3.1.;                             |
|   | <b>3b.</b> Mounting of sections to slides performed using starch additives; | <b>3b.</b> Starch additives used in mounting sections may exhibit immunoreactivity and should therefore be avoided;    |
|   | <b>3c.</b> Insufficient rinsing of slides;                                  | <b>3c.</b> Use fresh solution in buffer baths and wash bottles;  |
|   | <b>3d.</b> Drying of sections during staining procedure;                    | <b>3d.</b> Use humidity chamber. Wipe only three to four slides at a time before applying reagent;                     |
|   | <b>3e.</b> Inappropriate fixation method;                                   | <b>3e.</b> Use only fixative as recommended herein. Aberrantly fixed tissue may exhibit excessive background staining; |
|   | <b>3f.</b> Non-specific binding of reagents to tissue;                      | <b>3f.</b> Check fixation method of the specimen and presence of necrosis;   |
| <b>4. Tissue detaches from slides</b>             | <b>4a.</b> Use of inappropriate slides;                                     | <b>4a.</b> Adhere to recommendation herein and use SuperFrost® Plus slides;  |
| <b>5. Excessively strong specific staining</b>    | <b>5a.</b> Inappropriate fixation method;                                   | <b>5a.</b> Ensure proper fixative and fixation method;   |
|   | <b>5b.</b> Prolonged reagent incubation times;                              | <b>5b.</b> Review and adhere to staining protocol given in sections 2.2 / 3.2 above;                                   |
|   | <b>5c.</b> Inappropriate wash solution;                                     | <b>5c.</b> Use the Wash Buffer (10 x) (catalog number 10215364001).  |

## XII. Symbols

**Symbol:**

**Explanation:**



Catalog number



Batch code



Global Trade Item Number



Unique Device Identifier



In vitro diagnostic medical device



Manufacturer



Contains sufficient for <n> tests



Consult instructions for use



Use by



Temperature limitation



Date of manufacturing



Do not re-use



Contact for technical support (Telephone)

Contains materials of animal origin

Content

### XIII. Manufacturer

**Manufactured by:** Roche Diagnostics GmbH  
Sandhofer Strasse 116  
68305 Mannheim  
Germany

<https://navifyportal.roche.com>

**Contact for technical support (Telephone):** +800 5505 6606

The summary of safety and performance can be found here:

<https://ec.europa.eu/tools/eudamed>

### XIV. Revision Status

The current Instructions for Use represent Version 1.0 released June 2024.

Changes to previous version for product 9511 / 06594441001 (2.9, released December 2022):

- Manufacturer changed from Roche mtm laboratories AG to Roche Diagnostics GmbH, and thus new product number, 10213370001, issued
- Editorial changes

### XV. Intellectual Property

CINtec and E6H4 are trademarks of Roche.

All other trademarks are the property of their respective owners.

© 2024 Roche

## Annex 1 References

1. Klaes R, Friedrich T, Spitkovsky D, Ridder R, Rudy W, Petry U, Dallenbach-Hellweg G, Schmidt D, and von Knebel Doeberitz M. Overexpression of p 16<sup>INK4a</sup> as a specific marker for dysplasia and neoplastic epithelial cells of the cervix uteri. *Int J Cancer* 2001, 92(2):276-84
2. Sano T, Oyama T, Kashiwabara K, Fukuda T, and Nakajima T. Expression status of p16 protein is associated with human papillomavirus oncogenic potential in cervical and genital lesions. *Am J Pathol* 1998, 153(6):1741-8
3. von Knebel Doeberitz M. New markers for cervical dysplasia to visualise the genomic chaos created by aberrant oncogenic papillomavirus infections. *Eur J Cancer* 2002, 38(17):2229-42
4. Klaes R, Benner A, Friedrich T, Ridder R, Herrington S, Jenkins D, Kurman RJ, Schmidt D, Stoler M, and von Knebel Doeberitz M. p16<sup>INK4a</sup> immunohistochemistry improves interobserver agreement in the diagnosis of cervical intraepithelial neoplasia. *Am J Surgical Pathol* 2002, 26(11):1389-99
5. Wang SS, Trunk M, Schiffman M, Herrero R, Sherman ME, Burk RD Hildesheim A, Concepcion Bratti M, Wright TC, Rodriguez AC, Chen S, Reichert A, von Knebel Doeberitz C, Ridder R, and von Knebel Doeberitz M. Validation of p16<sup>INK4a</sup> as a marker of oncogenic human papillomavirus infection in cervical biopsies from a population-based cohort in Costa Rica. *Cancer Epidemiol Biomarkers Prev* 2004, 13(8):1355-60
6. Agoff SN, Lin P, Morihara J, Mao C, Kiviat NB, and Koutsky LA. p16<sup>INK4a</sup> expression correlates with degree of cervical neoplasia: a comparison with Ki-67 expression and detection of high-risk HPV types. *Modern Pathol* 2003, 16(7):665-73
7. Negri G, Egarter-Vigl E, Kasal A, Romano F, Haitel A, and Mian C. p16<sup>INK4a</sup> is a useful marker for the diagnosis of adenocarcinoma of the cervix uteri and its precursors: an immunohistochemical study with immunocytochemical correlations. *Am J Surg Pathol* 2003, 27(2):187-93
8. Negri G, Vittadello F, Romano F, Kasal A, Rivasi F, Girlando S, Mian C, and Egarter-Vigl E. p16<sup>INK4a</sup> expression and progression risk of low-grade intraepithelial neoplasia of the cervix uteri. *Virchows Arch* 2004, 445(6):616-20
9. Schorge JO, Lea JS, Elias KJ, Rajanbabu R, Coleman RL, Miller DS, and Ashfaq R. p16<sup>INK4a</sup> as molecular biomarker of cervical adenocarcinoma. *Am J Obstet Gynecol* 2004, 190(3):668-73
10. Christal JL, and Valente PT. The utility of p16 immunohistochemistry in the diagnosis of cervical intraepithelial neoplasia. *Pathol Case Rev* 2006, 11(3):117-20
11. Kong CS, Balzer BL, Troxell ML, Patterson BK, and Longacre TA. p16<sup>INK4a</sup> Immunohistochemistry is Superior to HPV In Situ Hybridization for the Detection of High-risk HPV in Atypical Squamous Metaplasia. *Am J Surg Pathol* 2007, 31(1):33-43
12. Nucci MR, and Crum CP. Redefining Early Cervical Neoplasia: Recent Progress. *Adv Anat Pathol* 2007, 14(1):1-10
13. Ishikawa M, Fujii T, Saito M, Nindl I, Ono A, Kubushiro K, Tsukazaki K, Mukai M, and Nozawa S. Overexpression of p16<sup>INK4a</sup> as an indicator for human papillomavirus oncogenic activity in cervical squamous neoplasia. *Int J Gynecol Cancer* 2006, 16(1):347-53

14. Kalof AN, Evans MF, Simmons-Arnold L, Beatty BG, and Cooper K. p16<sup>INK4a</sup> Immunoeexpression and HPV In Situ Hybridization Signal Patterns Potential Markers of High-Grade Cervical Intraepithelial Neoplasia. *Am J Surg Pathol* 2005, 29(5):674-9
15. Ansari-Lari MA, Staebler A, Zaino RJ, Shah KV, and Ronnett BM. Distinction of Endocervical and Endometrial Adenocarcinomas: Immunohistochemical p16 Expression Correlated With Human Papillomavirus (HPV) DNA Detection. *Am J Surg Pathol* 2004, 28(2):160-7
16. Regauer S, and Reich O. CK17 and p16 expression patterns distinguish (atypical) immature squamous metaplasia from high-grade cervical intraepithelial neoplasia (CIN III). *Histopathology* 2007, 50(5):629-35
17. Stoler MH, and Schiffman M. Interobserver reproducibility of cervical cytologic and histologic interpretations. *JAMA* 2001, 285(11):1500-5
18. McCluggage WG, Walsh MY, Thornton CM, Hamilton PW, Date A, Caughley LM, and Bharucha H. Inter- and intra-observer variation in the histopathological reporting of intraepithelial lesions using a modified Bethesda grading system. *Br J Obstet Gynaecol* 1998, 105(2):206-10
19. Ismail SM, Colclough AB, Dinnen JS, Eakins D, Evans DMD, Gradwell E, O'Sullivan JPO, Summerell JM, and Newcombe R. Reporting cervical intraepithelial neoplasia (CIN): intra- and interpathologist variation and factors associated with disagreement. *Histopathology* 1990, 16(4):371-6
20. De Vet HCW, Knipschild PG, Schouten HJ, Koudstaal J, Kwee W-S, Willebrand D, Sturmans F, and Arends JW. Interobserver variation in histopathological grading of cervical dysplasia. *J Clin Epidemiol* 1990, 43(12):1395-8
21. Bergeron C, Ordi J, Schmidt D, Trunk MJ, Keller T, and Ridder R, for the CINtec Histology Study Group. Conjunctive p16<sup>INK4a</sup> testing significantly increases accuracy in diagnosing high-grade cervical intraepithelial neoplasia. *Am J Clin Pathol* 2010, 133(3):395-406
22. Sayed K, Korourian S, Ellison DA, Kozlowski K, Talley L, Horn HV, et al. Diagnosing cervical biopsies in adolescents: the use of p16 immunohistochemistry to improve reliability and reproducibility. *J Low Genit Tract Dis* 2007, 11(3):141-6
23. Gurrola-Diaz CM, Suarez-Rincon AE, Vazquez-Camacho G, Buonocunto-Vazquez G, RosalesQuintana S, Wentzensen N, et al. p16<sup>INK4a</sup> immunohistochemistry improves the reproducibility of the histological diagnosis of cervical intraepithelial neoplasia in cone biopsies. *Gynecol Oncol* 2008, 111(1):120-4
24. Horn LC, Reichert A, Oster A, Arndal SF, Trunk MJ, Ridder R, et al. Immunostaining for p16<sup>INK4a</sup> used as a conjunctive tool improves interobserver agreement of the histologic diagnosis of cervical intraepithelial neoplasia. *Am J Surg Pathol* 2008, 32(4):502-12
25. Haidopoulos D, Partsinevelos GA, Vlachos GD, Rodolakis A, Markaki S, Voulgaris Z, et al. p16<sup>INK4A</sup> is a strong biomarker for cervical intraepithelial neoplasia and invasive cervical carcinoma: a reappraisal. *Reprod Sci* 2009, 16(7):685-93
26. Lesnikova I, Lidang M, Hamilton-Dutoit S and Koch J. p16 as a diagnostic marker of cervical neoplasia: a tissue microarray study of 796 archival specimens. *Diagn Pathol* 2009, 4:22
27. Ordi J, Garcia S, del Pino M, Landolfi S, Alonso I, Quinto L, et al. p16<sup>INK4a</sup> immunostaining identifies occult CIN lesions in HPV-positive women. *Int J Gynecol Pathol* 2009, 28(1):90-7

28. Galgano MT, Castle PE, Atkins KA, Brix WK, Nassau SR and Stoler MH. Using biomarkers as objective standards in the diagnosis of cervical biopsies. *Am J Surg Pathol* 2010, 34(8):1077-87
29. Reuschenbach M, Seiz M, von Knebel Doeberitz C, Vinokurova S, Duwe A, Ridder R, et al. Evaluation of cervical cone biopsies for coexpression of p16<sup>INK4a</sup> and Ki-67 in epithelial cells. *Int J Cancer* 2012, 130(2):388-94
30. Darragh TM, Colgan TJ, Cox JT, Heller DS, Henry MR, Luff RD, McCalmont T, Nayar R, Palefsky JM, Stoler MH, Wilkinson EJ, Zaino RJ, Wilbur DC; Members of LAST Project Work Groups. The Lower Anogenital Squamous Terminology Standardization Project for HPV-Associated Lesions: background and consensus recommendations from the College of American Pathologists and the American Society for Colposcopy and Cervical Pathology. *J Low Genit Tract Dis* 2012, 16(3):205-42. Erratum in: *J Low Genit Tract Dis*. 2013, 17(3):368
31. Liao GD, Sellors JW, Sun HK, Zhang X, Bao YP, Jeronimo J, et al. p16<sup>INK4A</sup> immunohistochemical staining and predictive value for progression of cervical intraepithelial neoplasia grade 1: a prospective study in China. *Int J Cancer* 2014, 134(7):1715-24
32. Reuschenbach M, Wentzensen N, Dijkstra MG, von Knebel Doeberitz M and Arbyn M. p16<sup>INK4a</sup> immunohistochemistry in cervical biopsy specimens: A systematic review and meta-analysis of the interobserver agreement. *Am J Clin Pathol* 2014, 142(6):767-72
33. Shah AA, Jeffus SK, Zhao Z, Stoler MH and Stelow EB. Adjunct p16(INK4a) immunohistochemistry aids the detection of high-grade squamous intraepithelial lesions in endocervical curettage specimens. *Am J Clin Pathol* 2014, 141(3):342-7
34. Bergeron C, Ronco G, Reuschenbach M, Wentzensen N, Arbyn M, Stoler M, et al. The clinical impact of using p16(INK4a) immunochemistry in cervical histopathology and cytology: an update of recent developments. *Int J Cancer* 2015, 136(12):2741-51
35. Clinton LK, Miyazaki K, Ayabe A, Davis J, Tauchi-Nishi P and Shimizu D. The LAST guidelines in clinical practice: implementing recommendations for p16 use. *Am J Clin Pathol* 2015, 144(6):844-9
36. Zhang G, Yang B and Abdul-Karim FW. p16 Immunohistochemistry is useful in confirming highgrade squamous intraepithelial lesions (HSIL) in women with negative HPV testing. *Int J Gynecol Pathol* 2015, 34(2):180-6
37. Darragh TM. The LAST Project and the diagnostic bottom line. *Cytopathology* 2015, 26(6):343-5
38. de Sam Lazaro S, Newbill CP, Berlin M and Morgan TK. p16 Staining of Cervical Biopsies May Decrease the Frequency of Unnecessary Loop Electrosurgical Excision Procedures. *J Low Genit Tract Dis* 2016, 20(3):201-6
39. Stoler MH, Wright TC Jr, Ferenczy A, Ranger-Moore J, Fang Q, Kapadia M, Ridder R. Routine Use of Adjunctive p16 Immunohistochemistry Improves Diagnostic Agreement of Cervical Biopsy Interpretation: Results From the CERTAIN Study. *Am J Surg Pathol*. 2018, 42(8):1001-9

## Annex 2

### Example for diffuse staining pattern




Fig. 1: CIN 3

### Example for focal staining pattern




Fig. 2: Squamous metaplasia, mature