

cobas u pack

| REF         | CONTENT | SYSTEM                     |
|-------------|---------|----------------------------|
| 06334601001 | ▽ 400   | cobas u 601 urine analyzer |

## English

## Caution

Do not open the inner bag prior to use.  
Immediately insert cassette into analyzer!

## Intended use

The **cobas u** pack is a cassette with teststrips for the in vitro qualitative or semi-quantitative determination of pH, leukocytes, nitrite, protein, glucose, ketones, urobilinogen, bilirubin, color and erythrocytes in urine with the **cobas u** 601 urine analyzer. These measurements are useful in the evaluation of renal, urinary, hepatic and metabolic disorders.

For professional use only.

## Test principle

**pH:** The test paper contains the indicators methyl red, o-cresolphthalein and bromothymol blue and reacts specifically with H<sup>+</sup>-ions.

**Leukocytes (LEU):** The test reveals the presence of granulocyte esterases. These esterases cleave an indoxyl ester, and the indoxyl so liberated reacts with a diazonium salt to produce a violet dye.

**Nitrite (NIT):** The test is based on the principle of the Griess test and is specific for nitrite. The reaction reveals the presence of nitrite and hence indirectly nitrite-forming bacteria in the urine by a pink-to-red coloration of the test patch. Even a slight pink coloration is indicative of significant bacteriuria.

**Protein (PRO):** The test is based on the principle of the protein error of a pH indicator. It is particularly sensitive to albumin.

**Glucose (GLU):** The glucose determination is based on the specific glucose-oxidase/peroxidase reaction (GOD/POD method).

**Ketones (KET):** This test is based on the principle of Legal's test and is more sensitive to acetoacetic acid than to acetone.

**Urobilinogen (UBG):** A stable diazonium salt reacts almost immediately with urobilinogen to give a red azo dye.

**Bilirubin (BIL):** The test is based on the coupling of bilirubin with a diazonium salt. Even the slightest pink coloration constitutes a positive, i.e. pathologic, result. Other urinary constituents produce a more or less intense yellow coloration.

**Blood (ERY/Hb):** The peroxidase-like action of hemoglobin and myoglobin specifically catalyzes the oxidation of the indicator by means of the organic hydroperoxide contained in the test paper to give a blue-green coloration.

**Compensation area (COMP):** This white area, which is not impregnated with reagents, allows instrumental compensation for the intrinsic color of the urine while testing leukocytes, nitrite, glucose, ketones, urobilinogen, bilirubin, erythrocytes; and determination of the urine color (COL).

## Reagents

Each test contains per 1 cm<sup>2</sup> test patch area the following:

**pH:** Bromothymol blue 13.9 µg; methyl red 1.1 µg; o-cresolphthalein 7.3 µg

**Leukocytes:** Indoxylcarbonic acid ester 15.5 µg; methoxymorpholinobenzene diazonium salt 5.5 µg

**Nitrite:** 3-hydroxy-1,2,3,4-tetrahydro-7,8-benzoquinoline 33.5 µg; sulfanilamide 29.1 µg

**Protein:** 3',3'',5',5''-tetrachlorophenol-3,4,5,6-tetrabromosulfophthalein 13.9 µg

**Glucose:** 3,3',5,5'-tetramethylbenzidine 103.5 µg; GOD 6 U, POD 35 U

**Ketones:** Sodium nitroprusside 96.5 µg

**Urobilinogen:** 4-methoxybenzene-diazonium-tetrafluoroborate 67.7 µg

**Bilirubin:** 2,6-dichlorobenzene-diazonium-tetrafluoroborate 16.7 µg

**Blood:** 3,3',5,5'-tetramethylbenzidine 52.8 µg; 2,5-dimethyl-2,5-dihydroperoxyhexane 297.2 µg

## Precautions and warnings

For in vitro diagnostic use for laboratory professionals. Exercise the normal precautions required for handling all laboratory reagents.

Infectious or microbial waste:

Warning: handle waste as potentially biohazardous material. Dispose of waste according to accepted laboratory instructions and procedures.

Environmental hazards:

Apply all relevant local disposal regulations to determine the safe disposal.

Safety data sheet available for professional user on request.

**Note:** If the cassette has been stored refrigerated, it must be left at room temperature for a minimum of one hour prior to use.

The cassette contains a non-toxic silicate-based desiccant which must not be removed. If ingested by accident, drink large quantities of water.

The stopper of the test strip vial contains a non-toxic silicate-based desiccant, which must not be removed. If ingested by accident, drink large quantities of water.

## Reagent handling

Ready for use.

## Storage and stability

Store the cassette at 2-30 °C.

After loading the cassette into the analyzer, the test strips are stable within the tightly closed cassette compartment for 14 days. After this period, the cassette has to be replaced by a new one.

Do not use the cassette after the specified expiry date.

## Specimen collection and preparation

Use only clean, well-rinsed vessels to collect urine.

Do not add preservatives to the urine.

Use fresh urine that has not been centrifuged.<sup>1</sup> The urine specimen should not stand for more than 2 hours before testing.<sup>1</sup> For specimen collection and preparation only use suitable tubes or collection containers, as false positive readings, particularly for glucose and protein, can result from residues of detergent or strongly oxidizing disinfectants in the specimen collection vessel.<sup>2</sup>

Using midstream urine is recommended to avoid contamination by commensal urethral flora in both sexes.<sup>2</sup> Do not expose urine specimens to sunlight as this induces oxidation of bilirubin and urobilinogen and hence leads to artificially low results for these 2 parameters.<sup>2</sup> Vaginal secretion or menstrual blood may contaminate urine from females.<sup>2</sup>

Diagnosis or therapy should never be based on one test result alone but should be established in the context of all other medical findings. In doubtful cases, it is therefore advisable to repeat the test and consider potential interferences.

In case of a positive result it is advisable to use a follow-up investigation.

## Materials provided

- [REF] 06334601001, Cassette with 400 test strips

## Materials required (but not provided)

- [REF] 06390498001, **cobas u** 601 urine analyzer
- [REF] 06390579001, **cobas u** calibration strip
- Controls as indicated below
- General laboratory equipment


## Assay

1. Unfold and cut open aluminum bag with scissors (illustration 1).
2. Remove test strip cassette from the packaging and remove the two protection pads (illustration 2).
3. Immediately place test strip cassette into the **cobas u** 601 urine analyzer (illustration 3).

Follow the instructions in the Operators Manual of the instrument for correct insertion and positioning. These instructions also contain information on further handling precautions of the cassette.

# cobas u pack

cobas u pack



Note: If the cassette is stored in the opened bag or exposed to air (humidity, nitrogen oxides) for more than 3 minutes, environmental conditions may cause a color change of the test patches and damage of the reagents.

This has to be avoided. Do not use the cassette if the packaging shows severe damages, or the test strip layers in the cassette are not correctly aligned, or the test strips show unusual coloring.

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate Operators Manual for analyzer-specific assay instructions.

## Calibration

**cobas u** calibration strips are used for the calibration of the photometer unit of the **cobas u 601** urine analyzer. For details see the Operators Manual of the analyzer.

## Quality control

For quality control, use commercially available urine controls, or other suitable control material.

Following quality controls from Bio-Rad are recommended to use:

- Bio-Rad qUAntify Plus Control
- Bio-Rad Liquichek Urinalysis Control

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

## Calculation

After the test strip has been accepted by the instrument, it is measured by means of reflectance photometry. The results are automatically calculated and printed on the report form in terms of "normal", "neg.", "pos." or as concentration values.

## Limitations - interference

Therapeutic drugs and endogenous substances were tested for a potential interference to the test parameters of the **cobas u** pack.

All parameters were tested with negative urine samples and samples spiked to the first positive concentration range.

Therapeutic drugs were tested at concentrations in urine occurring under medication with the therapeutic dosage and above.

There are no significant therapeutic drug interferences up to the concentrations as presented below:

| Parameter | Therapeutic drug | No interference up to | Effect above stated concentration |
|-----------|------------------|-----------------------|-----------------------------------|
| NIT       | Ascorbic acid    | 1500 mg/L             | false negative results            |
|           | Phenazopyridine  | 150 mg/L              | false positive results            |
| PRO       | Salicylic acid   | 4800 mg/L             | false negative results            |
| GLU       | Ascorbic acid    | 250 mg/L              | false normal results              |

| Parameter | Therapeutic drug | No interference up to | Effect above stated concentration                    |
|-----------|------------------|-----------------------|--|
| KET       | N-Acetylcysteine | 30 mg/L               | false positive results and elevated positive results |
|           | Levodopa         | 125 mg/L              | false positive results and elevated positive results |
|           | Methylodopa      | 100 mg/L              | false positive results and elevated positive results |
| UBG       | Gabapentin       | 2400 mg/L             | false normal results                                 |
|           | Phenazopyridine  | 150 mg/L              | false positive results and elevated positive results |
| BIL       | Amoxicillin      | 8000 mg/L             | false negative results                               |
|           | Ascorbic acid    | 600 mg/L              | false negative results                               |
| ERY       | Ascorbic acid    | 1000 mg/L             | false negative results                               |
|           | Gabapentin       | 7200 mg/L             | false negative results                               |
|           | Ibuprofen        | 750 mg/L              | false negative results                               |
|           | Levodopa         | 375 mg/L              | false positive results and elevated positive results |
|           | Methylodopa      | 800 mg/L              | false positive results and elevated positive results |
|           | Salicylic acid   | 2400 mg/L             | false negative results                               |

There are no significant endogenous substance interferences up to the concentrations as presented below:

| Parameter | Endogenous substance | No interference up to | Effect above stated concentration  |
|-----------|----------------------|-----------------------|--|
| LEU       | Bilirubin            | 200 mg/L              | false positive results and elevated positive results                         |
|           | Glucose              | 10000 mg/L            | false negative results   |
|           | Hemoglobin           | 200 mg/L              | false positive results and elevated positive results                         |
|           | Urobilinogen         | 150 mg/L              | false negative results, false positive results and elevated positive results |

| Parameter | Endogenous substance | No interference up to | Effect above stated concentration                    |
|-----------|----------------------|-----------------------|--|
| NIT       | Bilirubin            | 200 mg/L              | false positive results                               |
|           | Creatinine           | 9000 mg/L             | false negative results                               |
|           | Hemoglobin           | 400 mg/L              | false positive results                               |
|           | Urobilinogen         | 90 mg/L               | false negative results and false positive results    |
| PRO       | Ammonium             | 5000 mg/L             | false negative results                               |
|           | Creatinine           | 6000 mg/L             | elevated positive results                            |
|           | Hemoglobin           | 70 mg/L               | false positive results and elevated positive results |
|           | Urea                 | 75000 mg/L            | false positive results and elevated positive results |
|           | Urobilinogen         | 750 mg/L              | false positive results                               |
| GLU       | Ammonium             | 12500 mg/L            | false normal results                                 |
|           | Bilirubin            | 400 mg/L              | false normal results                                 |
|           | Urea                 | 90000 mg/L            | false normal results                                 |
|           | Urobilinogen         | 120 mg/L              | false positive results                               |
| KET       | Hemoglobin           | 600 mg/L              | false positive results and elevated positive results |
|           | Urobilinogen         | 1250 mg/L             | false negative results                               |
| UBG       | Bilirubin            | 400 mg/L              | false normal results                                 |
|           | Nitrite              | 10 mg/L               | false normal results                                 |
| BIL       | Urobilinogen         | 45 mg/L               | false positive results and elevated positive results |
|           | Nitrite              | 10 mg/L               | false negative results                               |
| ERY       | Nitrite              | 40 mg/L               | false negative results                               |
|           | Uric acid            | 800 mg/L              | false negative results                               |
|           | Urobilinogen         | 600 mg/L              | false positive results and elevated positive results |

**Common limitations:**

**NIT:** Prolonged urinary retention in the bladder (4-8 hours) is essential in order to obtain an accurate result.<sup>2</sup> More than 80 % of all bacteria responsible for urinary tract infections are Gram-negative rods (E.coli, Klebsiella, Enterobacter and Proteus species).<sup>3</sup> Gram-negative bacteria have the ability to reduce urinary nitrate to nitrite and can therefore be detected indirectly with the test strips.<sup>2</sup> Normal nutrition as a rule ensures a sufficiently high content of nitrate in the urine for the detection of bacteria.<sup>4</sup> Some common uropathogens, e.g. Enterococcus spp. and Staphylococcus spp. (5-15 % of bacteria responsible for urinary tract infections),<sup>3</sup> do not reduce urinary nitrate to nitrite and will therefore not be detected whatever their urinary concentration.<sup>2</sup> False-negative results may occur as a result of strong diuresis with frequent voiding of urine, insufficient intake or too short retention of urine in the bladder.<sup>2</sup>

Attention: Nitrogen oxides present in the atmosphere may have an influence on the stability of the nitrite test parameter.

**PRO:** False positive readings may be found after infusion of polyvinylpyrrolidone (blood substitute).<sup>2</sup>

**UBG:** Drugs that turn red in an acid environment (e.g. phenazopyridine) may produce false positive readings or reddish colorations on the test parameter for urobilinogen.<sup>5</sup>

**BIL:** Drugs that turn red in an acid environment (e.g. phenazopyridine) may produce false positive readings or reddish colorations on the test parameter for bilirubin.<sup>5</sup>

**Blood/ERY:** The result values refer to intact erythrocytes. At concentrations of about 5-50 Ery/ $\mu$ L, significant hemolysis (such as may occur on prolonged standing of the urine) leads to values which are higher than the corresponding concentrations given for intact erythrocytes. In women the test for blood may be falsified from 3 days before to 3 days after menstruation. It is therefore advisable not to perform the test during this time. After physical activity, e.g. strenuous jogging, raised values for erythrocytes and protein may occur without being signs of disease.<sup>6</sup>

**Note:** A selection of relevant commercially available drugs or their metabolites were tested. For questionable results, repeat the test after discontinuing a particular drug.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

**Expected values**

Based on literature. Current medical guidelines are leading.

| Parameter | Expected value  | Additional information   |
|-----------|---|--|
| pH        | 5 - 9 <sup>7</sup>  |  |
| LEU       | < 10 Leu/ $\mu$ L <sup>2</sup>  | 10 - 100 Leu/ $\mu$ L borderline <sup>2</sup>  |
| NIT       | < 1 $\mu$ mol/L (< 0.005 mg/dL) <sup>8</sup>                                      | A positive result is indicative of urinary tract infection, but a negative result does not rule out UTI. <sup>5</sup>          |
| PRO       | $\leq$ 30 mg/dL <sup>9</sup>  | > 30 mg/dL proteinuria <sup>9</sup>  |
| GLU       | < 25 mg/dL < 1.4 mmol/L <sup>10</sup>   | For daytime urine. Using semi-quantitative reagent strips, expected values in a healthy population are negative. <sup>11</sup> |
| KET       | $\leq$ 2 mg acetoacetic acid/dL <sup>12</sup>                                     | Borderline > 2 mg up to 50 mg acetoacetic acid/dL <sup>12</sup>  |
| UBG       | < 1 mg/dL <sup>a),4</sup>   | 1 - 4 mg/dL borderline (4 mg/dL corresponding to 2+, indicating liver damage) <sup>4</sup>                                     |
| BIL       | neg. <sup>12</sup>  | When this method is used, normal urine contains no detectable bilirubin.   |
| ERY/Hb    | < 18 Ery/ $\mu$ L (< 3 Ery/HPF) <sup>12</sup>                                     | Hematuria $\geq$ 18 Ery/ $\mu$ L ( $\geq$ 3 Ery/HPF) <sup>13,14</sup>  |
|           | Conversion factor 5.8 to translate chamber counting HPF into $\mu$ L <sup>2</sup> |  |

a) Values displayed by the instrument are rounded compared to conventional values.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

**Result values**

| Parameter | Result values                                     |
|-----------|---|
| pH        | 5, 6, 6.5, 7, 8, 9                                |
| LEU       | NEG, 25, 100, 500 Leu/ $\mu$ L<br>NEG, 1+, 2+, 3+ |
| NIT       | NEG, POS  |

| Parameter | Result values   |
|-----------|---|
| PRO       | NEG, 25, 75, 150, 500 mg/dL<br>NEG, 0.25, 0.75, 1.5, 5.0 g/L<br>NEG, 1+, 2+, 3+, 4+ |
| GLU       | NORM, 50, 100, 300, 1000 mg/dL<br>NORM, 3, 6, 17, 56 mmol/L<br>NORM, 1+, 2+, 3+, 4+ |
| KET       | NEG, 5, 15, 50, 150 mg/dL<br>NEG, 0.5, 1.5, 5, 15 mmol/L<br>NEG, 1+, 2+, 3+, 4+     |
| UBG       | NORM, 1, 4, 8, 12 mg/dL<br>NORM, 17, 68, 135, 203 µmol/L<br>NORM, 1+, 2+, 3+, 4+    |
| BIL       | NEG, 1, 3, 6 mg/dL<br>NEG, 17, 50, 100 µmol/L<br>NEG, 1+, 2+, 3+                    |
| ERY       | NEG, 10, 25, 50, 150, 250 Ery/µL<br>NEG, 1+, 2+, 3+, 4+, 5+                         |
| COL       | pale yellow, yellow, amber, brown, orange, red, green, others                       |

#### Specific performance data

Representative performance data are given below. Results obtained in individual laboratories may differ.

The values specified for the analytical sensitivity are defined as the concentration of the analyte which leads to a positive result in  $\geq 90\%$  of the examined urines. For pH, analytical sensitivity is not applicable (N.A.). The method comparison data of **cobas u 601** urine analyzer with **cobas u 411** urine analyzer and **cobas u 601** urine analyzer with **URISYS 2400** Analyzer (for COL) using at least 1348 clinical samples are presented below. Method comparison results for specific gravity and clarity are presented in the Performance Data Document of the **cobas u 601** urine analyzer, addendum to the Operators Manual.

| Parameter | Analytical sensitivity | Method comparison <sup>b)</sup>                                      |
|-----------|------------------------|--|
| pH        | N. A.                  | Ident.: 77 %<br>pH 5+6: 98 %<br>pH 8+9: 88 %                         |
| LEU       | 10 - 30 Leu/µL         | NEG: 91 %<br>POS: 97 %   |
| NIT       | 0.03 - 0.07 mg/dL      | NEG: 95 %<br>POS: 94 %   |
| PRO       | 10 - 18 mg/dL albumin  | NEG: 96 %<br>POS: 93 %   |
| GLU       | 20 - 40 mg/dL          | NEG: 98 %<br>POS: 100 %  |
| KET       | 3 - 7 mg/dL            | NEG: 94 %<br>POS: 96 %   |
| UBG       | 1.0 - 1.6 mg/dL        | NEG: 96 %<br>POS: 98 %   |
| BIL       | 0.4 - 0.6 mg/dL        | NEG: 93 %<br>POS: 95 %   |
| ERY       | 3 - 15 Ery/µL          | NEG: 95 %<br>POS: 96 %   |
| COL       | N. A.                  | Pale yellow+yellow: 94 %<br>Amber: 73 %<br>Brown: 92 %<br>Red: 100 % |

b) The values for neg and pos indicate the proportion of concordant negative or positive results.

#### Precision

Precision experiments comprised an assessment of repeatability (within-run precision) and intermediate precision.

Repeatability was checked in 2 separate runs with 21 measurements each for the tested controls. In total there were 42 measurements performed per used control.

Intermediate precision was assessed over 21 days with 2 runs per day and duplicate measurements per used control. In total there were 84 measurements performed per used control.

The following results were obtained:

| Repeatability |                       |            |                 |
|---------------|-----------------------|------------|-----------------|
| Parameter     | Control <sup>c)</sup> | Result     | Exact agreement |
| pH            | Level 1               | 6.5        | 100 %           |
|               | Level 2               | 7          | 100 %           |
| LEU           | Level 1               | NEG        | 100 %           |
|               | Level 2               | 500 Leu/µL | 100 %           |
| NIT           | Level 1               | NEG        | 100 %           |
|               | Level 2               | POS        | 100 %           |
| PRO           | Level 1               | NEG        | 100 %           |
|               | Level 2               | 150 mg/dL  | 100 %           |
| GLU           | Level 1               | NORM       | 100 %           |
|               | Level 2               | 1000 mg/dL | 100 %           |
| KET           | Level 1               | NEG        | 100 %           |
|               | Level 2               | 150 mg/dL  | 100 %           |
| UBG           | Level 1               | NORM       | 100 %           |
|               | Level 2               | 12 mg/dL   | 100 %           |
| BIL           | Level 1               | NEG        | 100 %           |
|               | Level 2               | 6 mg/dL    | 100 %           |
| ERY           | Level 1               | NEG        | 100 %           |
|               | Level 2               | 250 Ery/µL | 100 %           |
| COL           | Level 1               | Yellow     | 100 %           |
|               | Level 2               | Brown      | 100 %           |

c) BIO-RAD Liquichek

| Intermediate precision |                       |            |                 |
|------------------------|-----------------------|------------|-----------------|
| Parameter              | Control <sup>c)</sup> | Result     | Exact agreement |
| pH                     | Level 1               | 6.5        | 100 %           |
|                        | Level 2               | 7          | 100 %           |
| LEU                    | Level 1               | NEG        | 100 %           |
|                        | Level 2               | 500 Leu/µL | 100 %           |
| NIT                    | Level 1               | NEG        | 100 %           |
|                        | Level 2               | POS        | 100 %           |
| PRO                    | Level 1               | NEG        | 100 %           |
|                        | Level 2               | 150 mg/dL  | 100 %           |
| GLU                    | Level 1               | NORM       | 100 %           |
|                        | Level 2               | 1000 mg/dL | 100 %           |
| KET                    | Level 1               | NEG        | 100 %           |
|                        | Level 2               | 150 mg/dL  | 100 %           |
| UBG                    | Level 1               | NORM       | 100 %           |
|                        | Level 2               | 12 mg/dL   | 100 %           |
| BIL                    | Level 1               | NEG        | 100 %           |
|                        | Level 2               | 6 mg/dL    | 100 %           |

# cobas u pack

cobas u pack



| Intermediate precision |                      |                  |                 |
|------------------------|----------------------|------------------|-----------------|
| Parameter              | Control <sup>o</sup> | Result           | Exact agreement |
| ERY                    | Level 1              | NEG              | 100 %           |
|                        | Level 2              | 250 Ery/ $\mu$ L | 100 %           |
| COL                    | Level 1              | Yellow           | 100 %           |
|                        | Level 2              | Brown            | 100 %           |

## References

- GP16-A3 (Urinalysis; Approved Guideline - Third Edition).
- ECLM, European Urinalysis Guidelines. Scand J Clin Lab Invest, 2000. 60: p. 1-96.
- Ronald A. The Etiology of Urinary Tract Infection: Traditional and Emerging Pathogens, Dis Mon (2003);49:71-82.
- King Strasinger S, Schaub Di Lorenzo M. Urinalysis and Body Fluids (2008) Fifth Edition, ISBN 978-0-8036-1697.
- Simerville JA., Maxted WC and Pahira J.J., Urinalysis: a comprehensive review. Am Fam Physician, 2005. 71(6): p. 1153-62.
- Rao PV, Jones JS. How to evaluate 'dipstick hematuria': What to do before you refer. Cleveland Clinic Journal of Medicine (2008);75(3):227-233.
- Thomas L. ed. Clinical Laboratory Diagnostics. Electronic ed. (2020);TH-Books GmbH.
- Pannala AS., et al., The effect of dietary nitrate on salivary, plasma, and urinary nitrate metabolism in humans. Free Radic Biol Med, 2003. 34(5): p. 576-84.
- Johnson DW. Global Proteinuria Guidelines: Are We Nearly There Yet?; Clin Biochem Rev. (2011);Vol.32.
- Cowart SL, Stachura ME. Glucosuria, Clinical Methods: The History, Physical, and Laboratory Examinations. 3rd edition., Boston 1990, Chapter 139.
- Wu, A.H.B., Tietz Clinical Guide to Laboratory Tests; 4th Edition. 2006.
- McPherson RA, M.R.P., HENRY'S Clinical Diagnosis and Management by Laboratory Methods, 23rd edition. ISBN 9780323295680, 2017.
- Wollin T, Laroche B and Psooy K, Canadian guidelines for the management of asymptomatic microscopic hematuria in adults. Can Urol Assoc J, 2009. 3(1): p. 77-80.
- Nielsen M, Qaseem A; High Value Care Task Force of the American College of Physicians. Hematuria as a marker of occult urinary tract cancer: advice for high-value care from the American College of Physicians. Ann Intern Med. 2016;164(7):488-497.








For further information, please refer to the appropriate Operators Manual for the analyzer concerned, and the Method Sheets of all necessary components.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

## Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard:

|   |   |
|---|---|
|  | Contents of kit                                     |
|  | Analyzers/Instruments on which reagents can be used |
|  | Reagent   |
|  | Calibrator  |
|  | Volume for reconstitution                           |
|  | Global Trade Item Number                            |
|  | Unique Device Identifier                            |

Rx only

For USA: Caution: Federal law restricts this device to sale by or on the order of a physician.

COBAS, NAVIFY and URISYS 2400 are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2025, Roche Diagnostics




Roche Diagnostics GmbH  
Sandhofer Strasse 116  
68305 Mannheim, Germany  
www.roche.com

+800 5505 6606

