

For general laboratory use.



LightCycler® Multiplex RNA Virus Master

 **Version: 11**

Content Version: August 2025

Easy-to-use reaction mix for one-step RT-qPCR using the LightCycler® PRO, LightCycler® 480, or the LightCycler® 96 Real-Time PCR Systems.

Cat. No. 06 754 155 001 1 kit
200 reactions of 20 µL final volume each

Cat. No. 07 083 173 001 1 kit
1,000 reactions of 20 µL final volume each

Store the kit at –15 to –25°C.

1.	General Information	3
1.1.	Contents	3
1.2.	Storage and Stability	3
	Storage Conditions (Product)	3
	Storage Conditions (Working Solution).....	3
1.3.	Additional Equipment and Reagent required	4
1.4.	Application	4
1.5.	Preparation Time.....	4
	Typical Run Time.....	4
2.	How to Use this Product	5
2.1.	Before you Begin	5
	Sample Materials	5
	Control Reactions	5
	Primers.....	5
	Probe	5
	Mg ²⁺ Concentration.....	5
	General Considerations	5
	Precautions	5
2.2.	Protocols	6
	LightCycler [®] PRO, LightCycler [®] 480, and LightCycler [®] 96 System protocols.....	6
	Protocol for use with the LightCycler [®] PRO System (Multiwell Plate 96 or 384)	6
	Protocol for use with the LightCycler [®] 480 Instrument II (Multiwell Plate 96 or 384).....	7
	Protocol for use with the LightCycler [®] 96 Instrument.....	10
	Setup of the RT-qPCR reaction for the LightCycler [®] PRO, LightCycler [®] 480, and LightCycler [®] 96 Instruments	10
3.	Results	12
4.	Troubleshooting	15
5.	Additional Information on this Product	16
5.1.	Test Principle	16
	How this product works	16
5.2.	Quality Control.....	16
6.	Supplementary Information	16
6.1.	Conventions.....	16
6.2.	Changes to previous version	16
6.3.	Ordering Information.....	17
6.4.	Trademarks.....	18
6.5.	License Disclaimer	18
6.6.	Regulatory Disclaimer.....	18
6.7.	Safety Data Sheet	18
6.8.	Contact and Support.....	18

1. General Information

1.1. Contents

Vial / bottle	Cap	Label	Function / description	Catalog number	Content
1	blue	LightCycler® Multiplex RNA Virus Master, RT Enzyme Solution, 200x conc.	Contains Reverse Transcriptase.	06 754 155 001	1 vial, 28 µL
				07 083 173 001	5 vials, 28 µL each
2	red	LightCycler® Multiplex RNA Virus Master, RT-qPCR Reaction Mix, 5x conc.	Contains RT-qPCR Reaction Buffer, AptaTaq Polymerase, dATP, dCTP, dGTP, dUTP, MgCl ₂ , and proprietary additives.	06 754 155 001	1 vial, 880 µL
				07 083 173 001	5 vials, 880 µL each
3	colorless	LightCycler® Multiplex RNA Virus Master, Water, PCR Grade	To adjust the final reaction volume.	06 754 155 001	3 vials, 1 mL each
				07 083 173 001	15 vials, 1 mL each

1.2. Storage and Stability

Storage Conditions (Product)

The kit is shipped on dry ice.

When stored at –15 to –25°C, the kit is stable through the expiration date printed on the label. The kit is stable at +2 to +8°C for 4 weeks.

Vial / bottle	Cap	Label	Storage
1	blue	RT Enzyme Solution, 200x conc.	Store at –15 to –25°C or store at +2 to +8°C for a maximum of 4 weeks. ⚠ Close lid immediately after use.
2	red	RT-qPCR Reaction Mix, 5x conc.	Store at –15 to –25°C. ⚠ Avoid repeated freezing and thawing (more than 5 times). Aliquot Vial 2 and freeze or store Vial 2 at +2 to +8°C for a maximum of 4 weeks.
3	colorless	Water, PCR Grade	Store at –15 to –25°C or store at +2 to +8°C for a maximum of 4 weeks.

Storage Conditions (Working Solution)

Prepare the reagents right before use. The working solution, see Section, **Protocols, Setup of the RT-qPCR reaction**, is stable at +15 to +25°C for up to 4 hours, and is therefore ideal for use in automated workflows.

1.3. Additional Equipment and Reagent required

Standard laboratory equipment

- Nuclease-free pipette tips
- 1.5 mL RNase-free microcentrifuge tubes to prepare master mixes and dilutions
- To minimize risk of RNase contamination, autoclave all vessels

 *Wear gloves at all times.*

For RT-qPCR

- Real-Time PCR systems such as the LightCycler® PRO, LightCycler® 480, LightCycler® 96, or other Real-Time PCR Systems
- LightCycler® Multiwell Plate 96, white*
- LightCycler® Multiwell Plate 384, white*
- LightCycler® Multiwell Plate 96, white, 4 bar codes*
- LightCycler® Multiwell Plate 384, white, 4 bar codes*
- Sealing Foil Applicator*
- LightCycler® Sealing Foil*
- LightCycler® 8-Tube Strips*
- LightCycler® 8-Tube Strip Adapter Plate*
- Centrifuge with swinging-bucket rotor

1.4. Application

The LightCycler® Multiplex RNA Virus Master is designed for fast, highly sensitive and specific real-time one-step RT-qPCR analysis of viral RNA.

The 2-vial composition (separate RT enzyme and qPCR Mix) is ideally suited for use with RT-minus controls. The proprietary reaction buffer allows a fast and convenient hot start RT-qPCR without pre-activation of the Taq DNA Polymerase. The kit is optimized for hydrolysis probes and does not require optimization with MgCl₂.

1.5. Preparation Time

Typical Run Time

The LightCycler® Multiplex RNA Virus Master can be used for multiplex RT-qPCR protocols. For example, a septaplex protocol using 45 cycles requires less than 75 minutes when using the LightCycler® PRO System.

2. How to Use this Product

2.1. Before you Begin

Sample Materials

Use any viral template RNA suitable for RT-qPCR in terms of purity, concentration, and absence of RT-PCR inhibitors. For reproducible isolation of nucleic acids, use:

- Either a MagNA Pure System together with a dedicated reagent kit (for automated isolation),
- or a High Pure Nucleic Acid Isolation Kit (for manual isolation).

Control Reactions

Always run appropriate positive, negative, and RT-minus controls with the samples.

Negative control and RT-minus control reactions

Always run a negative control with the samples. To prepare a negative control:

- To check for the presence of contamination, prepare and include a negative control by replacing the template RNA with Water, PCR Grade (Vial 3).
- For the RT-minus control, omit the RT enzyme in the mix. This will verify whether your signal comes from the RNA target or from DNA contamination.

Primers

Suitable concentrations of PCR primers range from 0.2 to 0.5 μM (final concentration in RT-qPCR). The recommended starting concentration is 0.5 μM each.

Probe

Suitable concentrations of hydrolysis probes range from 0.05 to 0.5 μM (final concentration in PCR). The recommended starting concentration is 0.25 μM each.

- i** *The optimal probe concentration is the lowest concentration that results in the lowest quantification cycle value (C_q) and adequate fluorescence dynamics for a given target concentration.*
- i** *For a hydrolysis probe hybridization complex, the T_m of the hydrolysis probe has to be higher than the T_m of the primers.*

Mg²⁺ Concentration

The master mix of this kit is optimized with a fixed concentration of MgCl₂, which works with nearly all primer combinations. There is no need for adjustment.

General Considerations

Precautions

Always use RNase-free techniques. RNase contaminated reagents and reaction vessels will degrade template RNA. Please follow these guidelines to minimize risk of contamination:

- Wear disposable gloves and change them frequently.
- Avoid touching surfaces or materials that could cause RNase carryover.
- Use only reagents provided in this kit. Substitutions may introduce RNases.
- Clean and decontaminate work areas and instruments, including pipettes, with commercially available decontamination reagents.
- Use only new RNase-free aerosol-blocking pipette tips and microcentrifuge siliconized reaction tubes.
- Use a work area specifically designated for RNA work, and if possible, use reaction vessels and pipettes dedicated only for work with template RNA.

2.2. Protocols

LightCycler® PRO, LightCycler® 480, and LightCycler® 96 System protocols

The following procedure is optimized for use with the corresponding LightCycler® System you are using.

⚠ Program the LightCycler® Instrument before preparing the reaction mixes.

i For details on how to program the experimental protocol, see the *LightCycler® PRO System User Assistance*, *LightCycler® 480 Instrument II Operator's Manual*, or *LightCycler® 96 System Operator's Guide*.

Protocol for use with the LightCycler® PRO System (Multiwell Plate 96 or 384)

The LightCycler® PRO System protocol contains the following programs:

- **Reverse transcription** of the viral template RNA
- **Preincubation** for denaturation of the cDNA/RNA hybrid
- **Amplification** of the cDNA
- **Cooling** of the thermal cycler

Setup				
Thermal cycler type	Reaction volume [μL]			
96 (384)	10 – 100 (5 – 20)			
Detection format				
Select dyes used in your assays. If the dye is not predefined, use appropriate filter combination in the user-defined detection form				
Programs ⁽¹⁾				
Program name	Cycles			
Reverse transcription	1			
Preincubation	1			
Amplification	45 ⁽²⁾			
Cooling	1			
Temperature targets				
	Target [°C]	Acquisition mode	Duration [s]	Ramp rate [°C/s]
Reverse transcription	50 ⁽³⁾	None	600 ⁽⁴⁾	4.4 (4.8)
Preincubation	95	None	30	4.4 (4.8)
2-step amplification	95	None	5	4.4 (4.8)
	60 ⁽⁵⁾	Single	30	2.2 (2.5)
Cooling (automatically added)	40	None	30	2.2 (2.5)

⁽¹⁾ The LightCycler® Multiplex RNA Virus Master includes an aptamer for hot start. Anneal/Extend temperatures significantly below 60°C may not be well tolerated because aptamer will inhibit DNA polymerase activity at low temperatures.

⁽²⁾ Forty-five cycles is suitable for most assays. If the assay shows steep amplification curves and early crossing points, 40 cycles should be sufficient. Reducing the number of cycles will improve uniformity of the product and reduce the time required for the assay.

⁽³⁾ We recommend running the reverse transcription at +50°C. If necessary, it is also possible to run it at +55°C.

⁽⁴⁾ Ten minutes is the recommended default protocol. Depending on your assay, it should be possible to reduce the time for reverse transcription down to 4 minutes. For best results, run a test with 10, 8, 6, 4 minutes, etc.

⁽⁵⁾ Most available assays are designed for an annealing temperature of +60°C. If the T_m of the primers and probe indicate using different settings, select an annealing temperature +5°C below the calculated primer T_m for initial experiments. For assay optimization, choose the highest annealing temperature which still results in a high yield of amplicon with a low C_q and adequate fluorescence dynamics.

Color Compensation protocol for the LightCycler® PRO Instrument

The LightCycler® PRO Instrument does not require the creation of a color compensation object.

Protocol for use with the LightCycler® 480 Instrument II (Multiwell Plate 96 or 384)

The LightCycler® 480 Instrument II protocol contains the following programs:

- **Reverse transcription** of the viral template RNA
- **Preincubation** for denaturation of the cDNA/RNA hybrid
- **Amplification** of the cDNA
- **Cooling** of the thermal cycler

Setup					
Thermal block type			Reaction volume [µL]		
96 (384)			20 (10)		
Detection format		Excitation filter		Emission filter	
For example: Mono Color Hydrolysis Probe / UPL Probe					
FAM		465		510	
Programs ⁽¹⁾					
Program name		Cycles		Analysis mode	
Reverse transcription		1		None	
Preincubation		1		None	
Amplification		45 ⁽²⁾		Quantification	
Cooling		1		None	
Temperature targets					
	Target [°C]	Acquisition mode	Hold [hh:mm:ss]	Ramp rate [°C/s]	Acquisitions [per °C]
Reverse transcription	50 ⁽³⁾	None	00:10:00 ⁽⁴⁾	4.4 (4.8)	-
Preincubation	95	None	00:00:30	4.4 (4.8)	-
Amplification	95	None	00:00:05	4.4 (4.8)	-
	60 ⁽⁵⁾	Single	00:00:30	2.2 (2.5)	-
Cooling	40	None	00:00:30	2.2 (2.5)	-

⁽¹⁾ The LightCycler® Multiplex RNA Virus Master includes an aptamer for hot start. Anneal/Extend temperatures significantly below 60°C may not be well tolerated because aptamer will inhibit DNA polymerase activity at low temperatures.

⁽²⁾ Forty-five cycles is suitable for most assays. If the assay shows steep amplification curves and early crossing points, 40 cycles should be sufficient. Reducing the number of cycles will improve uniformity of the product and reduce the time required for the assay.

⁽³⁾ We recommend running the reverse transcription at +50°C. If necessary, it is also possible to run it at +55°C.

⁽⁴⁾ Ten minutes is the recommended default protocol. Depending on your assay, it should be possible to reduce the time for reverse transcription down to 4 minutes. For best results, run a test with 10, 8, 6, 4 minutes, etc.

⁽⁵⁾ Most available assays are designed for an annealing temperature of +60°C. If the T_m of the primers and probe indicate using different settings, select an annealing temperature +5°C below the calculated primer T_m for initial experiments. For assay optimization, choose the highest annealing temperature which still results in a high yield of amplicon with a low C_q and adequate fluorescence dynamics.

Color Compensation protocol for the LightCycler® 480 Instrument II

For a multicolor measurement, the application of a color compensation file is necessary to compensate for optical crosstalk between the different channels. For the LightCycler® 480 Instrument II, an instrument-specific color compensation file is mandatory for multicolor experiments, and a color compensation object can be generated as shown below.

2. How to Use this Product

The LightCycler® 480 Instrument II protocol contains the following programs:

- **Reverse transcription** of the viral template RNA
- **Preincubation** for denaturation of the cDNA/RNA hybrid
- **Amplification** of the cDNA
- **Temperature gradient step** to create the Color Compensation file
- **Cooling** of the thermal cyclers

i For details on how to program the experimental protocol, see the *LightCycler® 480 Instrument II Operator's Manual, Version 1.5*.

Programming a customized detection format for the LightCycler® 480 System filter combination selection

The Detection Format in the LightCycler® 480 Software, Version 1.5 setup must be customized for the applied multicolor hydrolysis probe format used in the RT-qPCR detection. In the **Tool** module, the **Detection Formats** option allows creating a new detection format specified by the user, including a **Detection Format** list, a **Filter Combination** selection area, and a **Selected Filter Combination List**. Different filter settings for the LightCycler® 480 Instrument II are defined.

The following table shows the RT-qPCR parameters that must be programmed for a LightCycler® 480 Instrument II Color Compensation file with a LightCycler® Multiwell Plate 96.

Setup					
Thermal cycler type			Reaction volume [µL]		
96			20		
Detection format		Excitation filter		Emission filter	
For example: 3 Color Hydrolysis Probe					
FAM		465		510	
LC Red 610		533		610	
Cy5		618		660	
For new customized hydrolysis probes detection formats, set for all selected filters in the "Selected Filter Combination List" (under Tools), the following values:					
Melt factor		1			
Quant factor		10			
Max integration time (sec)		2			
Programs					
Program Name		Cycles		Analysis mode	
Reverse transcription		1		None	
Preincubation		1		None	
Amplification		45		Quantification	
Temperature gradient step		1		Color compensation	
Cooling		1		None	
Temperature targets					
	Target [°C]	Acquisition mode	Hold [hh:mm:ss]	Ramp rate [°C/s]	Acquisitions [per °C]
Reverse transcription	50	None	00:10:00	4.4	-
Preincubation	95	None	00:00:30	4.4	-
Amplification	95	None	00:00:05	4.4	-
	60	Single	00:00:30	2.2	-
Temperature gradient step	95	None	00:00:10	4.4	-
	40	None	00:00:10	2.2	-
	95	Continuous	-	-	1 - 5
Cooling	40	None	00:00:30	2.2	-

Preparation of the Color Compensation run

Prepare the calibrator RT-qPCR mix for more than one reaction. Multiply the amount in the volume column by the number of reactions you need (minimum of 3 to 5 replicates) plus additional reactions since there will be a slight loss of liquid during the pipetting steps. In order to ensure accuracy, do not pipette volumes less than 1 μL . For each dye, set up the following reactions:

Reagent	1x buffer [μL]	1x for each dye [μL]
RT Enzyme Solution, 200x conc. (Vial 1)	0.1	0.1
RT-qPCR Reaction Mix, 5x conc. (Vial 2)	4.0	4.0
Detection mix for each dye	–	X (depending on the assay)
Water, PCR Grade (Vial 3)	15.9	Y (depending on the assay)
Template, such as viral RNA or positive samples eluates	–	5.0
Total volume	20.0	20.0

- 1 Pipette the three replicates of each different calibrator mix into a LightCycler® Multiwell Plate 96.

- 2 Seal the LightCycler® Multiwell Plate with a LightCycler® Sealing Foil.

- 3 Place the multiwell plate into the centrifuge with a swinging-bucket rotor and suitable adapter, and balance it with a suitable counterweight, such as another multiwell plate.
 - Centrifuge at $1,500 \times g$ for 2 minutes.

- 4 Load the multiwell plate into the LightCycler® 480 Instrument II and start the program.

Create Color Compensation object

When the experiment is finished, open the **Sample Editor**. Select workflow **Color Comp** and designate the used positions of the Multiwell Plate as “Water” for Buffer replicates, and the appropriate dyes respectively (e.g., FAM, LC Red 610, Cy5 for the example mentioned above).

Open the **Analysis** module Color Compensation, click on **Calculate**, and save the color compensation object in the database.

From this point on, you can apply this “CC Object” to each multicolor experiment performed with FAM, LC Red 610, and Cy5 on the same instrument.

Protocol for use with the LightCycler® 96 Instrument

Run editor				
Detection format	Reaction volume [µL]			
For example: Dyes 1: FAM	20			
Programs ⁽¹⁾				
	Temp. [°C]	Ramp [°C/s]	Duration [s]	Acquisition mode
Reverse transcription	50 ⁽²⁾	4.4	600 ⁽³⁾	None
Preincubation	95	4.4	30	None
2-step amplification	No. of cycles: 45 ⁽⁴⁾			
	95	4.4	5	None
	60 ⁽⁵⁾	2.2	30	Single

⁽¹⁾ The LightCycler® Multiplex RNA Virus Master includes an aptamer for hot start. Anneal/Extend temperatures significantly below 60°C may not be well tolerated because aptamer will inhibit DNA polymerase activity at low temperatures.

⁽²⁾ We recommend running the reverse transcription at +50°C. If necessary, it is also possible to run it at +55°C.

⁽³⁾ Ten minutes is the recommended default protocol. Depending on your assay, it should be possible to reduce the time for reverse transcription down to 4 minutes. For best results, run a test with 10, 8, 6, 4 minutes, etc.

⁽⁴⁾ Forty-five cycles is suitable for most assays. If the assay shows steep amplification curves and early crossing points, 40 cycles should be sufficient. Reducing the number of cycles will improve uniformity of the product and reduce the time required for the assay.

⁽⁵⁾ Most available assays are designed for an annealing temperature of +60°C. If the T_m of the primers and probe indicate using different settings, select an annealing temperature +5°C below the calculated primer T_m for initial experiments. For assay optimization, choose the highest annealing temperature which still results in a high yield of amplicon with a low C_q and adequate fluorescence dynamics.

Color Compensation protocol for the LightCycler® 96 Instrument

The LightCycler® 96 Instrument does not require the creation of a color compensation object.

Setup of the RT-qPCR reaction for the LightCycler® PRO, LightCycler® 480, and LightCycler® 96 Instruments

Follow the procedure below to prepare at least ten 20 µL standard reactions:

i Always wear gloves during handling.

1 Thaw the solutions and, to ensure recovery of all the contents, briefly spin vials in a microcentrifuge before opening.
– Mix carefully by pipetting up and down or vortex briefly.

2 Place samples on ice.

3 Prepare a 20x conc. solution of your primers and a 20x conc. solution of your probes.

- 4 In a 1.5 mL reaction tube, prepare the RT-qPCR Mix.
- For best results, prepare at least 10 reactions in order to reduce pipetting errors. To prepare more reactions, multiply the amount in the “Volume 1 reaction” column below by the number of reactions to be run, plus at least one additional reaction.

Reagent ⁽¹⁾	Volume 1 reaction [μL]		Volume 10 reactions [μL]		Final conc.
	96-well plate	384-well plate	96-well plate	384-well plate	
Water, PCR Grade (Vial 3)	8.9	4.45	89.0	44.5	–
RT-qPCR Reaction Mix, 5x conc. (Vial 2)	4.0	2.00	40.0	20.0	1x
Primer Mix, 20x conc. ⁽²⁾	1.0	0.50	10.0	5.0	1x
Probe Mix, 20x conc.	1.0	0.50	10.0	5.0	1x
RT Enzyme Solution, 200x conc. (Vial 1)	0.1	0.05	1.0	0.5	1x
Total volume	15.0	7.50	150.0	75.0	

⁽¹⁾ For eluates derived from stool samples, it is recommended to add 0.2 μg/μL (final) of molecular biology-grade Bovine Serum Albumin*. Adjust for BSA volume by subtracting from PCR Grade Water volume.

⁽²⁾ Due to possible primer/primer interactions that occur during storage, it may be necessary to preheat the PCR primer mix for 1 minute at +95°C before pipetting the full mixture. This extra step will ensure optimum sensitivity.

- 5 Mix carefully by pipetting up and down or vortex briefly.
- 6 Prepare reagents right before use; the working solution (everything combined except RNA template) is stable at +15 to +25°C up to 4 hours, and is therefore ideal for use in automated workflows.
- 7 Prepare sample concentration of the viral RNA and/or DNA.
- 8 Pipette 15 μL (7.5 μL) RT-qPCR Mix into a multiwell plate or LightCycler® 8-Tube Strip.
- Add 5 μL (2.5 μL) of the RNA and/or DNA template.
 - Seal multiwell plate with a LightCycler® Sealing Foil or seal the LightCycler® 8-Tube Strips using the corresponding lid.
- 9 Place the multiwell plate into the centrifuge with a swinging-bucket rotor and suitable adapter, and balance it with a suitable counterweight, such as another multiwell plate.
- Centrifuge at 1,500 × *g* for 2 minutes.
- 10 Load the reaction vessels into the LightCycler® PRO, LightCycler® 480, or LightCycler® 96 Instrument.
- 11 Start the PCR program described above.
- If you use reaction volumes other than 20 μL, be sure to adapt the right volume in the running protocol. To start, we recommend using the same hold times as for the 20 μL volume.

3. Results

Results on the LightCycler® PRO Instrument

The following results were obtained using the LightCycler® Multiplex RNA Virus Master on the LightCycler® PRO Instrument. A duplex reaction using primers and Universal ProbeLibrary probes specific for **G6PD** (FAM) and **β2M** (HEX) was run. Human total RNA served as the template. Total RNA from different tissues in a dilution series (50 ng, 5 ng, 500 pg, 50 pg, and 10 pg and a no template control [NTC]) was analyzed using real-time PCR.

All targets:




Fig. 1: Duplex RT-qPCR amplification reaction on the LightCycler® PRO Instrument within one picture.

The following two images show the individual amplification curves of the respective channels.

FAM Channel (494 - 523)




Fig. 2: The FAM channel shows the results for **G6PD**. Amplification curves shown were obtained from dilutions of 50 ng (far left), 5 ng, 500 pg, 50 pg, and 10 pg (far right) human total RNA per well, including a no template control (flat line). Duplex RT-qPCR was performed in a reaction volume of 20 µL per well in a LightCycler® Multiwell Plate 96.


HEX Channel (541 - 565)

Fig. 3: The HEX channel shows the results for $\beta 2M$. Amplification curves shown were obtained from dilutions of 50 ng (far left), 5 ng, 500 pg, 50 pg, and 10 pg (far right) human total RNA per well, including a no template control (flat line). Duplex RT-qPCR was performed in a reaction volume of 20 μ L per well in a LightCycler[®] Multiwell Plate 96.

3. Results

Results on the LightCycler® 480 Instrument II

The following results were obtained using the LightCycler® Multiplex RNA Virus Master on the LightCycler® 480 Instrument II. A duplex reaction using primers and UPL probes specific for **G6PD** (FAM) and **β2M** (Yellow 555) was run. Human total RNA served as the template. Total RNA from different tissues in a dilution series (50 ng, 5 ng, 500 pg, 50 pg, and 10 pg) and a no template control [NTC] was analyzed using real-time PCR.

FAM Channel (465 - 510)




Fig. 4: The FAM channel shows the results for **G6PD**. Amplification curves shown were obtained from dilutions of 50 ng (far left), 5 ng, 500 pg, 50 pg, and 10 pg (far right) human total RNA per well, including a no template control (flat line). Duplex RT-qPCR was performed in a reaction volume of 20 μ L per well in a LightCycler® Multiwell Plate 96.

Yellow 555 Channel (533 - 580)




Fig. 5: The Yellow 555 channel shows the results for **β2M**. Amplification curves shown were obtained from dilutions of 50 ng (far left), 5 ng, 500 pg, 50 pg, and 10 pg (far right) human total RNA per well, including a no template control (flat line). Duplex RT-qPCR was performed in a reaction volume of 20 μ L per well in a LightCycler® Multiwell Plate 96.

4. Troubleshooting

Observation	Possible cause	Recommendation
Increase specificity.		Some assays show higher specificity when using a higher reverse transcription temperature and/or higher annealing temperature.
Fluorescence intensity varies.	Some of the reagent is still in the upper part of the microwell/8-tube strip, or an air bubble is trapped in the microwell/8-tube strip.	Repeat centrifugation, but allow sufficient centrifugation time so all reagent is at the bottom of the microwell/8-tube strip and air bubbles are expelled.
	Skin oils or dirt are present on the lid or sealing foil.	Always wear gloves when handling the multiwell plate/8-tube strip.
Fluorescence intensity is very low.	Low concentration or deterioration of dyes in the reaction mixtures because dye was not stored properly.	Keep dye-labeled reagents away from light. Store the reagents at -15 to -25°C and avoid repeated freezing and thawing.
	Poor PCR efficiency (reaction conditions not optimized).	Primer concentration should be in the range of 0.2 to $0.5\ \mu\text{M}$; probe concentration should be in the range of 0.2 to $0.5\ \mu\text{M}$ and half of the primer concentration.
		Check annealing temperature of primers and probes.
		Check experimental protocol.
		Optimize annealing temperature in the reverse transcription step or in the PCR reaction.
	Always run a positive control along with your samples.	
	Chosen integration time is too low.	Choose the appropriate detection format in combination with “dynamic” detection mode, or Increase integration time when using “manual” detection mode. i For details, see the <i>LightCycler® 480 Instrument II Operator’s Manual</i> .
RT-qPCR primers and probes are not optimized.	Check sequence and location of the hydrolysis probe on the PCR product. Check RT-qPCR product on an agarose gel.	
PCR has not been optimized.	Check primer design (quality). Check RT-qPCR product on an agarose gel.	
RNA is degraded during isolation or improper storage.	If possible, check RNA quality on a gel. Check RNA with an established RT-qPCR primer when available.	
Pipetting errors and omitted reagents.	Check for missing reagents.	
	Check the pipetting procedure.	
Impure sample material inhibits reaction.	Dilute sample 1:10 and repeat the analysis.	
	Repurify the nucleic acids to ensure removal of inhibitory agents.	
Negative control sample gives a positive signal.	Contamination	Remake all critical reaction mixes. Be sure to use special RT-qPCR setup working areas.

5. Additional Information on this Product

5.1. Test Principle

How this product works

The LightCycler® Multiplex RNA Virus Master consists of 3 different vials:

- Vial 1: RT Enzyme Solution (Reverse Transcriptase)
- Vial 2: RT-qPCR Reaction Mix
- Vial 3: Water, PCR Grade

The separate vial of Reverse Transcriptase makes it possible to prepare a RT-qPCR reaction mix for running RT-minus controls (lacking reverse transcriptase). This is important for verifying whether the obtained results are derived from RNA transcripts or from residual (contaminating) genomic DNA. The kit also provides sufficient vials of Water, PCR Grade to ensure that fresh (unopened) vials can be used. This minimizes the risk of contamination of RT-qPCR reaction mixes by RNases and other substances.

The Reverse Transcriptase provided in this kit is a recombinant reverse transcriptase with higher stability than native reverse transcriptase. This feature allows for higher reverse transcription temperatures of up to +55°C. This Roche recombinant Reverse Transcriptase also has lower affinity for DNA than other commonly used reverse transcriptases. The resulting RT-qPCR reaction mix also has AptaTaq DNA Polymerase, nucleotides, and additives, ensuring a hot start amplification system with high specificity. The mix contains an optimized concentration of MgCl₂, eliminating the need for additional adjustments. For greater convenience, the RT-qPCR mixture, including Reverse Transcriptase, can be used for both RNA and/or DNA templates in parallel.



5.2. Quality Control

The LightCycler® Multiplex RNA Virus Master is function tested using the LightCycler® System.

6. Supplementary Information

6.1. Conventions

To make information consistent and easier to read, the following text conventions and symbols are used in this document to highlight important information:

Text convention and symbols	
 <i>Information Note: Additional information about the current topic or procedure.</i>	
 Important Note: Information critical to the success of the current procedure or use of the product.	
① ② ③ etc.	Stages in a process that usually occur in the order listed.
① ② ③ etc.	Steps in a procedure that must be performed in the order listed.
* (Asterisk)	The Asterisk denotes a product available from Roche Diagnostics.

6.2. Changes to previous version

Editorial changes.

Reagents and consumables have been updated regarding new product names.

Removal of LightCycler 2.0® Instrument-specific information due to product phase-out.

6.3. Ordering Information

Roche offers a large selection of reagents and systems for life science research. For a full overview of related products and manuals, please visit and bookmark our homepage lifescience.roche.com.

Product	Pack Size	Cat. No.
Accessories general (hardware)		
Sealing Foil Applicator	1 piece	10 018 607 001
LightCycler® 8-Tube Strip Adapter Plate	1 piece, The adapter plate can be used multiple times.	06 612 598 001
LightCycler® 480 Block Kit 96 Silver	1 block kit	05 015 219 001
Thermal Cycler Assembly 384	1 piece	09 742 581 001
Thermal Cycler Assembly 96	1 piece	09 742 565 001
LightCycler® 480 Block Kit 384 Silver	1 block kit	05 015 197 001
Consumables		
LightCycler® Multiwell Plate 96	5 x 10 plates without sealing foils	05 220 319 001
LightCycler® Sealing Foil	50 foils	04 729 757 001
LightCycler® Multiwell Plate 384	5 x 10 plates without sealing foils	05 217 555 001
MagNA Pure 96 Internal Control Tube	150 tubes	06 374 905 001
LightCycler® Multiwell Plate 96, white	5 x 10 plates	04 729 692 001
LightCycler® Multiwell Plate 384, white	5 x 10 plates	04 729 749 001
LightCycler® 8-Tube Strips (white)	10 x 12 white strips and clear caps	06 612 601 001
Instruments		
MagNA Pure 96 Instrument	1 instrument	06 541 089 001
LightCycler® 96 Instrument	1 instrument	05 815 916 001
MagNA Pure 24 Instrument	Instrument with built-in control unit, software and accessories	07 290 519 001
LightCycler® 480 Instrument II	1 instrument, 96-well version	05 015 278 001
	1 instrument, 384-well version	05 015 243 001
LightCycler® PRO system	1 instrument, 96-well version	09 541 713 001
	1 instrument, 384-well version	09 582 487 001
Reagents, kits		
MagNA Pure 96 DNA and Viral NA Large Volume Kit	Kit for up to 3 x 96 isolations	06 374 891 001
Bovine Serum Albumin	20 mg, 1 mL	10 711 454 001
MagNA Pure 24 Total NA Isolation Kit	Kit for up to 96 isolations (200 µL)	07 658 036 001
LightCycler® Uracil-DNA Glycosylase	50 µL, 100 U, (2 U/µL)	03 539 806 001
High Pure Viral Nucleic Acid Kit	1 kit, up to 100 isolations	11 858 874 001
MagNA Pure 96 DNA and Viral NA Small Volume Kit	For up to 3 x 192 isolations	06 543 588 001

6. Supplementary Information

6.4. Trademarks

MAGNA PURE, APTATAQ and LIGHTCYCLER are trademarks of Roche.
All other product names and trademarks are the property of their respective owners.

6.5. License Disclaimer

For additional documentation such as certificates and safety data sheets, please visit:
documentation.roche.com.

6.6. Regulatory Disclaimer

For general laboratory use.

6.7. Safety Data Sheet

Please follow the instructions in the Safety Data Sheet (SDS).

6.8. Contact and Support

If you have questions or experience problems with this or any Roche product for Life Science, please contact our Technical Support staff. Our scientists are committed to providing rapid and effective help.
Please also contact us if you have suggestions for enhancing Roche product performance or using our products in new or specialized ways. Such customer information has repeatedly proven invaluable to the research community worldwide.

To ask questions, solve problems, suggest enhancements or report new applications, please visit our **Online Technical Support Site**.

Visit documentation.roche.com, to download or request copies of the following Materials:

- Instructions for Use
- Safety Data Sheets
- Certificates of Analysis
- Information Material

To call, write, fax, or email us, visit lifescience.roche.com and select your home country to display country-specific contact information.

