

cobas[®] SARS-CoV-2 Duo

Qualitative and quantitative assay for use on the cobas[®] 6800/8800 Systems

For Research Use Only.
Not for use in diagnostic procedures.

cobas[®] SARS-CoV-2 Duo - 192T

P/N: 09626824190

cobas[®] SARS-CoV-2 Duo Control Kit

P/N: 09626832190

cobas[®] Buffer Negative Control Kit

P/N: 07002238190 or
P/N: 09051953190

Table of contents

| | |
|---|-----------|
| Proposed use | 4 |
| Summary and explanation of the test..... | 4 |
| Reagents and materials..... | 6 |
| cobas® SARS-CoV-2 Duo reagents and controls..... | 6 |
| cobas omni reagents for sample preparation..... | 8 |
| Reagent storage and handling requirements..... | 9 |
| Additional materials required..... | 10 |
| Instrumentation and software required..... | 11 |
| Precautions and handling requirements | 12 |
| Warnings and precautions | 12 |
| Reagent handling..... | 12 |
| Good laboratory practice..... | 13 |
| Sample collection, transport, and storage..... | 13 |
| Sample collection..... | 13 |
| Transport and storage | 13 |
| Instructions for use..... | 14 |
| Procedural notes | 14 |
| Running cobas® SARS-CoV-2 Duo | 14 |
| Specimens collected in UTM-RT or UVT..... | 14 |
| cobas® SARS-CoV-2 Duo procedure | 15 |
| Results..... | 16 |
| Quality control and validity of results | 16 |
| Interpretation of results | 16 |
| cobas® SARS-CoV-2 Duo for System Software v1.4 or higher..... | 16 |
| Procedural limitations..... | 18 |

Additional information 19

Key test features 19

Symbols 20

Technical support 21

Manufacturer, importer, and distributor 21

Trademarks and patents 21

Copyright 21

References 22

Document revision 23

Proposed use

The **cobas**® SARS-CoV-2 Duo for use on the **cobas**® 6800/8800 Systems (**cobas**® SARS-CoV-2 Duo) is an automated real-time RT-PCR assay for the in vitro qualitative and quantitative detection of SARS-CoV-2 RNA in nasal and nasopharyngeal swab specimens.

cobas® SARS-CoV-2 Duo is intended for **research use only and is not for use in diagnostic procedures**.

Summary and explanation of the test

Explanation of the test

cobas® SARS-CoV-2 Duo is an automated real-time RT-PCR assay for the in vitro qualitative and quantitative detection of SARS-CoV-2 RNA in collected nasal and nasopharyngeal swab specimens collected in Copan Universal Transport Medium System (UTM-RT) or BD™ Universal Viral Transport System (UVT) from individuals suspected of COVID-19. The viral load is quantified against a non-SARS-CoV-2 armored RNA quantitation standard (RNA-QS), which is introduced into each specimen during sample preparation. The RNA-QS also functions as an internal control to monitor the entire sample preparation and PCR amplification process. In addition, the test utilizes three external controls: a high titer positive, a low titer positive, and a negative control. The high positive and low positive external controls are manufactured by dilution from stock material with a titer traceable to the First WHO International Standard for SARS-CoV-2 RNA (NIBSC code: 20/146). Each Amplification/Detection kit lot is calibrated traceable to the First WHO International Standard for SARS-CoV-2 RNA (NIBSC code: 20/146).

Principles of the procedure

cobas® SARS-CoV-2 Duo is based on fully automated sample preparation (nucleic acid extraction and purification) followed by PCR amplification and detection. The **cobas**® 6800/8800 Systems consist of the sample supply module, the transfer module, the processing module, and the analytic module. Automated data management is performed by the **cobas**® 6800/8800 software, which assigns results for all tests. Results can be reviewed directly on the system screen, and printed as a report.

Nucleic acid from patient samples and added RNA QS molecules are simultaneously extracted. Nucleic acid is released by addition of proteinase and lysis reagent to the sample. The released nucleic acid binds to the silica surface of the added magnetic glass particles. Unbound substances and impurities, such as denatured protein, cellular debris and potential PCR inhibitors, are removed with subsequent wash steps and purified nucleic acid is eluted from the magnetic glass particles with elution buffer at elevated temperature.

External controls are processed in the same way with each **cobas**® SARS-CoV-2 Duo run. The test utilizes three external controls: a high titer positive, a low titer positive, and a negative control.

Selective amplification of SARS-CoV-2 target nucleic acid from the sample is achieved by the use of a dual target virus specific approach from highly-conserved regions of SARS-CoV-2 located in the ORF 1a and ORF 1a/b non-structural regions. Selective amplification of RNA QS is achieved by the use of non-competitive sequence specific forward and reverse primers which have no homology with the SARS-CoV-2 genome.

A thermostable DNA polymerase enzyme is used for amplification. The target and RNA QS sequences are amplified simultaneously utilizing a universal PCR amplification profile with pre-defined temperature steps and number of cycles.

The master mix includes deoxyuridine triphosphate (dUTP), instead of deoxythymidine triphosphate (dTTP), which is incorporated into the newly synthesized DNA (amplicon). Any contaminating amplicons from previous PCR runs are destroyed by the AmpErase enzyme [uracil-N-glycosylase], which is included in the PCR mix, when heated in the first thermal cycling step. However, newly formed amplicons are not destroyed since the AmpErase enzyme is inactivated once exposed to temperatures above 55°C.

Amplified target is detected by cleavage of fluorescently labeled oligonucleotide probes. The **cobas**[®] SARS-CoV-2 Duo master mix contains two detection probes specific for SARS-CoV-2 target sequences and one for the RNA QS. The probes are labeled with target-specific fluorescent reporter dyes allowing simultaneous detection of SARS-CoV-2 target and RNA QS in two different target channels. The fluorescent signal of the intact probe is suppressed by the quencher dye. During the PCR amplification step, hybridization of the probes to the specific single-stranded DNA template results in cleavage by the 5' to 3' nuclease activity of the DNA polymerase resulting in separation of the reporter and quencher dyes and the generation of a fluorescent signal. With each PCR cycle, increasing amounts of cleaved probes are generated and the cumulative signal of the reporter dye increases concomitantly. Real-time detection and discrimination of PCR products are accomplished by measuring the fluorescence of the released reporter dyes for the viral targets and RNA QS.

Reagents and materials

The materials provided for cobas® SARS-CoV-2 Duo can be found in Table 1. Materials required, but not provided can be found in Table 2, Table 3, Table 4, Table 7, Table 8 and Table 9.

Refer to the **Reagents and materials** section and **Precautions and handling requirements** section for the hazard information for the product.

cobas® SARS-CoV-2 Duo reagents and controls

All unopened reagents and controls shall be stored as recommended in Table 1 to Table 4.

Table 1 cobas® SARS-CoV-2 Duo

(SARS-CoV-2 Duo)

Store at 2-8°C

192 test cassette (P/N 09626824190)

| Kit components | Reagent ingredients | Quantity per kit 192 tests |
|--|---|-------------------------------|
| Proteinase Solution (PASE) | Tris buffer, < 0.05% EDTA, calcium chloride, calcium acetate, 8% proteinase, glycerol EUH210: Safety data sheet available on request. EUH208: Contains Subtilisin. May produce an allergic reaction. | 22.3 mL |
| RNA Quantitation Standard (RNA-QS) | Tris buffer, < 0.05% EDTA, < 0.001% non-Sarbecovirus related armored RNA construct containing primer and probe specific primer sequence regions (non-infectious RNA in MS2 bacteriophage), < 0.1% sodium azide | 21.2 mL |
| Elution Buffer (EB) | Tris buffer, 0.2% methyl-4 hydroxybenzoate | 21.2 mL |
| Master Mix Reagent 1 (MMX-R1) | Manganese acetate, potassium hydroxide, < 0.1% sodium azide | 7.5 mL |
| SARS-CoV-2 Duo Master Mix Reagent 2 (SARS-CoV-2 Duo MMX-R2) | Tricine buffer, potassium acetate, < 18% dimethyl sulfoxide, glycerol, < 0.1% Tween 20, EDTA, < 0.12% dATP, dCTP, dGTP, dUTPs, < 0.01% upstream and downstream SARS-CoV-2 and primers, < 0.01% Internal Control forward and reverse primers, < 0.01% fluorescent-labeled oligonucleotide probes specific for SARS-CoV-2, and the RNA Internal Control, < 0.01% oligonucleotide aptamer, < 0.1% Z05D DNA polymerase, < 0.10% AmpErase (uracil-N-glycosylase) enzyme (microbial), < 0.1% sodium azide | 9.7 mL |

Table 2 cobas® SARS-CoV-2 Duo Control Kit**(SARS-CoV-2 Duo CTL)**

Store at 2–8°C

(P/N 09626832190)

| Kit components | Reagent ingredients | Quantity per kit |
|---|---|-------------------------|
| SARS-CoV-2 Low Positive Control (SARS-CoV-2 L(+)C) | Tris buffer, < 0.05% Sodium azide, < 0.005% EDTA, < 0.003% Poly rA, < 0.01% Non-infectious plasmid DNA (microbial) containing SARS-CoV-2 sequence | 5.2 mL (8 x 0.65 mL) |
| SARS-CoV-2 High Positive Control (SARS-CoV-2 H(+)C) | Tris buffer, < 0.05% Sodium azide, < 0.005% EDTA, < 0.003% Poly rA, < 0.01% Non-infectious plasmid DNA (microbial) containing SARS-CoV-2 sequence | 5.2 mL (8 x 0.65 mL) |

Table 3 cobas® Buffer Negative Control Kit


Store at 2–8°C

(P/N 07002238190 or 09051953190)

| Kit components | Reagent ingredients | Quantity per kit |
|---|--|-------------------------|
| cobas® Buffer Negative Control (BUF (-) C) | Tris buffer, < 0.1% sodium azide, EDTA, < 0.002% Poly rA RNA (synthetic) | 16 mL (16 x 1 mL) |

cobas omni reagents for sample preparation

Table 4 cobas omni reagents for sample preparation*

| Reagents | Reagent ingredients | Quantity per kit | Safety symbol and warning** |
|--|---|------------------|---|
| cobas omni MGP Reagent (MGP) Store at 2–8°C (P/N 06997546190) | Magnetic glass particles, Tris buffer, 0.1% methyl-4 hydroxybenzoate, < 0.1% sodium azide | 480 tests | Not applicable |
| cobas omni Specimen Diluent (SPEC DIL) Store at 2–8°C (P/N 06997511190) | Tris buffer, 0.1% methyl-4 hydroxybenzoate, < 0.1% sodium azide | 4 x 875 mL | Not applicable |
| cobas omni Lysis Reagent (LYS) Store at 2–8°C (P/N 06997538190) | 43% (w/w) guanidine thiocyanate***, 5% (w/v) polydocanol***, 2% (w/v) dithiothreitol***, dihydro sodium citrate | 4 x 875 mL |  <p>DANGER</p> <p>H302 + H332: Harmful if swallowed or if inhaled. H314: Causes severe skin burns and eye damage. H411: Toxic to aquatic life with long lasting effects. EUH032: Contact with acids liberates very toxic gas. P261: Avoid breathing dust/fume/gas/mist/vapours/ spray. P273: Avoid release to the environment. P280: Wear protective gloves/ protective clothing/ eye protection/ face protection/hearing protection. P303 + P361 + P353: IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water. P304 + P340 + P310: IF INHALED: Remove person to fresh air and keep comfortable for breathing. Immediately call a POISON CENTER/doctor. P305 + P351 + P338 + P310: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/ doctor. P391: Collect spillage. 593-84-0 Guanidinium thiocyanate 9002-92-0 Polidocanol 3483-12-3 (R*,R*)-1,4-dimercaptobutane-2,3-diol</p> |
| cobas omni Wash Reagent (WASH) Store at 15–30°C (P/N 06997503190) | Sodium citrate dihydrate, 0.1% methyl-4 hydroxybenzoate | 4.2 L | Not applicable |

* These reagents are not included in the cobas® SARS-CoV-2 Duo test kits. See listing of additional materials required (Table 7).

** Product safety labeling primarily follows EU GHS guidance

***Hazardous substance

Reagent storage and handling requirements

Reagents shall be stored and will be handled as specified in Table 5 and Table 6.

When reagents are not loaded on the cobas® 6800/8800 Systems, store them at the corresponding temperature specified in Table 5.

Table 5 Reagent storage (when reagent is not on the system)

| Reagent | Storage temperature |
|------------------------------------|---------------------|
| cobas® SARS-CoV-2 Duo - 192T | 2–8°C |
| cobas® SARS-CoV-2 Duo Control Kit | 2–8°C |
| cobas® Buffer Negative Control Kit | 2–8°C |
| cobas omni Lysis Reagent | 2–8°C |
| cobas omni MGP Reagent | 2–8°C |
| cobas omni Specimen Diluent | 2–8°C |
| cobas omni Wash Reagent | 15–30°C |

Reagents loaded onto the cobas® 6800/8800 Systems are stored at appropriate temperatures and their expiration is monitored by the system. The cobas® 6800/8800 Systems allow reagents to be used only if all of the conditions shown in Table 6 are met. The system automatically prevents use of expired reagents. Table 6 allows the user to understand the reagent handling conditions enforced by the cobas® 6800/8800 Systems.

Table 6 Reagent expiry conditions enforced by the cobas® 6800/8800 Systems

| Reagent | Kit expiration date | Open-kit stability | Number of runs for which this kit can be used | On-board stability (cumulative time on board outside refrigerator) |
|------------------------------------|------------------------------|---------------------------------------|---|--|
| cobas® SARS-CoV-2 Duo - 192T | Date not passed [†] | 90 days from first usage [†] | Max 40 runs [†] | Max 40 hours [†] |
| cobas® SARS-CoV-2 Duo Control Kit | Date not passed [†] | Not applicable ^a | Not applicable | Max 8 hours [†] |
| cobas® Buffer Negative Control Kit | Date not passed | Not applicable ^a | Not applicable | Max 10 hours |
| cobas omni Lysis Reagent | Date not passed | 30 days from loading* | Not applicable | Not applicable |
| cobas omni MGP Reagent | Date not passed | 30 days from loading* | Not applicable | Not applicable |
| cobas omni Specimen Diluent | Date not passed | 30 days from loading* | Not applicable | Not applicable |
| cobas omni Wash Reagent | Date not passed | 30 days from loading* | Not applicable | Not applicable |

^a Single use reagents

*Time is measured from the first time that reagent is loaded onto the cobas® 6800/8800 Systems.

[†] The performance has not been established for suggested use cycles and time, but is based on similar reagents used on the same system.

Additional materials required

Table 7 Materials and consumables for use on **cobas®** 6800/8800 Systems

| Material | P/N |
|---|--|
| cobas omni Processing Plate | 05534917001 |
| cobas omni Amplification Plate | 05534941001 |
| cobas omni Pipette Tips | 05534925001 |
| cobas omni Liquid Waste Container | 07094388001 |
| cobas omni Lysis Reagent | 06997538190 |
| cobas omni MGP Reagent | 06997546190 |
| cobas omni Specimen Diluent | 06997511190 |
| cobas omni Wash Reagent | 06997503190 |
| Solid Waste Bag and Solid Waste Container or Solid Waste Bag With Insert and Kit Drawer | 07435967001 and 07094361001 or 08030073001 and 08387281001 |
| cobas omni Secondary Tubes 13x75 (optional) | 06438776001 |

Instrumentation and software required

The cobas® 6800/8800 software v1.4 and cobas® SARS-CoV-2 Duo analysis package must be installed on the instrument(s). The Instrument Gateway (IG) server will be provided with the system.

Table 8 Instrumentation

| Equipment | P/N |
|--------------------------------------|-----------------------------|
| cobas® 6800 System (Option Moveable) | 05524245001 and 06379672001 |
| cobas® 6800 System (Fix) | 05524245001 and 06379664001 |
| cobas® 8800 System | 05412722001 |
| Sample Supply Module | 06301037001 |
| Instrument Gateway | 06349595001 |

For additional information, please refer to the cobas® 6800/8800 Systems – User Assistance and/or User Guide.

Note: Contact your local Roche representative for a detailed order list for sample racks, racks for clotted tips and rack trays accepted on the instruments.

Precautions and handling requirements

Warnings and precautions

As with any test procedure, good laboratory practice is essential to the proper performance of this assay. Due to the high sensitivity of this test, care should be taken to keep reagents and amplification mixtures free of contamination.

- For Research Use Only. Not for use in diagnostic procedures.
- All patient samples should be handled as if infectious, using good laboratory procedures as outlined in Biosafety in Microbiological and Biomedical Laboratories and in the CLSI Document M29-A4.^{1,2} Only personnel proficient in handling infectious materials and the use of **cobas® SARS-CoV-2 Duo** and **cobas® 6800/8800 Systems** should perform this procedure.
- All human-sourced materials should be considered potentially infectious and should be handled with universal precautions. If spillage occurs, immediately disinfect with a freshly prepared solution of 0.5% sodium hypochlorite in distilled or deionized water (dilute household bleach 1:10) or follow appropriate site procedures.
- The use of sterile disposable pipettes and nuclease-free pipette tips is recommended. Use only supplied or specified required consumables to ensure optimal test performance.
- Safety Data Sheets (SDS) are available on request from your local Roche representative.
- Closely follow procedures and guidelines provided to ensure that the test is performed correctly. Any deviation from the procedures and guidelines may affect optimal test performance.
- False positive results may occur if carryover of samples is not adequately controlled during sample handling and processing.

Reagent handling

- Handle all reagents, controls, and samples according to good laboratory practice in order to prevent carryover of samples or controls.
- Before use, visually inspect each reagent cassette, diluent, lysis reagent, and wash reagent to ensure that there are no signs of leakage. If there is any evidence of leakage, do not use that material for testing.
- **cobas omni** Lysis Reagent contains guanidine thiocyanate, a potentially hazardous chemical. Avoid contact of reagents with the skin, eyes, or mucous membranes. If contact does occur, immediately wash with generous amounts of water; otherwise, burns can occur.
- **cobas® SARS-CoV-2 Duo** test kits, **cobas® SARS-CoV-2 Duo Control kit**, **cobas® Buffer Negative Control kit**, **cobas omni** MGP Reagent, and **cobas omni** Specimen Diluent contain sodium azide as a preservative. Avoid contact of reagents with the skin, eyes, or mucous membranes. If contact does occur, immediately wash with generous amounts of water; otherwise, burns can occur. If these reagents are spilled, dilute with water before wiping dry.
- Do not allow **cobas omni** Lysis Reagent, which contains guanidine thiocyanate, to contact sodium hypochlorite (bleach) solution. This mixture can produce a highly toxic gas.
- Dispose of all materials that have come in contact with samples and reagents in accordance with country, state, and local regulations.

Good laboratory practice

- Do not pipette by mouth.
- Do not eat, drink, or smoke in designated work areas.
- Wear laboratory gloves, laboratory coats, and eye protection when handling samples and reagents. Gloves must be changed between handling samples and cobas® SARS-CoV-2 Duo kits, cobas® SARS-CoV-2 Duo Control kit, cobas® Buffer Negative Control kit and cobas omni reagents to prevent contamination. Avoid contaminating gloves when handling samples and controls.
- Wash hands thoroughly after handling samples and kit reagents, and after removing the gloves.
- Thoroughly clean and disinfect all laboratory work surfaces with a freshly prepared solution of 0.5% sodium hypochlorite in distilled or deionized water (dilute household bleach 1:10). Follow by wiping the surface with 70% ethanol.
- If spills occur on the cobas® 6800/8800 instrument, follow the instructions in the cobas® 6800/8800 Systems – User Assistance and/or User Guide to properly clean and decontaminate the surface of instrument(s).

Sample collection, transport, and storage

Note: Handle all samples and controls as if they are capable of transmitting infectious agents.

Sample collection

Ensure that the correct collection device is used with the appropriate sample type by referring to the table below:

Table 9 Overview of collection devices and sample types

| Collection Device | Sample Type | Sample Type |
|--|----------------|-------------|
| | Nasopharyngeal | Nasal |
| Copan Universal Transport Media (UTM-RT) | √ | √ |
| BD™ Universal Viral Transport (UVT) | √ | √ |

- Collect nasal and nasopharyngeal specimens according to standard collection technique using flocked or polyester-tipped swabs and immediately place in 3 mL of Copan Universal Transport Medium (UTM-RT) or BD™ Universal Viral Transport (UVT).
- Refer to the Instructions for Use of the Collection Devices for hazard information.

Transport and storage

- Transportation of collected specimens must comply with all applicable regulations for the transport of etiologic agents.
- Sample stability when using cobas® SARS-CoV-2 Duo has not been established for suggested temperatures and time, but is based on viability data from testing similar viruses in the UTM-RT or UVT Systems as stated in the Copan UTM-RT System Instructions For Use and in the following bullets:
 - After collection, the specimen should be stored at 2-25°C and processed within 48 hours.
 - If delivery and processing exceed 48 hours, specimens should be transported in dry ice and once in laboratory frozen at -70°C or colder.

Instructions for use

Procedural notes

- Do not use cobas® SARS-CoV-2 Duo reagents, cobas® SARS-CoV-2 Duo Control Kit, cobas® Buffer Negative Control Kit, or cobas omni reagents after their expiry dates.
- Do not reuse consumables. They are for one-time use only.
- Refer to the cobas® 6800/8800 Systems – User Assistance and/or User Guide for proper maintenance of instruments.

Running cobas® SARS-CoV-2 Duo

cobas® SARS-CoV-2 Duo can be run with a minimum required sample volume of 0.6 mL in the cobas omni secondary tube for specimens collected in Copan Universal Transport Medium (UTM-RT) and BD™ Universal Viral Transport (UVT).

Specimens collected in UTM-RT or UVT

Specimens collected in tubes compatible with cobas® 6800/8800 Systems may be loaded directly onto the cobas® 6800/8800 Systems. The swab must be removed from the sample tube prior to direct loading onto the system. Specimens collected in Copan Universal Transport Medium (UTM-RT) or BD™ Universal Viral Transport (UVT) tubes which are not compatible with the cobas® 6800/8800 Systems must be transferred into a secondary tube prior to processing on the cobas® 6800/8800 Systems. The cobas omni Secondary Tube is the preferred option.

Additional tubes for testing cobas® SARS-CoV-2 Duo are available. Contact your local Roche representative for detailed testing instructions and an order list of primary tubes and secondary tubes compatible with the instruments.

Always use caution when transferring specimens from a primary collection tube to a secondary tube.

Use pipettes with aerosol-barrier or positive-displacement tips to handle specimens.

Always use a new pipette tip for each specimen.

Ensure samples are equilibrated to room temperature prior to transfer into a cobas omni Secondary Tube.

Follow the steps below to transfer patient sample from a primary collection tube into a cobas omni Secondary Tube:

- Unscrew the primary sample tube cap.
- Lift the cap and any attached swab to allow a pipette to be inserted into the sample tube.
- Transfer 0.6 mL into the prepared barcoded secondary tube.
- Transfer secondary tube to a rack. Close the primary sample tube cap.

Table 10 Sample type selection in the user interface of the cobas® SARS-CoV-2 Duo

| Collection kit/Matrix type | Minimum volume (mL) Processing tube | Process as Sample Type |
|---|--|------------------------|
| Copan Universal Transport Medium BD™ Universal Viral Transport | 0.6 mL cobas omni Secondary tube | VTM |

cobas® SARS-CoV-2 Duo procedure

The test procedure is described in detail in the **cobas® 6800/8800 Systems – User Assistance and/or User Guide**. Figure 1 below summarizes the procedure.

Figure 1 cobas® SARS-CoV-2 Duo procedure

- 1** Log onto the system
Press Start to prepare the system
Order tests
- 2** Refill reagents and consumables as prompted by the system
 - Load test specific reagent cassette
 - Load control cassettes
 - Load pipette tips
 - Load processing plates
 - Load MGP reagent
 - Load amplification plates
 - Refill specimen diluent
 - Refill lysis reagent
 - Refill wash reagent
- 3** Loading samples onto the system
 - Load sample racks and clotted tip racks onto the sample supply module
 - Confirm samples have been accepted into the transfer module
- 4** Start the run by choosing the Start manually button on the user interface or have it start automatically after 120 minutes or if the batch is full
- 5** Review and export results
- 6** Remove and cap any sample tubes meeting the minimum volume requirements if needed for future use
Clean up the instrument
 - Unload empty control cassettes
 - Empty amplification plate drawer
 - Empty liquid waste
 - Empty solid waste

Results

The **cobas**® 6800/8800 Systems automatically detect the SARS-CoV-2 RNA and determine the RNA concentration, for each sample and control. Individual target results for samples as well as test validity and overall results for controls are displayed on the user interface. The SARS-CoV-2 RNA concentration is expressed in International Units per milliliter (IU/mL).

Quality control and validity of results

- One **cobas**® Buffer Negative Control [(-) Ctrl] and two positive controls, a low positive control [SARS-CoV-2 L (+) C] and a high positive control [SARS-CoV-2 H (+) C] are processed with each batch.
- In the **cobas**® 6800/8800 software and/or report, check for flags and their associated results to ensure the batch validity.
- All flags are described in the **cobas**® 6800/8800 Systems User Guide.
- The batch is valid if no flags appear for any controls. If the batch is invalid, repeat testing of the entire batch.

Validation of results is performed automatically by the **cobas**® 6800/8800 software based on negative and positive control performance.

Interpretation of results

cobas® SARS-CoV-2 Duo for System Software v1.4 or higher

Display examples for **cobas**® SARS-CoV-2 Duo for System Software v1.4 or higher are shown in Figure 2.

Figure 2 Example of **cobas**® SARS-CoV-2 Duo results display for System Software v1.4 or higher

| Test | Sample ID | Valid | Flags | Sample Type | Overall Result | Target 1 (Quantitative) | Target 2 (Qualitative) |
|----------------|------------------------|-------|-------|--------------------|---------------------|-------------------------|------------------------|
| SARS-CoV-2-Duo | Sample_01 | Yes | | VTM | Target Not Detected | Target Not Detected | Negative |
| SARS-CoV-2-Duo | Sample_02 | No | Y40T | VTM | Invalid | Invalid | Invalid |
| SARS-CoV-2-Duo | Sample_03 | Yes | | VTM | Titer | 4.87e+007 IU/ml | Positive |
| SARS-CoV-2-Duo | Sample_04 | Yes | | VTM | > Titer max | > Titer max | Positive |
| SARS-CoV-2-Duo | Sample_05 | Yes | | VTM | < Titer min | < Titer min | Positive |
| SARS-CoV-2-Duo | C161420284090428828404 | Yes | | SARS-CoV-2 H (+) C | Titer | 17.75* | 17.75* |
| SARS-CoV-2-Duo | C161420284093009580264 | Yes | | SARS-CoV-2 L (+) C | Titer | 32.98* | 32.98* |
| SARS-CoV-2-Duo | C161420284093009554953 | Yes | | (-) Ctrl | Target Not Detected | - | - |

* Illustrative Ct values

For a valid batch, check each individual sample for flags in the **cobas**® 6800/8800 software and/or report. The result interpretation should be as follows:

A valid batch may include both valid and invalid sample results.

Table 11 Target results for individual target result interpretation

| Results (quantitative) | Results (qualitative) | Interpretation |
|--------------------------|-----------------------|--|
| Target Not Detected | Negative | SARS-CoV-2 RNA not detected. Report results as "SARS-CoV-2 not detected." |
| < Titer min | Positive | SARS-CoV-2 RNA is detected. Calculated titer is below the Lower Limit of Quantitation (LLoQ) of the assay. Report results as "SARS-CoV-2 detected, less than (Titer min)" Titer min = 100 IU/mL |
| Titer (IU/mL) | Positive | SARS-CoV-2 RNA is detected. Calculated titer is within the Linear Range of the assay – greater than or equal to Titer min and less than or equal to Titer max. Report results as "(Titer) of SARS-CoV-2 detected." |
| > Titer max ^a | Positive | SARS-CoV-2 RNA is detected. Calculated titer is above the Upper Limit of Quantitation (ULoQ) of the assay. Report results as "SARS-CoV-2 detected, greater than (Titer max)." Titer max = 1.00E+09 IU/mL |
| Invalid | Invalid | Results are invalid. Sample should be retested. If the result is still invalid, a new specimen should be obtained. |

^a Sample result > Titer max refers to SARS-CoV-2 positive samples detected with titers above the upper limit of quantitation (ULoQ). If a quantitative result is desired, the original sample should be diluted with SARS-CoV-2-negative transport media depending on the type of the original sample, and the test should be repeated. Multiply the reported result by the dilution factor.

Procedural limitations

- **cobas**® SARS-CoV-2 Duo has been evaluated only for use in combination with the **cobas**® SARS-CoV-2 Duo Control Kit, **cobas**® Buffer Negative Control Kit, **cobas omni** MGP Reagent, **cobas omni** Lysis Reagent, **cobas omni** Specimen Diluent, and **cobas omni** Wash Reagent for use on the **cobas**® 6800/8800 Systems.
- Reliable results depend on proper sample collection, storage and handling procedures.
- This test can be used for the detection of SARS-CoV-2 RNA in nasal and nasopharyngeal swab samples collected in a Copan UTM-RT System (UTM-RT) or BD™ Universal Viral Transport System (UVT). Testing of other sample types with **cobas**® SARS-CoV-2 Duo may result in inaccurate results.
- Detection of SARS-CoV-2 RNA may be affected by sample collection methods, patient factors (e.g., presence of symptoms), and/or stage of infection.
- As with any molecular test, mutations within the target regions of **cobas**® SARS-CoV-2 Duo could affect primer and/or probe binding resulting in failure to detect the presence of virus.
- Due to inherent differences between technologies, it is recommended that, prior to switching from one technology to the next, users perform method correlation studies in their laboratory to qualify technology differences. One hundred percent agreement between the results should not be expected due to aforementioned differences between technologies. Users should follow their own specific policies/procedures.
- False negative or invalid results may occur due to interference. The Internal Control is included in **cobas**® SARS-CoV-2 Duo to help identify the specimens containing substances that may interfere with nucleic acid isolation and PCR amplification.
- The addition of AmpErase enzyme into the **cobas**® SARS-CoV-2 Duo Master Mix reagent enables selective amplification of target RNA; however, good laboratory practices and careful adherence to the procedures specified in this Instructions For Use document are necessary to avoid contamination of reagents.

Additional information

Key test features





















































| | |
|--|---|
| Sample type | Nasopharyngeal swab samples collected in the Copan UTM-RT System or the BD™ UVT System Nasal swab samples collected in the Copan UTM-RT System or the BD™ UVT System |
| Minimum amount of sample required | 0.6 mL* |
| Sample processing volume | 0.4 mL |

*Dead volume of 0.2 mL should be considered for the **cobas omni** Secondary tubes. Other tubes compatible with **cobas**® 6800/8800 Systems (consult User Assistance Guide) may have different dead volume and require more or less minimum volume.

Symbols

The following symbols are used in labeling for Roche PCR diagnostic products.

Table 12 Symbols used in labeling for Roche PCR diagnostics products

| | | |
|---|--|---|
|  Age or Date of Birth |  Device not for near-patient testing |  QS IU/PCR QS IU per PCR reaction, use the QS International Units (IU) per PCR reaction in calculation of the results. |
|  Ancillary Software |  Device not for self-testing | |
|  Assigned Range [copies/mL] Assigned Range (copies/mL) |  Distributor <i>(Note: The applicable country/region may be designated beneath the symbol)</i> |  Serial number |
|  Assigned Range [IU/mL] Assigned Range (IU/mL) |  Do not re-use |  Site |
|  Authorized representative in the European Community |  Female |  Standard Procedure |
|  Barcode Data Sheet |  For IVD performance evaluation only |  Sterilized using ethylene oxide |
|  Batch code |  Global Trade Item Number |  Store in dark |
|  Biological risks |  Importer |  Temperature limit |
|  Catalogue number |  In vitro diagnostic medical device |  Test Definition File |
|  CE marking of conformity; this device is in conformity with the applicable requirements for CE marking of an in vitro diagnostic medical device |  Lower Limit of Assigned Range |  This way up |
| |  Male |  Ultrasensitive Procedure |
|  Collect date |  Manufacturer |  Unique Device Identifier |
|  Consult instructions for use |  Negative control |  Upper Limit of Assigned Range |
|  Contains sufficient for <n> tests |  Non-sterile |  Urine Fill Line |
|  Content of kit |  Patient Name |  US Only: Federal law restricts this device to sale by or on the order of a physician. |
|  Control |  Patient number |  Use-by date |
|  Date of manufacture |  Peel here | |
|  Device for near-patient testing |  Positive control | |
|  Device for self-testing |  QS copies / PCR QS copies per PCR reaction, use the QS copies per PCR reaction in calculation of the results. | |

Technical support

For technical support (assistance) please reach out to your local affiliate:

https://www.roche.com/about/business/roche_worldwide.htm

Manufacturer, importer, and distributor

Table 13 Manufacturer, importer, and distributor



Roche Molecular Systems, Inc.
1080 US Highway 202 South
Branchburg, NJ 08876 USA
www.roche.com

Made in USA



Roche Diagnostics GmbH
Sandhofer Strasse 116
68305 Mannheim, Germany

EC Importer¹

Distributed by

Roche Diagnostics
9115 Hague Road
Indianapolis, IN 46250-0457 USA
(For Technical Assistance call the
Roche Response Center
toll-free: 1-800-526-1247)

1 Symbol text required by USA only.

2 For USA only.

Trademarks and patents

See <http://www.roche-diagnostics.us/patents>

Copyright

©2022 Roche Molecular Systems, Inc.

References

1. Centers for Disease Control and Prevention. Biosafety in Microbiological and Biomedical Laboratories, 6th ed. U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, National Institutes of Health HHS Publication No. (CDC) 300859, revised June 2020.
2. Clinical and Laboratory Standards Institute (CLSI). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline, 4th ed. CLSI Document M29-A4. Wayne, PA: CLSI, 2014.

Document revision

| Document Revision Information | |
|-------------------------------|-------------------|
| Doc Rev. 1.0 04/2022 | First Publishing. |