
**Handheld Base Unit**  
*Operator's Manual*

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## Revision History

<b>Manual version</b>	<b>Revision date</b>	<b>Changes</b>
Version 1.0	2008-07	Initial Version
Version 2.0	2015-05	Added reference to CoaguChek XS Pro and CoaguChek Pro II systems, updated safety related information, updated section "Contact Roche", minor changes

# **Handheld Base Unit**

## **Operator's Manual**

**Version 2.0**

© 2008-2015 Roche Diagnostics. All rights reserved.

The contents of this document, including all graphics, are the property of Roche Diagnostics. No part of this document may be reproduced or transmitted in any form or by any means, electronic or mechanical, for any purpose, without the express written permission of Roche Diagnostics.

Roche Diagnostics has made every reasonable effort to ensure that all the information contained in this manual is correct at the time of printing. However, Roche Diagnostics reserves the right to make any changes necessary without notice as part of ongoing product development.

Please send questions or comments about this manual to your local Roche representative.

COAGUCHEK, COBAS and COBAS H are trademarks of Roche.

All other trademarks are the property of their respective owners.

On the packaging, or the Handheld Base Unit you may encounter the following symbols, shown here with their meaning

Caution, consult accompanying documents. Refer to safety-related notes in the instructions for use accompanying this product.

Temperature limitation (Store at)

Manufacturer

Catalog number

In vitro diagnostic medical device

This product fulfills the requirements of the European Directive 98/79/EC on *in vitro* diagnostic medical devices.

The system fulfills the Canadian and U.S. safety requirements in accordance with CAN/CSA C22.2 No. 61010-1 and UL 61010-1.

Power supply connection

This page intentionally left blank.

<b>1</b>	<b>Introduction</b>	9
	Purpose of document .....	9
	Contents of the pack .....	9
	If you need help .....	9
	Labeling of safety information .....	10
	Safety information.....	11
	Disposal of the system.....	12
	General care.....	13
	Electromagnetic Interference .....	13
	HBU operating conditions .....	13
<b>2</b>	<b>Overview of the Handheld Base Unit (HBU)</b>	14
	Connection ports.....	17
	Configuration switch.....	18
<b>3</b>	<b>Operating the device</b>	19
3.1	Connecting the CoaguChek XS Plus, CoaguChek XS Pro, and CoaguChek Pro II meter .....	20
3.2	Connecting the <b>cobas h 232</b> meter .....	21
	Important notes.....	24
	Self testing.....	24
<b>4</b>	<b>Cleaning</b>	25
	Introduction .....	25
	Cleaning agents.....	25
	Cleaning the Handheld Base Unit .....	26
<b>5</b>	<b>Troubleshooting</b>	27
<b>6</b>	<b>Product specification</b>	28
	Technical data and operating conditions .....	28
	Ordering.....	29
	Repairs .....	29
	Contact Roche.....	29
<b>7</b>	<b>Information about software licenses</b>	30

---

# 1 Introduction

## *Purpose of document*

This document describes the use of the Handheld Base Unit (also called the **HBU** in this document). The HBU is used for connectivity between a data management system (DMS) and the following Point of Care (POC) meters:

- CoaguChek XS Plus meter
- CoaguChek XS Pro meter
- CoaguChek Pro II meter
- **cobas h** 232 meter

## *Contents of the pack*

- Handheld Base Unit
- Operator's Manual
- USB data cable
- Power supply unit
- Extension piece for **cobas h** 232 meter
- Back cover for use with CoaguChek meters

## *If you need help*

For all questions about the HBU that are not answered in this manual, contact your Roche representative. In order to expedite troubleshooting, please have your device, its serial number, and this manual ready when you call.

## ***Labeling of safety information***

This section explains how safety-related messages and information related to the proper handling of the system are presented in the Handheld Base Unit manual. Please read these passages carefully.

The safety alert symbol alone (without a signal word) promotes awareness to hazards which are generic or directs the reader to related safety information.

These symbols and signal words are used for specific hazards:

### **WARNING**

Indicates a hazardous situation which, if not avoided, could result in death or serious injury.

### **CAUTION**

Indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.

### **NOTICE**

Indicates a hazardous situation which, if not avoided, may result in damage to the system.

Important information that is not safety relevant is presented against a colored background (without a symbol). Here you will find additional information on correct use of the Handheld Base Unit or useful tips.

---

## **Safety information**

---

### **Operator qualification**

Only trained healthcare professionals may operate the Handheld Base Unit. Operators must have received comprehensive instruction in the operation and care of the Handheld Base Unit.

---

### **WARNING**

#### **Protection against infection and blood-borne pathogens**

Healthcare professionals using the Handheld Base Unit must be aware that any object coming into contact with human blood is a potential source of infection. Operators need to adhere to Standard Precautions when handling or using the Handheld Base Unit. All parts of this system should be considered potentially infectious and are capable of transmitting blood-borne pathogens.

- Use gloves.
  - Wash hands thoroughly with soap and water before putting on a new pair of gloves.
  - Follow all health and safety regulations in force locally.
-

---

 **CAUTION****Allergy or injury caused by reagents and other working solutions**

Direct contact with reagents, detergents, cleaning/disinfection solutions, or other working solutions may cause skin irritation or inflammation.

- Always use protective gloves.
  - Observe the cautions given in the package inserts of the reagents and cleaning/disinfection solutions.
  - If a reagent, control, linearity, or cleaning/disinfection solution comes into contact with your skin, wash it off immediately with water.
  - Follow all health and safety regulations in force locally.
- 

***Disposal of the system*** **WARNING****Infection by a potentially biohazardous instrument**

The Handheld Base Unit or its components must be treated as potentially biohazardous waste. Decontamination (i.e., a combination of processes including cleaning, disinfection and/or sterilization) is required before reuse, recycling, or disposal.

Dispose of the system or its components according to the appropriate local regulations. For more information, contact your Roche representative.

---

---

## **General care**

### **NOTICE**

Clean the Handheld Base Unit only with the solutions recommended (see chapter 4). Using other solutions may result in incorrect operation and possible system failure. Make sure that the Handheld Base Unit is thoroughly dried after cleaning or disinfecting.

---

## **Electromagnetic Interference**

The Handheld Base Unit fulfills the IEC 61326-2-6 requirements for emitted interference and interference immunity.

---

Do not use the Handheld Base Unit near strong electromagnetic fields, which could interfere with the proper operation of the Handheld Base Unit.

Electrostatic discharges may cause malfunction of the Handheld Base Unit.

---

## **HBU operating conditions**

To ensure that the Handheld Base Unit functions properly, please observe the following guidelines:

- Only use the HBU at a room temperature between 18 °C and 32 °C (65 °F and 90 °F).
- Only use the device at a relative humidity between 10% and 85% (no condensation).
- When operating the device using the power adapter, use only a voltage of 100V to 240 V ( $\pm 10\%$ ), 50/60 Hz.

## 2 Overview of the Handheld Base Unit (HBU)

The HBU is used for:

- Charging the rechargeable meter battery pack when the meter is docked.
- Supporting communication with a data management system (DMS) via USB and Ethernet (RJ45) data ports.

The device has the following elements:

- 1 Charging contacts
- 2 Infrared window for communication with the meter
- 3 Status LED (lights up when power is connected), charge indicator
- 4 Extension piece for **cobas h 232** meter
- 5 Removable back cover for use with CoaguChek meters

*HBU status display*

- Lights up red: Power supply is connected, application is starting up
  - Lights up green: Ready
  - Flashes red (once per second): Error
  - Lights up blue: Configuration mode
-

**Connection ports**

The following connection ports are located at the back of the Handheld Base Unit:

- 6** Network connection (RJ45)
- 7** USB connection
- 8** Power input jack for the power supply unit
- 9** Removable cover for configuration switch

**Configuration switch**


The configuration switch on the back sets the HBU to one of three different modes of operation:

Pos	Mode
1	Configuration
2	Serial & Ethernet
3	Ethernet

For further information please refer to the *ROCHE HBU-BU-BUH Technical Note.pdf*, which can be found on the HBU (attached to a computer, acting as a removable disk). See page 19 for instructions on how to access the PDF file.

### 3 Operating the device

To help the System Administrator configure the HBU, technical documentation is provided. This technical note is stored on the HBU as a PDF file. You can access this file by connecting the HBU to a PC with the USB cable provided. Follow the directions below.

- 1 Connect the power cord to the power input jack (8). After a minimum of 30 seconds, the status display shows that the HBU is ready.
- 2 Remove the cover (9) to access the configuration switch of the HBU.
- 3 Make sure that the switch is at position 1 (to your left in the illustration, see page 18). This is the default position when shipped.
- 4 Connect the HBU to a PC using the USB cable. The HBU appears as a removable disk on the PC screen. The LED on the HBU lights up blue, indicating that the HBU and the PC are connected.
- 5 Copy the file *ROCHE HBU-BU-BUH Technical Note.pdf* to the desired location on your PC.
- 6 After changing the configuration settings, move the configuration switch to either position 2 or 3, depending on your local settings. Close the cover (9) again.
- 7 If you want the Handheld Base Unit to communicate
  - via USB: connect the USB cable connector to the appropriate port (7) or
  - via Ethernet: connect the RJ45 cable connector to the appropriate port (6).


You should plug in either a USB cable connector or an Ethernet cable connector, but not both at the same time.

Make sure the configuration switch is set to the correct position for the operating mode you selected.

### 3.1 Connecting the CoaguChek XS Plus, CoaguChek XS Pro, and CoaguChek Pro II meter

---


Place the CoaguChek XS Plus meter, the CoaguChek XS Pro meter or CoaguChek Pro II meter on the HBU, pushing it gently towards the charging contacts, until it snaps into position.


When the meter is correctly docked on the HBU and all cables are connected properly, battery charging and data transfer will start automatically.

## 3.2 Connecting the cobas h 232 meter

For using the HBU with a **cobas h 232** meter, first remove the back cover (if mounted) to install the required extension piece. To remove the back cover, use appropriate tools like, e.g., two slotted screwdrivers.


- Insert the screwdrivers into the two holes and push the latches holding the cover sideways as shown.
- Once the cover is unlocked, pull it upwards and off the HBU.

- Install the HBU extension piece as shown above.

- Place the **cobas h 232** meter on the HBU, pushing it gently towards the charging contacts, until it snaps into position.

When the meter is correctly docked on the HBU and all cables are connected properly, battery charging and data transfer will start automatically.

**Important notes****Always...**

- ... handle the device carefully; avoid dropping or knocking it
- ... operate the device at an ambient temperature according to its operating conditions (see page 13 and/or page 28).
- ... place the device on a level, stable surface
- ... keep the device clean and dry

**Never ...**

- ... spill liquids on the device
- ... immerse the device in any liquid
- ... touch the charging terminals with metallic items, such as a paper clip, when plugged in

**Self testing**

The Handheld Base Unit performs a self-test every 7 days. The self-test causes the Handheld Base Unit software to reboot. To avoid interrupting communication between the meter and Handheld Base Unit, the software will only reboot when communication has been idle for at least 10 minutes. The test procedure itself takes about 3 minutes — in this time no communication between the meter and the host computer is possible. After completion of the self-test, the HBU will return automatically to the status it was in before the self-test began.

As part of the self-test, the integrity of the Handheld Base Unit's internal memory is checked. In case of a failure, the HBU will stay in USB mode even if the configuration switch is set to a different mode, and no serial or Ethernet communication of the Handheld Base Unit is possible.

---

## 4 Cleaning

### ***Introduction***

Healthcare professionals should wear gloves and follow their institution's infection control procedures when handling blood testing equipment. Use the following procedures as warranted.

#### **NOTICE**

Do not use any other cleaning solutions (e.g. cleaners containing the chemicals ether, polyhexanide, or prepared solutions or wipes containing a mixture of bleach and detergent) on the Handheld Base Unit. Use of other cleaning solutions could result in damage to the Handheld Base Unit.

---

### ***Cleaning agents***

Before using a cleaner on the HBU, check the active ingredients. Acceptable active ingredients are:

- Ammonium chloride solution (2 %)
- Mild soapy water
- Clorox Healthcare® Bleach Germicidal Cleaner
- Citric acid (2.5 %)
- Hydrogen peroxide (0.5 %)
- Sodium hypochlorite solution (0.6 %)
- 70 % isopropyl alcohol

Use of cleaners containing other active ingredients could result in damage to the HBU.

### ***Cleaning the Handheld Base Unit***

- Unplug the HBU before cleaning.
- Wipe the surfaces with a soft cloth slightly dampened (**not wet**) with the cleaning agents as recommended above.
- Avoid wiping the electrical connections at the back of the base unit
- **Dry** the HBU thoroughly after cleaning. Visually verify that no solution is seen in the HBU power input jack or ports at the completion of cleaning. If cleaning fluid is allowed to collect in any port, severe damage can occur to the meter and the HBU.

#### **NOTICE**

Ensure the HBU including all openings at the back (jack, ports) are thoroughly dried after cleaning or disinfecting. A flashing LED (red) on the HBU indicates a fault condition.

- If cleaning solution does get into any opening (jack, ports), unplug the HBU, then **dry the components** with a dry cloth or gauze pad before returning the meter to the HBU.
- Plug in the HBU.

## 5 Troubleshooting

Display/symptom	Root cause	Possible solution
HBU light is not on	Not connected to power supply unit or power supply unit is defective, the HBU is defective, or mains power outlet not active.	Disconnect and connect the power supply again.
HBU light flashes	Charging contacts may be shorted.	Disconnect and connect the power supply again.
Unit stays in USB configuration and LED flashes	The Handheld Base Unit performs a self-check regularly. In case of a failure, the HBU will stay in USB mode even if the configuration switch is set to a different mode, and no serial or Ethernet communication of the Handheld Base Unit is possible.	Check network connections, then disconnect and connect the power supply again.

If the error persists, contact your Roche representative.

## 6 Product specification

### *Technical data and operating conditions*

Specification	Handheld Base Unit	Power supply unit
Height	79.3 mm / 3.12 in	34 mm / 1.34 in + AC plug (28-40 mm / 1.1-1.6 in)
Width	117 mm / 4.61 in	52 mm / 2.05 in
Length	226 mm / 8.9 in 320 mm / 12.6 in with extension	88 mm / 3.46 in
Weight	380 g/13.4 oz. 410 g/14.46 oz with extension	N/A
User interface	LED	LED
Interfaces	Charging contacts IR port RJ45 Ethernet USB type B	DC connector Replaceable AC input contacts
Operating temperature	18 to 32 °C 65 to 90 °F	0 to 40 °C 32 to 104 °F
Storage temperature	-25 to 70 °C -13 to 158 °F	-40 to 70 °C -40 to 158 °F
Humidity (operating)	10 to 85% (no condensation)	N/A
Relative humidity (storage)	< 95%	N/A
Input voltage	7.5 V DC	90 to 264 V AC
Input frequency	N/A	50 to 60 Hz
Input current	1.7 A (max)	350-150 mA (max) at 100-240 V AC / 50 Hz
Data transfer rate	IR: 9.6-115 Kbit/s Ethernet: 10 Mbit/s USB: 12 Mbit/s	N/A

## Ordering

Item	Description	REF
Handheld Base Unit		04805658001
Power Supply	Handheld base unit power supply	04805666001

## Repairs

Please note that repairs and other modifications to the meter may only be performed by persons authorized by Roche Diagnostics.

## Contact Roche

For all questions about the Handheld Base Unit that are not answered in this manual, contact your Roche representative. If you do not already have contact details:

- Visit our website at [www.roche.com](http://www.roche.com). Select “Roche in your country” at the top of the page and then select your country to find the appropriate local office contact information, or:
- Visit our website at [www.coaguChek.com](http://www.coaguChek.com). Locate the “CoaguChek Worldwide” box on the page and select your country.

The Handheld Base Unit is manufactured for and distributed by:

In Australia:

Roche Diagnostics Australia Pty Limited  
 ABN 29 003 001 205  
 31 Victoria Avenue  
 Castle Hill, NSW, 2154

## 7 Information about software licenses

This product incorporates software modules developed under open source licenses. The source code of this software can be requested on a standard data exchange medium from the manufacturer at the following address:

Roche Diagnostics GmbH  
Sandhofer Str. 116  
68305 Mannheim  
Germany

The complete license agreements are stored as a pdf file (file name "*license.pdf*") on the Handheld Base Unit. The file *license.pdf* is located in the same folder as the file *ROCHE HBU-BU-BUH Technical Note.pdf*.

You can access this file by connecting the Handheld Base Unit to a PC with the USB cable. For detailed instructions on how to do this, see the section "Operating the device" in this Operator's Manual.


Everyone is permitted to copy and distribute verbatim copies of this license document, but changing it is not allowed.

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