# aPTT Lupus

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## Letteratura

- Langdell RD, Wagner RH, Brinkhous KM. Effect of antihemophilic factor on one-stage clotting tests. *J Lab Clin Med* 1953;41:637-647.
- Ignjatovic V. Activated Partial Thromboplastin Time. Haemostasis: methods and protocols (ed. Monagle P). *Methods in molecular biology* 2013;992:111-120.
- Bowyer A, Kitchen S, Markris M. The responsiveness of different APTT reagents to mild factor VIII, IX and XI deficiencies. *Int Jnl Lab Hem* 2011;33:154-158.
- Hayward CPM, Moffat KA, Liu Y. Laboratory Investigations for Bleeding Disorders. *Semin Thromb Hemost* 2012;38:742-752.
- Tripodi A, Lippi G, Plebani M. How to report results of prothrombin and activated partial thromboplastin times. *Clin Chem Lab Med* 2016;54(2):215-222.

- 6 Greaves M. Limitations of the laboratory monitoring of heparin therapy. Scientific and Standardization Committee Communications: on behalf of the Control of Anticoagulation Subcommittee of the Scientific and Standardization Committee of the International Society of Thrombosis and Haemostasis. *Thromb Haemost* 2002 Jan;87(1):163-4.
- 7 Faust AC, Kanyer D, Wittkowsky AK. Managing transitions from oral factor Xa inhibitors to unfractionated heparin infusions. *Am J Health-Syst Pharm* 2016;73:2037-2041.
- 8 Aarab R, van Es J, de Pont ACJM et al. Monitoring of unfractionated heparin in critically ill patients. *Neth J Med* 2013;71(9):466-471.
- 9 Vandiver JW, Vondracek TG. Antifactor Xa Levels versus Activated Partial Thromboplastin Time for Monitoring Unfractionated Heparin. *Pharmacotherapy* 2012;32(6):546-558.
- 10 Mehta TP, Smythe MA, Mattson JC. Strategies for Managing Heparin Therapy in Patients with Antiphospholipid Antibody Syndrome. *Pharmacotherapy* 2011;31(12):1221-1231.
- 11 Van Cott EM, Roberts AJ, Dager WE. Laboratory Monitoring of Parenteral Direct Thrombin Inhibitors. *Semin Thromb Hemost* 2017; [Epub ahead of print]
- 12 Adcock DM, Gosselin R. Direct Oral Anticoagulants (DOACs) in the Laboratory: 2015 Review. *Thromb Res* 2015;136:7-12.
- 13 Samuelson BT, Cuker A, Siegal DM et al. Laboratory assessment of the anticoagulant activity of direct oral anticoagulants. A systematic review. *Chest* 2017; 151(1):127-138.
- 14 Ames PRJ, Graf M, Archer J et al. Prolonged Activated Partial Thromboplastin Time: Difficulties in Discriminating Coexistent Factor VIII Inhibitor and Lupus Anticoagulant. *Clin Appl Thromb Hemost* 2015;21(2):149-154.
- 15 Kershaw G, Favaloro EJ. Laboratory identification of factor inhibitors: an update. *Pathology* 2012;44(4):293-302.
- 16 Tripodi A. Liver Disease and Hemostatic (Dys)function. *Semin Thromb Hemost* 2015;41:462-467.
- 17 Ha NB, Regal RE. Anticoagulation in Patients With Cirrhosis: Caught Between a Rock-Liver and a Hard Place. *Ann Pharmacother* 2016, Vol. 50(5):402-409.
- 18 Boral BM, Williams DJ, Boral LI. Disseminated Intravascular Coagulation. *Am J Clin Pathol* 2016;146:670-680.
- 19 Tripodi A, Chantarangkul V, Asti D, et al. Activated partial thromboplastin time: results of a case-control study evaluating six commercial reagents in assessing the risk of venous thromboembolism. *J Thromb Haemost* 2006;4:1407-1409.
- 20 Hron G, Eichinger S, Weltermann A, Quehenberger P, Halbmayer WM, Kyrle PA. Prediction of recurrent venous thromboembolism by the activated partial thromboplastin time. *J Thromb Haemost* 2006;4:752-6.
- 21 Tripodi A, Chantarangkul V, Martinelli I, et al. A shortened activated partial thromboplastin time is associated with the risk of venous thromboembolism. *Blood* 2004;104:12:3631-3634.
- 22 Boekel E, Bartels Piet. Abnormally Short Activated Partial Thromboplastin Times Are Related to Elevated Plasma Levels of TAT, F1+2, D-Dimer and FVIII:C. *Pathophysiol Haemost Thromb* 2002;32:137-142.
- 23 Sørensen B, Ingerslev J. Dynamic APTT parameters: applications in thrombophilia. *J Thromb Haemost* 2012;10:244-50.
- 24 Lin CH, Kuo YW, Kuo CY, et al. Shortened Activated Partial Thromboplastin Time Is Associated With Acute Ischemic Stroke, Stroke Severity, and Neurological Worsening. *J Stroke Cerebrovasc Dis* 2015;24(10):2270-2276.
- 25 CLSI Document H21-A5, Vol.28, No.5, 2008. Collection, transport, and processing of blood specimens for testing plasma-based coagulation assays and molecular hemostasis assays; approved guideline, 5th edition.
- 26 CLSI Document H3-A6. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard - Sixth Edition, vol. 27, No. 26, 2007.
- 27 van den Besselaar AMHP, Meeuwisse-Braun J, Jansen-Gruter R, et al. Monitoring heparin therapy by the activated partial thromboplastin time: die effect of pre-analytical conditions. *Thromb Haenwst.* 1987;57:226-231.
- 28 Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. *Clin Chem* 1986;32:470-475.
- 29 Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". *Eur J Clin Chem Clin Biochem* 1996;34:385-386.
- 30 Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. *Ann Clin Biochem* 2001;38:376-385.
- 31 Devreese KMJ, Verfaillie CJ, De Bisschop F et al. Interference of C-reactive protein with clotting times. *Clin Chem Lab Med* 2015; 53(5): e141-e145.
- 32 Price EA, Jin J, Nguyen HM, et al. Discordant aPTT and Anti-Xa Values and Outcomes in Hospitalized Patients Treated with Intravenous Unfractionated Heparin. *Ann Pharmacother* 2013;47:151-8.
- 33 Vandiver JW, Vondracek TG. Antifactor Xa Levels versus Activated Partial Thromboplastin Time for Monitoring Unfractionated. *Pharmacotherapy* 2012;32(6):546-558.
- 34 Marlar RA, Clement B, Gausman J, et al. Activated Partial Thromboplastin Time Monitoring of Unfractionated Heparin Therapy: Issues and Recommendations. *Semin Thromb Hemost* 2017;43(3):253-260.
- 35 Cawkwell RD. Patient's age and activated partial thromboplastin time test. *Thromb Haemostasis* 1983;39:780-781.
- 36 Gallisti S, Muntean W, Leschnick B, et al. Longer aPTT Values in Healthy Children than in Adults: No Single Cause. *Thrombosis Research* 1997;88:355-359.
- 37 CLSI Document EP05-A3. Evaluation of Precision of Quantitative Measurement Procedures. Vol. 24, No. 25, 2014. Approved guideline, 3rd Edition.
- 38 Martin RF. General Deming Regression for Estimating Systematic Bias and its Confidence Interval in Method Comparison Studies. *Clinical Chemistry* 2000;46(1):100-104.

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© 2020, Roche Diagnostics
Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim