cobas h 232

Addendum to Operator's Manual V4.0

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Addendum du manuel d'utilisation V 4.0

Addendum al Manuale Operatore V 4.0

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Suplemento ao Manual do Operador, v. 4.0

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cobas[®]

Life needs answers

Document version	Revision date	Changes
Version 1.0	2016-04	New Document, created as an amendment to Operator's Manual Version 4.0

Purpose of document

This document provides corrections to information in the **cobas h 232** Operator's Manual version 4.0. This addendum is valid for and intended to be used in conjunction with the **cobas h 232** Operator's Manual version 4.0. It is not a replacement for the complete **cobas h 232** Operator's Manual version 4.0.

Be sure to read the complete **cobas h 232** Operator's Manual version 4.0.

© 2016, Roche Diagnostics. All rights reserved.

The contents of this document, including all graphics, are the property of Roche Diagnostics. No part of this document may be reproduced or transmitted in any form or by any means, electronic or mechanical, for any purpose, without the express written permission of Roche Diagnostics. Roche Diagnostics has made every reasonable effort to ensure that all the information contained in this manual is correct at the time of printing. However, Roche Diagnostics reserves the right to make any changes necessary without notice as part of ongoing product development.


Please send questions or comments about this document to your local Roche representative.

ROCHE CARDIAC, COBAS, COBAS H and IQC are trademarks of Roche.

cobas h 232 • Addendum to Operator's Manual

EN

The following text replaces the related text on TnT Elevated messages in chapter 6.2 *Performing a test* on page 107 of the Operator's Manual version 4.

For Troponin T only:

As soon as the meter detects a signal (this may take a few minutes), a qualitative message (*TnT Elevated*) appears under the hourglass icon (see illustration).

The final quantitative result is displayed when the test is completed. Decisions on further treatment of the patient should be based on the final result only. Refer to the Roche CARDIAC T Quantitative or Roche CARDIAC POC Troponin T package insert for more details.

The following warning has to be added to chapter 7 *Control Testing and Quality Control* on page 113 of the Operator's Manual version 4.

Potential risk of incorrect patient results

Failure to check the system regularly by performing quality control tests may lead to incorrect results and therefore to inappropriate therapy decisions. To ensure that the system is functioning properly adhere to the manufacturer recommendations on quality control testing given in the tables below.

Manufacturer recommended frequency of quality control tests:

Frequency	Material to test	Reason	Release criteria
Daily*	CARDIAC IQC High and CARDIAC IQC Low	Verification of the proper functioning of the optical system	Both levels are "Pass"

*only on days when the **cobas h 232** meter is used for measurements

Manufacturer recommended use:

Frequency	Material to test	Reason	Release criteria
Once a month		To verify adequate storage at customer site	
With each new shipment of test strips		To verify proper functioning after transportation to customer site	
When an implausible result is obtained	Roche CARDIAC Controls		Results are within the range provided by the manufacturer
When an operator uses the cobas h 232 system (meter with or without scanner) for the first time		To qualify the new operator; to verify proper handling and sample dosing	
After adverse events such as dropping the meter or other strong mechanical shocks.	Roche CARDIAC Controls and CARDIAC IQC High and CARDIAC IQC Low	To verify proper functioning of the meter	Results are within the range provided by the manufacturer

QC recommendations by national associations or local hospital guidelines take precedence over manufacturer recommendations. For more details on the procedure refer to the Operator's Manual and Instructions for Use.

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Roche Diagnostics GmbH