Tina-quant Lipoprotein (a) Gen.2

REF	i	CONTENT		Analizadores adecuados para el cobas c pack
08106126190	08106126500	Tina-quant Lipoprotein (a) Gen.2 (150 pruebas)	ID del sistema 2086 001	cobas c 303, cobas c 503, cobas c 703

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- Kostner KM, Kostner GM. Lipoprotein (a): a historical appraisal. *J Lipid Res* 2017 Jan;58(1):1-14.
- Schmidt K, Noureen A, Kronenberg F, et al. Structure, function, and genetics of lipoprotein (a). *J Lipid Res* 2016 Aug;57(8):1339-59. doi: 10.1194/jlr.R067314. Epub 2016 Apr 13. PMID: 27074913; PMCID: PMC4959873.
- Marcovina SM, Albers JJ. Lipoprotein (a) measurements for clinical application. *J Lipid Res* 2016 Apr;57(4):526-537.
- Maranhão RC, Carvalho PO, Strunz CC, et al. Lipoprotein (a): structure, pathophysiology and clinical implications. *Arq Bras Cardiol*. 2014 Jul;103(1):76-84.
- Nordestgaard BG, Langsted A. Lipoprotein (a) as a cause of cardiovascular disease: insights from epidemiology, genetics, and biology. *J Lipid Res*. 2016 Nov;57(11):1953-1975.
- Reyes-Soffer G, Ginsberg HN, Ramakrishnan R. The metabolism of lipoprotein (a): an ever-evolving story. *J Lipid Res* 2017 Sep;58(9):1756-1764.
- Coassin S, Kronenberg F. Lipoprotein(a) beyond the kringle IV repeat polymorphism: The complexity of genetic variation in the LPA gene. *Atherosclerosis*. 2022 May;349:17-35.
- Kronenberg F. Human Genetics and the Causal Role of Lipoprotein(a) for Various Diseases. *Cardiovasc Drugs Ther*. 2016 Feb;30(1):87-100.
- Langsted A, Kamstrup PR, Nordestgaard BG. Lipoprotein(a): fasting and nonfasting levels, inflammation, and cardiovascular risk. *Atherosclerosis*. 2014 May;234(1):95-101.
- Kronenberg F. Lipoprotein(a) in various conditions: to keep a sense of proportions. *Atherosclerosis*. 2014 May;234(1):249-251.
- Reyes-Soffer G. The impact of race and ethnicity on lipoprotein(a) levels and cardiovascular risk. *Curr Opin Lipidol*. 2021 Jun 1;32(3):163-166.
- Cobbaert CM, Althaus H, Begcevic Brkovic I, et al. IFCC Working Group for Standardization of Apolipoproteins by Mass Spectrometry. Towards an SI-Traceable Reference Measurement System for Seven Serum Apolipoproteins Using Bottom-Up Quantitative Proteomics: Conceptual Approach Enabled by Cross-Disciplinary/Cross-Sector Collaboration. *Clin Chem*. 2021 Mar 1;67(3):478-489.
- Ward NC, Watts GF, Bishop W, et al. Australian Atherosclerosis Society Position Statement on Lipoprotein(a): Clinical and Implementation Recommendations. *Heart Lung Circ*. 2023 Mar;32(3):287-296.
- Cegla J, France M, Marcovina SM, et al. Lp(a): When and how to measure it. *Ann Clin Biochem*. 2021 Jan;58(1):16-21.
- Kronenberg F, Mora S, Stroes ESG, et al. Lipoprotein(a) in atherosclerotic cardiovascular disease and aortic stenosis: a European Atherosclerosis Society consensus statement. *Eur Heart J* 2022 Oct 14;43(39):3925-3946.
- Forbes CA, Quek RG, Deshpande S, et al. The relationship between Lp(a) and CVD outcomes: a systematic review. *Lipids Health Dis*. 2016 May 17;15:95.
- Pearson K, Rodriguez F. Lipoprotein(a) and Cardiovascular Disease Prevention across Diverse Populations. *Cardiol Ther* 2020 Dec;9(2):275-292.
- Wilson DP, Jacobson TA, Jones PH, et al. Use of Lipoprotein(a) in clinical practice: A biomarker whose time has come. A scientific statement from the National Lipid Association. *J Clin Lipidol* 2019 May-Jun;13(3):374-392.
- Cegla J, Neely RDG, France M, et al. HEART UK Medical, Scientific and Research Committee. HEART UK consensus statement on Lipoprotein(a): A call to action. *Atherosclerosis*. 2019 Dec;291:62-70.
- Grundy SM, Stone NJ, Bailey AL, et al. 2018 AHA/ACC/AACVPR/AAPA/ABC/ACPM/ADA/AGS/APHA/ASPC/NLA/PCNA Guideline on the Management of Blood Cholesterol: A Report of the American College of Cardiology/American Heart Association Task Force on Clinical Practice Guidelines. *Circulation*. 2019 Jun 18;139(25):e1082-e1143.
- Mach F, Baigent C, Catapano AL, et al. ESC Scientific Document Group. 2019 ESC/EAS Guidelines for the management of dyslipidaemias: lipid modification to reduce cardiovascular risk. *Eur Heart J* 2020 Jan 1;41(1):111-188.
- Handelsman Y, Jellinger PS, Guerin CK, et al. Consensus Statement by the American Association of Clinical Endocrinologists and American College of Endocrinology on the Management of Dyslipidemia and Prevention of Cardiovascular Disease Algorithm - 2020 Executive Summary. *Endocr Pract*. 2020 Oct;26(10):1196-1224.
- Pearson GJ, Thanassoulis G, Anderson TJ, et al. 2021 Canadian Cardiovascular Society Guidelines for the Management of Dyslipidemia for the Prevention of Cardiovascular Disease in Adults. *Can J Cardiol* 2021 Aug;37(8):1129-1150.
- Kronenberg F, Mora S, Stroes ESG, et al. Frequent questions and responses on the 2022 lipoprotein(a) consensus statement of the European Atherosclerosis Society. *Atherosclerosis*. 2023 Jun;374:107-120.
- Simó JM, Camps J, Gómez F, et al. Evaluation of a Fully Automated Particle-enhanced Turbidimetric Immunoassay for the Measurement of Plasma Lipoprotein(a). Population-Based Reference Values in an Area with Low Incidence of Cardiovascular Disease. *Clin Biochem* 2003 Mar;36(2):129-134.
- National Committee for Clinical Laboratory Standards, Procedures for the handling and Processing of Blood Specimens, Approved Guideline, NCCLS publication H18-A, Villanova, 1990.
- Simó JM, Camps J, Vilella E, et al. Instability of Lipoprotein (a) in Plasma Stored at -70 °C: Effects of Concentration, Apolipoprotein (a) Genotype, and Donor Cardiovascular Disease. *Clin Chem* 2001 Sep;47(9):1673-1678.
- Sgoutas DS, Tuten T. Effect of Freezing and Thawing of Serum on the Immunoassay of Lipoprotein(a). *Clin Chem* 1992;38(9):1873-1877.
- Marcovina SM, Albers JJ, Scanu AM, et al. Use of a reference Material Proposed by the International federation of Clinical Chemistry and Laboratory Medicine to Evaluate Analytical methods for the Determination of Plasma Lipoprotein (a). *Clin Chem* 2000 Dec;46(12):1956-1967.
- Anne Langsted, Pia R. Kamstrup, Borge G. Nordestgaard. High lipoprotein(a) and high risk of mortality. *European Heart Journal* (2019) 40:2760-2770.
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. *Clin Chem* 1986;32:470-475.
- Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". *Eur J Clin Chem Clin Biochem* 1996;34:385-386.
- Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. *Ann Clin Biochem* 2001;38:376-385.
- Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. *Clin Chem Lab Med* 2007;45(9):1240-1243.
- Marcovina SM, Koschinsky ML. A Critical Evaluation of the Role of Lp(a) in Cardiovascular Disease: Can Lp(a) Be Useful in Risk Assessment? *Semin Vasc Med* 2002 Aug;2(3):335-344.
- Shai I, Rimm EB, Hankinson SE, et al. Lipoprotein (a) and Coronary Heart Disease Among Women: Beyond a Cholesterol Carrier? *Eur Heart J* 2005;26:1633-1639.
- Nordestgaard BG, Chapman MJ, Ray K, et al. Lipoprotein (a) as a cardiovascular risk factor: current status. *Eur Heart J* 2010 Dec;31(23):2844-2853.

- 38 Marcovina SM, Koschinsky ML, Albers JJ, et al. Report of the National Heart, Lung, and Blood Institute Workshop on Lipoprotein (a) and Cardiovascular Disease: Recent Advances and Future Directions. Clin Chem 2003 Nov;49(11):1785-1796.
- 39 Tsimikas S, Clopton P, Brilakis ES, et al. Relationship of Oxidized Phospholipids on Apolipoprotein B-100 Particles to Race/Ethnicity, Apolipoprotein (a) Isoform Size, and Cardiovascular Risk Factors: Results From the Dallas Heart Study, Circulation 2009 Apr;119(13):1711-1719.
- 40 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.

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© 2024, Roche Diagnostics
Roche Diagnostics GmbH