

Urea/BUN

Order information

REF		CONTENT		Analyzer(s) on which cobas c pack(s) can be used
04460715190	04460715500	Urea/BUN (500 tests)	System-ID 07 6303 9	cobas c 311 , cobas c 501/502 , COBAS INTEGRA 400 plus

Materials required (but not provided):

		cobas c 311 , cobas c 501/502	COBAS INTEGRA 400 plus
10759350190	Calibrator f.a.s. (12 x 3 mL)	Code 401	System-ID 07 3718 6
12149435122	Precinorm U plus (10 x 3 mL)	Code 300	System-ID 07 7999 7
12149443122	Precipath U plus (10 x 3 mL)	Code 301	System-ID 07 8000 6
05117003190	PreciControl ClinChem Multi 1 (20 x 5 mL)	Code 391	System-ID 07 7469 3
05947626190	PreciControl ClinChem Multi 1 (4 x 5 mL)	Code 391	System-ID 07 7469 3
05117216190	PreciControl ClinChem Multi 2 (20 x 5 mL)	Code 392	System-ID 07 7470 7
05947774190	PreciControl ClinChem Multi 2 (4 x 5 mL)	Code 392	System-ID 07 7470 7
04489357190	Diluent NaCl 9 % (50 mL)	System-ID 07 6869 3	n.a.

English

Intended use

In vitro test for the quantitative determination of urea/urea nitrogen in human serum, plasma and urine on **cobas c** and COBAS INTEGRA systems.

Summary

Measurements of urea/urea nitrogen in human serum, plasma and urine, performed with this assay, are used as screening tests and as an aid in diagnosis and monitoring of renal function.

Urea is the major end product of protein nitrogen metabolism. It is synthesized by the urea cycle in the liver from ammonia which is produced by amino acid deamination. Urea is excreted mostly by the kidneys but minimal amounts are also excreted in sweat and degraded in the intestines by bacterial action.¹ Serum urea mass concentration is either specified for the complete urea molecule or for nitrogen equivalents [blood urea nitrogen (BUN)].² Determination of blood urea nitrogen is primarily used as a screening test for renal function. When used in conjunction with serum creatinine determinations it can aid in the differential diagnosis of the three types of azotemia: prerenal, renal, and postrenal. The urea to creatinine ratio has been proposed as a crude discriminator between prerenal and intrinsic azotemia.¹ Elevations in blood urea nitrogen concentration are seen in inadequate renal perfusion, shock, diminished blood volume (prerenal causes), chronic nephritis, nephrosclerosis, tubular necrosis, glomerular-nephritis (renal causes), and urinary tract obstruction (postrenal causes). Transient elevations may also be seen during periods of high protein intake. Liver diseases may lead to unpredictable blood urea nitrogen concentrations, including abnormally low levels. Low blood urea nitrogen concentrations are not common, but can be found in cases such as malnutrition, lack of protein in the diet, or overhydration.^{1,3}


Test principle

Kinetic test with urease and glutamate dehydrogenase.^{4,5,6,7}

Urea is hydrolyzed by urease to form ammonium and carbonate.



In the second reaction 2-oxoglutarate reacts with ammonium in the presence of glutamate dehydrogenase (GLDH) and the coenzyme NADH to produce L-glutamate. In this reaction two moles of NADH are oxidized to NAD⁺ for each mole of urea hydrolyzed.



The rate of decrease in the NADH concentration is directly proportional to the urea concentration in the specimen and is measured photometrically.

Precautions and warnings

For in vitro diagnostic use for laboratory professionals. Exercise the normal precautions required for handling all laboratory reagents.

Infectious or microbial waste:

Warning: handle waste as potentially biohazardous material. Dispose of waste according to accepted laboratory instructions and procedures.

Environmental hazards:

Apply all relevant local disposal regulations to determine the safe disposal.

Safety data sheet available for professional user on request.

Reagent handling

Ready for use

Specimen collection and preparation

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum

Plasma: Li-heparin and K₂-EDTA plasma. Do not use ammonium heparin.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Urine

Bacterial growth in the specimen and high atmospheric ammonia concentrations as well as contamination by ammonium ions may cause erroneously elevated results.

Stability in <i>serum/plasma</i> . ⁸	7 days at 15-25 °C
	7 days at 2-8 °C
	1 year at (-15)-(-25) °C

Freeze only once.

Stability in <i>urine</i> . ⁸	2 days at 15-25 °C
	7 days at 2-8 °C
	1 month at (-15)-(-25) °C

Freeze only once.

Centrifuge samples containing precipitates before performing the assay.

See the limitations and interferences section for details about possible sample interferences.

Materials provided

See "Reagents – working solutions" section for reagents.

Materials required (but not provided)

See "Order information" section

General laboratory equipment

Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

The performance of applications not validated by Roche is not warranted and must be defined by the user.

Calculation

The systems automatically calculate the analyte concentration of each sample.

Conversion factors: mmol/L urea x 6.006 = mg/dL urea
 mmol/L urea x 0.06006 = g/L urea
 mmol/L urea nitrogen x 2.801 = mg/dL urea nitrogen
 mmol/L urea nitrogen x 0.02801 = g/L urea nitrogen
 mg/dL urea x 0.467 = mg/dL urea nitrogen

When 24-hour urine is used as the specimen, multiply the result by the 24-hour volume to obtain values in g or mmol/24 hours.

Expected values

Urea:

*Serum/plasma*⁹

Adults 2.76-8.07 mmol/L (16.6-48.5 mg/dL)

Urine

24-hour urine¹⁰ 428-714 mmol/24 h (25.7-42.9 g/24 h),
 corresponding to
 286-595 mmol/L (1.71-3.57 g/dL)^a

a) Based on average urine output of 1.2-1.5 L/24 h

Urea nitrogen (BUN):

*Serum/plasma*¹⁰

Adults (18-60 years) 2.14-7.14 mmol/L 6-20 mg/dL

Adults (60-90 years) 2.86-8.21 mmol/L 8-23 mg/dL

Infants (< 1 year) 1.43-6.78 mmol/L 4-19 mg/dL

Infants/children 1.79-6.43 mmol/L 5-18 mg/dL

Urine

24-hour urine¹⁰ 428-714 mmol/24 h (12-20 g/24 h),
 corresponding to
 286-595 mmol/L (801-1666 mg/dL)^b

b) Based on average urine output of 1.2-1.5 L/24 h

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

cobas c systems**System information**

For **cobas c 311** analyzer:

UREAL: ACN 418 (serum/plasma)

U-BUN: ACN 421 (serum/plasma)

URELU: ACN 417 (urine)

UBUNU: ACN 428 (urine)

SUREA: ACN 419 (STAT, reaction time: 4, serum/plasma)

SUBUN: ACN 427 (STAT, reaction time: 4, serum/plasma)

SUREU: ACN 420 (STAT, reaction time: 4, urine)

SBUNU: ACN 429 (STAT, reaction time: 4, urine)

For **cobas c 501** analyzer:

UREAL: ACN 418 (serum/plasma/urine)

U-BUN: ACN 421 (serum/plasma/urine)

SUREA: ACN 419 (STAT, reaction time: 4, serum/plasma/urine)

SUBUN: ACN 427 (STAT, reaction time: 4, serum/plasma/urine)

For **cobas c 502** analyzer:

UREAL: ACN 8418 (serum/plasma)

U-BUN: ACN 8421 (serum/plasma)

URELU: ACN 8417 (urine)

UBUNU: ACN 8428 (urine)

SUREA: ACN 8419 (STAT, reaction time: 4, serum/plasma)

SUBUN: ACN 8427 (STAT, reaction time: 4, serum/plasma)

SUREU: ACN 8420 (STAT, reaction time: 4, urine)

SBUNU: ACN 8429 (STAT, reaction time: 4, urine)

Reagents - working solutions

R1 NaCl 9 %

R2 TRIS buffer: 220 mmol/L, pH 8.6; 2-oxoglutarate: 73 mmol/L;
 NADH: 2.5 mmol/L; ADP: 6.5 mmol/L; urease (jack bean):
 ≥ 300 μkat/L; GLDH (bovine liver): ≥ 80 μkat/L; preservative;
 nonreactive stabilizers

R1 is in position C and R2 is in position B.

Storage and stability

Shelf life at 2-8 °C:

See expiration date
 on **cobas c** pack
 label.

On-board in use and refrigerated on the analyzer: 8 weeks

Application for serum and plasma**cobas c 311 test definition**

Assay type	Rate A		
Reaction time / Assay points	10 / 10-19 (STAT 4 / 10-19)		
Wavelength (sub/main)	700/340 nm		
Reaction direction	Decrease		
Units	mmol/L (mg/dL, g/L)		
Reagent pipetting	Diluent (H ₂ O)		
R1	10 μL	90 μL	
R2	38 μL	110 μL	
Sample volumes	Sample	Sample dilution	
		Sample	Diluent (NaCl)
Normal	2 μL	–	–
Decreased	6 μL	15 μL	120 μL
Increased	2 μL	–	–

cobas c 501 test definition

Assay type	Rate A		
Reaction time / Assay points	10 / 16-28 (STAT 4 / 16-28)		
Wavelength (sub/main)	700/340 nm		
Reaction direction	Decrease		
Units	mmol/L (mg/dL, g/L)		
Reagent pipetting	Diluent (H ₂ O)		
R1	10 μL	90 μL	
R2	38 μL	110 μL	
Sample volumes	Sample	Sample dilution	
		Sample	Diluent (NaCl)
Normal	2 μL	–	–
Decreased	6 μL	15 μL	120 μL

UREAL

Urea/BUN



Increased 2 µL – –

cobas c 502 test definition

Assay type	Rate A		
Reaction time / Assay points	10 / 16-28 (STAT 4 / 16-28)		
Wavelength (sub/main)	700/340 nm		
Reaction direction	Decrease		
Units	mmol/L (mg/dL, g/L)		
Reagent pipetting	Diluent (H ₂ O)		
R1	10 µL	90 µL	
R2	38 µL	110 µL	
Sample volumes	Sample	Sample dilution	
		Sample	Diluent (NaCl)
Normal	2 µL	–	–
Decreased	6 µL	15 µL	120 µL
Increased	4 µL	–	–

Application for urine

cobas c 311 test definition

Assay type	Rate A		
Reaction time / Assay points	10 / 10-19 (STAT 4 / 10-19)		
Wavelength (sub/main)	700/340 nm		
Reaction direction	Decrease		
Units	mmol/L (mg/dL, g/L)		
Reagent pipetting	Diluent (H ₂ O)		
R1	10 µL	90 µL	
R2	38 µL	110 µL	
Sample volumes	Sample	Sample dilution	
		Sample	Diluent (NaCl)
Normal	2 µL	3 µL	147 µL
Decreased	2 µL	2 µL	178 µL
Increased	2 µL	–	–

cobas c 501/502 test definition

Assay type	Rate A		
Reaction time / Assay points	10 / 16-28 (STAT 4 / 16-28)		
Wavelength (sub/main)	700/340 nm		
Reaction direction	Decrease		
Units	mmol/L (mg/dL, g/L)		
Reagent pipetting	Diluent (H ₂ O)		
R1	10 µL	90 µL	
R2	38 µL	110 µL	
Sample volumes	Sample	Sample dilution	
		Sample	Diluent (NaCl)
Normal	2 µL	3 µL	147 µL
Decreased	2 µL	2 µL	178 µL
Increased	2 µL	–	–

Calibration

Calibrators	S1: H ₂ O S2: C.f.a.s.
Calibration mode	Linear

Calibration frequency 2-point calibration

- every 4 weeks on board
- after reagent lot change
- as required following quality control procedures

Calibration interval may be extended based on acceptable verification of calibration by the laboratory.

Traceability: This method has been standardized against ID/MS.

Quality control

Serum/plasma

For quality control, use control materials as listed in the "Order information" section.

In addition, other suitable control material can be used.

Urine

Quantitative urine controls are recommended for routine quality control.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

Limitations - interference

In very rare cases, gammopathy, in particular type IgM (Waldenström's macroglobulinemia), may cause unreliable results.¹¹

Serum/plasma

Criterion: Recovery within ± 0.83 mmol/L of initial values of samples ≤ 8.3 mmol/L and within $\pm 10\%$ for samples > 8.3 mmol/L.

Icterus:¹² No significant interference up to an I index of 60 for conjugated and unconjugated bilirubin (approximate conjugated and unconjugated bilirubin concentration: 1026 µmol/L (60 mg/dL)).

Hemolysis:¹² No significant interference up to an H index of 1000 (approximate hemoglobin concentration: 621 µmol/L (1000 mg/dL)).

Lipemia (Intralipid):¹² No significant interference up to an L index of 1000. There is poor correlation between the L index (corresponds to turbidity) and triglycerides concentration.

Ammonium ions may cause erroneously elevated results.

Drugs: No interference was found at therapeutic concentrations using common drug panels.^{13,14}

Urine

Criterion: Recovery within ± 15 mmol/L of initial values of samples ≤ 150 mmol/L and within $\pm 10\%$ for samples > 150 mmol/L.

Drugs: No interference was found at therapeutic concentrations using common drug panels.¹⁴

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

ACTION REQUIRED

Special Wash Programming: The use of special wash steps is mandatory when certain test combinations are run together on **cobas c** systems. The latest version of the carry-over evasion list can be found with the NaOHD-SMS-SmpCln1+2-SCCS Method Sheets. For further instructions refer to the operator's manual. **cobas c 502** analyzer: All special wash programming necessary for avoiding carry-over is available via the **cobas** link, manual input is required in certain cases.

Where required, special wash/carry-over evasion programming must be implemented prior to reporting results with this test.

Limits and ranges

Measuring range

Serum/plasma

0.5-40 mmol/L (3.0-240 mg/dL urea, 1.4-112 mg/dL urea nitrogen)

Determine samples having higher concentrations via the rerun function. Dilution of samples via the rerun function is a 1:3 dilution. Results from samples diluted using the rerun function are automatically multiplied by a factor of 3.

Urine

1-2000 mmol/L (6-12000 mg/dL urea, 2.8-5600 mg/dL urea nitrogen)

Determine samples having higher concentrations via the rerun function. Dilution of samples via the rerun function is a 1:1.8 dilution. Results from samples diluted using the rerun function are automatically multiplied by a factor of 1.8.

Determine samples having concentrations lower than the technical limit of 40 mmol/L (240 mg/dL urea and 112 mg/dL urea nitrogen) via the rerun function. Samples are measured undiluted.

Lower limits of measurement

Lower detection limit of the test

Serum/plasma

0.5 mmol/L (3.0 mg/dL urea, 1.4 mg/dL urea nitrogen)

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying 3 standard deviations above that of the lowest standard (standard 1 + 3 SD, repeatability, n = 21).

Urine

1 mmol/L (6 mg/dL urea, 2.8 mg/dL urea nitrogen)

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying 3 standard deviations above that of the lowest standard (standard 1 + 3 SD, repeatability, n = 21).

Specific performance data

Representative performance data on the analyzers are given below. Results obtained in individual laboratories may differ.

Precision

Precision was determined using human samples and controls in an internal protocol with repeatability (n = 21) and intermediate precision (*serum/plasma*: 3 aliquots per run, 1 run per day, 21 days; *urine*: 3 aliquots per run, 1 run per day, 10 days). The following results were obtained on the **cobas c 501** analyzer:

Serum/plasma

Repeatability	Mean mmol/L (mg/dL urea)	SD mmol/L (mg/dL urea)	CV %
Precinorm U	6.74 (40.5)	0.07 (0.4)	1.0
Precipath U	23.4 (141)	0.2 (1)	0.9
Human serum 1	9.18 (55.1)	0.09 (0.5)	1.0
Human serum 2	15.1 (90.7)	0.1 (0.6)	0.9

Intermediate precision	Mean mmol/L (mg/dL urea)	SD mmol/L (mg/dL urea)	CV %
Precinorm U	6.66 (40.0)	0.08 (0.5)	1.2
Precipath U	23.2 (139)	0.3 (2)	1.1
Human serum 3	9.13 (54.8)	0.10 (0.6)	1.1
Human serum 4	14.9 (89.5)	0.2 (1.2)	1.3

Urine

Repeatability	Mean mmol/L (mg/dL urea)	SD mmol/L (mg/dL urea)	CV %
Control level 1	161 (967)	4 (24)	2.2
Control level 2	288 (1730)	3 (18)	1.2
Human urine 1	324 (1946)	4 (24)	1.3
Human urine 2	137 (823)	3 (18)	1.9

Intermediate precision	Mean mmol/L (mg/dL urea)	SD mmol/L (mg/dL urea)	CV %
Control level 1	154 (925)	4 (24)	2.7
Control level 2	280 (1682)	6 (36)	2.3
Human urine 3	316 (1898)	6 (36)	2.0
Human urine 4	133 (799)	3 (18)	2.4

The data obtained on **cobas c 501** analyzer(s) are representative for **cobas c 311** analyzer(s).

Method comparison

Urea values for human serum, plasma and urine samples obtained on a **cobas c 501** analyzer (y) were compared with those determined on Roche/Hitachi 917/MODULAR P analyzers (x), using the corresponding Roche/Hitachi reagent.

Serum/plasma

Sample size (n) = 175

Passing/Bablok ¹⁵	Linear regression
$y = 0.990x + 0.138$ mmol/L	$y = 0.976x + 0.303$ mmol/L
$\tau = 0.959$	$r = 0.998$

The sample concentrations were between 2.27 and 39.4 mmol/L (13.6 and 237 mg/dL urea).

Urine

Sample size (n) = 267

Passing/Bablok ¹⁵	Linear regression
$y = 1.006x - 6.50$ mmol/L	$y = 1.035x - 14.1$ mmol/L
$\tau = 0.949$	$r = 0.998$

The sample concentrations were between 39.0 and 1314 mmol/L (234 and 7892 mg/dL urea).

The data obtained on **cobas c 501** analyzer(s) are representative for **cobas c 311** analyzer(s).

COBAS INTEGRA systems

System information

UREL: Test ID 0-003 for serum and plasma

URELU: Test ID 0-004 for urine

Reagents - working solutions

R1 TRIS buffer: 220 mmol/L; 2-oxoglutarate: 73 mmol/L; NADH: 2.5 mmol/L; ADP: 6.5 mmol/L; urease (jack bean): ≥ 300 μ kat/L; GLDH (bovine): ≥ 80 μ kat/L; stabilizers; pH 8.6

R1 is in position B.

Storage and stability

Shelf life at 2-8 °C	See expiration date on cobas c pack label
On-board in use at 10-15 °C	8 weeks

Application for serum, plasma and urine

Test definition

Measuring mode	Absorbance
Abs. calculation mode	Kinetic
Reaction direction	Decrease
Wavelength A/B	340/409 nm
Calc. first/last	23/28
Unit	mmol/L

Serum, plasma

Reaction mode	R-S	
<i>Urine</i>		
Reaction mode	D-R-S	
Predilution factor	50	
Pipetting parameters		
<i>Serum/plasma/urine</i>	Diluent (H ₂ O)	
R1	50 µL	95 µL
Sample	2 µL	98 µL
Total volume	245 µL	
Calibration		
Calibrator	Calibrator f.a.s. Use deionized water as zero calibrator.	
Calibration mode	Linear regression	
Calibration replicate	Duplicate recommended	
Calibration interval	Each cobas c pack, every 4 weeks, and as required following quality control procedures	

Calibration interval may be extended based on acceptable verification of calibration by the laboratory.

Traceability: This method has been standardized against ID/MS.

Quality control

Quality control serum, plasma	Precinorm U plus or PreciControl ClinChem Multi 1 Precipath U plus or PreciControl ClinChem Multi 2
Quality control urine	Quantitative urine controls are recommended for routine quality control.
Control interval	24 hours recommended
Control sequence	User defined
Control after calibration	Recommended

For quality control, use control materials as listed in the "Order information" section. In addition, other suitable control material can be used.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

Limitations - interference

In very rare cases, gammopathy, in particular type IgM (Waldenström's macroglobulinemia), may cause unreliable results.¹¹

Serum/plasma

Criterion: Recovery within ± 0.83 mmol/L of initial values of samples ≤ 8.3 mmol/L and within $\pm 10\%$ for samples > 8.3 mmol/L.

Icterus:¹² No significant interference up to an I index of 60 for conjugated and unconjugated bilirubin (approximate conjugated and unconjugated bilirubin concentration: 1026 µmol/L or 60 mg/dL).

Hemolysis:¹² No significant interference up to an H index of 1000 (approximate hemoglobin concentration: 621 µmol/L or 1000 mg/dL). Hemolytic specimens may cause high absorbance flagging. Choose diluted sample treatment for automatic rerun.

Lipemia (Intralipid):¹² No significant interference up to an L index of 2000. There is poor correlation between the L index (corresponds to turbidity) and triglycerides concentration. Lipemic specimens may cause high absorbance flagging.

Anticoagulants: Do not use ammonium heparin as an anticoagulant.

Drugs: No interference was found at therapeutic concentrations using common drug panels.^{13,14}

Ammonium ions may cause erroneously elevated results.

Urine

Criterion: Recovery within ± 15 mmol/L of initial values of samples ≤ 150 mmol/L and within $\pm 10\%$ for samples > 150 mmol/L

Drugs: No interference was found at therapeutic concentrations using common drug panels.¹⁴

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

ACTION REQUIRED

Special Wash Programming: The use of special wash steps is mandatory when certain test combinations are run together on COBAS INTEGRA analyzers. Refer to the CLEAN Method Sheet for further instructions and for the latest version of the Extra wash cycle list.

Where required, special wash/carry-over evasion programming must be implemented prior to reporting results with this test.

Limits and ranges**Measuring range***Serum/plasma*

0.5-40 mmol/L (3.0-240 mg/dL urea, 1.4-112 mg/dL urea nitrogen)

Determine samples having higher concentrations via the rerun function. Dilution of samples via the rerun function is a 1:10 dilution. Results from samples diluted using the rerun function are automatically multiplied by a factor of 10.

Urine

1.0-2000 mmol/L (0.006-12 g/dL urea, 2.8-5600 mg/dL urea nitrogen)

Determine samples having higher concentrations via the rerun function. Dilution of samples via the rerun function is a 1:3 dilution. Results from samples diluted using the rerun function are automatically multiplied by a factor of 3.

Determine samples having lower concentrations via the rerun function. For samples with concentrations lower than 40 mmol/L, the rerun function reduces the sample predilution factor to 2 (final dilution 1 + 1). The results are automatically multiplied by the reduced predilution factor.

Lower limits of measurement*Serum/plasma*

Lower detection limit of the test:

0.5 mmol/L (3.0 mg/dL urea, 1.4 mg/dL urea nitrogen)

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying 3 standard deviations above that of a zero sample (zero sample + 3 SD, repeatability, n = 30).

Urine

Lower detection limit of the test:

1.0 mmol/L (0.006 g/dL urea, 2.8 mg/dL urea nitrogen)

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying 3 standard deviations above that of a zero sample (zero sample + 3 SD, repeatability, n = 30).

Specific performance data

Representative performance data on the analyzers are given below. Results obtained in individual laboratories may differ.

Precision

Precision was determined using human samples and controls in an internal protocol with repeatability and intermediate precision (2 aliquots per run, 2 runs per day, 20 days). The following results were obtained on the COBAS INTEGRA 700 analyzer:

Serum/plasma

	Level 1	Level 2
Mean	4.08 mmol/L (24.6 mg/dL)	31.0 mmol/L (186 mg/dL)

	Level 1	Level 2
CV repeatability	2.3 %	0.9 %
CV intermediate precision	3.9 %	2.8 %

Urine

	Level 1	Level 2
Mean	421 mmol/L (2.53 g/dL)	679 mmol/L (4.08 g/dL)
CV repeatability	1.3 %	1.2 %
CV intermediate precision	1.8 %	1.8 %

The data obtained on COBAS INTEGRA 700 analyzer(s) are representative for COBAS INTEGRA 400 analyzer(s).

Method comparison**Serum/plasma**

Urea values for human serum and plasma samples obtained on a COBAS INTEGRA 700 analyzer using the COBAS INTEGRA Urea/BUN reagent (y) were compared with those determined using commercially available reagents for urea on a COBAS INTEGRA 700 analyzer (x) and an alternative manufacturer's clinical chemistry system (x). Samples were measured in duplicate. Sample size (n) represents all replicates.

COBAS INTEGRA 700 analyzer

Sample size	(n)	236
Corr. coefficient	(r)	0.999
	(r _s)	0.999

Linear regression $y = 1.004x + 0.071$ mmol/L

Passing/Bablok¹⁵ $y = 1.001x + 0.014$ mmol/L

Alternative system

Sample size	(n)	236
Corr. coefficient	(r)	0.999
	(r _s)	0.999

Linear regression $y = 0.983x + 0.176$ mmol/L

Passing/Bablok¹⁵ $y = 0.995x + 0.041$ mmol/L

The sample concentrations were between 1.1 and 38.1 mmol/L (6.61 and 229 mg/dL).

Urine

Urea values for human urine samples obtained on a COBAS INTEGRA 700 analyzer using the COBAS INTEGRA Urea/BUN reagent (y) were compared with those determined using commercially available reagents for urea on a COBAS INTEGRA 700 analyzer (x).

COBAS INTEGRA 700 analyzer

Sample size	(n)	120
Corr. coefficient	(r)	0.999
	(r _s)	0.998

Linear regression $y = 1.000x + 1.30$ mmol/L

Passing/Bablok¹⁵ $y = 0.999x + 3.47$ mmol/L

The sample concentrations were between 56.6 and 796 mmol/L (0.340 and 4.78 g/dL).

The data obtained on COBAS INTEGRA 700 analyzer(s) are representative for COBAS INTEGRA 400 analyzer(s).

References

- Lamb EJ, Jones GRD. Kidney Function Tests. In: Rifai N, Chiu RWK, Young I, Burnham CAD, Wittwer CT, editors. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. 7th ed. 2023; p. 352e1-352e59.


- Vanholder R, Gryp T, Glorieux G. Urea and chronic kidney disease: the comeback of the century? (in uraemia research). Nephrol Dial Transplant. 2018 Jan 1;33(1):4-12. doi: 10.1093/ndt/gfx039.
- Rosenberg WMC, Badrick T, Lo SF, et al. Liver disease. In: Rifai N, Chiu RWK, Young I, Burnham CAD, Wittwer CT, editors. Tietz Textbook of Laboratory Medicine, Saunders Elsevier, Philadelphia, 7th edition, 2023, chapter 51, p. 701-763.e21.
- Richterich R, Colombo JP. Klinische Chemie. 4th ed. Basel: Karger S 1978:319-324.
- Talke H, Schubert GA. Enzymatische Harnstoffbestimmung in Blut und Serum im optischen Test nach Warburg. Klin Wochenschr 1965;43:174.
- Tiffany TO, Jansen JM, Burtis CA, et al. Enzymatic kinetic rate and end-point analyses of substrate, by use of a GeMSAEC Fast Analyzer. Clin Chem 1972;18:829-840.
- Sampson EJ, Baired MA, Burtis CA, et al. A coupled-enzyme equilibrium method for measuring urea in serum: Optimization and evaluation of the AACC study group on urea candidate reference method. Clin Chem 1980;26:816-826.
- WHO Publication: Use of anticoagulants in diagnostic laboratory investigations, WHO/DIL/LAB/99.1 Rev.2:Jan 2002.
- Löhr B, El-Samalouti V, Junge W, et al. Reference Range Study for Various Parameters on Roche Clinical Chemistry Analyzers. Clin Lab 2009;55:465-471.
- Wu AHB, ed. Tietz Clinical Guide to Laboratory Tests, 4th edition. St. Louis (MO): Saunders Elsevier 2006;1096.
- Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007;45(9):1240-1243.
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-475.
- Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". Eur J Clin Chem Clin Biochem 1996;34:385-386.
- Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.
- Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard:

	Contents of kit
	Volume for reconstitution
	Global Trade Item Number

Rx only For USA: Caution: Federal law restricts this device to sale by or on the order of a physician.

COBAS, NAVIFY, PRECINORM, PRECIPATH and PRECICONTROL are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2024, Roche Diagnostics

04460715500V16.0

UREAL

Urea/BUN

CE 0123



Roche Diagnostics GmbH
Sandhofer Strasse 116
68305 Mannheim, Germany
www.roche.com

+800 5505 6606



cobas®