


REF		CONTENT		Analyzer(s) on which <b>cobas c</b> pack(s) can be used
08105529190	08105529500	Tina-quant $\beta$ 2-Microglobulin (150 tests)	System-ID 2025 001	<b>cobas c 303, cobas c 503, cobas c 703</b>

Materials required (but not provided):

08047545190	Calibrator $\beta$ 2-Microglobulin (2 x 1 mL)	Code 20474	
08362785190	Control Set $\beta$ 2-Microglobulin - Level I (2 x 1 mL) - Level II (2 x 1 mL)	Code 20144 Code 20145	
08063494190	Diluent NaCl 9 % (123 mL)	System-ID 2906 001	

**English****System information****B2MG:** ACN 20250**B2MGU:** ACN 20251**Intended use**

Immunoturbidimetric assay for the quantitative in vitro determination of  $\beta$ 2-microglobulin (B2MG) in human serum, plasma and urine on **cobas c** systems.

**Summary**

$\beta$ 2-microglobulin (B2MG) measurement, performed with this assay in human serum, plasma and urine can be used to aid in diagnosis of kidney diseases (such as glomerulopathies, tubulopathies, renal failure, and kidney-related amyloidosis), acute kidney allograft rejections, active rheumatoid arthritis and hematological diseases (including malignancies such as lymphoma and myeloma).

B2MG is a small globular peptide made up of 100 amino acid residues and with a molecular weight of 11800 D. It non-covalently associates with the major histocompatibility complex-I (MHC-I)/human leukocyte antigen-I (HLA-I) on the surfaces of all nucleated cells (exception: trophoblasts) and it is present in most biological fluids including serum, urine and synovial fluid.<sup>1,2</sup> Lymphocytes are the major location of synthesis. B2MG is constantly released into the blood in small quantities. B2MG is normally cleared by glomerular filtration in the kidneys. Thereafter, up to 99.9 % is reabsorbed by the proximal tubules.<sup>1</sup> The appearance of notable quantities of this protein in the urine reflects failure of the tubular resorptive mechanisms. Levels of serum and urinary B2MG reflect the functions of glomeruli and proximal tubules.<sup>3</sup>

Patients with renal dysfunction show high plasma concentrations of B2MG due to the impaired kidneys excretion.<sup>2</sup> Patients with acute kidney injury have a significantly higher serum B2MG concentration compared to healthy people.<sup>4</sup> In addition, serum B2MG levels are known to increase in renal diseases which disrupt glomerular filtration and tubular reabsorption.<sup>1</sup> B2MG clearance during dialysis treatment is clinically significant since B2MG can deposit as amyloid, leading to systemic amyloidosis in many patients on long-term dialysis and resulting in injury of many tissues, including kidneys, brain, heart and peripheral nerves, depending on the site of amyloid deposition.<sup>2,5</sup>

In urine, increased levels of B2MG are usually indicative of renal tubular injury.<sup>3,6</sup> In healthy patients only small amounts of B2MG can be detected. Drug nephrotoxicity causing severe, acute changes in tubular reabsorption and progressive renal diseases causing irreversible structural tubular defects impair tubular reabsorption of numerous smaller proteins including B2MG. Therefore, urinary B2MG is discussed as a marker for the diagnosis and monitoring of tubulointerstitial renal damage.<sup>1,3</sup> Excretion of B2MG is increased in Fanconi syndrome, a generalized dysfunction of the proximal tubules.<sup>7</sup> Causes for acquired Fanconi syndrome include exposure to toxins and drugs.<sup>8</sup> Elevated B2MG values may identify patients at higher risk of glomerular filtration rate (GFR) decline in other kidney diseases such as membranous nephropathy.<sup>9</sup> Furthermore, there is evidence that B2MG excretion is associated with acute allograft rejection in renal transplant recipients.<sup>10</sup>

B2MG serum concentration reflects the activation of the cellular immune system.<sup>11</sup> Increased serum B2MG levels have been shown in autoimmune diseases such as rheumatoid arthritis.<sup>12,13,14,15</sup>

In oncology increased serum B2MG levels are observed in patients with several malignancies, including lymphoma and multiple myeloma. In

multiple myeloma, measurement of serum B2MG levels is considered essential for staging and recommended by the EHA-ESMO Clinical Practice Guidelines.<sup>16,17,18</sup> B2MG is also used as a prognostic factor for lymphoma malignancies and recommended as part of the initial workup in the ESMO guidelines (2013).<sup>17,19,20,21,22,23,24,25,26,27,28</sup> Furthermore, a prognostic index was created for patients with large B-cell lymphoma, comprising 4 factors: age, performance status, stage and B2MG. This was suggested as a useful risk model for categorising lymphoma patients into risk groups and predicting the various outcomes.<sup>28</sup>

Various assay methods are available for B2MG determination, such as radioimmunoassays (RIA), enzyme-linked immunosorbent assays (ELISA), nephelometric immunoassays, and turbidimetric methods.<sup>29</sup> The Roche B2MG assay is based on the principle of immunological agglutination with latex reaction enhancement.

**Test principle**

Immunoturbidimetric assay.

Latex-bound anti- $\beta$ 2-microglobulin antibodies react with antigen from the sample to form antigen/antibody complexes which are determined turbidimetrically after agglutination.<sup>30</sup>

**Reagents - working solutions**

- R1** TRIS/HCl buffer: 23 g/L, pH 8.2; NaCl: 19 g/L; EDTA: 1.3 g/L; preservative
- R3** Latex particles coated with polyclonal anti-human  $\beta$ 2-microglobulin antibody (rabbit); preservative

R1 is in position B and R3 is in position C.

**Precautions and warnings**

For in vitro diagnostic use for health care professionals. Exercise the normal precautions required for handling all laboratory reagents.

Infectious or microbial waste:

Warning: handle waste as potentially biohazardous material. Dispose of waste according to accepted laboratory instructions and procedures.

Environmental hazards:

Apply all relevant local disposal regulations to determine the safe disposal.

Safety data sheet available for professional user on request.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:



Warning

- H317 May cause an allergic skin reaction.
- H412 Harmful to aquatic life with long lasting effects.

**Prevention:**

- P261 Avoid breathing mist or vapours.
- P273 Avoid release to the environment.

**Tina-quant  $\beta$ 2-Microglobulin**

P280 Wear protective gloves.

**Response:**

P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.

P362 + P364 Take off contaminated clothing and wash it before reuse.

**Disposal:**

P501 Dispose of contents/container to an approved waste disposal plant.

Product safety labeling follows EU GHS guidance.

Contact phone: all countries: +49-621-7590

**Reagent handling**

Ready for use

Carefully invert reagent container several times prior to use to ensure that the reagent components are mixed.

**Storage and stability**

Shelf life at 2-8 °C: See expiration date on **cobas c** pack label.

On-board in use and refrigerated on the analyzer: 12 weeks

**Specimen collection and preparation**

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum

Plasma: Li-heparin, K<sub>2</sub>- and K<sub>3</sub>-EDTA plasma

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Urine

Each urine sample must be centrifuged (10 minutes at approximately 3000 × g) prior to testing.<sup>31</sup>

B2MG is unstable in acidic conditions and degradation occurs at pH < 6 within 2 hours.<sup>29,7</sup> Thus, pre-analytical conditions are highly important. Since the degradation also takes place in the bladder, collection of a spot urine sample should not be performed in the morning due to a lower urine pH.<sup>29</sup>

A strict control of the urine pH after collection is required: urine samples must be adjusted to pH 7-9 by the addition of 1 N NaOH as soon as possible after receipt.<sup>7</sup>

Stability in serum and Li-heparin, K<sub>2</sub>- and K<sub>3</sub>-EDTA plasma: 3 days at 15-25 °C  
6 months at -20 °C (± 5 °C)  
Freezing and thawing up to 2 times is allowed.

Stability in urine: 5 days at 15-25 °C  
14 days at 2-8 °C  
12 weeks at -20 °C (± 5 °C)  
Freezing and thawing up to 2 times is allowed.

Centrifuge samples containing precipitates before performing the assay. See the limitations and interferences section for details about possible sample interferences.

**Materials provided**

See "Reagents – working solutions" section for reagents.

**Materials required (but not provided)**

See "Order information" section

General laboratory equipment

**Assay**

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

The performance of applications not validated by Roche is not warranted and must be defined by the user.

**Application for serum, plasma and urine****Test definition**

Reporting time	10 min		
Wavelength (sub/main)	-/700 nm		
Reagent pipetting		Diluent (H <sub>2</sub> O)	
R1	63 $\mu$ L	-	
R3	63 $\mu$ L	-	
<i>Sample volumes</i>	<i>Sample</i>	<i>Sample dilution</i>	
		<i>Sample</i>	<i>Diluent (NaCl)</i>
Normal	1 $\mu$ L	-	-
Decreased	1 $\mu$ L	10 $\mu$ L	100 $\mu$ L
Increased	1 $\mu$ L	-	-

For further information about the assay test definitions refer to the application parameters setting screen of the corresponding analyzer and assay.

**Calibration***Application for serum/plasma (ACN 20250)*

Calibrators S1: H<sub>2</sub>O  
S2: Calibrator  $\beta$ 2-Microglobulin

Calibration mode Linear

Calibration frequency Full calibration  
- after reagent lot change  
- as required following quality control procedures

*Application for urine (ACN 20251)*

Transfer of calibration from serum/plasma application (ACN 20250)

Calibration interval may be extended based on acceptable verification of calibration by the laboratory.

Traceability: This method has been standardized against the WHO standard.

**Quality control**

For quality control, use control materials as listed in the "Order information" section. In addition, other suitable control material can be used.

The control intervals and limits should be adapted to each laboratory's individual requirements. It is recommended to perform quality control always after lot calibration and subsequently at least every 12 weeks. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

**Calculation**

**cobas c** systems automatically calculate the analyte concentration of each sample in the unit mg/L (nmol/L,  $\mu$ g/L).

Conversion factors: mg/L  $\times$  84.7 = nmol/L  
mg/L  $\times$  1000 =  $\mu$ g/L

**Limitations - interference***Serum/plasma*

Criterion: Recovery within  $\pm 0.22$  mg/L of initial values at a  $\beta$ 2-microglobulin concentration of  $\leq 2.2$  mg/L and within  $\pm 10\%$  for samples  $> 2.2$  mg/L.

Icterus:<sup>32</sup> No significant interference up to an I index of 60 for conjugated and unconjugated bilirubin (approximate conjugated and unconjugated bilirubin concentration: 1026  $\mu$ mol/L or 60 mg/dL).

Hemolysis:<sup>32</sup> No significant interference up to an H index of 1000 (approximate hemoglobin concentration: 621  $\mu$ mol/L or 1000 mg/dL).

Lipemia (Intralipid):<sup>32</sup> No significant interference up to an L index of 1000. There is poor correlation between the L index (corresponds to turbidity) and triglycerides concentration.

High-dose hook effect: No false result occurs up to a B2MG concentration of 240 mg/L.

Drugs: No interference was found at therapeutic concentrations using common drug panels.<sup>33,34</sup>

Rheumatoid factors: No significant interference from rheumatoid factors up to a concentration of 1200 IU/mL.

In very rare cases, gammopathy, in particular type IgM (Waldenström's macroglobulinemia), may cause unreliable results.<sup>35</sup>

*Urine*

Criterion: Recovery within  $\pm 0.1$  mg/L of initial values at a  $\beta$ 2-microglobulin concentration of  $\leq 1.0$  mg/L and within  $\pm 10\%$  for samples  $> 1.0$  mg/L.

Hemolysis: No significant interference up to an H index of 1100 (approximate hemoglobin concentration: 1100 mg/dL).<sup>32</sup>

High-dose hook effect: No false result occurs up to a B2MG concentration of 240 mg/L.

Drugs: No interference was found at therapeutic concentrations using common drug panels.<sup>34</sup>

In very rare cases, gammopathy, in particular type IgM (Waldenström's macroglobulinemia), may cause unreliable results.<sup>35</sup>

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

**ACTION REQUIRED**

**Special Wash Programming:** The use of special wash steps is mandatory when certain test combinations are run together on **cobas c** systems. All special wash programming necessary for avoiding carry-over is available via the **cobas** link. The latest version of the carry-over evasion list can be found with the NaOHD/SMS/SCCS Method Sheet. For further instructions, refer to the operator's manual.

**Limits and ranges****Measuring range***Serum/plasma*

0.2-8.0 mg/L (16.9-678 nmol/L)

*Urine*

0.2-5.8 mg/L (16.9-491 nmol/L)

Determine samples having higher concentrations via the rerun function. Dilution of samples via the rerun function is a 1:11 dilution. Results from samples diluted using the rerun function are automatically multiplied by a factor of 11.

**Lower limits of measurement***Limit of Blank, Limit of Detection and Limit of Quantitation*

Limit of Blank = 0.1 mg/L (8.5 nmol/L)

Limit of Detection = 0.15 mg/L (12.7 nmol/L)

Limit of Quantitation = 0.2 mg/L (16.9 nmol/L)

The Limit of Blank, Limit of Detection and Limit of Quantitation were determined in accordance with the CLSI (Clinical and Laboratory Standards Institute) EP17-A2 requirements.

The Limit of Blank is the 95<sup>th</sup> percentile value from  $n \geq 60$  measurements of analyte-free samples over several independent series. The Limit of Blank corresponds to the concentration below which analyte-free samples are found with a probability of 95 %.

The Limit of Detection is determined based on the Limit of Blank and the standard deviation of low concentration samples.

The Limit of Detection corresponds to the lowest analyte concentration which can be detected (value above the Limit of Blank with a probability of 95 %).

The Limit of Quantitation is the lowest analyte concentration that can be reproducibly measured with a total error of 20 %. It has been determined using low concentration human  $\beta$ 2-microglobulin samples.

**Expected values***Serum/plasma*

< 60 years 0.8-2.4 mg/L 68-203 nmol/L<sup>29</sup>

> 60 years  $\leq 3.0$  mg/L 254 nmol/L<sup>29</sup>

*Urine*

male:  $\leq 0.300$  mg/L (25.4 nmol/L)<sup>36</sup>

female:  $\leq 0.183$  mg/L (15.5 nmol/L)<sup>36</sup>

B2MG/creatinine (urine):  $\leq 2.46$  nmol/mmol<sup>36</sup>

24 h urine: 2.80-30.7 nmol<sup>7</sup>

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

**Specific performance data**

Representative performance data on the analyzers are given below. These data represent the performance of the analytical procedure itself.

Results obtained in individual laboratories may differ due to heterogenous sample materials, aging of analyzer components and mixture of reagents running on the analyzer.

**Precision**

Precision was determined using human samples and controls in accordance with the CLSI (Clinical and Laboratory Standards Institute) EP05-A3 requirements with repeatability ( $n = 84$ ) and intermediate precision (2 aliquots per run, 2 runs per day, 21 days). Results for repeatability and intermediate precision were obtained on the **cobas c** 503 analyzer.

*Serum/plasma*

<i>Repeatability</i>	<i>Mean</i>	<i>SD</i>	<i>CV</i>
	<i>mg/L</i>	<i>mg/L</i>	<i>%</i>
CS B2MG1 <sup>a)</sup>	2.53	0.0405	1.6
CS B2MG2 <sup>b)</sup>	4.53	0.0303	0.7
Human serum 1	0.243	0.0247	10.2
Human serum 2	1.74	0.0273	1.6
Human serum 3	2.80	0.0191	0.7
Human serum 4	3.94	0.0393	1.0
Human serum 5	7.06	0.0300	0.4

*Intermediate precision*

	<i>Mean</i>	<i>SD</i>	<i>CV</i>
	<i>mg/L</i>	<i>mg/L</i>	<i>%</i>
CS B2MG1 <sup>a)</sup>	2.53	0.0459	1.8
CS B2MG2 <sup>b)</sup>	4.51	0.0538	1.2
Human serum 1	0.243	0.0265	10.9
Human serum 2	1.74	0.0318	1.8
Human serum 3	2.81	0.0251	0.9
Human serum 4	3.94	0.0410	1.0
Human serum 5	7.06	0.0402	0.6

a) Control Set  $\beta$ 2-Microglobulin, Level I

b) Control Set  $\beta$ 2-Microglobulin, Level II

*Urine*

<i>Repeatability</i>	<i>Mean</i>	<i>SD</i>	<i>CV</i>
	<i>mg/L</i>	<i>mg/L</i>	<i>%</i>

**Tina-quant β2-Microglobulin**

CS B2MG1 <sup>a)</sup>	2.54	0.0295	1.2
CS B2MG2 <sup>b)</sup>	4.57	0.0355	0.8
Human urine 1	0.274	0.0373	13.6
Human urine 2	0.395	0.0385	9.8
Human urine 3	0.807	0.0457	5.7
Human urine 4	2.53	0.0425	1.7
Human urine 5	5.07	0.0389	0.8

<i>Intermediate precision</i>	<i>Mean</i>	<i>SD</i>	<i>CV</i>
	<i>mg/L</i>	<i>mg/L</i>	<i>%</i>
CS B2MG1 <sup>c)</sup>	2.54	0.0385	1.5
CS B2MG2 <sup>d)</sup>	4.57	0.0692	1.5
Human urine 1	0.274	0.0385	14.0
Human urine 2	0.395	0.0424	10.7
Human urine 3	0.807	0.0473	5.9
Human urine 4	2.53	0.0449	1.8
Human urine 5	5.06	0.0458	0.9

c) Control Set β2-Microglobulin Level I

d) Control Set β2-Microglobulin Level II

The data obtained on **cobas c** 503 analyzer(s) are representative for **cobas c** 303 analyzer(s) and **cobas c** 703 analyzer(s).

**Method comparison***Serum/plasma*

β2-Microglobulin values for human serum samples obtained on a **cobas c** 503 analyzer (y) were compared with those determined using the corresponding reagent on a **cobas c** 501 analyzer (x).

Sample size (n) = 64

Passing/Bablok <sup>37</sup>	Linear regression
$y = 1.000x + 0.0150 \text{ mg/L}$	$y = 1.001x + 0.00921 \text{ mg/L}$
$\tau = 0.991$	$r = 1.000$

The sample concentrations were between 0.320 and 7.82 mg/L.

β2-Microglobulin values for human serum samples obtained on a **cobas c** 303 analyzer (y) were compared with those determined using the corresponding reagent on a **cobas c** 501 analyzer (x).

Sample size (n) = 67

Passing/Bablok <sup>37</sup>	Linear regression
$y = 1.009x + 0.00991 \text{ mg/L}$	$y = 1.004x + 0.0380 \text{ mg/L}$
$\tau = 0.976$	$r = 0.999$

The sample concentrations were between 0.240 and 7.78 mg/L.

β2-Microglobulin values for human serum samples obtained on a **cobas c** 703 analyzer (y) were compared with those determined using the corresponding reagent on a **cobas c** 503 analyzer (x).

Sample size (n) = 71

Passing/Bablok <sup>37</sup>	Linear regression
$y = 1.019x - 0.0612 \text{ mg/L}$	$y = 1.015x - 0.0528 \text{ mg/L}$
$\tau = 0.981$	$r = 0.999$

The sample concentrations were between 0.484 and 7.65 mg/L.

*Urine*

β2-Microglobulin values for human urine samples obtained on a **cobas c** 503 analyzer (y) were compared with those determined using the corresponding reagent on a **cobas c** 501 analyzer (x).

Sample size (n) = 68

Passing/Bablok <sup>37</sup>	Linear regression
$y = 1.000x + 0.00300 \text{ mg/L}$	$y = 0.998x + 0.00986 \text{ mg/L}$

 $\tau = 0.983$   $r = 1.000$ 

The sample concentrations were between 0.260 and 5.52 mg/L.

β2-Microglobulin values for human urine samples obtained on a **cobas c** 303 analyzer (y) were compared with those determined using the corresponding reagent on a **cobas c** 501 analyzer (x).

Sample size (n) = 76

Passing/Bablok <sup>37</sup>	Linear regression
$y = 1.021x - 0.00731 \text{ mg/L}$	$y = 1.023x - 0.0105 \text{ mg/L}$
$\tau = 0.968$	$r = 0.999$

The sample concentrations were between 0.200 and 5.39 mg/L.

β2-Microglobulin values for human urine samples obtained on a **cobas c** 703 analyzer (y) were compared with those determined using the corresponding reagent on a **cobas c** 503 analyzer (x).

Sample size (n) = 73

Passing/Bablok <sup>37</sup>	Linear regression
$y = 0.966x + 0.0129 \text{ mg/L}$	$y = 0.969x + 0.00734 \text{ mg/L}$
$\tau = 0.975$	$r = 0.999$

The sample concentrations were between 0.221 and 5.74 mg/L.

**References**

- Argyropoulos CP, Chen SS, Ng YH, et al. Rediscovering Beta-2 Microglobulin As a Biomarker across the Spectrum of Kidney Diseases. *Front Med (Lausanne)*. 2017 Jun 15;4:73. doi: 10.3389/fmed.2017.00073.
- Hortin GL. Amino Acids, Peptides and Proteins. In: Burtis C, Ashwood E, Bruns D, editors. *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics*. Saunders Elsevier, Philadelphia, 5th edition, 2012, chapter 21, p. 509-564.
- Zeng X, Hossain D, Bostwick DG, et al. Urinary β2-Microglobulin Is a Good Indicator of Proximal Tubule Injury: A Correlative Study with Renal Biopsies. *J Biomark* 2014;2014:492838. doi: 10.1155/2014/492838.
- Barton KT, Kakajiwala A, Dietzen DJ, et al. Using the newer Kidney Disease: Improving Global Outcomes criteria, beta-2-microglobulin levels associate with severity of acute kidney injury. *Clin Kidney J* 2018 Dec;11(6):797-802. doi: 10.1093/ckj/sfy056.
- Mumtaz A, Anees M, Bilal M, et al. Beta-2 microglobulin levels in hemodialysis patients. *Saudi J Kidney Dis Transpl* 2010 Jul;21(4):701-706
- Ikeda A, Konta T, Takasaki S, et al. In a non-diabetic Japanese population, the combination of macroalbuminuria and increased urine beta 2-microglobulin predicts a decline of renal function: the Takahata study. *Nephrol Dial Transplant*. 2009 Mar;24(3):841-847. doi: 10.1093/ndt/gfn591.
- Schardijn G, Statius van Eps LW, Swaak AJ, et al. Urinary beta2-microglobulin in upper and lower urinary-tract infections. *The Lancet* 1979 Apr;1(8120):805-807.
- Solano A, Lew SQ, Ing TS. Dent-Wrong disease and other rare causes of the Fanconi syndrome. *Clin Kidney J* 2014 Aug;7(4):344-347.
- Branten AJ, du Buf-Vereijken PW, Klasen IS, et al. Urinary excretion of beta2-microglobulin and IgG predict prognosis in idiopathic membranous nephropathy: a validation study. *J Am Soc Nephrol* 2005 Jan;16(1):169-174. doi: 10.1681/ASN.2004040287.
- Oetting WS, Rogers TB, Krick TP, et al. Urinary beta2-microglobulin is associated with acute renal allograft rejection. *Am J Kidney Dis* 2006 May;47(5):898-904.
- Li L, Dong M, Wang XG. The Implication and Significance of Beta 2 Microglobulin: A Conservative Multifunctional Regulator. *Chin Med J (Engl)*. 2016 Feb 20;129(4):448-455. doi: 10.4103/0366-6999.176084.
- Gottenberg JE, Miceli-Richard C, Ducot B, et al. Markers of B-lymphocyte activation are elevated in patients with early rheumatoid arthritis and correlated with disease activity in the ESPOIR cohort. *Arthritis Res Ther*. 2009;11(4):R114. doi: 10.1186/ar2773. Epub 2009 Jul 23.

- 13 Talal N, Grey HM, Zvaifler N, et al. Elevated salivary and synovial fluid beta2-microglobulin in Sjogren's syndrome and rheumatoid arthritis. *Science*. 1975 Mar 28;187(4182):1196-1198. doi: 10.1126/science.46621.
- 14 Aho K, Heliövaara M, Knekt P, et al. Serum immunoglobulins and the risk of rheumatoid arthritis. *Ann Rheum Dis*. 1997 Jun;56(6):351-356. doi: 10.1136/ard.56.6.351.
- 15 Gottenberg JE, Aucouturier F, Goetz J, et al. Serum immunoglobulin free light chain assessment in rheumatoid arthritis and primary Sjogren's syndrome. *Ann Rheum Dis*. 2007 Jan;66(1):23-27. doi: 10.1136/ard.2006.052159.
- 16 Yoo C, Yoon DH, Suh C. Serum beta-2 microglobulin in malignant lymphomas: an old but powerful prognostic factor. *Blood Res*. 2014 Sep;49(3):148-153. doi: 10.5045/br.2014.49.3.148.
- 17 Rossi D, Fangazio M, De Paoli L, et al. Beta-2-microglobulin is an independent predictor of progression in asymptomatic multiple myeloma. *Cancer*. 2010 May 1;116(9):2188-2200. doi: 10.1002/cncr.24959.
- 18 Dimopoulos MA, Moreau P, Terpos E, et al. Multiple myeloma: EHA-ESMO Clinical Practice Guidelines for diagnosis, treatment and follow-up†. *Ann Oncol* 2021 Mar;32(3):309-322.
- 19 Cánovas A, Alonso JJ, Barreiro G, et al. Prognostic factors in follicular lymphoma: the importance of beta-2 microglobulin. *Tumori*. 2010 Jan-Feb;96(1):117-121. doi: 10.1177/030089161009600119.
- 20 Yoo C, Yoon DH, Jo JC, et al. Prognostic impact of beta-2 microglobulin in patients with extranodal natural killer/T cell lymphoma. *Ann Hematol*. 2014 Jun;93(6):995-1000. doi: 10.1007/s00277-014-2015-2.
- 21 Wu L, Wang T, Gui W, et al. Prognostic significance of serum beta-2 microglobulin in patients with non-Hodgkin lymphoma. *Oncology*. 2014;87(1):40-47. doi: 10.1159/000362670.
- 22 Nakajima Y, Tomita N, Watanabe R, et al. Prognostic significance of serum beta-2 microglobulin level in Hodgkin lymphoma treated with ABVD-based therapy. *Med Oncol*. 2014 Sep;31(9):185. doi: 10.1007/s12032-014-0185-3.
- 23 Atesoglu EB, Hacıhanefioglu A, Gulbas Z. Beta-2 microglobulin predicts the outcome after autologous stem cell transplantation in non-Hodgkin lymphoma. *Transfus Apher Sci*. 2015 Feb;52(1):65-71. doi: 10.1016/j.transci.2014.12.007.
- 24 Yoo C, Yoon DH, Kim S, et al. Serum beta-2 microglobulin as a prognostic biomarker in patients with mantle cell lymphoma. *Hematol Oncol*. 2016 Mar;34(1):22-27. doi: 10.1002/hon.2188.
- 25 Wang Q, Qin Y, Zhou S, et al. Prognostic value of pretreatment serum beta-2 microglobulin level in advanced classical Hodgkin lymphoma treated in the modern era. *Oncotarget*. 2016 Nov 1;7(44):72219-72228. doi: 10.18632/oncotarget.12663.
- 26 Seo S, Hong JY, Yoon S, et al. Prognostic significance of serum beta-2 microglobulin in patients with diffuse large B-cell lymphoma in the rituximab era. *Oncotarget*. 2016 Nov 22;7(47):76934-76943. doi: 10.18632/oncotarget.12734.
- 27 Kanemasa Y, Shimoyama T, Sasaki Y, et al. Beta-2 microglobulin as a significant prognostic factor and a new risk model for patients with diffuse large B-cell lymphoma. *Hematol Oncol*. 2017 Dec;35(4):440-446. doi: 10.1002/hon.2312.
- 28 Buske C, Hutchings M, Ladetto M, et al. ESMO Consensus Conference on malignant lymphoma: general perspectives and recommendations for the clinical management of the elderly patient with malignant lymphoma. *Ann Oncol*. 2018 Mar 1;29(3):544-562. doi: 10.1093/annonc/mdx413.
- 29 Thomas L. *Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory Results*. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft 1998;685-688.
- 30 Bernard AM, Vyskocil A, Lauwerys RR. Determination of b2-Microglobulin in Human Urine and Serum by Latex Immunoassay. *Clin Chem* 1981;27(6):832-837.
- 31 European Confederation of Laboratory Medicine (ECLM). *European Urinalysis Guidelines*. *Scand J Clin Lab Invest* 2000;60:1-96.
- 32 Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. *Clin Chem* 1986;32:470-475.
- 33 Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". *Eur J Clin Chem Clin Biochem* 1996;34:385-386.
- 34 Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. *Ann Clin Biochem* 2001;38:376-385.
- 35 Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. *Clin Chem Lab Med* 2007;45(9):1240-1243.
- 36 Chan PC, Kulasingam V, Lem-Ragosnig B. Validating urinary measurement of beta-2-microglobulin with a Roche reagent kit designed for serum measurements. *Clin Biochem* 2012 Nov;45(16-17):1533-1535.
- 37 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. *J Clin Chem Clin Biochem* 1988 Nov;26(11):783-790.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

#### Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard:

 CONTENT

Contents of kit



Volume for reconstitution

 GTIN

Global Trade Item Number

Rx only

For USA: Caution: Federal law restricts this device to sale by or on the order of a physician.

COBAS, NAVIFY and TINA-QUANT are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2024, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim  
www.roche.com

+800 5505 6606

