



Density Reference Standard Beads (DRSB)

Version 43

Content version: January 2022

Beads for one-point density calibration

Cat. No. 06 422 667 001

Batch B

Store Beads at +2 to +8°C

1.	Introduction	3
	Contents	3
	Storage and Stability	3
2.	How to Use this Product	4
3.	Protocol	5
3.1	Preparation of the DRSB solution	5
3.2	Checking the FlowFactor (FF) (see Figure 2)	6
4.	Lot Specific Data	7
5.	Adjusting the FlowFactor	18
5.1	How to Calculate and Change the FlowFactor	18
5.2	FlowFactor Calibration	19
6.	Supplementary Information	20
6.1	Conventions	20
6.1.1	Text Conventions	20
6.1.2	Symbols	20
6.2	Changes to Previous Version	21
6.3	Trademarks	21
6.4	Regulatory Disclaimer	21
6.5	Contact Support	21

1. Introduction

The Cedex HiRes Analyzer measures the cell density, also known as concentration, of a cellular suspension, along with its viability status and cell characteristics such as diameter and compactness. In order to check the correct calibration with regard to the density, use the Density Reference Standard Beads instead of an ordinary cell sample.

This product is traceable in the following aspects:

- Particle Size: NIST, USA
- Particle Concentration: Physikalisch-Technische Bundesanstalt, Berlin, Germany

The Density Reference Standard Beads are made to mimic cell behavior in flow dynamics. Due to their size and optical properties, they will appear as dead cells in the Cedex Software.

Contents

Content	Volume	Cat. No.
Beads for one-point density calibration	10 ml	06 422 667 001

Storage and Stability

Store Beads at +2 to +8°C.

The product is stable until the expiry date printed on the label, when handled as described in these Instructions for Use.

2. How to Use this Product

There is no general advice with regard to how often or how many counts should be done in order to ensure that your analyzer is working correctly. Roche Diagnostics has had good results using the DRSB on a monthly basis, carrying out 10 samples per run with the Cedex HiRes Analyzer.

Two factors influence the acceptance range for calibration with Density Reference Standard Beads in connection with the Cedex HiRes Analyzer.

- **Sample preparation:** Mixing, pipette operation, and pipette quality (precision, accuracy, service state) have been shown in the field to add approximately 1.5 – 2% to the variability in density measurements.
- **Measurement precision is based on the statistical nature of the measurement process.** It depends on the density of the DRSB used, Cell Type parameter settings, and the level of precision used for the measurement.

Sampling quality is essential for the evaluation of the status of the instrument. Special care should be taken to ascertain, for example, that among other factors:

- Beads were not frozen, but stored properly at +2 to +8°C.
- Beads were allowed to acclimate to a temperature of +23 to +27°C prior to use.
- Weight of the unopened bottle is correct (see label on the bottle).
- Ultrasonic bath is used for mixing.
- Bottle is rocked gently, including rocking upside down.
- No more than 2 samples are drawn from the bottle without intermediate remixing.
- Only calibrated pipettes are used.
- Only trained staff are performing the sample preparation.

3. Protocol

3.1 Preparation of the DRSB solution

- Verify that the beads have been stored correctly at temperatures of +2 to +8°C (BEADS CANNOT BE FROZEN).
 - Verify that the bottle was securely closed before use (check the weight of the unopened bottle; the correct value is given on the bottle).
 - Allow the beads to acclimate to a temperature of +23 to +27°C prior to use.
 - Use an ultrasonic bath at a temperature of +23 to +27°C and at the highest available intensity for 5 minutes (with cap slightly loosened but secured against falling over) to shake the beads.
 - Ensure that no beads are sticking to the base or side of the bottle before use.
- 🕒 The DRSB solution contains SDS, which may show signs of some coagulation or crystallization at low temperatures. Crystals and/or signs of coagulation can be removed by allowing the beads to acclimate, with occasional mixing, at +25°C until the coagulation disappears. Alternatively, the DRSB bottle can be gently rolled between the palms of the hands until the coagulation has disappeared. Note that as long as the DRSB solution has been allowed to acclimate to +23 to +27°C, and all steps in this Instructions for Use have been carried out, any remaining coagulation or crystallization will have no effect on the performance or quality of the DRSB solution when used in a Cedex HiRes Analyzer.

3.2 Checking the FlowFactor (FF) (see Figure 2)

-
- ① Pipet 1 sample of 0.3 ml DRSB into a Cedex Sample cup, and run the sample with factory settings for default Cell Type Std. Size immediately. Select the maximum possible setting for “precision”.
 - ② Mix the DRSB thoroughly, then pipet the next sample of 0.3 ml into a Cedex Sample cup, and run the sample immediately.
 - ③ Repeat this procedure until 10 samples are processed.
 - ④ Calculate the mean value of the Total Cell Density (TCD) of the 10 samples used.
 - ⑤ Calculate the relative standard deviation of the TCD values of the 10 samples used, and verify that the relative standard deviation is less than or equal to 5%. Otherwise, the Cedex HiRes Analyzer, the beads, or the handling have to be checked and the calibration must be repeated.
 - ⑥ Calculate the deviation of the mean TCD value of the 10 samples used from the actual value (given as Particle number/ml on the bottle of beads).
 - ⑦ Verify that the deviation of the mean TCD value is less than or equal to 5%, or as specified by your requirements, from the actual value given on the bottle for the beads. If the value falls outside of the acceptable range, skip to Step 9.
 - ⑧ Close bottle tightly and store beads at +2 to +8°C (BEADS CANNOT BE FROZEN). The current FF is correct and no change is necessary.
 - ⑨ If the mean value falls out of range, repeat Steps 1 to 5 using a second/different LOT (batch) of Density Reference Standard Beads. Continue with Step 10.
 - ⑩ Calculate the FlowFactor (FF) of each measurement series (see 5.1, “How to Calculate and Change the FlowFactor”), and the mean value of the two FFs.
 - ⑪ Verify that the deviation of the two FFs from the mean value of the FFs, are less than or equal to 5%. Otherwise, the Cedex HiRes Analyzer, the beads, or the handling have to be checked and the calibration must be repeated.
 - ⑫ Calculate the new FF (mean value of the FFs); (see 5.1, “How to Calculate and Change the FlowFactor”), or follow your company’s requirements.
-

4. Lot Specific Data

Cat. No. 06 422 667 001, Batch B

Valid for Lot. No. 57130074

In this chapter, you will find lot specific data about your product. The table below provides the following information for each bottle produced for this lot.

Column 1: Bottle No. for the bottle.

Column 2: Actual concentration expressed in particle number/ml for the bottle.

Column 3: Total weight of bottle, including bottle, contents, and label.

Column 4: Check Box for marking which bottle was received.

④ Use this table as follows:

- ① Print out the table.
- ② Find the bottle number on the bottle label as shown in Figure 1.
- ③ Place a check mark in the "Bottle Received" column to mark the specific bottle received for future reference.

④ This product is traceable in the following aspects:

- Particle Size: NIST, USA
- Particle Concentration: Physikalisch-Technische Bundesanstalt, Berlin, Germany

LOT 57130074/B30



10

Density Reference Standard Beads Batch B

06 422 667 001

Particle diameter 10 µm +/- 0,2

10 ml

Particle number/ml 10.06 x 10⁵

Store at +2 to +8°C

Total weight 29.991 g

Roche Diagnostics GmbH
Mannheim, Germany +49 621 759 0
Roche Diagnostics
Indianapolis, IN, USA 001 (800) 428 5433




Fig. 1: Example of how to find the bottle number on the bottle label. The bottle number is circled.

Lot Specific Data

ID-Nr. LOT 57130074	Concentration Particle number/ml (10 ⁶)	weight (g)	Bottle received
B1	9.98	29.576	
B2	10.03	30.300	
B3	10.00	30.181	
B4	10.05	30.104	
B5	10.02	30.631	
B6	10.05	29.538	
B7	9.96	29.660	
B8	10.01	29.644	
B9	9.91	29.671	
B10	9.98	29.666	
B11	9.95	29.717	
B12	9.95	29.790	
B13	9.97	29.508	
B14	10.08	29.762	
B15	9.99	29.914	
B16	10.04	29.987	
B17	9.99	30.328	
B18	10.00	29.865	
B19	10.01	29.952	
B20	10.00	30.010	
B21	10.05	29.572	
B22	9.99	29.606	
B23	9.94	29.716	
B24	9.94	29.547	
B25	10.05	29.595	
B26	10.03	29.999	
B27	9.92	29.562	
B28	9.96	29.765	
B29	10.05	29.940	

ID-Nr. LOT 57130074	Concentration Particle number/ml (10 ⁵)	weight (g)	Bottle received
B30	9.93	29.770	
B31	10.02	30.322	
B32	10.00	30.008	
B33	9.99	30.185	
B34	10.00	29.714	
B35	10.02	29.926	
B36	10.04	29.863	
B37	10.00	29.730	
B38	9.97	29.726	
B39	9.99	29.881	
B40	9.94	29.718	
B41	10.03	29.654	
B42	10.02	29.368	
B43	10.01	29.768	
B44	9.99	29.824	
B45	10.02	29.646	
B46	9.94	29.942	
B47	10.03	29.550	
B48	10.03	29.540	
B49	9.98	29.677	
B50	9.99	29.654	
B51	10.02	29.835	
B52	10.01	29.514	
B53	9.95	29.675	
B54	10.06	29.721	
B55	10.03	29.681	
B56	9.97	29.673	
B57	10.06	29.959	
B58	10.04	29.720	

ID-Nr. LOT 57130074	Concentration Particle number/ml (10 ⁵)	weight (g)	Bottle received
B59	9.98	29.594	
B60	9.96	29.689	
B61	9.97	29.778	
B62	10.00	29.873	
B63	9.93	29.547	
B64	10.02	29.816	
B65	10.03	30.077	
B66	10.04	29.995	
B67	9.91	29.744	
B68	9.95	29.770	
B69	9.99	29.646	
B70	10.06	29.522	
B71	9.96	29.687	
B72	10.03	29.673	
B73	9.97	29.784	
B74	10.01	29.652	
B75	10.05	29.698	
B76	9.98	30.122	
B77	9.95	29.389	
B78	9.97	29.645	
B79	9.99	30.120	
B80	9.94	29.567	
B81	10.04	29.626	
B82	9.99	29.884	
B83	10.00	29.599	
B84	10.02	29.541	
B85	10.02	29.722	
B86	9.95	29.676	
B87	10.05	29.787	

ID-Nr. LOT 57130074	Concentration Particle number/ml (10 ⁵)	weight (g)	Bottle received
B88	9.92	29.699	
B89	9.96	29.545	
B90	9.98	29.718	
B91	10.04	29.554	
B92	10.01	29.849	
B93	9.98	30.425	
B94	9.98	29.660	
B95	9.99	29.968	
B96	9.97	29.526	
B97	9.98	29.558	
B98	10.04	29.621	
B99	9.96	29.558	
B100	9.94	30.014	
B101	10.06	29.681	
B102	9.93	29.513	
B103	10.07	29.521	
B104	10.04	30.330	
B105	10.00	29.524	
B106	10.04	29.635	
B107	10.04	29.373	
B108	10.03	30.003	
B109	10.06	29.527	
B110	10.02	29.862	
B111	10.02	29.545	
B112	9.98	29.612	
B113	9.95	29.400	
B114	10.03	29.494	
B115	10.04	29.865	
B116	9.99	29.867	

Lot Specific Data

ID-Nr. LOT 57130074	Concentration Particle number/ml (10 ⁶)	weight (g)	Bottle received
B117	9.98	29.440	
B118	9.97	29.582	
B119	9.97	29.679	
B120	9.92	29.724	
B121	9.99	29.565	
B122	10.03	29.717	
B123	9.92	29.561	
B124	9.98	29.812	
B125	10.00	29.921	
B126	9.98	29.524	
B127	10.02	29.627	
B128	9.99	29.674	
B129	10.01	29.978	
B130	9.96	29.693	
B131	10.00	29.607	
B132	9.99	29.586	
B133	10.05	29.564	
B134	10.04	29.490	
B135	10.05	29.747	
B136	9.96	29.897	
B137	10.04	29.842	
B138	9.93	29.671	
B139	10.01	29.542	
B140	9.94	29.692	
B141	10.06	29.755	
B142	9.96	29.598	
B143	9.96	29.641	
B144	9.95	30.202	
B145	9.90	30.122	

ID-Nr. LOT 57130074	Concentration Particle number/ml (10 ⁵)	weight (g)	Bottle received
B146	9.96	30.076	
B147	10.14	30.085	
B148	9.78	29.983	
B149	9.86	30.720	
B150	9.91	29.731	
B151	10.01	29.822	
B152	9.91	30.348	
B153	9.94	30.703	
B154	10.12	30.222	
B155	9.74	30.879	
B156	10.07	30.700	
B157	10.03	29.426	
B158	10.00	29.558	
B159	10.11	29.517	
B160	10.04	29.453	
B161	10.03	29.527	
B162	10.20	29.416	
B163	10.01	29.651	
B164	10.07	29.472	
B165	9.88	29.580	
B166	9.89	29.423	
B167	9.87	29.464	
B168	10.07	29.527	
B169	9.94	29.550	
B170	9.90	29.549	
B171	9.82	29.725	
B172	9.81	29.622	
B173	9.72	29.703	
B174	10.10	29.639	

ID-Nr. LOT 57130074	Concentration Particle number/ml (10 ⁵)	weight (g)	Bottle received
B175	9.80	29.598	
B176	9.93	29.631	
B177	10.00	29.584	
B178	9.72	29.787	
B179	10.01	29.262	
B180	9.85	29.588	
B181	10.23	29.506	
B182	10.03	29.604	
B183	10.00	29.453	
B184	9.84	29.693	
B185	10.18	29.536	
B186	10.01	29.568	
B187	9.98	29.514	
B188	10.21	29.580	
B189	9.96	29.538	
B190	10.07	29.618	
B191	10.27	29.437	
B192	10.23	29.515	
B193	9.92	29.571	
B194	10.04	29.581	
B195	9.99	29.570	
B196	10.06	29.316	
B197	9.85	29.598	
B198	10.05	29.520	
B199	10.05	29.508	
B200	10.08	29.385	
B201	10.18	29.496	
B202	10.29	29.133	
B203	9.82	29.535	

ID-Nr. LOT 57130074	Concentration Particle number/ml (10 ⁵)	weight (g)	Bottle received
B204	10.13	29.652	
B205	9.91	29.431	
B206	10.06	29.584	
B207	9.89	29.605	
B208	9.97	29.444	
B209	9.89	29.649	
B210	10.00	29.575	
B211	10.30	29.691	
B212	10.13	29.402	
B213	9.94	29.603	
B214	10.02	29.668	
B215	9.94	29.625	
B216	10.06	29.448	
B217	9.92	29.384	
B218	10.14	29.782	
B219	10.08	29.706	
B220	10.08	29.712	
B221	9.74	29.656	
B222	9.98	29.709	
B223	10.02	29.615	
B224	9.89	29.730	
B225	9.97	29.593	
B226	9.94	29.687	
B227	9.90	29.566	
B228	9.87	29.661	
B229	9.86	29.786	
B230	10.02	29.745	
B231	9.86	29.465	
B232	9.73	29.488	

Lot Specific Data

ID-Nr. LOT 57130074	Concentration Particle number/ml (10 ⁶)	weight (g)	Bottle received
B233	9.97	29.603	
B234	10.20	29.515	
B235	10.10	29.545	
B236	10.15	29.630	
B237	10.17	29.623	
B238	10.12	29.550	
B239	9.99	29.688	
B240	10.03	29.557	
B241	9.87	29.536	
B242	9.92	29.603	
B243	9.91	29.310	
B244	10.05	29.497	
B245	9.93	29.572	
B246	9.96	29.476	
B247	10.04	29.558	
B248	10.08	29.621	
B249	9.93	29.569	
B250	10.00	29.546	
B251	10.26	29.625	
B252	10.09	29.520	
B253	10.09	29.702	
B254	9.81	29.602	
B255	9.96	29.709	
B256	10.04	29.697	
B257	10.07	29.802	
B258	9.94	29.599	
B259	9.93	29.616	
B260	10.13	29.479	
B261	9.85	29.733	

ID-Nr. LOT 57130074	Concentration Particle number/ml (10 ⁵)	weight (g)	Bottle received
B262	9.93	29.476	
B263	9.95	29.643	
B264	9.76	29.524	
B265	9.92	29.710	
B266	9.91	29.592	
B267	10.17	29.590	
B268	10.17	29.618	
B269	9.86	29.687	
B270	9.98	29.423	
B271	9.99	29.612	
B272	9.93	29.832	
B273	10.01	29.537	
B274	9.94	29.599	
B275	10.20	29.583	
B276	10.01	29.491	
B277	10.06	29.606	
B278	9.75	29.571	
B279	9.76	29.736	
B280	9.91	29.688	
B281	10.04	29.776	
B282	9.87	29.659	
B283	9.91	29.676	
B284	9.85	29.441	
B285	9.84	29.652	
B286	9.87	29.636	
B287	9.83	29.419	
B288	10.15	29.631	
B289	9.75	29.843	
B290	9.79	29.740	

ID-Nr. LOT 57130074	Concentration Particle number/ml (10 ⁵)	weight (g)	Bottle received
B291	9.99	29.582	
B292	10.00	29.640	
B293	9.92	29.689	
B294	9.86	29.696	
B295	9.78	29.636	
B296	9.93	29.468	
B297	9.92	29.719	
B298	9.75	29.598	
B299	9.88	29.819	
B300	10.17	29.551	
B301	10.07	29.823	
B302	9.95	29.632	
B303	10.03	29.896	
B304	9.88	29.800	
B305	9.85	29.708	
B306	9.92	29.787	
B307	9.90	29.671	
B308	9.78	29.945	
B309	9.90	29.852	
B310	9.86	29.634	
B311	9.85	29.665	
B312	9.98	29.929	
B313	9.80	29.743	
B314	10.07	29.570	
B315	9.82	29.654	
B316	9.93	29.682	
B317	10.17	29.584	
B318	10.04	29.562	
B319	10.16	29.600	

ID-Nr. LOT 57130074	Concentration Particle number/ml (10 ⁵)	weight (g)	Bottle received
B320	10.20	29.808	
B321	10.23	29.783	
B322	10.05	29.540	
B323	7.52	29.547	
B324	9.99	29.692	
B325	10.20	29.631	
B326	10.09	29.574	
B327	10.16	29.643	
B328	9.91	29.511	
B329	9.81	29.823	
B330	10.00	29.648	
B331	9.80	29.756	
B332	9.99	29.604	
B333	9.93	29.680	
B334	9.74	29.647	
B335	9.94	29.605	
B336	10.06	29.323	
B337	9.94	29.684	
B338	10.11	29.751	
B339	10.02	29.532	
B340	9.81	29.729	
B341	9.93	29.686	
B342	9.99	29.753	
B343	10.00	29.896	
B344	9.96	29.722	
B345	9.87	29.653	
B346	9.72	29.477	
B347	10.06	29.464	
B348	10.24	29.626	

Lot Specific Data

ID-Nr. LOT 57130074	Concentration Particle number/ml (10 ⁵)	weight (g)	Bottle received
B349	10.09	29.669	
B350	10.13	29.742	
B351	10.05	29.624	
B352	9.87	29.563	
B353	9.80	29.750	
B354	10.02	29.657	
B355	10.07	29.686	
B356	9.88	29.626	
B357	10.04	29.662	
B358	10.04	29.653	
B359	9.80	29.749	
B360	9.82	29.705	
B361	9.97	29.354	
B362	9.95	29.853	
B363	9.93	29.614	
B364	10.07	29.730	
B365	9.88	29.804	
B366	10.24	29.689	
B367	10.07	29.613	
B368	10.13	29.665	
B369	10.11	29.737	
B370	9.84	29.428	
B371	10.25	29.829	
B372	9.84	29.676	
B373	9.93	29.716	
B374	9.79	29.490	
B375	9.78	29.563	
B376	9.89	29.669	
B377	10.07	29.526	

ID-Nr. LOT 57130074	Concentration Particle number/ml (10 ⁵)	weight (g)	Bottle received
B378	10.11	29.554	
B379	10.00	29.815	
B380	9.94	29.838	
B381	9.94	29.586	
B382	9.81	29.678	
B383	9.85	29.608	
B384	9.84	30.789	
B385	9.73	29.643	
B386	10.04	29.798	
B387	9.97	29.942	
B388	9.88	29.715	
B389	9.79	29.848	
B390	9.87	29.561	
B391	10.02	29.772	
B392	9.99	29.633	
B393	9.88	29.604	
B394	9.76	29.859	
B395	9.86	29.707	
B396	9.91	29.565	
B397	9.83	29.785	
B398	10.00	29.651	
B399	9.94	29.699	
B400	9.98	29.745	
B401	10.11	29.647	
B402	9.94	29.864	
B403	10.08	29.764	
B404	9.79	29.576	
B405	10.18	29.712	
B406	9.95	29.740	

ID-Nr. LOT 57130074	Concentration Particle number/ml (10 ⁶)	weight (g)	Bottle received
B407	10.08	29.632	
B408	9.88	29.546	
B409	10.03	29.615	
B410	10.21	29.850	
B411	9.76	29.714	
B412	9.82	29.684	
B413	9.85	29.613	
B414	9.90	29.663	
B415	9.93	29.708	
B416	9.89	29.740	
B417	9.89	29.585	
B418	9.97	29.753	
B419	9.73	29.670	
B420	9.80	29.770	
B421	9.92	29.693	
B422	9.99	29.579	
B423	9.97	29.536	
B424	9.91	29.728	
B425	10.01	29.715	
B426	10.05	29.776	
B427	10.07	29.745	
B428	10.01	29.611	
B429	9.96	29.479	
B430	9.90	29.528	
B431	9.87	29.599	
B432	9.98	29.626	
B433	10.02	29.456	
B434	9.77	29.658	
B435	9.73	29.661	

ID-Nr. LOT 57130074	Concentration Particle number/ml (10 ⁵)	weight (g)	Bottle received
B436	9.92	29.631	
B437	9.88	29.681	
B438	9.84	29.719	
B439	9.81	29.604	
B440	9.71	29.723	
B441	9.90	29.628	
B442	9.97	29.724	
B443	9.86	29.638	
B444	9.81	29.774	
B445	9.95	29.681	
B446	9.71	29.482	
B447	9.75	29.690	
B448	9.98	29.735	
B449	9.97	29.560	
B450	10.25	29.599	
B451	10.01	29.486	
B452	9.86	29.503	
B453	10.11	29.568	
B454	9.94	29.507	
B455	9.99	29.558	
B456	9.97	29.504	
B457	9.94	29.494	
B458	9.77	29.405	
B459	9.80	29.627	
B460	9.82	29.584	
B461	10.03	29.470	
B462	9.97	29.399	
B463	9.84	29.388	
B464	9.82	29.752	

Lot Specific Data

ID-Nr. LOT 57130074	Concentration Particle number/ml (10 ⁶)	weight (g)	Bottle received
B465	9.95	29.578	
B466	10.14	29.684	
B467	9.97	29.882	
B468	9.81	29.713	
B469	9.84	29.719	
B470	9.92	29.835	
B471	10.01	29.677	
B472	9.87	29.780	
B473	9.82	29.889	
B474	10.24	29.703	
B475	9.93	29.777	
B476	9.83	29.565	
B477	9.76	29.818	
B478	10.17	29.769	
B479	9.88	29.752	
B480	10.08	29.590	
B481	9.99	29.722	
B482	9.77	29.706	
B483	9.84	29.717	
B484	9.93	29.641	
B485	10.28	29.770	
B486	10.03	29.538	
B487	10.06	29.627	
B488	9.95	29.856	
B489	9.79	29.748	
B490	9.90	29.776	
B491	10.07	29.757	
B492	9.99	29.811	
B493	9.77	29.712	

ID-Nr. LOT 57130074	Concentration Particle number/ml (10 ⁵)	weight (g)	Bottle received
B494	9.81	29.708	
B495	9.86	29.725	
B496	9.86	29.591	
B497	10.07	29.645	
B498	9.86	30.305	
B499	9.77	30.285	
B500	9.89	29.920	
B501	10.12	29.711	
B502	10.09	29.958	
B503	10.04	29.945	
B504	9.94	30.025	
B505	9.96	29.964	
B506	9.74	30.041	
B507	10.02	29.882	
B508	10.24	29.856	
B509	9.97	29.545	
B510	10.12	29.582	
B511	10.02	29.862	
B512	9.94	29.890	
B513	9.85	30.075	
B514	9.89	30.023	
B515	9.92	30.212	
B516	9.87	30.149	
B517	9.85	29.971	
B518	9.80	29.941	
B519	10.01	29.725	
B520	9.81	29.563	
B521	9.74	29.482	
B522	9.76	29.462	

ID-Nr. LOT 57130074	Concentration Particle number/ml (10 ⁵)	weight (g)	Bottle received
B523	9.86	29.579	
B524	9.95	29.464	
B525	9.97	29.627	
B526	9.90	29.471	
B527	9.84	29.586	
B528	9.88	29.661	
B529	10.03	29.750	
B530	10.08	29.630	
B531	9.97	29.536	
B532	9.99	29.689	
B533	10.15	29.746	
B534	9.88	29.581	
B535	10.02	29.676	
B536	10.01	29.569	
B537	10.02	29.520	
B538	9.91	29.849	
B539	10.17	29.681	
B540	10.24	29.556	
B541	10.10	29.520	
B542	10.17	29.617	
B543	10.15	29.597	
B544	10.02	29.552	
B545	10.26	29.547	
B546	10.08	29.554	
B547	10.31	29.472	
B548	9.99	29.604	
B549	10.14	29.608	
B550	10.26	29.605	
B551	10.19	29.816	

ID-Nr. LOT 57130074	Concentration Particle number/ml (10 ⁵)	weight (g)	Bottle received
B552	9.96	30.112	
B553	10.07	29.626	
B554	10.11	29.499	
B555	10.10	29.712	
B556	10.02	29.458	
B557	9.97	29.688	
B558	9.95	29.652	
B559	10.25	29.670	
B560	10.00	29.613	

5. Adjusting the FlowFactor

The FlowFactor (FF) is analyzer specific and part of the conversion factor that relates the number of objects detected in the Cedex HiRes Analyzer to the actual density in the analyzed sample.

The conversion factor is linearly dependent on the FlowFactor, thus allowing for the possibility to compute the appropriate setting for this parameter via the comparison of Cedex HiRes Analyzer results versus a known density of a sample (*e.g.*, Density Reference Standard Beads).

Refer to the relevant Cedex HiRes Operator's Manual for the location of the current FlowFactor. The location depends on the installed Software version.

5.1 How to Calculate and Change the FlowFactor

- 1 Write down the current FlowFactor (FFold) and calculate a new FlowFactor as follows:

$$\text{FF (new}_1\text{)} = \frac{\text{actual density (according to bottle label)}}{\text{mean value TCD of measurement series 1}} \times \text{FF (old)}$$

$$\text{FF (new}_2\text{)} = \frac{\text{actual density (according to bottle label)}}{\text{mean value TCD of measurement series 2}} \times \text{FF (old)}$$

("actual density" is specified as Particle number/ml on the label of the bottle of beads used for the calibration.)

$$\text{FF (new)} = \frac{\text{FF (new}_1\text{)} + \text{FF (new}_2\text{)}}{2}$$

- 2 Refer to the relevant Cedex HiRes Operator's Manual for information about the location of the FlowFactor. Update the FlowFactor in that location based on the result calculated in Step 1.

5.2 FlowFactor Calibration




Fig. 2: Calibration Cedex HiRes Analyzer

6. Supplementary Information

6.1 Conventions

6.1.1 Text Conventions

To make information consistent and easier to read, the following text conventions are used in this document:

Text Convention	Usage
Numbered stages labeled ①, ②, etc.	Stages in a process that usually occur in the order listed.
Numbered instructions labeled ❶, ❷, etc.	Steps in a procedure that must be performed in the order listed.
Asterisk *	Denotes a product available from Roche Diagnostics.

6.1.2 Symbols

In this document, the following symbols are used to highlight important information:

Symbol	Description
ⓘ	Information Note: Additional information about the current topic or procedure.
⚠	Important Note: Information critical to the success of the procedure or use of the product.

6.2 Changes to Previous Version

- Updated to include lot-specific data for new lot.

6.3 Trademarks

CEDEX is a trademark of Roche.

Other brands or product names are trademarks of their respective holders.

6.4 Regulatory Disclaimer

For use in quality control/manufacturing process only.

6.5 Contact Support

For additional documentation such as certificates and safety data sheets, please visit documentation.roche.com.

Published by

Roche Diagnostics GmbH
Sandhofer Straße 116
68305 Mannheim
Germany

© 2022 Roche Diagnostics.
All rights reserved.

06326056001 (43) 0122

For more information about this product please visit documentation.roche.com

For more documentation such as instructions for use and safety data sheets, please visit documentation.roche.com

Your Roche CustomBiotech Customer Service

**Europe, Middle East, Africa,
and Latin America**

Phone +49 621 759 8580
Fax +49 621 759 6385
mannheim.custombiotech@roche.com

United States

Phone +1 800 428 5433 (toll-free)
Fax +1 317 521 4065
custombiotech.ussales@roche.com

Canada

Phone +1 450 686 7050
Fax +1 450 686 7012
custombiotech.can@roche.com

Japan

Phone +81 3 6634 1046
Fax +81 3 5479 0585
japan.custombiotech@roche.com

Asia Pacific

Phone +65 6371 6638
Fax +65 6371 6601
apac.custombiotech@roche.com