

For general laboratory use.
This product is not available in all territories due to
different national regulations.
This document is not intended for use in the USA.



DNA Process Control Kit

 **Version: 04**

Content version: September 2016

Includes DNA Process Control, Control Assay, and multiplex master mix for monitoring the whole workflow, from nucleic acid extraction to qPCR.

Cat. No. 07 339 542 001 1 kit
6 x 96 DNA Control reactions and 600 qPCR reactions of 20 µl
final volume each

Cat. No. 07 339 666 001 1 kit
2 x 96 DNA Control reactions and 200 qPCR reactions of 20 µl
final volume each

Refill Kit for the DNA Process Control detection only.

DNA Process Control Detection Kit

Cat. No. 07 339 623 001 1 kit
400 qPCR reactions of 20 µl final volume each

Store the kits at -15 to -25°C.

1.	General Information	3
1.1.	Contents	3
	DNA Process Control Kit	3
	DNA Process Control Detection Kit	3
1.2.	Storage and Stability	4
	Storage Conditions (Product)	4
1.3.	Additional Equipment and Reagents Required	4
1.4.	Application	5
	Product Description	5
1.5.	Preparation Time	5
	Typical run time	5
2.	How to Use this Product	6
2.1.	Before you Begin	6
	Sample Materials	6
	Control Reactions	6
	Primers	6
	Probe	6
	General Considerations	7
	Precautions	7
	Safety Information	7
	Laboratory Procedures	7
	Working Solution	7
2.2.	Protocols	8
	Isolation of DNA	8
	Preparation of the qPCR	9
	Setting up the qPCR Reaction	10
	LightCycler® 480 Instrument II Protocol	10
	LightCycler® 96 Instrument Protocol	13
3.	Results	14
	Use with Different Sample Types	14
	Interpretation of Data	15
	Negative Control/ No Template Control (NTC)	15
	Positive Controls	15
	Sample Result	15
	DNA Process Control/Internal Control	15
4.	Troubleshooting	16
5.	Additional Information on this Product	17
5.1.	Test Principle	17
	Background Information	17
	How this Product Works	17
5.2.	Quality Control	18
6.	Supplementary Information	19
6.1.	Conventions	19
6.2.	Changes to previous version	19
6.3.	Ordering Information	20
6.4.	Trademarks	21
6.5.	License Disclaimer	21
6.6.	Regulatory Disclaimer	21
6.7.	Safety Data Sheet	21
6.8.	Contact and Support	21

1. General Information

1.1. Contents

DNA Process Control Kit

Vial / Bottle	Cap	Label	Function / Description	Catalog Number	Content
1	orange	DNA Process Control Kit, DNA Process Control, conc.	Nuclease-resistant DNA concentrate in storage buffer	07 339 542 001	6 vials, 400 µl each
				07 339 666 001	2 vials, 400 µl each
2	white	DNA Process Control Kit, DNA Process Control Diluent	DNA Process Control Diluent	07 339 542 001	3 bottles, 17 ml each
				07 339 666 001	1 bottle, 17 ml
3	red	DNA Process Control Kit, qPCR Reaction Mix, 5x conc.	LightCycler® DNA Multiplex Master	07 339 542 001	3 vials, 880 µl each
				07 339 666 001	1 vial, 880 µl
4	yellow	DNA Process Control Kit, DNA Process Control Detection Assay, 20x conc.	Primer/Probe Mix for detection of the DNA Process Control	07 339 542 001	6 vials, 120 µl each
				07 339 666 001	2 vials, 120 µl each
5	white	DNA Process Control Kit, Water, PCR Grade	To adjust the final reaction volume	07 339 542 001	12 vials, 1 ml each
				07 339 666 001	4 vials, 1 ml each

DNA Process Control Detection Kit

Vial / Bottle	Cap	Label	Function / Description	Catalog Number	Content
3	red	DNA Process Control Detection Kit, qPCR Reaction Mix, 5x conc.	LightCycler® DNA Multiplex Master	07 339 623 001	2 vials, 880 µl each
4	yellow	DNA Process Control Detection Kit, DNA Process Control Detection Assay, 20x conc.	Primer/Probe Mix for detection of the DNA Process Control	07 339 623 001	4 vials, 120 µl each
5	white	DNA Process Control Detection Kit, Water, PCR Grade	To adjust the final reaction volume	07 339 623 001	8 vials, 1 ml each

i The DNA Process Control Detection Kit can only be used in combination with the DNA Process Control that is contained in the DNA Process Control Kit.

1.2. Storage and Stability

Storage Conditions (Product)

The kits are shipped on dry ice.

When stored at –15 to –25°C, the kits are stable through the expiration date printed on the label.

⚠ Store the kits protected from light.

Vial / Bottle	Cap	Label	Storage
1	orange	DNA Process Control, conc.	Store at –15 to –25°C for up to 12 months. ⚠ Avoid repeated freezing and thawing.
2	white	DNA Process Control Diluent	Store at –15 to –25°C for up to 12 months or +2 to +8°C for a maximum of 4 weeks.
3	red	qPCR Reaction Mix, 5x conc.	Store at –15 to –25°C. ⚠ Avoid repeated freezing and thawing (more than 5 times). Aliquot Vial 3 and store at –15 to –25°C or store Vial 3 at +2 to +8°C for a maximum of 4 weeks.
4	yellow	DNA Process Control Detection Assay, 20x conc.	Store at –15 to –25°C for up to 12 months. ⚠ Avoid repeated freezing and thawing.
5	white	Water, PCR Grade	Store at +2 to +8°C for up to 12 months.

1.3. Additional Equipment and Reagents Required

Standard Laboratory Equipment

- Nuclease-free pipette tips
- 1.5 ml RNase/DNase-free microcentrifuge tubes to prepare master mixes and dilutions.
- To minimize risk of nuclease contamination, autoclave all vessels.
- Wear gloves at all times.

For qPCR

- LightCycler® 480 Instrument II* or the LightCycler® 96 Instrument*
- LightCycler® 480 Multiwell Plate 96, white*
- LightCycler® 480 Multiwell Plate 384, white*
- Standard swinging-bucket centrifuge with rotor for multiwell plates

For DNA Purification

- MagNA Pure 96 Instrument* including disposables
- MagNA Pure 96 Internal Control Tube*, optional
- MagNA Pure 96 DNA and Viral NA Kit, Large Volume* or
- MagNA Pure 96 DNA and Viral NA Kit, Small Volume*

Alternatively, use other nucleic acid purification instruments and manual sample preparation products such as:

- MagNA Pure LC 2.0 Instrument* with MagNA Pure LC Total Nucleic Acid Isolation Kit*
- MagNA Pure Compact Instrument* with MagNA Pure Compact Nucleic Acid Isolation Kit I*
- High Pure Viral Nucleic Acid Kit*

For Color Compensation with the LightCycler® 480 Instrument II

- LightCycler® Multiplex DNA Master*

1.4. Application

The DNA Process Control Kit is a tool to be used to control for potential failures of sample preparation, amplification, detection, and handling errors. The product is intended for use with a variety of sample materials (*e.g.*, blood, serum, stool, urine) and a variety of DNA targets.

The DNA Process Control Kit is intended for general laboratory use. Any use of the product for *in vitro* diagnostic tests is the sole responsibility of the operator and must be validated by the operator following the relevant national regulations.

Product Description

The DNA Process Control Kit contains a nuclease-resistant DNA that can be added at various steps during the nucleic acid purification of DNA detection workflows. Successful detection of this control DNA serves as a positive control, verifying the functionality of the DNA purification as well as the detection using quantitative PCR (qPCR). It also contains the LightCycler® Multiplex DNA Master optimized for multiplex reactions, ensuring performance for up to 4 targets plus the internal process control.

1.5. Preparation Time

Typical run time

The DNA Process Control purification step has various run times, from 20 minutes (High Pure) up to 50 – 60 minutes for a typical run on the MagNA Pure 96 Instrument with 96 samples, depending on the protocol.

The DNA Process Control detection can be used with a fast qPCR protocol with run times of less than 60 minutes using the LightCycler® Multiplex DNA Master on the LightCycler® 480 Instrument II or LightCycler® 96 Instrument.

2. How to Use this Product

2.1. Before you Begin

Sample Materials

Use any template DNA suitable for qPCR in terms of purity, concentration, and absence of PCR inhibitors. For reproducible isolation of nucleic acids, we recommend:

- Either a MagNA Pure System together with a dedicated reagent kit (for automated isolation),
- or a High Pure Nucleic Acid Isolation Kit (for manual isolation).

Options for Nucleic Acid Isolation and qPCR Using the DNA Process Control Kit	
Nucleic Acid Isolation	qPCR
MagNA Pure 96 Instrument	
or	
MagNA Pure LC Instrument	LightCycler® 480 Instrument II
or	or
MagNA Pure Compact Instrument	LightCycler® 96 Instrument
or	
High Pure Viral Nucleic Acid Kit	

Control Reactions

Control Samples and No Template Controls (NTCs)

Always run appropriate negative and positive controls for each of your target parameters.

- In addition to the DNA Process Control, we highly recommend testing known negative and known positive samples as controls in each run to check the complete procedure, including sample preparation and qPCR.
- If the DNA Process Control is added manually to the sample material, we recommend having one known positive and one known negative sample without the DNA Process Control.
- If the DNA Process Control is added automatically to all samples, we recommend performing a no template control (NTC) PCR using Water, PCR Grade (Vial 5) instead of sample eluate in the qPCR. The negative control is necessary for determining potential contamination issues.

Primers

Suitable concentrations of qPCR primers range from 0.2 to 0.5 μM (final concentration in qPCR). The recommended starting concentration is 0.5 μM each.

Probe

Suitable concentrations of hydrolysis probes range from 0.05 to 0.5 μM (final concentration in qPCR). The recommended starting concentration is 0.25 μM each.

- i** *The optimal probe concentration is the lowest concentration that results in the lowest quantification cycle value (Cq) and adequate fluorescence dynamics for a given target concentration.*
- i** *To ensure efficient probe cleavage, the T_m of the hydrolysis probe has to be higher than the T_m of the primers.*

General Considerations

Precautions

Use nuclease-free techniques. Nuclease-contaminated reagents and reaction vessels may degrade template DNA. Please follow these guidelines to minimize risk of contamination:

- Wear disposable gloves and change them frequently.
- Avoid touching surfaces or materials that could cause nuclease carryover.
- Use only reagents provided in this kit. Substitutions may introduce nucleases.
- Clean and decontaminate work areas and instruments, including pipettes, with commercially available decontamination reagents.
- Use only new nuclease-free aerosol-blocking pipette tips.
- Use a work area specifically designated for nucleic acid work, and if possible, use reaction vessels and pipettes dedicated only for work with template nucleic acid.

This product is for use by experienced personnel who have training in standardized molecular testing procedures and expertise in viral research, in laboratories with appropriate biosafety equipment and containment procedures.

Safety Information

Laboratory Procedures

All sample material and resulting waste should be considered potentially infectious. Thoroughly clean and disinfect all work surfaces with disinfectants recommended by the local authorities.

- Do not eat, drink, or smoke in the laboratory work area.
- Do not pipette by mouth.
- Wear protective disposable gloves, laboratory coats, and eye protection when handling samples and kit reagents.
- Avoid microbial and nuclease contamination of reagents when removing aliquots from reagent bottles. Use sterile disposable pipette tips.
- Wash hands thoroughly after handling samples and kit reagents.

In addition, to minimize the risk of carryover contamination which may result in false-positive results, follow the guidelines listed below:

- Use a PCR hood.
- Wipe and UV-illuminate PCR workstations and biosafety cabinets before use.
- Have separate areas for sample preparation, PCR reaction setup, and PCR amplification.
- Discard pipette tips in sealed containers to prevent airborne contamination.
- The DNA Process Control concentrates and the working solution must be handled with care; open and prepare the solutions in a separate location.
- Avoid opening LightCycler® 480 Multiwell Plates containing amplification products.

Working Solution

Preparation of the DNA Process Control Working Solution

The procedure below will result in a consistent DNA Process Control concentration that can be added to the sample material. The process control C_q values may vary depending on the type of sample material, the applied purification protocol, and workflow. It is recommended to use the protocol as outlined below and add 20 µl of the working solution to the sample material. If lower or higher C_qs from the DNA Process Control assay are required, the dilutions can be adjusted accordingly.

- 1 Thaw one aliquot of the DNA Process Control, conc. (Vial 1) and one bottle of DNA Process Control Diluent (Bottle 2).

- 2 Pipette 400 µl from Vial 1 and add to 3.6 ml DNA Process Control Diluent (Bottle 2). This will serve as the DNA Process Control working solution.

- 3 Vortex briefly.

2.2. Protocols

The DNA Process Control is added to the sample material directly or with the lysis buffer at the beginning of the procedure. Purification can be done using a variety of automated systems as indicated above. To control the workflow for a particular sample, the DNA Process Control can be used with target-specific assays using either a (one well) multicolor-multiplex format, or in two independent monoplex reactions.

For very low-copy targets, it is best to perform monoplex assays, as they generally produce higher sensitivity. In combination with appropriate positive target samples and valid NTCs, a positive DNA Process Control assay ensures correct sample processing during the nucleic acid purification and verifies the functionality of the qPCR detection reaction. The DNA Process Control assay is best used with the qPCR procedure described below.

The workflow for the MagNA Pure 96 nucleic acid purification and subsequent qPCR detection is described in detail below.

In case of automated primary sample handling and automated qPCR setup workflows, it can be beneficial to increase the pipetting volume of the individual reaction components. For this purpose, a 4x DNA Process Control Detection Assay working solution can be prepared.

Isolation of DNA

Example of a MagNA Pure 96 Instrument Run Using an Internal Control Position

1 Prepare the MagNA Pure 96 Instrument according to the Operator's Manual.

2 Depending on the kit used, select one of the following protocols:

For large volume kit with 500 µl sample volume:

- Pathogen Universal 500
- Viral NA Universal LV
- Viral NA Plasma LV

For small volume kit with 200 µl sample volume:

- Pathogen Universal 200
 - Viral NA Universal SV
 - Viral NA Plasma SV
-

3 Select elution volume:

- 50 µl or 100 µl
-

4 Add the required volume of **DNA Process Control working solution** to a MagNA Pure 96 Internal Control Tube and proceed according to the MagNA Pure 96 User Training Guide.

Other Nucleic Acid Purification Methods

Add the **DNA Process Control working solution** volume shown in the table below, directly into the sample material just prior to starting the purification process.

Method	Volume of DNA Process Control Working Solution [µl]
MagNA Pure LC	20.0
MagNA Pure Compact	20.0
High Pure Viral Nucleic Acid Kit	20.0

Follow the procedure below (Preparation of the qPCR) to prepare 96, 20 µl standard reactions with the LightCycler® Multiplex DNA Master, using either the LightCycler® 480 System or the LightCycler® 96 System and the LightCycler® 480 Multiwell Plate 96, white.

⚠ Do not touch the surface of the LightCycler® 480 Multiwell Plate 96 or 384.

Preparation of the qPCR

Preparation of the Master Mix for Manual qPCR Setup

- Keep all reagents on ice.
- Prepare master mix on ice as shown below for 95 samples plus 1 NTC (20 µl each, 10 µl amplifications to be used for 384-well plate setups).

Vial	Reagent(1)	Reagent Conc.	Master Mix for qPCR			
			1 Reaction [µl]		100 Reactions [µl]	
			96-well plate	384-well plate	96-well plate	384-well plate
3	qPCR Reaction Mix	5x	4.0	2.0	400	200
4	DNA Process Control Detection Assay	20x	1.0	0.5	100	50
X	Gene-Specific Primers		2.0	1.0	200	100
6	Water, PCR Grade		8.0	4.0	800	400
Master Mix Volume			15.0	7.5	1,500	750

⁽¹⁾ For eluates derived from stool samples, it is recommended to add 0.2 µg/µl (final) of molecular biology-grade Bovine Serum Albumin. Adjust for BSA volume by subtracting from PCR water volume.

Preparation of the Master Mix for Automated qPCR Setup with a DNA Process Control Detection Assay Working Solution

For automated qPCR setup environments that require larger pipetting volumes, the DNA Process Control Detection Assay, 20x conc. can be used to create a 4x working solution.

- Thaw one vial of the DNA Process Control Detection Assay, 20x conc. (Vial 4) and one vial of Water, PCR Grade (Vial 5).
- Add 480 µl water to Vial 4 (DNA Process Control Detection Assay, 20x conc.) from Step 1.
- Vortex briefly.
- Add 5 µl of this 4x working solution to a 20 µl final qPCR reaction volume to achieve a final concentration of 1x.
- Keep all reagents on ice.
- Prepare master mix on ice as shown below for 95 samples plus 1 NTC.

Vial	Reagent(1)	Reagent Conc.	Master Mix for qPCR	
			1 Reaction [µl]	100 Reactions [µl]
			3	qPCR Reaction Mix
4	DNA Process Control detection Assay	4x	5.0	500
X	Gene-Specific Primers		2.0	200
5	Water, PCR Grade		4.0	400
Master Mix Volume			15.0	1,500

⁽¹⁾ For eluates derived from stool samples, it is recommended to add 0.2 µg/µl (final) of molecular biology-grade Bovine Serum Albumin. Adjust for BSA volume by subtracting from PCR Water volume.

2. How to Use this Product

Setting up the qPCR Reaction

- 1 Place your samples on ice. Prepare the qPCR master mix as described in the pipetting protocols above (**Preparation of the Master Mix for Manual qPCR Setup** or **Preparation of the Master Mix for Automated qPCR Setup with a DNA Process Control Detection Assay Working Solution**), and place on ice.
- 2 Dispense 15 µl of the reaction mixture qPCR master mix (Step 1) and the sample eluates to the respective wells of the precooled LightCycler® 480 Multiwell Plate 96. Always start with the negative control, followed by the prepared samples. Positive controls are pipetted last.
 - Add 5 µl Water, PCR Grade (Vial 5) into the NTC position.
 - Add 5 µl sample eluate to all sample positions.
- 3 Seal the LightCycler® 480 Multiwell Plate 96 with a LightCycler® 480 Sealing Foil.
- 4 Place the multiwell plate in a standard swinging-bucket centrifuge containing a rotor for multiwell plates with suitable adapters. Balance it with a suitable counterweight (e.g., another multiwell plate), and centrifuge for 2 minutes at 1,500 × g.
- 5 Proceed immediately with the qPCR.
 - ⚠ **Avoid leaving the plate at +15 to +25°C.**

LightCycler® 480 Instrument II Protocol

The DNA Process Control Kit is compatible with the LightCycler® 480 Instrument II, 96-well or 384-well. The following procedure is optimized for use with the corresponding LightCycler® System that you are using.

⚠ **Program the LightCycler® Instrument before preparing the reaction mixes.**

A LightCycler® Instrument protocol that uses the LightCycler® Multiplex DNA Master contains the following programs:

- **Pre-Incubation** for activation of DNA polymerase and denaturation of the DNA
- **Amplification** of the DNA
- **Cooling** of the thermal block

i For details on how to program the experimental protocol, see the current LightCycler® 480 Instrument Operator's Manual.

Programming a Customized Detection Format for the LightCycler® System Filter Combination Selection

The detection format in the LightCycler® 480 Instrument II Software, Version 1.5 setup needs to be customized for the applied dual-color or triple-color-hydrolysis format used in the qPCR detection.

- Dual-color measurement and the simultaneous analysis of one additional parameter (e.g., FAM, Yellow 555, or Red 610 channel), and the DNA Process Control (Cy5).
- Triple-color measurement and simultaneous analysis of multiple parameters (e.g., FAM and Yellow 555 or FAM and Red 610), and the DNA Process Control (Cy5).

In the **Tool** module, the **Detection Formats** option allows creating new detection formats specified by the user, including a **Detection Format** list, a **Filter Combination** selection area, and **Selected Filter Combination List**. Different filter settings for the LightCycler® 480 Instrument II are defined.

Protocol for use with the LightCycler® 480 Instrument II (Multiwell Plate 96 or 384)

Setup					
Block Type		Reaction Volume [µl]			
96 (384)		20 (10)			
Detection Format	Excitation Filter		Emission Filter		
For example: 3 Color Hydrolysis Probe					
FAM	465		510		
Yellow 555 (or Red 610)	533 (533)		580 (610)		
Cy5	618		660		
For new customized hydrolysis probes detection formats, set for all selected filters in the “Selected Filter Combination List” (under Tools), the following values:					
Melt Factor	1				
Quant Factor	10				
Max Integration Time (Sec)	2				
Programs					
Program Name	Cycles	Analysis Mode			
Pre-Incubation	1	None			
Amplification	45 ⁽¹⁾	Quantification			
Cooling	1	None			
Temperature Targets					
	Target [°C]	Acquisition Mode	Hold [hh:mm:ss]	Ramp Rate [°C/s]	Acquisitions [per °C]
Pre-Incubation	95	None	00:00:30	4.4 (4.8)	–
Amplification	95	None	00:00:05	4.4 (4.8)	–
	60 ⁽²⁾	Single	00:00:30	2.2 (2.5)	–
Cooling	40	None	00:00:30	2.2 (2.5)	–

⁽¹⁾ 45 cycles is suitable for most assays. If the assay shows steep amplification curves and early crossing points, 40 cycles should be sufficient. Reducing the number of cycles will improve uniformity of the product and reduce the time required for the assay.

⁽²⁾ The LightCycler® Multiplex DNA Master includes an aptamer for hot start. Anneal/Extend temperatures significantly below 60°C may not be well tolerated because aptamers will inhibit DNA polymerase activity at low temperatures.

Color Compensation Protocol for the LightCycler® 480 Instrument II

For a multicolor, multiplex amplification, the application of a color compensation file is necessary to compensate for optical crosstalk between the different channels.

For the LightCycler® 480 Instrument II, an instrument-specific color compensation file is mandatory for multicolor experiments, and a color compensation object can be generated by performing the following experiment.

The LightCycler® 480 Instrument II protocol contains the following programs:

- **Pre-Incubation** for activation of DNA polymerase and denaturation of the DNA
- **Amplification** of the DNA
- **Temperature Gradient Step** to create the Color Compensation file
- **Cooling** of the thermal block

i For details on how to program the experimental protocol, see the *LightCycler® 480 Software Operator's Manual, Version 1.5*.

2. How to Use this Product

The following table shows the qPCR parameters that must be programmed for a LightCycler® 480 System Color Compensation file run with a LightCycler® 480 Multiwell Plate 96, white.

Setup					
Block Type		Reaction Volume [µl]			
96		20			
Detection Format					
Customized (see previous section)					
Programs					
Program Name	Cycles	Analysis Mode			
Pre-Incubation	1	None			
Amplification	45	Quantification			
Temperature Gradient Step	1	Color Compensation			
Cooling	1	None			
Temperature Targets					
	Target [°C]	Acquisition Mode	Hold [hh:mm:ss]	Ramp Rate [°C/s]	Acquisitions [per °C]
Pre-Incubation	95	None	00:00:30	4.4	-
Amplification	95	None	00:00:05	4.4	-
	60	Single	00:00:30	2.2	-
Temperature Gradient Step	95	None	00:00:10	4.4	-
	40	None	00:00:10	2.2	-
	95	Continuous	-	-	5
Cooling	40	None	00:00:30	2.2	-

Preparation of the Color Compensation Run

Prepare the calibrator qPCR mix for more than one reaction; multiply the amount in the columns below by the number of reactions (minimum of 3 to 5 replicates). Note that there will be a slight loss of liquid during the pipetting steps. Please calculate extra volume of the qPCR mix by adding at least one additional reaction volume. In order to ensure accuracy, we do not recommend pipetting volumes below 1 µl when adding the individual reagents.

Reagent	1x Buffer [µl]	1x for Each Dye [µl]	1x DNA Process Control [µl]
qPCR Reaction Mix, 5x conc. (Vial 3)	4.0	4.0	4.0
DNA Process Control Detection Assay, 20x (Vial 4)	-	-	1.0
Detection mix for each dye (one dye per well)	-	X (depending on the assay)	-
Water, PCR Grade (Vial 5)	16.0 µl	Y (depending on the assay)	10.0
Template, such as DNA or positive samples eluates	-	5.0	-
DNA Process Control eluate	-	-	5.0
Total Volume	20.0	20.0	20.0

- 1 Pipette the replicates of each different calibrator mix into a precooled LightCycler® 480 Multiwell Plate 96.
- 2 Seal the LightCycler® 480 Multiwell Plate using a sealing foil.
- 3 Place the multiwell plate in a standard swinging-bucket centrifuge containing a rotor for multiwell plates with suitable adaptors. Balance it with a suitable counterweight (e.g., another plate), and centrifuge for 2 minutes at $1,500 \times g$.
- 4 Load the multiwell plate into the LightCycler® 480 Instrument II and start the program.

Create Color Compensation Object

When the experiment is finished, open the **Sample Editor**. Select workflow **Color Comp** and designate the used positions of the Multiwell Plate as “Water” for Buffer replicates, and the appropriate dyes respectively (e.g., FAM, Yellow 555 or Red 610, Cy5 for the example mentioned above).

Open the **Analysis** module “Color Compensation”, click on **Calculate**, and save the color compensation object in the database.

From this point on, you can apply this “CC Object” to each dual-, triple- or quadruple-color experiment performed with FAM, Yellow 555, Red 610, and Cy5 on the same instrument.

LightCycler® 96 Instrument Protocol

The following procedure is optimized for use with the LightCycler® 96 System. Program the LightCycler® 96 Instrument before preparing the reaction mixes.

i For details on how to program the experimental protocol, see the *LightCycler® 96 Operator's Manual*.

Run Editor				
Detection Format		Reaction Volume (μl)		
Select Cy5 and any other dye used in your assays.		20		
Programs				
	Temp. [°C]	Ramp [°C/s]	Duration [s]	Acquisition Mode
Pre-Incubation	95	4.4	30	None
2-Step Amplification	No. of Cycles: 45 ⁽¹⁾			
	95	4.4	5	None
	60 ⁽²⁾	2.2	30	Single

⁽¹⁾ 45 cycles is suitable for most assays. If the assay shows steep amplification curves and early crossing points, 40 cycles should be sufficient. Reducing the number of cycles will improve uniformity of the product and reduce the time required for the assay.

⁽²⁾ The LightCycler® Multiplex DNA Master includes an aptamer for hot start. Anneal/Extend temperatures significantly below 60°C may not be well tolerated because aptamers will inhibit DNA polymerase activity at low temperatures.

Color Compensation Protocol for the LightCycler® 96 Instrument

The LightCycler® 96 Instrument does not require the creation of a color compensation object.

3. Results

The DNA Process Control Kit was used according to the protocols described above. Plasma samples (180 µl) were processed on the MagNA Pure 96 Instrument with ± 20 µl of the DNA Process Control working solution. Heat-inactivated CMV virus was also processed on the MagNA Pure 96 Instrument, then amplified in either monoplex or duplex qPCR amplifications. As demonstrated below, Cqs generated from monoplex and duplex qPCR amplifications yielded very similar results, indicating that the presence and amplification of the DNA Process Control does not affect target amplification performance.




Fig. 1: Duplex/monoplex performance. Amplification curves of a CMV target and the DNA Process Control assay. The DNA Process Control (DPC) was added (±) to plasma samples and processed with the MagNA Pure 96 DNA and Viral NA Small Volume Kit using the Pathogen Universal 200 V3.0 protocol.

Use with Different Sample Types

The DNA Process Control Kit was used in a panel of 5 different sample materials (EDTA-plasma, serum, whole blood, swab, and urine) according to the same standard protocol as described above (MagNA Pure 96 DNA and Viral NA Small Volume Kit using the Pathogen Universal 200 protocol with 100 µl eluate). Eluates were post-spiked with Human Reference cDNA and qPCR amplified in a duplex with the DNA Process Control. All NTCs are clean and the DNA Process Control assay does not produce false positives in any of the sample materials. The DNA Process Control assay performs comparably in all different sample types.




Fig. 2: Amplification curves of B2M2 target and the DNA Process Control (DPC) in various sample types. Twenty microliters DPC added to 180 µl of various sample types, then purified on the MagNA Pure 96 with the Viral NA Small Volume Kit using the Pathogen Universal 200 V3.0 protocol.

Interpretation of Data

The robust workflow described above ensures that a constant amount of DNA Process Control is added to the sample material. However, the detection of Cq values may vary depending on the overall workflow setup (used sample material, purification instrument and protocol, detection instrument, multiplexing parameters, researcher, *etc.*). In general, expect the Cq value to vary in a range of $\sim \pm 2 - 3$ Cqs for one particular setup.

Using the standard procedures with different nucleic acid purification protocols on different instruments with different sample materials, the DNA Process Control resulted in Cq values of $\sim 30 - 35$ (100 μ l eluate and 5 μ l eluate in a 20 μ l qPCR detection reaction with the LightCycler® Multiplex DNA Master on a LightCycler® 480 Instrument II). For one such particular workflow, the following table shows an exemplary interpretation guideline.

Target (Channel 465 – 510)	DNA Process Control (Channel 618 – 660)	Results/Target Status
No Cq	Valid Cq ($\sim 30 - 35$)	Valid/target negative
Cq <40	Valid Cq or no Cq	Valid/target positive
No Cq	No valid Cq <ul style="list-style-type: none"> ▪ For example, outside of $\sim 30 - 35$ ▪ No sigmoid curve 	Invalid/target unknown

⚠ For any result assessments, do not solely rely on Cq callings. Ensure that the amplification curves are sigmoid in shape and that the Cq values are reasonable. The qPCR should be repeated if Cq values show an uncertain Character Code and Description of the status in the result screen (LightCycler® 480 Software, Version 1.5).

⚠ For a valid test batch (including negative controls with a negative result and positive controls with a positive result), verify each individual sample for its result, including the DNA Process Control result.

Negative Control/ No Template Control (NTC)

The assay result for a negative control, measured in the appropriate channel, always needs to be negative. If the result is “positive”, all samples results controlled by the corresponding negative control are invalid because of potential contamination of the reagents. In the case of a “positive” result for the negative control, sample preparation and qPCR with the appropriate detector must be repeated.

Positive Controls

The assay result for the positive controls, measured in the appropriate channel, always needs to be positive. If the result is “negative”, all samples results controlled by the corresponding positive control are invalid, and the respective qPCR with the appropriate detector must be repeated.

Sample Result

Check if the results of the negative and positive controls in the run are valid, and interpret the sample results for each target and DNA Process Control according to the table above (Interpretation of Data).

DNA Process Control/Internal Control

If samples results are “negative” for a target, the simultaneous DNA Process Control measured in channel 618 – 660 always needs to be positive to prove that there actually was sample material in the reaction and that the qPCR was not inhibited. For positive samples with a high amount of target DNA, the extraction control may be “negative” because of the competition of the two reactions.

4. Troubleshooting

Observation	Possible cause	Recommendation
Fluorescence intensity varies.	Some of the reagent is still in the upper part of the microwell or an air bubble is trapped in the microwell.	Repeat centrifugation, but allow sufficient centrifugation time so that all reagent is at the bottom of the microwell and air bubbles are expelled.
	Skin oils or dirt on the surface of the microwell plate.	Always wear gloves when handling the multiwell plate.
Fluorescence intensity is very low.	Low concentration or deterioration of dyes in the reaction mixtures because dye was not stored properly.	Keep dye-labeled reagents protected from light. Store the reagents at -15 to -25°C and avoid repeated freezing and thawing.
	Poor qPCR efficiency (reaction conditions not optimized).	Check annealing temperature of primers and probes.
		Check experimental protocol.
		Optimize annealing temperature in the qPCR reaction.
	Always run a positive control along with your samples.	
	DNA is degraded during isolation or improper storage.	If possible, check DNA quality. Check DNA with an established qPCR primer when available.
Pipetting errors and omitted reagents.	Check for missing reagents. Check the pipetting procedure.	
Impure sample material inhibits reaction.	Dilute sample 1:10 and repeat the analysis. Repurify the nucleic acids to ensure removal of inhibitory agents.	
Negative control sample gives a positive signal.	Contamination	Remake all critical reaction mixes. Use separate qPCR setup working areas.

5. Additional Information on this Product

5.1. Test Principle

Background Information

Two major outcomes can occur during a typical detection workflow for DNA targets. The sample can be either positive or negative. Whereas a positive detection also verifies the functionality of workflow components, a negative test outcome could be due either to a true negative sample or to a failure of critical workflow components and thus be a false negative. To rule out the possibility of workflow failures, a positive control that is similar to the target material but inert to the target detection is supplied in this kit and can be used in the workflow.

In order to prevent any cross-reaction with sample-derived nucleic acids or target-specific detection systems, the DNA Process Control Assay amplicon has been designed to have no significant homologies to any other known sequence. The DNA Process Control Detection Assay primers and probe are added in a low concentration to further lower any possible competition effects in multiplex qPCR amplifications.

The DNA Process Control concentrate is formulated to achieve a robust C_q value within one specified workflow. Different sample materials and workflows may require adjusted dilution of the DNA Process Control concentrate.

The DNA Process Control Kit is designed for the detection of an endogenous heterologous control (DNA Process Control) during a sample purification procedure when using:

- either a MagNA Pure System together with a dedicated reagent kit (for automated isolation),
- or a High Pure Nucleic Acid Isolation Kit (for manual isolation).

How this Product Works

The DNA Process Control Kit provides a non-competitive internal control to monitor nucleic acid purification and detection processes in order to prevent false-negative results. The DNA Process Control working solution contains a constant amount of nuclease-resistant DNA⁽¹⁾ that can be added to a large variety of sample materials without impairing the purification of other sample intrinsic nucleic acid⁽²⁾, such as from DNA viruses.

⁽¹⁾ The strictly monitored production and quality processes in conjunction with the easy-to-use workflow ensure that constant amounts of the control DNA are added to each and every sample.

⁽²⁾ The DNA Process Control has a unique and completely artificial sequence that makes primer competition with any other target parameters less likely. The primer and the probe of the DNA Process Control are designed to only amplify and detect the unique sequence. The Cy5-labeled probe ensures that other target-specific assays with common dyes such as FAM, Yellow 555, or Red 610 can be combined in multiplex setups.

Test Principle

The control is added to the sample material and co-purified with all other sample endogenous nucleic acids. After purification, the DNA Process Control as well as any target-specific parameter is detected in a qPCR reaction.

- In monoplex reactions in single wells, the DNA Process Control must be positive.
- In multiplex reactions, the DNA Process Control detection must be positive in all target-negative samples. In target-positive samples, the DNA Process Control may be out-competed by the target-specific assay and a negative control result is allowed.

Any failure of the DNA Process Control Detection Assay in target-negative samples indicates an erroneous purification/detection workflow and the sample must be retested.

5. Additional Information on this Product

The key steps in the process are:

- ① DNA Process Control working solution is added to the sample material.

- ② The protective coat of the internal control is lysed during the extraction process, thus enabling co-purification of the control DNA with other sample endogenous nucleic acids.

- ③ The sample eluates are used for qPCR reactions.
 - The DNA Process Control Detection Assay specifically detects the internal control.
 - Laboratory-developed assays assess the status of other targets in the sample material.Since the DNA Process Control is added to all samples, its successful detection proves the correct sample processing as well as the functionality of the generic detection reagents.



5.2. Quality Control

Each component of the three DNA Process Control Kits is function tested using the MagNA Pure 96 Instrument for the nucleic acid purification and the LightCycler® 480 Instrument II for qPCR.

6. Supplementary Information

6.1. Conventions

To make information consistent and easier to read, the following text conventions and symbols are used in this document to highlight important information:

Text convention and symbols	
	<i>Information Note: Additional information about the current topic or procedure.</i>
	Important Note: Information critical to the success of the current procedure or use of the product.
① ② ③ etc.	Stages in a process that usually occur in the order listed.
❶ ❷ ❸ etc.	Steps in a procedure that must be performed in the order listed.
* (Asterisk)	The Asterisk denotes a product available from Roche Diagnostics.

6.2. Changes to previous version

Layout changes.
Editorial changes.

6.3. Ordering Information

Roche offers a large selection of reagents and systems for life science research. For a complete overview of related products and manuals, please visit and bookmark our homepage lifescience.roche.com.

Product	Pack Size	Cat. No.
Accessories general (hardware)		
LightCycler® 480 Block Kit 96 Silver	1 block kit	05 015 219 001
LightCycler® 480 Block Kit 384 Silver	1 block kit	05 015 197 001
Accessories software		
LightCycler® 480 Software, Version 1.5	1 software package	04 994 884 001
Consumables		
LightCycler® 480 Multiwell Plate 96, white	5 x 10 plates	04 729 692 001
LightCycler® 480 Multiwell Plate 384, white	5 x 10 plates	04 729 749 001
LightCycler® 480 Sealing Foil	50 foils	04 729 757 001
MagNA Pure 96 Internal Control Tube		06 374 905 001
Instruments		
LightCycler® 480 Instrument II	1 instrument	05 015 278 001
	1 instrument	05 015 243 001
MagNA Pure LC 2.0 Instrument	1 instrument	05 197 686 001
MagNA Pure Compact Instrument	1 instrument with integrated PC, touchscreen monitor and bar-code reader	03 731 146 001
MagNA Pure 96 Instrument		06 541 089 001
LightCycler® 96 Instrument	1 instrument	05 815 916 001
Reagents , kits		
Bovine Serum Albumin	1 ml, 20 mg/ml	10 711 454 001
High Pure Viral Nucleic Acid Kit	1 kit, up to 100 isolations	11 858 874 001
MagNA Pure LC Total Nucleic Acid Isolation Kit	1 kit, up to 192 isolations	03 038 505 001
MagNA Pure Compact Nucleic Acid Isolation Kit I	1 kit, 4 x 8 sealed cartridges, 32 isolations	03 730 964 001
MagNA Pure 96 DNA and Viral NA Small Volume Kit		06 543 588 001
MagNA Pure 96 DNA and Viral NA Large Volume Kit		06 374 891 001
LightCycler® Multiplex DNA Master	1 kit, 1,000 reactions of 20 µl final volume each	07 339 577 001
	1 kit, 200 reactions of 20 µl final volume each	07 339 585 001

6.4. Trademarks

HIGH PURE, LIGHTCYCLER and MAGNA PURE are trademarks of Roche.
All third party product names and trademarks are the property of their respective owners.

6.5. License Disclaimer

For patent license limitations for individual products please refer to: <http://technical-support.roche.com>.

6.6. Regulatory Disclaimer

For general laboratory use.
This product is not available in all territories due to different national regulations.
This document is not intended for use in the USA.

6.7. Safety Data Sheet

Please follow the instructions in the Safety Data Sheet (SDS).

6.8. Contact and Support

If you have questions or experience problems with this or any Roche product for Life Science, please contact our Technical Support staff. Our scientists are committed to providing rapid and effective help.
Please also contact us if you have suggestions for enhancing Roche product performance or using our products in new or specialized ways. Such customer information has repeatedly proven invaluable to the research community worldwide.

To ask questions, solve problems, suggest enhancements or report new applications, please visit our **Online Technical Support** Site.

Visit **lifescience.roche.com**, to download or request copies of the following **Materials**:

- Instructions for Use
- Safety Data Sheets
- Certificates of Analysis
- Information Material

To call, write, fax, or email us, visit **lifescience.roche.com** and select your home country to display country-specific contact information.

