

## **cobas<sup>®</sup> HEV**

---

For in vitro diagnostic use



<b>cobas<sup>®</sup> HEV – 192</b>	P/N: 09040986190
<b>cobas<sup>®</sup> HEV Control Kit</b>	P/N: 09040889190
<b>cobas<sup>®</sup> NHP Negative Control Kit</b>	P/N: 09051554190
<b>cobas<sup>®</sup> omni MGP Reagent</b>	P/N: 06997546190
<b>cobas<sup>®</sup> omni Specimen Diluent</b>	P/N: 06997511190
<b>cobas<sup>®</sup> omni Lysis Reagent</b>	P/N: 06997538190
<b>cobas<sup>®</sup> omni Wash Reagent</b>	P/N: 06997503190

# Table of contents

<b>Intended use .....</b>	<b>4</b>
<b>Summary and explanation of the test.....</b>	<b>4</b>
<b>Reagents and materials.....</b>	<b>7</b>
cobas® HEV reagents and controls.....	7
cobas® omni reagents for sample preparation .....	10
Reagent storage and handling requirements.....	11
Reagent handling requirements for the cobas® 5800 System.....	11
Reagent handling requirements for the cobas® 6800/8800 Systems .....	12
Additional materials required for the cobas® 5800 System .....	13
Additional materials required for the cobas® 6800/8800 Systems.....	14
Instrumentation and software required .....	14
<b>Precautions and handling requirements .....</b>	<b>15</b>
Warnings and precautions .....	15
Reagent handling .....	15
Good laboratory practice.....	16
<b>Sample collection, transport, storage, and pooling .....</b>	<b>16</b>
Living donor samples.....	16
<b>Instructions for use.....</b>	<b>20</b>
Automated sample pipetting and pooling (optional).....	20
Procedural notes .....	20
Running cobas® HEV on the cobas® 5800 System .....	20
Running cobas® HEV on the cobas® 6800/8800 Systems.....	21
<b>Results .....</b>	<b>22</b>
Quality control and validity of results on the cobas® 5800 System.....	22
Control results on the cobas® 5800 System .....	22
Quality control and validity of results on the cobas® 6800/8800 Systems.....	22
Control flags on the cobas® 6800/8800 Systems .....	23
Interpretation of results .....	23

Additional information for interpretation of results on the cobas® 5800 System .....	23
Repeat testing of individual sample(s) .....	24
Procedural limitations.....	24
System equivalency/system comparison .....	24
<b>Non-clinical performance evaluation performed on the cobas® 6800/8800 Systems ....</b>	<b>25</b>
Key performance characteristics - Living donor samples .....	25
Limit of Detection (LoD) .....	25
Reproducibility.....	26
Genotype verification .....	27
Analytical specificity .....	28
Analytical specificity – interfering substances.....	29
Correlation.....	30
Whole system failure .....	30
<b>Additional information .....</b>	<b>31</b>
Key test features .....	31
Symbols.....	32
Technical support.....	33
Manufacturer and importer .....	33
Trademarks and patents .....	33
Copyright.....	33
References.....	34
Document revision.....	37

## Intended use

The cobas® HEV test for use on cobas® 5800/6800/8800 Systems is a qualitative in vitro nucleic acid amplification test for the direct detection of hepatitis E virus (HEV) RNA (genotypes 1-4) in human plasma.

This test is intended for use to screen donor samples for HEV RNA in plasma samples from individual human donors, including donors of whole blood, blood components (red cells, platelets, and plasma), and other living donors. Plasma from all donors may be screened as individual samples. For donations of whole blood and blood components, plasma samples may be tested individually or plasma may be tested in pools comprised of aliquots of individual samples.

This test is not intended for use on samples of cord blood.

This test is not intended for use as an aid in diagnosis for HEV.

## Summary and explanation of the test

### Background: Screening of blood for transfusion-transmitted viral infections

Hepatitis E virus (HEV), a small, non-enveloped, RNA virus belonging to the Hepevirus genus (family Hepeviridae), is a human pathogen with a worldwide distribution.<sup>1</sup> The virus consists of an icosahedral particle that encloses a positive-sense, single-stranded RNA genome of 7.2kb.<sup>2</sup> Four major HEV genotypes, representing a single serotype, have been identified in humans and animals, including domestic pigs, wild boar, deer, and rodents.<sup>1,3,4</sup>

Molecular characterization of various HEV strains circulating among humans and animals has led to the recognition of four major genotypes.<sup>1</sup> Genotype 1, which occurs mainly in Asia, and genotype 2, which occurs in Africa and Mexico, are restricted to humans and transmitted via contaminated water in developing countries.<sup>1,5</sup> Genotypes 3 and 4 infect humans, pigs, and other mammalian species and cause sporadic cases of autochthonous HEV in both developing and developed countries.<sup>6</sup> Genotype 3 is the only genotype currently identified as the cause of autochthonous infection in the U.S.<sup>7</sup> and is the cause of the vast majority of infections in Europe, New Zealand, and North America.<sup>1,8-12</sup> Genotypes 3 and 4 are both present in Japan.<sup>1</sup> HEV genotypes 1, 3, and 4 are endemic to China.<sup>1</sup> Acute Hepatitis E is more common than Hepatitis A in China, France, the U.K., and Japan.<sup>1</sup>

The main mode of HEV transmission is the fecal-oral route through contaminated drinking water,<sup>1</sup> although food-borne transmission from consumption of undercooked or raw pork, organ meat, or shellfish, as well as zoonotic transmission as a result of contact with infected swine, domestic, or wild animals have been reported.<sup>7,13-16</sup> The full range of reservoirs for HEV is unknown.<sup>1</sup>

HEV infection usually causes a mild or subclinical infection with a self-limiting illness that lasts 4 to 6 weeks.<sup>1,8-11,17,18</sup> The symptoms are very similar to those of other forms of viral hepatitis infection, particularly Hepatitis A, with fatigue, jaundice, fever, malaise, nausea, vomiting, anorexia, and abdominal pain.<sup>1</sup> Patients often present with elevated alanine transaminase levels (~1,500 IU/L) and many present with jaundice.<sup>1</sup> HEV infections may occasionally be more severe and result in fulminant hepatic failure, particularly in pregnant woman, where the mortality rate can reach 10% to 25%; infants and children under 2 years of age, individuals with underlying liver disease (e.g., cirrhosis), and immunocompromised persons.<sup>1,19-23</sup> HEV infection causes more than 3 million symptomatic cases of acute Hepatitis E worldwide each year, which results in approximately 70,000 deaths annually.<sup>24</sup> The overall mortality rate ranges from 0.2% to 4.0%.<sup>1</sup> Most deaths from genotype 3 (HEV3) infection result from acute or subacute liver failure in patients with pre-existing liver disease, such as alcohol-related liver disease.<sup>1,22,25,26</sup>

HEV3 causes chronic infection, including up to 60% of infected immunosuppressed individuals, about 10% of whom develop cirrhosis.<sup>1</sup> Chronic infection is defined as persistent HEV RNA in serum or stool for 6 months or more.<sup>1</sup> Most

cases occur in solid-organ transplant recipients, although infection in individuals with hematologic disorders receiving transfusions and chemotherapy and a few individuals with HIV have also been reported.<sup>27-33</sup> Chronic infection has not been reported with HEV1 or HEV2.<sup>1</sup>

HEV infection has also been associated with neurologic syndromes, including Guillain-Barré syndrome, Bell's palsy, acute transverse myelitis, acute meningoencephalitis, ataxia, and encephalitis; the neurologic symptoms typically resolve in patients who clear the virus.<sup>1</sup> Membranoproliferative and membranous glomerulonephritis, acute pancreatitis, and severe thrombocytopenia have been reported during acute HEV infection, although the pathophysiologic mechanisms and causal relation, if any, have not been established.<sup>1</sup> Ribavirin therapy has been shown to be an effective treatment for acute severe HEV3 infection, and transplant recipients with chronic HEV infection are typically treated with reduction of immunosuppression (especially drugs that target T cells), interferon- $\alpha$ , and ribavirin.<sup>1</sup>

### Rationale for NAT testing

Like other hepatitises, HEV can be transmitted via transfusion of blood or blood products. Post-transfusion Hepatitis E has been reported in many countries.<sup>1,33-38</sup> The seropositive rate for HEV among the world's blood donors has been reported to vary from 0.4%–20.6%.<sup>39-47</sup> Asymptomatic HEV infections occur at a high rate around the world and, due to the prevalence of the virus, many blood donors may be infected and transmit the virus to recipients of their blood products. For example, a recent study of British blood donors showed 11% of donor sera was HEV IgG reactive, indicating past infection, and 0.7% of donor sera was IgM reactive, indicating acute infection.<sup>45</sup> In addition, 0.7% of plasma minipools from English donors contained HEV RNA.<sup>46</sup> A study of Chinese blood donors revealed similar results: 32.6% of donor sera was IgG reactive; 0.94% of donor sera was IgM reactive; and 0.07% of donations demonstrated HEV viremia.<sup>47</sup> A global investigation of plasma fractionation pools reported 10% of pools tested were HEV-RNA positive.<sup>47,48</sup>

### Explanation of the test

The **cobas**® HEV test is a qualitative PCR test for the detection of HEV RNA that is run on the **cobas**® 5800/6800/8800 Systems. The **cobas**® HEV test enables the simultaneous detection of HEV RNA and the internal control in a single test of an infected, individual donation or pooled plasma from individual donations.

### Principles of the procedure

The **cobas**® HEV test is based on real time PCR technology on a fully automated sample preparation (nucleic acid extraction and purification) followed by PCR amplification and detection system. The **cobas**® 5800 System consists of a single, integrated instrument. The **cobas**® 6800/8800 Systems consist of the sample supply module, the transfer module, the processing module, and the analytic module. Automated data management is performed by the **cobas**® 5800 or 6800/8800 System software which assigns test results for all tests as non-reactive, reactive, or invalid. When using the **cobas**® 5800 System, the **cobas**® **Synergy** software is recommended for reviewing results and printing reports, and it is required to send results to a Laboratory Information Management System (LIMS) or other result management system. When using the **cobas**® 6800/8800 System, results can be reviewed directly on the system screen, and printed as a report, or sent to a Laboratory Information Management System (LIMS) or other result management system.

Samples can either be tested individually or tested in pools consisting of multiple samples.

When testing on **cobas**® 5800 System, use of the **cobas**® **Synergy** software is required even if pooling is not performed. If pooling in a pre-analytical step is performed, **cobas**® **Synergy** software with the Hamilton MICROLAB® STAR/STARlet IVD may be used.

When testing on the **cobas**® 6800/8800 Systems, the **cobas** p 680 instrument, or **cobas**® **Synergy** software with the Hamilton

MICROLAB® STAR/STARlet IVD, may optionally be used in a pre-analytical step if pooling is to be performed.

Nucleic acid from the sample and added armored RNA internal control (IC) molecules (which serve as the sample preparation and amplification/detection process control) is simultaneously extracted. In addition the test utilizes two kit controls: a positive and a negative control. Viral nucleic acid is released by addition of proteinase and lysis reagent to the sample. The released nucleic acid binds to the silica surface of the added magnetic glass particles. Unbound substances and impurities, such as denatured protein, cellular debris, and potential PCR inhibitors (such as hemoglobin) are removed with subsequent wash reagent steps and purified nucleic acid is eluted from the glass particles with elution buffer at elevated temperature.

Selective amplification of target nucleic acid from the donor sample is achieved by the use of virus-specific forward and reverse primers which are selected from highly conserved regions of the viral nucleic acid. A thermostable DNA polymerase enzyme is used for both reverse-transcription and amplification. The master mix includes deoxyuridine triphosphate (dUTP), instead of deoxythymidine triphosphate (dTTP), which is incorporated into the newly synthesized DNA (amplicon).<sup>49-51</sup> Any contaminating amplicon from previous PCR run is destroyed by the AmpErase enzyme [uracil-N-glycosylase], which is included in the PCR master mix, when heated in the first thermal cycling step. However, newly formed amplicon is not destroyed since the AmpErase enzyme is inactivated once exposed to temperatures above 55°C.

The **cobas**® HEV master mix contains detection probes which are specific for HEV and IC nucleic acid. The specific HEV and IC detection probes are each labeled with one of two unique fluorescent dyes which act as a reporter. Each probe also has a second dye which acts as a quencher. The two reporter dyes are measured at defined wavelengths, thus permitting simultaneous detection and discrimination of the amplified HEV target and the IC.<sup>52,53</sup> The fluorescent signal of the intact probe is suppressed by the quencher dye. During the PCR amplification step, hybridization of the probes to the specific single-stranded DNA template results in cleavage by the 5' to 3' nuclease activity of the DNA polymerase resulting in separation of the reporter and quencher dyes and the generation of a fluorescent signal. With each PCR cycle, increasing amounts of cleaved probes are generated and the cumulative signal of the reporter dye is concomitantly increased. Since the two specific reporter dyes are measured at defined wavelengths, simultaneous detection and discrimination of the amplified HEV target and the IC are possible.

# Reagents and materials

## cobas® HEV reagents and controls

All unopened reagents and controls shall be stored as recommended in Table 1 to Table 4.

**Table 1** cobas® HEV test  
(HEV)



Store at 2–8°C

192 test cassette (P/N 09040986190)

Kit components	Reagent ingredients	Quantity per kit 192 tests
<b>Proteinase Solution (PASE)</b>	Tris buffer, < 0.05% EDTA, calcium chloride, calcium acetate, 8% (w/v) proteinase, glycerol  EUH210: Safety data sheets available on request. EUH208: Contains subtilisin from <i>Bacillus subtilis</i> . May produce an allergic reaction.	22.3 mL
<b>Internal Control (IC)</b>	Tris buffer, < 0.05% EDTA, < 0.001% internal control armored RNA construct (non-infectious RNA encapsulated in MS2 bacteriophage), < 0.002% Poly rA RNA (synthetic), < 0.1% sodium azide	21.2 mL
<b>Elution Buffer (EB)</b>	Tris buffer, 0.2% methyl-4 hydroxybenzoate	21.2 mL
<b>Master Mix Reagent 1 (MMX-R1)</b>	Manganese acetate, potassium hydroxide, < 0.1% sodium azide	7.5 mL
<b>HEV Master Mix Reagent 2 (HEV MMX-R2)</b>	Tricine buffer, potassium acetate, glycerol, 18% dimethyl sulfoxide, Tween 20, EDTA, < 0.06% dATP, dGTP, dCTP, < 0.14% dUTP, < 0.01% upstream and downstream HEV and internal control primers, < 0.01% fluorescent-labeled HEV probes, < 0.01% fluorescent-labeled internal control probe, < 0.01% oligonucleotide aptamer, < 0.01% Z05D DNA polymerase, < 0.01% AmpErase (uracil-N-glycosylase) enzyme (microbial), < 0.1% sodium azide	9.7 mL

**Table 2** cobas® HEV Control Kit  
(HEV (+) C)

Store at 2-8°C  
(P/N 09040889190)

Kit components	Reagent ingredients	Quantity per kit	Safety symbol and warning*
<b>HEV Positive Control (HEV (+) C)</b>	< 0.001% synthetic (armored) HEV RNA encapsulated in MS2 bacteriophage coat protein, normal human plasma; HEV RNA not detectable by PCR methods. 0.1% ProClin® 300 preservative**	16 mL (16 x 1 mL)	  <p><b>WARNING</b></p> <p>H317: May cause an allergic skin reaction. H412: Harmful to aquatic life with long lasting effects. P261: Avoid breathing mist or vapours. P273: Avoid release to the environment. P280: Wear protective gloves. P333 + P313: If skin irritation or rash occurs: Get medical advice/attention. P362 + P364: Take off contaminated clothing and wash it before reuse. P501: Dispose of contents/container to an approved waste disposal plant. 55965-84-9 Reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1)</p>

\* Product safety labeling primarily follows EU GHS guidance



\*\*Hazardous substance



**Table 3** cobas® NHP Negative Control Kit  
(NHP-NC)

Store at 2-8°C

(P/N 09051554190)


Kit components	Reagent ingredients	Quantity per kit	Safety symbol and warning*
<b>Normal Human Plasma Negative Control (NHP-NC)</b>	Normal human plasma; HEV RNA not detectable by PCR methods. < 0.1% ProClin® 300 preservative**	16 mL (16 x 1 mL)	  <b>WARNING</b> H317: May cause an allergic skin reaction. P261: Avoid breathing mists or vapours. P272: Contaminated work clothing should not be allowed out of the workplace. P280: Wear protective gloves. P333 + P313: If skin irritation or rash occurs: Get medical advice/attention. P362 + P364: Take off contaminated clothing and wash it before reuse. P501: Dispose of contents/container to an approved waste disposal plant. 55965-84-9 Reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1)

\* Product safety labeling primarily follows EU GHS guidance

\*\*Hazardous substance

## cobas omni reagents for sample preparation

**Table 4** cobas® omni reagents for sample preparation\*

Reagents	Reagent ingredients	Quantity per kit	Safety symbol and warning**
<b>cobas® omni MGP Reagent (MGP)</b> Store at 2–8°C (P/N 06997546190)	Magnetic glass particles, Tris buffer, 0.1% methyl-4 hydroxybenzoate, < 0.1% sodium azide	480 tests	Not applicable
<b>cobas® omni Specimen Diluent (SPEC DIL)</b> Store at 2–8°C (P/N 06997511190)	Tris buffer, 0.1% methyl-4 hydroxybenzoate, < 0.1% sodium azide	4 x 875 mL	Not applicable
<b>cobas® omni Lysis Reagent (LYS)</b> Store at 2–8°C (P/N 06997538190)	42.56% (w/w) guanidine thiocyanate***, 5% (w/v) polydocanol***, 2% (w/v) dithiothreitol, dihydro sodium citrate	4 x 875mL	 <p><b>DANGER</b></p> <p>H302: Harmful if swallowed.  H314: Causes severe skin burns and eye damage.  H411: Toxic to aquatic life with long lasting effects.  EUH032: Contact with acids liberates very toxic gas.  EUH071: Corrosive to the respiratory tract.  P273: Avoid release to the environment.  P280: Wear protective gloves/ protective clothing/ eye protection/ face protection/ hearing protection.  P303 + P361 + P353: IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.  P304 + P340 + P310: IF INHALED: Remove person to fresh air and keep comfortable for breathing. Immediately call a POISON CENTER/ doctor.  P305 + P351 + P338 + P310: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/ doctor.  P391: Collect spillage.  593-84-0 Guanidinium thiocyanate  9002-92-0 Polidocanol  3483-12-3 (R*,R*)-1,4-dimercaptobutane-2,3-diol</p>
<b>cobas® omni Wash Reagent (WASH)</b> Store at 15–30°C (P/N 06997503190)	Sodium citrate dihydrate, 0.1% methyl-4 hydroxybenzoate	4.2 L	Not applicable

\* These reagents are not included in the cobas® HEV test kit. See listing of additional materials required (Table 8 and Table 9).

\*\* Product safety labeling primarily follows EU GHS guidance

\*\*\*Hazardous substance

## Reagent storage and handling requirements

Reagents shall be stored and will be handled as specified in Table 5, Table 6 and Table 7.

When reagents are not loaded on the cobas® 5800/6800/8800 Systems, store them at the corresponding temperature specified in Table 5.

**Table 5** Reagent storage (when reagent is not on the system)

Reagent	Storage temperature
cobas® HEV – 192	2–8°C
cobas® HEV Control Kit	2–8°C
cobas® NHP Negative Control Kit	2–8°C
cobas® omni Lysis Reagent	2–8°C
cobas® omni MGP Reagent	2–8°C
cobas® omni Specimen Diluent	2–8°C
cobas® omni Wash Reagent	15–30°C

## Reagent handling requirements for the cobas® 5800 System

Reagents loaded onto the cobas® 5800 System are stored at appropriate temperatures and their expiration is monitored by the system. The system allows reagents to be used only if all of the conditions shown in Table 6 are met. The system automatically prevents use of expired reagents. Table 6 allows the user to understand the reagent handling conditions enforced by the cobas® 5800 System.

**Table 6** Reagent expiry conditions enforced by the cobas® 5800 System

Reagent	Kit expiration date	Open-kit stability	Number of runs for which this kit can be used	On-board stability (cumulative time on board outside refrigerator)
cobas® HEV – 192	Date not passed	90 days from first usage	Max 40 runs	Max 36 days
cobas® HEV Control Kit	Date not passed	Not applicable*	Not applicable	Max 10 hours
cobas® NHP Negative Control Kit	Date not passed	Not applicable*	Not applicable	Max 10 hours
cobas® omni Lysis Reagent	Date not passed	30 days from loading**	Not applicable	Not applicable
cobas® omni MGP Reagent	Date not passed	30 days from loading**	Not applicable	Not applicable
cobas® omni Specimen Diluent	Date not passed	30 days from loading**	Not applicable	Not applicable
cobas® omni Wash Reagent	Date not passed	30 days from loading**	Not applicable	Not applicable

\* Single use reagents

\*\*Time is measured from the first time that reagent is loaded onto the cobas® 5800 System.

## Reagent handling requirements for the cobas® 6800/8800 Systems

Reagents loaded onto the cobas® 6800/8800 Systems are stored at appropriate temperatures and their expiration is monitored by the system. The system allows reagents to be used only if all of the conditions shown in Table 7 are met. The system automatically prevents use of expired reagents. Table 7 allows the user to understand the reagent handling conditions enforced by the cobas® 6800/8800 Systems.

**Table 7** Reagent expiry conditions enforced by the cobas® 6800/8800 Systems

Reagent	Kit expiration date	Open-kit stability	Number of runs for which this kit can be used	On-board stability (cumulative time on board outside refrigerator)
cobas® HEV – 192	Date not passed	90 days from first usage	Max 40 runs	Max 40 hours
cobas® HEV Control Kit	Date not passed	Not applicable*	Not applicable	Max 10 hours
cobas® NHP Negative Control Kit	Date not passed	Not applicable*	Not applicable	Max 10 hours
cobas® omni Lysis Reagent	Date not passed	30 days from loading**	Not applicable	Not applicable
cobas® omni MGP Reagent	Date not passed	30 days from loading**	Not applicable	Not applicable
cobas® omni Specimen Diluent	Date not passed	30 days from loading**	Not applicable	Not applicable
cobas® omni Wash Reagent	Date not passed	30 days from loading**	Not applicable	Not applicable

\* Single use reagents

\*\*Time is measured from the first time that reagent is loaded onto the cobas® 6800/8800 Systems.

## Additional materials required for the cobas® 5800 System

**Table 8** Material and consumables for use on cobas® 5800 System

Material	P/N
<b>cobas® omni</b> Processing Plate 24	08413975001
<b>cobas® omni</b> Amplification Plate 24	08499853001
<b>cobas® omni</b> Liquid Waste Plate 24	08413983001
CORE Tips with Filter, 1 mL	04639642001
CORE Tips with Filter, 300 µL	07345607001
<b>cobas® omni</b> Liquid Waste Container	07094388001
<b>cobas® omni</b> Lysis Reagent	06997538190
<b>cobas® omni</b> MGP Reagent	06997546190
<b>cobas® omni</b> Specimen Diluent	06997511190
<b>cobas® omni</b> Wash Reagent	06997503190
Solid Waste Bag	07435967001
or	or
Solid Waste Bag with Insert	08030073001

## Additional materials required for the cobas® 6800/8800 Systems

**Table 9** Material and consumables for use on cobas® 6800/8800 Systems

Material	P/N
<b>cobas® omni</b> Processing Plate	05534917001
<b>cobas® omni</b> Amplification Plate	05534941001
<b>cobas® omni</b> Pipette Tips	05534925001
<b>cobas® omni</b> Liquid Waste Container	07094388001
<b>cobas® omni</b> Lysis Reagent	06997538190
<b>cobas® omni</b> MGP Reagent	06997546190
<b>cobas® omni</b> Specimen Diluent	06997511190
<b>cobas® omni</b> Wash Reagent	06997503190
Solid Waste Bag	07435967001
Solid Waste Container or Solid Waste Bag with Insert	07094361001 or 08030073001

## Instrumentation and software required

The cobas® HEV analysis package for the cobas® 5800 System shall be installed on the cobas® 5800 System. The x800 Data Manager software for the cobas® 5800 System will be provided with the system. The cobas® Synergy software shall be installed.

The cobas® 6800/8800 software and cobas® HEV analysis package shall be installed on the instrument(s). The Instrument Gateway (IG) server will be provided with the system. The cobas® Synergy software shall be installed, if applicable.

**Table 10** Instrumentation

<b>cobas® 6800/8800 Systems</b>	<b>P/N</b>
<b>cobas®</b> 6800 System (Option Moveable)	05524245001 and 06379672001
<b>cobas®</b> 6800 System (Fix)	05524245001 and 06379664001
<b>cobas®</b> 8800 System	05412722001
Sample Supply Module for <b>cobas®</b> 6800/8800 Systems	06301037001
<b>cobas® Synergy</b> software electronic license (for <b>cobas®</b> 5800 System only)	09311246001
<b>Options for pipetting and pooling</b>	<b>P/N</b>
<b>cobas® p</b> 680 instrument	06570577001
<b>cobas® Synergy</b> software electronic license ( <b>cobas®</b> 6800/8800 Systems) (Optional)	09311238001
Hamilton MICROLAB® STAR IVD	04640535001
Hamilton MICROLAB® STAR/Iet IVD	04872649001

Refer to the cobas® 5800 System User Assistance or the cobas® 6800/8800 Systems User Assistance for additional information. Refer to the cobas® p 680 instrument User Assistance, or to the cobas® Synergy software User Assistance, for additional information about primary and secondary sample tubes accepted on the instruments.

Note: Contact your local Roche representative for a detailed order list for sample racks, racks for clotted tips and rack trays accepted on the instruments.

# Precautions and handling requirements

## Warnings and precautions

As with any test procedure, good laboratory practice is essential to the proper performance of this assay. Due to the high sensitivity of this test, care should be taken to keep reagents and amplification mixtures free of contamination.

- For in vitro diagnostic use only.
- All samples should be handled as if infectious, using good laboratory procedures as outlined in Biosafety in Microbiological and Biomedical Laboratories and in the CLSI Document M29-A4.<sup>54,55</sup> Only personnel proficient in handling infectious materials and the use of the **cobas®** HEV test, the **cobas®** 5800/6800/8800 Systems, the **cobas® p 680** instrument (for **cobas®** 6800/8800 Systems) or the Hamilton MICROLAB® STAR/STARlet IVD with **cobas® Synergy** software should perform this procedure.
- All human-sourced materials should be considered potentially infectious and should be handled with universal precautions. If spillage occurs, immediately disinfect with a freshly prepared solution of 0.6% sodium hypochlorite in distilled or deionized water or follow appropriate site procedures.
- **cobas®** HEV Control Kit and **cobas®** NHP Negative Control Kit contain plasma derived from human blood. Testing of normal human plasma by PCR methods showed no detectable HEV RNA. No known test method can offer complete assurance that products derived from human blood will not transmit infectious agents.
- Do not freeze whole blood.
- The use of sterile disposable pipettes and nuclease-free pipette tips is recommended. Use only supplied or specified required consumables to ensure optimal test performance.
- Closely follow procedures and guidelines provided to ensure that the test is performed correctly. Any deviation from the procedures and guidelines may affect optimal test performance.
- Disruption of the cell-plasma interface or diffusion of material post-centrifugation may result in higher invalid rates.
- False positive results may occur if carryover of samples is not adequately controlled during sample handling and processing.
- Inform your local competent authority and manufacturer about any serious incidents which may occur when using this assay.

## Reagent handling

- Handle all reagents, controls, and samples according to good laboratory practice in order to prevent carryover of samples or controls.
- Before use, visually inspect each reagent cassette, diluent, lysis reagent, and wash reagent to ensure that there are no signs of leakage. If there is any evidence of leakage, do not use that material for testing.
- **cobas® omni** Lysis Reagent contains guanidine thiocyanate, a potentially hazardous chemical. Avoid contact of reagents with the skin, eyes, or mucous membranes. If contact does occur, immediately wash with generous amounts of water; otherwise, burns can occur.
- **cobas®** HEV test kits, **cobas® omni** MGP Reagent, and **cobas® omni** Specimen Diluent contain sodium azide as a preservative. Avoid contact of reagents with the skin, eyes, or mucous membranes. If contact does occur, immediately wash with generous amounts of water; otherwise, burns can occur. If these reagents are spilled, dilute with water before wiping dry.
- Do not allow **cobas® omni** Lysis Reagent, which contains guanidine thiocyanate, to contact sodium hypochlorite (bleach) solution. This mixture can produce a highly toxic gas.

- Safety Data Sheets (SDS) are available on request from your local Roche representative.
- Dispose of all materials that have come in contact with samples and reagents in accordance with country, state, and local regulations.

## Good laboratory practice

- Do not pipette by mouth.
- Do not eat, drink, or smoke in designated work areas.
- Wear laboratory gloves, laboratory coats, and eye protection when handling samples and reagents. Gloves must be changed between handling samples and cobas® HEV test kits and cobas® omni reagents to prevent contamination. Avoid contaminating gloves when handling samples and controls.
- Wash hands thoroughly after handling samples and kit reagents, and after removing the gloves.
- If spills occur on the cobas® 5800 instrument, follow the instructions in the cobas® 5800 System User Assistance to properly clean and decontaminate the surfaces of the instrument(s).
- Thoroughly clean and disinfect all laboratory work surfaces with a freshly prepared solution of 0.6% sodium hypochlorite in distilled or deionized water. Follow by wiping the surface with 70% ethanol.
- If spills occur on the cobas® 6800/8800 instruments, follow the instructions in the cobas® 6800/8800 Systems User Assistance to properly clean and decontaminate the surface of instrument(s).

## Sample collection, transport, storage, and pooling

**Note:** Handle all samples and controls as if they are capable of transmitting infectious agents.

Store all donor samples at specified temperatures.


Sample stability is affected by elevated temperatures.

## Living donor samples

- Plasma collected in EDTA, CPD, CPDA1, CP2D and 4% Sodium Citrate anticoagulant may be used with the cobas® HEV test. Follow the sample collection tube/bag manufacturer instructions for handling and centrifugation.
- Blood collected in EDTA anticoagulant, Becton-Dickinson EDTA Plasma Preparation Tubes (BD PPT™) or Greiner Vacuette® K2EDTA Plasma Gel Tubes may undergo additional centrifugation at 600 x g for 5 minutes prior to loading, optional pooling or retesting.
- Blood collected in EDTA anticoagulant may be stored for up to 12 days with the following conditions:
  - Samples must be centrifuged within 72 hours of draw.
  - For storage above 8°C, samples may be stored for 72 hours at up to 25°C, and up to 30°C for 24 hours during the 72 hours.


Other than noted above, samples are stored at 2-8°C. In addition, plasma separated from the cells may be stored for up to 12 months at ≤-18°C with three freeze/thaw cycles. Refer to Figure 1.



**Figure 1** Sample storage conditions in EDTA anticoagulant


- Blood collected in CPD, CPDA1, CP2D anticoagulant, Becton-Dickinson EDTA Plasma Preparation Tubes (BD PPT™) or Greiner Vacuette® K2EDTA Plasma Gel Tubes may be stored for up to 12 days with the following conditions:
  - Samples must be centrifuged within 72 hours of draw.
  - For storage above 8°C, samples may be stored for 72 hours at up to 25°C, and up to 30°C for 24 hours during the 72 hours.

Other than noted above, samples are stored at 2-8°C. In addition, plasma separated from the cells may be stored for up to 30 days at ≤ -18°C with three freeze/thaw cycles. Refer to Figure 2.

**Figure 2** Sample storage conditions

- Plasma collected in 4% Sodium Citrate anticoagulant may be stored for up to 30 days at 2-8°C.
  - For storage above 8°C, samples may be stored for 72 hours at up to 25°C, and up to 30°C for 24 hours during the 72 hours.

Other than noted above, samples are stored at 2-8°C. In addition, plasma separated from the cells may be stored for up to 12 months at ≤ -18°C with three freeze/thaw cycles. Refer to Figure 3.

**Figure 3** Sample storage conditions in 4% sodium citrate anticoagulant

- If samples are to be shipped, they should be packaged and labeled in compliance with applicable country and/or international regulations covering the transport of samples and etiologic agents.

# Instructions for use

## Automated sample pipetting and pooling (optional)

Either the **cobas® p 680** instrument, or **cobas® Synergy** software with the Hamilton MICROLAB® STAR/STARlet can be used as an optional instrument with the **cobas® 6800/8800** Systems for automated pipetting and pooling of aliquots of multiple primary samples into one pooled sample.

**cobas® Synergy** software with the Hamilton MICROLAB® STAR/STARlet IVD shall be used as an accessory to the **cobas® 5800** System for automated pipetting and pooling of aliquots of multiple primary samples into one pooled sample.

Refer to the **cobas® p 680** instrument User Assistance or to the **cobas® Synergy** software User Assistance for more information.

## Procedural notes

- Do not use **cobas®** HEV test reagents, **cobas®** HEV Control Kit, **cobas®** NHP Negative Control Kit or **cobas® omni** reagents after their expiry dates.
- Do not reuse consumables. They are for one-time use only.
- Refer to the **cobas® 5800** System User Assistance for proper maintenance of instruments.
- Refer to the **cobas® 6800/8800** Systems User Assistance for proper maintenance of instruments.

## Running cobas® HEV on the cobas® 5800 System

The test procedure is described in detail in the **cobas® 5800** System User Assistance. Figure 4 below summarizes the procedure. Refer to the **cobas® Synergy** software User Assistance as applicable for details on optional pooling procedures.


**Figure 4** **cobas®** HEV test procedure on the **cobas® 5800** System

1	Pipetting and pooling
2	Loading sample racks onto the system <ul style="list-style-type: none"> <li>• Load sample racks onto the system</li> <li>• Order tests manually if no LIS orders are available</li> </ul>
3	Refill reagent and consumables as prompted by the system <ul style="list-style-type: none"> <li>• Load test specific reagent cassette(s)</li> <li>• Load control mini racks</li> <li>• Load processing tips</li> <li>• Load elution tips</li> <li>• Load processing plates</li> <li>• Load liquid waste plates</li> <li>• Load amplification plates</li> <li>• Load MGP cassette</li> <li>• Refill specimen diluent</li> <li>• Refill lysis reagent</li> <li>• Refill wash reagent</li> </ul>
4	Start the run by choosing the Start button manually on the user interface. All subsequent runs will start automatically if not manually postponed.
5	Review results
6	Remove any sample tubes Clean up the instrument <ul style="list-style-type: none"> <li>• Empty reagent cassettes</li> <li>• Empty Control mini racks</li> <li>• Empty amplification plate drawer</li> <li>• Empty liquid waste</li> <li>• Empty solid waste</li> </ul>

## Running cobas® HEV on the cobas® 6800/8800 Systems

The test procedure is described in detail in the **cobas® 6800/8800 Systems User Assistance**; refer to the **cobas® p 680 instrument User Assistance** or to the **cobas® Synergy software User Assistance** as applicable for details on optional pooling procedures. Figure 5 below summarizes the procedure.

**Figure 5** cobas® HEV test procedure on the cobas® 6800/8800 Systems



## Results

The cobas® 5800 and cobas® 6800/8800 Systems automatically detect HEV RNA simultaneously for the samples and controls.

### Quality control and validity of results on the cobas® 5800 System

The cobas® 5800 System will be delivered with the default setting of controls (positive and negative) scheduled with every run, but can be configured to a less frequent control schedule, by a Roche service engineer or by contacting Roche customer technical support, based on laboratory procedures and/or local regulations.

- One negative control [(-) C] and one positive control, [HEV (+) C] are processed at least every 24 hours or with every new kit lot.
- In the cobas® 5800 System and/or report, check for flags and their associated results to ensure control validity.
- The associated samples are valid if no flags appear for all controls.

Invalidation of results is performed automatically by the cobas® 5800 System based on negative or positive control failures.

### Control results on the cobas® 5800 System

The results of the controls are shown in the cobas® 5800 software in the “Controls” app.

- Controls are marked with “Valid” in the column “Control result” if all Targets of the control are reported valid. Controls are marked with ‘Invalid’ in the column “Control result” if all or one Target of the control are reported invalid.
- Controls marked with ‘Invalid’ show a flag in the “Flags” column. More information on why the control is reported invalid including flag information will be shown in the detail view.
- If the positive control is invalid, repeat testing of the positive controls and all associated samples. If the negative control is invalid, repeat testing of all controls and all associated samples.

**Table 11** Control flags for negative and positive controls on the cobas® 5800 System

Negative Control	Flag	Control Result	Interpretation
(-) C	A flag is shown	Invalid	The entire batch is assigned invalid if the result for the (-) C is invalid.
Positive Control	Flag	Control Result	Interpretation
HEV (+) C	A flag is shown	Invalid	The entire batch is assigned invalid if the result for the HEV (+) C is invalid.

### Quality control and validity of results on the cobas® 6800/8800 Systems

- One negative control [(-) C] and one positive control [HEV (+) C] are processed with each batch.
- In the cobas® 6800/8800 software and/or report, check for flags and their associated results to ensure the batch validity.
- The batch is valid if no flags appear for both controls.

Invalidation of results is performed automatically by the cobas® 6800/8800 software based on negative and positive control failures.

## Control flags on the cobas® 6800/8800 Systems

**Table 12** Control flags for negative and positive controls

Negative Control	Flag	Result	Interpretation
(-) C	Q02	Invalid	The entire batch is assigned invalid if the result for the (-) C is invalid.
Positive Control	Flag	Result	Interpretation
HEV (+) C	Q02	Invalid	The entire batch is assigned invalid if the result for the HEV (+) C is invalid.

If the batch is invalid, repeat testing of the entire batch including samples and controls.

## Interpretation of results

For a valid batch, check each individual sample for flags in the **cobas®** 5800/6800/8800 Systems software and/or report. The result interpretation should be as follows:

- A valid batch may include both valid and invalid donor sample results dependent on flags obtained for the individual samples.
- Sample results are valid only if the respective positive control and the negative control of the corresponding batch are valid.

Two parameters are measured simultaneously for each sample: HEV and the internal control. Final sample results for the **cobas®** HEV test are reported by the software. In addition to the overall results, individual target results will be displayed in the **cobas®** 5800/6800/8800 Systems software and should be interpreted as follows:

**Table 13** Target results for individual target result interpretation

Target results	Interpretation
HEV Non-Reactive	No target signal detected for HEV and IC signal detected.
HEV Reactive	Target signal detected for HEV and IC signal may be or may not be detected.
Invalid	Target and internal control signal not detected.

If using the **cobas®** Synergy software, review of the final result calculation should be performed through the **cobas®** Synergy software.

## Additional information for interpretation of results on the cobas® 5800 System

The results of the samples are shown in the **cobas®** 5800 System. It is recommended to review results in the **cobas®** Synergy software.

- Samples associated with a valid control batch (as defined by your system control configuration) are shown as 'Valid' in the "Control result" column. Samples associated with a failed control batch are shown as 'Invalid' in the "Control result" column.
- If the associated controls of a sample result are invalid, a specific flag will be added to the sample result as follows:
  - Q05D : Result validation failure because of an invalid positive control
  - Q06D : Result validation failure because of an invalid negative control
- The values in "Results" column for individual sample target result should be interpreted as shown in Table 13 above.
  - The **cobas®** 5800 System will display individual target results. The overall result will be shown only in the result view of the **cobas®** Synergy software.
  - For more detailed information on sample results and flags refer to the **cobas®** 5800 System User Assistance.

## Repeat testing of individual sample(s)

Sample tubes with a final result of Invalid for the target require repeat testing.

An additional centrifugation at 600 x g for 5 minutes may help to reduce repeat invalid results for blood collected in EDTA anticoagulant, Becton-Dickinson EDTA Plasma Preparation Tubes (BD PPT™) or Greiner Vacuette® K2EDTA Plasma Gel Tubes.

## Procedural limitations

- The **cobas®** HEV test has been evaluated only for use in combination with the **cobas®** HEV Control Kit, **cobas®** NHP Negative Control Kit, **cobas®** **omni** MGP Reagent, **cobas®** **omni** Lysis Reagent, **cobas®** **omni** Specimen Diluent, and **cobas®** **omni** Wash Reagent for use on the **cobas®** 6800/8800 Systems.
- Reliable results depend on proper sample collection, storage and handling procedures.
- Do not use heparinized plasma with this test because heparin has been shown to inhibit PCR.
- Detection of HEV RNA is dependent on the number of virus particles present in the sample and may be affected by sample collection, storage and handling, patient factors (i.e., age, presence of symptoms), and/or stage of infection and pool size.
- Though rare, mutations within the highly conserved regions of a viral genome covered by the **cobas®** HEV test, may affect primers and/or probe binding resulting in the failure to detect presence of virus.
- Due to inherent differences between technologies, it is recommended that, prior to switching from one technology to the next, users perform method correlation studies in their laboratory to qualify technology differences. Users should follow their own specific policies/procedures.

## System equivalency/system comparison

System equivalency of the **cobas®** 5800 System with the **cobas®** 6800/8800 Systems was demonstrated via equivalency studies.

The results presented in these Instructions for Use are based on equivalent performance for all systems.



# Non-clinical performance evaluation performed on the cobas® 6800/8800 Systems

## Key performance characteristics - Living donor samples

### Limit of Detection (LoD)

#### WHO International Standard

The limit of detection (LoD) of the cobas® HEV test for HEV RNA was determined using the WHO International Standard for HEV (PEI code 6329/10).

For the WHO International standard, 3 independent dilution series of the viral standard were prepared with normal, virus-negative (HEV) human EDTA-plasma. Each dilution series was tested using 3 different lots of the cobas® HEV test kits with approximately 63 replicates per lot, for a total of approximately 189 replicates per concentration. For HEV virus, 95% PROBIT analysis (Table 14) and 50% PROBIT analysis (Table 15) on the data combined across dilution series and reagent lots was used to estimate the LoD, along with the lower and upper limit of 95% confidence intervals. The reactivity rates observed in the LoD studies for HEV are summarized in Table 16.

**Table 14** Results of 95% PROBIT analysis on LoD data collected with viral standard in EDTA plasma

Analyte	Measuring units	LoD	Lower 95% confidence limit	Upper 95% confidence limit
HEV	IU/mL	18.6	15.9	22.6

**Table 15** Results of 50% PROBIT analysis on LoD data collected with viral standard in EDTA plasma

Analyte	Measuring units	LoD	Lower 95% confidence limit	Upper 95% confidence limit
HEV	IU/mL	3.9	3.4	4.3

**Table 16** Reactivity rates summary for HEV in EDTA plasma

HEV RNA concentration (IU/mL)	Number reactive	Number of valid replicates	% Reactive	95% Lower confidence bound (one-sided)
40	187	187	100.0%	98.4%
20	179	188	95.2%	91.8%
10	165	189	87.3%	82.6%
6	113	187	60.4%	54.2%
2	52	189	27.5%	22.2%

## Reproducibility

The reproducibility of the **cobas**® HEV test on the **cobas**® 6800/8800 Systems was determined using the WHO International Standard for HEV (PEI code 6329/10). This study consisted of testing 3 panels of HEV at concentrations of approximately 0.5 x, 1 x and 2 x the LoD of the **cobas**® HEV test. Testing was performed for the following variability components:

- day-to-day variability over 3 days
- lot-to-lot variability using 3 different reagent lots of the **cobas**® HEV test
- instrument-to-instrument variability using 3 different **cobas**® 8800 Systems

Approximately 21 replicates were tested with each of the 3 panels for total of 63 replicates with each reagent lot. All valid reproducibility data were evaluated by calculating the percentage of reactive test results for each concentration level across all variable components.

The limits of two-sided 95% Confidence Intervals for each Reactive Rate were calculated for each of the three levels of HEV tested across 3 days, 3 reagent lots, and 3 **cobas**® 8800 Systems. The **cobas**® HEV test is reproducible over multiple days, reagent lots and multiple instruments. The results from reagent lot-to-lot variability are summarized in Table 17.

**Table 17** **cobas**® HEV test reagent lot-to-lot reproducibility summary

Analyte	Concentration	Reagent lot	% Reactive (reactive/valid replicates)	Lower limit of 95% confidence interval	Upper limit of 95% confidence interval
HEV	2 x LoD	1	100.0% (61/61)	94.1%	100.0%
HEV	2 x LoD	2	100.0% (63/63)	94.3%	100.0%
HEV	2 x LoD	3	100.0% (63/63)	94.3%	100.0%
HEV	1 x LoD	1	88.9% (56/63)	78.4%	95.4%
HEV	1 x LoD	2	96.8% (60/62)	88.8%	99.6%
HEV	1 x LoD	3	100.0% (63/63)	94.3%	100.0%
HEV	0.5 x LoD	1	82.5% (52/63)	70.9%	90.9%
HEV	0.5 x LoD	2	95.2% (60/63)	86.7%	99.0%
HEV	0.5 x LoD	3	84.1% (53/63)	72.7%	92.1%

## Genotype verification

The performance of the **cobas**® HEV test to detect 4 genotypes of HEV was determined by testing a total of 16 unique clinical samples and 7 HEV cultured isolates with known genotypes. All samples were quantified traceable to the HEV WHO standard. All 16 clinical samples were tested after dilution with normal, virus-negative (HEV) human EDTA-plasma to 5 x LoD of the **cobas**® HEV test, of which 10 samples were also tested neat. All 7 cultured isolates were tested after dilution with normal, virus-negative (HEV) human EDTA-plasma to 5 x LoD of the **cobas**® HEV test. All clinical samples and cultured isolates were detected at neat and/or at 5 x LoD (Table 18).

**Table 18** HEV clinical samples and cultured isolates

Genotype	Clinical samples	Clinical samples	Cultured isolates
	% Reactive (reactive/samples tested) neat	% Reactive (reactive/samples tested) diluted to 5 x LoD	% Reactive (reactive/samples tested) diluted to 5 x LoD
1	Not tested*	Not tested*	100.0% (3/3)
2	Not tested*	Not tested*	100.0% (1/1)
3	100.0% (10/10)	100.0% (10/10)	Not tested*
4	Not tested*	100.0% (6/6)	100.0% (3/3)

\*Insufficient volume to test at neat/diluted

## Analytical specificity

The analytical specificity of the **cobas®** HEV test was evaluated for cross-reactivity with 28 microorganisms at  $10^6$  particles, copies, or PFU/mL, which included 21 viral isolates, 6 bacterial strains and 1 yeast isolate (Table 19). The microorganisms were added to normal, virus-negative human EDTA-plasma and tested with and without HEV added to a concentration of approximately 3 x LoD of the **cobas®** HEV test. The tested microorganisms do not cross-react or interfere with the **cobas®** HEV test.

**Table 19** Microorganisms tested for analytical specificity

Viruses	Flavivirus	Bacteria	Yeast
Adenovirus 5	West Nile Virus	<i>Escherichia coli</i>	<i>Candida albicans</i>
Cytomegalovirus	Dengue Virus type 1	<i>Propionibacterium acnes</i>	-
Epstein-Barr Virus	Usutu Virus	<i>Staphylococcus aureus</i>	-
Herpes Simplex Virus type 1	-	<i>Staphylococcus epidermidis</i>	-
Herpes Simplex Virus type 2	-	<i>Streptococcus viridans</i>	-
Hepatitis A Virus	-	<i>Staphylococcus haemolyticus</i>	-
Hepatitis B Virus	-	-	-
Hepatitis C Virus	-	-	-
Hepatitis G Virus	-	-	-
Human Immunodeficiency Virus (HIV-1 Group M)	-	-	-
Human Immunodeficiency Virus (HIV-2)	-	-	-
Human T-cell lymphotropic Virus type I	-	-	-
Human T-cell lymphotropic Virus type II	-	-	-
Human Herpes Virus 6	-	-	-
Influenza Virus A	-	-	-
Parvovirus B19	-	-	-
Chikungunya Virus	-	-	-
Varicella Zoster Virus	-	-	-

Plasma samples from each of the disease states (Table 20) were tested without HEV and with HEV added to a concentration of approximately 3 x LoD of the **cobas®** HEV test. These disease states do not cross-react or interfere with the **cobas®** HEV test.

**Table 20** Disease states samples tested for analytical specificity

Disease state	Disease state	Disease state
Adenovirus type 5	Hepatitis B Virus	Human T-cell lymphotropic Virus type I
Cytomegalovirus	Hepatitis C Virus	Human T-cell lymphotropic Virus type II
Dengue Virus	Herpes Simplex Virus type 1	Parvovirus B19
Epstein-Barr Virus	Herpes Simplex Virus type 2	West Nile Virus
Hepatitis A Virus	Human Immunodeficiency Virus (HIV-1)	-

## Analytical specificity – interfering substances

### Endogenous interference substances

Plasma samples with abnormally high levels of triglycerides (up to 33.2 g/L), hemoglobin (up to 4.7 g/L), unconjugated bilirubin (up to 0.28 g/L), albumin (up to 60 g/L), and human DNA (up to 0.004 g/L) were tested with and without HEV added to a concentration of approximately 3 x LoD of the **cobas®** HEV test. Samples containing these endogenous substances did not interfere with the sensitivity or specificity of **cobas®** HEV test.

### Exogenous interference substances

Normal, virus-negative (HEV) human EDTA-plasma samples containing abnormally high concentrations of drugs (Table 21) were tested with and without HEV added to a concentration of 3 x LoD of the **cobas®** HEV test. These exogenous substances did not interfere with the sensitivity or specificity of the **cobas®** HEV test.

**Table 21** Clinical samples tested with drugs

Name of drug tested	Concentration
Acetaminophen	1324 µmol/L
Acetylsalicylic Acid	3620 µmol/L
Ascorbic Acid	342 µmol/L
Atorvastatin	600 µg Eq/L
Fluoxetine	11.2 µmol/L
Ibuprofen	2425 µmol/L
Loratadine	0.78 µmol/L
Nadolol	3.88 µmol/L
Naproxen	2170 µmol/L
Paroxetine	3.04 µmol/L
Phenylephrine HCL	491 µmol/L
Sertraline	1.96 µmol/L

## Correlation

### Performance evaluation of the cobas® HEV test compared to the Realstar® HEV RT-PCR Kit 1.0 test

The performance of the cobas® HEV test and the Realstar® HEV RT-PCR Kit 1.0 test (Altona Diagnostics) was compared using 100 individual HEV NAT-positive plasma samples. One hundred positive samples were tested neat and 67 positive samples were tested diluted 1:6. In addition, 100 HEV negative plasma samples were tested neat with both methods.

The seronegative samples demonstrated 100% specificity by generating 100 out of 100 non-reactive results with both methods.

For positive samples, both methods were in agreement based on the McNemar's test indicating that the performance of cobas® HEV test and Realstar® HEV RT-PCR Kit 1.0 test are equivalent (Table 22).

**Table 22** Correlation of positive samples (neat)

Methods	Methods	HEV results	HEV results
Realstar® HEV RT-PCR Kit 1.0 test	cobas® HEV	Neat	Diluted 1:6
Non-reactive	Non-reactive	0	3
Reactive	Non-Reactive	0	3
Non-reactive	Reactive	1	9
Reactive	Reactive	99	52
Total	Total	100	67
McNemar's Test, p-value (two-sided, $\alpha=0.05$ )	McNemar's Test, p-value (two-sided, $\alpha=0.05$ )	1.00	0.09

## Whole system failure

The whole system failure rate for the cobas® HEV test was determined by testing 100 replicates of EDTA plasma spiked with HEV. These samples were tested at a target concentration of approximately 3 x LoD and were run in pools of 1 (undiluted). The study was performed using the cobas® 8800 System with cobas® p 680 instrument (pipetting and pooling).

The results of this study determined that all replicates were reactive for HEV, resulting in a whole system failure rate of 0%. The two-sided 95% exact confidence interval was 0% for the lower bound and 3.62% for the upper bound [0%: 3.62%].

## Additional information

### Key test features





















































<b>Sample type</b>	Plasma
<b>Amount of sample required</b>	1000 µL*
<b>Amount of sample processed</b>	850 µL

\*Tubes used for testing may have different dead volumes and require more or less minimum volume. Contact your local Roche service representative for further information.

## Symbols

The following symbols are used in labeling for Roche PCR diagnostic products.

**Table 23** Symbols used in labeling for Roche PCR diagnostics products

 Age/DOB	 Device not for near-patient testing	 QS IU/PCR
 Ancillary Software	 Device not for self-testing	QS IU per PCR reaction, use the QS International Units (IU) per PCR reaction in calculation of the results.
 Assigned Range [copies/mL]	 Distributor (Note: The applicable country/region may be designated beneath the symbol)	 SN
 Assigned Range [IU/mL]	 Do not re-use	 Site
 EC REP	 Female	 Procedure Standard
 Barcode Data Sheet	 For IVD performance evaluation only	 STERILE EO
 LOT	 GTIN	 Store in dark
 Biological risks	 Importer	 Temperature limit
 REF	 IVD	 Test Definition File
 CE marking of conformity; this device is in conformity with the applicable requirements for CE marking of an in vitro diagnostic medical device	 LLR	 This way up
 Collect Date	 Male	 Procedure UltraSensitive
 Consult instructions for use	 Manufacturer	 UDI
 Contains sufficient for <n> tests	 CONTROL -	 ULR
 CONTENT	 Non-sterile	 Urine Fill Line
 CONTROL	 Patient Name	 Rx Only
 Date of manufacture	 Patient number	US Only: Federal law restricts this device to sale by or on the order of a physician.
 Device for near-patient testing	 Peel here	 Use-by date
 Device for self-testing	 CONTROL +	
	 QS copies / PCR	QS copies per PCR reaction, use the QS copies per PCR reaction in calculation of the results.



## Technical support

For technical support (assistance) please reach out to your local affiliate:  
[https://www.roche.com/about/business/roche\\_worldwide.htm](https://www.roche.com/about/business/roche_worldwide.htm)

## Manufacturer and importer

**Table 24** Manufacturer and importer



Roche Molecular Systems, Inc.  
1080 US Highway 202 South  
Branchburg, NJ 08876 USA  
[www.roche.com](http://www.roche.com)

Made in USA



Roche Diagnostics GmbH  
Sandhofer Strasse 116  
68305 Mannheim, Germany

## Trademarks and patents

See <https://diagnostics.roche.com/us/en/about-us/patents>

## Copyright

©2023 Roche Molecular Systems, Inc.



Roche Diagnostics GmbH  
Sandhofer Str. 116  
68305 Mannheim  
Germany



## References

1. Kamar N, Bendall R, Legrand-Abravanel F, et al. Hepatitis E. *Lancet*. 2012;379:2477-88.
2. Ahmad I, Holla RP, Jameel S. Molecular virology of hepatitis E virus. *Virus Res*. 2011;161:47-58.
3. Meng XJ, Purcell RH, Halbur PG, et al. A novel virus in swine is closely related to the human hepatitis E virus. *Proc Natl Acad Sci U S A*. 1997;94:9860-5.
4. Yugo DM, Meng XJ. Hepatitis E virus: foodborne, waterborne and zoonotic transmission. *Int J Environ Res Public Health*. 2013;10:4507-33.
5. Balayan MS, Andjaparidze AG, Savinskaya SS, et al. Evidence for a virus in non-A, non-B hepatitis transmitted via the fecal-oral route. *Intervirology*. 1983;20:23-31.
6. Dalton HR, Bendall R, Ijaz S, Banks M. Hepatitis E: an emerging infection in developed countries. *Lancet Infect Dis*. 2008;8:698-709.
7. Meng XJ. Hepatitis E as a zoonotic disease in the United States. Oral presentation at: National Institutes of Health Research Workshop: Hepatitis E in the United States; 7-8 March 2012; Bethesda, Maryland.
8. Mansuy JM, Peron JM, Abravanel F, et al. Hepatitis E in the south west of France in individuals who have never visited an endemic area. *J Med Virol*. 2004;74:419-24.
9. Dalton HR, Stableforth W, Thurairajah P, et al. Autochthonous hepatitis E in Southwest England: natural history, complications and seasonal variation, and hepatitis E virus IgG seroprevalence in blood donors, the elderly and patients with chronic liver disease. *Eur J Gastroenterol Hepatol*. 2008;20:784-90.
10. Wichmann O, Schimanski S, Koch J, et al. Phylogenetic and case-control study on hepatitis E virus infection in Germany. *J Infect Dis*. 2008;198:1732-41.
11. Dalton HR, Fellows HJ, Gane EJ, et al. Hepatitis E in New Zealand. *J Gastroenterol Hepatol*. 2007;22:1236-40.
12. Tsang TH, Denison EK, Williams HV, et al. Acute hepatitis E infection acquired in California. *Clin Infect Dis*. 2000;30:618-9.
13. Teo CG. Much meat, much malady: changing perceptions of the epidemiology of hepatitis E. *Clin Microbiol Infect*. 2010;16:24-32.
14. Colson P, Borentain P, Queyriaux B, et al. Pig liver sausage as a source of hepatitis E virus transmission to humans. *J Infect Dis*. 2010;202:825-34.
15. Matsuda H, Okada K, Takahashi K, Mishiro S. Severe hepatitis E virus infection after ingestion of uncooked liver from a wild boar. *J Infect Dis*. 2003;188:944.
16. Tei S, Kitajima N, Takahashi K, Mishiro S. Zoonotic transmission of hepatitis E virus from deer to human beings. *Lancet*. 2003;362:371-3.
17. Gotanda Y, Iwata A, Ohnuma H, et al. Ongoing subclinical infection of hepatitis E virus among blood donors with an elevated alanine aminotransferase level in Japan. *J Med Virol*. 2007;79:734-42.
18. Mitsui T, Tsukamoto Y, Suzuki S, et al. Serological and molecular studies on subclinical hepatitis E virus infection using periodic serum samples obtained from healthy individuals. *J Med Virol*. 2005;76:526-33.
19. Nelson KE, Kmush B, Labrique AB. The epidemiology of hepatitis E virus infections in developed countries and among immunocompromised patients. *Expert Rev Anti Infect Ther*. 2011;9:1133-48.

20. Navaneethan U, Al Mohajer M, Shata MT. Hepatitis E and pregnancy: understanding the pathogenesis. *Liver Int.* 2008;28:1190-9.
21. Jilani N, Das BC, Husain SA, et al. Hepatitis E virus infection and fulminant hepatic failure during pregnancy. *J Gastroenterol Hepatol.* 2007;22:676-82.
22. Peron JM, Bureau C, Poirson H, et al. Fulminant liver failure from acute autochthonous hepatitis E in France: description of seven patients with acute hepatitis E and encephalopathy. *J Viral Hepat.* 2007;14:298-303.
23. Ramachandran J, Eapen CE, Kang G, et al. Hepatitis E superinfection produces severe decompensation in patients with chronic liver disease. *J Gastroenterol Hepatol.* 2004;19:134-8.
24. Rein DB, Stevens GA, Theaker J, Wittenborn JS, Wiersma ST. The global burden of hepatitis E virus genotypes 1 and 2 in 2005. *Hepatology.* 2012;55:988-97.
25. Dalton HR, Bendall RP, Rashid M, et al. Host risk factors and autochthonous hepatitis E infection. *Eur J Gastroenterol Hepatol.* 2011;23:1200-5.
26. Dalton HR, Hazeldine S, Banks M, Ijaz S, Bendall R. Locally acquired hepatitis E in chronic liver disease. *Lancet.* 2007;369:1260.
27. Kamar N, Selves J, Mansuy JM, et al. Hepatitis E virus and chronic hepatitis in organ-transplant recipients. *N Engl J Med.* 2008;358:811-7.
28. Gerolami R, Moal V, Colson P. Chronic hepatitis E with cirrhosis in a kidney-transplant recipient. *N Engl J Med.* 2008;358:859-60.
29. Haagsma EB, van den Berg AP, Porte RJ, et al. Chronic hepatitis E virus infection in liver transplant recipients. *Liver Transpl.* 2008;14:547-53.
30. Kamar N, Garrouste C, Haagsma EB, et al. Factors associated with chronic hepatitis in patients with hepatitis E virus infection who have received solid organ transplants. *Gastroenterology.* 2011;140:1481-9.
31. Colson P, Kaba M, Moreau J, Brouqui P. Hepatitis E in an HIV-infected patient. *J Clin Virol.* 2009;45:269-71.
32. Peron JM, Mansuy JM, Recher C, et al. Prolonged hepatitis E in an immunocompromised patient. *J Gastroenterol Hepatol.* 2006;21:1223-4.
33. Tamura A, Shimizu YK, Tanaka T, et al. Persistent infection of hepatitis E virus transmitted by blood transfusion in a patient with T-cell lymphoma. *Hepatol Res.* 2007;37:113-20.
34. Matsubayashi K, Nagaoka Y, Sakata H, et al. Transfusion-transmitted hepatitis E caused by apparently indigenous hepatitis E virus strain in Hokkaido, Japan. *Transfusion.* 2004;44:934-40.
35. Mitsui T, Tsukamoto Y, Yamazaki C, et al. Prevalence of hepatitis E virus infection among hemodialysis patients in Japan: evidence for infection with a genotype 3 HEV by blood transfusion. *J Med Virol.* 2004;74:563-72.
36. Boxall E, Herborn A, Kochethu G, et al. Transfusion-transmitted hepatitis E in a 'nonhyperendemic' country. *Transfus Med.* 2006;16:79-83.
37. Colson P, Coze C, Gallian P, et al. Transfusion-associated hepatitis E, France. *Emerg Infect Dis.* 2007;13:648-9.
38. Khuroo MS, Kamili S, Yattoo GN. Hepatitis E virus infection may be transmitted through blood transfusions in an endemic area. *J Gastroenterol Hepatol.* 2004;19:778-84.
39. Zaaijer HL, Kok M, Lelie PN, et al. Hepatitis E in the Netherlands: imported and endemic. *Lancet.* 1993;341:826.

40. Christensen PB, Engle RE, Hjort C, et al. Time trend of the prevalence of hepatitis E antibodies among farmers and blood donors: a potential zoonosis in Denmark. *Clin Infect Dis*. 2008;47:1026-31.
41. Boutrouille A, Bakkali-Kassimi L, Cruciere C, Pavio N. Prevalence of anti-hepatitis E virus antibodies in French blood donors. *J Clin Microbiol*. 2007;45:2009-10.
42. Kaufmann A, Kenfak-Foguena A, Andre C, et al. Hepatitis E virus seroprevalence among blood donors in southwest Switzerland. *PLoS One*. 2011;6:e21150.
43. Vollmer T, Diekmann J, Johne R, et al. Novel approach for detection of hepatitis E virus infection in German blood donors. *J Clin Microbiol*. 2012;50:2708-13.
44. Takeda H, Matsubayashi K, Sakata H, et al. A nationwide survey for prevalence of hepatitis E virus antibody in qualified blood donors in Japan. *Vox Sang*. 2010;99:307-13.
45. Beale MA, Tettmar K, Szypulska R, Tedder RS, Ijaz S. Is there evidence of recent hepatitis E virus infection in English and North Welsh blood donors? *Vox Sang*. 2011;100:340-2.
46. Ijaz S, Szypulska R, Tettmar KI, Kitchen A, Tedder RS. Detection of hepatitis E virus RNA in plasma mini-pools from blood donors in England. *Vox Sang*. 2012;102:272.
47. Guo QS, Yan Q, Xiong JH, et al. Prevalence of hepatitis E virus in Chinese blood donors. *J Clin Microbiol*. 2010;48:317-8.
48. Baylis SA, Hanschmann KM, Blumel J, Nubling CM, HEV Collaborative Study Group. Standardization of hepatitis E virus (HEV) nucleic acid amplification technique-based assays: an initial study to evaluate a panel of HEV strains and investigate laboratory performance. *J Clin Microbiol*. 2011;49:1234-9.
49. Longo MC, Berninger MS, Hartley JL. Use of uracil DNA glycosylase to control carry-over contamination in polymerase chain reactions. *Gene*. 1990;93:125-8.
50. Savva R, McAuley-Hecht K, Brown T, Pearl L. The structural basis of specific base-excision repair by uracil-DNA glycosylase. *Nature*. 1995;373:487-93.
51. Mol CD, Arvai AS, Slupphaug G, et al. Crystal structure and mutational analysis of human uracil-DNA glycosylase: structural basis for specificity and catalysis. *Cell*. 1995;80:869-78.
52. Higuchi R, Dollinger G, Walsh PS, Griffith R. Simultaneous amplification and detection of specific DNA sequences. *Biotechnology (N Y)*. 1992;10:413-7.
53. Heid CA, Stevens J, Livak KJ, Williams PM. Real time quantitative PCR. *Genome Res*. 1996;6:986-94.
54. Chosewood LC, Wilson DE, eds. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. HHS Publication No. (CDC) 21-1112. US Department of Health and Human Services; 2009.
55. Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections*. 4th ed. M29-A4. Clinical and Laboratory Standards Institute: Wayne, PA; 2014.

## Document revision

Document Revision Information	
Doc Rev. 2.0 09/2023	<p>Addition of <b>cobas®</b> 5800 instrument throughout.</p> <p>Updated <b>Principles of the procedure</b> section.</p> <p>Added <b>Reagent handling requirements for cobas® 5800 System</b> section.</p> <p>Added <b>Additional materials required for the cobas® 5800</b> section.</p> <p>Updated <b>Instrumentation and software required</b> section.</p> <p>Updated <b>Precautions and handling requirements</b> section.</p> <p>Updated <b>Instructions for use</b> section.</p> <p>Updated <b>Results</b> section and added <b>System equivalency/system comparison</b> section.</p> <p>Updated <b>cobas®</b> branding.</p> <p>Please contact your local Roche Representative if you have any questions.</p>

The summary of safety and performance report can be found using the following link: <https://ec.europa.eu/tools/eudamed>