

**ONLINE DAT Cannabinoids II****Order information**

| REF         | CONTENT                                       |                    | Analyzer(s) on which <b>cobas c</b> pack(s) can be used |
|-------------|---|--------------------|---|
| 08058644190 | ONLINE DAT Cannabinoids II (850 tests)        | System-ID 2107 001 | <b>cobas c 503</b>                                      |
| 08771669190 | ONLINE DAT Cannabinoids II (150 tests)        | System-ID 2107 002 | <b>cobas c 503</b>                                      |
| 03304671190 | Preciset DAT Plus I calibrator CAL 5          | Code 20435         |   |
| 07978766190 | Serum DAT Control Low (ACQ Partner Channel*)  |                    |   |
| 07978740190 | Serum DAT Control High (ACQ Partner Channel*) |                    |   |
| 08063494190 | NaCl Diluent 9 % (123 mL)                     | System-ID 2906 001 |   |

\*Roche does not hold the product registration for Partner Channels. The legal manufacturer indicated on the kit is solely responsible for all of the design, legal, and regulatory aspects of the product.

**English****System information**

**THQ5S:** ACN 21077 (Serum/plasma): for qualitative assay, 50 ng/mL

**Intended use**

Cannabinoids II (THCII) is an in vitro diagnostic test for the qualitative detection of cannabinoids in human serum and plasma on Roche/Hitachi **cobas c** systems at a cutoff concentration of 50 ng/mL.

**Cannabinoids II provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC-MS) or liquid chromatography coupled with tandem mass spectrometry (LC-MS/MS) is the preferred confirmatory method.<sup>1</sup> Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.**

**Summary**

The principal psychoactive component of the hemp plant, *Cannabis sativa*, is generally accepted to be  $\Delta^9$  tetrahydrocannabinol ( $\Delta^9$  THC), although other cannabinoids may contribute to the psychological and physiological actions of marijuana. The acute effects of marijuana use, concomitant with the desired "high", are memory impairment, time confusion, interference with learning, impaired motor skills and depersonalization.<sup>2,3,4</sup> These effects are also manifested in chronic users in addition to cardiovascular, pulmonary, and reproductive effects. Marijuana is usually smoked, but may be ingested, either incorporated into food or as a liquid extract (tea). It is rapidly absorbed from the lungs into the blood with rapid onset of effects; the onset is slower but prolonged when ingested. The natural cannabinoids and their metabolic products are fat soluble and are stored in the body's fatty tissues, including brain tissue, for prolonged periods after use.<sup>5</sup>

Cannabinoid metabolites are found in blood, bile, feces, and urine and may be detected in urine within hours of exposure. Because of their fat solubility, they also remain in the body's fatty tissues with slow release and subsequent urinary excretion for days, weeks, and even months after the last exposure, depending on the intensity and frequency of use.<sup>1</sup> The prominent  $\Delta^9$  THC metabolite, 11-nor- $\Delta^9$  THC-9-carboxylic acid ( $\Delta^9$  COOH-THC), is the primary urinary marker for detecting marijuana use.

**Test principle**

The assay is based on the kinetic interaction of microparticles in a solution (KIMS)<sup>6,7</sup> as measured by changes in light transmission. In the absence of sample drug, soluble drug conjugates bind to antibody-bound microparticles, causing the formation of particle aggregates. As the aggregation reaction proceeds in the absence of sample drug, the absorbance increases.

When a serum sample contains the drug in question, this drug competes with the drug derivative conjugate for microparticle-bound antibody. Antibody bound to sample drug is no longer available to promote particle aggregation, and subsequent particle lattice formation is inhibited. The presence of sample drug diminishes the increasing absorbance in proportion to the concentration of drug in the sample. Sample drug content is determined relative to the value obtained for a known cutoff concentration of drug.<sup>8</sup>

**Reagents - working solutions**

**R1** Conjugated cannabinoid derivative; buffer; bovine serum albumin; 0.09 % sodium azide

**R3** Microparticles attached to cannabinoid antibody (mouse monoclonal); buffer; bovine serum albumin; 0.09 % sodium azide

R1 is in position B and R3 is in position C.

**Precautions and warnings**

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines. Safety data sheet available for professional user on request.

**Reagent handling**

Ready for use

Carefully invert reagent container several times prior to use to ensure that the reagent components are mixed.

**Storage and stability**

Shelf life at 2-8 °C: See expiration date on **cobas c** pack label

On-board in use and refrigerated on the analyzer: 8 weeks

**Do not freeze.****Specimen collection and preparation**

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum: Serum tubes with and without separating gel.

Plasma: K<sub>2</sub>- or K<sub>3</sub>-EDTA, lithium heparin.

Stability: 5 days capped at 15-25 °C  
14 days capped at 2-8 °C  
6 months capped at -20 °C

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

See the limitations and interferences section for details about possible sample interferences.

Specimens can be repeatedly frozen and thawed up to 3 times.

Invert thawed specimens several times prior to testing.

**CAUTION:** Specimen dilutions should only be used to interpret results of Calc.? and Samp.? alarms, or when estimating concentration in preparation for GC-MS. Dilution results are not intended for patient values. Dilution procedures, when used, should be validated.

**Materials provided**

See "Reagents – working solutions" section for reagents.

**Materials required (but not provided)**

See "Order information" section

General laboratory equipment

**Assay**

THC and its derivatives may adsorb onto plastics.<sup>9</sup> To minimize the potential for lowering the drug concentration of any sample containing THC, the following is recommended:

1. Dispense > 0.5 mL of each sample (calibrators, controls and patient specimens) into separate analyzer sample cups by pouring over from the primary container or by dispensing with a glass pipette.
2. Avoid the use of plastic pipettes and/or tips due to the potential for adsorbance and possible decrease of THC concentration.
3. Assay the samples within 2 hours of dispensing into the sample cup.
4. Do not return any unused material back into the original sample container.

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

The performance of applications not validated by Roche is not warranted and must be defined by the user.

**Application for serum and plasma****Test definition**

|                       |               |                            |                       |
|-----------------------|---------------|----------------------------|-----------------------|
|                       | Qualitative   |                            |                       |
| Reporting time        | 10 min        |                            |                       |
| Wavelength (sub/main) | - /570 nm     |                            |                       |
| Reagent pipetting     |               | Diluent (H <sub>2</sub> O) |                       |
| R1                    | 64 µL         | -                          |                       |
| R3                    | 28 µL         | -                          |                       |
| <i>Sample volumes</i> | <i>Sample</i> | <i>Sample dilution</i>     |                       |
|                       |               | <i>Sample</i>              | <i>Diluent (NaCl)</i> |
| Normal                | 1.8 µL        | -                          | -                     |
| Decreased             | 1.8 µL        | -                          | -                     |
| Increased             | 1.8 µL        | -                          | -                     |

For further information about the assay test definitions refer to the application parameters setting screen of the corresponding analyzer and assay.

**Calibration**

|                       |  |
|-----------------------|--|
| Calibrators           | <i>Qualitative application</i><br><i>50 ng/mL cutoff assay</i><br>S1: Preciset DAT Plus I calibrator - CAL 5,<br>200 ng/mL with automatic pre-dilution<br>The drug concentration of the calibrator has been verified by GC-MS. |
| Calibration K factor  | For the qualitative application a K factor of -1000 is predefined in the application settings.   |
| Calibration mode      | <i>Qualitative application</i><br>Linear   |
| Calibration frequency | Full calibration<br>- after reagent lot change<br>- as required following quality control procedures   |

For the cutoff calibrator a value of "0" is encoded in the e-barcode in order to ensure flagging of positive samples with >Test and negative absorbance values for negative samples.

Calibration interval may be extended based on acceptable verification of calibration by the laboratory.

Traceability: This method has been standardized against a primary reference method (GC-MS).

**Quality control**

For quality control, use control materials as listed in the "Order information" section. In addition, other suitable control material can be used.

Drug concentrations of the high and low controls have been verified by GC-MS.

The control intervals and limits should be adapted to each laboratory's individual requirements. It is recommended to perform quality control always after lot calibration and subsequently at least every 8 weeks.

Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

**Results**

The cutoff calibrator is used as a reference in distinguishing between preliminary positive and negative samples. Samples producing a positive or "0" absorbance value are considered preliminary positive. Preliminary positive samples are flagged with >Test. Samples producing a negative absorbance value are considered negative. Negative samples are preceded by a minus sign.

As with any sensitive test for drugs of abuse on automated clinical chemistry analyzers, the possibility exists for analyte carry-over from a sample with an extremely high concentration to a normal (negative) sample which immediately follows it.

Preliminary positive results should be confirmed by another method.

**Limitations - interference**

Criterion: No cross-over at initial values of samples of 25 ng/mL and 75 ng/mL (control levels).

See the "Specific performance data" section of this document for information on substances tested with this assay. There is the possibility that other substances and/or factors may interfere with the test and cause erroneous results (e.g., technical or procedural errors).

A preliminary positive result with this assay indicates the presence of cannabinoids and/or cannabinoid metabolites in serum. It does not measure the level of intoxication.

Icterus:<sup>10</sup> No significant interference up to an I index of 60 for conjugated and unconjugated bilirubin (approximate conjugated and unconjugated bilirubin concentration: 1026 µmol/L or 60 mg/dL).

Hemolysis:<sup>10</sup> No significant interference up to an H index of 1000 (approximate hemoglobin concentration: 622 µmol/L or 1000 mg/dL).

Lipemia (Intralipid):<sup>10</sup> No significant interference up to an L index of 900. There is poor correlation between the L index (corresponds to turbidity) and triglycerides concentration.

Rheumatoid factors: No significant interference from rheumatoid factors up to a concentration of 1200 IU/mL.

Immunoglobulins: No significant interference from immunoglobulins up to a concentration of 16 g/L (simulated by human immunoglobulin A), up to a concentration of 70 g/L (simulated by human immunoglobulin G) and up to a concentration of 10 g/L (simulated by human immunoglobulin M).

Albumin: No significant interference from human serum albumin up to a concentration of 70 g/L.

As with any assay employing mouse antibodies, the possibility exists for interference by human anti-mouse antibodies (HAMA) in the sample, which could cause falsely lowered results.

In very rare cases, gammopathy, in particular type IgM (Waldenström's macroglobulinemia), may cause unreliable results.<sup>11</sup>

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

**ACTION REQUIRED**

**Special Wash Programming:** The use of special wash steps is mandatory when certain test combinations are run together on **cobas c** systems. All special wash programming necessary for avoiding carry-over is available via the **cobas** link. The latest version of the carry-over evasion list can be found with the NaOHD/SMS/SCCS Method Sheet for information. For further instructions refer to the operator's manual.

**Expected values***Qualitative assay*

Results of this assay distinguish preliminary positive (≥ 50 ng/mL) from negative samples only. The amount of drug detected in a preliminary positive sample cannot be estimated.

### Specific performance data

Representative performance data on the analyzers are given below. These data represent the performance of the analytical procedure itself.

Results obtained in individual laboratories may differ due to heterogenous sample materials, aging of analyzer components and mixture of reagents running on the analyzer.

### Precision

Precision was determined using human samples and controls in accordance with the CLSI (Clinical and Laboratory Standards Institute) EP05-A3 requirements with repeatability (n = 84) and intermediate precision (2 aliquots per run, 2 runs per day, 21 days). The following results were obtained:

### Qualitative precision - 50 ng/mL

| Cutoff (50)  | Number tested | Correct results | Confidence level        |
|--------------|---------------|-----------------|-------------------------|
| Serum -75 %  | 84            | 84              | > 95 % negative reading |
| ACQ-L        | 84            | 84              | > 95 % negative reading |
| Cutoff serum | 84            | n.a.*           | n.a.*                   |
| ACQ-H        | 84            | 84              | > 95 % positive reading |
| Serum +75 %  | 84            | 84              | > 95 % positive reading |

\*n.a. = not applicable

### Accuracy

110 serum samples, screened negative for cannabinoids on a **cobas c 501** analyzer were evaluated with the Cannabinoids II assay on a **cobas c 503** analyzer. 100 % of these normal serums were negative for all cutoffs with the Cannabinoids II assay on a **cobas c 503** analyzer. 54 serum samples screened positive for cannabinoids relative to the 50 ng/mL cutoff on a **cobas c 501** analyzer were evaluated with the Cannabinoids II assay on a **cobas c 503** analyzer. At the 50 ng/mL cutoff, 100 % of the samples were positive on both the **cobas c 501** analyzer and the **cobas c 503** analyzer.

| Cannabinoids II correlation (cutoff = 50 ng/mL) |   |                      |     |
|---|---|----------------------|-----|
|   |   | cobas c 501 analyzer |     |
|   |   | +                    | -   |
| cobas c 503 analyzer                            | + | 54                   | 0   |
|   | - | 0                    | 110 |

### Analytical specificity

The specificity of this assay for various cannabinoids and cannabinoid metabolites was determined by generating inhibition curves for each of the compounds listed and determining the approximate quantity of each compound that is equivalent in assay reactivity to the 50 ng/mL  $\Delta^9$  COOH-THC assay cutoff. The following results were obtained on a **cobas c 501** analyzer.

| Compound                                 | ng/mL                  | Approximate % cross-reactivity |
|--|------------------------|--------------------------------|
|  | Equivalent to 50 ng/mL |                                |
| <b>11-nor-9-carboxy-THC</b>              |                        |                                |
| 8- $\beta$ -11-Dihydroxy- $\Delta^9$ THC | 94.9                   | 52.7                           |
| 8- $\alpha$ -Hydroxy- $\Delta^9$ THC     | 91.8                   | 54.5                           |
| 11-Hydroxy- $\Delta^9$ THC               | 327                    | 15.3                           |

### Drug interference

Interfering substances were added to serum containing 11-nor-9-carboxy-THC at -50 % and +50 % of the cutoff level at the concentration listed below. Samples were tested and the following results were obtained on a **cobas c 501** analyzer.

| Compound       | Comp. conc. mg/L | Neg. level | Pos. level |
|----------------|------------------|------------|------------|
| Acetaminophen  | 200              | neg        | pos        |
| Acetylcysteine | 1660             | neg        | pos        |

|                                   |          |     |     |
|-----------------------------------|----------|-----|-----|
| Acetylsalicylic acid              | 1000     | neg | pos |
| Amitriptyline                     | 1.00     | neg | pos |
| Ampicillin-Na                     | 1000     | neg | pos |
| Ascorbic acid                     | 300      | neg | pos |
| Caffeine                          | 59.8     | neg | pos |
| Cefoxitin                         | 2500     | neg | pos |
| Cyclosporine                      | 5.00     | neg | pos |
| <i>d</i> -Amphetamine             | 1.36     | neg | pos |
| Doxycycline                       | 50.0     | neg | pos |
| Erythromycin                      | 59.9     | neg | pos |
| Fenoprofen                        | 195      | neg | pos |
| Furosemide                        | 59.9     | neg | pos |
| Gentisic acid                     | 18.0     | neg | pos |
| Heparin                           | 5000 U/L | neg | pos |
| Hydrochlorothiazide               | 6.02     | neg | pos |
| <i>l</i> -Amphetamine             | 1.00     | neg | pos |
| Ibuprofen                         | 50.0     | neg | pos |
| Imipramine                        | 0.70     | neg | pos |
| Ketamine                          | 10.0     | neg | pos |
| Levodopa                          | 20.0     | neg | pos |
| Lidocaine                         | 12.0     | neg | pos |
| Methyldopa + 1.5 H <sub>2</sub> O | 20.0     | neg | pos |
| Metronidazole                     | 200      | neg | pos |
| Naproxen                          | 499      | neg | pos |
| Phenylbutazone                    | 400      | neg | pos |
| Procaine                          | 20.0     | neg | pos |
| Promethazine                      | 1.20     | neg | pos |
| Quinidine                         | 12.0     | neg | pos |
| Quinine                           | 48.0     | neg | pos |
| Rifampicin                        | 60.0     | neg | pos |
| Tetracycline                      | 15.1     | neg | pos |
| Theophylline                      | 100      | neg | pos |
| Trifluoperazine                   | 1.00     | neg | pos |

### References

- Karch SB, ed. Drug Abuse Handbook. Boca Raton, FL: CRC Press LLC 1998.
- Tinklenberg JR, Darley CF. Psychological and cognitive effects of cannabis. In: Connell H, Dorn N, eds. Cannabis and Man: Proceedings of Third International Cannabinoids Conference, London: Churchill Livingstone 1975.
- Klonoff H. Marijuana and driving in real-life situations. Science 1974;186:317-324.
- Melges FT, Tinklenberg JR, Hollister LE, et al. Temporal disintegration and depersonalization during marijuana intoxication. Arch Gen Psychiatry 1970;23:204-210.
- Lemberger L, Tamarkin NR, Axelrod J, et al. Delta-9-tetrahydrocannabinol: metabolism and disposition in long-term marijuana smokers. Science 1971 Jul 2;173(3991):72-74.
- Armbruster DA, Schwarzhoff RH, Pierce BL, et al. Method comparison of EMIT II and ONLINE with RIA for drug screening. J Forensic Sci 1993;38:1326-1341.
- Armbruster DA, Schwarzhoff RH, Hubster EC, et al. Enzyme immunoassay, kinetic microparticle immunoassay, radioimmunoassay, and fluorescence polarization immunoassay compared for drugs-of-abuse screening. Clin Chem 1993;39:2137-2146.




**ONLINE DAT Cannabinoids II**

- 8 Antonian E, McNally AJ, Ng C, et al. An Abuscreen immunoassay for THC metabolites in urine on the Olympus AU5000 Series Clinical Analyzers. In: American Academy of Forensic Sciences. Program: The Forensic Sciences and Government. Abstract 1991;177.
- 9 Decker WJ. Laboratory support of drug abuse control programs: An overview. *Clinical Toxicology* 1977;10(1):28.
- 10 Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. *Clin Chem* 1986;32:470-475.
- 11 Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. *Clin Chem Lab Med* 2007;45(9):1240-1243.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

**Symbols**

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard (for USA: see [dialog.roche.com](http://dialog.roche.com) for definition of symbols used):

|   |                                       |
|---|---------------------------------------|
|  | Contents of kit                       |
|  | Volume after reconstitution or mixing |
|  | Global Trade Item Number              |

ABUSCREEN, COBAS, COBAS C, ONLINE DAT and PRECISET are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2020, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim  
[www.roche.com](http://www.roche.com)

