

Elecsys Anti-TPO

REF			SYSTEM
06368590190	06368590500	100	cobas e 411 cobas e 601 cobas e 602

English

System information

For **cobas e 411** analyzer: test number 720
 For **cobas e 601** and **cobas e 602** analyzers: Application Code Number 137

Intended use

Immunoassay for the in vitro quantitative determination of antibodies to thyroid peroxidase in human serum and plasma. The anti-TPO determination is used as an aid in the diagnosis of autoimmune thyroid diseases.

The electrochemiluminescence immunoassay "ECLIA" is intended for use on **cobas e** immunoassay analyzers.

Summary

Thyroid-specific peroxidase (TPO) is synthesized in the endoplasmic reticulum, where it is folded to its native state and undergoes core glycosylation, before being transported to the apical plasma membrane of thyrocytes.^{1,2}

In synergy with thyroglobulin (Tg) this enzyme has an essential function in the iodination of L-tyrosine and the chemical coupling of the resulting mono- and di-iodotyrosine to form the thyroid hormones T4, T3, and rT3.³

TPO is a potential autoantigen. Elevated serum titers of antibodies to TPO are found in several forms of thyroiditis caused by autoimmunity.^{4,5} TPO was identified as the causative antigen in 1985 when studies demonstrated that human antisera reacting to "microsomal antigen" precipitated TPO prepared from Graves' disease thyroid tissue.^{6,7} Clinically the two terms anti-TPO and microsomal antibody can be used synonymously; there are differences, however, with regard to the test methods.

High anti-TPO titers are found in up to 90 % of patients with chronic Hashimoto's thyroiditis. In Graves' disease, 70 % of the patients have an elevated titer.^{4,8,9} Although the sensitivity of the procedure can be increased by simultaneously determining other thyroid antibodies (anti-Tg, TSH-receptor-antibody - TRAb), a negative finding does not rule out the possibility of an autoimmune disease. The magnitude of the antibody titer does not correlate with the clinical activity of the disease.^{8,9,10} Initially elevated titers can become negative after lengthy periods of illness or during remission. If antibodies reappear following remission, then a relapse is probable.¹¹

Whereas the usual microsomal antibody tests employ unpurified microsomes as an antigen preparation, the anti-TPO tests use a purified peroxidase. The two procedures are of comparable performance in terms of clinical sensitivity, but better lot-to-lot consistency and higher clinical specificity can be expected from anti-TPO tests due to the higher quality of the antigen used.

Recombinant antigen and polyclonal anti-TPO antibodies are used in the Elecsys Anti-TPO assay.

Test principle

Competition principle. Total duration of assay: 18 minutes.

- 1st incubation: 20 µL of sample are incubated with anti-TPO-antibodies labeled with a ruthenium complex^{a)}.
- 2nd incubation: After addition of biotinylated TPO and streptavidin-coated microparticles, the anti-TPO antibodies in the sample compete with the ruthenium-labeled anti-TPO antibodies for the biotinylated TPO antigen. The entire complex becomes bound to the solid phase via interaction of biotin and streptavidin.
- The reaction mixture is aspirated into the measuring cell where the microparticles are magnetically captured onto the surface of the electrode. Unbound substances are then removed with ProCell/ProCell M. Application of a voltage to the electrode then induces chemiluminescent emission which is measured by a photomultiplier.
- Results are determined via a calibration curve which is instrument-specifically generated by 2-point calibration and a master curve provided via the reagent barcode or e-barcode.

a) Tris(2,2'-bipyridyl)ruthenium(II)-complex (Ru(bpy)₃²⁺)

Reagents - working solutions

The reagent rackpack is labeled as A-TPO.

- M Streptavidin-coated microparticles (transparent cap), 1 bottle, 6.5 mL:
 Streptavidin-coated microparticles 0.72 mg/mL; preservative.
- R1 Anti-TPO-Ab~Ru(bpy)₃²⁺ (gray cap), 1 bottle, 9 mL:
 Polyclonal anti-TPO antibody (sheep) labeled with ruthenium complex 1.0 mg/L; TRIS buffer 100 mmol/L, pH 7.2; preservative.
- R2 TPO-biotin (black cap), 1 bottle, 9 mL:
 Biotinylated TPO (recombinant) 0.15 mg/L; TRIS buffer 30 mmol/L, pH 7.0; preservative.

Precautions and warnings

For in vitro diagnostic use for laboratory professionals. Exercise the normal precautions required for handling all laboratory reagents.

Infectious or microbial waste:

Warning: handle waste as potentially biohazardous material. Dispose of waste according to accepted laboratory instructions and procedures.

Environmental hazards:

Apply all relevant local disposal regulations to determine the safe disposal. Safety data sheet available for professional user on request.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:



Warning

H317 May cause an allergic skin reaction.

Prevention:

P261 Avoid breathing mist or vapours.

P272 Contaminated work clothing should not be allowed out of the workplace.

P280 Wear protective gloves.

Response:

P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.

P362 + P364 Take off contaminated clothing and wash it before reuse.

Disposal:

P501 Dispose of contents/container to an approved waste disposal plant.

Hazardous components:

- 2-methyl-2H-isothiazol-3-one hydrochloride
- Product safety labeling follows EU GHS guidance.

Contact phone: all countries: +49-621-7590

Avoid foam formation in all reagents and sample types (specimens, calibrators and controls).

Elecsys Anti-TPO

Reagent handling

The reagents in the kit have been assembled into a ready-for-use unit that cannot be separated.

All information required for correct operation is read in from the respective reagent barcodes.

Storage and stability

Store at 2-8 °C.

Do not freeze.

Store the Elecsys reagent kit **upright** in order to ensure complete availability of the microparticles during automatic mixing prior to use.

Stability:	
unopened at 2-8 °C	up to the stated expiration date
after opening at 2-8 °C	6 weeks
on the analyzers	2 weeks

Specimen collection and preparation

Only the specimens listed below were tested and found acceptable.

Serum collected using standard sampling tubes or tubes containing separating gel.

Li-heparin plasma.

Criterion: Slope 0.9-1.1 + intercept within $\pm 2x$ analytical sensitivity (LDL) + coefficient of correlation ≥ 0.95 .

Stable for 3 days at 2-8 °C, 1 month at -20 °C. Freeze only once.¹²

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

Do not use heat-inactivated samples.

Do not use samples and controls stabilized with azide.

Ensure the samples, calibrators and controls are at 20-25 °C prior to measurement.

Due to possible evaporation effects, samples, calibrators and controls on the analyzers should be analyzed/measured within 2 hours.

Materials provided

See "Reagents – working solutions" section for reagents.

Materials required (but not provided)

- [REF] 06472931190, Anti-TPO CalSet, for 4 x 1.5 mL
- [REF] 05042666191, PreciControl ThyroAB, for 4 x 2.0 mL
- [REF] 11732277122, Diluent Universal, 2 x 16 mL sample diluent
- [REF] 03183971122, Diluent Universal, 2 x 36 mL sample diluent
- General laboratory equipment

▪ **cobas e** analyzer

Additional materials for the **cobas e 411** analyzer:

- [REF] 11662988122, ProCell, 6 x 380 mL system buffer
- [REF] 11662970122, CleanCell, 6 x 380 mL measuring cell cleaning solution
- [REF] 11930346122, Elecsys SysWash, 1 x 500 mL washwater additive
- [REF] 11933159001, Adapter for SysClean
- [REF] 11706802001, AssayCup, 60 x 60 reaction cups
- [REF] 11706799001, AssayTip, 30 x 120 pipette tips
- [REF] 11800507001, Clean-Liner

Additional materials for **cobas e 601** and **cobas e 602** analyzers:

- [REF] 04880340190, ProCell M, 2 x 2 L system buffer
- [REF] 04880293190, CleanCell M, 2 x 2 L measuring cell cleaning solution

- [REF] 03023141001, PC/CC-Cups, 12 cups to prewarm ProCell M and CleanCell M before use
- [REF] 03005712190, ProbeWash M, 12 x 70 mL cleaning solution for run finalization and rinsing during reagent change
- [REF] 12102137001, AssayTip/AssayCup, 48 magazines x 84 reaction cups or pipette tips, waste bags
- [REF] 03023150001, WasteLiner, waste bags
- [REF] 03027651001, SysClean Adapter M

Additional materials for all analyzers:

- [REF] 11298500316, ISE Cleaning Solution/Elecsys SysClean, 5 x 100 mL system cleaning solution

Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

Resuspension of the microparticles takes place automatically prior to use. Read in the test-specific parameters via the reagent barcode. If in exceptional cases the barcode cannot be read, enter the 15-digit sequence of numbers.

Bring the cooled reagents to approximately 20 °C and place on the reagent disk (20 °C) of the analyzer. Avoid foam formation. The system automatically regulates the temperature of the reagents and the opening/closing of the bottles.

Calibration

Traceability: This method has been standardized against the NIBSC (National Institute for Biological Standards and Control) 66/387 Standard.

Every Elecsys reagent set has a barcoded label containing specific information for calibration of the particular reagent lot. The predefined master curve is adapted to the analyzer using the relevant CalSet.

Calibration frequency: Perform calibration on all analyzers as follows:

- with every reagent kit
- Renewed calibration on all analyzers:
- daily: when using the same reagent kit on the analyzers
 - as required: e.g. quality control findings outside the defined limits

Quality control

Use PreciControl ThyroAB or other suitable controls for routine quality control procedures.

Controls for the various concentration ranges should be run individually at least once every 24 hours when the test is in use, once per reagent kit, and following each calibration.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

If necessary, repeat the measurement of the samples concerned.

Follow the applicable government regulations and local guidelines for quality control.

Calculation

The analyzer automatically calculates the analyte concentration of each sample (either in IU/mL or kIU/L).

Limitations - interference

The assay is unaffected by icterus (bilirubin $\leq 1129 \mu\text{mol/L}$ or $\leq 66 \text{ mg/dL}$), hemolysis (Hb $\leq 0.15 \text{ mmol/L}$ or $\leq 0.24 \text{ g/dL}$), lipemia (triglycerides $\leq 23.9 \text{ mmol/L}$ or $\leq 2100 \text{ mg/dL}$) and biotin ($\leq 40.9 \text{ nmol/L}$ or $\leq 10 \text{ ng/mL}$).

Criterion: Recovery within $\pm 10 \%$ of initial value.

Samples should not be taken from patients receiving therapy with high biotin doses (i.e. $> 5 \text{ mg/day}$) until at least 8 hours following the last biotin administration.

No interference was observed from rheumatoid factors up to a concentration of 450 IU/mL.

Pharmaceutical substances

In vitro tests were performed on 16 commonly used pharmaceuticals. No interference with the assay was found.

Elecsys Anti-TPO

Analytical specificity

A 0.3 % cross-reactivity with human autoantibodies to thyroglobulin (4000 IU/mL) were found, tested with anti-TPO concentrations of approximately 50 IU/mL and 250 IU/mL.

References

- Fayadat L, Niccoli-Sire P, Lanet J, et al. Human thyroperoxidase is largely retained and rapidly degraded in the endoplasmic reticulum. Its N-glycans are required for folding and intracellular trafficking. *Endocrinology* 1998;139(10):4277-4285.
- Kuliawat R, Ramos-Castañeda J, Liu Y, et al. Intracellular trafficking of thyroid peroxidase to the cell surface. *J Biol Chem* 2005;280(30):27713-27718.
- Suzuki K, Kawashima A, Yoshihara A, et al. Role of thyroglobulin on negative feedback autoregulation of thyroid follicular function and growth. *J Endocrinol* 2011;209:169-174.
- Effraimidis G, Wiersinga WM. Autoimmune thyroid disease: old and new players. *Eur J Endocrinol* 2014;170(6):241-252.
- McIntosh RS, Asghar MS, Weetman AP. The antibody response in human autoimmune thyroid disease. *Clin Sci* 1997;(92)6:529-541.
- Czarnocka B, Ruf J, Ferrand M, et al. Purification of the human thyroid peroxidase and its identification as the microsomal antigen involved in autoimmune thyroid diseases. *FEBS Letters* 1985;190:147-152.
- Portmann L, Hamada N, Heinrich G, et al. Antithyroid peroxidase antibody in patients with autoimmune thyroid disease: possible identity with anti-microsomal antibody. *J Clin Endocrinol Metab* 1985;61:1001-1003.
- Volpé R. Rational Use of Thyroid Function Tests. *Crit Rev Clin Lab Sci* 1997;34(5):405-438.
- Feldt-Rasmussen U. Analytical and clinical performance goals for testing autoantibodies to thyroperoxidase, thyroglobulin, and thyrotropin receptor. *Clin Chem* 1996;42(1):160-163.
- Utiger RD. The pathogenesis of autoimmune thyroid disease. *N Eng J Med* 1991;325:278-279.
- Schott M, Eckstein A, Willenberg HS, et al. Improved prediction of relapse of Graves' thyrotoxicosis by combined determination of TSH receptor and thyroperoxidase antibodies. *Horm Metab Res* 2007;39(1):56-61.
- Greiling H, Gressner AM. *Lehrbuch der Klinischen Chemie und Pathobiochemie*. 3rd edition, Stuttgart; New York: Schattauer 1995:1012.
- Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. *J Clin Chem Clin Biochem* 1988 Nov;26(11):783-790.

For further information, please refer to the appropriate user guide or operator's manual for the analyzer concerned, the respective application sheets and the Method Sheets of all necessary components (if available in your country).

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard:

	Contents of kit
	Analyzers/Instruments on which reagents can be used
	Reagent
	Calibrator
	Volume for reconstitution
	Global Trade Item Number

Rx only

For USA: Caution: Federal law restricts this device to sale by or on the order of a physician.

COBAS, NAVIFY, ELECSYS and PRECICONTROL are trademarks of Roche.
All other product names and trademarks are the property of their respective owners.
Additions, deletions or changes are indicated by a change bar in the margin.
© 2024, Roche Diagnostics

0123

Roche Diagnostics GmbH
Sandhofer Strasse 116
68305 Mannheim, Germany
www.roche.com
+800 5505 6606

