

# **cobas<sup>®</sup> HEV**

---

## **Nucleic acid test for use on the cobas<sup>®</sup> 6800/8800 Systems**

*For in vitro diagnostic use*

**cobas<sup>®</sup> HEV – 96**

P/N: 07001045190

**cobas<sup>®</sup> HEV Control Kit**

P/N: 07001100190

**cobas<sup>®</sup> NHP Negative Control Kit**

P/N: 07002220190

# Table of contents

<b>Intended use .....</b>	<b>4</b>
<b>Summary and explanation of the test.....</b>	<b>4</b>
<b>Reagents and materials.....</b>	<b>6</b>
cobas® HEV reagents and controls.....	6
cobas omni reagents for sample preparation.....	9
Reagent storage and handling requirements.....	10
Additional materials required .....	11
Instrumentation and software required .....	11
<b>Precautions and handling requirements .....</b>	<b>12</b>
Warnings and precautions.....	12
Reagent handling .....	12
Good laboratory practice .....	13
<b>Sample collection, transport, storage, and pooling .....</b>	<b>13</b>
Living donor samples .....	13
<b>Instructions for use.....</b>	<b>15</b>
Automated sample pipetting and pooling (optional) .....	15
Procedural notes .....	15
Running the cobas® HEV test.....	15
<b>Results.....</b>	<b>16</b>
Quality control and validity of results.....	16
Interpretation of results .....	17
Repeat testing of individual sample(s) .....	17
Procedural limitations .....	17
<b>Non-clinical performance evaluation .....</b>	<b>18</b>
Key performance characteristics .....	18
Living donor samples.....	18
Limit of Detection (LoD) .....	18
Reproducibility .....	19
Genotype verification .....	20
Analytical specificity .....	21

---

Analytical specificity – interfering substances .....	22
Correlation.....	23
Whole System Failure.....	23
<b>Additional information .....</b>	<b>24</b>
Key test features .....	24
Manufacturer and distributors .....	26
Trademarks and patents.....	26
Copyright .....	26
References .....	27
Document revision .....	30

## Intended use

The cobas® HEV test is a qualitative *in vitro* nucleic acid amplification test for the direct detection of hepatitis E virus (HEV) RNA (genotypes 1-4) in human plasma.

This test is intended for use to screen donor samples for HEV RNA in plasma samples from individual human donors, including donors of whole blood, blood components (red cells, platelets, and plasma), and other living donors. Plasma from all donors may be screened as individual samples. For donations of whole blood and blood components, plasma samples may be tested individually or plasma may be tested in pools comprised of aliquots of individual samples.

This test is not intended for use on samples of cord blood.

This test is not intended for use as an aid in diagnosis for HEV.

## Summary and explanation of the test

### Background: Screening of blood for transfusion-transmitted viral infections

Hepatitis E virus (HEV), a small, non-enveloped, RNA virus belonging to the Hepevirus genus (family Hepeviridae), is a human pathogen with a worldwide distribution.<sup>1</sup> The virus consists of an icosahedral particle that encloses a positive-sense, single-stranded RNA genome of 7.2kb.<sup>2</sup> Four major HEV genotypes, representing a single serotype, have been identified in humans and animals, including domestic pigs, wild boar, deer, and rodents.<sup>1,3,4</sup>

Molecular characterization of various HEV strains circulating among humans and animals has led to the recognition of four major genotypes.<sup>1</sup> Genotype 1, which occurs mainly in Asia, and genotype 2, which occurs in Africa and Mexico, are restricted to humans and transmitted via contaminated water in developing countries.<sup>1,5</sup> Genotypes 3 and 4 infect humans, pigs, and other mammalian species and cause sporadic cases of autochthonous HEV in both developing and developed countries.<sup>6</sup> Genotype 3 is the only genotype currently identified as the cause of autochthonous infection in the U.S.<sup>7</sup> and is the cause of the vast majority of infections in Europe, New Zealand, and North America.<sup>1,8-12</sup> Genotypes 3 and 4 are both present in Japan.<sup>1</sup> HEV genotypes 1, 3, and 4 are endemic to China.<sup>1</sup> Acute Hepatitis E is more common than Hepatitis A in China, France, the U.K., and Japan.<sup>1</sup>

The main mode of HEV transmission is the fecal-oral route through contaminated drinking water,<sup>1</sup> although food-borne transmission from consumption of undercooked or raw pork, organ meat, or shellfish, as well as zoonotic transmission as a result of contact with infected swine, domestic, or wild animals have been reported.<sup>7,13</sup> The full range of reservoirs for HEV is unknown.<sup>1</sup>

HEV infection usually causes a mild or subclinical infection with a self-limiting illness that lasts 4 to 6 weeks.<sup>1,8-11,17,18</sup> The symptoms are very similar to those of other forms of viral hepatitis infection, particularly Hepatitis A, with fatigue, jaundice, fever, malaise, nausea, vomiting, anorexia, and abdominal pain.<sup>1</sup> Patients often present with elevated alanine transaminase levels (~1,500 IU/L) and many present with jaundice.<sup>1</sup> HEV infections may occasionally be more severe and result in fulminant hepatic failure, particularly in pregnant woman, where the mortality rate can reach 10% to 25%; infants and children under 2 years of age, individuals with underlying liver disease (e.g., cirrhosis), and immunocompromised persons.<sup>1,19-23</sup> HEV infection causes more than 3 million symptomatic cases of acute Hepatitis E worldwide each year, which results in approximately 70,000 deaths annually.<sup>24</sup> The overall mortality rate ranges from 0.2% to 4.0%.<sup>1,25</sup> Most deaths from genotype 3 (HEV3) infection result from acute or subacute liver failure in patients with pre-existing liver disease, such as alcohol-related liver disease.<sup>1,22,26,27</sup>

HEV3 causes chronic infection, including up to 60% of infected immunosuppressed individuals, about 10% of whom develop cirrhosis.<sup>1</sup> Chronic infection is defined as persistent HEV RNA in serum or stool for 6 months or more.<sup>1</sup> Most cases occur in solid-organ transplant recipients, although infection in individuals with hematologic disorders receiving transfusions and chemotherapy and a few individuals with HIV have also been reported.<sup>28-34</sup> Chronic infection has not been reported with HEV1 or HEV2.<sup>1</sup>

HEV infection has also been associated with neurologic syndromes, including Guillain-Barré syndrome, Bell's palsy, acute transverse myelitis, acute meningoencephalitis, ataxia, and encephalitis; the neurologic symptoms typically resolve in patients who clear the virus.<sup>1</sup> Membranoproliferative and membranous glomerulonephritis, acute pancreatitis, and severe thrombocytopenia have been reported during acute HEV infection, although the pathophysiologic mechanisms and causal

relation, if any, have not been established.<sup>1</sup> Ribavirin therapy has been shown to be an effective treatment for acute severe HEV3 infection, and transplant recipients with chronic HEV infection are typically treated with reduction of immunosuppression (especially drugs that target T cells), interferon- $\alpha$ , and ribavirin.<sup>1</sup>

## Rationale for NAT testing

Like other hepatitises, HEV can be transmitted via transfusion of blood or blood products. Post-transfusion Hepatitis E has been reported in many countries.<sup>1,34-39</sup> The seropositive rate for HEV among the world's blood donors has been reported to vary from 0.4%–20.6%.<sup>40-48</sup> Asymptomatic HEV infections occur at a high rate around the world and, due to the prevalence of the virus, many blood donors may be infected and transmit the virus to recipients of their blood products. For example, a recent study of British blood donors showed 11% of donor sera was HEV IgG reactive, indicating past infection, and 0.7% of donor sera was IgM reactive, indicating acute infection.<sup>46</sup> In addition, 0.7% of plasma minipools from English donors contained HEV RNA.<sup>47</sup> A study of Chinese blood donors revealed similar results: 32.6% of donor sera was IgG reactive; 0.94% of donor sera was IgM reactive; and 0.07% of donations demonstrated HEV viremia.<sup>48</sup> A global investigation of plasma fractionation pools reported 10% of pools tested were HEV-RNA positive.<sup>48,49</sup>

## Explanation of the test

The cobas® HEV test is a qualitative PCR test for the detection of HEV RNA that is run on the cobas® 6800 System and cobas® 8800 System. The cobas® HEV test enables the simultaneous detection of HEV RNA and the internal control in a single test of an infected, individual donation or pooled plasma from individual donations.

## Principles of the procedure

The cobas® HEV test is based on fully automated sample preparation (nucleic acid extraction and purification) followed by PCR amplification and detection. The cobas® 6800/8800 Systems consist of the sample supply module, the transfer module, the processing module, and the analytic module. Automated data management is performed by the cobas® 6800/8800 software which assigns test results for all tests as non-reactive, reactive, or invalid. Results can be reviewed directly on the system screen, and printed as a report.

Samples can either be tested individually or tested in pools consisting of multiple samples. The cobas p 680 instrument may optionally be used in a pre-analytical step if pooling is to be performed.

Nucleic acids from the sample and added armored RNA internal control (IC) molecules (which serve as the sample preparation and amplification/detection process control) are simultaneously extracted. In addition the test utilizes two external controls: a positive and a negative control. Viral nucleic acids are released by addition of proteinase and lysis reagent to the sample. The released nucleic acids bind to the silica surface of the added magnetic glass particles. Unbound substances and impurities, such as denatured proteins, cellular debris, and potential PCR inhibitors (such as hemoglobin) are removed with subsequent wash reagent steps and purified nucleic acids are eluted from the glass particles with elution buffer at elevated temperature.

Selective amplification of target nucleic acid from the donor sample is achieved by the use of virus-specific forward and reverse primers which are selected from highly conserved regions of the viral nucleic acid. A thermostable DNA polymerase enzyme is used for both reverse-transcription and amplification. The master mix includes deoxyuridine triphosphate (dUTP), instead of deoxythymidine triphosphate (dTTP), which is incorporated into the newly synthesized DNA (amplicon).<sup>50-52</sup> Any contaminating amplicons from previous PCR runs are destroyed by the AmpErase enzyme [uracil-N-glycosylase], which is included in the PCR mix, when heated in the first thermal cycling step. However, newly formed amplicons are not destroyed since the AmpErase enzyme is inactivated once exposed to temperatures above 55°C.

The cobas® HEV master mix contains detection probes which are specific for HEV and IC nucleic acid. The specific HEV and IC detection probes are each labeled with one of two unique fluorescent dyes which act as a reporter. Each probe also has a second dye which acts as a quencher. The two reporter dyes are measured at defined wavelengths, thus permitting simultaneous detection and discrimination of the amplified HEV target and the IC.<sup>53,54</sup> The fluorescent signals of the intact probes are suppressed by the quencher dye. During the PCR amplification step, hybridization of the probes to the specific single-stranded DNA template results in cleavage by the 5' to 3' nuclease activity of the DNA polymerase resulting in separation of the reporter and quencher dyes and the generation of a fluorescent signal. With each PCR cycle, increasing amounts of cleaved probes are generated and the cumulative signal of the reporter dye is concomitantly increased. Since the two specific reporter dyes are measured at defined wavelengths, simultaneous detection and discrimination of the amplified HEV target and the IC are possible.

# Reagents and materials

## cobas® HEV reagents and controls

All unopened reagents and controls shall be stored as recommended in Table 1 to Table 4.



**Table 1** cobas® HEV test

cobas® HEV test Store at 2-8°C 96 test cassette (P/N 07001045190)		
Kit components	Reagent ingredients	Quantity per kit 96 tests
<b>Proteinase Solution (PASE)</b>	Tris buffer, < 0.05% EDTA, calcium chloride, calcium acetate, 8% (w/v) proteinase  EUH210: Safety data sheets available on request. EUH208: May produce an allergic reaction. Contains: Subtilisin, 9014-01-1	13 mL
<b>Internal Control (IC)</b>	Tris buffer, < 0.05% EDTA, < 0.001% internal control armored RNA construct (non-infectious RNA encapsulated in MS2 bacteriophage), < 0.002% Poly rA RNA (synthetic), < 0.1% sodium azide	13 mL
<b>Elution Buffer (EB)</b>	Tris buffer, 0.2% methyl-4 hydroxybenzoate	13 mL
<b>Master Mix Reagent 1 (MMX-R1)</b>	Manganese acetate, potassium hydroxide, < 0.1% sodium azide	5.5 mL
<b>HEV Master Mix Reagent 2 (HEV MMX-R2)</b>	Tricine buffer, potassium acetate, glycerol, 18% dimethyl sulfoxide, Tween 20, EDTA, < 0.06% dATP, dGTP, dCTP, < 0.14% dUTP, < 0.01% upstream and downstream HEV and internal control primers, < 0.01% fluorescent-labeled HEV probes, < 0.01% fluorescent-labeled internal control probe, < 0.01% oligonucleotide aptamer, < 0.01% ZO5D DNA polymerase, < 0.01% AmpErase (uracil-N-glycosylase) enzyme (microbial), < 0.1% sodium azide	6 mL

**Table 2** cobas® HEV Control Kit**cobas® HEV Control Kit**

Store at 2-8°C



(P/N 07001100190)

Kit components	Reagent ingredients	Quantity per kit	Safety symbol and warning
<b>HEV Positive Control (HEV (+) C)</b>	<p>&lt; 0.001% synthetic (armored) HEV RNA encapsulated in MS2 bacteriophage coat protein, normal human plasma, non-reactive by licensed tests for antibody to HCV, antibody to HIV-1/2, HBsAg, antibody to HBc; HIV-1 RNA, HIV-2 RNA, HCV RNA, HBV DNA, HEV RNA, WNV RNA, and CMV DNA not detectable by PCR methods.</p> <p>0.1% ProClin® 300 preservative</p>	16 mL (16 x 1 mL)	  <p>Warning</p> <p>H317: May cause an allergic skin reaction.</p> <p>P261: Avoid breathing dust/fumes/gas/mist/ vapors/spray.</p> <p>P272: Contaminated work clothing should not be allowed out of the workplace.</p> <p>P280: Wear protective gloves.</p> <p>P302+ P352: IF ON SKIN wash with plenty of soap and water.</p> <p>P333 + P313: If skin irritation or rash occurs: Get medical advice/attention.</p> <p>P363: Wash contaminated clothing before reuse.</p>

**Table 3** cobas® NHP Negative Control Kit**cobas® NHP Negative Control Kit**

Store at 2-8°C



(P/N 07002220190)

Kit components	Reagent ingredients	Quantity per kit	Safety symbol and warning
<b>Normal Human Plasma Negative Control (NHP-NC)</b>	<p>Normal human plasma, non-reactive by licensed tests for antibody to HCV, antibody to HIV-1/2, HBsAg, antibody to HBc; HIV-1 RNA, HIV-2 RNA, HCV RNA, HBV DNA, HEV RNA, WNV RNA, and CMV DNA not detectable by PCR methods.</p> <p>&lt; 0.1% ProClin® 300 preservative</p>	16 mL (16 x 1 mL)	<div>   </div> <p>Warning</p> <p>H317: May cause an allergic skin reaction.</p> <p>P261: Avoid breathing dust/fumes/gas/mist/vapors/spray.</p> <p>P272: Contaminated work clothing should not be allowed out of the workplace.</p> <p>P280: Wear protective gloves.</p> <p>P302+ P352: IF ON SKIN wash with plenty of soap and water.</p> <p>P333 + P313: If skin irritation or rash occurs: Get medical advice/attention.</p> <p>P363: Wash contaminated clothing before reuse.</p>



## cobas omni reagents for sample preparation

**Table 4** cobas omni reagents for sample preparation\*

Reagents	Reagent ingredients	Quantity per kit	Safety symbol and warning
<b>cobas omni MGP Reagent (MGP)</b> Store at 2–8°C (P/N 06997546190)	Magnetic glass particles, Tris buffer, 0.1% methyl-4 hydroxybenzoate, < 0.1% sodium azide	480 tests	Not applicable
<b>cobas omni Specimen Diluent (SPEC DIL)</b> Store at 2–8°C (P/N 06997511190)	Tris buffer, 0.1% methyl-4 hydroxybenzoate, < 0.1% sodium azide	4 x 875 mL	Not applicable
<b>cobas omni Lysis Reagent (LYS)</b> Store at 2–8°C (P/N 06997538190)	42.56% (w/w) guanidine thiocyanate, 5% (w/v) polydocanol, 2% (w/v) dithiothreitol, dihydro sodium citrate	4 x 875mL	  Danger H302: Harmful if swallowed. H318: Causes serious eye damage. H412: Harmful to aquatic life with long lasting effects. EUH032: Contact with acids liberates very toxic gas. P301 + P312: IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. P264: Wash skin thoroughly after handling. P270: Do not eat, drink or smoke when using this product. P273: Avoid release to the environment. P280: Wear protective gloves/eye protection/face protection. P305+P351+P338: IF IN EYES Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P310: Immediately call a POISON CENTER or doctor/Physician if you feel unwell. P330: Rinse mouth
<b>cobas omni Wash Reagent (WASH)</b> Store at 15–30°C (P/N 06997503190)	Sodium citrate dihydrate, 0.1% methyl-4 hydroxybenzoate	4.2L	Not applicable

\* These reagents are not included in the cobas® HEV test kit. See listing of additional materials required (Table 7).

## Reagent storage and handling requirements

Reagents shall be stored and will be handled as specified in Table 5 and Table 6.

When reagents are not loaded on the cobas® 6800/8800 Systems, store them at the corresponding temperature specified in Table 5.

**Table 5** Reagent storage (when reagent is not on the system)

Reagent	Storage temperature
cobas® HEV – Test 96	2–8°C
cobas® HEV Control Kit	2–8°C
cobas® NHP Negative Control Kit	2–8°C
cobas omni Lysis Reagent	2–8°C
cobas omni MGP Reagent	2–8°C
cobas omni Specimen Diluent	2–8°C
cobas omni Wash Reagent	15–30°C

Reagents loaded onto the cobas® 6800/8800 Systems are stored at appropriate temperatures and their expiration is monitored by the system. The system allows reagents to be used only if all of the conditions shown in Table 6 are met. The system automatically prevents use of expired reagents. Table 6 allows the user to understand the reagent handling conditions enforced by the cobas® 6800/8800 Systems.

**Table 6** Reagent expiry conditions enforced by the cobas® 6800/8800 Systems

Reagent	Kit expiration date	Open-kit stability	Number of runs for which this kit can be used	On-board stability (cumulative time on board outside refrigerator)
cobas® HEV – Test 96	Date not passed	30 days since first usage	Max 10 runs	Max 8 hours
cobas® HEV Control Kit	Date not passed	Not applicable	Not applicable	Max 8 hours
cobas® NHP Negative Control Kit	Date not passed	Not applicable	Not applicable	Max 8 hours
cobas omni Lysis Reagent	Date not passed	30 days since loading*	Not applicable	Not applicable
cobas omni MGP Reagent	Date not passed	30 days since loading*	Not applicable	Not applicable
cobas omni Specimen Diluent	Date not passed	30 days since loading*	Not applicable	Not applicable
cobas omni Wash Reagent	Date not passed	30 days since loading*	Not applicable	Not applicable

\* Time is measured from the first time that reagent is loaded onto the cobas® 6800/8800 Systems.

## Additional materials required

**Table 7** Material and consumables for use on **cobas®** 6800/8800 Systems

Material	P/N
<b>cobas omni</b> Processing Plate	05534917001
<b>cobas omni</b> Amplification Plate	05534941001
<b>cobas omni</b> Pipette Tips	05534925001
<b>cobas omni</b> Liquid Waste Container	07094388001
<b>cobas omni</b> Lysis Reagent	06997538190
<b>cobas omni</b> MGP Reagent	06997546190
<b>cobas omni</b> Specimen Diluent	06997511190
<b>cobas omni</b> Wash Reagent	06997503190
Solid Waste Bag	07435967001
Solid Waste Container	07094361001

## Instrumentation and software required

The **cobas®** 6800/8800 software and **cobas®** HEV analysis package shall be installed on the instrument(s). The Instrument Gateway (IG) server will be provided with the system.

**Table 8** Instrumentation

Equipment	P/N
<b>cobas®</b> 6800 System (Option Moveable)	05524245001 and 06379672001
<b>cobas®</b> 6800 System (Fix)	05524245001 and 06379664001
<b>cobas®</b> 8800 System	05412722001
<b>cobas p</b> 680 Instrument (optional for pipetting and pooling)	06578624001
Sample Supply Module	06301037001

Refer to the **cobas®** 6800/8800 Systems Operator's Manual and **cobas p** 680 instrument Operator's Manual for additional information for primary and secondary sample tubes accepted on the instruments.

Note: Contact your local Roche representative for a detailed order list for sample racks, racks for clotted tips and rack trays accepted on the instruments.

# Precautions and handling requirements

## Warnings and precautions

As with any test procedure, good laboratory practice is essential to the proper performance of this assay. Due to the high sensitivity of this test, care should be taken to keep reagents and amplification mixtures free of contamination.

- For in vitro diagnostic use only.
- All samples should be handled as if infectious, using good laboratory procedures as outlined in Biosafety in Microbiological and Biomedical Laboratories and in the CLSI Document M29-A3.<sup>38,39</sup> Only personnel proficient in handling infectious materials and the use of the cobas® HEV test, cobas® 6800/8800 Systems and cobas p 680 instrument (optional) should perform this procedure.
- All human-sourced materials should be considered potentially infectious and should be handled with universal precautions. If spillage occurs, immediately disinfect with a freshly prepared solution of 0.5% sodium hypochlorite in distilled or deionized water (dilute household bleach 1:10) or follow appropriate site procedures.
- cobas® HEV Control Kit and cobas® NHP Negative Control Kit contain plasma derived from human blood. The source material has been tested by licensed antibody tests and found non-reactive for the presence of antibody to HCV, antibody to HIV-1/2, HBsAg, and antibody to HBc. Testing of normal human plasma by PCR methods also showed no detectable HEV RNA, HIV-1 (Groups M and O) RNA, HIV-2 RNA, HCV RNA, HBV DNA, WNV RNA, or CMV DNA. No known test method can offer complete assurance that products derived from human blood will not transmit infectious agents.
- Do not freeze whole blood.
- The use of sterile disposable pipettes and nuclease-free pipette tips is recommended. Use only supplied or specified required consumables to ensure optimal test performance.
- Safety Data Sheets (SDS) are available on request from your local Roche representative.
- Closely follow procedures and guidelines provided to ensure that the test is performed correctly. Any deviation from the procedures and guidelines may affect optimal test performance.
- False positive results may occur if carryover of samples is not adequately controlled during sample handling and processing.

## Reagent handling

- Handle all reagents, controls, and samples according to good laboratory practice in order to prevent carryover of samples or controls.
- Before use, visually inspect each reagent cassette, diluent, lysis reagent, and wash reagent to ensure that there are no signs of leakage. If there is any evidence of leakage, do not use that material for testing.
- cobas omni Lysis Reagent contains guanidine thiocyanate, a potentially hazardous chemical. Avoid contact of reagents with the skin, eyes, or mucous membranes. If contact does occur, immediately wash with generous amounts of water; otherwise, burns can occur.
- cobas® HEV test kits, cobas omni MGP Reagent, and cobas omni Specimen Diluent contain sodium azide as a preservative. Avoid contact of reagents with the skin, eyes, or mucous membranes. If contact does occur, immediately wash with generous amounts of water; otherwise, burns can occur. If these reagents are spilled, dilute with water before wiping dry.
- Do not allow cobas omni Lysis Reagent, which contains guanidine thiocyanate, to contact sodium hypochlorite (bleach) solution. This mixture can produce a highly toxic gas.
- Dispose of all materials that have come in contact with samples and reagents in accordance with country, state, and local regulations.

## Good laboratory practice

- Do not pipette by mouth.
- Do not eat, drink, or smoke in designated work areas.
- Wear laboratory gloves, laboratory coats, and eye protection when handling samples and reagents. Gloves must be changed between handling samples and cobas® HEV test kits and cobas omni reagents to prevent contamination. Avoid contaminating gloves when handling samples and controls.
- Wash hands thoroughly after handling samples and kit reagents, and after removing the gloves.
- Thoroughly clean and disinfect all laboratory work surfaces with a freshly prepared solution of 0.5% sodium hypochlorite in distilled or deionized water (dilute household bleach 1:10). Follow by wiping the surface with 70% ethanol.
- If spills occur on the cobas® 6800/8800 instruments, follow the instructions in the cobas® 6800/8800 Systems Operator's Manual to properly clean and decontaminate the surface of instrument(s).

## Sample collection, transport, storage, and pooling

**Note:** Handle all samples and controls as if they are capable of transmitting infectious agents.


Store all donor samples at specified temperatures.

Sample stability is affected by elevated temperatures.

## Living donor samples

- Plasma collected in EDTA, CPD, CPDA1, CP2D and 4% Sodium Citrate (Apheresis plasma) anticoagulant may be used with the cobas® HEV test. Follow the sample collection tube/bag manufacturer instructions for handling and centrifugation.
- Blood collected in EDTA, CPD, CPDA1, CP2D anticoagulant or Becton-Dickinson EDTA Plasma Preparation Tubes (BD PPT™) may be stored for up to 12 days with the following conditions:
  - Samples must be centrifuged within 72 hours of draw.
  - For storage above 8°C, samples may be stored for 72 hours at up to 25°C, and up to 30°C for 24 hours during the 72 hours.

Other than noted above, samples are stored at 2-8°C. In addition, plasma separated from the cells may be stored for up to 30 days at  $\leq -18^{\circ}\text{C}$  with three freeze/thaw cycles. Refer to Figure 1.

**Figure 1** Sample storage conditions

- Apheresis plasma collected in 4% Sodium Citrate anticoagulant may be stored for up to 30 days at 2-8°C.
  - For storage above 8°C, samples may be stored for 72 hours at up to 25°C, and up to 30°C for 24 hours during the 72 hours.
- If samples are to be shipped, they should be packaged and labeled in compliance with applicable country and/or international regulations covering the transport of samples and etiologic agents.

# Instructions for use

## Automated sample pipetting and pooling (optional)

The cobas p 680 instrument is an optional component of the cobas® 6800/8800 Systems used for automated pipetting and pooling of aliquots of multiple primary samples into one pooled sample. Refer to the cobas p 680 instrument Operator's Manual for more information.


## Procedural notes

- Do not use cobas® HEV test reagents, cobas® HEV Control Kit, cobas® NHP Negative Control Kit or cobas omni reagents after their expiry dates.
- Do not reuse consumables. They are for one-time use only.
- Refer to the cobas® 6800/8800 Systems Operator's Manual for proper maintenance of instruments.

## Running the cobas® HEV test

The test procedure is described in detail in the cobas® 6800/8800 Systems Operator's Manual and the cobas p 680 instrument Operator's Manual. Figure 2 below summarizes the procedure.

**Figure 2** cobas® HEV test procedure



## Results

The cobas® 6800/8800 Systems automatically detect HEV RNA simultaneously for the samples and controls.

### Quality control and validity of results

- One negative control [(-) C] and one positive control [HEV (+) C] is processed with each batch.
- In the cobas® 6800/8800 software and/or report, check for flags and their associated results to ensure the batch validity.
- The batch is valid if no flags appear for both controls.

Invalidation of results is performed automatically by the cobas® 6800/8800 software based on negative and positive control failures.

### Control flags

**Table 9** Control flags for negative and positive controls

Negative Control	Flag	Result	Interpretation
(-) C	Q02	Invalid	The entire batch is assigned invalid if the result for the (-) C is invalid.
Positive Control	Flag	Result	Interpretation
HEV (+) C	Q02	Invalid	The entire batch is assigned invalid if the result for the HEV (+) C is invalid.

If the batch is invalid, repeat testing of the entire batch including samples and controls.



## Interpretation of results

For a valid batch, check each individual sample for flags in the cobas® 6800/8800 software and/or report. The result interpretation should be as follows:

- A valid batch may include both valid and invalid donor sample results dependent on flags obtained for the individual samples.
- Sample results are valid only if the respective positive controls and the negative control of the corresponding batch are valid.

Two parameters are measured simultaneously for each sample: HEV and the internal control. Final sample results for the cobas® HEV test are reported by the software. In addition to the overall results, individual target results will be displayed in the cobas® 6800/8800 software and should be interpreted as follows:

**Table 10** Target results for individual target result interpretation

Target results	Interpretation
HEV Non-Reactive	No target signal detected for HEV and IC signal detected.
HEV Reactive	Target signal detected for HEV and IC signal may be or may not be detected.
Invalid	Target and internal control signal not detected.

## Repeat testing of individual sample(s)

Sample tubes with a final result of Invalid for the target require repeat testing

## Procedural limitations

- The cobas® HEV test has been evaluated only for use in combination with the cobas® HEV Control Kit, cobas® NHP Negative Control Kit, cobas omni MGP Reagent, cobas omni Lysis Reagent, cobas omni Specimen Diluent, and cobas omni Wash Reagent for use on the cobas® 6800/8800 Systems.
- Reliable results depend on proper sample collection, storage and handling procedures.
- Do not use heparinized plasma with this test because heparin has been shown to inhibit PCR.
- Detection of HEV RNA is dependent on the number of virus particles present in the sample and may be affected by sample collection, storage and handling, patient factors (i.e., age, presence of symptoms), and/or stage of infection and pool size.
- Though rare, mutations within the highly conserved regions of a viral genome covered by the cobas® HEV test, may affect primers and/or probe binding resulting in the failure to detect presence of virus.
- Due to inherent differences between technologies, it is recommended that, prior to switching from one technology to the next, users perform method correlation studies in their laboratory to qualify technology differences. Users should follow their own specific policies/procedures.

# Non-clinical performance evaluation

## Key performance characteristics

### Living donor samples

### Limit of Detection (LoD)

#### WHO International Standard

The limit of detection (LoD) of the cobas® HEV test for HEV RNA was determined using the WHO International Standard for HEV (PEI code 6329/10).

For the WHO International standard, 3 independent dilution series of the viral standard were prepared with normal, virus-negative (HEV) human EDTA-plasma. Each dilution series was tested using 3 different lots of the cobas® HEVtest kits with approximately 63 replicates per lot, for a total of approximately 189 replicates per concentration. For HEV virus, PROBIT analysis on the data combined across dilution series and reagent lots was used to estimate the LoD, along with the lower and upper limit of 95% confidence interval (Table 11). The reactivity rates observed in the LoD studies for HEV are summarized in Table 12.

**Table 11** Results of PROBIT analysis on LoD data collected with viral standard in EDTA plasma

Analyte	Measuring units	LoD	Lower 95% confidence limit	Upper 95% confidence limit
HEV	IU/mL	18.6	15.9	22.6

**Table 12** Reactivity rates summary for HEV in EDTA plasma

HEV RNA concentration (IU/mL)	Number reactive	Number of valid replicates	% Reactive	95% Lower confidence bound (one-sided)
40	187	187	100.0%	98.4%
20	179	188	95.2%	91.8%
10	165	189	87.3%	82.6%
6	113	187	60.4%	54.2%
2	52	189	27.5%	22.2%

## Reproducibility

The reproducibility of the cobas® HEV test on the cobas® 6800/8800 Systems was determined using the WHO International Standard for HEV (PEI code 6329/10). This study consisted of testing 3 panels of HEV at concentrations of approximately 0.5 x, 1 x and 2 x the LoD of the cobas® HEV test. Testing was performed for the following variability components:

- day-to-day variability over 3 days
- lot-to-lot variability using 3 different reagent lots of the cobas® HEV test
- instrument-to-instrument variability using 3 different cobas® 8800 Systems

Approximately 21 replicates were tested with each of the 3 panels for total of 63 replicates with each reagent lot. All valid reproducibility data were evaluated by calculating the percentage of reactive test results for each concentration level across all variable components.

The limits of two-sided 95% Confidence Intervals for each Reactive Rate were calculated for each of the three levels of HEV tested across 3 days, 3 reagent lots, and 3 cobas® 8800 Systems. The cobas® HEV test is reproducible over multiple days, reagent lots and multiple instruments. The results from reagent lot-to-lot variability are summarized in Table 13.

**Table 13** cobas® HEV test reagent lot-to-lot reproducibility summary

Analyte	Concentration	Reagent lot	% Reactive (reactive/valid replicates)	Lower limit of 95% confidence interval	Upper limit of 95% confidence interval
HEV	2 x LoD	1	100.0% (61/61)	94.1%	100.0%
		2	100.0% (63/63)	94.3%	100.0%
		3	100.0% (63/63)	94.3%	100.0%
	1 x LoD	1	88.9% (56/63)	78.4%	95.4%
		2	96.8% (60/62)	88.8%	99.6%
		3	100.0% (63/63)	94.3%	100.0%
	0.5 x LoD	1	82.5% (52/63)	70.9%	90.9%
		2	95.2% (60/63)	86.7%	99.0%
		3	84.1% (53/63)	72.7%	92.1%

## Genotype verification

The performance of the cobas® HEV test to detect 4 genotypes of HEV was determined by testing a total of 16 unique clinical samples and 7 HEV cultured isolates with known genotypes. All samples were quantified traceable to the HEV WHO standard. All 16 clinical samples were tested after dilution with normal, virus-negative (HEV) human EDTA-plasma to 5 x LoD of the cobas® HEV test, of which 10 samples were also tested neat. All 7 cultured isolates were tested after dilution with normal, virus-negative (HEV) human EDTA-plasma to 5 x LoD of the cobas® HEV test. All clinical samples and cultured isolates were detected at neat and/or at 5 x LoD (Table 14).

**Table 14** HEV clinical samples and cultured isolates

Genotype	Clinical samples		Cultured isolates
	% Reactive (reactive/samples tested) neat	% Reactive (reactive/samples tested) diluted to 5 x LoD	% Reactive (reactive/samples tested) diluted to 5 x LoD
1	Not tested*	Not tested*	100.0% (3/3)
2	Not tested*	Not tested*	100.0% (1/1)
3	100.0% (10/10)	100.0% (10/10)	Not tested*
4	Not tested*	100.0% (6/6)	100.0% (3/3)

\*Insufficient volume to test at neat/diluted

## Analytical specificity

The analytical specificity of the cobas® HEV test was evaluated for cross-reactivity with 28 microorganisms at  $10^6$  particles, copies, or PFU/mL, which included 21 viral isolates, 6 bacterial strains and 1 yeast isolate (Table 15). The microorganisms were added to normal, virus-negative human EDTA-plasma and tested with and without HEV added to a concentration of approximately 3 x LoD of the cobas® HEV test for each virus. The tested microorganisms do not cross-react or interfere with the cobas® HEV test.

**Table 15** Microorganisms tested for analytical specificity

Viruses	Flavivirus	Bacteria	Yeast
Adenovirus 5	West Nile Virus	<i>Escherichia coli</i>	<i>Candida albicans</i>
Cytomegalovirus	Dengue Virus type 1	<i>Propionibacterium acnes</i>	
Epstein-Barr Virus	Usutu Virus	<i>Staphylococcus aureus</i>	
Herpes Simplex Virus type 1		<i>Staphylococcus epidermidis</i>	
Herpes Simplex Virus type 2		<i>Streptococcus viridans</i>	
Hepatitis A Virus		<i>Staphylococcus haemolyticus</i>	
Hepatitis B Virus			
Hepatitis C Virus			
Hepatitis G Virus			
Human Immunodeficiency Virus (HIV-1 Group M)			
Human Immunodeficiency Virus (HIV-2)			
Human T-cell lymphotropic Virus type I			
Human T-cell lymphotropic Virus type II			
Human Herpes Virus 6			
Influenza Virus A			
Parvovirus B19			
Chikungunya Virus			
Varicella Zoster Virus			

Plasma samples from each of the disease states (Table 16) were tested without HEV and with HEV added to a concentration of approximately 3 x LoD of the cobas® HEV test. These disease states do not cross-react or interfere with the cobas® HEV test.

**Table 16** Disease states samples tested for analytical specificity

Disease state		
Adenovirus type 5	Hepatitis B Virus	Human T-cell lymphotropic Virus type I
Cytomegalovirus	Hepatitis C Virus	Human T-cell lymphotropic Virus type II
Dengue Virus	Herpes Simplex Virus type 1	Parvovirus B19
Epstein-Barr Virus	Herpes Simplex Virus type 2	West Nile Virus
Hepatitis A Virus	Human Immunodeficiency Virus (HIV-1)	

## Analytical specificity – interfering substances

### Endogenous interference substances

Plasma samples with abnormally high levels of triglycerides (up to 33.2 g/L), hemoglobin (up to 4.7 g/L), unconjugated bilirubin (up to 0.28 g/L), albumin (up to 60 g/L), and human DNA (up to 0.004 g/L) were tested with and without HEV added to a concentration of approximately 3 x LoD of the cobas® HEV test. Samples containing these endogenous substances did not interfere with the sensitivity or specificity of cobas® HEV test.

### Exogenous interference substances

Normal, virus-negative (HEV) human EDTA-plasma samples containing abnormally high concentrations of drugs (Table 17) were tested with and without HEV added to a concentration of 3 x LoD of the cobas® HEV test. These exogenous substances did not interfere with the sensitivity or specificity of the cobas® HEV test.

**Table 17** Clinical samples tested with drugs

Name of drug tested	Concentration
Acetaminophen	1324 µmol/L
Acetylsalicylic Acid	3620 µmol /L
Ascorbic Acid	342 µmol/L
Atorvastatin	600 µg Eq/L
Fluoxetine	11.2 µmol/L
Ibuprofen	2425 µmol/L
Loratadine	0.78 µmol/L
Nadolol	3.88 µmol/L
Naproxen	2170 µmol/L
Paroxetine	3.04 µmol/L
Phenylephrine HCL	491 µmol/L
Sertraline	1.96 µmol/L

## Correlation

### Performance evaluation of the cobas® HEV test compared to the Realstar HEV RT-PCR Kit 1.0 test

The performance of the cobas® HEV test and the Realstar® HEV RT-PCR Kit 1.0 test (Altona Diagnostics) was compared using 100 individual HEV NAT-positive plasma samples. 100 positive samples were tested neat and 67 positive samples were tested diluted 1:6. In addition, 100 HEV negative plasma samples were tested neat with both methods.

The seronegative samples demonstrated 100% specificity by generating 100 out of 100 non-reactive results with both methods.

For positive samples, both methods were in agreement based on the McNemars's test indicating that the performance of cobas® HEV test and Realstar® HEV RT-PCR Kit 1.0 test are equivalent (Table 18).

**Table 18** Correlation of positive samples (neat)

Methods		HEV results	
Realstar® HEV RT-PCR Kit 1.0 test	cobas® HEV	Neat	Diluted 1:6
Non-reactive	Non-reactive	0	3
Reactive	Non-Reactive	0	3
Non-reactive	Reactive	1	9
Reactive	Reactive	99	52
Total		100	67
McNemar's Test, p-value (two-sided, $\alpha=0.05$ )		1.00	0.09

## Whole System Failure

The Whole System Failure rate for the cobas® HEV test was determined by testing 100 replicates of EDTA plasma spiked with HEV. These samples were tested at a target concentration of approximately 3 x LoD and were run in pools of 1 (undiluted). The study was performed using the cobas® 8800 System with cobas p 680 instrument (pipetting and pooling).

The results of this study determined that all replicates were reactive for HEV, resulting in a Whole System Failure rate of 0%. The two-sided 95% exact confidence interval was 0% for the lower bound and 3.62% for the upper bound [0%: 3.62%].

## Additional information

### Key test features

<b>Sample type</b>	Plasma
<b>Amount of sample required</b>	1000 µL
<b>Amount of sample processed</b>	850 µL
<b>Test duration</b>	Results are available within less than 3.5 hours after loading the sample on the system.



The following symbols are used in labeling for Roche PCR diagnostic products.

**Table 19** Symbols used in labeling for Roche PCR diagnostics products



Ancillary Software



*In-Vitro*-Diagnostic Medical Device



Authorized Representative  
in the European community



Lower Limit of Assigned Range



Barcode Data Sheet



Manufacturer



Batch code



Store in the dark



Biological Risks



Contains Sufficient for <*n*> tests



Catalogue number



Temperature Limit



Consult instructions for use



Test Definition File



Contents of kit



Upper Limit of Assigned Range



Distributed by



Use By



For IVD Performance Evaluation Only



This product fulfills the requirements of the European Directive 98/79 EC for *in vitro* diagnostic medical devices.

US Customer Technical Support 1-800-526-1247

## Manufacturer and distributors

**Table 20** Manufacturer and distributors

	<p>Manufactured in the United States</p> <p>Roche Diagnostics GmbH Sandhofer Straße 116 68305 Mannheim, Germany</p>	
	<p>Roche Diagnostics (Schweiz) AG Industriestrasse 7 6343 Rotkreuz, Switzerland</p> <p>Roche Diagnostics GmbH Sandhofer Straße 116 68305 Mannheim, Germany</p> <p>Roche Diagnostics, SL Avda. Generalitat, 171-173 E-08174 Sant Cugat del Vallès Barcelona, Spain</p> <p>Roche Diagnostica Brasil Ltda. Av. Engenheiro Billings, 1729 Jaguaré, Building 10 05321-010 São Paulo, SP Brazil</p>	<p>Roche Diagnostics 201, boulevard Armand-Frappier H7V 4A2 Laval, Québec, Canada (For Technical Assistance call: Pour toute assistance technique, appeler le: 1-877 273-3433)</p> <p>Roche Diagnostics 2, Avenue du Vercors 38240 Meylan, France</p> <p>Distributore in Italia: Roche Diagnostics S.p.A. Viale G. B. Stucchi 110 20052 Monza, Milano, Italy</p> <p>Distribuidor em Portugal: Roche Sistemas de Diagnósticos Lda. Estrada Nacional, 249-1 2720-413 Amadora, Portugal</p>

## Trademarks and patents

See <http://www.roche-diagnostics.us/patents>

## Copyright

©2014, Roche Molecular Systems, Inc.



## References

1. Kamar N, Bendall R, Legrand-Abravanel F, et al. Hepatitis E. *Lancet*. 2012;379:2477-2488.
2. Ahmad I, Holla RP, Jameel S. Molecular virology of hepatitis E virus. *Virus Res*. 2011;161:47-58.
3. Meng XJ, Purcell RH, Halbur PG, et al. A novel virus in swine is closely related to the human hepatitis E virus. *Proc Natl Acad Sci USA*. 1997;94:9860-9865.
4. Yugo DM, Meng XJ. Hepatitis E virus: foodborne, waterborne, and zoonotic transmission. *Int J Environ Res Public Health*. 2013;10:4507-4533.
5. Balayan MS, Andjaparidze AG, Savinskaya SS, et al. Evidence for a virus in non-A, non-B hepatitis transmitted via the fecal-oral route. *Intervirology*. 1983;20:23-31.
6. Dalton HR, Bendall R, Ijaz S, Banks M. Hepatitis E: an emerging infection in developed countries. *Lancet Infect Dis*. 2008;8:698-709.
7. Meng XJ. Hepatitis E as a zoonotic disease in the United States. Oral presentation by X.J. Meng at the NIH Research Workshop: Hepatitis E in the United States, Bethesda, Maryland. March 2012;7-8.  
<http://www3.niddk.nih.gov/fund/other/HepE2012/HepatitisE2012ProgramBook.pdf>
8. Mansuy JM, Péron JM, Abravanel F, et al. Hepatitis E in the south west of France in individuals who have never visited an endemic area. *J Med Virol*. 2004;74:419-424.
9. Dalton HR, Stableforth W, Thuraiajah P, et al. Autochthonous hepatitis E in Southwest England: natural history, complications and seasonal variation, and hepatitis E virus IgG seroprevalence in blood donors, the elderly and patients with chronic liver disease. *Eur J Gastroenterol Hepatol*. 2008;20:784-790.
10. Wichmann O, Schimanski S, Koch J, et al. Phylogenetic and case-control study on hepatitis E virus in Germany. *J Infect Dis*. 2008;198:1732-1741.
11. Dalton HR, Fellows HJ, Gane EJ, et al. Hepatitis E in New Zealand. *J Gastroenterol Hepatol*. 2007;22:1236-1240.
12. Tsang TH, Denison EK, Williams HV, Venczel LV, Ginsberg MM, Vugia DJ. Acute hepatitis E infection acquired in California. *Clin Infect Dis*. 2000;30:618-619.
13. Teo CG. Much meat, much malady: changing perceptions of the epidemiology of hepatitis E. *Clin Microbiol Infect*. 2010;16:24-32.
14. Colson P, Borentain P, Queyriaux B, et al. Pig liver sausage as a source of hepatitis E virus transmission to humans. *J Infect Dis*. 2010;202:825-834.
15. Matsuda H, Okada K, Takahashi K, Mishoro S. Severe hepatitis E virus infection after ingestion of uncooked liver from a wild boar. *J Infect Dis*. 2003;188:944.
16. Tei S, Kitajima N, Takahashi K, Mishoro S. Zoonotic transmission of hepatitis E virus from deer to human beings. *Lancet* 2003;362:371-373.
17. Gotanda Y, Iwata A, Ohnuma H, et al. Ongoing subclinical infection of hepatitis E virus among blood donors with an elevated alanine aminotransferase level in Japan. *J Med Virol*. 2007;79:734-742.
18. Mitsui T, Tsukamoto Y, Suzuki S, et al. Serological and molecular studies on subclinical hepatitis E virus infection using periodic serum samples obtained from healthy individuals. *J Med Virol* 2005;76:526-533.
19. Nelson KE, Kmush B, Labrique AB. The epidemiology of hepatitis E virus infections in developed countries and among immunocompromised patients. *Expert Rev Anti Infect Ther*. 2011;9:1133-1148.
20. Navaneethan U, Al Mohajer M, Shata MT. Hepatitis E and pregnancy: understanding the pathogenesis. *Liver Int*. 2008;28:1190-1199.
21. Jilani N, Das BC, Husain SA, et al. Hepatitis E virus infection and fulminant hepatic failure during pregnancy. *J Gastroenterol Hepatol*. 2007;22:676-682.
22. Péron JM, Beureau C, Poirson H, et al. Fulminant liver failure from acute autochthonous hepatitis E in France: description of seven patients with acute hepatitis E and encephalopathy. *J Viral Hepat*. 2007;14:298-303.

23. Ramachandran J, Eapen CE, Kang G, et al. Hepatitis E superinfection produces severe decompensation in patients with chronic liver disease. *J Gastroenterol Hepatol*. 2004;19:134-138.
24. Rein DB, Stevens GA, Theaker J, et al. The global burden of hepatitis E virus genotypes 1 and 2 in 2005. *Hepatology*. 2012;55:988-997.
25. World Health Organization, Department of Communicable Disease Surveillance and Response. Hepatitis E. [http://www.who.int/scr/disease/hepatitis/Hepatitis\\_E\\_whoedscrdc2001\\_12.pdf](http://www.who.int/scr/disease/hepatitis/Hepatitis_E_whoedscrdc2001_12.pdf) (reviewed 13 Feb 2014).
26. Dalton HR, Bendall RP, Rashid M, et al. Host risk factors and autochthonous hepatitis E infection. *Eur J Gastroenterol Hepatol*. 2011;23:1200-1205.
27. Dalton HR, Hazeldine S, Banks M, Ijaz S, Bendall R. Locally acquired hepatitis E in chronic liver disease. *Lancet* 2007;369:1260.
28. Kamar N, Selves J, Mansuy JM, et al. Hepatitis E virus and chronic hepatitis in organ-transplant recipients. *N Engl J Med*. 2008;358:811-817.
29. Gérolami R, Moal V, Colson P, et al. Chronic hepatitis E with cirrhosis in a kidney-transplant patient. *N Engl J Med*. 2008;358:859-860.
30. Haagsma EB, van den Berg AP, Porte RJ, et al. Chronic hepatitis E virus infection in liver transplant recipients. *Liver Transpl*. 2008;14:547-553.
31. Kamar N, Garrouste C, Haagsma EB, et al. Factors associated with chronic hepatitis in patients with hepatitis E virus infection who have received solid organ transplants. *Gastroenterology*. 2011;140:1481-1489.
32. Colson P, Kaba M, Moreau J, et al. Hepatitis E in an HIV-infected patient. *J Clin Virol*. 2009;45:269-271.
33. Péron JM, Mansuy JM, Récher C, et al. Prolonged hepatitis E in an immunocompromised patient. *J Gastroenterol Hepatol*. 2006;21:1223-1224.
34. Tamura A, Shimizu YK, Tanaka T, et al. Persistent infection of hepatitis E virus transmitted by blood transfusion in a patient with T-cell lymphoma. *Hepatol Res*. 2007;37:113-120.
35. Matsubayashi K, Nagaoka Y, Sakata H, et al. Transfusion-transmitted hepatitis E caused by apparently indigenous hepatitis E virus strain in Hokkaido, Japan. *Transfusion*. 2004;44:934-940.
36. Mitsui T, Tsukamoto Y, Yamazaki C, et al. Prevalence of hepatitis E virus infection among hemodialysis patients in Japan: evidence for infection with a genotype 3 HEV by blood transfusion. *J Med Virol*. 2004;74:563-572.
37. Boxall E, Herborn A, Kochethu G, et al. Transfusion-transmitted hepatitis E in a “nonhyperendemic” country. *Transfus Med*. 2006;16:79-83.
38. Colson P, Coze C, Gallian P, et al. Transfusion-associated hepatitis E, France. *Emerg Infect Dis*. 2007;13:648-649.
39. Khuroo MS, Kamili S, Yattoo GN. Hepatitis E virus infection may be transmitted through blood transfusions in an endemic area. *J Gastroenterol Hepatol*. 2004;19:778-784.
40. Zaaijer HL, Kok M, Lelie PN, Timmerman RJ, Chau K, van der Pal HJ. Hepatitis E in The Netherlands: imported and endemic. *Lancet*. 1993;341:826.
41. Christensen PB, Engle RE, Hjort C, et al. Time trend of the prevalence of hepatitis E antibodies among farmers and blood donors: a potential zoonosis in Denmark. *Clin Infect Dis*. 2008;47:1026-1031.
42. Boutrouille A, Bakkali-Kassimi L, Crucièrè C, et al. Prevalence of anti-hepatitis E virus antibodies in French blood donors. *J Clin Microbiol*. 2007;45:2009-2010.
43. Kaufmann A, Kenfak-Foguena A, André C, et al. Hepatitis E virus seroprevalence among blood donors in southwest Switzerland. *PLoS One*. 2011;6:e21150.
44. Vollmer T, Diekmann J, Johne R, et al. Novel approach for detection of hepatitis E virus infection in German blood donors. *J Clin Microbiol*. 2012;50:2708-2713.
45. Takeda H, Matsubayashi K, Sakata H, et al. A nationwide survey for prevalence of hepatitis E virus antibody in qualified blood donors in Japan. *Vox Sang*. 2010;99:307-313.

46. Beale MA, Tettmar K, Szypulska R, Tedder RS, Ijaz S. Is there evidence of recent hepatitis E virus infection in English and North Welsh blood donors? *Vox Sang*. 2011;100:340-342.
47. Ijaz S, Szypulska R, Tettmar KI, Kitchen A, Teddler RS. Detection of hepatitis E virus RNA in plasma mini-pools from blood donors in England. *Vox Sang*. 2012;102:272.
48. Guo QS, Yan Q, Xiong JH, et al. Prevalence of hepatitis E virus in Chinese blood donors. *J Clin Microbiol*. 2010;48:317-318.
49. Baylis SA, Hanschmann KM, Blümel J, Nubling CM; HEV Collaborative Study Group. Standardization of hepatitis E virus (HEV) nucleic acid amplification technique-based assays: an initial study to evaluate a panel of HEV strains and investigate laboratory performance. *J Clin Microbiol*. 2011;49:1234-1239.
50. Longo MC, Berninger MS, Hartley JL. Use of uracil DNA glycosylase to control carry-over contamination in polymerase chain reactions. *Gene*. 1990;93:125-128.
51. Savva R, McAuley-Hecht K, Brown T, Pearl L. The structural basis of specific base-excision repair by uracil-DNA glycosylase. *Nature*. 1995;373:487-493.
52. Mol CD, Arvai AS, Slupphaug G, et al. Crystal structure and mutational analysis of human uracil-DNA glycosylase: structural basis for specificity and catalysis. *Cell*. 1995;80:869-878.
53. Higuchi R, Dollinger G, Walsh PS, Griffith R. Simultaneous amplification and detection of specific DNA sequences. *Biotechnology (NY)*. 1992;10:413-417.
54. Heid CA, Stevens J, Livak JK, Williams PM. Real-time quantitative PCR. *Genome Res* 1996;986-994.
55. Center for Disease Control and Prevention. Biosafety in Microbiological and Biomedical Laboratories, 5th ed. U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, National Institutes of Health HHS Publication No. (CDC) 21-1112, revised December 2009.
56. Clinical and Laboratory Standards Institute (CLSI). Protection of Laboratory Workers from Occupationally Acquired Infections. Approved Guideline-Third Edition. CLSI Document M29-A3 Villanova, PA: CLSI, 2005.

## Document revision

Document Revision Information	
Doc Rev. 1.0 07/2014	First Publishing.