

For general laboratory use. Not for use in diagnostic procedures. FOR *IN VITRO* USE ONLY.



*Roche Applied Science*

---

# LightCycler<sup>®</sup> TaqMan<sup>®</sup> Master

---

Version March 2006

Ready-to-use hot start reaction mix for PCR on the LightCycler<sup>®</sup> 2.0 System\*  
using Hydrolysis (TaqMan<sup>®</sup>) Probes

**Cat. No. 04 535 286 001**

Kit for 96 reactions

**Store the kit at –15 to –25°C**

Ⓢ \*LightCycler<sup>®</sup> 2.0 Instrument,  
LightCycler<sup>®</sup> 1.5 Instrument or an  
instrument version below!

# Table of Contents

---

<b>1. What this Product Does</b> .....	<b>3</b>
Number of Tests	3
Kit Contents	3
Storage and Stability	3
Additional Equipment and Reagents Required	4
Application	4
Assay Time	5
<b>2. How to Use this Product</b> .....	<b>6</b>
2.1 Before You Begin	6
Sample Material	6
Negative Control	6
Primers	6
Hydrolysis Probes	6
MgCl <sub>2</sub>	6
2.2 Procedure	7
LightCycler® 2.0 System Protocol	7
Fluorescence and Run Setup Parameters	8
Preparation of the Master Mix	9
Preparation of the PCR Mix	10
2.3 Related Procedures	11
Prevention of Carry-Over Contamination	11
Two-step RT-PCR	11
<b>3. Results</b> .....	<b>12</b>
<b>4. Troubleshooting</b> .....	<b>13</b>
<b>5. Additional Information on this Product</b> .....	<b>15</b>
5.1 How this Product Works	15
Test principle	15
References	15
Quality Control	15
<b>6. Supplementary Information</b> .....	<b>16</b>
6.1 Conventions	16
Text Conventions	16
Symbols	16
6.2 Changes to Previous Version	16
6.3 Ordering Information	17
6.4 License Disclaimer	18
6.5 Trademarks	18

# 1. What this Product Does

---

**Number of Tests** The kit is designed for 96 reactions with a final reaction volume of 20  $\mu$ l each.

## Kit Contents

Vial/Cap	Label	Contents/Function
1a white cap	Enzyme	<ul style="list-style-type: none"><li>• 1 <math>\times</math> vial 1a, 3 <math>\times</math> vial 1b, for 3 <math>\times</math> 128 <math>\mu</math>l Master Mix (5 <math>\times</math> conc.)</li><li>• Ready-to-use hot start PCR reaction mix (after pipetting 10 <math>\mu</math>l from vial 1a into one vial 1b)</li><li>• Contains FastStart Taq DNA Polymerase, reaction buffer, MgCl<sub>2</sub> and dNTP mix (with dUTP instead of dTTP)</li></ul>
1b red cap	Reaction Mix	
2 colorless cap	Water, PCR grade	<ul style="list-style-type: none"><li>• 2 <math>\times</math> 1 ml</li><li>• to adjust the final reaction volume</li></ul>

## Storage and Stability

- Store the kit at  $-15$  to  $-25^{\circ}\text{C}$  through the expiration date printed on the label.
- The kit is shipped on dry ice.
- Once the kit is opened, store the kit components as described in the following table:

Vial	Label	Storage
1a white cap	Enzyme	<ul style="list-style-type: none"><li>• Store at <math>-15</math> to <math>-25^{\circ}\text{C}</math>.</li><li>• <b>Avoid repeated freezing and thawing!</b></li><li>• <b>Keep vial 1b away from light!</b></li></ul>
1b red cap	Reaction Mix	
1 red cap (after addition of 1a to 1b)	Master Mix	<ul style="list-style-type: none"><li>• Store at <math>-15^{\circ}\text{C}</math> to <math>-25^{\circ}\text{C}</math> for a maximum of 3 months.</li><li>• After thawing, store at <math>+2</math> to <math>+8^{\circ}\text{C}</math> for a maximum of one week.</li><li>• <b>Avoid repeated freezing and thawing!</b></li></ul>
2 colorless cap	Water, PCR-grade	<ul style="list-style-type: none"><li>• Store at <math>-15</math> to <math>-25^{\circ}\text{C}</math>.</li></ul>

## 1. What this Product Does, continued

---

### Additional Equipment and Reagents Required

Refer to the list below for additional reagents and equipment required to perform PCR reactions with the LightCycler<sup>®</sup> TaqMan<sup>®</sup> Master using the LightCycler<sup>®</sup> 2.0 System:

- LightCycler<sup>®</sup> 2.0 System\* (LightCycler<sup>®</sup> 2.0 Instrument\*, LightCycler<sup>®</sup> 1.5 Instrument\*, or an instrument version below)
- LightCycler<sup>®</sup> Capillaries\* (20  $\mu$ l)
- Standard benchtop microcentrifuge containing a rotor for 2.0 ml reaction tubes
  - Ⓢ The LightCycler<sup>®</sup> 2.0 System provides adapters that allow LightCycler<sup>®</sup> Capillaries to be centrifuged in a standard microcentrifuge rotor.

or

- Ⓢ LC Carousel Centrifuge 2.0\* for use with the LightCycler<sup>®</sup> 2.0 Sample Carousel (optional)
  - Ⓢ If you use a LightCycler<sup>®</sup> Instrument version below 2.0, you need in addition the LC Carousel Centrifuge 2.0 Bucket 2.1\*. To adapt the LightCycler<sup>®</sup> 2.0 Sample Carousel to the former LC Carousel Centrifuge, you need the LC Carousel Centrifuge 2.0 Rotor Set\*.
- PCR template (genomic DNA or cDNA)
- PCR primers
- Hydrolysis probe
- LightCycler<sup>®</sup> Uracil-DNA Glycosylase\* (optional<sup>†</sup>)
- Nuclease-free, aerosol-resistant pipette tips
- Pipettes with disposable, positive-displacement tips
- Sterile reaction tubes for preparing master mixes and dilutions
- Ⓢ <sup>#</sup> If you want to perform color compensation when using the LightCycler<sup>®</sup> Red 640 and 705-labeled HybProbe pairs in dual color experiments in the same capillary, see Related Procedures section for details.
- Ⓢ <sup>†</sup>For prevention of carry-over contamination; see Related Procedures section for details.

*\*available from Roche Applied Science; see Ordering Information for details*

## 1. What this Product Does, continued

---

### Application

The LightCycler® TaqMan® Master is designed for research studies. In combination with the LightCycler® 2.0 System, the kit enables high sensitive detection and quantification of defined DNA sequences (with suitable PCR primers and detection probes). It can also be used to detect and quantify defined RNA sequences in a two-step RT-PCR (with additional reagents for reverse transcription).

To prevent false-positives due to contamination with amplified DNA, the kit can be used with LightCycler® Uracil-DNA Glycosylase.

In principle, the kit can be used for the amplification and detection of any DNA or cDNA target. However, you would have to adapt each amplification protocol to the reaction conditions of the LightCycler® 2.0 System and design a pair of specific PCR primers and a Hydrolysis probe for each target.

- ⚠ The amplicon size should not exceed 1 kb in length. For optimum results, select a product length of 700 bp or less.
- ⚠ The performance of the kit described in this Instruction Manual is warranted only when it is used with the LightCycler® 2.0 System.

### Assay Time

<b>Procedure</b>	<b>Time</b>
PCR Set-up	20 min
LightCycler® 2.0 System PCR run	40 min
<b>Total assay time</b>	<b>60 min</b>

## 2. How to Use this Product

---

### 2.1 Before You Begin

**Sample Material** Use any template DNA (*e.g.*, genomic or plasmid DNA, cDNA) suitable for PCR in terms of purity, concentration, and absence of inhibitors. For reproducible isolation of nucleic acids use:

- either the MagNA Pure LC Instrument and a dedicated MagNA Pure LC reagent kit (for automated isolation) or
- a High Pure Nucleic Acid Isolation Kit (for manual isolation).

For details see the Roche Applied Science catalogue or the website: [www.roche-applied-science.com](http://www.roche-applied-science.com)

Ⓢ If you are using a non-purified cDNA sample from reverse transcription, especially if it contains high background concentrations of RNA and oligonucleotides, you can improve your results by using 2  $\mu$ l or less of sample in the reaction.

**Negative Control** Always run a negative control with the samples. To prepare a negative control, replace the template DNA with PCR-grade water (vial 2, colorless cap).

**Primers** Use PCR primers at a final concentration of 0.1 – 1  $\mu$ M. The recommended starting concentration is 0.5  $\mu$ M each.

**Hydrolysis Probes** Use Hydrolysis probes at a final concentration of 0.05 – 0.1  $\mu$ M. In some cases the fluorescence signal can be increased with a higher probe concentration up to 0.2  $\mu$ M.

**MgCl<sub>2</sub>** All components in the reaction mix of the LightCycler® TaqMan® Master are optimized for almost all primer combinations. You do not need to add additional MgCl<sub>2</sub> to the mix to get efficient and specific PCR.

### LightCycler® 2.0 System Protocol

The following procedure is optimized for use with the LightCycler® 2.0 System.

⚠ Program the LightCycler® Instrument before preparing the reaction mixes. A LightCycler® 2.0 System protocol that uses LightCycler® TaqMan® Master contains the following programs:

- **Pre-Incubation** (activation of FastStart DNA polymerase and denaturation of the DNA)
- **Amplification** of the target DNA
- **Cooling** the rotor and thermal chamber

For details on how to program the experimental protocol, see the LightCycler® 2.0 System Operator's Manual.

⚠ Set all other protocol parameters not listed in the tables below to '0'.

The following table shows the PCR parameters that must be programmed for a LightCycler® 2.0 System PCR run with the LightCycler® TaqMan® Master.

Analysis Mode	Cycles	Segment	Target Temperature <sup>1)</sup>	Hold Time	Acquisition Mode
<b>Pre-Incubation</b>					
None	1		95°C	10 min <sup>2)</sup>	none
<b>Amplification</b>					
Quantification	45 <sup>3)</sup>	Denaturation	95°C	10 s	none
		Annealing	primer dependent <sup>4)</sup>	20 – 40 s <sup>5)</sup>	none
		Extension	72°C	1 s	single
<b>Cooling</b>					
None	1		40°C	30 s	none

<sup>1)</sup> Temperature Transition Rate/Slope is 20°C/sec.

<sup>2)</sup> If high polymerase activity is needed in early cycles, you can sometimes improve results by extending the pre-incubation to 15 min.

<sup>3)</sup> 45 cycles are suitable for most of the assays. Nevertheless, for a well optimized assay with steep amplification curves and early crossing points, even with low target concentrations, 40 cycles should be sufficient (resulting in a reduced assay time!).

<sup>4)</sup> For initial experiments, set the target temperature 5°C below the calculated primer melting temperatures.

<sup>5)</sup> For higher precision in target quantification experiments, in some cases it can be advantageous to choose relative long annealing times, up to 45 s when amplifying sequences > 500 bp.

**Fluorescence and  
Run Setup  
Parameters**

Parameter	Setting	
<b>All LightCycler® Software Versions</b>		
Seek Temperature	30°C	
LightCycler® Software prior to Version 3.5		
Display Mode	<ul style="list-style-type: none"> <li>• during the run      • fluorescence channel F1*</li> <li>• for analysis        • fluorescence channel F1*/F3</li> </ul> <p style="margin-left: 100px;">* for using a fluorescein-labeled Hydrolysis probe only</p>	
Fluorescence Gains	Fluorimeter	Gain Value
	Channel 1 (F1)	5
	Channel 2 (F2)	1
	Channel 3 (F3)	1
LightCycler® Software Version 3.5		
Display Mode	<ul style="list-style-type: none"> <li>• during the run      • fluorescence channel F1*</li> <li>• for analysis        • fluorescence channel F1*/F3</li> </ul>	
Fluorescence Gains	Not required	
	<p>Ⓢ In data created with LightCycler® Software Version 3.5, all fluorescence values are normalized to a fluorescence gain of “1”. This produces a different scale on the Y-axis than that obtained with previous LightCycler® software versions. This difference does not affect the crossing points or any calculated concentrations obtained.</p>	

Parameter	Setting
LightCycler® Software Version 4.0 (or higher)	
Default Channel	<ul style="list-style-type: none"> <li>during run</li> <li>for analysis</li> </ul>
	<ul style="list-style-type: none"> <li>Depending on the fluorescent reporter dye used for labeling the Hydrolysis probe, choose Channel 530 or 560.</li> <li>530/705 or 560/705</li> <li>⚠ Channel 560 is available on a LightCycler® 2.0 Instrument only.</li> </ul>
Fluorescence Gains	Not required
“Max. Seek Pos”	Enter the number of sample positions the instrument should look for.
“Instrument Type”	<ul style="list-style-type: none"> <li>“6 Ch.”: for LightCycler® 2.0 Instrument (selected by default)</li> <li>“3 Ch.”: for LightCycler® 1.5 Instrument and instrument versions below.</li> </ul>
“Capillary Size”	Select “20 $\mu$ l” as the capillary size for the experiment. ⚠ For the “6 Ch.” instrument type only.

### Preparation of the Master Mix

Prepare the 5 $\times$  conc. Master Mix as described below.

- 1 Thaw one vial of “Reaction Mix” (vial 1b, red cap).
  - 2 Briefly centrifuge one vial 1a (“Enzyme”, white cap) and one thawed vial 1b (“Reaction Mix”, red cap, from Step 1).
  - 3 Pipet 10  $\mu$ l from vial 1a (white cap) into vial 1b (red cap).  
 Ⓞ Each vial 1a contains enough enzyme for three vials of Reaction Mix.
  - 4 Mix gently by pipetting up and down. Do not vortex.
  - 5 Re-label vial 1b (red cap) with the new labels (vial 1: Master Mix) provided with the kit. Place one on the top of the cap and one on the side of the vial.
- Ⓞ The volume of the resulting Master Mix (5 $\times$  conc.) is sufficient for 32 reactions with a final volume of 20  $\mu$ l each.

**Preparation of the PCR Mix**

Proceed as described below for a 20  $\mu\text{l}$  standard reaction.

Ⓞ The protocol is designed for a final reaction volume of 20  $\mu\text{l}$ . For volumes < 20  $\mu\text{l}$ , the reaction and cycle conditions must be optimized.

⚠ Do not touch the surface of the capillaries. Always wear gloves when handling the capillaries.

① Depending on the total number of reactions, place the required number of LightCycler® Capillaries in precooled centrifuge adapters or in a LightCycler® Sample Carousel in a precooled LightCycler® Centrifuge Bucket.

② Prepare a 10x conc. solution that contains PCR primers and Hydrolysis probe.

③ Thaw the LightCycler® TaqMan® Master, 5 $\times$  conc. (vial 1, red cap), mix gently and store on ice.

④ In a 1.5 ml reaction tube on ice, prepare one reaction by adding the following components in the order listed below, then mix gently by pipetting up and down:

Component	Volume
Water, PCR-grade (vial 2, colorless cap)	9 $\mu\text{l}$
Primers/Probe, 10 $\times$ conc.	2 $\mu\text{l}$
Master Mix, 5 $\times$ conc. (vial 1)	4 $\mu\text{l}$
<b>Final volume</b>	<b>15 <math>\mu\text{l}</math></b>

Ⓞ To prepare the PCR Mix for more than one reaction, multiply the amount in the "Volume" column above by the number of reactions to be run + one additional reaction.

⑤ • Pipet 15  $\mu\text{l}$  PCR mix into each precooled LightCycler® Capillary.  
• Add 5  $\mu\text{l}$  of the DNA template.

Ⓞ Use up to 500 ng complex genomic DNA or  $10^1 - 10^{10}$  copies plasmid DNA.

⑥ • Seal each capillary with a stopper and place the adapters (containing the capillaries) in a standard benchtop microcentrifuge.

⚠ Place the centrifuge adapters in a balanced arrangement in the centrifuge.

• Centrifuge at  $700 \times g$  for 5 s (3,000 rpm in a standard benchtop microcentrifuge). Alternatively, use the LC Carousel Centrifuge for spinning the capillaries.

⑦ Transfer the capillaries to the sample carousel of the LightCycler® Instrument.

⑧ Cycle the samples as described above.

## 2.3 Related Procedures

---

### Prevention of Carry-Over Contamination

Uracil DNA N-Glycosylase (UNG) is suitable for preventing carry-over contamination in PCR. This carry-over prevention technique involves incorporating deoxyuridine triphosphate (dUTP, a component of the Master Mix in this kit) into amplification products, then pretreating later PCR mixtures with UNG. If a dUTP-containing contaminant is present in the later PCRs, it will be cleaved by a combination of the UNG and the high temperatures of the initial denaturation step; it will not serve as a PCR template.

⚠ Since your target DNA template contains thymidine rather than uridine it is not affected by this procedure.

⚠ Refer to the pack insert of LightCycler® Uracil-DNA Glycosylase for details on application.

### Two-step RT-PCR

The LightCycler® TaqMan® Master can also be used to perform two-step RT-PCR.

In two-step RT-PCR, the reverse transcription of RNA into cDNA is separated from the other reaction steps and is performed outside the LightCycler® 2.0 System. Subsequent amplification and online monitoring is performed according to the LightCycler® 2.0 System standard procedure, using cDNA as starting sample material.

One of the following reagents is required for reverse transcription of RNA into cDNA (see Ordering Information for details):

- Transcriptor Reverse Transcriptase
- Transcriptor First Strand cDNA Synthesis Kit
- First Strand cDNA Synthesis Kit for RT-PCR (AMV)


Synthesis of cDNA is performed according to the detailed instructions provided with the cDNA synthesis reagent.

⚠ Do not use more than 8 µl of undiluted cDNA template per 20 µl final reaction volume because greater amounts may inhibit PCR. For initial experiments, we recommend running undiluted, 1:10 diluted, and 1:100 diluted cDNA template in parallel to determine the optimum template amount.

### 3. Results

---

The following amplification curves were obtained with dilutions of human genomic DNA in 1:10 steps and a set of primers and a FAM/TAMRA-labeled probe, specific for a fragment of 442 bp of the cyclophilin A gene.





**Fig. 1:** Cyclophilin A amplification curves in the analysis module for absolute quantification of LightCycler<sup>®</sup> Software 4.0.

## 4. Troubleshooting

	Possible cause	Recommendation
<b>Fluorescence curves reach maximum long before cycling is finalized.</b>	Very high starting amount of nucleic acid.	Cycling can be finished by clicking the <b>End Program</b> button. The next program part will continue automatically.
	Number of cycles is too high.	Reduce the number of cycles in the protocol.
<b>Log-linear phase of amplification just starts when cycle program finishes.</b>	Number of cycles is too low.	<ul style="list-style-type: none"> <li>• While cycling is still ongoing use the <b>Add 10 Cycles</b> button to increase the number of cycles.</li> <li>• Increase the number of cycles in the protocol.</li> <li>• Use higher amount of starting material.</li> <li>• Optimize PCR conditions (primer/probe design, protocol).</li> </ul>
	<b>No amplification detectable.</b>	Wrong channel has been chosen to display amplification on screen.
<b>No amplification detectable.</b>	Inhibitory effects of the sample material due to insufficient purification.	<ul style="list-style-type: none"> <li>• Try a 1 : 10 dilution of your sample.</li> <li>• Isolate the nucleic acids of your sample material to ensure removal of inhibitory agents.</li> </ul>
	FastStart DNA Polymerase is not enough activated.	<ul style="list-style-type: none"> <li>• Check whether PCR was started with 10 min pre-incubation step at 95°C.</li> <li>• Check whether denaturation time in cycles was about 10 sec.</li> </ul>
	Pipetting errors or omitted reagents.	<ul style="list-style-type: none"> <li>• Check for missing or defective reagents.</li> </ul>
	Amplicon length is > 1 kb.	Do not use primers for amplicons > 1 kb. Amplification efficiency is reduced with amplicons of this size. Optimal results are obtained for amplicons < 700 bp.
	Difficult template <i>e.g.</i> , unusual GC-rich sequence.	<ul style="list-style-type: none"> <li>• Optimize temperatures and times in the amplification cycles.</li> <li>• Optimize primer/probe sequences.</li> <li>• Repeat PCR with increasing content of DMSO, up to 10 % final concentration.</li> </ul>
<b>Fluorescence intensity varies.</b>	Some of the reagent is still in the upper vessel of the capillary, or an air bubble is trapped in the capillary tip.	Ensure the centrifugation step was properly performed.
	Skin oils or dirt on the surface of the capillary.	Always wear gloves when handling the capillaries.

#### 4. Troubleshooting, continued

	Possible cause	Recommendation
<b>Fluorescence intensity is very low.</b>	Low concentration or deterioration of dyes in the reaction mixtures due to unsuitable storage conditions.	<ul style="list-style-type: none"><li>• Protect dye-labeled reagents from light.</li><li>• Store the reagents at <math>-15</math> to <math>-25^{\circ}\text{C}</math> and avoid repeated freezing and thawing.</li><li>• After thawing, store the Master Mix at <math>2 - 8^{\circ}\text{C}</math> for a maximum of 1 week.</li></ul>
	Poor PCR efficiency due to non-optimized reaction conditions.	Check concentrations of reagents and LightCycler 2.0 System protocol.
	Chosen gain settings are too low. (Only with LightCycler® Software version < 3.3)	Use the <b>Real Time Fluorimeter</b> to find suitable gain settings. Change gain settings in the protocol and repeat the run.  With LightCycler® software version > 3.3 no gain setting is required.
<b>Overflow of fluorescence intensity. (Only with LightCycler® Software version 3.3)</b>	Unsuitable gain settings.	Gain settings cannot be corrected during or after a run. Before repeating the run, use the <b>Real Time Fluorimeter</b> to find suitable gain settings.  With LightCycler® software version > 3.3 no gain setting is required.
	<b>Negative control is positive.</b>	Contamination. <ul style="list-style-type: none"><li>• Exchange all critical solutions.</li><li>• Pipet on a clean bench.</li><li>• Close lid of the negative control reaction tube immediately after pipetting.</li><li>• Use UNG for prevention of contamination by amplicons.</li></ul>
<b>High background</b>	Very low fluorescence signals, therefore the background seems to be relatively high.	Follow general optimization strategies for LightCycler® 2.0 System PCR.
	Insufficient quality of probe.	Prepare a new solution of probe.

## 5. Additional Information on this Product

---

### 5.1 How this Product Works

LightCycler® TaqMan® Master is a ready-to-use reaction mix designed for the TaqMan® detection format using the LightCycler® 2.0 System.

It contains FastStart polymerase for a "Hot Start" PCR, which has been shown to significantly improve the specificity and sensitivity of PCR by minimizing the formation of non-specific amplification products [1, 2, 3, 4]. FastStart Taq DNA Polymerase is a chemically modified form of thermostable recombinant Taq DNA polymerase that is inactive at room temperature and below. The enzyme is active only at high temperatures, where primers no longer bind non-specifically. The enzyme is completely activated (by removal of blocking groups) in a single pre-incubation step (95°C, 10 minutes) before cycling begins. Activation does not require the extra handling steps typical of other "Hot Start" techniques.

The LightCycler® TaqMan® Master provides convenience, excellent performance, reproducibility, and minimal contamination risk. All you have to supply is PCR primers, a detection probe and your template DNA.

- Ⓢ The reaction mix in this kit is optimized for a fixed MgCl<sub>2</sub> concentration, which works with nearly all primer combinations. You do not need to adjust the MgCl<sub>2</sub> concentration to amplify different sequences.
- ⚠ The performance of the kit described in this Instruction Manual is guaranteed only when it is used with the LightCycler® 2.0 System.

### Test principle

Hydrolysis probe assays, also called TaqMan® assays, use a single probe containing two labels, a fluorescent reporter dye and a fluorescent quencher. While the probe is intact, the quencher is close to the reporter dye and suppresses the reporter fluorescence via fluorescence resonance energy transfer (FRET). When the probe is hybridized to the target sequence, the 5'-nuclease activity of the polymerase can cleave the hydrolysis probe, separating reporter and quencher. With a rising amount of target sequence during PCR, more probe is cleaved, and the fluorescence signal of the unquenched reporter dye increases.

- ⚠ As the principle of Hydrolysis probe assays is probe cleavage during PCR, Hydrolysis probes cannot be used to perform a melting curve analysis. In contrast, HybProbe probes which consist of two specially designed, sequence-specific oligonucleotide probes labeled with different dyes, are still intact at the end of amplification, and thus may be used in a subsequent melting curve experiment (*e.g.*, for mutation detection or SNP analysis).

### References

- 1 Chou, Q. *et al* (1992) *Nucleic Acid Res.* **20**,1717-1723.
- 2 Kellogg, D.E. *et al* (1994) *BioTechniques* **16**,1134-1137.
- 3 Birch, D.E. *et al* (1996) *Nature* **381**,445-446.
- 4 PCR Manual, Roche Diagnostics (1999) 2nd edition (1999) **2**, 52-58.

### Quality Control

The LightCycler® TaqMan® Master is function tested using the LightCycler® 2.0 System.

## 6. Supplementary Information

---

### 6.1 Conventions

**Text Conventions** To make information consistent and memorable, the following text conventions are used in this Instruction Manual:

Text Convention	Usage
Numbered instructions labeled ❶, ❷, etc.	Steps in a procedure that must be performed in the order listed.
Asterisk *	Denotes a product available from Roche Applied Science.

### Symbols

In this Instruction Manual, the following symbols are used to highlight important information:

Symbol	Description
ⓘ	Information Note: Additional information about the current topic or procedure.
⚠	Important Note: Information critical to the success of the procedure or use of the product.

### 6.2 Changes to Previous Version

- Wrong gain settings corrected.  
Correct setting is Channel F1: 5, Channel F2: 1, Channel F3: 1.
- ⓘ Only relevant when working with LightCycler<sup>®</sup> Software prior to Version 3.5.

## 6.3 Ordering Information

Roche Applied Science offers a large selection of reagents and systems for life science research. For a complete overview of related products and manuals, please visit and bookmark our home page, [www.roche-applied-science.com](http://www.roche-applied-science.com), and our Special Interest Sites including:

- The LightCycler® 2.0 System family for real-time, online PCR:  
<http://www.lightcycler-online.com>
- The MagNA Pure family for automated nucleic acid isolation:  
<http://www.roche-applied-science.com/sis/magnapure/>
- Manual Nucleic Acid Isolation and Purification:  
<http://www.roche-applied-science.com/napure/>

	Product	Pack Size	Cat. No.	
<b>Instrument and Accessories</b>	LightCycler® 2.0 Instrument	1 instrument plus accessories	03 351 414 001	
	LightCycler® 1.5 Instrument	1 instrument plus accessories	04 484 495 001	
	LightCycler® Capillaries (20 µl)	1 pack (8 boxes, each with 96 capillaries and stoppers)	11 909 339 001	
	LC Carousel Centrifuge 2.0	1 rotor + 2 buckets	03 724 697 001	
	LC Carousel Centrifuge 2.0 Bucket 2.1	1 bucket	03 724 689 001	
	LC Carousel Centrifuge 2.0	1 centrifuge plus rotor and bucket	03 709 507 001 (115 V) 03 709 582 001 (230 V)	
<b>Software</b>	LightCycler® Software 4.05	1 software package	04 717 392 001	
	LightCycler® Probe Design Software 2.0	1 software package	04 342 054 001	
<b>LightCycler® Kits for PCR</b>	LightCycler® DNA Master Hyb-Probe	1 kit (96 reactions) 1 kit (480 reactions)	12 015 102 001 12 158 825 001	
	LightCycler® FastStart DNA Master HybProbe	1 kit (96 reactions) 1 kit (480 reactions)	03 003 248 001 12 239 272 001	
	LightCycler® FastStart DNA Master <sup>PLUS</sup> HybProbe	1 kit (96 reactions) 1 kit (480 reactions)	03 515 575 001 03 515 567 001	
	LightCycler® FastStart DNA Master <sup>PLUS</sup> HybProbe, 100 µl Reactions	1 kit (384 reactions)	03 752 178 001	
	LightCycler® DNA Master SYBR Green I	1 kit (96 reactions) 1 kit (480 reactions)	12 015 099 001 12 158 817 001	
	LightCycler® FastStart DNA Master SYBR Green I	1 kit (96 reactions) 1 kit (480 reactions)	03 003 230 001 12 239 264 001	
	LightCycler® FastStart DNA Master <sup>PLUS</sup> SYBR Green I	1 kit (96 reactions) 1 kit (480 reactions)	03 515 869 001 03 515 885 001	
	LightCycler® FastStart DNA Master <sup>PLUS</sup> SYBR Green I, 100 µl Reactions	1 kit (384 reactions)	03 752 186 001	
	<b>Associated Kits and Reagents</b>	LightCycler® Uracil-DNA Glycosylase	100 U (50 µl)	03 539 86 001

## 6.4 License Disclaimer

---

LIMITED LICENSE: A license under the non-U.S. counterparts of U.S. Patents Nos. 4,683,202, 4,683,195 and 4,965,188 owned by F. Hoffmann-La Roche Ltd. (Roche), for use in research and development, has an up-front fee component and a running-royalty component. The purchase price of the product includes a limited, non-transferable license under the running-royalty component to use only this amount of product to practice the Polymerase Chain Reaction (PCR) and related processes described in said patents solely for the research and development activities of the purchaser when this product is used in conjunction with a thermal cycler whose use is covered by the up-front fee component. Rights to the up-front fee component must be obtained by the end user in order to have a complete license. These rights under the up-front fee component may be purchased from Applied Biosystems or obtained by purchasing an authorized thermal cycler. This product is also an authorized reagent and may also be used under service sublicenses from Applied Biosystems under the foregoing patents.

A license to perform the 5' nuclease process for research requires the use of a Licensed 5' Nuclease Kit (containing Licensed Probe), or the combination of an Authorized Core Kit plus Licensed Probe, or license rights that may be purchased from Applied Biosystems. This product is an Authorized Core Kit without Licensed Probe. Its purchase price includes a limited, non-transferable immunity from suit under U.S. Patents Nos. 5,210,015, 5,487,972, 5,476,774, and 5,219,727, and corresponding patent claims outside the United States, owned by Roche Molecular Systems, Inc. or F. Hoffmann-La Roche Ltd ("Roche"), for using only this amount of the product in the practice of the 5' nuclease process solely for the purchaser's own internal research and development activities. This product is also an Authorized Core Kit for use with service sublicenses available from Applied Biosystems. This product conveys no rights under U.S. Patents Nos. 5,804,375, 6,214,979, 5,538,848, 5,723,591, 5,876,930, 6,030,787, or 6,258,569, or corresponding patent claims outside the United States, expressly, by implication, or by estoppel.

No right under any other patent claims (such as apparatus or system claims in U.S. Patent No. 6,814,934) and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is hereby granted expressly, by implication, or by estoppel. This product is for research purposes only. Diagnostic uses require a separate license from Roche.

Further information on purchasing licenses to practice real-time PCR processes may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

The technology used for the LightCycler<sup>®</sup> System is licensed from Idaho Technology Inc., Salt Lake City, UT, USA.

## 6.5 Trademarks

LIGHTCYCLER, LC, MAGNA PURE, HIGH PURE, FASTSTART, HYBPROBE, and TAQMAN are Trademarks of Roche.



---

## Contact and Support

If you have questions or experience problems with this or any Roche Applied Science (RAS) product, please contact our Technical Support staff. Our scientists commit themselves to providing rapid and effective help.

We also want you to contact us if you have suggestions for enhancing RAS product performance or using our products in new or specialized ways. Such customer information has repeatedly proven invaluable to RAS and the worldwide research community.

To ask questions, solve problems, suggest enhancements or report new applications, please visit our **Online Technical Support** Site at:

**[www.roche-applied-science.com/support](http://www.roche-applied-science.com/support)**

To call, write, fax, or email us, visit the Roche Applied Science home page, [www.roche-applied-science.com](http://www.roche-applied-science.com), and select your home country. Country-specific contact information will be displayed.

On the Roche Applied Science home page select **Printed Materials** to find:

- in-depth Technical Manuals
- Lab FAQs: Protocols and references for life science research
- our quarterly Biochemica Newsletter
- Material Safety Data Sheets
- Pack Inserts and Product Instructions

or to request hard copies of printed materials.

---



## Diagnostics

Roche Diagnostics GmbH  
Roche Applied Science  
68298 Mannheim  
Germany