

# Glucose Bio

Cedex Bio

**REF 06 343 732 001**

2 × 100 tests

Applications GLC2B (932), GLC2D (901)

**Version 06**  
 Content version: January 2015


Store at +2 to +8°C

## Intended Use

This product is intended for quantitative determination of glucose in aqueous solutions using the Cedex Bio Instrument.

## Test Principle

Glucose is phosphorylated by ATP in presence of hexokinase (HK) to glucose-6-phosphate (G-6-P), which is oxidized by NADH in presence of glucose-6-phosphat dehydrogenase (G-6-PDH). The rate of NADPH formation is measured UV-photometrically and is directly proportional to the glucose concentration



## Contents

Vial/Cap	Content	Function / Composition
<b>1</b> R1 white cap	2 vials, 100 tests each	<b>Reagent 1:</b> TRIS, 100 mmol/L; Mg <sup>2+</sup> , 4 mmol/L; ATP, 1.7 mmol/L; NADP, ≥ 1 mmol/L; preservative; pH 7.8
<b>2</b> SR black cap	2 vials, 100 tests each	<b>Start Reagent:</b> HEPES, 30 mmol/L; Mg <sup>2+</sup> , 4 mmol/L; HK (yeast), ≥ 130 µkat/L; G-6-PDH ( <i>E. coli</i> ), ≥ 250 µkat/L; preservative; pH 7.0

## Storage and Stability

Store at +2 to +8°C and at +2 to +10°C for on-board use.

The kit is stable at +2 to +8°C until the expiration date printed on the label when stored unopened in original vials and kept free of contamination.

On-board stability: After first use on the instrument, the reagents are stable for up to 4 weeks.

## Additional Materials Required

- Cedex Bio Instrument, Cat. No. 06 395 554 001, with general accessories and disposables
- Standard laboratory equipment
- Calibrator A Bio, Cat. No. 06 682 189 001
- Control A Level 1 Bio, Cat. No. 06 682 197 001
- Control A Level 2 Bio, Cat. No. 06 682 227 001
- Control A Level 3 Bio, Cat. No. 06 682 545 001

## Test Protocol

GLC2B is the standard protocol for glucose determination. If the glucose concentration of samples is expected higher than 7.5 g/L (41.6 mmol/L), then protocol GLC2D shall be used (same protocol with preset 1:10 sample predilution).

Results higher than the valid measuring range will be flagged and a rerun with an automatic 1:10 predilution will be recommended by the software.

## Cedex Bio test definition:

Measuring mode	Absorbance
Abs. calculation mode	Endpoint
Reaction direction	Increase
Wavelength A/B	340/409 nm
Calc. first/last	16/37
Unit	mmol/L, mg/L
Reaction mode	R1-S-SR

## Pipetting parameters:

R1	150 µl
Sample	2 µl + 20 µl Diluent (H <sub>2</sub> O)
SR	30 µl
Total volume	202 µl

## Calibration:

Calibrator	Calibrator A Bio (CAL A, 599)
Calibration mode	Linear regression
Calibration intervall	Each new lot, and if recalibration is required due to QC results

## Quality control:

Controls	Control A Level 1 Bio (CONA1, 111) Control A Level 2 Bio (CONA2, 112) Control A Level 3 Bio (CONA3, 113)
----------	----------------------------------------------------------------------------------------------------------------

Use the recommended control material. Other suitable control material can be used in addition.

Intervals / limits Control intervals and acceptance limits should be adapted to each laboratory's individual requirements.

If values do not fall into the defined limits, corrective measures and recalibration are required

## Measuring Range

Glucose concentration can be determined in following ranges:

- GLC2B: 0.02 - 7.50 g/L (0.111 - 41.6 mmol/L)
- GLC2D: 0.2 - 75.0 g/L (1.11 - 416.0 mmol/L), up to maximal solubility with postdilution

## Conversion Factors

Conversion factors for glucose concentration:

- 1 mmol/L = 180.2 mg/L
- 1 g/L = 5.550 mmol/L

## Traceability

This method has been standardized against Isotope dilution mass spectrometry (ID-MS).

## Precision

Representative performance data on Cedex Bio analyzers are shown. Results obtained in individual laboratories may differ.

Precision was determined in samples of three concentration levels. Coefficients of variation (CV) were calculated for in-run precision (n=21) and inter-run precision (on 10 days).


	Level 1	Level 2	Level 3
Mean	2.52 mmol/L (1.14 mg/L)	9.89 mmol/L (1.75 mg/L)	14.1 mmol/L (4.86 mg/L)
CV in-run	1.1%	0.8%	0.5%
CV inter-run	0.5%	0.6%	0.5%

## Changes to Previous Version


- Editorial Changes
- Correction of "Calc. first/last" Value

## Test Protocol Barcodes

Glucose, GLC2B (ACN 932) v02



Glucose with 1:10 predilution, GLC2D (ACN 901) v02



## Trademarks

CEDEX is a trademark of Roche.

## Regulatory Disclaimer

For use in quality control/ manufacturing process only.

## Online Technical Support

Please visit our Online Technical Support site for additional information about this product: [www.technical-support.roche.com](http://www.technical-support.roche.com)

**For more information about this product, as well as documentation such as Instructions for Use and Material Safety Data Sheets, please visit [custombiotech.roche.com](http://custombiotech.roche.com)**

© 2015 Roche Diagnostics GmbH. All rights reserved.

**For more information, contact your local Roche representative.**

Europe, Middle East, Africa,  
and Latin America:  
Phone +49 621 759 8580  
Fax +49 621 759 6385

Asia Pacific:  
Phone +65 6371 6638  
Fax +65 6371 6601

Japan:  
Phone +81 3 5443 5285  
Fax +81 3 5443 7934

Canada:  
Phone +1 450 686 7050  
Fax +1 450 686 7012

United States:  
Phone +1 800 428 5433,  
ext. 14649 (toll-free)  
Fax +1 317 521 4065

Roche Diagnostics GmbH  
Sandhofer Strasse 116  
68305 Mannheim  
Germany

