



Certificate of Analysis

Product name

Glycogen , MB Grade

Material No.

10901393001

Lot No.

64438320

Date of Manufacturing

2022-09-20

Expiry Date

2025-09

Parameter**Result****Unit**

Appearance

clear,colourless solution

Concentration

20

mg/ml

Endonucleases up to 200 mcg w.lambda-DNA

not detectable

Nicking act. up to 200 mcg w.pBR 322-DNA

not detectable

Ribonucleases up to 200 mcg w. MS II RNA

not detectable

Nucleic acid up to 200mcg Glycog.applic.

not detectable

Proteases up to 200 mcg Glycogen applic.

not detectable

This lot meets our specification

This certificate is a computer printout and has therefore not been signed by hand.

Issue date 2022-09-20

Page 1 of 1

Roche Diagnostics GmbH

Sandhofer Straße 116; D#68305 Mannheim; Telefon +49#621#759#0; Telefax +49#621#759#2890

Sitz der Gesellschaft: Mannheim - Registergericht: AG Mannheim HRB 3962 - Geschäftsführung: Claus Haberda;

Clemens Schmid - Aufsichtsratsvorsitzender: Dr. Thomas Schinecker