




## Certificate of Analysis

**Product Name:** Cytokeratin 5 (SP27) Rabbit Monoclonal Primary Antibody  
**Catalog Number(s):** 760-4935  
**Roche Part Number:** 07292821001  
**Lot Number:** V0000854  
**Date of Manufacture:** 2017-08-08  
**Expiration Date:** 2020-05-31

Test	Acceptance Criteria	Results
Functional testing on the following tissue(s): Mesothelioma	Intensity difference between sample and a reference must be less than or equal to 0.5	<input checked="" type="checkbox"/> PASS <input type="checkbox"/> FAIL
	Background difference between sample and a reference must be less than or equal to 0.5	<input checked="" type="checkbox"/> PASS <input type="checkbox"/> FAIL
Appearance	Transparent liquid and free of visible foreign material.	<input checked="" type="checkbox"/> PASS <input type="checkbox"/> FAIL
Bioburden	Visible bacterial colonies must be less than 100CFU/mL.	<input checked="" type="checkbox"/> PASS <input type="checkbox"/> FAIL
ELISA	0.01-1.0 µg/mL	0.17 µg/mL
Source	N/A	Supernatant

Additional information: N/A

 **CELL MARQUE**  
6600 Sierra College Blvd.  
Rocklin, California 95677  
(916) 746-8900 Fax: (916) 746-8989

Authorized Signature: 

Date: 10/31/17

Quality Assurance Representative

Issued: 6/3/16 Rev 1.0

Document #QF0901-008-02 [CM Legacy #]

Manufactured by: Cell Marque Corporation  
6600 Sierra College Boulevard • Rocklin, California 95677 • USA • Telephone (916) 746-8900 • Fax (916) 746-8976