

**Certificate of Analysis**

**Product Name:** GCDFP-15 (EP1582Y) Rabbit Monoclonal Antibody  
**Catalog Number(s):** 760-4386  
**Roche Part Number:** 05463530001  
**Lot Number:** V0000929  
**Date of Manufacture:** 2017-08-29  
**Expiration Date:** 2020-07-31

Test	Acceptance Criteria	Results
Functional testing on the following tissue(s): Breast	Intensity difference between sample and a reference must be less than or equal to 0.5	<input checked="" type="checkbox"/> PASS <input type="checkbox"/> FAIL
	Background difference between sample and a reference must be less than or equal to 0.5	<input checked="" type="checkbox"/> PASS <input type="checkbox"/> FAIL
Appearance	Transparent liquid and free of visible foreign material.	<input checked="" type="checkbox"/> PASS <input type="checkbox"/> FAIL
Bioburden	Visible bacterial colonies must be less than 100CFU/mL.	<input checked="" type="checkbox"/> PASS <input type="checkbox"/> FAIL
ELISA	0.01-1.0 µg/mL	0.39 µg/mL
Source	N/A	Supernatant

Additional information: N/A



Authorized Signature:  Date: 10/30/17  
Quality Assurance Representative

Issued: 6/3/16 Rev 1.0

Document #QF0901-008-02 [CM Legacy #]