



## Certificate of Analysis

**Product name**

KIT COBAS 6800/8800 MPX 96T CE-IVD

**Material No.**

06997708190

**Lot No.**

242567

**Date of Manufacturing**

2017-05-03

**Expiry Date**

2019-01

Parameter	Result	Unit
RFID Verification	Pass	
The RFID label matches: a. The kit material number b. The corresponding human readable label c. The batch and date of expiration in SAP		
HIV 2.5X LOD Solution	100	%
HBV 2.5X LOD Solution	100	%
HCV 2.5X LOD Solution	100	%
HIV-1O 2.5X LOD Solution	100	%
HIV-2 2.5X LOD Solution	100	%
HIV 0.5X LOD Solution	95	%
HBV 0.5X LOD Solution	90	%
HCV 0.5X LOD Solution	90	%
HIV-1O 0.5X LOD Solution	90	%

**Issue date** 2017-06-23

Page 1 of 2



## Certificate of Analysis

**Product name**

KIT COBAS 6800/8800 MPX 96T CE-IVD

**Material No.**

06997708190

**Lot No.**

242567

**Date of Manufacturing**

2017-05-03

**Expiry Date**

2019-01

**Parameter****Result****Unit**

HIV-2 0.5X LOD Solution

76

%

NHP

100

%

Component Check

1

This lot meets our specification

**Disposition by**

Martina Gjini

This certificate is a computer printout and has therefore not been signed by hand.

**Issue date** 2017-06-23

Page 2 of 2

Roche Diagnostics International Ltd  
Forrenstrasse 2  
CH-6343 Rotkreuz  
Switzerland  
Tel: +41-41- 799 22 44