



a member of the Roche group

CERTIFICATE OF ANALYSIS

Product: CONFIRM anti-estrogen Receptor (ER) (SP1) Rabbit Monoclonal Primary Antibody
Ventana REF: 790-4325
Roche Number: 05278414001
Batch Number: D08925

Test	Acceptance Criteria	Result
Functional Staining on test and reference slides	Nuclear staining must be within a 0.5 intensity score when compared to the reference slide on a scale of 0-4+	PASS
	Approximate percent 0-50% (negative) tumor cell staining	PASS
	Approximate percent 51-100% (positive) tumor cell staining	PASS
	Intensity difference between sample and a reference must be less than or equal to 0.5+	PASS
	Background staining reading must score 0.5+ or less	PASS
Visual Inspection	Reagent must be free of foreign material and substantially free of intrinsic particulate matter	PASS
	Reagent must be clear and free of turbidity	PASS
Bioburden	Visible bacterial colonies must be less than or equal to 100 CFU/mL	PASS

This lot met the specifications.

Expiry date: 2015-10-28

Manufacture date: 07-Nov-2013

Laura A. Jimenez

12-NOV-2013

Quality Department

Date