



VENTANA®

a member of the Roche group

CERTIFICATE OF ANALYSIS

Product Name: PAS Staining Kit
Ventana [REF]: 860-014
Roche Number: 05279291001
Batch Number: G04868

Test	Acceptance Criteria	Result
Functional Staining on at least two of the following: kidney, small intestine, liver, skin, and aorta tissue	Glycogen should stain bright magenta	PASS
	Nuclei should stain purple	PASS
	Negative Control – cells where only the cell nuclei stain purple	PASS

This lot met the specifications.

Expiry Date: 2017-11-28

Manufacture Date : 28-JUL-2016

Ken Eblen

K Eblen *Aug 6/2016*

Quality Department

Date