



Certificate of Analysis

Product name

Hematoxylin II

Material No.

05277965001

Lot No.

N25784

Date of Manufacturing

2026-01-28

Expiry Date

2027-07-20

Parameter**Result****Unit**

Functional Staining

Pass

Tissue type (liver, kidney, colon, skin, tonsil)

Nucleus of cell is stained dark purple/blue.

Cytoplasm and stromal areas are stained light purple/blue.

Slide does not exhibit excess precipitate that interferes with stain interpretation

Staining intensity must fall between the acceptable light and dark reference slides

Visual Inspection

Pass

Reagent must be clear of foreign material and particulate matter, and free of turbidity

This lot meets our specification

Disposition by

Lupe Villarreal

This certificate is a computer printout and has therefore not been signed by hand.

Issue date 2026-02-09

Page 1 of 1

Ventana Medical Systems, Inc.
1910 Innovation Park Drive
Tucson, Arizona 85755
Phone: + 1 800 227 2155
<http://www.ventana.com>