



Certificate of Analysis

Product name

CONFIRM anti-Melanosome (HMB45) Mouse M

Material No.

05479282001

Lot No.

H03140

Date of Manufacturing

2021-02-26

Expiry Date

2023-02-03

Parameter**Result****Unit**

Final Acceptance Testing

Pass

Melanoma

The test batch must score 3.5+ or greater on melanoma tissue

The specific stain intensity must also be within 0.5 pt. of the reference batch for a passing result

The background intensity must be 0.5+ or less

The intensity of the negative reagent control slide must be 0.5 or less

Bioburden

Pass

Visible bacterial colonies must be 100 or less CFU/mL

Visual Inspection

Pass

Reagent must be clear, and free of turbidity.

Reagent must be free of foreign material and substantially free of intrinsic particulate matter.

Issue date 2021-03-16

Page 1 of 2

Ventana Medical Systems, Inc.
1910 Innovation Park Drive
Tucson, Arizona 85755
Phone: + 1 520 887 2155
<http://www.ventana.com>



Certificate of Analysis

Product name

CONFIRM anti-Melanosome (HMB45) Mouse M

Material No.

05479282001

Lot No.

H03140

Date of Manufacturing

2021-02-26

Expiry Date

2023-02-03

This lot meets our specification

Disposition by

Mia Cazares

This certificate is a computer printout and has therefore not been signed by hand.

Issue date 2021-03-16

Page 2 of 2

Ventana Medical Systems, Inc.
1910 Innovation Park Drive
Tucson, Arizona 85755
Phone: + 1 520 887 2155
<http://www.ventana.com>