



# Certificate of Analysis

**Product name**

GIEMSA STAINING KIT

**Material No.**

05279224001

**Lot No.**

K28837

**Date of Manufacturing**

2024-02-08

**Expiry Date**

2025-09-30

**Parameter****Result****Unit**

FINAL ACCEPTANCE TESTING

Pass

on tonsil and tissue infected with H. pylori. Monocytes and lymphocytes stain purple, erythrocytes stain pink, mast cells stain blue. H. pylori stain blue. There is no negative control for this stain.

This lot meets our specification

**Disposition by**

Joshua Burgette

This certificate is a computer printout and has therefore not been signed by hand.

**Issue date** 2024-02-26

Page 1 of 1

Ventana Medical Systems, Inc.  
1910 Innovation Park Drive  
Tucson, Arizona 85755  
Phone: + 1 800 227 2155  
<http://www.ventana.com>