



## Certificate of Analysis

**Product name**

BRAV V600E (VE1) Mouse Monoclonal PAB US

**Material No.**

07862270001

**Lot No.**

F09199

**Date of Manufacturing**

2019-05-01

**Expiry Date**

2021-04-30

**Parameter****Result****Unit**

Final Acceptance Testing

Pass

Positive Colon Carcinoma

The test batch of VENTANA anti-BRAF V600E (VE1) Mouse Monoclonal Primary Antibody must score 1.0 point or greater on colon CRC tissue with acceptable background.

The specific stain intensity must also be within 0.5 pt. of the reference batch for a passing result.

The background intensity must be equal to or less than 0.5 point on a 0-3 scale.

The intensity of the negative reagent control slide must be 0.5 point or less on a 0-3 scale.

Bioburden

Pass

Visible bacterial colonies must be 100 or less CFU/mL

Visual Inspection

Pass

**Issue date** 2019-05-21

Page 1 of 2



## Certificate of Analysis

**Product name**

BRAF V600E (VE1) Mouse Monoclonal PAB US

**Material No.**

07862270001

**Lot No.**

F09199

**Date of Manufacturing**

2019-05-01

**Expiry Date**

2021-04-30

**Parameter****Result****Unit**

Reagent must be clear, and free of turbidity.

Reagent must be free of foreign material and substantially free of intrinsic particulate matter.

This lot meets our specification

**Disposition by**

Lupe Villarreal

This certificate is a computer printout and has therefore not been signed by hand.

**Issue date** 2019-05-21

Page 2 of 2

Ventana Medical Systems, Inc.  
1910 Innovation Park Drive  
Tucson, Arizona 85755  
Phone: + 1 520 887 2155  
<http://www.ventana.com>