

**CERTIFICATE OF ANALYSIS**


---

Product: FITC Anti-Fibrinogen Primary Antibody  
 Ventana REF: 760-2685  
 Roche Number: 05267960001  
 Batch Number: G06214

Test	Acceptance Criteria	Result
Functional Staining on kidney or skin tissue	Kidney or skin slide staining intensity must be equivalent or stronger to the slide stained with the reference lot	PASS
Visual Inspection	Reagent must be free of foreign material and substantially free of intrinsic particulate matter	PASS
	Reagent must be clear and free of turbidity	PASS
Bioburden	Visible bacterial colonies must be less than or equal 100 CFU/mL	PASS


---

This lot met the specifications.

Expiry Date: 2018/08/28

Manufacture Date : 28/AUG/2016

Ken Eblen



AUG/2016

---

Quality Department

Date